Brazilian Journal of Microbiology (2008) 39:357-361 ISSN 1517-8382

# EVALUATION OF CULTURE MEDIA FOR COUNTS OF *BIFIDOBACTERIUM ANIMALIS* SUBSP. *LACTIS* BB 12 IN YOGHURT AFTER REFRIGERATED STORAGE

#### Luciano Fachin; Juliana Moryia; Ana Lourdes Neves Gândara<sup>1</sup>; Walkiria Hanada Viotto<sup>1\*</sup>

<sup>1</sup>Departamento de Tecnologia de Alimentos, Faculdade de Engenharia de Alimentos, Universidade Estadual de Campinas, Campinas, SP, Brasil

Submitted: August 23, 2006; Approved: February 20, 2008.

## ABSTRACT

The agar RCPB pH5 has been considered a good alternative for counts of *Bifidobacterium* in yoghurt. However, during the refrigerated storage of yoghurt it is extremely difficult to count this microorganism due to the size of the colonies, which are so small they require the aid of a stereoscope to count them. Another agar, MRS-LP, has been also recommended for counts of *Bifidobacterium* in the presence of yoghurt bacteria. This study evaluated the supplementation of RCPB pH5 agar with dehydrated liver extract and the salts KH<sub>2</sub>PO<sub>4</sub>, K<sub>2</sub>HPO<sub>4</sub>, FeSO<sub>4</sub>7H<sub>2</sub>O, MnSO<sub>4</sub>H<sub>2</sub>O and MgSO<sub>4</sub>7H<sub>2</sub>O, aiming at improving the differentiation of *Bifidobacterium* in yoghurt using the agar MRS-LP. The agar MRS-LP presented the same cell recovery as non-fortified RCPB pH5 agar, used as a standard medium, thus being considered a good option for counts of *Bifidobacterium* in yoghurt. The fortified RCPB pH5 also presented the same recovery as the standard RCPB pH5 medium, however, the addition of dehydrated liver extract to the RCPB pH5 agar considerably increased the size of the *Bifidobacterium* colonies after refrigerated storage, making differentiation of the colonies much easier and reliable when compared to the standard non-fortified RPCP pH5. The addition of the salts (KH<sub>2</sub>PO<sub>4</sub>, K<sub>2</sub>HPO<sub>4</sub>, FeSO<sub>4</sub>7H<sub>2</sub>O, MnSO<sub>4</sub>H<sub>2</sub>O and MgSO<sub>4</sub>7H<sub>2</sub>O) had no influence on the performance of the RCPB pH5 agar.

Key-words: count, Bifidobacterium, yoghurt, refrigerated storage.

## **INTRODUCTION**

Currently much attention is given to the use of probiotic yoghurts due to their potential benefits to human health (2). It is recommended that each product should indicate the minimum daily amount required for it to confer health benefit(s) based on *in vitro* (animal) and human studies (3). However, the yoghurt-manufacturing conditions are far from favourable for the survival of *Bifidobacterium*, resulting in a low number of viable cells during refrigerated storage. In addition, the lack of a standard method for the count of *Bifidobacterium* in the presence of yoghurt bacteria reduces confidence in the technology of production of these yoghurts. Few methods are efficient to correctly count *Bifidobacterium* in the presence of yoghurt cultures, causing difficulties in the quality control (11).

Onggo and Fleet (8) used RCPB agar to isolate and differentiate *Bifidobacterium* spp., *L. bulgaricus* and *S. thermophilus* in yoghurts. Rybka and Kailasapath (10) proposed adjustment of the pH of the agar to 5, to avoid development of the non-lactic flora. Ghodussi and Robinson (4) obtained good results with RCPB agar for the differential count of *Bifidobacterium* in yoghurts. On the other hand, Payne, Morris and Beers (9) did not get good differentiation between the *Bifidobacterium* colonies and those of the yoghurt in this agar, the differential count of *Bifidobacterium* not being possible. In laboratory tests by the present authors, RCPB pH5 agar gave excellent differentiation and recovery of *Bifidobacterium* in

<sup>\*</sup>Corresponding Author. Mailing address: Department of Food Technology, CP 6121, Faculty of Food Engineering, State University of Campinas, Brazil. Tel.: +55-019-3521-3988 Fax: +55 019-3289-3617. E-mail: walkiria@fea.unicamp.br

relation to the yoghurt cultures when the count was made soon after fermentation of the yoghurt, and also when inoculated by spreading the individual microorganisms, grown separately, on the surface of the agar (7). However, when this agar was used to count *Bifidobacterium* in yoghurts during refrigerated storage, the diameter of the colonies was so small, they could only be visualised with the aid of a stereoscope, magnifying from 15 to 20 times. Thus the capacity of *Bifidobacterium* to develop in RCPB pH5 agar under conditions of cell stress, as occurs during refrigerated storage, was considerably lower, indicating that this agar is not a good option for the count of *Bifidobacterium* in yoghurts during the shelf life.

Beehrens (1) described a medium containing a combination of low pH and propionic acid which makes the medium both elective and selective for all species of *Bifidobacterium*. When propionic acid was added to Columbia agar at pH 5.0 it enhanced the growth of bifidobacteria. Colonies sometimes reached 2 mm in diameter compared with only 0.2-0.5 mm on the control medium at pH 5.0 without propionic acid. Lapierre, Undeland and Cox (5) proposed the use of lithium chloride and sodium propionate to inhibit yoghurt bacteria in association with lactose-cystine-liver agar, in order to count *Bifidobacterium* in fermented milk products. Vinderola and Reinheimer (11) used these inhibitors associated with MRS agar, to count *Bifidobacterium* in yoghurts.

Based on the fortifications proposed by Teraguchi, Uehara, Ogasa and Mitsuoka (1978) in their study on the selective count of *Bifidobacterium*, as described in Laroia and Martin (6), the objective of this research was to improve the differentiation of RCPB pH5 agar and evaluate the use of MRS-LP agar (11) for counts of *Bifidobacterium* in yoghurt after refrigerated storage.

## MATERIAL AND METHODS

## Cultures

Pure freeze dried cultures of *Bifidobacterium animalis* subsp. *lactis* Bb 12 (*B. animalis* Bb 12, Chr. Hansen Indústria e Comércio Ltda, Valinhos, SP), *Streptococcus thermophilus* A (*S. thermophilus*, old brand Wiesby Starter Cultures & Media, current Danisco Brasil Ltda, Cotia, SP), *Lactobacillus delbrueckii subsp. bulgaricus* 18 (*L. bulgaricus*, old brand Wiesby Starter Cultures & Media, current Danisco Brasil Ltda, Cotia, SP) and *S. thermophilus* TA 040 and *L. bulgaricus* LB 340, both provided by Rhodia Food (current Danisco Brasil Ltda, Cotia, SP) were used.

## Culture Media MRS agar

MRS agar (DIFCO or MERCK) was used as the standard agar for the pure cultures of *B. animalis* Bb 12. The plates were maintained at 37°C for 72 hours under anaerobic conditions (Anaerogen<sup>®</sup> - anaerobiosis generator from Oxoid).

## RCPB pH5 agar

The RCA (OXOID) agar and the Prussian Blue dye (SIGMA-ALDRICH) were separately sterilised at 121°C for 15 minutes, then cooled to 50°C and mixed aseptically. An amount of the Prussian Blue dye equivalent to 0.03% by weight of the final agar was dissolved in distilled water (12% by volume of the agar), the rest of the water (88% by volume) being used to dissolve the commercial RCA agar ingredients, which was subsequently adjusted to pH5. The dye solution was not adjusted to pH5.

## **MRS-LP** agar

The MRS agar (OXOID or MERCK) was dissolved in an amount of water equivalent to 93% of the agar final volume and sterilised at 121°C for 15 minutes. The inhibitory agents (0.2 g/ L lithium chloride and 0.3 g/L sodium propionate) were dissolved in the remaining 7% of the water and this solution sterilised using a 0.22 mm filter. The filtered solution of antimicrobial agents was carefully mixed with the MRS agar at 50°C in a laminar flow chamber, avoiding the incorporation of air.

### Fortified RCPB pH5 agar

The RCPB pH5 agar was fortified with the following salts: KH<sub>2</sub>PO<sub>4</sub> (1 g/L), K<sub>2</sub>HPO<sub>4</sub> (1 g/L), FeSO<sub>4</sub>7H2O (0.01 g/L), MnSO<sub>4</sub>H<sub>2</sub>O (0.01 g/L) and MgSO<sub>4</sub>7H<sub>2</sub>O (0.2 g/L) and also with the dehydrated liver infusion (150 mL/L). The liver infusion was prepared according to Laroia and Martin (6) and sterilised separately at 121°C/15 minutes (both the precipitate and the supernatant of the extraction were sterilised in the same flask). For the subsequent addition to the agar, the sterilised extract was left to decant at room temperature, and the supernatant (150 mL/L) added to the RCA agar, already containing the salts and the water, in a laminar flow chamber. The pH of the RCA agar containing the salts, liver extract and water was adjusted to 5. The volume of water used to dissolve the Prussian Blue dye was 12% of the total volume of the RCA agar. For example: to produce 300 mL of fortified RCPB pH5 agar, the salts and the RCA agar (weighed according to the manufacturer's instructions) were dissolved in 219 mL Milli-Q water (Millipore®) and mixed with 45 mL liver infusion, before adjusting the pH to 5 and sterilising at 121°C/15 minutes. Simultaneously, the Prussian Blue dye (0.03%) was dissolved in 36 mL water and also sterilised at 121°C/15 minutes. The mixture of the dye with the agar was carried out in a laminar flow chamber, in the same way as the RCPB pH5 agar. In a second experiment, the RCPB pH5 agar was only fortified with the liver extract, without the addition of the salts (RCPB pH5-Liver).

## **Inoculum** activation

The inoculum of *B. animalis* Bb 12 was prepared in milk powder, reconstituted to 12% solids, supplemented with 1% yeast extract and sterilised at 115°C for 10 minutes. The freezedried culture of *B. animalis* Bb 12 was incubated at 37°C for 6 hours. The inoculum of the yoghurt cultures was prepared with skimmed milk powder, reconstituted to 11% solids and sterilised in an autoclave at 115°C for 10 minutes. The freeze-dried cultures of *S. thermophilus* and *L. bulgaricus* were activated at 37°C until coagulation of the milk (6-7 hours for *S. thermophilus* and 14-15 hours for *L. bulgaricus*).

# Evaluation of the agars RCPB pH5, MRS-LP and fortified RCPB pH5 for counts made after the refrigerated storage of yoghurt

The media RCPB pH5, MRS-LP and fortified RCPB pH5 were tested as following:

- a) pure culture test: the activated inoculum of *B. animalis* Bb 12, S. thermophilus TA 040 and L. bulgaricus LB 340 were enumerated on RCPB pH5, MRS-LP and fortified RCPB pH5. The B. animalis Bb 12 cultures were also enumerated on MRS agar used as standard medium. For enumeration of each pure culture, each inoculum was decimally diluted in 0.1% peptone water to 10<sup>-6</sup> and a 0.1 mL aliquot of each dilution being plated on the agar under study by the surface spread technique. The 10<sup>-4</sup> to 10<sup>-6</sup> dilutions were plated on the RCPB pH5 agar and in MRS standard medium for Bifidobacterium. The 10<sup>-1</sup> to 10<sup>-6</sup> dilutions were plated on the selective MRS-LP agar and the fortified RCPB pH5 agar. The same inoculum of each microrganism were used in each medium studied and the plates of each medium were incubated at the same conditions (37°C for 72 hours) in the same anaerobic jars using Anaerogen® anaerobiosis generator from Oxoid.
- b) yoghurt test: a probiotic yoghurt were produced incubating *B. animalis* Bb 12 (1% inoculum), *S. thermophilus* TA 040 (1% inoculum) and *L. bulgaricus* LB 340 (1% inoculum) at 37°C until pH 4,6 was reached and stored at 7°C for 30 days. After 30 days of storage time, *B. animalis* Bb 12 were enumerated on RCPB pH5, MRS-LP and fortified RCPB pH5 to evaluate the performance of each medium for *Bifidobacterium* counting after refrigerated storage of yoghurt. RCPB pH5 was used as standard medium for counting *B. animalis* Bb12 and *L. bulgarucius* in yoghurt because this medium was already described by Rybka and Kailasapath (10) for this purpose.

After incubation under anaerobiosis, the RCPB pH5 and fortified RCPB pH5 plates were exposed to air for a few hours to allow for the formation of a blue halo around the *L. bulgaricus* colonies. After incubation, the media were analysed with respect to differentiated or selective development and with respect to cell recovery of *B. animalis* Bb 12, being compared with the standard MRS agar in the case of the pure culture and with the RCPB pH5 agar for the probiotic yoghurt (Fig. 1).

#### Effect of adding salts in the fortification of RCPB pH5 agar

The following experiment was designed in order to evaluate the addition of salts ( $KH_2PO_4$ ,  $K_2HPO_4$ ,  $FeSO_47H_2O$ ,

MnSO<sub>4</sub>H<sub>2</sub>O and MgSO<sub>4</sub>7H<sub>2</sub>O) on the performance of fortified RCPB pH5 agar: a probiotic yoghurt was produced incubating 1% of each activated inoculum of *B. animalis* Bb 12, *S. thermophilus* A and *L. bulgaricus* 18 until pH 4,6 was reached. The count of *Bifidobacterium* cells was carried out after 30 days of refrigerated storage of yoghurt at 7°C in RCPB pH5 agar fortified with salts and infusion of dehydrated liver and in RCPB pH5 fortified with only the dehydrated liver infusion (RCPB pH5-Liver).

## Statistical analysis

The Tukey test was applied to check for differences in microbial recovery in the media studied using the Software Statistica 5.0. All trials were made in triplicate.

## **RESULTS AND DISCUSSION**

Evaluation of the agars RCPB pH5, MRS-LP and fortified RCPB pH5 for counts made after the refrigerated storage of yoghurt

Table 1 shows the recovery of *B. animalis* Bb 12, *S. thermophilus* and *L. bulgaricus* in each medium studied when in pure culture and the recovery of *B. animalis* Bb12 counts in MRS medium used as standard for pure culture test and *B. animalis* Bb12 counts in yoghurt after 30 days of refrigerated storage. The *B. animalis* Bb 12 count in all media under study in pure culture test was not significantly different from that in the standard MRS agar (p<0.05). Pure cultures of *S. thermophilus* were completely inhibited in the RCPB pH5 and fortified RCPB pH5 agars, probably due to the low pH value, as described by Rybka and Kailasapathy (10). In previous studies carried out in the author's laboratory, it was seen that *S. thermophilus* usually developed to about 8-9 log CFU/mL in

**Table 1.** *S. thermophilus, L. bulgaricus* and *B. animalis* Bb 12 counts (log CFU/mL) in the various media studied in pure cultures (Bb 12, St and Lb) and *B. animalis* Bb 12 in yoghurt (Yoghurt) after 30 days refrigerated storage.

Cultures	MRS	RCPB pH5	MRS-LP	Fortified RCPB pH5
Bb12 <sup>a</sup>	9,47±0,12	9.54±0,25*	9.49±0,21*	9.53±0,22*
St <sup>a</sup>	nd	ng	5.19±0,15	ng
Lb <sup>a</sup>	nd	8.55±0,15	7.60+0,20	6.50±0,18
Yoghurt <sup>b</sup>	nd	6.40±0,07	6.41±0,18**	$6.60\pm0,12^{**}$


<sup>a</sup> Bb12 - *B. animalis* Bb 12 count in pure culture, St – *S. thermophilus* TA 040 count in pure culture and Lb – *L. bulgaricus* LB 340 count in pure culture. <sup>b</sup>*B. animalis* Bb 12 count in yoghurt after 30 days storage.\* The counts did not differ from the count in the control MRS agar at p < 0.05.\*\* The counts did not differ significantly from the count in RCPB pH5 at p < 0.05.ng – no growth.nd – not determined.

yoghurt at the conditions used in this study, so a partial inhibition of *S. thermophilus* in MRS-LP agar can be seen at Table 1 by the low count obtained in this agar (5.19 log CFU/mL).

The MRS-LP agar caused a reduction of approximately 1 log in the L. bulgaricus LB 340 count when compared to the recovery of this organism in RCPB pH5 agar, which is not an inhibitory media for L. bulgaricus and can be used for the differential count of this microorganism in yoghurt (10). Although the level of reduction was satisfactory for the selective count of Bifidobacterium in this study, it is known that the sensitivity to antimicrobial agents depends on the strain of microorganism used, and thus a reduction of only 1 logarithmic cycle in the number of Lactobacillus cells makes the use of MRS-LP agar risky, since this agar is selective and not differential and there is no difference between the colonies of Bifidobacterium and Lactobacillus strains in this medium. Surprisingly, the fortified RCPB pH5 agar also presented a reduction in the L. bulgaricus count, of about 2 logarithmic cycles, when this was inoculated as a pure culture (6.50 log CFU/mL) when compared to the RCPB pH5 agar (8.55 log CFU/mL).

In yoghurt, the *B. animalis* colonies developed presenting a mean diameter of about 1 to 2 mm on the fortified RCPB pH5 agar after refrigerated storage, although on the RCPB pH5 agar they could only be seen with the aid of a stereoscope being magnified at 15 to 20 times. On the MRS-LP agar, the *B. animalis* colonies presented a diameter of about 0.1 mm.

After refrigerated storage, the *B. animalis* counts in RCPB pH5, MRS-LP and in fortified RCPB pH5 agar was not significantly different at p<0.05 (Table 1). Although the RCPB pH5 agar showed good recovery of the cells after the storage period of the yoghurt, the *Bifidobacterium* colonies could only be visualised with the aid of a stereoscope, making the count difficult and prone to errors. Payne *et al.* (9) did not manage to



**Figure 1.** Petri dishes for counting *B. animalis* Bb 12 (narrow arrows, white colonies) after 30 days of storage at 7°C. White colonies of *L. bulgaricus* surrounded by blue halo (wide arrows). Left–RCPBpH5. Right–RCPB pH5-Liver: RCPB pH5 fortified with 150 mL/L of liver infusion.

distinguish the species of *Bifidobacterium* in relation to the yoghurt cultures, when mixed cultures were inoculated onto RCPB agar. The fortification of RCPB pH5 in this study with salts and a liver infusion overcame this problem, allowing the *Bifidobacterium* colonies to develop to a diameter of about 1 to 2 mm, making the counting and differentiation of the colonies much easier. The MRS-LP agar also presented selectivity and could be used for the count of *Bifidobacterium* during the refrigerated storage of yoghurt, although the final diameter of the colonies is much smaller than that found using the fortified RCPB pH5 (1 to 2 mm against 0.1 mm in the MRS-LP medium)

# Effect of the addition of salts in the fortification of RCPB pH5 agar

Table 2 shows the effect of the addition of salts (KH<sub>2</sub>PO<sub>4</sub>, K<sub>2</sub>HPO<sub>4</sub>, FeSO<sub>4</sub>7H<sub>2</sub>O, MnSO<sub>4</sub>H<sub>2</sub>O, MgSO<sub>4</sub>7H<sub>2</sub>O) on the performance of the fortified RCPB pH5 agar. There was no difference (p<0.05) in the *Bifidobacterium* counts between the two types of fortification (with and without salts addition), showing that the presence of the salts did not affect the performance of the fortified agar. Thus the addition of dehydrated liver extract alone was sufficient for the RCPB pH5 agar to considerably increase its power to differentiate the *Bifidobacterium* during refrigerated storage of the yoghurt, when these cells were more stressed due to the long time of exposure at acid pH of yoghurt and low temperatures of the refrigerated storage.

Concluding the three culture media evaluated (RCPB pH5, MRS-LP and fortified RCPB pH5) presented the capacity to count *B. animalis* Bb 12 in yoghurt after refrigerated storage. However, in RCPB pH5 agar, the *Bifidobacterium* colonies were extremely small and could only be counted with the aid of a stereoscope, making the count difficult and prone to errors. The MRS-LP agar showed good selectivity for *Bifidobacterium* and could be used to count this microorganism in yoghurt after refrigerated storage; however, the count of *Bifidobacterium* on this medium could be risky, as this medium only inhibited about 1 logarithmic cycle of *L. bulgaricus* cells when in pure culture. Fortified RCPB pH5 agar also presented the same cell recovery

**Table 2.** Bifidobacterium animalis subsp. lactis Bb 12 counts(log CFU/mL) in RCPB pH5 agar with diverse fortifications.

RCPB pH5	RCPB pH5 agar fortified with liver extract and salts*	RCPB pH5 agar fortified with liver extract (RCPB pH5- Liver agar)
$7.30 \pm 0.40$	$7.24 \pm 0.51$	$7.41 \pm 0.59$

\* KH<sub>2</sub>PO<sub>4</sub> (1 g/L), K<sub>2</sub>HPO<sub>4</sub> (1 g/L), FeSO<sub>4</sub>7H2O (0.01 g/L), MnSO<sub>4</sub>H<sub>2</sub>O (0.01 g/L) and MgSO<sub>4</sub>7H<sub>2</sub>O (0.2 g/L).

of *Bifidobacterium* as RCPB pH5 agar and its fortification with 150mL/L of dehydrated liver extract resulted in a significant increase in the diameter of the *Bifidobacterium* colonies, making the count easier and reliable after the refrigerated storage of yoghurt, when *Bifidobacterium* cells are more stressed.

## ACKNOWLEDGEMENTS

To Chr. Hansen Indústria and Comércio Ltda. for the donation of *Bifidobacterium animalis* Bb 12. To Rhodia Food (current Danisco Brasil Ltda, Cotia, SP) for donation of *Streptococcus thermophilus* TA 040 and *Lactobacillus delbrueckii subsp. bulgaricus* LB 340. To Wiesby Starter Cultures & Media (current Danisco Brasil Ltda, Cotia, SP) for donation of *S. thermophilus* A and *Lactobacillus delbrueckii subsp. bulgaricus* 18. To FAPESP for financial support.

#### **RESUMO**

# Avaliação de meios de cultura para contagem de Bifidobacterium animalis subsp. lactis Bb 12 em iogurte após a estocagem refrigerada

O meio RCPB pH5 tem sido considerado uma boa opção para a contagem de Bifidobacterium em iogurte. Entretanto, durante a estocagem refrigerada do iogurte é extremante difícil a contagem deste microrganismo devido ao pequeno diâmetro desenvolvido pelas colônias de Bifidobacterium neste meio, sendo que a sua contagem somente se torna possível com o auxílio de um estereoscópio. Outro meio, MRS-LP, também tem sido recomendado para a contagem de Bifidobaterium em iogurte. Este estudo avaliou a suplementação do meio RCPB pH5 com extrato de figado desidratado e com os sais KH<sub>2</sub>PO<sub>4</sub>, K<sub>2</sub>HPO<sub>4</sub>, FeSO<sub>4</sub>7H<sub>2</sub>O, MnSO<sub>4</sub>H<sub>2</sub>O e MgSO<sub>4</sub>7H<sub>2</sub>O, visando melhorar a diferenciação de Bifidobacterium em iogurte durante a estocagem refrigerada e também avaliou a contagem seletiva de Bifidobacterium em iogurte usando o meio MRS-LP. O meio MRS-LP apresentou a mesma recuperação de células que o meio RCPB pH5, usado como padrão, após 30 dias de estocagem refrigerada do iogurte, sendo considerado uma boa opção para a contagem de Bifidobacterium em iogurtes durante a estocagem

refrigerada. O meio RCPB pH5 fortificado também apresentou a mesma recuperação de células de *Bifidobacterium* que o meio padrão RCPB pH5; entretanto, a adição de extrato de figado desidratado aumentou consideravelmente o diâmetro das colônias de *Bifidobacterium*, tornando a diferenciação destas bastante fácil e confiável quando comparadas à sua diferenciação no meio RCPB pH5 sem a fortificação. A adição dos sais (KH<sub>2</sub>PO<sub>4</sub>, K<sub>2</sub>HPO<sub>4</sub>, FeSO<sub>4</sub>7H<sub>2</sub>O, MnSO<sub>4</sub>H<sub>2</sub>O e MgSO<sub>4</sub>7H<sub>2</sub>O) não exerceu influência no desempenho do meio RCPB pH5.

Palavras-chave: contagem, *Bifidobacterium*, iogurte, estocagem refrigerada

## REFERENCES

- Beehrens, H. (1990). An elective and selective isolation medium for Bifidobacterium spp. Lett. Appl. Microbiol., 11 (3), 155-157.
- Bomba, A.; Nemcová, R.; Mudronová, D.; Guba, P. (2002). The possibilities of potentiating the efficacy of probiotics. *Trends Food Sci. Technol.*, 13 (4), 121-126.
- FAO/WHO. Report of a Joint FAO/WHO Expert Consultation on Evaluation of Health and Nutritional Properties of Probiotics in Food Including Powder Milk with Live Lactic Acid Bacteria, Amerian Córdoba Park Hotel, Córdoba, Argentina, 1-4 October, 2001.
- 4. Ghodussi, H.B.; Robinson, R.K. (1996). Enumeration of starter cultures in fermented milks. J. Dairy Res., 63 (1), 151-158.
- 5. Lapierre, L.; Undeland, P.; Cox, L. J. (1992). Lithium chloridesodium propionate agar for the enumeration of Bifidobacteria in fermented dairy products. *J. Dairy Sci.*, 75 (5), 1192-1196.
- 6. Laroia, S.; Martin, J.H. (1991). Methods for enumerating and propagating *Bifidobactéria*. *Cultured Dairy Prod. J.*, 26 (2), 32-33.
- Moriya, J.; Fachin, L.; Gândara, A.L.N.; Viotto, W.H. (2006). Evaluation of culture media for counts of *Bifidobacterium animalis* in the presence of yoghurt bacteria. *Braz. J. Microbiol.*, 37, 510-514.
- Onggo, I.; Fleet, G.H. (1993). Media for the isolation and enumeration of lactic acid bacteria from yoghurts. *Aust. J. Dairy Technol.*, 48 (2), 89-92.
- Payne, J.F.; Morris, A.E.J.; Beers, P. (1999). Evaluation of selective media for the enumeration of *Bifidobacterium* sp. in milk. *J. Appl. Microbiol.*, 86 (2), 353-358.
- Rybka, S.; Kailasapathy, K. (1996). Media for the enumeration of yoghurt Bacteria. *Int. Dairy J.*, 6 (8-9), 839-850.
- Vinderola, C.G.; Reinheimer, J.A. (1999). Culture media for the enumeration of *Bifidobacterium bifidum* and *Lactobacillus acidophilus* in the presence of yoghurt bacteria. *Int. Dairy J.*, 9 (8), 497-505.