

## Monoclonal Antibodies

### 1B2.1E12 anti-TNP

#### ANTIGEN USED FOR IMMUNIZATION

OVA-TNP conjugate: Type VI ovalbumin (OVA) and 2,4,6-trinitrobenzenesulfonic acid (TNBS), both from Sigma Chemical Co. (St. Louis, MO), were conjugated as described.<sup>(1)</sup>

#### METHOD OF IMMUNIZATION

BALB/c mice were immunized intraperitoneally with 100  $\mu$ g of TNP-OVA conjugate and complete Freund's adjuvant (1:1, v/v). On Days 20 and 30, the animals were boosted with the same dose of antigen diluted in 0.15 M saline. Three days later, immunized mice were sacrificed in a CO<sub>2</sub> chamber, and their spleen cells were aseptically collected for fusion experiment.

#### PARENTAL CELL LINE USED FOR FUSION

SP2/0 Ag-14

#### SELECTION AND CLONING PROCEDURES

Spleen cells of three immunized mice were pooled and then fused with SP2/0 Ag-14 myeloma cells (4:1 ratio) by mixing them with polyethylene glycol:dimethylsulfoxide (PEG: DMSO; 50%:10%) as described.<sup>(2)</sup> After fusion, the cells were suspended in HAT medium, dispensed into each well of two 24-well flat bottomed tissue culture plates containing macrophage feeder layers, and then cultured at 37°C in a humid atmosphere of 5% CO<sub>2</sub> in air. The hybridoma cells producing anti-TNP antibodies, as detected by ELISA using TNP-BSA, were cloned by limiting dilution in 96-well microplates (Corning) containing macrophage feeder layers.

#### HEAVY AND LIGHT CHAINS OF IMMUNOGLOBULIN

IgG<sub>1k</sub>

#### SPECIFICITY

The supernatants and ascites obtained from the clone 1B2.1E12 were assayed for their reactivity with TNP-OVA, TNP-BSA, OVA, and BSA. These materials gave positive reactions only when TNP-BSA was used as antigen in ELISA tests.

#### SPECIFIC ANTIGEN IDENTIFIED

TNP

#### AVAILABILITY

Tissue culture supernatant	Yes ✓	No
Ascitic fluid	Yes ✓	No
Hybridoma cells	Yes ✓	No

#### ADDRESS CORRESPONDENCE TO:

Wirla M.S.C. Tamashiro  
Departamento de Microbiologia e Imunologia  
Instituto de Biologia  
Universidade Estadual de Campinas  
CP 6109, 13081-970 Campinas, SP  
Brazil  
Telephone: +55-19-7887911  
Fax: +55-19-7888190  
E-mail: wirlatam@unicamp.br

#### REFERENCES

1. Rittenberg MB, and Amkraut AA: Immunogenicity of trinitrophenyl-hemocyanin: Production of primary and secondary anti-haptenprecipitins. *J Immunol* 1966;97:421-430.
2. de St. Groth FS, and Scheidegger D: Production of monoclonal antibodies: Strategy and tactics. *J Immunol Methods* 1980;35:1-21.