## Repository of the Max Delbrück Center for Molecular Medicine (MDC) in the Helmholtz Association

http://edoc.mdc-berlin.de/16056

## Efficient CRISPR-mediated mutagenesis in primary immune cells using CrispRGold and a C57BL/6 Cas9 transgenic mouse line

Chu, V.T. and Graf, R. and Wirtz, T. and Weber, T. and Favret, J. and Li, X. and Petsch, K. and Tran, N.T. and Sieweke, M.H. and Berek, C. and Kuehn, R. and Rajewsky, K.

This is the final version of the accepted manuscript. The original article has been published in final edited form in:

Proceedings of the National Academy of Sciences of the United States of America 2016 NOV 01 ; 113(44): 12514-12519 2016 OCT 11 (first published online) doi: 10.1073/pnas.1613884113

Publisher: National Academy of Science

© 2016 The Author(s)

Biological sciences: Immunology and inflammation

# Efficient CRISPR-mediated mutagenesis in primary immune cells using CrispRGold and a new Cas9 transgenic mouse line

#### Efficient CRISPR-based mutagenesis in primary B cells

Van Trung Chu<sup>a,1,2</sup>, Robin Graf<sup>a,1,2</sup>, Tristan Wirtz<sup>a</sup>, Timm Weber<sup>a</sup>, Jeremy Favret<sup>a,b</sup>, Xun Li<sup>a</sup>, Kerstin Petsch<sup>a</sup>, Ngoc Tung Tran<sup>a</sup>, Michael H. Sieweke<sup>a,b,c,d</sup>, Claudia Berek<sup>e</sup>, Ralf Kuehn<sup>a,f</sup> and Klaus Rajewsky<sup>a,2</sup>.

<sup>a</sup>Max-Delbrück-Center for Molecular Medicine, 13125 Berlin, Germany
<sup>b</sup>Centre d'Immunologie de Marseille-Luminy, Université Aix-Marseille, UM2, Campus de Luminy, Case 906, 13288 Marseille Cedex 09, France.
<sup>c</sup>Institut National de la Santé et de la Recherche Médicale (INSERM), U1104, Marseille, France
<sup>d</sup>Centre National de la Recherche Scientifique (CNRS), UMR7280, Marseille, France
<sup>e</sup>Deutsches Rheuma-Forschungszentrum (DRFZ), 10117 Berlin, Germany
<sup>f</sup>Berlin Institute of Health, Kapelle-Ufer 2, 10117 Berlin, Germany
<sup>1</sup>These authors contributed equally to this work
<sup>2</sup>Correspondence should be addressed to K.R (Klaus.Rajewsky@mdcberlin.de), V.T.C (Vantrung.Chu@mdc-berlin.de) and R.G (Robin.Graf@mdcberlin.de)

Keywords: CRISPR/Cas9, sgRNA design, primary cells, B cells, knockout efficiency

#### Abstract

Applying CRISPR/Cas9-mediated mutagenesis to primary mouse immune cells, we used high-fidelity sgRNAs designed with a novel sgRNA design tool (CrispRGold) to target genes in primary B cells, T cells and macrophages isolated from a new Cas9 transgenic mouse line. Using this system, we achieved an average knockout efficiency of 80% in B cells. On this basis, we established a robust small-scale CRISPRmediated screen in these cells and identified new genes essential for B cell activation and plasma cell differentiation. This screening system does not require deep sequencing and may serve as a precedent for the application of CRISPR/Cas9 to primary mouse cells.

#### Significance

Efficient CRISPR/Cas9-mediated mutagenesis is necessary for robust genetic screening in primary cells and requires sufficiently high levels of Cas9 and reliable sgRNAs. We provide a new sgRNA design tool that selects high-fidelity sgRNAs and a new Cas9 transgenic mouse line that expresses Cas9 at high levels in all cells analyzed. Using this system, we achieved an average knockout efficiency of 80% in primary B cells and established a robust CRSPR/Cas9-mediated genetic screening system that can be adapted to other primary cell types. The CRISPR/Cas9 technology is a powerful tool for gene editing (1-3). Application of this technology to primary cells is still in its initial phase, but promises to accelerate research. In many instances, such as the primary hematopoietic cells studied here, high knockout efficiencies are essential for studying gene function, as the limited life span of the cells often does not allow for selection and expansion of knockout clones prior to in vitro experiments. Moreover, screens that are based on readouts that do not affect survival require high knockout efficiencies to identify true positive hits. The two most decisive factors for high CRISPR/Cas9 knockout efficiencies are sufficient nuclear Cas9 levels and the selection of potent sgRNAs (4, 5). To overcome the requirement of Cas9 delivery to primary cells, transgenic mice expressing Cas9 and GFP from the Rosa26 locus have been generated (6; a related approach, targeting the hematopoietic system, has been lentiviral Cas9 transduction into hematopoietic stem cells, ref. 7). While primary dendritic cells of these mice have been successfully used for functional screening in vitro (6, 8), their lymphocytes express only low levels of GFP, and their use in CRISPR-mediated screens has not been reported so far. In addition, these mice were not generated on a pure C57BL/6 background. We thus generated a new C57BL/6 Cas9-transgenic line with an improved Cas9 expression cassette (9). As we were interested in small-scale screening, a key consideration was to use only few, but reliable sgRNAs without the necessity of sgRNA-testing experiments prior to screening. Previous studies have shown large variations in sgRNA efficiencies (8, 10-12). Multiple sgRNA design programs have been developed, but none of them combines all the criteria that we developed in-lab to reach optimal sgRNA activity and specificity, including on-target efficiency, genome-wide off-target predictions taking into account the impact of mismatch-type and -position, targeting of all isoforms of genes, and conservation of functional secondary structure of the sgRNA (13). We thus developed a new sgRNA design tool that selects sgRNAs having high specificity, high on-target efficiency and intact secondary structure important for the recognition by Cas9. Combining this new program with our Cas9 transgenic C57BL/6 mice, we established an easy and robust screening method to study gene function in primary hematopoietic (and likely other) cells with high signal-to-noise ratio.

#### Results

#### Design of specific and efficient sgRNAs

Coding genes typically contain hundreds of protospacer adjacent motif (PAM) sequences, allowing for a wide choice among possible sgRNAs. To design the most specific and efficient sgRNAs to target a gene or set of genes, we implemented our sgRNA design rules in a new tool that we termed CrispRGold. CrispRGold first projects the coding sequences (CDS) of the different isoforms of a gene to a single 'minimal CDS', in order to identify sgRNAs that would target most or all isoforms of that gene. For each sgRNA candidate, the position in the minimal CDS is stored and the sgRNA is aligned to the genome using the Burrows-Wheeler Aligner (BWA) (14) to identify offtarget sites. Based on the quantitative measurement of off-target site activity (15), we defined a type- and position-dependent mismatch-penalty matrix, such that the sum of mismatch-penalties (or mutation distance) for the various off-target sites would inversely correlate with the frequency of detected offtarget site activities (Fig. 1A and B). This mismatch penalty matrix especially weighs the increased impact of mismatches close to the 3'-end of the sgRNA and the decreased relevance of mutations causing wobble base pairing. CrispRGold uses this matrix to calculate the risk of each individual predicted off-target site genome-wide for a given sgRNA, sorts these sites according to their risk, and attributes a global off-target risk (or specificity score) to each sgRNA based on the three off-target sites with the highest risk. In addition, CrispRGold computes sgRNA-intrinsic properties, such as the GC-content and the folding energy between the targeting sequence and the scaffold-RNA of the sgRNA (16). A high binding energy between the targeting sequence and the scaffold may alter the secondary structure of the sgRNA and potentially disturb recognition of the sgRNA by Cas9 or binding of the targeting sequence to the DNA (Fig. 1C). Finally, CrispRGold goes through iterative loops of decreasing stringency to find the required number of sgRNAs per gene or sequence. In each loop, sgRNAs are prioritized according to their specificity (Fig. 1D).

#### Characterization of a new Cas9 transgenic mouse line

In the mice created by Platt et al. (6) (hereafter called R26-Cas9p2aGFP), the Cas9 protein is fused to eGFP via a self-cleaving P2A peptide. In our Rosa26-LSL-Cas9iGFP mice, generated by CRISPR-mediated editing of C57BL/6 zygotes, Cas9 is also expressed from the Rosa26 locus in a Cre-dependent manner, but with the Cas9 coding sequence linked to eGFP via an internal ribosomal entry site (IRES) (9). To apply the CRISPR/Cas9 technology to primary immune cells, we generated Cas9-transgenic mice with ubiguitous expression of Cas9 by crossing Rosa26-LSL-Cas9iGFP with Cre-deleter mice (Fig. 1E). We then evaluated the expression of Cas9 in the resulting R26-Cas9iGFP/+ animals by monitoring GFP reporter levels by FACS. All hematopoietic sub-populations analyzed exhibited high GFP levels, including hematopoietic stem cells (HSC) in the bone marrow (BM), and various hematopoietic lineages, such as T, B and myeloid cells, independent of their anatomical location (Fig. 1F,G and Fig. S1A,B and C). The percentages of HSCs, B cells, T cells and myeloid cells were similar between R26-Cas9iGFP/+ and control mice, indicating that the expression of Cas9iGFP is not toxic to any of these cell types (Fig. 1H). GFP levels in R26-Cas9iGFP/+ mice were higher than in R26-Cas9p2aGFP/+ mice (6), likely due to the presence of the IRES instead of the P2A peptide between Cas9 and GFP (Fig. S2A and B). To confirm that Cas9 is indeed present in the GFP<sup>+</sup> cells, we performed western blotting on splenic B cells that were activated in vitro and on BM-derived macrophages (BMDM). In both cases, we detected high levels of Cas9 protein (Fig. S3A). Moreover, Cas9 was localized in the nucleus in the BMDM (Fig. S3B). Thus, primary immune cells of R26-Cas9iGFP/+ animals express high levels of GFP and nuclear Cas9 to nontoxic levels, allowing for genome editing and easy discrimination of Cas9<sup>+</sup> and Cas9<sup>-</sup> cells in subsequent experiments.

#### Efficient CRISPR/Cas9-mediated mutagenesis in primary immune cells

To assess the knockout efficiency of the sgRNAs designed by CrispRGold in primary immune cells, we first selected 12 genes encoding B cell surface markers and designed 2 sgRNAs per gene (Fig. S4A). These sgRNAs were individually cloned into retroviruses bearing a BFP reporter and a Puromycin resistance gene. We then isolated B cells from wildtype and R26-Cas9iGFP/+ mice and mixed them in a 1:4 ratio, such that we could use the Cas9<sup>-</sup> (GFP<sup>-</sup>) cells as internal controls. These mixed B cells were activated with anti-CD40 antibodies and IL-4 for 2 days and then transduced with retroviral particles. The transduced cells were further cultured on 40LB feeder cells (17) (NIH3/T3 cells stably expressing CD40 ligand and BAFF) in the presence of IL-21 and Puromycin (Fig. 2A). Four days after transduction, all sgRNAs had led to cell surface marker knockouts, mostly with very high efficiency and an average knockout frequency above 70% (Fig. 2B,C and Fig. S4B). Only one sgRNA led to significantly lower knockout efficiencies due to a previously unknown sequence peculiarity (R.G., unpublished data), now implemented in CrispRGold. To confirm that these high knockout efficiencies are not unique to B cells, we performed a similar experiment using primary T cells that we activated with anti-CD3 and anti-CD28 (Fig. S5A). Four days after transduction with retroviral particles encoding sgRNAs against CD44, about 50% of the cells showed a complete knockout of CD44 (Fig. 2D). Finally, to determine if we reach similarly high knockout frequencies in quiescent primary myeloid cells, we differentiated bone marrow cells into resting BMDMs using M-CSF and transduced them with lentiviral vectors expressing sgRNAs targeting CD64 and CD14 (Fig. S5B). Four days after transduction, 40-80% of the cells showed a knockout of the targeted genes (Fig. 2E and Fig. S5C). Thus, primary immune cells of R26-Cas9iGFP/+ mice are suitable for CRISPR/Cas9-mediated genome editing, and CrispRGold designs reliably highly efficient sgRNAs. In addition to the intrinsic properties of sgRNAs, the on-target dosage of sgRNA/Cas9 complexes is essential for efficient mutagenesis. The lower GFP levels in lymphocytes of R26-Cas9p2aGFP/+ mice (6) suggested that they have lower levels of Cas9 than the R26-Cas9iGFP/+ mice here described. To compare the impact of Cas9-dosage on knockout efficiency, we thus repeated the experiment depicted in Fig. 2A,

using the same sgRNAs and B cells of both Cas9-transgenic lines side by side. Indeed, all sgRNAs tested led to significantly lower knockout efficiencies in B cells from R26-Cas9p2aGFP/+ mice compared to R26-Cas9iGFP/+ mice (Fig. 2*F*,*G* and Fig. S3*C*), consistent with lower Cas9 levels in these cells (Fig. S3*D*). Thus, the combination of optimal sgRNAs and high levels of Cas9 in primary cells of R26-Cas9iGFP/+ mice led to the high knockout efficiencies observed in primary cells isolated from these animals.

## Robust inactivation of transcription factors known to be important for B cell differentiation

Extending the analysis to transcription factors (TFs), we selected 5 TFs, Prdm1, Irf4, Myc, Xbp1, and Pou2af1, that are known to be important for B cell survival, proliferation and differentiation (18-24) (Fig. S6A). We designed 3 sgRNAs per TF (Fig. S6B) and performed retroviral transduction experiments as before, but this time analyzing B cell survival and plasma cell differentiation (Fig. 3A). Expression of sgRNAs targeting Xbp1, Irf4, Pou2af1 and Myc led to a strong survival disadvantage of Cas9-expressing cells, indicated by the decreased percentage of  $GFP^+$  cells (Fig. 3B,C and Fig. S6C). As for plasma cell differentiation, monitored by CD138 expression, targeting of all 5 TFs led to a strong block of differentiation (Fig. 3D, E and Fig. S6D). These findings are consistent with previous reports, showing that Xbp1, Irf4, Pou2af1 and Myc are important for both B cell survival, proliferation and terminal differentiation (18-22), whereas Prdm1 is only important for the latter (23, 24). Moreover, in all cases, the 3 individual sgRNAs showed a strong and consistent effect on the biological readout, further demonstrating that sgRNAs designed by CrispRGold work with high efficiency and consistency.

## Identification of novel genes important for B cell activation and plasma cell differentiation

Benefitting from the robustness of our system and the consistently high knockout frequencies, we performed a small-scale screen to identify novel genes important for B cell proliferation and terminal differentiation. Based on a microarray analysis of B cells and plasma cells isolated from immunized mice, we selected 83 candidate genes up-regulated during plasma cell

differentiation in vivo (Fig. S7). We designed one sgRNA against each of these genes and 13 control genes. We then performed the same experiment as before in a 96-well system, using the percentage of GFP<sup>+</sup> and CD138<sup>+</sup> cells as readout for survival/proliferation and plasma cell differentiation, respectively (Fig. 3A and Fig. S8A). Based on two experiments, which led to highly consistent results, we identified 22 genes to have a strong impact on survival/proliferation and 8 genes to selectively affect plasma cell differentiation (Fig. 3F,G and Fig. S8B). The 22 genes important for B cell survival/proliferation contained genes whose function in this context was previously unknown, such as Egr2, Dis3, Ost4, Preb and Pomp, which we validated in further experiments (Fig. S9A). Dis3 is potentially involved in immunoglobulin class switch recombination via targeting AID (25), whereas Pomp might be involved in plasma cell differentiation (26). Furthermore, we identified Arf4, Creld2 and Zfp36 among the novel genes enhancing or blocking plasma cell differentiation (Fig. 3H and Fig. S9B). Of note, mice deficient for Zfp36 have been shown earlier to develop autoimmune disease, a finding that could connect to our observation of enhanced plasma cell differentiation in its absence (27). These results show that the screening system as described here leads to clear and consistent functional results, permitting small-scale screens in primary mouse cells without the need of high numbers of sgRNAs per gene or deep sequencing.

#### Discussion

We provide a new and flexible sgRNA design tool that reliably identifies highly specific and efficient sgRNAs for a given set of genes or DNA sequences. CrispRGold can also be used for other organism, and will be available on a web platform. Primary cells of R26-Cas9iGFP/+ mice express high levels of GFP and Cas9 to non-toxic levels, allowing for efficient CRISPR-mediated gene editing and the use of GFP<sup>-</sup> internal controls. We show that these tools can be combined in an easy and robust screening assay of CRISPR mutagenesis in primary hematopoietic cells, and identify novel genes important for B cell activation and terminal differentiation. The same approach can be easily applied to other types of primary cells.

#### Methods

#### Mice

Cas9 knock-in mice were reported previously (9). Briefly, we knocked a Cas9-IRES-GFP cassette preceded by a LoxP-flanked Stop cassette into the mouse Rosa26 locus. The LSL-Cas9 mice were bred with Cre-deleter mice to generate Cas9 transgenic mice (available on Jackson Laboratory, stock No. 028555). The C57BL/6 mice were obtained by local breeding. The heterozygous R26-Cas9p2aGFP/+ mice were purchased from the Jackson Laboratory (Stock No: 026179). BALB/c mice were from Charles River Laboratories. For the microarray of *in vivo* generated plasma cells, BALB/c mice were immunized intraperitoneally with 100µg phOx coupled to chicken serum albumin (CSA), precipitated with aluminum hydroxide. After 6–8 weeks, mice were boosted intravenously with 100µg soluble antigen. Plasma cells were sorted from spleen and the BM at day 6 and day 60 after secondary immunization. Animal experiments were approved by the Institution Animal Care and use Committee.

#### CrispRGold

Minimal CDS (minCDS) are generated based on the refseq tables (mm9). sgRNA candidates are defined by the presence of NGG and cutting site inside the minCDS. sgRNA candidates are mapped to the mouse genome (mm9) using BWA (14) (version 0.7.12) to identify potential off-target sites, ignoring the first base of the targeting sequence and allowing for 4 mismatches. The folding energy between the targeting sequence and the scaffold RNA is calculated using RNAduplex (Vienna Package) (16). CrispRGold is written in Perl, as all scripts that were implemented in CrispRGold if not otherwise stated. Briefly, sgRNA candidates are processed and scored as described in the main text. The first loop considers sgRNAs within the first 45% of the minCDS, with  $T_m <= 60^{\circ}$ C, the lowest off-target risk score > 11, a scaffold-folding energy <= 20 kcal/mol, targeting the max N isoforms, without low-efficiency features and distance to the CDS-start >= 50 nts. The second loop considers sgRNAs as the first loop, but within the first 60% and with the lowest off-target risk score > 6. The third loop considers sgRNAs as the

second loop, but with  $T_m \le 65^{\circ}C$  and distance to CDS-start  $\ge 10$  nts. The fourth loop considers sgRNAs as the third loop, but with distance to the CDS-start  $\ge 1$  nt and neglecting  $T_m$ , scaffold-folding energy and low-efficiency features. The last loop considers sgRNAs as the fourth loop, but extending the search space to 90% off the minCDS.

#### 96-well cloning approach

The MSCV hU6 CcdB PGK Puro T2A BFP vector was generated by cloning the PCR-amplified hU6-BbsI-CcdB-BbsI-gRNA fragment into the Sall and Xhol sites of the MSCV vector. The PGK-Puromycin-T2A-BFP fragment was amplified by overlapping PCR and cloned into the Mlul site of the MSCVhU6-BbsI-CcdB-BbsI-gRNA vector. For generating the mini-library, forward and reverse oligos were separately ordered in 96 deep-well plates. Each forward and reverse oligo was mixed and phosphorylated individually. Then annealed oligo duplexes were cloned into the Bbsl sites of the MSCV U6 CcdB PGK Puro T2A BFP vector. The plasmids were transformed into DH5 $\alpha$  bacteria using a heat-shock 96-well system. After 30 min pre-culture at 37°C, the transformed bacteria were transferred into 96 deep-well plates containing 1.5ml LB liquid medium and sealed with PCR seals (Thermo Scientific). These plates were cultured for 12 h then split into 2 new 96 deep-well plates and further cultured for 10-12h. Bacteria were collected by centrifugation at 4000 rpm for 1 min and plasmids were isolated using the NucleoSpin® 96 plasmid core kit (Macherey-Nagel).

#### Cell culture

Retroviral Plat-E packaging cells were maintained in DMEM (Gibco) supplied with 10% FCS (Gibco), 2mM L-Glutamine (Gibco), and 2mM Sodium Pyruvate (Gibco). 40LB feeder cells, producing BAFF and CD40L, were previously generated by Nojima et al(17) and maintained in completed DMEM medium. To prepare the feeder layer, 40LB feeder cells were irradiated with  $12G\gamma$  and plated at  $5x10^4$  cells/cm. Naïve B cells were isolated from the spleen of R26-Cas9iGFP/+, R26-Cas9p2aGFP/+ or C57BL/6 mice by depletion of CD43<sup>+</sup> cells using CD43 microbeads (Miltenyi Biotec). Resting B cells were plated at

 $10^6$  cells/ml in DMEM (Gibco) supplied with 10% FCS (Gibco), 2mM L-Glutamine, 2mM Sodium Pyruvate, 2mM HEPES (Gibco), 1x NAA (Gibco), beta-mercaptoethanol (Sigma), and  $10\mu$ g/ml Gentamicin (Lonza). Cas9 (GFP<sup>+</sup>)-expressing B cells were mixed wildtype (GFP<sup>-</sup>) B cells and stimulated with anti-CD40 antibody (Biolegend) and IL-4 (Perprotech) for 48h before retroviral transduction.

#### High-throughput retroviral production

Plat-E cells were plated at  $2.2 \times 10^4$  cells/100µl medium in each well of 96-well plates at day 0. 1 day later, the cells were transfected with 100ng DNA using Fugene®HD transfection reagents (Promega). 12 h after transfection, the supernatant was discarded and the cells were cultured with 150µl new complete medium. The transfected cells were cultured in a 32°C incubator for 24 h. The retroviral supernatant was collected 48 h and 72 h after transfection into 96 deep-well plates. The 96-deep-well plates with supernatant were centrifuged at 1300 rpm for 5 min and the supernatant was transferred into new 96 deep-well plates and stored at -80°C.

#### High-throughput retroviral transduction

48 h after stimulation,  $GFP^+/GFP^-$  mixed B cells were harvested and resuspended at a density of  $2x10^5$  cells/ml in completed B cells medium supplied with  $16\mu$ g/ml Polybrene (Sigma).  $100\mu$ l cell suspension were transferred into each well of 96-well plates and  $100\mu$ l of retroviral supernatant were added, then activated B cells were spin-transduced at 2500rpm for 90min. The 96-well plates were placed in a 37 °C incubator and on then next day the transduced  $GFP^+/GFP^-$  B cells were transferred into 96-well 40LB feeder cell plates in the presence of 20ng/ml of IL-21 (Perprotech) and were selected with  $1.25\mu$ g/ml Puromycin (Sigma). The transduced  $GFP^+/GFP^-$  B cells were harvested and analysed at different time points by using BD Fortessa with a 96-well HTS unit.

#### T cell activation and transduction

Total T cells were isolated from R26-Cas9iGFP/+ and C57BL/6 mice using Pan T Cell Isolation Kit II, mouse (Miltenyi Biotec) and subsequently CD4 and CD8 T cells were sorted by FACS, the Cas9 (GFP<sup>+</sup>)-expressing T cells were mixed with wild-type (GFP<sup>-</sup>) T cells and activated with plate-bound anti-CD3 (145-2C11, Biolegend) and anti-CD28 (37.51, Biolegend) in the presence of 25 ng/ml rmIL-2 (Peprotech) for 2 days. The activated T cells were spin-transduced with retroviral particles expressing sgRNAs. The transduced T cells were selected with Puromycin for 4 days. After Puromycin selection, the transduced T cells were analysed by flow cytometry or sorted.

#### **BMDM** generation and lentiviral transduction

BM-derived macrophages were differentiated from total mouse bone marrow cells for 14 days in macrophage growth medium (Dubelco's modified medium (DMEM) with 10% FCS, 20% L-cell supernatant, 1% sodium pyruvate, 1% L-Glutamine). sgRNAs targeting CD64 and CD14 surface makers were cloned into BbsI sites of pLKV-U6gRNA(BbsI)-PGKPuroP2ABFP vector (Addgene, #50946). Lentiviral particles were produced at the Centre International de Recherche en Infectiologie (CIRI) U1111/UMR5308 Inserm-CNRS-UCBL-ENS de Lyon (Lyon, France). BMDMs were seeded at a density of  $2x10^5$  cells per well of 24 well plates and transduced with lentiviral supernatant for 4 hours. The medium was then replaced and 2 days later transduced cells were selected with 2 µg/ml Puromycin for 4 days. The transduced BMDMs were analysed by FACS.

#### FACS analysis and sorting

For FACS analysis, single cell suspensions were prepared from the bone marrow, spleen, peritoneal cavity and mLN from R26-Cas9iGFP/+, R26-Cas9p2aGFP/+ or C57BL/6 mice and cells were blocked with FcγR antibody (Biolegend) for 10 min. The surface antigens were stained with fluorescent-conjugated antibodies for 15 min. The cells were washed with FACS buffer (PBS/1%BSA) and analysed by BD Fortessa. The data was analysed by FlowJo software.

For cell sorting, the BFP<sup>+</sup> (sgRNA) GFP<sup>+</sup> (Cas9) cells were sorted into 15ml Falcon tubes and centrifuged prior to DNA isolation for further analaysis. For single cell sequencing analysis, single BFP<sup>+</sup>GFP<sup>+</sup> cells were sorted into 96-well plates containing  $5\mu$ l of DNA quick extraction buffer (Epicenter). The plates were briefly centrifuged and DNA was denatured and PCR was performed.

#### Immunofluorescent staining and western blotting

1x10<sup>5</sup> BMDMs were transferred into lysine-coated coverslip and further cultured for 24 hours. The cells were washed with cold PBS and fixed with 4% of Paraformaldehyde for 20 min at room temperature (RT). The cells were washed and permeabilized with PBS/0.25% TritonX-100 for 10 min at RT. The cells were washed and stained overnight with anti-Flag tag M2 antibody at 4°C. Anti-Flag tag antibody was developed with rat anti-mouse IgG1 PE (Biolegend) for 1h at room temperature. The nuclei were counterstained with DAPI (Sigma). The images were analyzed by Keyence fluorescent microscopy. To detect Cas9 protein, activated B cell and BMDM lysates were run on SDS-PAGE gel and transferred into PVDF membrane (Healthcare). The Cas9 protein was developed with anti-Flag tag and anti-Cas9 antibodies;  $\beta$ -actin was used for loading controls.

#### Acknowledgements

We thank H.P. Rahn for excellent FACS-related support. This work was supported by the European Research Council (ERC Advanced Grant 268921, to K.R.), the German Ministry of Education and Research within the VIP program (TAL-CUT 03V0261, to R.K.) and a BIH-Einstein fellowship to M.H.S. M.H.S is INSERM-Helmholtz group leader and received grants from the Agence Nationale de la Recherche (ANR-11-BSV3-026-01) and InCA (13-10/405/AB-LC-HS).

#### **Author contributions**

V.T.C. R.G. and K.R. designed the project, V.T.C. R.G. T.W. J.F. T.D.W. X.L. K.P. and N.T.T performed experiments and acquired the data, V.T.C. R.G. and T.W. analysed and interpreted the data, R.G. wrote CrispRGold and did the computational work, M.H.S. C.B. and R.K. provided materials and advice, V.T.C. R.G. T.W. T.D.W and K.R. wrote the paper.

#### **Competing financial interests**

The authors declare no competing financial interests.

#### References

- 1. Mali P, *et al.* (2013) RNA-guided human genome engineering via Cas9. *Science* 339(6121):823-826.
- 2. Cong L, *et al.* (2013) Multiplex genome engineering using CRISPR/Cas systems. *Science* 339(6121):819-823.
- 3. Hsu PD, *et al.* (2013) DNA targeting specificity of RNA-guided Cas9 nucleases. *Nature biotechnology* 31(9):827-832.
- Doench JG, *et al.* (2014) Rational design of highly active sgRNAs for CRISPR-Cas9-mediated gene inactivation. *Nature biotechnology* 32(12):1262-1267.
- 5. Osakabe Y, *et al.* (2016) Optimization of CRISPR/Cas9 genome editing to modify abiotic stress responses in plants. *Sci Rep* 6:26685.
- 6. Platt RJ, *et al.* (2014) CRISPR-Cas9 knockin mice for genome editing and cancer modeling. *Cell* 159(2):440-455.
- 7. Aubrey BJ, *et al.* (2015) An inducible lentiviral guide RNA platform enables the identification of tumor-essential genes and tumor-promoting mutations in vivo. *Cell Rep* 10(8):1422-1432.
- 8. Parnas O, *et al.* (2015) A Genome-wide CRISPR Screen in Primary Immune Cells to Dissect Regulatory Networks. *Cell* 162(3):675-686.
- 9. Chu VT, *et al.* (2016) Efficient generation of Rosa26 knock-in mice using CRISPR/Cas9 in C57BL/6 zygotes. *BMC Biotechnol* 16(1):4.
- Koike-Yusa H, Li Y, Tan EP, Velasco-Herrera Mdel C, & Yusa K (2014) Genome-wide recessive genetic screening in mammalian cells with a lentiviral CRISPR-guide RNA library. *Nature biotechnology* 32(3):267-273.
- 11. Shalem O, *et al.* (2014) Genome-scale CRISPR-Cas9 knockout screening in human cells. *Science* 343(6166):84-87.
- 12. Wang T, Wei JJ, Sabatini DM, & Lander ES (2014) Genetic screens in human cells using the CRISPR-Cas9 system. *Science* 343(6166):80-84.
- 13. Graham DB & Root DE (2015) Resources for the design of CRISPR gene editing experiments. *Genome biology* 16:260.
- 14. Li H & Durbin R (2009) Fast and accurate short read alignment with Burrows-Wheeler transform. *Bioinformatics* 25(14):1754-1760.

- Tsai SQ, et al. (2015) GUIDE-seq enables genome-wide profiling of offtarget cleavage by CRISPR-Cas nucleases. *Nature biotechnology* 33(2):187-197.
- 16. Lorenz R, et al. (2011) ViennaRNA Package 2.0. Algorithms Mol Biol 6:26.
- 17. Nojima T, *et al.* (2011) In-vitro derived germinal centre B cells differentially generate memory B or plasma cells in vivo. *Nature communications* 2:465.
- Schubart DB, Rolink A, Kosco-Vilbois MH, Botteri F, & Matthias P (1996)
   B-cell-specific coactivator OBF-1/OCA-B/Bob1 required for immune response and germinal centre formation. *Nature* 383(6600):538-542.
- Lin KI, Lin Y, & Calame K (2000) Repression of c-myc is necessary but not sufficient for terminal differentiation of B lymphocytes in vitro. *Mol Cell Biol* 20(23):8684-8695.
- 20. Reimold AM, *et al.* (2001) Plasma cell differentiation requires the transcription factor XBP-1. *Nature* 412(6844):300-307.
- 21. Klein U, *et al.* (2006) Transcription factor IRF4 controls plasma cell differentiation and class-switch recombination. *Nature immunology* 7(7):773-782.
- 22. Corcoran LM, *et al.* (2005) Differential requirement for OBF-1 during antibody-secreting cell differentiation. *The Journal of experimental medicine* 201(9):1385-1396.
- Shapiro-Shelef M, et al. (2003) Blimp-1 is required for the formation of immunoglobulin secreting plasma cells and pre-plasma memory B cells. *Immunity* 19(4):607-620.
- Lin Y, Wong K, & Calame K (1997) Repression of c-myc transcription by Blimp-1, an inducer of terminal B cell differentiation. *Science* 276(5312):596-599.
- Basu U, et al. (2011) The RNA exosome targets the AID cytidine deaminase to both strands of transcribed duplex DNA substrates. *Cell* 144(3):353-363.
- Zhang X, et al. (2015) MicroRNA-101 Suppresses Tumor Cell Proliferation by Acting as an Endogenous Proteasome Inhibitor via Targeting the Proteasome Assembly Factor POMP. *Mol Cell* 59(2):243-257.

27. Taylor GA, *et al.* (1996) A pathogenetic role for TNF alpha in the syndrome of cachexia, arthritis, and autoimmunity resulting from tristetraprolin (TTP) deficiency. *Immunity* 4(5):445-454.

#### **Figure legends**

#### Fig. 1. Design of efficient sgRNAs and characterization of new Cas9transgenic animals.

(*A*) Heatmap showing mismatch penalties depending on the position within the targeting sequence of the sgRNA and the type of mismatch. Scores were chosen to fit the off-target frequencies detected by Tsai et al.

(*B*) Dotplot of the predicted inverse offtarget risk of each off-target site of the VEGFA\_3 sgRNA versus the reported GUIDEseq read numbers. Read numbers were extracted from Tsai et al.

(*C*) Predicted secondary structure of two sgRNAs targeting mouse Btk. sgRNAs with low binding energy between the targeting sequence and the scaffold RNA do not disturb the secondary structure of the scaffold RNA (left), in contrast to sgRNAs with high binding energies between the targeting sequence and the scaffold RNA (right).

(D) Scheme of the CrispRGold sgRNA design tool.

(*E*) Scheme of the R26-LSL-Cas9iGFP allele before (top) and after crossing to Cre-deleter mice (bottom). Cre-mediated removal of the LSL cassette leads to ubiquitous expression of Cas9 and GFP in the entire mouse.

(*F*) FACS analysis of GFP expression in hematopoietic stem cells (HSC), defined as the Lin<sup>-</sup>,CD45<sup>+</sup>,Sca1<sup>+</sup>cKit<sup>+</sup> LSK population. LSK cells from wildtype C57BL/6 mice (red) and R26-Cas9iGFP/+ mice (blue) are shown on the same FACS plot (right).

(*G*) FACS analysis of GFP expression in B cells (CD19<sup>+</sup>), monocytes/dendritic cells (CD11b<sup>+</sup>), granulocytes (Gr-1<sup>+</sup>) and T cells (CD3<sup>+</sup>). Cells were isolated from the BM (top) or spleen (bottom). Cells from wildtype C57BL/6 mice (red) and R26-Cas9iGFP/+ mice (blue) are shown in the same FACS plots.

(*H*) Percentages of the indicated cell populations within all cells in the BM (left) or the spleen (right) of wildtype animals (red) or R26-Cas9iGFP/+ animals (blue). Each dot represents one mouse.

## Fig. 2. Efficient CRISPR-mediated gene inactivation in primary immune cells.

(*A*) Scheme of CRISPR-mediated gene knockout in Cas9-expressing B cells. Splenic B cells isolated from wildtype and Cas9 transgenic mice were mixed at a 1:4 ratio. These cells were activated with anti-CD40 and IL-4 for 2 days and transduced with retroviral particles expressing specific sgRNAs and a Puromycin resistance gene. One day after transduction, the cells were transferred into 40LB feeder cell plates and cultured in the presence of IL-21 and Puromycin.

(*B*) FACS plots of CD19<sup>+</sup>BFP<sup>+</sup> cells targeted with sgRNAs against B220 (left), CD24 (middle) and CD44 (right) 4 days after transduction. Knockout efficiencies were calculated based on the GFP<sup>+</sup> cells that lost the surface marker (sub-gates of the GFP<sup>+</sup> population).

(*C*) Summary of the knockout efficiencies for the 12 indicated surface marker genes with 2 sgRNAs/gene. Knockout efficiencies were determined by FACS 4 days after transduction of Cas9-expressing primary B cells with the respective sgRNAs.

(*D*) FACS plots showing the CD44 knockout efficiency in Cas9-expressing primary T cells. Mixed Cas9 (GFP<sup>+</sup>) and wildtype (GFP<sup>-</sup>) splenic T cells were activated and transduced with retroviral particles expressing sgCD44-1 and sgCD44-2. The gates indicate the knockout efficiency in GFP<sup>+</sup> T cells.

(*E*) FACS plots showing the knockout efficiency of CD64 and CD14 in nondividing primary BMDMs. BMDMs from C57BL/6 (top) and R26-Cas9iGFP/+ (bottom) mice were transduced with lentiviruses expressing sgCD64-2 and sgCD14-2.

(*F*) FACS plots of knockout efficiencies in primary B cells from R26-Cas9iGFP/+ (upper panel) and R26-Cas9p2aGFP/+ (lower panel) mice, using the same setup as in (a) and sgRNAs against B220, CD24 and CD44.

(*G*) Summary of knockout efficiencies using primary B cells from R26-Cas9iGFP/+ (blue) and R26-Cas9p2aGFP/+ (red) mice.

## Fig. 3. Identification of novel genes involved in B cell activation and differentiation using robust CRISPR-mediated screening.

(*A*) Scheme of CRISPR-mediated screening in primary B cells. GFP<sup>+</sup> splenic B cells from R26-Cas9iGFP/+ mice were mixed with GFP<sup>-</sup> wildtype B cells and treated as indicated in the 96-well system. The cells were analysed by flow cytometry at the indicated time points.

(*B*) FACS analysis of mixed  $GFP^+/GFP^-B$  cells 4 or 6 days after transduction with empty vector (control), sgPrdm1-1, sgXbp1-1, sgPou2af1-1, sgIrf4-1 or sgMyc-1. The gates show the percentage of  $GFP^+$  and  $GFP^-B$  cells after gating on  $BFP^+$  cells.

(*C*) Heatmap of the percentages of GFP<sup>+</sup> B cells at day 4, 6 or 8 normalized to day 2. One control sgRNA and three different sgRNAs per TF were used to knockout *Prdm1*, *Xbp1*, *Pou2af1*, *Irf4* or *Myc*.

(*D*) FACS analysis of plasma cell differentiation, indicated by the expression of CD138, in the same experiment as in (b). Gates are set on CD138<sup>+</sup> plasmablasts in GFP<sup>-</sup> (top) and GFP<sup>+</sup> (bottom) cells.

(*E*) Heatmap of the percentages of CD138<sup>+</sup> plasmablasts within the GFP+ and GFP- at the indicated time points.

(*F*) Fold changes of the percentages of  $GFP^+$  B cells at day 8 after transduction with sgRNAs against the candidate genes in two independent experiments. The repressor genes are indicated in green and the enhancer genes are shown in magenta. Fold changes were normalized to the initial percentage of  $GFP^+$  cells.

(*G*) Fold changes of the percentages of CD138<sup>+</sup>GFP<sup>+</sup> plasmablasts at day 8 in the same experiment as in (F). The repressor genes are indicated in blue and the enhancer genes are shown in red.

(*H*) Validation of the newly identified genes important in plasma cell differentiation. Histogram showing the percentages of GFP<sup>+</sup>CD138<sup>+</sup> and GFP<sup>-</sup>CD138<sup>+</sup> plasmablasts 6 days after transduction with sgRNAs targeting *Rosa26*, *Prdm1*, *Ost4*, *Arf4*, *Creld2*, *Zfp36* and *Edem1*.







### **Supporting information**

Chu et al.

#### Antibodies

The following anti-mouse antibodies were purchased from Biolegend: anti-B220 (RA3-6B2), anti-CD19 (6D5), anti-CD138 (281-2), anti-CD22 (OX-97), anti-CD24 (M1/69), anti-CD44 (IM7), anti-CD80 (16-10A1), anti-CD81 (Eat-2), anti-CD83 (Michel-19), anti-CD86 (GL-1), anti-CD45 (30-F11), anti-CD11b (M1/70), anti-Gr-1 (RB6-8C5), anti-CD3 $\epsilon$  (145-2C11), anti-Ter-119 (Ter-119), anti-Ckit (2B8), anti-Sca1 (D7), and anti-CD64 (clone X54-5/7.1). Anti-CD95 (15A7), anti-CD14 (clone Sa2-8) and anti-Tnfrsf13b (ebio8F10-3) antibodies were ordered from eBioscience. Anti-Cas9 (7A9-3A3) antibody was purchased from Cell Signaling Technology, anti-Flag (M2) antibody and anti- $\beta$ -actin from Sigma. HRP-conjugated secondary anti-mouse Ig antibody was ordered from Rockland Immunochemicals.



**Fig. S1.** Characterization of heterozygous R26-Cas9iGFP transgenic mice. (*A*) FACS plots of T cells (CD3<sup>+</sup>), myeloid cells (CD11b<sup>+</sup>) and B cells (CD19<sup>+</sup>) isolated from the mesenteric lymph nodes (mLN) of Cas9 transgenic or wildtype animals. (*B*) FACS plots of B cells (CD19<sup>+</sup>) isolated from the peritoneal cavity and spleen of Cas9 transgenic or wildtype animals. Gates were set on B1 and B2 cells in the spleen and B2, B1a and B1b in the peritoneal cavity based on CD19, B220, CD43, CD5 and IgM markers as indicated. (*C*) FACS plots of B220<sup>+</sup>CD19<sup>+</sup>Fas<sup>hi</sup>CD38<sup>low</sup> germinal center (GC) B cells isolated from the spleen and mLN. Data are representative for four independent experiments. Cells from wildtype C57BL/6 animals (red) and R26-Cas9iGFP/+ animals (blue) are shown in the same FACS plots.



**Fig. S2.** Comparison of GFP expression levels in hematopoietic cells of R26-Cas9iGFP and R26-Cas9p2aGFP mice. (*A*) FACS plots of hematopoietic stem cells (defined as Lin<sup>-</sup>Sca1<sup>+</sup>c-Kit<sup>+</sup> or LSK), isolated from the bone marrow (BM) of R26-Cas9iGFP/+ (blue), R26-Cas9p2aGFP/+ (red) and C57BL/6 (black) mice. (*B*) FACS plots of B cells (CD19<sup>+</sup>), T cells (CD3<sup>+</sup>), granulocytes (Gr1<sup>+</sup>CD11b<sup>+</sup>) and monocytes/dendritic cells (CD11b<sup>+</sup>Gr1<sup>-</sup>) isolated from the BM and spleen of R26-Cas9iGFP (blue), R26-Cas9p2aGFP (red) and wild-type (black) mice. The data are representative for 2 independent experiments.



Fig. S3. Cas9 levels and localization in primary cells of R26-Cas9iGFP/+ mice and Cas9-dosedependent knockout efficiency. (A) Western blots of Cas9 proteins (detected by anti-Cas9 and anti-Flag antibodies) in anti-CD40/IL-4-activated B cells and BM-derived macrophages (BMDMs) isolated from R26-Cas9iGFP/+ and wildtype animals. Beta-Actin was used as a loading control. (B) Immunostaining of Cas9 in BMDMs derived from R26-Cas9iGFP/+ and wildtype animals. Cas9 was stained using anti-Flag (M2) antibodies (red). DAPI was used for counterstaining of nuclei (green). (C) Knockout efficiencies in B cells from R26-Cas9iGFP/+ (top) and R26-Cas9p2aGFP/+ (bottom) mice that were activated and transduced with retroviruses expressing sgCD38-1, sgCD38-2, sgFas-1 and sgFas-2 side by side. The gates in the FACS plots show the percentage of cells that lost the surface marker 4 days after transduction. (D) Western blots of Cas9 proteins in anti-CD40/IL-4-activated B cells isolated R26-Cas9iGFP/+, R26-Cas9p2aGFP/+ from and wildtype animals. The data are representative for 2 independent experiments.



**Fig. S4.** Knockout efficiencies in Cas9-expressing B cells. (*A*) Selection of sgRNAs designed with CrispRGold targeting 12 B cell surface makers. Targeting sequences are shown 5'-3'. (*B*) FACS plots showing the knockout efficiencies in Cas9-expressing B cells 4 days after transduction with sgRNAs against the indicated surface markers. Note that some surface markers are only induced upon activation (CD83 and Tnfrsf13b). The sgRNAs used are indicated in the title of each FACS plot. The data are representative for 3 independent experiments.



**Fig. S5.** CRISPR-mediated gene knockout in primary T cells and BMDMs. (*A*) Cas9expressing T cells were mixed with wild-type T cells. The cells were activated with anti-CD3 and anti-CD28 antibodies for 2 days. The activated T cells were transduced with retroviral particles expressing sgRNAs. The transduced T cells were selected with Puromycin for 4 days before assessment of knockout efficiencies. (*B*) Scheme of CRISPR-mediated gene knockout in primary BMDMs. BM cells from R26-Cas9iGFP/+ mice and wild-type mice were differentiated for 2 weeks in the presence of M-CSF. BMDMs were transduced with lentiviruses expressing sgRNAs targeting mouse CD64 and CD14 genes. The transduced cells were selected with Puromycin for 4 days and the gene knockout efficiency was analysed by flow cytometry. (*C*) FACS plots of the transduced BMDMs derived from control animals (top) and R26-Cas9iGFP/+ animals (bottom) stained with antibodies against the targeted surface markers. The summary of the experiments is shown on the right side. The data are representative for 3 individual Cas9 transgenic mice (each representing a data point in the summary).



**Fig. S6.** Robust detection of transcription factors important for B cell activation and differentiation. (*A*) Scheme of B cell activation and differentiation. Splenic B cells were activated with anti-CD40+IL4 and subsequently differentiated into plasmablasts using IL21 and CD40-ligand expressing feeder cells. (*B*) sgRNAs used to target *Prdm1*, *Xbp1*, *Pou2af1*, *Irf4* and *Myc*. (*C*) Analysis of the percentage of GFP<sup>+</sup> Cas9-expressing B cells and GFP<sup>-</sup> wild-type B cells. These cells were mixed, activated and transduced (individually) with 3 sgRNAs per gene targeting *Prdm1*, *Xbp1*, *Pou2af1*, *Irf4* and *Myc*. The FACS plots show the GFP<sup>+</sup> Cas9-expressing B cells 4 and 8 days after transduction. (*D*) Analysis of the induction of CD138 in GFP<sup>+</sup> Cas9-expressing and GFP<sup>-</sup> wildtype B cells 6 days after transduction. The data are representative for 3 independent experiments.



**Fig. S7.** Gene set used for the small-scale screen. Total RNA was isolated from follicular B, GC and plasma cells that were isolated from the spleen and BM of immunized animals. Microarrays were performed and data were normalized prior to analysis. The heat-map shows the expression levels of the selected genes with differential expression in the plasma cell populations.

Α



**Fig. S8.** Small scale-CRISPR-mediated screening to detect novel genes important for B cell activation and plasma cell differentiation. (*A*) Scheme of the high-throughput 96 well-based screening approach. One sgRNA was designed for each gene of interest. Each sgRNA was separately cloned into a retroviral vector. The cloning, retrovirus production and transduction of B cells was performed in 96 well plates. The transduced cells were selected with Puromycin. Survival/proliferation was analyzed based on monitoring the percentage of GFP<sup>+</sup> Cas9-expressing B cells at different time points. The fold changes were calculated by comparing to the percentage of GFP<sup>+</sup> cells at day 2, allowing for the identification of repressor or enhancer genes. The plasmablast assay was based on the normalization of GFP<sup>+</sup>-derived CD138<sup>+</sup> plasmablasts. (*B*) Diagram showing the screen set and identification of important genes for survival and proliferation (middle, top) and plasma cell differentiation (middle, bottom). 22 genes were crucial for activated B cell survival and proliferation, and 8 genes were essential for plasma cell generation.10 genes were essential for both survival/proliferation and differentiation.



**Fig. S9.** Validation of novel candidates important for B cell survival/proliferation and plasma cell differentiation. (*A*) FACS data showing the percentage of GFP<sup>+</sup> Cas9 expressing and GFP<sup>-</sup> wildtype B cells 2, 4, 6 and 8 days after transduction with the indicated sgRNAs. As before, Cas9-expressing B cells were mixed with wild-type B cells, activated and transduced with retroviral particles expressing sgRNA targeting *Rosa26 (as control), Myc, Egr2, Dis3, Ost4, Preb,* and *Pomp. (B)* Analysis of the induction of CD138 in GFP<sup>+</sup> Cas9-expressing and GFP<sup>-</sup> wildtype B cells. Cas9-expressing B cells were mixed with wild-type B cells, activated and transduced and transduced with retroviral particles expressing sgRNAs targeting *Rosa26 (as control), Myc, Egr2, Dis3, Ost4, Preb,* and *Pomp. (B)* Analysis of the induction of CD138 in GFP<sup>+</sup> Cas9-expressing and GFP<sup>-</sup> wildtype B cells. Cas9-expressing B cells were mixed with wild-type B cells, activated and transduced and transduced with retroviral particles expressing sgRNAs targeting *Rosa26 (as control), Prdm1, Ost4, Arf4, Creld2, Zfp36,* and *Edem1.*