


Study of the carriage and characterization of a common bacterium, *Staphylococcus aureus*, on the campus of Concordia University St. Paul

Hannah Pedersen, Dr. Mandy Brosnahan, and Dr. Taylor Mach

BACKGROUND: *Staphylococcus aureus* is a dual role bacterium and is able to live commensally in some patients, but can cause disease in others. Methicillin resistant *S. aureus* (MRSA) causes life threatening disease in patients and can be very hard to treat since it is resistant to many antibiotics. In order for *S. aureus* to wreak havoc on the body, it must be able to have specific genes expressed to secrete toxins. These toxins are what causes the patients to get a wide variety of symptoms like boils, scalded skin, toxic shock syndrome, pneumonia or sepsis. This study sets out to find the amount of people who carry *S. aureus* in their nasal cavity and if they do, is is MRSA and does it have the specific genes expressed to secrete toxins. As of 02/26/20, there have been 1065 swabs collected from the community of Concordia St. Paul's campus and 880 of those swabs have been processed through culture tests. Out of those 880 swabs, 234 have been positive for *S. aureus* and this is a 26.6% carriage rate.

- TEST SIGNIFICANCES
- Mannitol Salt Agar (MSA) plate contains a high salt concentration and mimics what it would be like for a bacterium to grow on human skin. *S. aureus* can grow on this plate because it can thrive under high salt conditions. When it grows, it causes the fermentation of mannitol, which gives it the yellow color change on the plate (2).
 - Deoxyribonuclease (DNase) plate looks for the presence of DNase enzymes. The plate's medium contains DNA in it that can be broken down by a *S. aureus* enzyme called deoxyribonuclease, which you can see by the halo color change. The deoxyribonuclease catalyzes the hydrolytic cleavage of phosphodiester linkages in DNA backbone. This gives the bacteria nutrients and helps to support its growth (3).
 - Coagulase test is done to see if the bacterium has the coagulase enzyme, which will cause aggregation of red particles on the blue background. The coagulase enzyme converts fibrinogen to fibrin which gives *S.aureus* the ability to clot blood/plasma.Fibrin also can coat the bacteria and prevent it from getting phagocytosed (4).
 - Columbia Nalidixic Acid Blood agar (CNA) plate contains two antibiotics that doesn't allow gram negative organisms to grow on it. It also contains blood which allows for beta hemolysis to be seen. In order for the yellow glow to occur, the bacterium must possess hemolysins that cause blood cell lysis. *S.aureus* use their hemolysins to lyse RBC, which causes hemostatic disturbances (1).
 - Catalase test is looking for the presence of a catalase enzyme, and when it comes into contact with hydrogen peroxide, will cause O2 gas to be released (bubbles). *S.aureus* uses catalase enzyme to protect itself from oxidative damage (5).
 - Gram stain is done in order to look at the bacterium under the microscope to make sure that it looks like a *S.aureus* bacterium. Typical *S.aureus* bacterium are Gram Positive (Purple color) and in cocci shape and grape like orientations.



Methods of Testing



RESULTS

Isolate Tests							
Isolate #'s	MSA	DNASE	Coagulase	CNA	Catalase	Gram Stain	Final Results
D2D0561	Negative	n/a	n/a	n/a	n/a	n/a	Negative
D2D0563	Positive	Negative	Positive	Positive	Positive	Gram +, purple cocci/cluster	Negative
D2D0296	Positive	Positive	Positive	Positive	Positive	Gram +, purple cocci/cluster	Positive
D2D0480	Negative	n/a	n/a	n/a	n/a	n/a	Negative
D2D0550	Positive	Positive	Positive	Positive	Positive	Gram +, purple cocci/cluster	Positive
D2D0467	Negative	n/a	n/a	n/a	n/a	n/a	Negative
D2D0482	Positive	Negative	Positive	Positive	Positive	Gram +, purple cocci/cluster	Negative
D2D0554	Negative	n/a	n/a	n/a	n/a	n/a	Negative
D2D0439	Positive	Negative	Positive	Positive	Positive	Gram +, purple cocci/cluster	Negative
D2D0567	Positive	Negative	Positive	Positive	Positive	Gram +, purple cocci/cluster	Negative
D2D0569	Negative	n/a	n/a	n/a	n/a	n/a	Negative

There were eleven isolates that were tested in this study and out of the eleven, there were two that tested positive for being *S. Aureus*. This comes to a 18% carriage rate for these eleven isolates. For the study as the whole, as of 02/26/20 the carriage rate for positive isolates is 26.6%.



References:

1. A.Aksimentiev, K. Schulten, Bacterial toxin alpha-hemolysin. *Theoretical & Computational Biophysics Group*. **88**.3745-3761, (2005)
2. S. Aryal, Mannitol salt agar for the isolation for *staphylococcus aureus*. *Microbio. Info.com*. **1** (2019) [MicroBiology Information.com]
3. A.Tankeshwar, Deoxyribonuclease test: principle, procedures and results. *Microbe online*. **1** (2014) [Microbe Online].
4. M.McAdow, D.M. Missiakas, O.Schneewind, *Staphylococcus aureus* secretes coagulase & von willebrand factor binding protein to modify the coagulation cascade and establish host infections. *J. Innate Immun.* **4**(2), 141–148 (2012).
5. K.Reiner, Catalase test protocol. *American Society for Microbio*. **1** (2010) [American Society for Microbiology].

Acknowledgements:

Special thanks to Dr. Patrick Schlievert (University of Iowa) for helpful conversations. This research was partially funded by eight CSP Faculty Development Grants. This work has IRB approval from CSP (studies 2016_42 & 2018_37). Special thanks to Dr.Brosnahan and Dr.Mach