The Coronary Intravascular Lithotripsy System

Maria Natalia Tovar Forero and Joost Daemen

Thoraxcenter, Erasmus University Medical Center, Rotterdam, the Netherlands

Abstract

Calcified lesions often mean percutaneous intervention results are suboptimal and increase the risk of procedural complications and future adverse events. Available plaque-modifying devices rely on tissue compression or debulking, with the intention of fracturing calcium and facilitating optimal stent deployment. In contrast, coronary intravascular lithotripsy delivers unfocused, circumferential, pulsatile mechanical energy to safely disrupt the calcium within the target lesion. The present review summarises the evidence available so far on this therapy and includes a practical description of the components and function of the Shockwave Intravascular Lithotripsy System (Shockwave Medical).

Keywords

Shockwave Intravascular Lithotripsy System, coronary calcification, coronary plaque modification, coronary lithotripsy.

Disclosure: JD has received institutional research support from Pie Medical, Acist Medical, PulseCath, Medtronic, Boston Scientific and Abbott Vascular, and speaker and consultancy fees from PulseCath, Medtronic, ReCor Medical and Acist Medical. MNTF has no conflict of interest to declare. Acknowledgements: The authors thank Jason Hokama and Suzanne Wallace from Shockwave Medical for their critical review of the manuscript and supporting material.

Received: 1 July 2019 Accepted: 8 August 2019 Citation: Interventional Cardiology Review 2019;14(3):174–81. DOI: https://doi.org/10.15420/icr.2019.18.R1 Correspondence: Joost Daemen, Department of Cardiology, Room Rg-628, Erasmus University Medical Center, PO Box 2040, 3000 CA Rotterdam, the Netherlands. E: j.daemen@erasmusmc.nl

Open Access: This work is open access under the CC-BY-NC 4.0 License which allows users to copy, redistribute and make derivative works for non-commercial purposes, provided the original work is cited correctly.

The amount of coronary artery calcification increases with age and the presence of cardiovascular risk factors and comorbidities.^{1,2} Up to 20% of percutaneous coronary intervention (PCI) procedures are challenged by severe calcifications, and coronary calcifications have been shown to be an independent predictor of PCI failure and future adverse cardiac events.^{3,4} Lesion calcification increases procedural complexity and time. More specifically, calcium localisation (superficial or deep), distribution (focal, circumferential and longitudinal extension) and thickness influence procedural success, stent delivery and deployment.⁵

Several techniques to treat calcified lesions in native coronary arteries are available, including high-pressure and super-high pressure noncompliant balloons, cutting/scoring balloons, atherectomy devices, both rotational and orbital, and excimer lasers.⁶ These devices rely on tissue compression and or tissue debulking, and have higher rates of procedural complications, such as dissections, perforations and distal embolisation. Moreover, their success rate is reduced when deep, thick or eccentric calcifications are present, and the induced tissue injury might accelerate uncontrolled neointimal growth and restenosis.^{6,7} So far, neither specialty balloons nor atherectomy devices have been proved to be superior to high-pressure non-compliant balloons in improving clinical outcomes.^{8,9}

Recently, an alternative way to disrupt calcium has been developed that is based on the lithotripsy concept used to treat kidney and ureteral stones. The Intravascular Lithotripsy (IVL) System (Shockwave Medical) transforms electrical energy into mechanical energy during low-pressure balloon inflation.¹⁰ The technology does

not rely on direct vascular tissue injury for plaque modification but on sonic waves, which travel from the balloon-based catheter to the surrounding tissue with the intention of safely and selectively breaking both superficial and deep calcium deposits with minimal soft tissue impairment, while improving vessel compliance. In contrast to debulking techniques, the calcium fragments resulting from the IVL therapy remain in situ, reducing the likelihood of distal embolisation. The safety and efficacy of the IVL with minimal vessel injury was demonstrated in prospective single-arm studies, where moderate to severely calcified lesions in peripheral artery disease were tackled; these were the Shockwave Lithoplasty Disrupt Trial for PAD 1 and 2 (Disrupt PAD 1 and 2) and Safety and Feasibility of the Shockwave Lithoplasty® System for the Treatment of Peripheral Vascular Stenosis (BTK) trials.¹¹⁻¹⁴

The Coronary Intravascular Lithotripsy System

The Shockwave Coronary Rx Lithoplasty® Study (Disrupt CAD I), a multicentre, prospective, single-arm study conducted in seven centres in Europe and Australia, was the first to assess the safety and efficacy of the Shockwave Coronary IVL in 60 patients with severe (100%) calcified lesions in native coronary arteries before drug-eluting stent implantation. The balloon was successfully delivered in 59 patients and stent implantation was successful in all cases, with no major procedural complications, such as slow flow, no reflow, distal embolisation or perforation reported. Device delivery was facilitated by pre-dilatation with a small balloon in 37% of patients. Clinical success was 95%, defined as residual diameter stenosis <50% without in-hospital major cardiac adverse events (MACE; a composite of death, MI and target-

vessel revascularisation) at 30 days; MACE at 6 months increased to 8.5% with three non-Q wave MI events within the first 30 days and two cardiac deaths. $^{\rm 15}$

An optical coherence tomography sub-study in 31 patients in the Disrupt CAD I trial confirmed the effects of IVL on the vessel wall.¹⁶ Calcium fractures were evident in 42.9% of the lesions and multiple circumferential cracks in the same cross-sectional area were witnessed in more than 25% of the cases, with a higher incidence of fractures in heavier calcified plaques (highest tertile versus lowest tertile; p=0.009). The presence of calcium ruptures allowed to increase acute lumen gain (mean acute area gain = 2.1 mm^2) independently of the degree of calcification, enabling successful stent implantation with uniform expansion. Coronary dissections (type B or greater) occurred in four cases during angioplasty and were successfully treated with stent implantation; no other procedural complications were reported. This data provided the first results encouraging the use of the Shockwave Coronary IVL for the treatment of calcified lesions in the coronary vasculature. Following these findings, the device received European CE mark approval in May 2017 for commercial use.

The findings of Disrupt CAD I were further confirmed in a small (n=26) real-world study including patients with both stable and unstable disease, either as an upfront calcium modification technique or as a bailout after suboptimal results with standard balloon dilatation. Angiographic success occurred in all cases (<20% residual stenosis) without procedural complications.¹⁷

At present, the results of the Disrupt CAD II (post-market prospective, multicentre, single-arm study), which enrolled 120 patients across Europe, are pending. More recently, the Disrupt CAD III study was announced. The study will aim to enrol 392 patients in 50 centres across the US and Europe with the intention to obtain Food and Drug Administration (FDA) approval.

The available studies and case reports on the use of coronary IVL are summarised in *Table 1*.

The System

The Coronary IVL System consists of a portable, rechargeable generator, a connector cable with a push button to allow manually controlled delivery of electric pulses, and a 6 Fr compatible, rapid-exchange, semi-compliant balloon catheter to be used following standard angioplasty practice over a 0.014" guidewire (*Figure 1*).

The semi-compliant balloon integrates two radiopaque lithotripsy emitters 6 mm apart and two conventional markers at the proximal and distal edges of the balloon. These emitters receive electrical pulses from the generator vaporising the fluid (a standard mixture of 50% NaCl 0.9% and 50% radiopaque contrast) within the balloon and creating a rapidly expanding and collapsing bubble. This bubble can transmit unfocused circumferential pulsatile mechanical energy into the vessel wall, in the form of sonic pressure waves equivalent to approximately 50 atmospheres (atm). The balloons are available in diameters ranging from 2.5 mm to 4.0 mm with a standard length of 12 mm (*Table 2*); their crossing profiles range from 0.043" to 0.046" (*Figure 2*).

The IVL therapy consists of a series of 10 pulses (1 cycle) or 10 seconds (1 pulse per second). The number of therapies needed per

lesion will depend on lesion resistance; however, the maximum number of pulses to be delivered by each individual catheter is limited to 80 pulses (eight cycles).

Intravascular Lithotripsy Procedure

The IVL procedure does not require high-level additional training for interventional cardiologists. The shockwave balloon must be sized accordingly to the reference vessel diameter (ratio 1:1), placed in the target calcified lesion and inflated up to 4 atm to ensure apposition to the vessel wall; the lithotripsy emitters are then activated to deliver the acoustic pulses by pushing the button on the connector cable. Once a cycle of 10 pulses has been delivered, the balloon can be inflated up to 6 atm (nominal pressure) to increase balloon compliance and to assess symmetrical expansion, confirming calcium modification. Next, the balloon is deflated carefully to allow small air bubbles to escape. The previous steps must be repeated for each intended IVL cycle and at least two IVL cycles are recommended to treat the target area. For the treatment of lesions longer than 12 mm, the catheter needs to be repositioned and overlapping treatment areas might occur. See *Supplementary Videos 1–4*.

Due to the slightly higher profile of the shockwave catheter, predilatation with standard balloons might be necessary in some cases to facilitate deliverability and positioning, especially when lumen reduction is severe. Notwithstanding this, the balloon allows the use of guide-catheter extenders and buddy-wire support.¹⁸ Furthermore, although the system is labelled 6 Fr compatible, it could be used with a 5 Fr guiding catheter where the radial artery is small.¹⁹ The use of dilatation with non-compliant balloons after IVL, although not mandatory, could be considered to expand the lumen further. Moreover, aggressive plaque modification devices such as cutting/ scoring balloons or atherectomy could be used as adjuvant therapy in challenging lesions to improve results.²⁰

Potential Uses

At present, the instructions for use of the Shockwave Coronary IVL System restrict its use to lesion preparation in native coronary arteries (*Figure 3*). Given it is assumed to be safer than previous approaches, the number of case reports and small case series reporting on its use in more challenging scenarios is increasing.

Acute Coronary Syndromes

Calcified lesions in culprit vessels are common in patients presenting with acute coronary syndromes (moderate calcification is found 26.1% of these patients and severe calcification in 5.9%), and their presence is a strong predictor of definite stent thrombosis (HR 1.62; 95% CI [1.14–2.30]; p=0.007) and target lesion revascularisation (HR 1.44; 95% CI [1.17–1.78]; p<0.001).¹ The Disrupt CAD study included only patients with stable and unstable angina. Although there is not enough evidence to support the use of the IVL during primary PCI, early experience has shown favourable results.^{17,21}

Unprotected Left Main Calcified Stenosis

PCI has become an option for the treatment of left main (LM) disease, with a class IA recommendation for patients with a SYNTAX score \leq 22 and class IIA for a SYNTAX score 23–32.²² Calcification increases procedural complexity and therefore the risk of complications. Although atherectomy has previously been proposed as a valid and feasible option, the high-risk profile of those patients means safer approaches would be welcome.^{23,24}

Reference	Country	Study	Study Design	Sample Size	Sex	Age (Years)	NSTEMI/STEMI	Lesion type	Imaging	Follow-up (Months)	Outcome
Ali et al. 2017 ¹⁶	International (five countries)	Disrupt CAD OCT substudy)	Prospective, multicentre, single arm	31	80% men	71 ± 10	%0	<i>De novo</i> lesions with severe calcification in 87.1%	OCT (100%)	NA	Acute area gain 2.08 \pm 1.65 mm ²
2017 ¹⁵	International (five countries)	Disrupt CAD	Prospective, multicentre, single arm	09	80% men	72 (66-79)	%0	De novo lesions with severe calcification in 100%	Ч	0	Acute gain: mean 1.7 mm (range 1.3–2.1 mm) Angiographic complications: 0% Clinical success (residual stenosis <50% without in-hospital MACE: 95% MACE 30 days: 5% MACE 6 months: 8.3%
De Silva et al. 2017⁴⁰	Ж	Case report	NA	~	Man	69	%0	<i>De novo</i> lesion Severe calcification	OCT	NA	Angiographic result and OCT result: good No complications
Watking et al. 2018 ³⁸		Case report	AA	←	Man	67	0%	Stent underexpansion Severe calcification	OCT	AN	Angiographic and OCT result: good No complications
Ali et al. 2018%	NS	Case report	ИА	~	Man	73	%0	Stent underexpansion Severe calcification	OCT	АЛ	Angiographic result and OCT result: good; acute gain 5.38 mm ² No complications
Tovar Forero et al. 2018 ²²	The Netherlands	Case report	AN	-	Man	74	%0	stent underexpansion Severe calcification	OCT	۲Z	Angiographic and OCT result: good; acute gain: 4.26 mm ^a No complications
Morabito et al. 2018³	Italy	Case report	ИА	~	Woman	77	%0	Stent underexpansion Severe calcification	SUVI	AN	Angiographic and IVUS result: good; acute gain: 3.83 mm ² No complications
Chen et al. 2018 ⁴⁵	Germany	Case report	ИА	~	Man	61	100%	ISR Severe calcification between stent layers	OCT	AA	Angiographic and OCT result: good No complications
Rodriguez Costoya et al. 2018 ⁵⁰	Spain	Case series	NA	m	66.6% men	63, 73 and 81	%0	De novo lesion Severe calcification	OCT	NA	Angiographic and OCT result: good No complications
Venuti et al. 201	944 Italy	Case report	NA	~	Men	67	%0	Undilatatble lesion Severe calcification	SUVI	NA	Angiographic and IVUS result: good result No complications
Soriano et al. 2019⁴8	Italy	Case report	AA	-	Man	47	0%	De novo lesion Severe calcification	OCT	NA	Complication: balloon rupture and vessel dissection
Tassone et al. 2019 ³⁷	Italy	Case report	NA	-	Man	60	%0	Stent underexpansion Severe calcification	SUVI	NA	Angiographic and IVUS result: good Acute area gain 4.19 mm ² No complications
ISR = in-stent reste	nosis; NUS = intra-vasc	ular ultrasound; IV	/L = intra-vascular lithot	tripsy; MACE= I	major adverse cal	diac events; NA = r	not applicable; OCT = op	otical coherence tomography.			(Continued)

Coronary Calcified Lesions

Reference	Country	Study	Study Design	Sample Size	Sex	Age (Years)	NSTEMI/STEMI	Lesion Type	Imaging	Follow-up (Months)	Outcome
López-Lluva et al. 201947	Spain	Case report	NA	-	Woman	73	0%	<i>De novo</i> lesion Severe calcification	None	NA	Complication: Balloon rupture and vessel dissection
Salazar et al. 2019 ³	Spain	Case report	NA		Man	71	100%	Stent underexpansion Severe calcification	OCT	NA	Angiographic and OCT result: good Acute gain: 6.27 mm ² No complications
Legutko et al. 201951	Poland	Case report	NA	-	Woman	79	100%	<i>De novo</i> Severe calcification.	IVUS/OCT	NA	Angiographic and OCT result: good No complications
Vainer et al. 2019 ⁴⁶	The Netherlands	Case report	NA	_	Woman	70	0%	<i>De novo</i> lesion Severe calcification	OCT	NA	Angiographic and OCT result: good No complications
Sgueglia et al. 2019 ¹⁹	Italy	Case report	NA	_	No info	74	100%	<i>De novo</i> lesion Severe calcification	SUM	NA	Angiorgaphic and IVUS result: good No complications
Cicovic et al. 2019⁴¹	New Zeland	Case report	NA		Woman	73	100%	<i>De novo</i> lesion Severe calcification	None	NA	Cardiac pacing secondary to IVL
Warisawa et al. 2019 ¹⁸	UK	Case report	NA	-	Man	74	0%	<i>De novo</i> lesion Severe calcification.	None	NA	Angiographic result: good No complications
Azzalini et al. 2019 ²⁶	Italy	Case series	NA	2	Man	62 and 71	0%	CTO Severe calcification	None	NA	Angiographic result: good No complications.
Yeoh et al. 2019 ²⁷	UK	Case report	NA	-	Woman	81	0%	CTO Severe calcification	SUM	NA	Angiographic and IVUS result: good No complications
Urbano Carillo et al. 2019 ³⁴	Spain	Case report	NA	_	Male	71	100%	Stent underexpansion Severe calcification	OCT	NA	Angiographic and OCT result: good No complications
Wilson et al. 2019 ⁴²	UK	NA	Retrospective, single centre	54	No capture: 66.7%	No capture: 77.2 ± 6.7	51.8%	No information	None	NA	IVL provoked ventricular capture in 77.8% of the cases
					Capture: 81%	Capture: 71.6 ± 9.3					Heart rate <65 bpm increases risk of ventricular pacing (OR 16.3; (95 CI% [2.4–110.8], p=0.004)
Alfonso et al. 2019 ³⁵	Spain	Case report	NA	_	Man	63	100%	Stent underexpansion Severe calcification	OCT	NA	Angiographic and OCT result: good No complications
Wong et al. 2019 ²¹	New Zealand	Case series	NA	ω	100%	61, 70 and 76	100%	<i>De novo</i> lesions (66.6%) Stent underexpansion (33.3%) Severe calcification	None	NA	Angiographic result: good No complications
Wong et al. 2019 ²⁵	New Zealand	Case series	NA	ω	66.6% men	60, 64 and 96	66.6%	<i>De novo</i> lesions in left main coronary artery Severe calcification	None	NA	Angiographic result: good No complications
Wong et al. 2019 ¹⁷ ISR = in-stent reste	New Zealand	NA nular ultrasound: 1	Prospective, single centre, observational	26 httinsv: MACE	69% men	72 ± 8	54%	De novo lesions	None	NA	Procedural and clinical success: 100% No complication
ISR = in-stent reste	nosis; IVUS = intra-vas	cular ultrasound; I	VL = intra-vascular lithc	tripsy; MACE=	major adverse ca	rdiac events; NA = r	not applicable; OCT = op	tical coherence tomography.			

INTERVENTIONAL CARDIOLOGY REVIEW

Coronary Intravascular Lithotripsy System

Table 1: Cont.

Table 2: Intravascular Lithotripsy Balloon Sizes and Compliance for Coronary Intervention

Pressure (atr	n)			Balloon size			
	2.5 × 12 mm	2.75 × 12 mm	3.0 × 12 mm	3.25 × 12 mm	3.5 × 12 mm	3.75 × 12 mm	4.0 × 12 mm
4 (IVL)	2.4	2.6	2.9	3.2	3.4	3.7	3.8
5	2.4	2.7	2.9	3.2	3.4	3.7	3.8
6 (NP)	2.5	2.7	2.9	3.2	3.4	3.7	3.9
7	2.5	2.7	3.0	3.3	3.5	3.8	3.9
8	2.5	2.7	3.0	3.3	3.5	3.8	4.0
9	2.5	2.8	3.1	3.3	3.5	3.9	4.1
10 (RBP)	2.5	2.8	3.1	3.4	3.6	3.9	4.1

atm = atmospheres; IVL = intravascular lithotripsy; NP = nominal pressure; RBP = rated burst pressure.

Figure 1: Shockwave Coronary Intravascular Lithoplasty System



Source: Adapted with permission from Shockwave Medical

Figure 2: Intravascular Lithoplasty Therapy



Source: Adapted with permission from Shockwave Medical.

The Coronary IVL System, with controlled pulses delivered under low pressure, might potentially improve plaque modification with a lower risk of vessel closure, perforation or embolisation.²⁵

Chronic Total Occlusions

Moderate to severe calcification is frequently found in chronic total occlusions, and debulking devices are usually avoided because the

procedure is difficult and has a high risk of complications; where standard balloons fail, the Coronary IVL System might be useful in facilitating lumen dilatation and communication with the subintimal space.^{26,27}

Stent Underexpansion due to Underlying Calcification

Although the technique has been developed to treat calcified lesions in native coronary arteries before stenting, patients with severe stent underexpansion because of heavy calcification are at a higher risk of stent failure and future adverse events. Until now, undilatable lesions in previously stented segments have been courageously approached with debulking devices such as cutting balloons and atherectomy, with unpredictable results and risks of procedural complications and stent damage.²⁸⁻³¹

Of note, the effectiveness of those techniques is limited by the presence of metallic struts, and deeper calcifications therefore remain unaffected. The circumferential sonic waves of the Coronary IVL System, conversely, have the advantage of extending beyond strut layers and fracture deeper calcium deposits (*Figure 4*). Several case reports have supported the use of the technology for optimising stent expansion without complications.³²⁻³⁹

Of note, the efficacy of the system in segments with multiple layers of stents has not been demonstrated and its impact on stent backbone/ polymer integrity and drug elution is still unknown. Nonetheless, at present, there are no alternative percutaneous options for patients left with underexpanded stents due to heavy calcification.

Effects on Cardiac Rhythm

Electric signals similar to pacing spikes on the electrocardiogram (ECG) tracing during pulse-delivery have been described.^{17,40-42} These so-called 'shocktopics' and asynchronous cardiac pacing have been reported in up to 77.8% of the cases, with a 16-fold increased risk in patients with a heart rate <65 BPM. Cardiac pacing has not been linked to any specific number of IVL cycles or coronary artery anatomy, although its frequency is higher when either the left anterior descending artery or the right coronary artery are treated.⁴²

The exact mechanism behind this phenomenon is still unclear. A potential explanation might be that the transformed mechanical energy reaches and couples with the cardiac conduction system producing ectopic atrial and or ventricular captures.⁴³ Although no relevant clinical events have been reported, it warrants paying special attention to ECG and aortic pressure waveforms changes during IVL administration; the resulting VOO pacing mode is theoretically

Figure 3: Intravascular Lithoplasty Therapy: Effect on Heavy Calcified Coronary Lesions



 $LA = 2.55 \text{ mm}^2$

 $LA = 6.15 \text{ mm}^2$

 $LA = 8.32 \text{ mm}^2$

Angiography left anterior descending coronary artery: a: before intravascular lithotripsy (IVL); b: after IVL; c: after stenting. Red arrow shows severe calcification. Optical coherence tomography: d: before IVL; e: after IVL; f: after stenting. Asterisks show calcium deposits in the vessel wall. Arrows show calcium fractures. LA = lumen area.

pro-arrhythmic (potential R on T phenomenon) and, until further data become available, pacemaker carriers should be assessed for inappropriate device sensing during the IVL cycles and to ensure correct pacing function post-procedure. Further investigation in the matter will be provided by the substudies in the Disrupt CAD III trial.

Challenging Lesions

Some lesions might respond better to IVL therapy than to other plaque modification approaches. Reports have pointed out the usefulness of the IVL therapy on creating calcium fractures as assessed with intravascular ultrasound, and achieving optimal stent expansion in undilatable lesions that have been resistant to specialty balloons and rotational atherectomy.^{44–46}

In contrast, some lesions might not be suitable for IVL treatment or may remain resistant after the application of all 80 pulses. Severe tortuosity or angulation, critical lumen reduction, plaque indentation into the lumen and a very low vessel expansion compliance (small vessels and multiple stent layers present), could impact balloon deliverability and positioning. Up to 46% of the lesions might also require dedicated lesion predilatation and/or post-dilatation with non-compliant balloons or could benefit from adjuvant lesion preparation with conventional devices such as specialty balloons or atherectomy to either facilitate balloon delivery or increase calcium compliance after lithotripsy therapy.^{17,45}

Figure 4: Intravascular Lithotripsy Therapy Effect on Stent Underexpansion due to Calcified Lesion



Fluoroscopy: A: before intravascular lithotripsy (IVL); B: after IVL. Red arrow shows 'dog bone' effect. Optical coherence tomography: C: before IVL d: after IVL. Asterisks show calcium deposits in the vessel wall. Arrows show calcium fractures. LA = lumen area.

While balloon rupture is uncommon, it can cause vessel complications. Case reports have described sudden IVL balloon burst during lithotripsy therapy with important vessel dissection; however, it is fair to highlight that either critical stenosis or severe vessel tortuosity were present, which suggests that, in some anatomies, the IVL system might not be suitable or should be used with caution.^{47,48}

Furthermore, vessels with a diameter >4 mm (maximum shockwave balloon size) or important plaque eccentricity preclude appropriate IVL balloon apposition to the vessel wall, and may reduce the efficacy of the therapy. More data are needed on the specific efficacy of IVL in concentric versus eccentric lesions. In the Disrupt CAD study, 22% of the patients had eccentric plaques; nevertheless, overall device success was 98%.¹⁵

Moreover, performing intracoronary imaging where important coronary calcification is suspected during the angiographic assessment could help to accurately assess calcium distribution, localisation and thickness. The use of intracoronary imaging before and after lithotripsy therapy could not only assist the appropriate selection of IVL balloon sizes but also potentially help to identify IVL responders and identify patients who might need adjuvant therapy from other plaque-modification devices .

Further clinical data are needed to assess the implications of patient, vessel and lesion characteristics on the effectiveness of this therapy,

- Genereux P, Redfors B, Witzenbichler B, et al. Two-year outcomes after percutaneous coronary intervention of calcified lesions with drug-eluting stents. *Int J Cardiol* 2017;231:61–7. https://doi.org/10.1016/j.ijcard.2016.12.150; PMID: 28040289.
- Madhavan MV, Tarigopula M, Mintz GS, Maehara A, et al. Coronary artery calcification: pathogenesis and prognostic implications. J Am Coll Cardiol 2014;63:1703–14. https://doi. org/10.1016/j.jacc.2014.01.017; PMID: 24530667.
- Bourantas CV, Zhang YJ, Garg S, et al. Prognostic implications of coronary calcification in patients with obstructive coronary artery disease treated by percutaneous coronary intervention: a patient-level pooled analysis of 7 contemporary stent trials. *Heart* 2014;100:1158–64. https://doi. org/10.1136/heartjnl-2013-305180; PMID: 24846971.
 Lee MS, Shah N. The impact and pathophysiologic
- Lee MS, Shah N. The impact and pathophysiologic consequences of coronary artery calcium deposition in percutaneous coronary interventions. J Invasive Cardiol 2016;28:160–7. PMID: 26301561.
- Kobayashi Y, Okura H, Kume T, et al. Impact of target lesion coronary calcification on stent expansion. *Circ J* 2014;78:2209–14. https://doi.org/10.1253/circj.CJ-14-0108; PMID: 25017740.
- Barbato E, Shlofmitz E, Milkas A, et al. State of the art: evolving concepts in the treatment of heavily calcified and undilatable coronary stenoses – from debulking to plaque modification, a 40-year-long journey. *EuroIntervention* 2017;13:696–705. https://doi.org/10.4244/EIJ-D-17-00473; PMID: 28844031.
- Reifart N, Vandormael M, Krajcar M, et al. Randomized comparison of angioplasty of complex coronary lesions at a single center. Excimer Laser, Rotational Atherectomy, and Balloon Angioplasty Comparison (ERBAC) Study. *Circulation* 1997;96:91–8. https://doi.org/10.1161/01.CIR.96.1.91; PMID: 9236422.
- de Waha S, Allali A, Buttner HJ, et al. Rotational atherectomy before paclitaxel-eluting stent implantation in complex calcified coronary lesions: two-year clinical outcome of the randomized ROTAXUS trial. *Catheter Cardiovasc. Interv* 2016;87:691–700. https:// doi.org/10.1002/ccd.26290; PMID: 26525804.
- Mauri L, Bonan R, Weiner BH, et al. Cutting balloon angioplasty for the prevention of restenosis: results of the Cutting Balloon Global Randomized Trial. *Am J Cardiol* 2002;90:1079–83. https://doi.org/10.1016/S0002-9149(02)02773-X. PMID: 12423707.
- Serruys PW, Katagiri Y, Onuma Y. Shaking and breaking calcified plaque: lithoplasty, a breakthrough in interventional armamentarium? *JACC Cardiovasc Imaging* 2017;10:907–11. https://doi.org/10.1016/i.icmg.2017.05.011: PMID: 28797413.
- https://doi.org/10.1016/j.jcrng.2017.05.011; PMID: 28797413.
 Brodmann M, Werner M, Brinton TJ, et al. Safety and performance of lithoplasty for treatment of calcified peripheral artery lesions. *J Am Coll Cardiol* 2017;70:908–10.
 https://doi.org/10.1016/j.jacc.2017.06.022; PMID: 28797363.
- Tripolino C, Grillo P, Tassone EJ, et al. A case of critical calcified innominate artery stenosis successfully

treated with the shockwave lithoplasty. *Clin Med Insights Case Rep* 2019;12:1179547619828707. https://doi. org/10.1177/1179547619828707. PMID: 30792580

- org/10.1177/1179547619828707; PMID: 30792580. 13. Brodmann M, Holden A, Zeller T. Safety And feasibility of intravascular lithotripsy for treatment of below-the-knee arterial stenoses. *J Endovasc Ther* 2018;25:499–503. https://doi. org/10.1177/1526602818783989; PMID: 29911480.
- Brodmann M, Werner M, Holden A, et al. Primary outcomes and mechanism of action of intravascular lithotripsy in calcified, femoropopliteal lesions: results of Disrupt PAD II. *Catheter Cardiovasc. Interv* 2019;93:335–42. https://doi. org/10.1002/ccd.27943; PMID: 30474206.
- Brinton TJ, Ali ZA, Hill JM, et al. Feasibility of shockwave coronary intravascular lithotripsy for the treatment of calcified coronary stenoses. *Circulation* 2019;139:834–6. https://doi.org/10.1161/CIRCULATIONAHA.118.036531; PMID: 30715944.
- Ali ZA, Brinton TJ, Hill JM, et al. Optical coherence tomography characterization of coronary lithoplasty for treatment of calcified lesions: first description. *JACC Cardiovasc Imaging* 2017;10:897–906. https://doi.org/10.1016/j.jcmg.2017.05.012; PMID: 28797412.
- Wong B, El-Jack S, Newcombe R, et al. Shockwave intravascular lithotripsy for calcified coronary lesions: first real-world experience. *J Invasive Cardiol* 2019;31:46–8. https:// doi.org/10.1016/j.hlc.2019.05.022: PMID: 30765621.
- doi.org/10.1016/j.hlc.2019.05.022; PMID: 30765621.
 18. Warisawa T, Goto S, Salazar CH, et al. Safety and feasibility of coronary lithotripsy supported by guide extension catheter for the treatment of calcified lesion in angulated vessel. *Cardiovasc Revasc Med* 2019. https://doi.org/10.1016/j. carrev.2019.02.014: PMID: 30786960; eoub ahead of press.
- Sgueglia GA, Gioffre G, Piccioni F, Gaspardone A. Slender distal radial five French coronary shockwave lithotripsy. *Catheter Cardiovasc Interv* 2019;94:395–8. https://doi.org/10.1002/ ccd.28296; PMID: 31001899.
- Kassimis G, Raina T, Kontogiannis N, et al. How should we treat heavily calcified coronary artery disease in contemporary practice? From atherectomy to intravascular lithotripsy. Cardiovasc Revasc Med 2019. https://doi.org/10.1016/j. carrev.2019.01.010; PMID: 3071477; epub ahead of press.
 Wong B, El-Jack S, Newcombe R, et al. Shockwave
- Wong B, El-Jack S, Newcombe R, et al. Shockwave intravascular lithotripsy of calcified coronary lesions in ST-elevation myocardial infarction: first-in-man experience. *J Invasive Cardiol* 2019;31:E73–5. PMID: 31034437.
- Neumann FJ, Sousa-Uva M, Ahlsson A, et al. 2018 ESC/EACTS guidelines on myocardial revascularization. *EuroIntervention* 2019;14:1435–534. https://doi.org/10.4244/EJJY19M01_01; PMID: 30667361.
- Lee MS, Shlofmitz E, Kaplan B, Shlofmitz R. Percutaneous coronary intervention in severely calcified unprotected left main coronary artery disease: initial experience with orbital atherectomy. J Invasive Cardiol 2016;28:147–50. https://doi. org/10.1016/j.jcin.2015.12.054; PMID: 27031936.
- Lee MS, Shlofmitz E, Park KW, et al. Orbital atherectomy of severely calcified unprotected left main coronary artery

the need for further adjuvant therapy with conventional devices and potential complications.

Conclusion

The Coronary IVL System is a promising new treatment modality to tackle moderate to severe calcified coronary lesions, with a high rate of success and a low risk of complications. Larger studies and longerterm clinical data are needed to confirm the safety and efficacy of this technique with special attention to the effects on cardiac conduction and vessel healing response. Randomised controlled clinical trials are required to evaluate its superiority against currently available calciummodifying devices.

Supplementary data

Supplementary Video 1: Shockwave Coronary Intravascular Lithotripsy System. Reproduced with permission from Shockwave Medical.

Supplementary Video 2: IVL Mechanism of Action (*In-vitro* Demonstration). Reproduced with permission from Shockwave Medical.

Supplementary Video 3: Hard on Hard Tissue (*In-vitro* Demonstration). Reproduced with permission from Shockwave Medical.

Supplementary Video 4: Soft on Soft Tissue (*In-vitro* Demonstration). Reproduced with permission from Shockwave Medical. ■

> disease: one-year outcomes. J Invasive Cardiol 2018;30:270–4. https://doi.org/10.1111/joic.12362; PMID: 28116818.
> 25. Wong B, El-Jack S, Khan A, et al. Treatment of heavily calcified

- Wong B, El-Jack S, Khan A, et al. Treatment of heavily calcified unprotected left main disease with lithotripsy: the first case series. J Invasive Cardiol 2019;31:E143–7. PMID: 31158810.
- Azzalini L, Bellini B, Montorfano M, Carlino M. Intravascular lithotripsy in chronic total occlusion percutaneous coronary intervention. *EuroIntervention*. 2019. https://doi.org/10.4244/EIJ-D-19-00175; PMID: 31012852; epub ahead of press.
- Yeoh J, Hill J, Spratt JC. Intravascular lithotripsy assisted chronic total occlusion revascularization with reverse controlled antegrade retrograde tracking. *Catheter Cardiovasc Interv* 2019;93:1295–7. https://doi.org/10.1002/ccd.28165; PMID: 30838746.
- Whiteside HL, Nagabandi A, Kapoor D. Safety and efficacy of stentablation with rotational atherectomy for the management of underexpanded and undilatable coronary stents. *Cardiovasc Revasc Med* 2019. https://doi.org/10.1016/j. carrev.2019.01.013, PMID: 30685339; epub ahead of press.
- Latib A, Takagi K, Chizzola G, et al. Excimer Laser LEsion modification to expand non-dilatable stents: the ELLEMENT registry. Cardiovasc Revasc Med 2014;15:8–12. https://doi. org/10.1016/j.carrev.2013.10.005; PMID: 24290659.
- Wilson A, Ardehali R, Brinton TJ, et al. Cutting balloon inflation for drug-eluting stent underexpansion due to unrecognized coronary arterial calcification. *Cardiovasc Revasc Med* 2006;7:185–8. https://doi.org/10.1016/j.carrev.2006.05.005; PMID: 16945827.
- Edes IF, Ruzsa Z, Szabo G, et al. Rotational atherectomy of undilatable coronary stents: stentablation, a clinical perspective and recommendation. *EuroIntervention* 2016;12:e632–5. https://doi.org/10.4244/EUV12I5A103; PMID: 27497363.
- Tovar Forero MN, Wilschut J, Van Mieghem NM, Daemen J. Coronary lithoplasty: a novel treatment for stent underexpansion. *Eur Heart J* 2019;40:221. https://doi. org/10.1093/eurheartj/ehy593; PMID: 30289452.
- Morabito G, Tripolino C, Tassone EJ, et al. A case of stent under-expansion due to calcified plaque treated with shockwave lithoplasty. *Cardiology* 2018;141:75–7. https://doi. org/10.1159/000493747; PMID: 30408797.
- Urbano Carrillo CA, Cano Garcia M, Munoz Jimenez LD. Coronary lithoplasty in the treatment of stent underexpansion. *Rev Esp Cardiol (Engl Ed)* 2019. https://doi. org/10.1016/j.rec.2019.04.009; PMID: 31160266; epub ahead of press.
- Alfonso F, Bastante T, Antuna P, et al. Coronary lithoplasty for the treatment of undilatable calcified de novo and in-stent restenosis lesions. *JACC Cardiovasc Interv* 2019;12:497–9. https:// doi.org/10.1016/j.jcin.2018.12.025; PMID: 30772288.
- Salazar C, Escaned J, Tirado G, Gonzalo N. Undilatable calcific coronary stenosis causing stent underexpansion and late stent thrombosis: a complex scenario successfully managed with intravascular lithotripsy. *JACC Cardiovasc Interv* 2019;12:1510–2. https://doi.org/10.1016/j.jcin.2019.02.010;

PMID: 30981572

- 37. Tassone EJ, Tripolino C, Morabito G, et al. When calcium gets tough, the tough cardiologist starts to play. Cardiology 2018;141:167-71. https://doi.org/10.1159/000495177; PMID: 30630180.
- Watkins S, Good R, Hill J, et al. Intravascular lithotripsy to treat a severely underexpanded coronary stent. *EuroIntervention* 38. 2019;15:124-5. https://doi.org/10.4244/EIJ-D-18-00780; PMID: 30295290.
- Ali ZA, McEntegart M, Hill JM, Spratt JC. Intravascular lithotripsy 39. for treatment of stent underexpansion secondary to severe coronary calcification. Eur Heart J 2018. https://doi.org/10.1093/ eurheartj/ehy747; PMID: 30462174; epub ahead of press. 40. McQuillan C, Alkhalil M, Johnston PW. A paced heart without
- a pacemaker. Eur Heart J 2019;40:819a. https://doi.org/10.1093/ eurhearti/ehy749; PMID: 30428000.
- 41. Cicovic A, Cicovic S, Wong B, et al. A quicker pace: shockwave lithotripsy pacing with electromechanical capture. JACC Cardiovasc Interv 2019;12:1739–40. https://doi.org/10.1016/j. jcin.2019.04.024; PMID: 31153842. Wilson SJ, Spratt JC, Hill J, et al. Coronary intravascular
- 42. lithotripsy is associated with a high incidence of

"shocktopics" and asynchronous cardiac pacing. EuroIntervention 2019. https://doi.org/10.4244/EIJ-D-19-00484; PMID: 31130523; epub ahead of press.

- 43 Quinn TA, Kohl P. Comparing maximum rate and sustainability of pacing by mechanical vs. electrical stimulation in the Langendorff-perfused rabbit heart. *Europace* 2016;18:iv85–93. https://doi.org/10.1093/europace/euw354; PMID: 28011835.
- Venuti G, D'Agosta G, Tamburino C, La Manna A. Coronary lithotripsy for failed rotational atherectomy, cutting balloon, 44. scoring balloon, and ultra-high-pressure non-compliant balloon. Catheter Cardiovasc Interv 2019;94:E111–5. https://doi. org/10.1002/ccd.28287; PMID: 31020765s. Chen G, Zrenner B, Pyxaras SA. Combined rotational
- 45. atherectomy and intravascular lithotripsy for the treatment of severely calcified in-stent neoatherosclerosis: a mini-review Cardiovasc Revasc Med 2018;20:819–21. https://doi.org/10.1016/j. carrev.2018.10.007; PMID: 30409500.
- Vainer J, Lux A, Ilhan M, et al. Smart solution for hard times: 46. successful lithoplasty of an undilatable lesion. *Neth Heart J* 2019;27:216–7. https://doi.org/10.1007/s12471-019-1261-2; PMID: 30868549.

- 47. López-Lluva MT, Jurado-Román A, Sanchez-Pérez I, et al. Shockwave: useful but potentially dangerous. JACC Cardiovasc Interv 2019;12:500–1. https://doi.org/10.1016/j icin.2018.12.035; PMID: 30772291.
- Soriano F, Veas N, Piccinelli E, Oreglia J. Coronary dissection 48. due to intravascular lithoplasty balloon rupture. *EuroIntervention* 2019;15:e558–9. https://doi.org/10.4244/EIJ-D-19-00383; PMID: 31036536.
- De Silva K, Roy J, Webb I, et al. A calcific, undilatable stenosis: lithoplasty, a new tool in the box? *JACC Cardiovasc Interv* 49. 2017;10:304-6. https://doi.org/10.1016/j.jcin.2016.11.048; PMID: 28183471.
- Rodriguez Costoya I, Tizon Marcos H, Vaquerizo Montilla B, et al. Coronary lithoplasty: initial experience 50. in coronary calcified lesions. Rev Esp Cardiol (Engl Ed) 2019;72:788-90. https://doi.org/10.1016/j.rec.2018.11.017; PMID: 30782450.
- 51. Legutko J, Niewiara Ł, Tomala M, et al. Successful Shockwave intravascular lithotripsy for severely calcified, undilatable lesion of the left anterior descending coronary artery in patient with recurrent myocardial infarction. Kardiol Pol 2019;77:723–5. https:// doi.org/10.33963/KP.14859; PMID: 31171761.