

Evaluation of a semi-automated Seegene PCR workflow with MTB, MDR, and NTM detection for rapid screening of tuberculosis in a low-prevalence setting

BRUNO VINCENT LUUKINEN,^{1,2}  RISTO VUENTO¹ and JARI JUHANI HIRVONEN^{1,3}

¹Department of Clinical Microbiology, Fimlab Laboratories, Tampere; ²Faculty of Medicine, University of Helsinki, Helsinki; and ³Mobidiag, Espoo, Finland

Luukinen BV, Vuento R, Hirvonen JJ. Evaluation of a semi-automated Seegene PCR workflow with MTB, MDR, and NTM detection for rapid screening of tuberculosis in a low-prevalence setting. APMIS 2020; 128: 406–413.

In areas of low tuberculosis (TB) prevalence, laboratory diagnosis of TB may essentially cover non-tuberculous mycobacteria (NTM) in addition to *Mycobacterium tuberculosis* (MTB). In this study, a semi-automated PCR workflow distinguishing MTB and NTM (Anyplex™ MTB/NTMe, Seegene) and subsequently detecting MTB isoniazid/rifampicin resistance (Allplex™ MTB/MDRe, Seegene) was evaluated for replacing smear microscopy of acid-fast bacilli as the rapid screening method for TB. With 279 clinical samples, 47 cultures positive for MTB and 76 for NTM, the Anyplex™ MTB/NTMe assay and smear microscopy showed equal sensitivities (49.6% vs 50.8%, respectively) but Anyplex™ MTB/NTMe was more sensitive for MTB (63.8% vs 25.6%) than for NTM (40.8% vs 64.5%). Allplex™ MTB/MDRe showed a slightly higher sensitivity of 68.1% for MTB (32/47 positive, n = 222). Antibiotic resistance profiles were correctly identified for all MTB isolates (one MDR isolate). Specificity was 100% for both assays. Anyplex™ MTB/NTMe detected all the 18 NTM species present in the study. The analytical performance of the evaluated high-throughput workflow was relatively weak compared to culture but potentially adequate as a rapid screening method analogous to smear microscopy with additional differentiation between TB, MDR-TB, and NTM.

Key words: Allplex; Anyplex; molecular diagnostics; tuberculosis; clinical microbiology; mycobacteriology.

Bruno Vincent Luukinen, Department of Clinical Microbiology, Fimlab Laboratories, PL 66, Tampere 331013, Finland.
e-mail: bruno.luukinen@fimlab.fi

The detection of *Mycobacterium tuberculosis* (MTB), the causative agent of tuberculosis (TB), is undeniably the number one priority when considering the laboratory diagnosis of mycobacteria. However, especially clinical settings of low TB-prevalence, routine laboratory diagnosis may include non-tuberculous mycobacteria (NTM) in addition to MTB.

It is well known that the clinical significance of isolating environmental NTM in clinical samples is highly unknown compared to the obligatory parasitic MTB (1). Indeed, clinical symptoms and radiological findings are considered a prerequisite for the diagnosis of NTM pulmonary disease (NTMPD) (2). As NTMPD and TB may be clinically indistinguishable but differ greatly in antimicrobial therapy, a significant added value of

detecting NTM lies in the exclusion of TB (3,4). Several cases of MDR-TB misdiagnosis due to NTM infection have been described in areas where NTM are not routinely analyzed (5,6). In countries where they are routinely analyzed, NTM make up a significant share of all mycobacterial laboratory findings and rates are systematically increasing (7,8).

Smear microscopy is still a widely used initial screening method of TB (9). Although inexpensive, the method does not identify different acid-fast bacilli. Here, in a low TB-prevalence setting in Finland, we evaluated the performance of two PCR assays, Anyplex™ MTB/NTMe (Seegene, Seoul, Korea) and Allplex™ MTB/MDRe (Seegene) using the automated NIMBUS IVD nucleic acid extraction and PCR setup system (Hamilton Company, Reno, NV, USA). The system was also surveyed as a combined workflow of rapid TB screening

Received 22 January 2020. Accepted 2 March 2020

analogous to smear microscopy but with additional analysis of identification of TB, MDR-TB, and NTM.

MATERIALS AND METHODS

Clinical material

A total of 279 clinical NALC-NaOH sedimented samples (stored at -70°C) and 35 culture-enriched (Löwenstein-Jensen agar or BD BACTEC™ MGIT™ growth indicator tubes (Becton, Dickinson and Company, Franklin Lakes, NJ, USA)) mycobacterial isolates were collected for this study during the year 2018 with 230 of pulmonary origin (215 sputum or tracheal aspirate, 15 bronchoalveolar lavage) and 49 of extrapulmonary origin (23 soft tissue, 15 urine, 10 pus, and 1 bone). The samples were collected, and the study was performed retrospectively in the Clinical Microbiology Laboratory of Fimlab Laboratories, Tampere, Finland.

Data for smear microscopy (auramine staining) and culture (Löwenstein-Jensen agar and BD BACTEC™ MGIT™ growth indicator tubes) were available for all but 15 urine samples with culture results only. Additionally, mycobacterial species identification (GenoType® Mycobacterium CM VER 2.0 assay, Hain Lifescience GmbH, Nehren, Germany) and antimicrobial susceptibility results (determined by National Institute for Health and Welfare, Finland) were available for cultured mycobacterial isolates.

Sample pretreatment

To ensure safe working outside a biosafety cabinet, and to reduce sample viscosity, a separate protocol inactivation was employed, previously described by Qi et al. (10) for the RealTime MTB assay (Abbott Molecular, Des Plaines, IL, USA). Briefly, samples were mixed 1:3 with inactivation reagent (0.6% sodium hydroxide [wt/vol], 60% isopropanol [vol/vol], and 1.8% Tween 20 [vol/vol]) and vortexed twice during a 1–24 h incubation at room temperature. To assess whether the inactivation protocol had any adverse effects to the workflow, a subset of samples was analyzed without the inactivation pretreatment (non-paired comparison to inactivated samples) and another smaller subset with and without the inactivation pretreatment (paired comparison). Samples for which no inactivation pretreatment was employed were known to be negative for MTB.

DNA extraction, amplification, and result analysis

Two different PCR assays were used, Anyplex™ MTB/NTMe (Seegene) and Allplex™ MTB/MDRe (Seegene). DNA was extracted with NIMBUS IVD (Hamilton Company, Reno, NV, USA) system using the STARMag 96 × 4 Universal Cartridge Kit (Seegene) and amplified with CFX96™ Real-Time PCR Detection System (Bio-Rad Laboratories, Hercules, CA, USA). Extraction and PCR setup were controlled with Seegene Launcher IVD (Seegene) and result analysis with Seegene Viewer IVD software (Seegene). All protocols were preset by the assay manufacturer. DNA concentrations were measured

using NanoDrop™ 2000 Spectrophotometer (Thermo-Fisher Scientific, Waltham, MA, USA).

Workflow analysis

To assess the turnaround time of the workflow, including Anyplex™ MTB/NTMe and subsequent testing of MTB-positive samples with Allplex™ MTB/MDRe, durations of different steps with realized sample batches were timed. To assess the range of assay durations, the instrument software was tested for analysis with a minimal and maximal number of samples, while still allowing for single batch analysis.

Statistical analysis

Statistical tests were employed, when appropriate, using the Prism software (GraphPad Software, San Diego, CA, USA). Mann–Whitney *U* test and Wilcoxon signed-rank test were used to compare non-paired and paired data of DNA concentrations, respectively. Fisher's exact test was used for comparing sample clotting during DNA extraction and PCR results with different sample types. Sensitivity, specificity, and positive (PPV) and negative predictive values (NPV) with 95% confidence intervals (CI) and medians with interquartile range were used to describe qualitative and quantitative results, respectively.

Ethics declaration

Leftover samples were analyzed retrospectively in this study, and no patient data were collected. The study design required no ethical committee approval. Informed consent was not required for this study.

RESULTS

Clinical material

One hundred twenty-three of the 279 analyzed clinical samples were culture-positive for mycobacteria, 47 with MTB (including 3 MDR and 2 *M. africanum* isolates) and 76 with NTM (consisting of 9 different NTM species). With the culture-enriched isolates included, a total of 18 NTM species were analyzed in the study.

Performance of the Anyplex™ MTB/NTMe assay

A total of 279 clinical samples were analyzed with the Anyplex™ MTB/NTMe assay. The assay showed overall sensitivity of 49.6% (61/123) with 67.2% for smear-positive (SP) and 30.5% for smear-negative (SN) samples. Separately, sensitivity was 63.8% (30/47) for MTB detection (SP: 91.7%, 11/12; SN: 53.1%, 17/32) and 40.8% (31/76) for NTM detection (SP: 61.2%, 30/49; SN: 3.7%, 1/27) (Table 1 and 2). The overall sensitivity was equal to that of smear microscopy (50.8%, 61/120) with reciprocal figures for smear microscopy's theoretical

Table 1. Agreement of Anyplex MTB/NTMe and Allplex MTB/MDRe with mycobacterial culture and smear microscopy results

	Smear microscopy		Anyplex MTB/NTMe						Allplex MTB/MDRe	
	AFB detection		MTB detection		NTM detection		MTB/NTM detection		MTB detection	
	Positive	Negative	Positive	Negative	Positive	Negative	Positive	Negative	Positive	Negative
Culture (n = 279)										
Positive	61	59	30	17	31	45	61	62	32	15
Negative	5	139	3 ¹	229	0	203	3 ¹	153	4 ¹	171
Smear microscopy ²										
Positive			11	1	30	19	41	20	11	1
Negative			17	15	1	26	18	41	19	13

AFB, acid-fast bacilli.

¹All cases recently culture-positive for MTBC.

²All culture-positive.

sensitivities for MTB (25.6%, 12/47) and NTM (64.5%, 49/76) detection whatsoever. True assay specificity and PPV were 100%, as the three observed culture-negative, PCR-positive cases were all recently culture-positive for MTB. NPV for MTB/NTM, MTB, and NTM detection were 71.2%, 93.1% (229/246), and 81.9% (203/248), respectively.

Of the nine NTM species present in clinical samples, *M. avium*, *M. chelonae*, and *M. simiae* were detected by AnyplexTM MTB/NTMe (Table 3). *M. xenopi* was the only NTM species not identified in any smear-positive clinical sample. Examining PCR data of culture-enriched isolates indeed showed that the detection of *M. xenopi* was relatively weak, showing late PCR cycles of 35.19 and 36.06 compared to the median of the whole culture-enriched sample set (23.30, 21.62–26.21). All the 18 NTM species present in the study were detectable in culture-enriched samples (Table 4).

Performance of the AllplexTM MTB/MDRe assay

Two hundred twenty-two of the available 279 clinical samples were analyzed with the AllplexTM MTB/MDRe assay. The AllplexTM MTB/MDRe assay showed sensitivity of 68.1% (32/47) (SP: 91.7%, 11/12; SN: 59.4%, 19/32) and NPV of 91.9% (171/186) for MTB detection (Table 1 and 2). As with AnyplexTM MTB/NTMe, true specificity and PPV were 100%. 29/30 (96.7%) of the samples analyzed as MTB-positive with AnyplexTM MTB/NTMe were also analyzed as positive with AllplexTM MTB/MDRe.

There were three MDR-MTB-positive samples in the study, all culture-positive, smear-negative, none detected by the AllplexTM MTB/MDRe assay (MTB negative). However, both MTB and MDR were detected from a culture-enriched sample. No false-positive MDR results were detected.

Table 2. Performance data of AnyplexTM MTB/NTMe and AllplexTM MTB/MDRe assays

Method	Sensitivity (CI), %			Specificity (CI), %	PPV (CI), %	NPV (CI), %
	Smear+	Smear–	Total			
Smear microscopy						
AFB detection			50.8 (42.0–59.6)	96.5 (91.9–98.7)	92.4 (83.1–97.1)	70.2 (63.5–76.2)
Anyplex MTB/NTMe						
MTB detection	91.7 (62.5–100)	53.1 (36.5–69.1)	63.8 (49.5–76.1)	98.7 ¹ (96.1–99.7)	90.9 ¹ (75.7–97.6)	93.1 (89.2–95.7)
NTM detection	61.2 (47.2–73.6)	3.7 (0–19.8)	40.8 (30.4–52.0)	100 (97.8–100)	100 (86.9–100)	81.9 (76.6–86.2)
MTB/NTM detection	67.2 (54.7–77.7)	30.5 (20.2–43.2)	49.6 (40.9–58.3)	98.1 ¹ (94.2–99.6)	95.3 ¹ (86.6–98.9)	71.2 (64.8–76.8)
Allplex MTB/MDRe						
MTB detection	91.7 (62.5–100)	59.4 (42.2–74.5)	68.1 (53.8–79.7)	97.7 ¹ (94.1–99.3)	88.9 ¹ (74.1–96.2)	91.9 (87.0–95.1)

AFB, acid-fast bacilli.

¹True specificity and PPV were 100% as all culture-negative, PCR-positive cases had been recently culture-positive for MTBC.

Table 3. Array of non-tuberculous mycobacteria present in the study: clinical samples

Species	Smear+			Smear–			Total		
	Culture+	PCR+	%	Culture+	PCR+	%	Culture+	PCR+	%
<i>M. avium</i>	40	25	62.5	15	1	6.7	55	26	47.3
<i>M. intracellulare</i>				5	0	0	5	0	0
<i>M. xenopi</i>	3	0	0	1	0	0	4	0	0
<i>M. chelonae</i>	3	3	100				3	3	100
<i>M. simiae</i>	3	2	66.7				3	2	66.7
<i>M. gordonae</i>				2	0	0	2	0	0
<i>M. szulgai</i>				2	0	0	2	0	0
<i>M. fortuitum</i>				1	0	0	1	0	0
<i>M. mucogenicum</i>				1	0	0	1	0	0
Total	49	30	61.2	27	1	3.7	76	31	40.8

Table 4. Array of mycobacteria present in the study: culture-enriched isolates

Cultured isolates Species	n	Anyplex MTB/NTMe		Other species with positive identification tested by manufacturer
		MTB	NTM	
<i>M. tuberculosis</i>	2	2		<i>M. bovis</i> BCG
<i>M. avium</i>	3		3	<i>M. asiaticum</i>
<i>M. abscessus</i>	2		2	<i>M. fallax</i>
<i>M. chelonae</i>	2		2	<i>M. massiliense</i>
<i>M. fortuitum</i>	2		2	<i>M. neoaurum</i>
<i>M. gordonae</i>	2		2	<i>M. paraseoulense</i>
<i>M. interjectum</i> ¹	2		2	<i>M. phlei</i>
<i>M. intracellulare</i>	2		2	<i>M. seoulense</i>
<i>M. kansasii</i>	2		2	<i>M. smegmatis</i>
<i>M. lentiflavum</i> ¹	2		2	<i>M. vaccae</i>
<i>M. szulgai</i>	2		2	
<i>M. terrae</i>	2		2	
<i>M. xenopi</i> ¹	2		2	
<i>M. celatum</i> ¹	1		1	
<i>M. malmoense</i> ¹	1		1	
<i>M. marinum</i> ¹	1		1	
<i>M. mucogenicum</i>	1		1	
<i>M. simiae</i> ¹	1		1	
<i>M. schimoidei</i> ¹	1		1	
NTM total	31		31	
Other AFB	2		1	

AFB, acid-fast bacilli (non-mycobacterial).

¹Assay specificity not tested by manufacturer.

Effect of the inactivation protocol in DNA extraction and PCR results

The inactivation protocol had a significant effect on DNA extraction yields ($p = 0.003$) when compared between a subset of samples with the inactivation protocol employed (11.8, 9.5–16 ng/ μ L; $n = 96$) and a different subset of samples without the inactivation protocol employed (13.9, 10–32 ng/ μ L; $n = 96$) (non-paired subsets). The observed difference was strengthened by a smaller subset of samples analyzed both with and without inactivation pretreatment ($p = 0.0001$) (paired subset) (Fig. 1A). However, the analysis of paired samples showed a mean 2.2-fold (1.6–

2.6, $n = 14$) shift in DNA concentration, whereas the inactivation protocol diluted the sample four-fold. A similar comparison of the same subset of paired samples was done with the AnyplexTM MTB/NTMe assay. Here, only a slight shift in PCR cycle levels was seen with no significant difference ($p = 0.14$) (Fig. 1B). No PCR abortion occurred in the study. However, sample clotting events in the automated DNA extraction were significantly more frequent with the inactivated samples (13.6%, 25/184) compared to the non-inactivated samples (3.2%, 3/95) ($p = 0.006$). This increased the need for manual pipetting during the DNA extraction protocol.




Fig. 1. Effect of the inactivation protocol on sample quality and results on DNA extraction yield (A) and PCR cycles (B). *Statistically significant difference.

Semi-automated workflow

The workflow of rapid TB screening with Anyplex™ MTB/NTME and subsequent testing of MTB-positive samples with Allplex™ MTB/MDRe is illustrated in Fig. 2. Assuming a 96-well plate format, the workflow allowed analysis of 72 samples with 45 MTB-positive samples in one cycle. The highest walkaway time was achieved with batches of no more than 40 clinical samples (not to mix with analysis of already extracted DNA). More samples resulted in two separate interruptions of the DNA extraction protocol, as the instrument required more pipette tips and an empty waste bag. Occasional sample clots further reduced walkaway time and increased extraction time. There was a separate interrupt in the extraction protocol for manually transferring samples with clots. Without these interruptions, the turnaround time of a 40-sample batch was approximately 6 h from decontaminated samples to MTB/NTM (inactivation: 1.5 h, DNA extraction and PCR setup: 2.25 h, amplification and result analysis: 2.25 h, and hands-on time included). The turnaround time for subsequent MDR-analysis of five MTB-positive samples took approximately 3.5 h (PCR setup: 0.25 h, amplification: 3 h, result analysis, and hands-on time: 0.25 h), making up a total turnaround time of 9.5 h.

DISCUSSION

Rapid screening of TB is still widely based on smear microscopy of acid-fast bacilli – a method that could be replaced by more automated and objective technology such as PCR. It is clear that a PCR-based method is more expensive per test than




Fig. 2. Illustration of the suggested semi-automated PCR workflow. Vertical arrows show workflow to result with live (dark) and dead (light) mycobacteria. Hands-on time excluded.

a nonspecific staining method. Additionally, a multi-stage PCR workflow such as the one evaluated in this study requires a high level of infrastructure and therefore does not suit small laboratories or low-income settings. Compared to smear microscopy, however, the system greatly increases the sensitivity of rapid screening of TB as was shown in this study. It also reduces the need for a high-expertise-demanding method, and increases assay

reproducibility with respect to workflow and result interpretation. With MTB, smear microscopy has the status as a determinant of TB infectivity, although there are also studies where smear-negative cases have been a significant source of infection (11,12). No similar practice has been described for NTM (2). We previously showed that in this specific clinical setting of low TB-prevalence and low population, an average of one sample was positive for acid-fast bacilli per day (13). In such a clinical setting, it would be easy to perform a subsequent smear microscopy analysis for PCR-positive samples. Thus, the status as an infectivity determinant does not essentially impair the applicability of PCR screening.

The evaluated semi-automated system is not exclusive for TB. The proposed TB screening workflow can be employed to complement a uniform and comprehensive molecular diagnostic system. This also allows the possibility of assigning automated work steps for other laboratory personnel not working in a TB laboratory. Here, the evaluated workflow also enabled a cost-efficient reflex strategy of antibiotic resistance testing where MTB-positive samples could be specifically selected for subsequent antibiotic resistance analysis without the need for overlapping analysis of all samples (13). As a follow-up test for MTB-positive samples with Anyplex™ MTB/NTMe assay, the Allplex™ MTB/MDRe assay provided reliable results, although more resistant MTB isolates would have been needed to fully evaluate the assay's MDR feature with 7 isoniazid and 18 rifampicin resistance determining mutation targets.

Moreover, the evaluated semi-automated system is not exclusive for TB. As Seegene holds an array of different molecular diagnostic assays available for its automated systems, the proposed TB screening workflow can be employed to complement a uniform and comprehensive molecular diagnostic system. This also allows the possibility of assigning automated work steps for other laboratory personnel not working in a TB laboratory. Here, the evaluated workflow enabled a cost-efficient reflex strategy of antibiotic resistance testing where MTB-positive samples could be specifically selected for subsequent antibiotic resistance analysis without the need for overlapping analysis of all samples. The system is analogous to the RealTime MTB system employed by Abbott, but, as previously shown, the Seegene TB assay selection may currently better suit the needs of low TB-prevalence, low population settings (13). The use of the Seegene Launcher software allowed easy and fluent control of samples, as a single sample worklist could be used from extraction to result analysis and follow-up testing

of selected samples. As a follow-up test for MTB-positive samples with Anyplex™ MTB/NTMe assay, the Allplex™ MTB/MDRe assay provided reliable results, although more resistant MTB isolates would have been needed to fully evaluate the assay's MDR feature with 7 isoniazid and 18 rifampicin resistance determining mutation targets.

Considering analytical performance, the two evaluated PCR assays showed overall sensitivities close to smear microscopy. While such figures are far from definitive laboratory diagnosis and therefore do not compete with culture, the relatively low sensitivity may be adequate for replacing smear microscopy as the initial screening method. As studies even for similar PCR assays have been shown to have higher performance figures (14,15), one must consider the possible adverse effects of the inactivation protocol or sample storage on the quality of the results. In addition, performance evaluation studies of Seegene TB assays are generally performed using a manual extraction protocol, and to our knowledge, there are no previous studies assessing the performance of the NIMBUS IVD in molecular diagnosis of TB. The inactivation protocol used in this study seemed, in fact, not to have any apparent effect on assay results despite sample dilution. However, as sample clotting occurred significantly more often with inactivated samples, this step seemed to be redundant in terms of reducing sample viscosity. In terms of laboratory safety, though, the process requires a separate MTB killing step as a sample-handling instrument may not prevent pathogen exposure if not situated inside a biosafety cabinet or a fume hood with an accessory HEPA filter. This is exceptionally important when handling MTB as laboratory transmission due to prolonged exposure may be very difficult to identify due to the slow nature of the MTB course of infection (16).

The Anyplex™ MTB/NTMe assay showed poorest performance on smear-negative NTM-positive samples but also specifically on *M. xenopi*, even with smear-positive samples. Depending on geographical location, *M. xenopi* is one of the species most often reported to cause NTM pulmonary disease, in addition to *M. abscessus*, *M. avium* complex, *M. fortuitum*, *M. kansasii*, and *M. malmoense* (4). Even though the regional incidence of NTM has a great diversity, *M. xenopi* is globally the third most common NTM identified (8%) after the *M. avium* complex (47%) and *M. gordonae* (11%) (17-19). Closer inspection of the test insert revealed that the manufacturer had in fact not tested specificity to *M. xenopi*. Other relevant species assayed in this study, but not by the manufacturer, were *M. celatum*, *M. interjectum*, *M. lentiflavum*,

M. malmoense, *M. marinum*, *M. simiae*, and *M. schimoidei*. Except for *M. marinum*, which is the causative agent of a superficial infection known as fish tank granuloma, all of the mentioned NTM species have been reported to cause pulmonary disease (20). Within tuberculous mycobacteria, information on specificity to *M. africanum* also appeared to be lacking. Regionally important in West Africa, testing assay specificity for *M. africanum* was highly relevant even though different MTB complex species have been recently proposed to comprise a single species of MTB (21). Although only with individual isolates, specificity of the Anyplex™ MTB/NTMe assay for these eight NTM species and one MTB variant was shown in this study with, as far as we are aware, no previous demonstration elsewhere.

Most importantly, the Anyplex™ MTB/NTMe assay showed reliable resolution between tuberculous and non-tuberculous mycobacteria with no observed false-positive results. This was based on the relatively higher threshold for the generic mycobacterial target (i.e., lower sensitivity for NTM) than for the MTB-specific target, so that MTB could not be falsely reported as NTM. However, this configuration did not compromise the sensitivity of MTB detection, which was superior to that of smear microscopy. The assays amplify the known gene targets IS6110 and Mpb64 for MTB detection, and a panmycobacterial 16S rRNA gene target for NTM detection.

The Anyplex™ MTB/NTMe and Allplex™ MTB/MDRe assays accompanied with the NIMBUS IVD system showed relatively low but adequate overall sensitivity to be used as an initial screening method for mycobacteria analogous to smear microscopy. However, future research assessing the cost efficiency of PCR compared to smear microscopy is vital.

We thank Seegene Inc. and Triolab Ltd. for providing the NIMBUS IVD instrument and accompanying labware, Allplex™ MTB/MDRe kits, and product support by Aliisa Mäkinen. The study was funded by Tampere Tuberculosis Foundation and Finnish Anti-Tuberculosis Foundation.

CONFLICTS OF INTEREST

None to declare.

REFERENCES

- Piersimoni C, Scarparo C. Pulmonary infections associated with non-tuberculous mycobacteria in immunocompetent patients. *Lancet Infect Dis* 2008;8:323–34.
- Griffith DE, Aksamit T, Brown-Elliott BA, Catanzaro A, Daley C, Gordin F, et al. An official ATS/IDSA statement: diagnosis, treatment, and prevention of nontuberculous mycobacterial diseases. *Am J Respir Crit Care Med* 2007;175:367–416.
- Tortoli E. Clinical manifestations of nontuberculous mycobacteria infections. *Clin Microbiol Infect* 2009;15:906–10.
- Prevots DR, Marras TK. Epidemiology of human pulmonary infection with nontuberculous mycobacteria a review. *Clin Chest Med* 2015;36:13–34.
- Maiga M, Siddiqui S, Diallo S, Diarra B, Traoré B, Shea YR, et al. Failure to recognize nontuberculous mycobacteria leads to misdiagnosis of chronic pulmonary tuberculosis. *PLoS ONE* 2012;7:e36902.
- Shahraki AH, Heidarieh P, Bostanabad SZ, Khosravi AD, Hashemzadeh M, Khandan S, et al. "Multidrug-resistant tuberculosis" may be nontuberculous mycobacteria. *Eur J Intern Med* 2015;26:279–84.
- Bar-On O, Mussaffi H, Mei-Zahav M, Prais D, Steuer G, Staffler P, et al. Increasing nontuberculous mycobacteria infection in cystic fibrosis. *J Cyst Fibros* 2015;14:53–62.
- Donohue MJ. Increasing nontuberculous mycobacteria reporting rates and species diversity identified in clinical laboratory reports. *BMC Infect Dis* 2018;18:163.
- World Health Organization. Global tuberculosis report 2019. Geneva: World Health Organization, 2019.
- Qi C, Wallis C, Pahalawatta V, Frank A, Ramdin N, Viana R, et al. Evaluation of the efficiency of the sample inactivation reagent in the Abbott RealTime MTB assay for inactivation of *Mycobacterium tuberculosis*. *J Clin Microbiol* 2015;53:3001–2.
- Behr M, Warren S, Salamon H, Hopewell P, de Leon A, Daley C, et al. Transmission of *Mycobacterium tuberculosis* from patients smear-negative for acid-fast bacilli. *Lancet* 1999;353:444–9.
- Nolan C, Blumberg H, Taylor Z, Bernardo J, Brennan P, Dunlap N, et al. American thoracic society/centers for disease control and prevention/infectious diseases society of America: controlling tuberculosis in the United States. *Am J Respir Crit Care Med* 2005;172:1169–227.
- Luukinen BV, Vuento R, Hirvonen JJ. Evaluation of two tuberculosis PCR assays for routine use in a clinical setting of low population and low tuberculosis prevalence. *APMIS* 2019;127:462–7.
- Lim J, Kim J, Kim JW, Ihm C, Sohn Y, Cho H, et al. Multicenter evaluation of Seegene Anyplex TB PCR for the detection of *Mycobacterium tuberculosis* in respiratory specimens. *J Microbiol Biotechnol* 2014;24:1004–7.
- Sali M, De Maio F, Caccuri F, Campilongo F, Sanguinetti M, Fiorentini S, et al. Multicenter evaluation of anyplex plus MTB/NTM MDR-TB assay for rapid detection of *Mycobacterium tuberculosis* complex and multidrug-resistant isolates in pulmonary and extrapulmonary specimens. *J Clin Microbiol* 2016;54:59–63.
- Singh K. Laboratory-acquired infections. *Clin Infect Dis* 2009;49:142–7.
- Hoefsloot W, Van Ingen J, Andrejak C, Angeby K, Bauriaud R, Bemer P, et al. The geographic diversity

- of nontuberculous mycobacteria isolated from pulmonary samples an NTM-NET collaborative study. *Eur Respir J* 2013;42:1604–13.
18. van der Werf MJ, Kodmon C, Katalinic-Jankovic V, Kummik T, Soini H, Richter E, et al. Inventory study of non-tuberculous mycobacteria in the European Union. *BMC Infect Dis* 2014;6(14):62.
 19. Yu X, Liu P, Liu G, Zhao L, Hu Y, Wei G, et al. The prevalence of non-tuberculous mycobacterial infections in mainland China: systematic review and meta-analysis. *J Infect* 2016;73:558–67.
 20. Daley CL, Griffith DE. Pulmonary non-tuberculous mycobacterial infections. *Int J Tuberc Lung Dis* 2010;14:665–71.
 21. Riojas MA, McGough KJ, Rider-Riojas CJ, Rastogi N, Hazbon MH. Phylogenomic analysis of the species of the *Mycobacterium tuberculosis* complex demonstrates that *Mycobacterium africanum*, *Mycobacterium bovis*, *Mycobacterium caprae*, *Mycobacterium microti* and *Mycobacterium pinnipedii* are later heterotypic synonyms of *Mycobacterium tuberculosis*. *Int J Syst Evol Microbiol* 2018;68:324–32.