



Published in final edited form as:

Brachytherapy. 2012 ; 11(1): 47–52. doi:10.1016/j.brachy.2011.07.002.

The American Brachytherapy Society Treatment Recommendations for Locally Advanced Carcinoma of the Cervix Part II: High Dose-Rate Brachytherapy

Akila N. Viswanathan, M.D., MPH^{*}, Sushil Beriwal, M.D.^{‡‡}, Jennifer De Los Santos, MD^{††}, D. Jeffrey Demanes, M.D.[‡], David Gaffney, M.D.[§], Jorgen Hansen, MS^{*}, Ellen Jones, MD^{§§}, Christian Kirisits, Sc.D.[†], Bruce Thomadsen, Ph.D.[#], and Beth Erickson, M.D.^{**}

^{*}Brigham and Women's Hospital and Dana-Farber Cancer Institute, Harvard Medical School, Boston, Massachusetts, Boston, MA, USA

^{‡‡}University of Pittsburgh School of Medicine, Pittsburgh, PA, USA

^{††}University of Alabama, Birmingham, AL, USA

[‡]David Geffen School of Medicine UCLA, Los Angeles, CA, USA

^{**}Medical College of Wisconsin, Milwaukee, WI, USA

[§]Huntsman Cancer Center, University of Utah, Salt Lake City, UT, USA

^{§§}University of North Carolina, Chapel Hill, NC, USA

[†]Medical University of Vienna, Vienna, Austria

[#]University of Wisconsin, Madison, WI, USA

Abstract

Purpose—This report presents the 2011 update to the American Brachytherapy Society (ABS) high-dose-rate (HDR) brachytherapy guidelines for locally advanced cervical cancer.

Methods—Members of the American Brachytherapy Society (ABS) with expertise in cervical cancer brachytherapy formulated updated guidelines for HDR brachytherapy using tandem and ring, ovoids, cylinder or interstitial applicators for locally advanced cervical cancer were revised based on medical evidence in the literature and input of clinical experts in gynecologic brachytherapy.

Results—The Cervical Cancer Committee for Guideline Development affirms the essential curative role of tandem-based brachytherapy in the management of locally advanced cervical cancer. Proper applicator selection, insertion, and imaging are fundamental aspects of the procedure. Three-dimensional imaging with magnetic resonance or computed tomography or radiographic imaging may be used for treatment planning. Dosimetry must be performed after each insertion prior to treatment delivery. Applicator placement, dose specification and dose fractionation must be documented, quality assurance measures must be performed, and follow-up information must be obtained. A variety of dose/fractionation schedules and methods for integrating brachytherapy with external-beam radiation exist. The recommended tumor dose in 2 Gray (Gy) per fraction radiobiologic equivalence (EQD2) is 80–90 Gy, depending on tumor size at the time of brachytherapy. Dose limits for normal tissues are discussed.

Address correspondence to: Akila N. Viswanathan, MD, MPH, Brigham and Women's Hospital/Dana-Farber Cancer Institute, Department of Radiation Oncology, 75 Francis Street, ASB 1, L2, Boston, MA 02115, Tel: 617-732-6331, Fax: 617-278-6988, aviswanathan@iroc.harvard.edu.

Conclusion—These guidelines update those of 2000 and provide a comprehensive description of HDR cervical cancer brachytherapy in 2011.

Introduction

Brachytherapy is an important component in the curative management of carcinoma of the cervix, and significantly improves survival.^{1, 2} High-dose-rate (HDR) and low dose-rate (LDR) brachytherapy appear to be relatively equivalent treatments in terms of survival outcomes based on existing retrospective and prospective studies.^{3–11} Advantages of HDR brachytherapy include opportunities for outpatient treatment, avoidance of exposure to staff from the radiation source, consistent and reproducible applicator positioning, and dose optimization attained with a variable dwell-time stepping source.³ Virtually all modern clinical trials for cervical cancer allow either HDR or LDR brachytherapy.

The use of HDR brachytherapy for cervical cancer has substantially increased over the past 10 years in the U.S. and internationally. The most recent Quality Research in Radiation Oncology (QRRO, formerly Patterns of Care) survey from 2007–2009 shows that 62% of surveyed facilities use HDR compared to 13% in the 1996–1999 survey.¹² A total of 85% of respondents to surveys in the U.S.¹³ and internationally¹⁴ use HDR brachytherapy. Nevertheless, with HDR brachytherapy, there is significant variation of the total tumor dose, the dose delivered per fraction and the proportion of tumor dose delivered with external-beam radiotherapy (EBRT) versus brachytherapy.¹⁴

Given the potential for short- and long-term injury to normal tissues from large HDR doses per treatment, the radiation oncologist must carefully assess and minimize normal-tissue doses administered per fraction, and must calculate the summative total dose of EBRT and brachytherapy. In order to assess the normal tissue doses per fraction accurately, computer-assisted tomography (CT) or magnetic resonance imaging (MRI) with the brachytherapy apparatus in place is recommended.

This article will present current concepts in HDR brachytherapy for cervical cancer including three-dimensional (3D) image-based dose-specification methods and review standard practice recommendations.

Methods

Gynecologic radiation oncology experts in the U.S. were surveyed regarding their willingness to serve as authors for these guidelines. Those responding affirmatively reviewed and updated the 2000 guidelines of the American Brachytherapy Society (ABS).¹⁵ These authors evaluated the relevant literature, identified established and controversial topics via conference calls, and supplemented this information with their clinical experience in order to formulate the current guidelines. A consensus decision was made to integrate strategies utilizing 3D image-guidance when possible.

This report was reviewed and approved by the Board of Directors of the ABS.

Results

1.0 Treatment Issues with HDR Brachytherapy

1.1 External-beam radiation therapy issues related to HDR brachytherapy—

Treatment with EBRT and brachytherapy should be completed in less than 8 weeks, as better local tumor control and survival can be expected with relatively shorter treatment courses.^{16, 17} The HDR brachytherapy may be interdigitated with EBRT to shorten the total treatment duration, with the latter typically given in 1.8-Gy fractions to 45 Gy. Many

institutions administer as much EBRT as possible first to minimize the amount of residual disease, ensure that the lymph node regions of the pelvis receive 5 days of EBRT per week for as long as possible, administer concurrent chemotherapy for a minimum of 5 consecutive weeks, and improve brachytherapy geometry due to tumor shrinkage increasing the distance between the tumor and the organs at risk (OAR). Others facilities elect to administer the first brachytherapy fraction early in the course of EBRT and treat one fraction per week, with brachytherapy not given on the same day as EBRT, in order to minimize treatment duration. For patients with large bulky tumors, commencing the treatment too early and specifying the dose to point A may underdose the tumor volume leading to poor local control.¹⁰ In the United States, the most common HDR intracavitary regimen prescribes 2 fractions per week for a total of 5 fractions.¹⁴ The ABS recommends that additional radiation to the parametria/nodes via a boost may be administered on non-brachytherapy days.

1.2 Chemotherapy issues unique to HDR brachytherapy—The ABS recommends the use of concurrent cisplatin based chemotherapy for patients with adequate renal function. When administering weekly cisplatin, the 5th and 6th dose of chemotherapy may fall during weeks when HDR brachytherapy commences. Though no data support an increase in toxicity,³ given the large fraction sizes utilized with HDR, the ABS recommends that chemotherapy not be administered on a brachytherapy day but rather on an EBRT day, due to the potential for increased complications due to normal-tissue sensitization.

2.0 Treatment Planning

2.1 Optimization issues specific to HDR intracavitary brachytherapy—Adequate geometry of the implant is imperative regardless of the simulation method. Incorrect placement of the applicator will negatively impact disease-free survival, increasing rates of local recurrence and often toxicity.¹⁸ Optimization of brachytherapy will not compensate for poor applicator placement.

A treatment plan should be generated by a qualified physicist or trained brachytherapy dosimetrist in collaboration with the treating radiation oncologist. The term *optimization* refers to the sophisticated process of achieving certain dose values at points or volumes within the implant; it is not the simple generation of a standard dose distribution by using fixed dose points located around the applicator. With conventional LDR brachytherapy, the shape of the dose distribution is hard to customize because of the few sources used (usually three in the tandem and one each in ovoids) and the limited number of source strengths. HDR brachytherapy allows more precise shaping of the dose distribution to the extent desired by the radiation oncologist. Some institutions use a squared distribution conforming to the cervix while others use a narrow tapered distribution that extends further into the uterus. Still others attempt to match the physical distribution of the LDR brachytherapy applicators, even though that produces a very different biological dose distribution.

Achieving an acceptable dose distribution with HDR brachytherapy requires both proper insertion of the appliance and a good optimization process. With 3D dosimetry, matching the dose distribution to the high-risk clinical target volume (HR-CTV) while simultaneously avoiding the OAR can be challenging. Two factors complicate the physical aspect of this challenge: throughout the history of cervical brachytherapy, the dose to the tumor, as defined by the HR-CTV, was unknown; and, increasing the weight of a source pushes the dose in all directions, towards OAR as well as the target. Optimization should be performed with caution by observing changes in the dose, dose/volume parameters and the spatial dose distribution that results from the modified loading pattern. The exclusive use of dose-volume-histogram (DVH)-based parameters to select a source loading is not recommended because substantial and perhaps undesirable changes in the spatial dose distribution may

occur. Hot or cold spots in the target region and in non-contoured OAR, such as the vagina, connective tissue, nerves, vessels or the ureters, may result. Importantly, in 3D imaging, the spatial dose distribution should be analyzed carefully for the location of cold and hot spots within the HR-CTV. Displaying isodose lines higher than the 100% isodose line may be important in recognizing and altering regions of high dose.

Optimization in the 2000 ABS Guidelines referred to setting lateral dose points adjacent to the applicator based on radiographic localization. With 3D-imaging, optimization refers to starting with a customary loading of the full length of the tandem and the vaginal applicator (ovoids, ring or cylinder), then modifying the dwell positions to reduce the dose to the OAR and ensure maximal tumor coverage; this results in differences in specification and reporting. For example, a dose of 5.5 Gy may be specified to a 3D-imaging-contoured target of 50-mm width at the level of point A. In order to fully cover the target, one approach is to define two dose points 25 mm from the tandem and normalize the 100% isodose line to these points. In this case a dose of 5.5 Gy is specified to the target while the dose at point A will be greater than 5.5 Gy. In daily clinical practice, the planning aims sometimes cannot be achieved due to the dose limits for the OAR. In such cases, the initially planned dose values should be decreased and an optimal compromise reached between tumor and OAR goals.

For the tandem applicator with needles (Figure 1), evaluation of the spatial dose distribution through the whole implant, including each needle, in addition to DVH values, becomes even more important. The balance of dose delivered through each needle should also be evaluated in order to avoid undesired high-dose regions in the adjacent tissues, such as the vagina, ureters, connective tissues and the OAR.¹⁹ A reproducible and safe approach is to first optimize dwell-time for the intracavitary part of the implant taking into account OAR primarily, without activating the needle positions. The missing coverage of the CTV is compensated for in a second step by fine-tuning the overall dose distribution with activation and direct adjustment of the dwell times in the needles. With inverse or graphical optimization, the dwell times of the intracavitary and interstitial parts should be controlled by the physicist, since most optimization algorithms do not take into account the spatial dose distribution. In general, approximately 10–20% of the total dwell time is linked to source positions in the needles, and most of the dose should be delivered through the tandem/ring or tandem/ovoid.

In interstitial brachytherapy, the target volume is typically larger than with intracavitary. The desired dose distribution to the central core of an interstitial implant, where needles may lie in close proximity to the tandem sources and the cervical and paracervical tumor, also differs from an intracavitary implant. In contrast, at the periphery of the implant the needles are in close proximity to the OAR and dose is necessarily reduced. During the optimization process, dwell positions and dwell times will be determined to deliver the intended dose. As with all volume implants, one point dose or fraction size cannot adequately describe the implant.

2.2 Dose Calculations for HDR brachytherapy—The radiobiology of HDR brachytherapy and the use of the linear-quadratic model to convert HDR to LDR doses were discussed in detail in the 2000 ABS recommendations and in recent studies. A worksheet is available for download from the ABS website to facilitate conversion of HDR fractionations into biologically equivalent doses in 2-Gy fractions – normalized therapy doses (NTD) or EQD2. At the time of this publication, the website is www.americanbrachytherapy.org/guidelines.html. These worksheets, however, are for theoretical guidance and should not replace the empirical observations or judgment of physicians experienced with HDR brachytherapy.

3.0 Dose recommendations for HDR brachytherapy

Recommendations for dose depend on the methodology followed for treatment planning. In the United States, the most commonly used regimens are 45 Gy EBRT to the pelvis (possibly with a sidewall boost) with concurrent cisplatin-based chemotherapy and either 5.5 Gy per fraction for 5 fractions (for patients treated with concurrent chemotherapy who have had either a complete response or have <4 cm of residual disease) or 6 Gy for 5 fractions (for patients with tumors >4 cm after EBRT). Over the past decade, the most common HDR fraction size used in the United States for all stages of cervix cancer has been 6 Gy for 5 fractions, but concerns have been raised about potential toxicity to the sigmoid colon and rectum in patients treated with chemo-radiation.²⁰ As a result, recent clinical trials have included a range of lower fractional doses, such as 5.5 Gy for 5 fractions. Other fractionation regimens are listed in Table 1.

Many institutions use cross-sectional imaging to visualize the cervix and involved regions. In these cases, though the dose to point A should be recorded, the goal should be good coverage (i.e., a D₉₀) of the involved region with EQD2 = 80 Gy for patients with either a complete response or a partial response with residual disease less than 4 cm. For non-responders or those with tumors larger than 4 cm at the time of brachytherapy, tumor dose escalation to an EQD2 of 85–90 Gy is recommended in order to maximize local control.^{21, 22} Other fractionation regimens with EQD2 in the range of 80–85 Gy are acceptable, although the larger the fraction size, the higher the risk for normal-tissue toxicity. For the normal tissues, it is recommended that for each fraction of brachytherapy, the DVH values are calculated and the final dose to the bladder, rectum and sigmoid calculated. Dose limits for the normal tissues are listed in Table 2. The EQD2 limit for the rectum and sigmoid is 70–75 Gy and for the bladder is @90 Gy.²³

Careful consideration should be given to the potential need to boost residual parametrial or lymphnode disease to higher doses. In HDR brachytherapy, the per-fraction dose to the sidewall may be substantial and therefore patients with small tumors or a complete response with no pelvic-sidewall or lymph-node spread of disease do not require a sidewall boost whereas those with enlarged lymph nodes should receive a boost with EBRT.²⁴ With each fraction of brachytherapy, the tumor dose is kept relatively constant, though variations in the normal-tissue doses are to be expected with each fraction. The tumor will likely regress over the course of brachytherapy, and therefore, for point A-specified patients, the OAR doses may increase. If treatment to point A results in normal tissues at or beyond the recommended tolerance doses, consideration should be given to 3D target planning. Another option may be to change to an interstitial implant. In some circumstances, it may be necessary to exceed the usual normal-tissue doses to adequately treat the tumor.

HDR interstitial brachytherapy may be delivered by a variety of alternative fractionation schemes (Table 3). There is a paucity of published experience, and the number of implant procedures and the fractions per implant session are not standardized. The HDR fractionation schedules noted in the literature or used by panelists are presented in Table 3. The dose distribution obtained with the combination of intracavitary and interstitial implants is different from that of an intracavitary implant alone, and may require lower EQD2 doses to the HR-CTV than typically delivered with intracavitary brachytherapy. With all cervical brachytherapy, the central tandem delivers a higher central tumor dose compared to the periphery of the target volume and should be placed when a uterus is present, even when needles are used, to prevent a cold spot.

4.0. Quality management issues for HDR

The large fraction sizes used for HDR brachytherapy require careful monitoring and quality management (QM), given the potential for toxicity and misadministration. Protocol consistency within an institution will help to avoid errors. Institutions should routinely document insertion, planning parameters including normal-tissue dose, treatment, and follow up. A 1998 report from the American Association of Physicists in Medicine addresses QM methods for HDR brachytherapy.²⁵ The recommendations from this report should guide the procedures for any brachytherapy program. QM issues common for all brachytherapy modalities, including treatment planning, treatment delivery systems, applicator commissioning and periodic checks, will not be addressed in this document. Some aspects of quality assurance directed at preventing errors in treatment planning and delivery that are specific to cervical cancer brachytherapy are summarized below.

4.1. Verification of treatment plan—The plan should be verified independently by a qualified brachytherapy physicist not involved in the generation of the plan. This verification should at least include the following items:

1. the dose information matches the prescription
2. the treatment unit, applicator and isotope match the prescription
3. the applicator and dwell positions are correctly located in the patient (consistent with the imaging modalities used)
4. the reference distance from the treatment device to the most distal dwell position is consistent with the applicator in use
5. the individual dwell times and total treatment time are consistent with plans of similar type taking into account the decay of the isotope in use. This can be accomplished by performing an independent calculation to a chosen point in the plan, the use of indices or atlases.

4.2. Pre-treatment verification—Before any treatment is delivered, the pre-treatment information should be verified by a qualified physicist and should include the following items:

1. the correct patient information has been entered into the treatment device
2. the per-fraction dose is consistent with the prescription
3. the dwell times (compensated for isotope decay) and step size programmed into the treatment device are consistent with the treatment plan
4. the channel numbers connected via transfer tubes to the applicator are consistent with the catheter numbers on the plan.

The ABS recommends that radiation oncologists and medical physicists at a facility starting an HDR brachytherapy program for the treatment of patients with cancer of the cervix should attend courses designed to review HDR practice and QM and spend time learning the procedure at a facility with extensive experience in the treatment modality.

Conclusions

The updated ABS 2011 Guideline recommends that 3D imaging with ultrasound, CT or MRI be performed when feasible to estimate the cervical tumor dimensions and ensure adequate coverage of the tumor. Normal-tissue dosimetry using 3D parameters results in a

more accurate reflection of doses administered and may provide more reliable indicators of the risk of toxicity.

References

1. Lanciano RM, Won M, Coia LR, et al. Pretreatment and treatment factors associated with improved outcome in squamous cell carcinoma of the uterine cervix: a final report of the 1973 and 1978 patterns of care studies. *Int J Radiat Oncol Biol Phys.* 1991; 20:667–76. [PubMed: 2004942]
2. Montana GS, Hanlon AL, Brickner TJ, et al. Carcinoma of the cervix: patterns of care studies: review of 1978, 1983, and 1988–1989 surveys. *Int J Radiat Oncol Biol Phys.* 1995; 32:1481–6. [PubMed: 7635793]
3. Stewart AJ, Viswanathan AN. Current controversies in high-dose-rate versus low-dose-rate brachytherapy for cervical cancer. *Cancer.* 2006; 107:908–15. [PubMed: 16874815]
4. Patel FD, Sharma SC, Negi PS, et al. Low dose rate vs. high dose rate brachytherapy in the treatment of carcinoma of the uterine cervix: a clinical trial. *Int J Radiat Oncol Biol Phys.* 1994; 28:335–41. [PubMed: 8276647]
5. Fu KK, Phillips TL. High-dose-rate versus low-dose-rate intracavitary brachytherapy for carcinoma of the cervix. *Int J Radiat Oncol Biol Phys.* 1990; 19:791–6. [PubMed: 2145249]
6. Hareyama M, Sakata K, Oouchi A, et al. High-dose-rate versus low-dose-rate intracavitary therapy for carcinoma of the uterine cervix: a randomized trial. *Cancer.* 2002; 94:117–24. [PubMed: 11815967]
7. Inoue T. The trail of the development of high-dose-rate brachytherapy for cervical cancer in Japan. *Jpn J Clin Oncol.* 2003; 33:327–30. [PubMed: 12949057]
8. Teshima T, Inoue T, Ikeda H, et al. High-dose rate and low-dose rate intracavitary therapy for carcinoma of the uterine cervix. Final results of Osaka University Hospital. *Cancer.* 1993; 72:2409–14. [PubMed: 8402457]
9. Wang X, Liu R, Ma B, et al. High dose rate versus low dose rate intracavity brachytherapy for locally advanced uterine cervix cancer. *Cochrane Database Syst Rev.* 2010:7.
10. Petereit DG, Sarkaria JN, Potter DM, et al. High-dose-rate versus low-dose-rate brachytherapy in the treatment of cervical cancer: analysis of tumor recurrence--the University of Wisconsin experience. *Int J Radiat Oncol Biol Phys.* 1999; 45:1267–74. [PubMed: 10613322]
11. Lertsanguansinchai P, Lertbutsayanukul C, Shotelersuk K, et al. Phase III randomized trial comparing LDR and HDR brachytherapy in treatment of cervical carcinoma. *Int J Radiat Oncol Biol Phys.* 2004; 59:1424–31. [PubMed: 15275728]
12. Eifel P. Patterns of Radiotherapy Practice for Patients Treated for Intact Cervical Cancer in 2005–2007: A QRRO Study. *Int J Radiat Biol Phys.* 2010; 78:S119.
13. Viswanathan AN, Erickson BA. Three-dimensional imaging in gynecologic brachytherapy: a survey of the American Brachytherapy Society. *Int J Radiat Oncol Biol Phys.* 2010; 76:104–9. [PubMed: 19619956]
14. Viswanathan AN, Creutzberg C, Craighead P, et al. International Brachytherapy Practice Patterns: A Survey of the Gynecologic Cancer Intergroup. *Int J Radiat Biol Phys.* Epub 2010 Dec 22.
15. Nag S, Erickson B, Thomadsen B, et al. The American Brachytherapy Society recommendations for high-dose-rate brachytherapy for carcinoma of the cervix. *Int J Radiat Oncol Biol Phys.* 2000; 48:201–11. [PubMed: 10924990]
16. Petereit DG, Sarkaria JN, Chappell R, et al. The adverse effect of treatment prolongation in cervical carcinoma. *Int J Radiat Oncol Biol Phys.* 1995; 32:1301–7. [PubMed: 7635769]
17. Lanciano RM, Pajak TF, Martz K, et al. The influence of treatment time on outcome for squamous cell cancer of the uterine cervix treated with radiation: a Patterns-of-Care study. *Int J Radiat Oncol Biol Phys.* 1993; 25:391–7. [PubMed: 8436516]
18. Viswanathan AN, Moughan J, Small W, et al. The Quality of Cervical Cancer Brachytherapy Implantation and Impact on Local Recurrence and Disease-Free Survival in RTOG Prospective Trials 0116 and 0128. *Int J Radiat Biol Phys.* 2010; 75:S86.

19. Trnkova P, Potter R, Baltas D, et al. New inverse planning technology for image-guided cervical cancer brachytherapy: description and evaluation within a clinical frame. *Radiother Oncol.* 2009; 93:331–40. [PubMed: 19846230]
20. Forrest JL, Ackerman I, Barbera L, et al. Patient outcome study of concurrent chemoradiation, external beam radiotherapy, and high-dose rate brachytherapy in locally advanced carcinoma of the cervix. *Int J Gynecol Cancer.* 2010; 20:1074–8. [PubMed: 20683420]
21. Dimopoulos JC, Lang S, Kirisits C, et al. Dose-volume histogram parameters and local tumor control in magnetic resonance image-guided cervical cancer brachytherapy. *Int J Radiat Oncol Biol Phys.* 2009; 75:56–63. [PubMed: 19289267]
22. Dimopoulos JC, Potter R, Lang S, et al. Dose-effect relationship for local control of cervical cancer by magnetic resonance image-guided brachytherapy. *Radiother Oncol.* 2009; 93:311–5. [PubMed: 19679365]
23. Georg P, Lang S, Dimopoulos JC, et al. Dose-volume histogram parameters and late side effects in magnetic resonance image-guided adaptive cervical cancer brachytherapy. *Int J Radiat Oncol Biol Phys.* Epub 2011 Feb 22.
24. Lee LJ, Sadow CA, Russell A, et al. Correlation of point B and lymph node dose in 3D-planned high-dose-rate cervical cancer brachytherapy. *Int J Radiat Oncol Biol Phys.* 2009; 75:803–9. [PubMed: 19286328]
25. Glasgow GO, Pethel TD. AAPM Radiation Therapy Committee Task Group No 59 K, HD. High dose rate brachytherapy treatment delivery. *Med Phys.* 1998; 25:375–403. [PubMed: 9571605]



\$watermark-text

\$watermark-text

\$watermark-text



Figure 1.

(a) A tandem and ovoid with interstitial catheters (Utrecht applicator, Nucletron B.V., Veenendaal, Netherlands). The applicator uses interstitial catheters that extend above the ovoids and cover a greater width of the cervix higher up than standard ovoids. (b) A tandem and ring applicator with interstitial catheters inserted (Vienna applicator, Nucletron B.V., Veenendaal, Netherlands). Catheters similar to the Utrecht applicator extend the dose to a greater height and width than without. Care must be taken to dose only approximately 10% through the needles, allowing the majority of the dose contribution to be from the tandem and ovoids.

Table 1

Examples of regimens frequently used in the U.S. for tandem and ovoid or tandem and ring brachytherapy.

EBRT, dose to ICRU 52 point or median dose in case of IMRT	Fractionation to point A (Gy)	EQD2 (Gy) to the tumor (point A dose with $\alpha/\beta=10$ Gy)*	EQD2 (Gy) with 90% of the target dose to the OAR using $\alpha/\beta =3$ Gy)	EQD2 (Gy) with 70% of the target dose to the OAR using $\alpha/\beta =3$ Gy)
25 × 1.8 Gy	4 × 7 Gy	83.9	90.1	74.2
25 × 1.8 Gy	5 × 6 Gy	84.3	88.6	73.4
25 × 1.8 Gy	6 × 5 Gy	81.8	83.7	70.5
25 × 1.8 Gy	5 × 5.5 Gy	79.8	82.6	69.6

* For institutions that use radiographic imaging for treatment planning, these doses (e.g. D90) are recorded at point A. For institutions that utilize CT or MRI, these doses are recorded covering the target volume or HR-CTV.

Table 2

Dose limits to the target and to the organs at risk.

	Radiographs	3D imaging
Point A	5×5.25 Gy	variable
D ₉₀		>80 Gy EQD2 or >85 Gy EQD2
ICRU point bladder	$5 \times < 3.7$ Gy	
D _{2cc} bladder		< 90 Gy EQD2
ICRU point rectum	$5 \times < 3.7$ Gy	
D _{2cc} rectum		< 75 Gy EQD2
D _{2cc} sigmoid		< 75 Gy EQD2

Table 3

Suggested doses for template-based HDR interstitial brachytherapy after 45–50.4 Gy of external beam.

Dose of EB radiotherapy	Brachytherapy dose*	EQD2 (Gy) to CTV
45 Gy/25 fractions	3.5 Gy × 9	79.7
	4.25 Gy × 7	79.6
	5 Gy × 5	75.5
50.4 Gy/28 fractions	3 Gy × 9	78.8
	4.5 Gy × 5	76.7

* Twice a day treatments with approximately 6 hours between fractions (based on general radiobiological principles) over one week. The 9-fraction regimen is given over 4.5 days in one week with one insertion. Other regimens using other doses of external beam are also acceptable.