


# “Lab-on-a-chip” nanobiosensor platforms for real-time and early diagnostic of diseases



**Prof. Laura M. Lechuga**

Nanobiosensors and Molecular Nanobiophysics Group  
Research Center on Nanoscience and Nanotechnology **(CIN2: CSIC-ICN)**  
CIBER on Bioengineering, Biomaterials and Nanomedicine **(CIBER-BBN)**  
**Barcelona (Spain)**

## CIN2 (CSIC-ICN)

(Campus UAB,  
Barcelona, Spain)



● Building will be ready in summer 2009. Fully running at the end '09

● 13 Research Groups running, more than 150 people. Centre for 250 persons. 6000 m<sup>2</sup>.

● Access to the largest Clean Room in Spain (CNM-CSIC) with nanotech area.

### MULTIDISCIPLINARY RESEARCH

- Nanoionics, nanooptics, nanomagnetism
- Synthesis of nanoparticles
- Surface nanostructuring
- Nanostructured hybrid materials
- Transport in nanoentities (CNT)
- Nanobiosensors, Nanobiophysics
- Nanotechnology for energy
- Toxicity of Nanotechnology
- Common research with ALBA synchrotron

### TECHNOLOGICAL TRANSFER

SENSIA, SL

ENDOR, SL



More information in  
[www.cin2.es](http://www.cin2.es)




# Outline

- Introduction: Nanomedicine, Nanodiagnostics, Biosensors
- Nanobiosensors for diagnostics
  - Plasmonic sensors
  - Integrated Nanophotonic sensors
  - Nanomechanical sensors
- Lab-on-a-chip microsystems
- Applications: Early detection of cancer
- Future perspectives

# MATERIALS





CNTs, NPs  
Nanowires  
Nanoparticles  
Graphene  
Nanocomposites  
Qdots





## PROSPECTS IN NANOTECHNOLOGY


# ELECTRONICS



# MEDICINE



# ENERGY



Boeing, 2008

Fuel cells  
Solar cells  
Bateries  
Ultracapacitors

Molecular electronics  
Nanochips  
Wereable electronics  
MP3, GPS



# What is Nanomedicine?


Application of materials, devices and process from Nanotechnology to develop nano-sized tools for the **diagnosis, prevention and treatment** of diseases, mainly at the initial stage.

## Nanomedicine

**Nanodiagnostics:**  
Nanosensors, imaging

**Targeted Drug Delivery**


**Regenerative Medicine:**  
Gene Therapy  
Cell Therapy  
Tissue Engineering




**To improve the prevention, diagnosis and therapy in Human Health**

**“Find, fight and follow”**

**The ultimate goal is to improve the quality of life**



# BIOSENSOR DEVICE



## ELECTROCHEMICAL

- amperométricos
- potenciométricos
- conductométricos

- ✓ Sensibles
- ✓ Miniaturización y multiplexado
- ✗ Mediadores para amplificar la señal

## OPTICAL

- absorption
- fluorescence
- **evanescent wave**

- ✓ Altamente sensible
- ✓ Técnica robusta y reproducible
- ✓ Alta capacidad de multiplexado y miniaturización


## MECHANICAL

- piezoelectric
- **nanomechanical**


- ✓ Sensibles
- ✓ Miniaturización y multiplexado
- ✗ Principio de funcionamiento poco establecido

# NANOBIOSENSORS


The key is the bioengineering of the sensor interface




## Nanobiostuctures



Oriented, stable, optimal surface densidad superficial óptima, no adsorciones inespecíficas



# “Lab-on-a-chip” nanobiosensor microsystem



- Nanobiosensores
- Microfluidica
- Electrónica y procesamiento de datos
- fuentes y detectores



## “POCT” point-of care testing

- precision
- fast
- sensitivity
- stability
- selectivity
- no pre-treatment
- label-free



- Instant Diagnostic
- In any place at any time
- Can operate inside the human body


Some applications	Aim
Drug development	Targeted and rapid design of new active ingredients
Cancer research	Identification of cancer genes and individual
Microbiology	Classification of microorganisms
Food industry	Analysis of various disease pathogens in food



# **PLASMONIC NANOBIOSENSORS**



# Plasmonic Biosensor (SPR)




electronics


Portable device, two channels




Wireless transmission module



Software



(old version, new one will be release at the end '08)



Immobilization and regeneration






Flow delivery System



# SPR Biosensor: applications

## Environmental toxic pollutants in real samples (water safety)

### INHIBITION IMMUNOASSAY




- Direct detection in real water samples
- Monoclonal antibodies
- No pre-treatment, fast
- Cycle of analysis 20 min
- More than 200 regenerable
- Validated results
- up to 6 analytes same analysis!!


Anal. Bioanal. Chem. 387, 1449 (2007)  
 Biosens. & Bioelec. 22, 1410 (2007)  
 Anal. Bioanal. Chem. 388, 207 (2007)  
 Sens. Actu. B 118, 399 (2006)  
 Talanta 69 (2), 359 (2006)  
 Anal. Chim. Acta 561, 40 (2006)  
 Biosens. & Bioelec. 21, 2129 (2006)




### Human Growth hormone



### Luteinizing hormone (hLH), Folicle stimulating hormone (FSH)



**OPEN THE WAY TO DIRECT DETECTION OF CANCER  
BIOMARKERS IN HUMAN FLUIDS**







- Direct detection in real human samples
- No pre-treatment
- Detection at physiological levels
- More than 100 regenerable
- Validated results
- Fast
- Low sample volume

# SPR Biosensor: Early cancer diagnostics

## EARLY DETECTION OF INHERITED BREAST CANCER

Detection of point mutations at BRCA-1




185delAG = two bases (AG) deletion

916delTT = two bases (TT) deletion

R1443X = C to T transition


5382insC = C insertion



Covalent attachment of thiol-derivatized DNA probes


### DNA self-assembled monolayers

(Optimisation of buffer concentration, stringency conditions, additives, lateral and vertical spacers,  $T^a$ .)




# SPR Biosensor: Early cancer diagnostics


## Hybridisation detection




## Detection of mutations



## Non stringent conditions




## High stringent conditions




# Hybridization with PCR oligos-like: Tetra-analyte detection format

**185delAG** = two bases (AG) deletion  
**916delTT** = two bases (TT) deletion

**R1443X** = C to T transition  
**5382insC** = one base (C) insertion




**Gene BRCA-1**



**Sensor chip**


Time scale: less than 1 hour



**DNA Extraction**



**PCR**




**To the sensor for analysis**





# BUT with SPR biosensors.....

- No sensitivity for small analytes and very low concentrations (pM-femtoM, single-molecule detection)
- Although portable device, size and weight are not optimal
- Multianalyte capabilities for high throughput




**Higher sensitivity is required**  
**Microsystems platforms**  
**High throughput/multiplexing capabilities**


- 
- 1) Integrated photonic Mach-Zehnder nanointerferometers**
  - 2) Nanomechanical biosensors (standard and optical microcantilevers)**
  - 3) Localised Surface Plasmon Resonance in nanoparticles (LSPR) and magnetoplasmonic biosensors**

# Localised plasmon in Nanoparticles


**Au nanospheres**




EM strongly localized






**Au/SiO<sub>2</sub>**




**Au nanorods**



**Au decahedra**



# LSPR: Dark-field Microscopy





**OPEN THE POSSIBILITY OF HIGHLY DENSE NANOBIOSENSING  
ARRAYS FOR THOUSANDS OF ANALYTES**

# **PHOTONIC NANOBIOSENSORS**



# Photonic Nanosensor: Mach-Zehnder Interferometer




$$I = I_o \cdot \delta^2 \cdot [1 + V \cdot \cos(\Delta\Phi)]$$

$$\Delta\Phi = \frac{2\pi}{\lambda} \cdot \Delta N \cdot L$$


## Design of the system

- *Single mode behaviour*
- *High surface sensitivity*







# Photonic chip

## Optical waveguides




**Standard fabrication at Clean Room facilities.  
Robust and reproducible technology.**




# MZI Performance

## Sensitivity evaluation






$$\Delta n_{o,\min} = 1.1 \times 10^{-7} \quad \Delta N_{\text{eff},\min} = 2.0 \times 10^{-8}$$

The detection limit value corresponds to an average growth layer of  $1.10^{-4} \text{ nm}$  ( $3.10^5 \text{ mol.cm}^2$ ) and a Surface Sensitivity around  $2.10^{-4} \text{ nm}^{-1}$ . Direct detection in the picomolar range ( $10^{-12} \text{ M}$ ) is possible ( $60 \text{ fg/mm}^2$ )






NANOTECNOLOGIA  
CAMPUS UAB, BELLATERRA, BARCELONA




# MZI Sensor: Clinical diagnosis

## EARLY DETECTION OF INHERITED BREAST CANCER

Detection of point mutations at BRCA-1



185delAG = two bases (AG) deletion  
 916delTT = two bases (TT) deletion  
 R1443X = C to T transition  
 5382insC = one base insertion



### BRCA-1 Gene


**DNA probe** (28 mer): 5'-SH-(CH<sub>2</sub>)<sub>6</sub>-(15T)-**GTT CTG TCA AAC** T-3'

**DNA target** (58 mer): 5' - TGC CAC ATG GCT CCA CAT GCA **AGT TTG ACA GAA** CTA CCC TGA TTT TCT GCA C - 3'


**DNA (mutations):** 5' - TGC CAC ATG GCT CCA CAT GCA **AGT TTG AAA CA GA**CTAC CCT GAT ACT TTT CTG GAT GCC -3'

MPTMS  
SILANIZACIÓN


10 % T=25°C t=24 h



SH-ADN




# Nanophotonic chip: Results







# From *chip on a lab* to “*lab on a chip*”




Results using a macro-flow cell:  
only one sensor, high volume of  
reagents and samples.

**Integration of micro/nanobiosensors in  
platforms with microfluidics, biological  
reagents, excitation and signal acquisition  
and processing.**








# “Lab-on-a-chip” Biosensor microsystem: Our approach







Fully integration of optical functions within microfluidic chips is in its infancy  
**on-chip approach**

# Nanosensor/microfluidics integration



Channels of 20-115  $\mu\text{m}$  height  
and 50-150  $\mu\text{m}$  width




3D-microfluidic network

Microfluidic inlet and the  
embedded microchannel on top  
the MZI devices


# Nanosensor/microfluidics integration

## Connection macro/micro



- Modular technology using microfabricated PMMA housing external module
- Easy replacement of the external connections and packaging

- **Steady-state flow rates: 1-1000  $\mu\text{L}/\text{min}$**
- **Up to 10 Bar without any liquid leakage**




- Reduction of 100 times the chemical reagents consumption!!!
- SU-8 microfluidics is resistant to acid and organic solvents (cleaning and regeneration)



# Development of an implantable biosensor for continuous care and monitoring of diabetic patients (P. CEZANNE)




The through-the-skin version



The biosensor will be implanted in the human body during 6 months-1 year, working in a continuous way and with a wireless transmitting data capabilities.

CLALIT HEALTH	ISRAEL
PROTECH AF LTD	UK
LABMAN AUTOMATION	UK
MICROTECH S.R.L.	ITALY
AFCON IND. LTD	ISRAEL
FIMI S.R.L.	ITALY
SIVCO ROMANIA SA	ROMANIA
ROBERT BOSCH GMBH	GE
CNM-CSIC	SPAIN
FRAUNHOFER E.V.	GE






# **NANOMECHANICAL BIOSENSORS**


# What is a microcantilever sensor?

Si microcantilever




A new class of highly sensitive, label-free and direct biosensor which transduces the molecular recognition of biomolecules into a nanomechanical motion (SNPs-DNA and femtomolar-proteins)

## 1) Microcantilever functionalization





## 2) Biomolecular Recognition produces a bending of the cantilever due to change of surface stress



**Cantilever bending ~ nm**


$$z \approx 4 \frac{1-\nu}{Et^2} L^2 (\Delta\sigma_t - \Delta\sigma_b)$$

**OPTICAL READ-OUT:** High sensitivity





Array of cantilevers

# Nanomechanical Biosensors




Si cantilevers, 2560 cantilevers, Array of 20-cantilevers,  $0.334 \text{ mm}$  thick,  $k = 5.5 \cdot 10^{-3} \text{ N/m}$   
Six time more sensitive for biosensing than commercial ones  
minimum detectable deflection  $0.1 \text{ nm}$   
sensitivity expressed in  $\text{fm/Hz}^{1/2}$ : 100



# Integration in Microsystems

For measuring arrays of microcantilevers.....



Array of sensors



## ■ Advantages

- ✓ High sensitivity
- ✓ Sub-angstrom resolution

## ■ Disadvantages

- ✓ Complex integration (difficult alignment): Difficult read-out of arrays
- ✓ Sequential switching of the laser due to overlapping of the reflection beams
  - ✓ Complex detection by the PD chip
  - ✓ Microlens array onto laser chip
- ✓ Immobilisation in each cantilever using the 20-flow cell




Flow cell-20-channels

- Optics: 20-VSCEL + Drivers
- Photodetector array. CMOS circuitry
- Microarrays of 20 microcantilevers
- Microfluidics: flow cell with 20 channels


# Novel optical waveguide microcantilever sensor

- No alignment (except for light coupling)
- More integrated approach
- No further adjustment
- Conventional photodetectors






# OWC: Design and fabrication



## Tested materials:


LPCVD  $\text{Si}_3\text{N}_4$ , PECVD  $\text{Si}_3\text{N}_4$ ,  
PECVD  $\text{SiO}_2$ , Thermal  $\text{SiO}_2$

**Silicon oxide cantilever, 600 nm**  
**Silicon nitride receptor, 120 nm**  
**Wavelength 632.8 nm**  
**Gap: 3  $\mu\text{m}$**   
**2 propagation modes**




**Coupling efficiency at the junction: up to 70 %**


# OWC: Devices



**The chips on wafer**




**A separated chip**



**The input waveguides**





**The cavity and cantilevers**



**The output waveguides**



- 2500 cantilevers per wafer
- 20 cantilevers per chip
- Cantilevers are  $200 \times 40 \times 0.6 \mu\text{m}$
- direct incoupling
- Spring constant  $0.05 \text{ N/m}$

# OWC: Performance




- ✓ OWC chips tested in air
- ✓ operated in dynamic mode

Light intensity distribution



Applying an excitation V of variable frequency and amplitude





The cantilever displacement can be detected with resolution of **0.04 nm**



**HIGH SENSITIVITY, GOOD PERFORMANCES FOR BIOSENSING**


# Future opportunities

## Actual DNA/Protein Biochips



- Need of high amounts of sample: (time consuming, expensive)
- Non-specific interactions
- Indirect read-out: fluorescent labeling. Labels restrict assay types
- Low sensitivity for single mismatch detection
- Inconsistent activity of immobilised proteins

## “Lab-on-a-chip” Nanobiosensors as an advanced platform





### Array of microcantilevers

- Drastic reduction of sample amount
- **Direct read-out:** real time analysis
- High sensitivity
- Micro/nanotechnology: mass production with low cost



# In the future.....

## Nanobiosensores en urgencias



### Tecnologías

Nanofotónica  
Nanotubos,  
nanopartículas  
Nanomecánica, NEMS





### Beneficios

Datos en tiempo real e in-situ  
Imagen a nivel celular  
Herramientas quirúrgicas de precisión  
guiadas por sensores



## Nanobiosensores en la consulta



### Tecnologías

Biochips  
Nanoarrays de alta densidad



### Beneficios

Análisis completo en minutos  
Diagnósticos rápidos y precisos  
Tratamientos específicos y personalizados



## Nanobiosensores en casa



### Tecnologías

Wireless  
Dispositivos portátiles con batería  
Displays de alta resolución



### Beneficios

Auto-Pruebas diagnósticas simples  
Transmisión automática de datos a historial  
clínico





# Acknowledgements

## • Nanobiosensors Group (CIN2-CSIC)

Physicist, chemist, electronics engineers, molecular biologist, biophysics, experts in optoelectronics technology,

## • Financial funding from:

European Union (V and VI FP)

Spanish Ministry of Science

Spanish Ministry of Health (CIBER)

Spanish National Research Council



**M. Botín Foundation**






Barcelona, Spain

# Nanobiosensors Group: our activity

- SPR, Magneto-SPR and LSPR
- Integrated optics Nanophotonic biosensor (MZI)
- $\mu$ -Fluidics integration
- Nanomechanical biosensors (standard and optical)
- Carbon Nanotubes Biosensors
- Biofuncionalization
- Lab-on-a-chip technological platforms (in-vitro and in-vivo)
- Applications: environmental control, clinical diagnostics
- Single-molecule biophysics using Magnetic tweezers and AFM



**MORE INFO IN**  
**[www.cin2.es/biosensores](http://www.cin2.es/biosensores)**



# Sensia

## biosensingSolutions

*Sensia is a company of  Group  
with the participation of MONDRAGON  
COOPERATION (MCC)*

**NEW SPR VERSION AVAILABLE ON JUNE '08**



[www.sensia.es](http://www.sensia.es)