



Title	Array CGH testing in prenatal diagnosis: a promising new service with improved diagnostic yield and shortened reporting time
Author(s)	Kan, ASY; Lau, ET; Tang, WF; Chan, SY; Ding, CK; Lee, CP; Hui, PW; Leung, KY; Ma, T; Leung, WC; Chan, YK; Tang, MHY
Citation	The 2012 Hospital Authority Convention, Hong Kong, 7-8 May 2012.
Issued Date	2012
URL	http://hdl.handle.net/10722/160617
Rights	Creative Commons: Attribution 3.0 Hong Kong License



Service Priorities and Programmes
Electronic Presentations

Convention ID: 904

Submitting author:

Post title: Associate Consultant, Queen Mary Hospital, HKWC

Array CGH testing in Prenatal Diagnosis: A Promising New Service with Improved Diagnostic Yield and Shortened Reporting Time

Kan ASY(1), Lau ET(2)(3), Tang WF (3), Chan SY(2), Ding CK(2), Lee CP(1), Hui PW(1), Leung KY(4), Ma T(4), Leung WC(5), Chan YK(2), Tang MHY(2)
(1)Department of Obstetrics and Gynaecology, Queen Mary Hospital, (2)Prenatal Diagnostic and Counselling Department, Tsan Yuk Hospital, (3)Department of Obstetrics and Gynaecology, The University of Hong Kong, (4)Department of Obstetrics and Gynaecology,

Keywords:

Comparative Genomic Hybridization

Prenatal diagnosis

Introduction

Array Comparative Genomic Hybridization (aCGH) with genome-wide coverage has proved to be valuable for postnatal and prenatal studies. Traditionally, prenatal diagnosis of chromosomal abnormalities has relied on conventional cytogenetics which required cell culture and chromosome analysis under microscope.

Objectives

To evaluate the use of aCGH for prenatal diagnosis in Hong Kong by comparing the (1) concordance of results and (2) reporting time against that of conventional cytogenetics.

Methodology

Two hundred and fifty-two prenatal samples (mostly amniotic fluid or chorionic villi) were analysed by both NimbleGen whole genome oligo array CGX-12 with 135K oligo probes/array as well as conventional karyotyping from January 2011 to December 2011 in Tsan Yuk Hospital Prenatal Diagnostic Laboratory. The proportion of cases detected with clinically significant copy number variations (CNVs) by aCGH were compared to abnormal karyotypes by conventional cytogenetics.

Result

Clinically significant CNVs were detected in 42 (16.7%) out of 252 samples. 12 abnormalities (4.8%) were not detected by conventional karyotyping. Among 106 samples with ultrasound abnormalities, clinically significant CNVs were detected in 36 (34%). 12 abnormalities (11.3%) were not detected by conventional karyotyping. These additionally diagnosed cases include complex unbalanced rearrangements and microdeletions which are clinically important and relevant for counselling on future reproductive risks. Array CGH also enhanced the characterization of the origin and genetic content of supernumerary marker chromosomes. The reporting time has shortened from 14 days for karyotyping to around 5 days for aCGH initial reporting. Conclusion: Array CGH is a promising prenatal diagnostic test with improved

diagnostic yield and shortened reporting time by 64%.