

CORE

Title	Nutrient supplemented serum-free medium increases cardiomyogenesis efficiency of human pluripotent stem cells.
Author(s)	Ting, S; Lecina, M; Chan, YC; Tse, HF; Reuveny, S; Oh, SKW
Citation	World Journal of Stem Cells, 2013, v. 5 n. 3, p. 86-97
Issued Date	2013
URL	http://hdl.handle.net/10722/203106
Rights	Creative Commons: Attribution 3.0 Hong Kong License



Online Submissions: http://www.wjgnet.com/esps/ wjsc@wjgnet.com doi:10.4252/wjsc.v5.i3.86 World J Stem Cells 2013 July 26; 5(3): 86-97 ISSN 1948-0210 (online) © 2013 Baishideng. All rights reserved.

ORIGINAL ARTICLE

Nutrient supplemented serum-free medium increases cardiomyogenesis efficiency of human pluripotent stem cells

Sherwin Ting, Marti Lecina, Yau-Chi Chan, Hung Fat Tse, Shaul Reuveny, Steve KW Oh

Sherwin Ting, Marti Lecina, Shaul Reuveny, Steve KW Oh, Bioprocessing Technology Institute, A*STAR (Agency for Science, Technology and Research), Singapore 138668, Singapore Yau-Chi Chan, Hung Fat Tse, Cardiology Division, Department

of Medicine, The University of Hong Kong, Queen Mary Hospital, Hong Kong, China Author contributions: Ting S and Lecina M contributed to

Author contributions: Ting S and Lecina M contributed to conception and design; Ting S, Lecina M, and Chan YC contributed to collection and assembly of data; Ting S, Reuveny S and Reuveny S contributed to manuscript writing; all authors contributed to data analysis and interpretation, and final approval of manuscript.

Supported by Science Technology and Research (A*STAR); Hong Hong Kong Research Grant Council Collaborative Research Fund (HKU8/CRF/09) and Theme-based Research Scheme (T12-705/11); Tse HT and Oh SKW contributed to financial support

Correspondence to: Dr. Steve KW Oh, Principal Scientist, Bioprocessing Technology Institute, A*STAR (Agency for Science, Technology and Research), 20 Biopolis Way No. 06-01, Singapore 138668, Singapore. steve_oh@bti.a-star.edu.sg Telephone: +65-6407-0855 Fax: +65-6478-9561 Received: March 26, 2013 Revised: June 3, 2013

Accepted: June 8, 2013

Published online: July 26, 2013

Abstract

AIM: To develop an improved p38 MAPK inhibitor-based serum-free medium for embryoid body cardiomyocyte differentiation of human pluripotent stem cells.

METHODS: Human embryonic stem cells (hESC) differentiated to cardiomyocytes (CM) using a p38 MAPK inhibitor (SB203580) based serum-free medium (SB media). Nutrient supplements known to increase cell viability were added to SB medium. The ability of these supplements to improve cardiomyogenesis was evaluated by measurements of cell viability, total cell count, and the expression of cardiac markers *via* flow cytometry. An improved medium containing Soy hydrolysate (HySoy) and bovine serum albumin (BSA) (SupSB media) was developed and tested on 2 additional cell lines (H1 and Siu-hiPSC). Characterization of the cardiomyocytes was done by immunohistochemistry, electrophysiology and quantitative real-time reverse transcriptionpolymerase chain reaction.

RESULTS: hESC cell line, HES-3, differentiating in SB medium for 16 d resulted in a cardiomyocyte yield of 0.07 ± 0.03 CM/hESC. A new medium (SupSB media) was developed with the addition of HySoy and BSA to SB medium. This medium resulted in 2.6 fold increase in cardiomyocyte yield (0.21 ± 0.08 CM/hESC). The robustness of SupSB medium was further demonstrated using two additional pluripotent cell lines (H1, hESC and Siu1, hiPSC), showing a 15 and 9 fold increase in cardiomyocyte yield respectively. The age (passage number) of the pluripotent cells did not affect the cardiomyocyte yields. Embryoid body (EB) cardiomyocytes formed in SupSB medium expressed canonical cardiac markers (sarcomeric α -actinin, myosin heavy chain and troponin-T) and demonstrated all three major phenotypes: nodal-, atrial- and ventricular-like. Electrophysiological characteristics (maximum diastolic potentials and action potential durations) of cardiomyocytes derived from SB and SupSB media were similar.

CONCLUSION: The nutrient supplementation (HySoy and BSA) leads to increase in cell viability, cell yield and cardiac marker expression during cardiomyocyte differentiation, translating to an overall increase in cardiomyocyte yield.

© 2013 Baishideng. All rights reserved.

Key words: Soy hydrolysate; Bovine serum albumin; Differentiation; Cardiomyocyte; Human embryonic stem cells; Human induced pluripotent stem cells

Core tip: Nutrient supplements were screened for improving cell survival during the cardiomyocyte differentiation process of human embryonic stem cells (hESC) and human induced pluripotent stem cells (hiPSC) in



a serum-free medium based on the inhibition of p38 MAPK (SB media). Soy hydrolysate and bovine serum albumin supplementation was found to improve cell viability and therefore increased the yield of HES-3 cardiomyocytes by 2.6-fold over non-supplemented SB medium. The enhancing effect of this medium was demonstrated in an additional hESC line (H1) and Siu1-hiPSC cell line (15 and 9 fold respectively). The cardiomyocytes formed expressed canonical cardiac markers (sarcomeric α -actinin, myosin heavy chain and troponin-T) and demonstrated all three major phenotypes: nodal-, atrial- and ventricular-like.

Ting S, Lecina M, Chan YC, Tse HF, Reuveny S, Oh SKW. Nutrient supplemented serum-free medium increases cardiomyogenesis efficiency of human pluripotent stem cells. *World J Stem Cells* 2013; 5(3): 86-97 Available from: URL: http://www. wjgnet.com/1948-0210/full/v5/i3/86.htm DOI: http://dx.doi. org/10.4252/wjsc.v5.i3.86

INTRODUCTION

Heart disease is one of the most common causes of mortality in the world and accounts for more than 800000 deaths per year on average in the United States alone^[1]. After an episode of major cardiac insult (such as a myocardial infarction), the heart can lose up to 2 billion cardiomyocytes. Being an organ that cannot auto-regenerate, progressive heart failure develops. Cardiomyocyte cell therapy can thus be a potential cure for heart failure after myocardial infarction, and it is suggested that at least one billion cardiomyocytes will be required per patient for such treatment^[2].

Human embryonic stem cells (hESCs) present an attractive cell source for generating large amounts of cardiomyocytes, due to their pluripotency and ability to proliferate for multiple passages^[3]. Several studies have reported techniques to efficiently differentiate hESCs to cardiomyocytes via growth factors and small molecule inhibitors^[4-8]. Current cardiomyocyte differentiation protocols can be divided into two groups: differentiating hESC in 2D monolayers culture or in 3D suspended embryoid bodies (EBs) cultures^[9]. Although monolayer differentiation protocols have achieved high yields of cardiomyocytes [6,10], the scalability of these methods is problematic and they have limited capability in generating the amounts of cardiomyocytes needed for cell therapy. On the other hand, methods that involve EBs formation, which have better potential for scale up, have lower yields of cardiomyocyte, and requires extensive use of expensive growth factors like BMP4 and activin A at multiple specific time points during differentiation^[7,11-13]. In addition, the growth factors have to be optimized for different cell lines, growth platforms or passage numbers^[14,15]. As such, there is a lack of protocols for cardiomyocyte differentiation in EB suspended cultures that are cost-effective, scalable and most importantly robust. Previously, we have developed a simple scalable methodology to differentiate hESCs to cardiomyocytes using a serum-free differentiation medium^[16,17] containing a small molecule p38 MAP kinase inhibitor SB203580 (SB media)^[18,19]. The enhancing effect of SB203580 on cardiomyogenesis of hESC has been correlated to the expected inhibition of the p38 pathway as well as the activation of JNK^[20]. This suggests a regulatory interlink between the JNK and p38 pathways during cardiomyogenesis. Compared to protocols based on growth factors, small molecules are less costly and more amenable for good manufacturing practice (GMP) manufacturing of cells^[21]. However, the SB medium is essentially protein-free and lacks nutrients (e.g., lipids) and growth factors. From the low cell viability and yield observed, we hypothesized that SB medium has nutritional deficiencies that limit cardiomyogensis, especially in the initial stages of the differentiation process.

In this study, we sought to improve the survival of cells in SB medium and thereby enhance cardiomyogenesis using the embryoid body method of differentiation. For successful growth and maintenance of metabolic functions of differentiated human cells *in vitro*, appropriate culture conditions are required to mimic the physiological conditions *in vivo*^[22,23]. The culture medium is one of the most important factors in maintaining cell and tissue culture as it provides nutrients and salts, hormones and growth factors, buffering elements and oxygen supply^[22-24]. While media supplements have been developed for a variety of cell types, none have been performed for stem cell differentiation.

In this study, nutritional components were screened for an improvement in overall yields of cardiomyocytes of HES-3 cells using serum-free and insulin-free SB medium. Two supplements, Soy-hydrolysate (HySoy) and bovine serum albumin (BSA), resulted in improved cardiomyocyte differentiation efficiency. The concentrations of both supplements were optimized, resulting in an increase in cell growth and differentiation. The robustness of this new medium (SupSB media) was evaluated using the H1 and HES-3 hESC, as well as the Siu1-hiPSC cell lines (2.6 to 15.0 fold increase in cardiomyocyte yield). Similar cell yields and differentiation efficiency were obtained using 9 different batches of BSA and HySoy indicating that batch variability is not a major concern. Cardiomyocytes formed in SupSB media expressed canonical cardiac markers and electrophysiological studies demonstrated successful cardiac differentiation to give all the three major phenotypes of cardiomyocytes: nodal-, atrial- and ventricular-like. In summary, we have developed a cost-effective, scalable and robust protocol for cardiomyocyte differentiation by improving the p38 MAP kinase protocol with the addition of BSA and HySoy. SupSB media increases cell yield and cardiac expression markers during cardiomyocyte differentiation, translating to an overall increase in cardiomyocyte yield.

MATERIALS AND METHODS

Culture of hESCs and hiPSCs

HES-3 [(46, XX); ES Cell International], H1 [(46, XY);



Ting S et al. Cardiomyogenesis of human pluripotent stem cells




Figure 1 Dot plot of anti-myosin heavy chain and anti-sarcomeric α -actinin of HES-3 differentiated cardiomyocytes. A: Myosin heavy chain (MHC); B: Sarcomeric α -actinin (SA).

WiCell], and Siu1-hiPSC (Professor Tse HF, The University of Hong Kong) with normal karyotypes were cultured in KNOCKOUT medium on inactivated immortalized mouse feeders. The medium was refreshed daily and the cells were passaged weekly. Cultures were kept at 37 °C with 5% CO₂.

Differentiation to cardiomyocytes

hESC cultures were washed using phosphate buffered saline (PBS) (Invitrogen), cut into small clumps (EZpassage tool; Invitrogen), and seeded at 1.33 \times 10⁶ cells/ mL in ultra-low attachment 12-well plates (Nunc). The plates were agitated for 1 h and then cultured in static conditions at 37°C in a humidified atmosphere with 5% CO2. The differentiation medium (SB media) comprised Dulbecco's modified Eagle's medium (Invitrogen) supplemented with 2 mmol/L L-glutamine (Invitrogen), 0.182 mmol/L sodium pyruvate (Invitrogen), 1% nonessential amino acids (Invitrogen), 0.1 mmol/L 2-mercaptoethanol, 5.6 mg/L transferrin (Invitrogen), and 20 mg/L sodium selenite (Sigma). A solution of 5 mmol/L of p38-MAPK inhibitor (SB203580; Sigma), dissolved in dimethylsulfoxide (Sigma), was added to the medium at a final concentration of 5 µmol/L. The medium was refreshed every 2 d following the previously described protocol^[16].

Supplements for cardiomyocyte differentiation

Various supplements (Table 1) were added to the SB medium and tested for cell yields and cardiomyogenesis efficiency using the above differentiation protocol.

Cardiomyocytes harvesting and quantification

Following the protocol established by Lecina *et al*^[15] and Ting *et al*^[16], beating aggregates collected from 1 well were washed with 5 mL of PBS without Ca²/ Mg² (PBS-) (Invitrogen), and incubated for 30 mins in 500 µL PBSsolution containing 1.6 mg/mL collagenase (Sigma) and 20% fetal bovine serum (Hyclone). Thereafter, 500 µL of 0.25% trypsin EDTA (Invitrogen) was added to generate single-cell suspensions. Cell suspensions were finally filtered through a 40-µm cell strainer (Becton Dickinson), fixed, and permeabilized (Caltag Laboratories). For FACS analysis, cells were incubated with anti-myosin heavy chain (anti-MHC; MF20, dilution 1:200; Developmental Studies Hybridoma Bank), anti-sarcomeric α -actinin (anti-SA; 1:100; Sigma) and anti-Troponin T (anti-cTnT; 1:200; Thermo Scientific). Fluorescein isothiocyanateconjugated antimouse antibody (1:500; DAKO) was used as a secondary antibody. Dot plots for both anti-MHC and anti-SA are shown in Figure 1. All incubations were carried out at room temperature for 30 min. FACS measurements were done using Guava (Millipore).

Cell counting and calculations of cardiomyocyte/hESC yield and normalized yield

Cell concentration was determined by the nuclei count method using NucleoCounter (Chemometec). Apoptosis and cell proliferation were evaluated using Annexin V (Invitrogen) and Ki-67 (BD Biosciences) respectively via FACS. Cardiomyocyte yields were calculated following published protocols^[16,17], namely: CM hESC iield = (cell countDay 16/cell countDay 0: Initial hESC seeded) × (% positive cellsMHC, SA/100); normalised yield = (CM/hESC yieldcondition)/(CM/hESC yieldcontrol).

qRT PCR

Total RNA was isolated from the cells ($< 5 \times 10^{\circ}$) using RNeasy Mini Kit (Qiagen) following the supplier's protocol. Reverse transcription was carried out with 1 µg total RNA using SuperScript III (Invitrogen). Real-time PCR was performed by applying a standard two-step amplification protocol on an ABI 7500 system (Applied Biosystem) to detect mRNA expression. Normalization of the results was done using a house-keeping gene, GAPDH. Primer sequences for Nanog, OCT 4, T-bra, Mesp 1, Nkx 2.5, MHC and GAPDH are provided in Table 2.

Immunocytochemistry

Cell aggregates were harvested, washed with PBS-, mechanically dissociated by pipetting, and plated in 24-well plates for 2 d in the SupSB medium at 37°C in a humidified atmosphere with 5% CO₂. The cells were then washed twice with PBS- and fixed with 4% paraformaldehyde (2 mL for 20 minutes at room temperature). After washing twice with PBS-, permeabilization and blocking was done using 0.1% Triton X-100 and 10% goat serum



Table 1 List of nutritional supplements								
Supplements	Concentration/addition	Source						
MEM vitamin solution (× 100)	30 µL/mL	Gibco 11120						
Yeastolate ultrafiltrate (× 50)	20 µL/mL	Gibco 18200						
Vitamin E	70 μg/mL	Sigma T3376						
Synthechol (cholesterol) (× 500)	12 μL/mL	Sigma S5442						
Yeast extract	50 μg/mL	Sigma Y1625						
Bovine serum albumin	50 µg/mL	Gibco A10008-01						
Soy hydrolysate	50 μg/mL	Kerry Bio-science						
		5X59022						

respectively. The following antibodies were used: anti-MHC (MF20; Developmental Studies Hybridoma Bank), anti-SA (Sigma), and anti-troponin-T (Thermo Scientific). Nuclear staining was done using SlowFade Glow with DAPI (4'6-diamidino-2-phenylindole) (Invitrogen). The fluorescence was observed using an Olympus IX71 fluorescence microscope (Olympus) coupled with Olympus imaging software Cell P.

Electrophysiology

Standard whole-cell patch-clamp recordings were performed at 37 \pm 0.5 °C to record the action potential phenotypes (HEKA Instruments Inc. Southboro, MA, United States) of beating cardiomyocyte aggregates as previously described^[25,26] . Patch pipettes were prepared from 1.5-mm thin-walled borosilicate glass tubes using a Sutter micropipette puller P-97 and had typical resistances of 3-4 M Ω when filled with an internal solution containing (mmol/L): 110 K⁺ aspartate, 20 KCl, 1 MgCl₂, 0.1 Na-GTP, 5 Mg-ATP, 5 Na2-phosphocreatine, 5 EGTA, 10 HEPES, and pH adjusted to 7.3 with KOH. The external Tyrode's bath solution consisted of (mmol/L): 140 NaCl, 5 KCl, 1 MgCl2, 0.4 KH2PO4, 1.8 CaCl2, 10 Glucose, 5 HEPES, with pH adjusted to 7.4 with NaOH. Twenty consecutive action potentials from spontaneously firing HES3-derived cardiomyocytes were recorded per cell to ensure stable waveforms for analysis. For the electrically quiescent cardiomyocytes, a stimulation of 0.1-1 nA for 5 ms was given to elicit an action potential. The sampling frequency was 2.00 kHz and data were corrected for the liquid junction potentials of +15.9 mV. Maximal diastolic potential as well as action potential duration at 90% (APD90) and 50% repolarization level (APD50) were measured.

Statistical analysis

Experiments were done using three independent replicates. The significance of the results was calculated by Student's *t*-test or by one-way ANOVA (P < 0.05, P < 0.01).

RESULTS

Identifying nutritional supplements that can improve cardiomyocyte differentiation

In order to identify elements that can increase cardiomyocyte differentiation, we selected a panel of defined and




Figure 2 Effect of nutritional supplements added to SB medium on HES-3 differentiation to cardiomyocytes. Differentiation efficiency is based on two FACS markers: myosin heavy chain (MHC) and sarcomeric α -actinin (SA). Cell and cardiomyocyte (CM) yields were measured at day 16 of culture. Results were normalized against the yields obtained in SB medium (control). (n = 3, ${}^{a}P < 0.05$, ${}^{b}P < 0.01$ vs cell yield). BSA: Bovine serum albumin; HySoy: Soy-hydrolysate.

non-defined nutritional supplements (Table 1) which are known to support cell growth of a variety of cell lines^[27]. The supplements were used at concentrations typically reported for animal cell culture.

HES-3 cells were seeded at a concentration of 1.33 $\times 10^{\circ}$ cells/mL and differentiated in SB medium supplemented with the various supplements described in Table 1. Control cultures were differentiated in non-supplemented SB medium. After 16 d, cells were harvested and measured for cell yield and expression of the cardiac markers sarcomeric α -actinin (SA) and myosin heavy chain (MHC). A third antibody, cTnT which detects cardiac troponin was also used to verify the results (data not shown). The differentiation efficiency of the supplement was evaluated by dividing the yields of cardiomyocytes produced per seeded hESC in the supplemented culture with the one in the control culture (normalized cardiomyocyte yield). This parameter considers both the differentiation efficiency (percentage of cardiomyocytes) and final number of total cells, making it meaningful for process analysis and evaluation^[16,17].

The nutritional supplements did not significantly affect cell counts. However, supplementing the media with BSA and HySoy improved the metabolic state of the cells, leading to an increase in normalized cardiomyocyte yield with improvements of 1.47 ± 0.20 and 2.45 ± 0.33 , fold on average respectively, compared to the control (Figure 2). Addition of vitamins had no effect on the differentiation process. Yeastolate, Yeast extract, and Vitamin E reduced cardiomyocyte yields showing a negative effect on cardiac marker expression (data not shown), while cholesterol had a negative effect on cell yields.


In order to further increase cardiomyocyte yields, the effect of HySoy and BSA concentration on cardiomyocyte yield was further investigated (Figure 3). At a BSA




Ting S et al. Cardiomyogenesis of human pluripotent stem cells

Table 2 Quantitative real-time reverse transcription-polymerase chain reaction primers								
Primers for qPCR	Forward (5'→3')	Reverse (5'→3')						
Nanog	GAAAAACAACTGGCCGAAGAAT	GGTGCTGAGGCCTTCTGC						
4-Oct	AACGACCATCTGCCGCTTT	GGCCGCAGCTTACACATGTT						
T-bra	AATTTGGTCCAGCCTTGGAAT	CGTTGCTCACAGACCACAG						
Mesp 1	GACGTGCTGGCTCTGTTG	TGTCACTTGGGCTCCTCAG						
Nkx 2.5	CAAGTGTGCGTCTGCCTTT	TTGTCCGCCTCTGTCTTCTC						
MHC	ATTGCTGAAACCGAGAATGG	CGCTCCTTGAGGTTGAAAAG						
GAPDH	GTCGGAGTCAACGGATTTGG	AAAAGCAGCCCTGGTGACC						

qPCR: Quantitative polymerase chain reaction; MHC: Myosin heavy chain; GAPDH: Glyceraldehyde 3-phosphate dehydrogenase.








Figure 3 Additive and dose effect of Soy hydrolysate and bovine serum albumin supplementation on HES-3 differentiation to cardiomyocytes. HES-3 cells differentiated in SB medium containing varying concentrations of bovine serum albumin (BSA) and Soy hydrolysate (HySoy), as well as a combination of both supplements (SupSB media), were harvested on day 16 and evaluated for total cell count (A) and cardiac specific markers, myosin heavy chain (MHC) and sarcomeric α -actinin (SA) (B); C: The results were summarized by calculating the normalized yield which takes into account the ratio of cardiomyocyte at day 16 compared to the initial human embryonic stem cell seeded. Normalized yields are based on both MHC and SA markers. Dot plots of both MHC and SA are shown in Figure 1. (n = 4. ^bP < 0.01 vs SB medium; ^cP < 0.05; ^dP < 0.01 vs SupSB media).

concentration of 0.5% v/v, an increase in cell growth compared to the control was observed, achieving an average of 2.54×10^6 cells/mL (Figure 3A). Similarly, cultures supplemented with 0.25% w/v HySoy achieved a maximum cell density of 2.51×10^6 cells/mL. The

supplements also improved cell viability, showing > 85% viability as compared to 65% in the controls (SB media). This increase in cell growth and viability resulted in a 2-fold increase in normalized cardiomyocyte yield for both BSA and HySoy (Figure 3C). Although there was no

WJSC | www.wjgnet.com




Figure 4 Cell growth and apoptosis in HES-3 cells during cardiomyocyte differentiation. A: Cell growth kinetics of SB and SupSB media cultures; B: Kinetics of apoptosis in SupSB media cultures (measured by annexin V antibody) (n = 3, ${}^{b}P < 0.01$ vs in ten days).

statistical difference in the level of expression of cardiac markers (SA and MHC) between the cultures supplemented with different concentrations of BSA and HySoy as compared to the control, a general trend of increased level of expression was observed in cultures supplemented with 0.25% w/v HySoy and 0.5% v/v BSA (Figure 3B). Embryoid bodies of all cultures started beating at day 9, indicating that the supplements did not affect the process of cardiac differentiation temporally.

Additive effect of HySoy and BSA supplements further increases cardiomyocyte yields

Next, we explored the effects of combining both supplements on cell yield and cardiac marker expression (SA and MHC). HES-3 cells were seeded in medium containing both 0.25% w/v HySoy and 0.5% v/v BSA (SupSB media) and compared to control cultures supplemented by HySoy and BSA alone or not supplemented at all.

Results show a further increase in normalized cardiomyocyte yield (2.6-fold) when HySoy and BSA were combined, compared to only a 2-fold increase with the individual supplements alone (P < 0.01) (Figure 3C). This increase in cardiomyocyte yield can be attributed to a higher percentage of differentiated cardiomyocytes as indicated by a higher percentage of MHC-expressing cells (16.5% of cells) cultured in the SupSB medium vs. the individual supplements (10%-12%) (Figure 3B). In addition, there was also an increase in cell yield of cultures in SupSB media resulting in $2.79 \times 10^{\circ}$ cells/mL (Figures 3A) and Figure 4), that was statistically significant (P < 0.01). The variability between different batches of BSA and Hy-Soy were also tested. Differentiation experiments using 9 different batches of BSA and HySoy were conducted and little variability was observed. A high level of apoptotic cells (60% of total cells), measured by annexin V, was observed during the first day of differentiation, which correlates to cell death probably as a result of mechanical manipulation of the culture during seeding. Thereafter, gradual increase in cell yields and down regulation of

apoptotic markers were observed (Figure 4). In summary, the additive effect of HySoy and BSA increased cardiomyocyte yield further over individual supplements alone, through increasing cell yield and expression of cardiac markers.

Kinetics of cell growth and marker expression during HES-3 differentiation

Next, we investigated the underlying factors that resulted in the overall increase in cardiomyocyte yields using the optimal concentrations of both supplements. Differentiating HES-3 cultures supplemented with optimal concentrations of HySoy and BSA were monitored for 10-16 d. Every 2-3 d cells were harvested and analyzed for cell yield, expression of pluripotent markers (Nanog and OCT4), mesoderm markers (Mesp 1 and T-bra), the early cardiac marker Nkx 2.5 and the late cardiac marker MHC *via* qRT-PCR (Figure 5).

More than 30% of cell death was observed during the first day in both the supplemented as well as the control cultures. This phenomenon can be attributed to the mechanical cutting of the HES-3 cultures and the adaptation of the cells to the new differentiation media. Thereafter, BSA and HySoy supplemented cultures showed 34% and 39% increase in cell yields respectively by day 3 and achieved cell densities of 2.16×10^6 and 2.32×10^6 cells/mL on day 10, as compared to 1.65×10^6 and 1.44×10^6 cells/mL in the control cultures respectively (Figure 5A and B). At the end of the culture, addition of BSA and HySoy resulted in a rise in cell yield by 23% and 19% respectively as compared to the controls.

qRT-PCR analysis show that both supplemented cultures and SB control culture displayed similar trends of down regulation of pluripotent markers (Nanog and OCT4) as well as up-regulation of mesoderm markers (Mesp 1 and T-bra) (Figure 5C and D). However, cultures with supplements showed a significantly higher expression of cardiac progenitor gene Nkx 2.5 as well as a higher general trend of late cardiomyocyte gene MHC






Figure 5 Cell growth and gene expression kinetics of HES-3 cells grown in cultures with and without supplements. Cell densities at multiple times points over a period of 10 to 16 d were recorded with cultures with and without supplements bovine serum albumin (BSA) (A); Soy hydrolysate (Hysoy) (B). Gene expression profile recorded via quantitative real-time reverse transcription-polymerase chain reaction at multiple time points over 16 d during the course of cardiomyocytes (CM) differentiation. 6 genes: hESC pluripotency (OCT4, Nanog), Mesoderm (Mesp1, T-bra), cardiac progenitors (Nkx2.5), and cardiomyocytes (MHC) were monitored for BSA (C); Hysoy (D). n = 3, ${}^{\circ}P < 0.05$, ${}^{b}P < 0.01$ vs SB medium.

expression during the later stages as compared to SB control cultures (Figure 5C and D). In summary, the addition of supplements, BSA and HySoy, not only enhanced cardiomyogenesis *via* improvement in cell yield, but also increased the expressions of both Nkx 2.5 and MHC.

Robustness of SupSB medium: Increased cardiomyogenesis with other pluripotent cell lines

In order to demonstrate the universality and robustness of SupSB media for increased cardiomyogenesis induction efficiency, we cultured two additional cell lines, H1 hESC and Siu1-hiPSC. These cells were differentiated in SupSB medium for 16 d and tested for cell yields and cardiac marker expression, with cells differentiated in SB medium used as controls. Results show that culturing cells in SupSB medium led to increased cell density from 1.10×10^6 to 3.12×10^6 cells/mL for H1 and from 0.67×10^6 to 1.22×10^6 cells/mL for Siu1-hiPSC as compared to SB media (Table 3). Moreover, significantly higher expression of the cardiac specific marker MHC was observed in cultures with SupSB medium (Table 3). H1 cultures differentiated in SupSB medium showed a significant increase in normalized cardiomyocyte yield of 15-fold over the control, while Siu1-hiPSC showed a 9-fold increase (Table 3). Cardiomyocyte yield per hESC seeded was 0.59 cardiomyocyte/hESC for H1 cultures and 0.04 cardiomyocyte/hESC for Siu1-hiPSC cultures.

Another factor that can influence the yield of cardiomyocyte differentiation is the age or passage number of the cells. As such, HES-3 cells at early (P12 to P13) and late (P21 to P39) passages were differentiated using SupSB medium and cardiomyocyte yields were measured at day 16. Results in Figure 6 show that cells from early or late passage showed similar cardiomyocyte yields. Table 3 Summary of results: Robustness of SupSB vs SB media for cell lines HES-3, H1 and Siu1-human induced pluripotent stem cells

Conditions		Cell expa	nsion	Differentiation	efficiency	CM yield
Cell type	Medium	Cell yield (10 ⁶ cells/mL)	Cell expansion fold	Sarcomeric actinin	мнс	Cardiomyocyte/hESC
HES-3	SupSB	2.27 ± 0.58	1.65 ± 0.47	14.70% ± 6.30%	13.59% ± 3.40%	0.21 ± 0.08
	SB	1.66 ± 0.48	1.25 ± 0.36	9.25% ± 6.19%	6.17% ± 2.95%	0.07 ± 0.03
H1	SupSB	3.14 ± 0.24	2.35 ± 0.18	26.70% ± 2.09%	$25.20\% \pm 5.41\%$	0.59 ± 0.01
	SB	1.10 ± 0.19	0.82 ± 0.14	$20.23\% \pm 1.87\%$	$4.76\% \pm 1.62\%$	0.04 ± 0.01
Siu1-hiPSC	SupSB	1.22 ± 0.15	0.92 ± 0.11	$20.76\% \pm 4.74\%$	$2.11\% \pm 0.73\%$	0.04 ± 0.04
	SB	0.67 ± 0.26	0.51 ± 0.20	$16.90\% \pm 1.95\%$	$1.01\% \pm 0.20\%$	0.004 ± 0.001

MHC: Myosin heavy chain. Cultures with SupSB medium were compared against SB medium for cell expansion and differentiation efficiency at the end of culture (day 16). The results were summarized by calculating the cardiomyocyte yield [cardiomyocytes (CM) yield] which takes into account the ratio of cardiomyocytes at day 16 compared to the initial human embryonic stem cells (hESC) seeded (n = 3 to 9).




Figure 6 Robustness of SupSB medium over different cell ages (passage). HES-3 cells from different passages (P12 to P39) were differentiated using the supplements, Soy hydrolysate (HySoy) and bovine serum albumin (BSA). The passages were divided into two groups: Early passage (P12 to P13) and Late Passage (P21 to P39). There is no significant difference in cardiomyocyte yield at different passages. Student *t*-test and one-way ANOVA test were used to calculate the statistical significance of the data.

Characterization of cardiomyocytes obtained from SupSB differentiation medium

HES-3, H1 and Siu1-hiPSC cells differentiated in SupSB medium all showed increases in cell yield, expression of cardiac markers and cardiomyocyte yields (Table 3). The increase in cell yield can be attributed to better metabolic conditions manifested by an increase of 20%-30% in cell viability over the control and an increased expression of proliferation marker Ki-67^[28]. Specifically, at day 2 of differentiation, 12% increase in Ki-67 expression was observed in cultures differentiated in SupSB medium compared to SB medium. By day 4 onward, Ki-67 expression was similar in both cultures (data not shown). The percentage of beating aggregates in cultures with SupSB medium was also higher than that in SB medium (Figure 7).

Cardiomyocytes differentiated in SupSB medium were also characterized via immunostaining of cardiac proteins. 16 d old EBs were harvested, mechanically dissociated, and plated on 0.1% gelatin-coated 24-well plates before staining with a set of antibodies against cardiacspecific markers (Figure 8A-C). Cells stained positive for cytoskeleton structural proteins (SA) and contractile functional proteins responsible for motility (MHC, and troponin-T)^[29].Electrophysiological study *via* whole-cell patch-clamping was performed. Action potential durations at 90% and 50% repolarization levels as well as maximal diastolic potential were recorded using HES-3 cells differentiated for 23 d. Cells differentiated in both SB and SupSB medium gave similar results. All the three major phenotypes of cardiomyocytes: nodal-, atrial- and ventricular-like were obtained (Figure 8D). The APD and maximal diastolic potential of the cardiomyocytes were similar for cells cultured in both SB and SupSB media (Figure 8E-G) and are comparable to those of the control hESC cell line H7 (data not shown).

Furthermore, a normal karyotype was observed in 20 d old cardiomyocytes differentiated in SupSB medium (Figure 9).

DISCUSSION

In this paper, we have demonstrated that the media supplements, HySoy and BSA, can improve the differentiation efficiency of hESCs into cardiomyocytes when added to SB media containing the p38 MAP kinase inhibitor, SB203580. In particular, when combined at optimal concentrations (0.25%HySoy and 0.50%BSA), the cardiomyocyte yield is enhanced by 2.6-fold compared to the SB medium without supplements for the HES-3 cell line. This increase is due to the improvement in cell growth, as well as differentiation efficiency. Cell lines such as, H1 and Siu1-hiPSC, which initially did not differentiate efficiently to cardiomyocytes in the SB media, showed significant improvements in yields of cardiomyocyte/hESC or cardiomyocyte/hiPSC (15- and 9-fold respectively) when these supplements were added. The variability in cardiomyocyte differentiation yields between cell lines is not surprising and it likely reflects genetic and epigenetic differences between pluripotent stem cell lines^[30-32] that influences their cardiac differentiation efficiencies^[14,15,33] Moreover, the level of differentiation obtained in SupSB medium was consistent and not affected by cell age (passages 12 to 39). We have observed varying ratios of expression levels for MHC and SA, ranging from 1:1 to 1:16 (MHC: SA) (Table 3). These differences can be attributed to the stages in the differentiation process which the cells

Ting S et al. Cardiomyogenesis of human pluripotent stem cells




Figure 7 Phase contrast microscopy of cardiomyocytes obtained from differentiation of HES-3 cells. Pictures of aggregates formed from human embryonic stem cells differentiated in SB media (A) and SupSB media (B) were taken after 16 d. Beating areas are indicated by red circles. Total number of beating aggregates was higher in cultures differentiated in SupSB mediau(about 80% of total aggregates) compared to in SB mediau (about 60% of total aggregates). Scale bar = 200 µm.




Figure 8 Characterization of cardiomyocytes: Immunocytochemistry and electrophysiology. HES-3 differentiated cardiomyocyte aggregates at the end of differentiation (day 16) were mechanically dissociated and plated on 24-well plates and strained for markers. A: Sarcomeric α-actinin (Structural protein); B: Troponin-T (contractile function protein); C: Myosin heavy chain (MHC) (contractile function protein). Nuclei stained with DAPI (blue). Scale bar = 20 µm. MHC, myosin heavy chain; DAPI: 4'6-diamidino-2-phenylindole. D: Whole-cell patch-clamp recording was performed on beating cardiomyocyte aggregates at day 23 of differentiation. The recording shows the successful derivation of all three cardiac phenotypes as well as no difference between cells grown in SB media (control) and SupSB media (supplements); E: Maximal diastolic potential; F: Action potential duration at 90% repolarization (ADP90); G: Action potential duration at 50% repolarization (ADP50).

WJSC | www.wjgnet.com




Figure 9 Karyotype of cardiomyocytes differentiated from HES-3 in SupSB media.

are at. SA is expressed earlier than MHC during cardiac differentiation^[34]. Cultures expressing a higher ratio of SA to MHC (16:1) indicate an earlier stage of differentiation in comparison to cultures with lower ratios.

The results from qRT-PCR analysis showed no acceleration in the differentiation process when cells cultured with supplements (BSA and HySoy) were compared to SB media, with expression of pluripotent and mesoderm markers showing similar trends. However, there was a significantly higher expression of Nkx 2.5, indicating a higher population of cardiac progenitor cells in the cultures with supplements. MHC expression showed a higher general trend, but not a significant increase, when comparing cells cultured with individual supplements to those cultured in SB media only. However, when the supplements were combined, their additive effect produced a significantly higher expression of MHC, indicating that SupSB media improves the differentiation efficiency of cardiomyocytes.

Albumin, as a major serum protein with a typical concentration of 50 mg/mL in blood, accounts for approximately 60% of the total protein in mammals^[35]. Although there is such an abundance of albumin, its physiological actions and molecular mechanisms are not fully understood^[36,37]. The main functions of albumin have been summarized to include (1) maintenance of blood oncotic pressure and pH; (2) binding and transport of physiologically important nutrients, including lipids, metal ions, amino acids and other factors; and (3) antioxidant functions but mainly from the perspective of its role in blood circulation^[35].

Protein hydrolysates are known as a potential source of metabolizable materials including amino acids, oligopeptides, iron salts, lipids and trace elements. Their beneficial effects to the growth of animal cell culture have been known for more than two decades^[38] and are generally thought to act as a concentrated balanced nutrient mixture that may partly or fully replace serum^[39,40]. In recent years suppliers have improved their production processes to generate more consistent products with less batch-to-batch variability for the biopharma industry^[41]. Multiple reports have shown the effect of protein hydrolysates on growth in a variety of cell lines including both animal and rodent cells^[40,42]. This study further validates HySoy as a suitable supplement for cardiomyocyte differentiation.

The increase in cardiomyocyte efficiency was attributed to two reasons: (1) increased cell growth; and (2) increased differentiation efficiency, probably due to the improved metabolic state of the cells. Both BSA and Hy-Soy have been repeatedly reported to increase cell densities in various different cell lines, and thus, the increase in cell growth observed during cardiomyocyte differentiation is not surprising. BSA and HySoy are also more defined in comparison to Fetal Bovine Serum, which is widely used in current EB and monolayer differentiation protocols^[4,10]. We assume that the addition of albumin and HySoy to the differentiation medium improves cardiomyocyte cell yields and viability by their antioxidant activity, conferring needed metabolized nutrients and the ability to transport nutrients (e.g., lipids) to the cells. Previously it has been shown that insulin inhibits the process of cardiomyogenesis, therefore it was removed from the SB media^[43]. BSA and HySoy supplements were able to restore the proliferative capacity of cardiomyocytes in this insulin-free SupSB media. Moreover, it was reported that BSA can also help in the transport of small molecules such as the p38 MAPK inhibitor to the cells^[35]. Maturation of cardiomyocytes into ventricular, atrial and nodal phenotypes occurred within 16 d of differentiation as shown by the electrophysiology characterization. Ventricular phenotype maturation was much faster in comparison to other works which indicated a time frame of $60 d^{[8]}$.

Currently, there are multiple companies offering complete medium for hESC cardiac differentiation, these media use expensive ingredients and thus they are considerably more expensive then SupSB media. This work is the first step in developing an inexpensive and efficient cardiomyocyte differentiation protocol that can be used by commercial companies to produce cardiomyocytes. To the best of our knowledge, this is the first report of the benefits of BSA and HySoy on improving cardiomyogenesis of pluripotent human stem cells. Future works into cardiomyocyte purification and increase cardiomyocyte yield is needed.

In conclusion, we have created a simple, robust and cost effective media that significantly improves cardiomyocyte differentiation over many passages for multiple pluripotent cell lines that will be useful for research and cell therapy applications.

ACKNOWLEDGEMENTS

The authors want to acknowledge Dr. Allen Chen and Dr. Hia Hui Ching for their advice as well as Dr. Filip Laco for reviewing this article.

COMMENTS

Background

Heart disease is one of the most common causes of mortality in the world



and accounts for more than 800000 deaths per year on average in the United States alone. After an episode of major cardiac insult (such as a myocardial infarction), the heart can lose up to 2 billion cardiomyocytes. Being an organ that cannot auto-regenerate, progressive heart failure develops. Cardiomyocyte cell therapy can thus be a potential cure for heart failure after myocardial infarction, and it is suggested that at least one billion cardiomyocytes will be required per patient for such treatment

Research frontiers

Human embryonic stem cells (hESCs) present an attractive cell source for generating large amounts of cardiomyocytes, due to their pluripotency and ability to proliferate for multiple passages. As such, several studies have reported techniques to efficiently differentiate hESCs to cardiomyocytes via growth factors and small molecule inhibitors. Specifically, a protocol based on the inhibition of p38 MAPK with small molecules displayed the ability to be scaled up in a cost efficient manner. However, low cell yield and viability were inherent to this protocol, reducing the output of cardiomyocytes.

Innovations and breakthroughs

Soy hydrolysate (HySoy) and bovine serum albumin (BSA) were found to improve cell viability during cardiomyogenesis of human embryonic and induced pluripotent stem cells in serum free medium. The addition of both supplements leads to an increase in cell viability, cell yield and cardiac marker expression during cardiomyocyte differentiation, translating to an overall increase in cardiomyocyte yield (2.6 fold increase over non supplemented medium).

Applications

Authors have created a simple, robust and cost effective serum free medium (SupSB media) that significantly improves cardiomyocyte differentiation of human pluripotent stem cell in suspended embryoid bodies culture. The new medium can be used for research and cell therapy applications using human embryonic and induced pluripotent cultures.

Terminology

p38 MAPK are a class of protein kinases involved in pathways that deal with differentiation, apoptosis and autophagy. p38 MAPK was inhibited with small molecule, SB203580.

Peer review

The authors have carried out an interesting study in which they developed an inexpensive and efficient cardiomyocyte differentiation protocol that can be used to produce cardiomyocytes. They demonstrated the benefits of BSA and HySoy on improving cardiomyogenesis of pluripotent human stem cells.

REFERENCES

- 1 Roger VL, Go AS, Lloyd-Jones DM, Benjamin EJ, Berry JD, Borden WB, Bravata DM, Dai S, Ford ES, Fox CS, Fullerton HJ, Gillespie C, Hailpern SM, Heit JA, Howard VJ, Kissela BM, Kittner SJ, Lackland DT, Lichtman JH, Lisabeth LD, Makuc DM, Marcus GM, Marelli A, Matchar DB, Moy CS, Mozaffarian D, Mussolino ME, Nichol G, Paynter NP, Soliman EZ, Sorlie PD, Sotoodehnia N, Turan TN, Virani SS, Wong ND, Woo D, Turner MB. Executive summary: heart disease and stroke statistics--2012 update: a report from the American Heart Association. *Circulation* 2012; 125: 188-197 [PMID: 22215894 DOI: 10.1161/CIR.0b013e3182456d46]
- 2 Murry CE, Reinecke H, Pabon LM. Regeneration gaps: observations on stem cells and cardiac repair. J Am Coll Cardiol 2006; 47: 1777-1785 [PMID: 16682301 DOI: 10.1016/ j.jacc.2006.02.002]
- 3 Zweigerdt R. The art of cobbling a running pump--will human embryonic stem cells mend broken hearts? *Semin Cell Dev Biol* 2007; 18: 794-804 [PMID: 18006339 DOI: 10.1016/j.semcdb.2007.09.014]
- 4 Burridge PW, Anderson D, Priddle H, Barbadillo Muñoz MD, Chamberlain S, Allegrucci C, Young LE, Denning C. Improved human embryonic stem cell embryoid body homogeneity and cardiomyocyte differentiation from a novel V-96 plate aggregation system highlights interline variability. *Stem Cells* 2007; 25: 929-938 [PMID: 17185609 DOI: 10.1634/stemcells.2006-0598]
- 5 Burridge PW, Thompson S, Millrod MA, Weinberg S, Yuan

X, Peters A, Mahairaki V, Koliatsos VE, Tung L, Zambidis ET. A universal system for highly efficient cardiac differentiation of human induced pluripotent stem cells that eliminates interline variability. *PLoS One* 2011; **6**: e18293 [PMID: 21494607 DOI: 10.1371/journal.pone.0018293]

- 6 Laflamme MA, Chen KY, Naumova AV, Muskheli V, Fugate JA, Dupras SK, Reinecke H, Xu C, Hassanipour M, Police S, O'Sullivan C, Collins L, Chen Y, Minami E, Gill EA, Ueno S, Yuan C, Gold J, Murry CE. Cardiomyocytes derived from human embryonic stem cells in pro-survival factors enhance function of infarcted rat hearts. *Nat Biotechnol* 2007; 25: 1015-1024 [PMID: 17721512 DOI: 10.1038/nbt1327]
- 7 Yang L, Soonpaa MH, Adler ED, Roepke TK, Kattman SJ, Kennedy M, Henckaerts E, Bonham K, Abbott GW, Linden RM, Field LJ, Keller GM. Human cardiovascular progenitor cells develop from a KDR+ embryonic-stem-cell-derived population. *Nature* 2008; 453: 524-528 [PMID: 18432194 DOI: 10.1038/nature06894]
- 8 Lian X, Hsiao C, Wilson G, Zhu K, Hazeltine LB, Azarin SM, Raval KK, Zhang J, Kamp TJ, Palecek SP. Robust cardiomyocyte differentiation from human pluripotent stem cells via temporal modulation of canonical Wnt signaling. *Proc Natl Acad Sci USA* 2012; 109: E1848-E1857 [PMID: 22645348 DOI: 10.1073/pnas.1200250109]
- 9 Burridge PW, Keller G, Gold JD,Wu JC. Production of de novo cardiomyocytes: human pluripotent stem cell differentiation and direct reprogramming. *Cell Stem Cell* 2012; 10: 16-28 [PMID: 22226352 DOI: 10.1016/j.stem.2011.12.013]
- 10 Tohyama S, Hattori F, Sano M, Hishiki T, Nagahata Y, Matsuura T, Hashimoto H, Suzuki T, Yamashita H, Satoh Y, Egashira T, Seki T, Muraoka N, Yamakawa H, Ohgino Y, Tanaka T, Yoichi M, Yuasa S, Murata M, Suematsu M,Fukuda K. Distinct metabolic flow enables large-scale purification of mouse and human pluripotent stem cell-derived cardiomyocytes. *Cell Stem Cell* 2013; **12**: 127-137 [PMID: 23168164 DOI: 10.1016/j.stem.2012.09.013]
- 11 Kattman SJ, Witty AD, Gagliardi M, Dubois NC, Niapour M, Hotta A, Ellis J, Keller G. Stage-specific optimization of activin/nodal and BMP signaling promotes cardiac differentiation of mouse and human pluripotent stem cell lines. *Cell Stem Cell* 2011; 8: 228-240 [PMID: 21295278 DOI: 10.1016/j.stem.2010.12.008]
- 12 Willems E, Spiering S, Davidovics H, Lanier M, Xia Z, Dawson M, Cashman J, Mercola M. Small-molecule inhibitors of the Wnt pathway potently promote cardiomyocytes from human embryonic stem cell-derived mesoderm. *Circ Res* 2011; **109**: 360-364 [PMID: 21737789 DOI: 10.1161/CIRCRE-SAHA.111.249540]
- 13 Shiba Y, Fernandes S, Zhu WZ, Filice D, Muskheli V, Kim J, Palpant NJ, Gantz J, Moyes KW, Reinecke H, Van Biber B, Dardas T, Mignone JL, Izawa A, Hanna R, Viswanathan M, Gold JD, Kotlikoff MI, Sarvazyan N, Kay MW, Murry CE, Laflamme MA. Human ES-cell-derived cardiomyocytes electrically couple and suppress arrhythmias in injured hearts. *Nature* 2012; **489**: 322-325 [PMID: 22864415 DOI: 10.1038/nature11317]
- 14 Moore JC, Fu J, Chan YC, Lin D, Tran H, Tse HF, Li RA. Distinct cardiogenic preferences of two human embryonic stem cell (hESC) lines are imprinted in their proteomes in the pluripotent state. *Biochem Biophys Res Commun* 2008; 372: 553-558 [PMID: 18503758 DOI: 10.1016/j.bbrc.2008.05.076]
- 15 Pekkanen-Mattila M, Kerkelä E, Tanskanen JM, Pietilä M, Pelto-Huikko M, Hyttinen J, Skottman H, Suuronen R, Aalto-Setälä K. Substantial variation in the cardiac differentiation of human embryonic stem cell lines derived and propagated under the same conditions--a comparison of multiple cell lines. *Ann Med* 2009; **41**: 360-370 [PMID: 19165643 DOI: 10.10 80/07853890802609542]
- 16 **Lecina M**, Ting S, Choo A, Reuveny S, Oh S. Scalable platform for human embryonic stem cell differentiation to car-



WJSC www.wjgnet.com

diomyocytes in suspended microcarrier cultures. *Tissue Eng Part C Methods* 2010; **16**: 1609-1619 [PMID: 20590381 DOI: 10.1089/ten.TEC.2010.0104]

- 17 Ting S, Lecina M, Reuveny S, Oh S. Differentiation of human embryonic stem cells to cardiomyocytes on microcarrier cultures. *Curr Protoc Stem Cell Biol* 2012; Chapter 1: Unit1D.7 [PMID: 22605644 DOI: 10.1002/9780470151808.sc01d07s21]
- 18 Graichen R, Xu X, Braam SR, Balakrishnan T, Norfiza S, Sieh S, Soo SY, Tham SC, Mummery C, Colman A, Zweigerdt R, Davidson BP. Enhanced cardiomyogenesis of human embry-onic stem cells by a small molecular inhibitor of p38 MAPK. *Differentiation* 2008; **76**: 357-370 [PMID: 18021257 DOI: 10.1111/j.1432-0436.2007.00236.x]
- 19 Xu XQ, Graichen R, Soo SY, Balakrishnan T, Rahmat SN, Sieh S, Tham SC, Freund C, Moore J, Mummery C, Colman A, Zweigerdt R, Davidson BP. Chemically defined medium supporting cardiomyocyte differentiation of human embryonic stem cells. *Differentiation* 2008; **76**: 958-970 [PMID: 18557764 DOI: 10.1111/j.1432-0436.2008.00284.x]
- 20 Kempf H, Lecina M, Ting S, Zweigerdt R, Oh S. Distinct regulation of mitogen-activated protein kinase activities is coupled with enhanced cardiac differentiation of human embryonic stem cells. *Stem Cell Res* 2011; 7: 198-209 [PMID: 21907163 DOI: 10.1016/j.scr.2011.06.001]
- 21 **Kirouac DC**, Zandstra PW. The systematic production of cells for cell therapies. *Cell Stem Cell* 2008; **3**: 369-381 [PMID: 18940729 DOI: 10.1016/j.stem.2008.09.001]
- 22 Davis JM. Basic cell culture: a practical approach. 2nd. Oxford: Oxford University Press, 2002
- 23 **Masters JRW.** Animal cell culture: a practical approach. 3rd. Oxford: Oxford University Press, 2000
- 24 **Freshney RI.** Culture of animal cells: a manual of basic technique. 4th. New York, Chichester: Wiley-Liss, 2005
- 25 Ng KM, Chan YC, Lee YK, Lai WH, Au KW, Fung ML, Siu CW, Li RA, Tse HF. Cobalt chloride pretreatment promotes cardiac differentiation of human embryonic stem cells under atmospheric oxygen level. *Cell Reprogram* 2011; 13: 527-537 [PMID: 22029419 DOI: 10.1089/cell.2011.0038]
- 26 Li W, Wang X, Fan W, Zhao P, Chan YC, Chen S, Zhang S, Guo X, Zhang Y, Li Y, Cai J, Qin D, Li X, Yang J, Peng T, Zychlinski D, Hoffmann D, Zhang R, Deng K, Ng KM, Menten B, Zhong M, Wu J, Li Z, Chen Y, Schambach A, Tse HF, Pei D, Esteban MA. Modeling abnormal early development with induced pluripotent stem cells from aneuploid syndromes. *Hum Mol Genet* 2012; **21**: 32-45 [PMID: 21949351 DOI: 10.1093/hmg/ddr435]
- 27 Brunner D, Frank J, Appl H, Schöffl H, Pfaller W, Gstraunthaler G. Serum-free cell culture: the serum-free media interactive online database. *ALTEX* 2010; 27: 53-62 [PMID: 20390239]
- 28 Gerdes J, Lemke H, Baisch H, Wacker HH, Schwab U, Stein H. Cell cycle analysis of a cell proliferation-associated human nuclear antigen defined by the monoclonal antibody Ki-67. J Immunol 1984; 133: 1710-1715 [PMID: 6206131]
- 29 Ebashi S. Ca2+ and the contractile proteins. *J Mol Cell Cardiol* 1984; **16**: 129-136 [PMID: 6232393]
- 30 Allegrucci C, Young LE. Differences between human embryonic stem cell lines. *Hum Reprod Update* 2007; **13**: 103-120 [PMID: 16936306 DOI: 10.1093/humupd/dml041]
- 31 Adewumi O, Aflatoonian B, Ahrlund-Richter L, Amit M, Andrews PW, Beighton G, Bello PA, Benvenisty N, Berry LS, Bevan S, Blum B, Brooking J, Chen KG, Choo AB, Churchill GA, Corbel M, Damjanov I, Draper JS, Dvorak P, Emanuelsson K, Fleck RA,

Ford A, Gertow K, Gertsenstein M, Gokhale PJ, Hamilton RS, Hampl A, Healy LE, Hovatta O, Hyllner J, Imreh MP, Itskovitz-Eldor J, Jackson J, Johnson JL, Jones M, Kee K, King BL, Knowles BB, Lako M, Lebrin F, Mallon BS, Manning D, Mayshar Y, McKay RD, Michalska AE, Mikkola M, Mileikovsky M, Minger SL, Moore HD, Mummery CL, Nagy A, Nakatsuji N, O'Brien CM, Oh SK, Olsson C, Otonkoski T, Park KY, Passier R, Patel H, Patel M, Pedersen R, Pera MF, Piekarczyk MS, Pera RA, Reubinoff BE, Robins AJ, Rossant J, Rugg-Gunn P, Schulz TC, Semb H, Sherrer ES, Siemen H, Stacey GN, Stojkovic M, Suemori H, Szatkiewicz J, Turetsky T, Tuuri T, van den Brink S, Vintersten K, Vuoristo S, Ward D, Weaver TA, Young LA, Zhang W. Characterization of human embryonic stem cell lines by the International Stem Cell Initiative. *Nat Biotechnol* 2007; **25**: 803-816 [PMID: 17572666 DOI: 10.1038/nbt1318]

- 32 Skottman H, Mikkola M, Lundin K, Olsson C, Strömberg AM, Tuuri T, Otonkoski T, Hovatta O, Lahesmaa R. Gene expression signatures of seven individual human embryonic stem cell lines. *Stem Cells* 2005; 23: 1343-1356 [PMID: 16081666 DOI: 10.1634/stemcells.2004-0341]
- 33 Osafune K, Caron L, Borowiak M, Martinez RJ, Fitz-Gerald CS, Sato Y, Cowan CA, Chien KR, Melton DA. Marked differences in differentiation propensity among human embryonic stem cell lines. *Nat Biotechnol* 2008; 26: 313-315 [PMID: 18278034 DOI: 10.1038/nbt1383]
- 34 **Gregorio CC**, Antin PB. To the heart of myofibril assembly. *Trends Cell Biol* 2000; **10**: 355-362 [PMID: 10932092]
- Francis GL. Albumin and mammalian cell culture: implications for biotechnology applications. *Cytotechnology* 2010; 62: 1-16 [PMID: 20373019 DOI: 10.1007/s10616-010-9263-3]
- 36 Quinlan GJ, Martin GS, Evans TW. Albumin: biochemical properties and therapeutic potential. *Hepatology* 2005; 41: 1211-1219 [PMID: 15915465 DOI: 10.1002/hep.20720]
- 37 Ahn SM, Byun K, Cho K, Kim JY, Yoo JS, Kim D, Paek SH, Kim SU, Simpson RJ, Lee B. Human microglial cells synthesize albumin in brain. *PLoS One* 2008; 3: e2829 [PMID: 18665237 DOI: 10.1371/journal.pone.0002829]
- 38 Mizrahi A. Primatone RL in mammalian cell culture media. Biotechnol Bioeng 1977; 19: 1557-1561 [PMID: 901933 DOI: 10.1002/bit.260191013]
- 39 Jan DC, Jones SJ, Emery AN, al-Rubeai M. Peptone, a lowcost growth-promoting nutrient for intensive animal cell culture. *Cytotechnology* 1994; 16: 17-26 [PMID: 7765786]
- 40 Schlaeger EJ. The protein hydrolysate, Primatone RL, is a cost-effective multiple growth promoter of mammalian cell culture in serum-containing and serum-free media and displays anti-apoptosis properties. *J Immunol Methods* 1996; 194: 191-199 [PMID: 8765172]
- 41 Pasupuleti VK, Braun S. State of the Art Manufacturing of Protein Hydrolysates. In: Protein Hydrolysates in Biotechnology. New York: Springer, 2010: 11-33
- 42 Franek F, Hohenwarter O, Katinger H. Plant protein hydrolysates: preparation of defined peptide fractions promoting growth and production in animal cells cultures. *Biotechnol Prog* 2000; 16: 688-692 [PMID: 11027157 DOI: 10.1021/ bp0001011]
- 43 Freund C, Ward-van Oostwaard D, Monshouwer-Kloots J, van den Brink S, van Rooijen M, Xu X, Zweigerdt R, Mummery C, Passier R. Insulin redirects differentiation from cardiogenic mesoderm and endoderm to neuroectoderm in differentiating human embryonic stem cells. *Stem Cells* 2008; 26: 724-733 [PMID: 18096723 DOI: 10.1634/stemcells.2007-0617]

P- Reviewers Houria B, Placantonakis DG, Resende R S- Editor Wen LL L- Editor A E- Editor Lu YJ



