

antibody interfering with ABO grouping and antibody screening. *Transfusion* 1975;14:16-22.

2. Beattie KM, Ferguson SJ, Burnie KL, Barr RM, Urbaniak ST, Atherton PH. Chloramphenicol antibody causing interference in antibody detection and identification tests. *Transfusion* 1976;16:174-7.
3. Gillund TD, Howard PL, Isham B. A serum agglutinating human red cells exposed to EDTA. *Vox Sang* 1972;23:369-70.
4. Reid ME, Battenfield LK, Toy PTCY, Ellisor SS, Hart CA. Agglutination of an EDTA blood sample caused by an EDTA-dependent panagglutinin. *Am J Clin Pathol* 1985;83:534-5.
5. Widmann FK, Ed. Technical manual, 9th ed. Arlington, VA: American Association of Blood Banks, 1985.

Monica Tobar, MT(ASCP), Senior Technologist, Immunohematology Reference Laboratory, American Red Cross Blood Services, Southeastern Michigan Region, PO Box 33351, Detroit, MI 48232.

TECH HINTS

MICROFILTRATION FOR CAPILLARY TUBE USE

Mary N. Crawford


Capillary tube testing requires reagents free of particulate matter and excessive lipids. It is difficult to filter, by the usual methods, the small amounts that are suitable for this technique. A microfilter that can

be assembled rapidly (2 min) and has proved useful for this purpose, is described. A Tampax^R tampon (Tam Brands, Inc., Lake Success, NY) should be pulled apart and enough fibers torn off to pack down to a 0.6 cm (1/4 inch) amount in a small 4.7 cm pipette filter tip (Electra Tips 620/650, Medical Laboratory Automation, Inc., Mt. Vernon, NY) from which the last 0.6 cm has been cut off to enhance the flow. Approximately 1–2 cm of Seitz fibers should be tamped down on top of the Tampax^R fibers. These fibers can be taken from the leftover Seitz material surrounding punched out pads in commercial filtration kits (11 or 22 mm filter discs made for the 2 or 25 mL sized Boerner Centrifugal Filters distributed by Arthur H. Thomas Co.). The microfilter is then inserted into a 400 μ L polyethylene microcentrifuge tube (MC-1 from Analytical Lab Accessories), and two or more drops of serum are put into the microfilter pipette tip. This whole assembly is put into a 1 mL Fisher plastic tube for centrifugation. A hard spin, 7000 rpm for 5 min, is necessary for passage of reagents through the combination of fibers, therefore a Fisher type centrifuge is helpful. For virtually total recovery of reagents, two drops of saline or 6 percent bovine albumin may first be spun through the microfilter and then completely removed from the microcentrifuge tube before centrifugation of the serum.

This procedure is depicted in Figure 1.

Mary N. Crawford, MD, Pearson C. Cummin Memorial Laboratory, 125 Ashwood Road, Villanova, PA 19085.

Figure 1. Improved microfilter for clarification of sera for capillary tube testing*



* With this apparatus, it is possible to filter 2–3 drops of serum. For almost total recovery of serum, 1–2 drops of saline or 6% albumin should be spun through first and then removed.

COMMUNICATIONS

To the Editor:

Hooray for John Judd! By condemning¹ the description of exceptionally powerful examples of anti-N as anti-N', he has struck a blow for clarity in communication between scientists.

If we are to succeed in understanding each other, it is important that we do not allow ourselves to lapse into jargon, or adopt loose terminology that departs from established conventions. The antibodies described in the reports by both Guizzo and Meadows² and Kosanke and Behzad³ were plainly examples of anti-N possessing a sufficient reserve of potency to give detectable agglutination with cells lacking N but possessing 'N'. To call these antibodies anti-N' is surely similar to calling anti-A, anti-A₂ because it agglutinates A₂ (as well as A₁) cells, or anti-D, anti-D^u because it is capable of detecting D^u (as well as D). Are we entitled to call anti-A, B anti-A_x? How about anti-i for an example of anti-I powerful enough that it requires dilution to show no agglutination of umbilical cord cells? Can we think of anti-Fy^b as anti-Fy^x if it detects the weak Fy^b supposedly produced by Fy^x?