



Vista
Analytical Laboratory


Extraction and Analytical Challenges for PFAS in Biosolids

Oliver Cawdell


Introduction

- ▶ PFAS sources in the environment.
- ▶ Why might PFAS be in biosolids, and why does it matter?
- ▶ Analytical challenges for biosolid and sludge matrices.
- ▶ Methods available for biosolids.
- ▶ Comparison of different extraction techniques.
- ▶ Summary of method performance

PFAS sources and pathways within the environment




- ▶ Biosolids are separated from wastewater during WWTP processing.
- ▶ Spreading of biosolids on farmland is common practice in North America.
- ▶ Land applied biosolids can cause PFAS to reenter the food chain creating a positive cycle.



Currently no published EPA method for biosolids

Highly variable
physical properties
and chemistry
between samples


Interferences Affecting Instrument Response: High TOC, Lipids and Industrial contaminants



Analytes with Different Chemical Properties

Hydrophobic analytes at the interface of biphasal samples

Surface sorption, including glassware





Isotope Dilution: Crucial for difficult matrices

Results are automatically recovery-corrected for effects of the analytical process

Mitigates the effects of ion enhancement / suppression caused by interferences

Improved analyte identification when matrix affects retention times

Results are highly accurate and precise



Selecting a method for biosolid analysis

- ▶ No official solid matrix currently available.
 - ▶ 537.1 – method for potable waters only, limited to 18 analytes.
 - ▶ 533 – isotope dilution method for potable waters. Not for solids.
- ▶ Wet or dry solid processing
- ▶ Pre extractions to remove the solid component:
 - ▶ Methanol extraction. Simple extraction to aqueous phase.
 - ▶ Ion Pair. Matrix digestion before solvent extraction. (Zhang et al. 2018)
- ▶ Cleanup via SPE, GCB or Lipid cleanup column.



Method Comparison

Two different variables were tested in this study:


- ▶ Sample processing comparison: Air dried versus centrifuge and separate processing.
- ▶ Solid Matrix Pre-extraction comparison: Ion par extraction (Zhang et al. 2018) versus methanol solid/liquid extraction.

For consistency, these were kept constant:


- ▶ Internal standards were added before extraction (isotope dilution).
- ▶ 1g dry weight equivalent of each sample was extracted.
- ▶ All extracts went through SPE and GCB cleanup before analysis by UPLC-MS/MS.




Internal standard recovery from a lab sand matrix



Internal standard recovery from a biosolid matrix



Internal standard recovery in wet and dried biosolids





Summary

- ▶ Contaminated biosolids are a significant source of PFAS to the environment.
- ▶ Isotope Dilution must be used for confidence in the results for complex matrices such as biosolids.
- ▶ Wet sample extractions yield better recovery of PFAS.
- ▶ Aggressive digestion techniques improve long chain recovery.
- ▶ Clean up should be sample/site specific to remove known contaminants and interferences.



thank you

ocawdell@vista-analytical.com

vista-analytical.com



Vista
Analytical Laboratory