

1 **Performance characteristics of five antigen-detecting rapid diagnostic test (Ag-**
2 **RDT) for SARS-CoV-2 asymptomatic infection: a head-to-head benchmark**
3 **comparison**

4

5 **Authors:**

6 Bàrbara Baro¹ PhD, Pau Rodo² MSc, Dan Ouchi² MSc, Antoni E. Bordoy³ PhD, Emilio
7 N. Saya Amaro³ Bsc, Sergi V. Salsench³ Bsc, Sònia Molinos³ PhD, Andrea Alemany^{2,3,4}
8 MB, Maria Ubals^{2,3,4} MB, Marc Corbacho-Monné^{2,4,5} MB, Pere Millat-Martinez¹ MB,
9 Michael Marks^{6,7} PhD, Bonaventura Clotet^{2,3,8,9} PhD, Nuria Prat¹⁰ MSc, Jordi Ara^{3,10}
10 PhD, Martí Vall-Mayans^{2,3} PhD, Camila G-Beiras² PhD, Quique Bassat^{1,11,12,13,14} PhD,
11 Ignacio Blanco^{3,10} PhD, Oriol Mitjà^{2,3,15} PhD.

12

13 ¹ISGlobal, Hospital Clínic, Universitat de Barcelona, Barcelona, Spain

14 ²Fight AIDS and Infectious Diseases Foundation, Badalona, Spain

15 ³Hospital Universitari Germans Trias i Pujol, Badalona, Spain

16 ⁴Faculty of Medicine and Health Sciences, University of Barcelona, Barcelona, Spain

17 ⁵Parc Taulí Hospital Universitari, Institut d'Investigació i Innovació Parc Taulí I3PT, Sabadell,
18 Spain

19 ⁶Clinical Research Department, Faculty of Infectious & Tropical Diseases, London School of
20 Hygiene & Tropical Medicine, London, UK

21 ⁷Hospital for Tropical Diseases, University College London Hospital, London, UK

22 ⁸IrsiCaixa AIDS Research Institute, Germans Trias i Pujol Research Institute (IGTP), Badalona,
23 Spain

24 ⁹Universitat de Vic-Universitat Central de Catalunya (UVIC-UCC), Vic, Spain

25 ¹⁰Gerència Territorial Metropolitana Nord, Institut Català de la Salut, Barcelona, Spain

26 ¹¹Centro de Investigação em Saúde de Manhiça (CISM), Maputo, Mozambique

27 ¹²ICREA, Pg. Lluís Companys 23, 08010 Barcelona, Spain

28 ¹³Pediatrics Department, Hospital Sant Joan de Déu (University of Barcelona), Barcelona, Spain

29 ¹⁴Consortio de Investigación Biomédica en Red de Epidemiología y Salud Pública (CIBERESP),
30 Madrid, Spain

31 ¹⁵Lihir Medical Centre - InternationalSOS, Lihir Island, Papua New Guinea

32

33 **Running title:** Benchmarking of Ag-RDTs for SARS-CoV-2

34

35 Corresponding author: barbara.baro@isglobal.org

36 Bàrbara Baro, Barcelona Institute for Global Health, Rosselló, 132, 08036 Barcelona, Spain

37

38

39 **Abstract (197 words)**

40 **Background:** Mass testing for early identification and isolation of infectious COVID-19
41 individuals is efficacious for reducing disease spread. Antigen-detecting rapid diagnostic tests
42 (Ag-RDT) may be suitable for testing strategies; however, benchmark comparisons are scarce.

43 **Methods:** We used 286 nasopharyngeal specimens from unexposed asymptomatic individuals
44 collected between December 2020 and January 2021 to assess five Ag-RDTs marketed by Abbott,
45 Siemens, Roche Diagnostics, Lepu Medical, and Surescreen.

46 **Results:** For the overall sample, the performance parameters of Ag-RDTs were as follows: Abbott
47 assay, sensitivity 38.6% (95%CI 29.1–48.8) and specificity 99.5% (97–100%); Siemens,
48 sensitivity 51.5% (41.3–61.6) and specificity 98.4% (95.3–99.6); Roche, sensitivity 43.6% (33.7–
49 53.8) and specificity 96.2% (92.4–98.5); Lepu, sensitivity 45.5% (35.6–55.8) and specificity
50 89.2% (83.8–93.3%); Surescreen, sensitivity 28.8% (20.2–38.6) and specificity 97.8% (94.5–
51 99.4%). For specimens with cycle threshold (Ct) <30 in RT-qPCR, all Ag-RDT achieved a
52 sensitivity $\geq 70\%$. The modelled negative- and positive-predictive value for 1% prevalence were
53 $>99\%$ and $<50\%$, respectively.

54 **Conclusions:** When screening unexposed asymptomatic individuals, two Ag-RDTs achieved
55 sensitivity $\geq 80\%$ for specimens with Ct<30 and specificity $\geq 96\%$. The estimated negative
56 predictive value suggests the suitability of Ag-RDTs for mass screenings of SARS-CoV-2
57 infection in the general population.

58

59 **Key words:** SARS-CoV-2, antigen-detecting rapid diagnostic test, head-to-head comparison,
60 mass screening

61

62 **Introduction**

63 Mass testing for early identification and isolation of individuals infected with the severe acute
64 respiratory syndrome coronavirus 2 (SARS-CoV-2), irrespective of symptoms, is potentially an
65 efficacious strategy to reduce disease transmission [1]. Recent advances on the validation of
66 Antigen-detecting Rapid Diagnostic Tests (Ag-RDTs) show promise to replace central laboratory
67 techniques for epidemiological control of the SARS-CoV-2 through mass testing.

68 Reverse transcription-polymerase chain reaction (RT-qPCR) is the current gold standard for
69 identifying the presence of the SARS-CoV-2 in respiratory specimens [2]. More recently,
70 transcription-mediated amplification (TMA) of the SARS-CoV-2 genome has been added to the
71 repertoire of nucleic acid amplification tests (NAAT) for SARS-CoV-2 detection [3]. Despite
72 their high sensitivity, NAATs are associated with drawbacks that limit their use for community-
73 based testing strategies, including the need for laboratory-processing, high cost, and long
74 turnaround from sampling to results release. Furthermore, there is cumulative evidence indicating
75 that the period of NAAT positivity in infected individuals largely exceeds the time window in
76 which infectious viral particles can be isolated from the respiratory tract, raising doubts about the
77 epidemiological meaning of a NAAT positive result [4].

78 Ag-RDTs, commonly used in diagnosing other infectious diseases, have emerged as an alternative
79 tool that meets the requirements for frequent testing at the point-of-care: rapid turnaround time,
80 low cost, and ease-of-use [5]. Overall, Ag-RDTs have lower sensitivity than NAATs; however,
81 clinical validation studies have consistently reported increasing sensitivities in specimens with
82 higher viral loads. These findings, along with the growing body of evidence on the lack of
83 infectivity of cases with low viral load [6–9] and the potential long-tail of positivity when using
84 highly sensitive methods such as PCR, suggest that frequent testing with Ag-RDTs—even those
85 with low sensitivity—may be more effective than less frequent testing with RT-qPCR or TMA
86 for mass screening campaigns to improve SARS-CoV-2 control [9,10].

87 The performance parameters of Ag-RDTs are mostly based on testing respiratory specimens from
88 clinically suspected cases [11–14] and contacts after exposure to a positive case [15–18].
89 However, the sensitivity bias associated with the viral load leads to high heterogeneity in the
90 reported performance parameters, which strongly depend on the disease status and potential
91 exposure (e.g., symptomatic vs. asymptomatic, contact vs. unexposed) of tested individuals. This
92 heterogeneity precludes comparative analyses between tests assessed in different studies and
93 challenges benchmarking of Ag-RDTs. Furthermore, head-to-head comparisons are scarce,
94 particularly in samples from asymptomatic individuals, the target population of community-based
95 screening strategies [19,20]. In this study, we used fresh nasopharyngeal samples collected in
96 routine mass screening campaigns of unexposed asymptomatic individuals to perform a head-to-
97 head comparison of five Ag-RDTs.

98 **Methods**

99 *Study design*

100 As part of the surveillance program for pandemic control in Catalonia (North-East Spain), the
101 local government launched NAAT-based systematic screenings in areas at high risk of an
102 outbreak. The University Hospital Germans Trias i Pujol processed nasopharyngeal specimens
103 collected in a healthcare area in North-East Spain (i.e., Metropolità Nord) with a catchment
104 population of ~1,400,000 people. These samples enabled us to assess the Ag-RDTs in line with
105 The Foundation for Innovative New Diagnostics (FIND) target product profile for lateral flow
106 assays that directly detect antigens of SARS-CoV-2 antigen assays [21], which recommends at
107 least 100 known negative samples and 100 known positive samples with a documented RT-PCR
108 result. In this study, we used samples collected between December 2020 and January 2021 (i.e.,
109 during the third wave of the epidemic in Spain) with RT-qPCR results available (i.e., data on
110 cycle threshold [Ct]) to perform a head-to-head assessment of five Ag-RDTs. Samples with
111 invalid results in any of the assessed Ag-RDTs were excluded from the analysis.

112 All samples used in this analysis had been collected in the setting of a public health surveillance
113 program, and data were handled according to the General Data Protection Regulation 2016/679
114 on data protection and privacy for all individuals within the European Union and the local
115 regulatory framework regarding data protection. The study protocol was approved by the ethics
116 committee of Hospital Germans Trias i Pujol (Badalona, Spain).

117 *Procedures*

118 Samples consisted of nasopharyngeal swabs collected by health care workers during mass testing
119 of unexposed asymptomatic individuals living in areas at high risk of an outbreak. Swab
120 specimens were placed into sterile tubes containing viral transport media (DeltaSwab Virus,
121 Deltalab; or UTM Universal Transport Medium, Copan). The reference test (i.e., RT-qPCR) was
122 performed on fresh samples stored at 2 – 8 °C for up to 24 hours; samples were then stored up to
123 12h at 2-8 °C until their use for the five Ag-RDTs.

124 RNA for RT-qPCR tests were extracted from fresh samples using the viral RNA/Pathogen
125 Nucleic Acid Isolation kit for the Microlab Starlet or Nimbus platforms (Hamilton, USA),
126 according to the manufacturer's instructions. PCR amplification was conducted according to the
127 recommendations of the 2019-nCoV RT-qPCR Diagnostic Panel of the Centers for Disease
128 Control and Prevention (CDC) (REF), using the Allplex™ 2019-nCoV assay (Seegene, South
129 Korea) on the CFX96 (Bio-Rad, USA) in line with manufacturer's instruction. Briefly, a 25 µL
130 PCR reaction mix was prepared that contained 8 µL of each sample's nucleic acids, 2019-nCoV
131 positive and negative controls, 5 µL of 2019-nCoV MOM (primer and probe mix) and 2 µL of
132 real-time one-step Enzyme. Thermal cycling was performed at the following conditions: 20 min
133 at 50 °C for reverse transcription, followed by 15 min at 95°C, and then 45 cycles of 15 sec at
134 94°C and 30 sec at 58°C. An RT-qPCR was considered positive according to the manufacturer's
135 instructions [22].

136 Index tests included the following Ag-RDTs: PanBio™ COVID-19 Ag Rapid test (Abbott),
137 CLINITEST® Rapid COVID-19 Antigen Test (Siemens), SARS-CoV-2 Rapid Antigen Test

138 (Roche Diagnostics), SARS-CoV-2 Antigen Rapid Test Kit (Lepu Medical), and COVID-19
139 Coronavirus Rapid Antigen Test Cassette (Surescreen). Supplementary Table 1 provides further
140 details regarding the specifications of each test. All Ag-RDT determinations were performed in
141 parallel by two blinded technicians, who used approximately 100 μ L of 1:2 mix of each kit buffer
142 and the sample previously homogenized. Samples were applied directly to the test cassette and
143 incubated for 15 minutes at room temperature before reading results at the naked eye, according
144 to the manufacturer instructions (i.e., the presence of any test line (T), no matter how faint,
145 indicates a positive result).

146 *Outcomes and statistical analysis*

147 We calculated that a sample size of at least 73 positive specimens and 165 negative specimens
148 would give 80% power to estimate overall sensitivity and specificity of Ag-RDT assays in our
149 study. We based our calculation on the expected sensitivity and specificity in asymptomatic
150 population of 65% and 96% [17,23], respectively, fixed precision of the point estimate of 2.5%,
151 and confidence level of 95%. The calculation was in line with FIND recommendations for
152 assessing Ag-RDTs that retrospective assessments should include a minimum of 100 samples per
153 RT-PCR result [21].

154 The primary analysis of the head-to-head comparison was the sensitivity and specificity of each
155 Ag-RDT. Sensitivity and specificity were calculated as defined by Altman et al. [24], and reported
156 as a percentage and the exact binomial 95% confidence interval (CI). Sensitivity was also
157 analysed in a subset of samples with $C_t < 30$, considered at high risk of transmission.

158 Secondary analyses were done assessing discordance between results obtained in each Ag-RDTs.
159 Positive and negative-predictive values for each Ag-RDT at population prevalence between 1%
160 and 15% for SARS-CoV-2 infection were modelled [25] and plotted with the exact binomial 95%
161 CI [26]. All analyses and plots were performed using R version 3.6 [27].

162 *Role of the funding source*

163 The funders of the study had no role in the study conception, design, conduct, data analysis, or
164 writing of the report. All authors had full access to all the data in the study and had final
165 responsibility for the decision to submit for publication.

166

167 **Results**

168 Our sample collection included 316 fresh nasopharyngeal swabs from unexposed asymptomatic
169 individuals who had a RT-qPCR result available. Of these, 30 were excluded because of lack of
170 documented Ct value (n=25), incomplete results due to limited sample volume (n=1), or invalid
171 results in any of the Ag-RDTs (n=4, all of them in the Lepu assay), resulting in a study set of 286
172 samples: 101 (35.3%) with positive RT-qPCR result and 185 (64.7%) with negative RT-qPCR
173 result (Figure 1).

174 The Ct value of samples with positive RT-qPCR result was <30 in 30 (29.7%) samples, 30-to-35
175 in 46 (45.5%), and >35 in 25 (24.8%). The overall sensitivity and specificity of the analysed Ag-
176 RDTs ranged from 28.7% to 51.5% and 89.2% to 99.5%, respectively (Table 1). When
177 considering only RT-qPCR positive samples with Ct <30 (i.e., indicates a high concentration of
178 viral genetic material which is typically associated with a higher risk of infectivity) [28], the
179 sensitivity of Ag-RDTs increased to 76.7% (95% CI 57.7 – 90.7) for the Abbott assay; 86.7%
180 (69.3–96.3) for the Siemens Assay; 83.3% (65.3 – 94.4) for the Roche assay; 83.3% (65.3–94.4)
181 for the Lepu assay; and to 70% (50.6–85.3%) for the Surescreen assay (Figure 2).

182 Of the 286 samples analysed by Ag-RDTs, 222 (77.6%) had concordant results across all Ag-
183 RDT assessed. The 29 samples with concordant positive results across Ag-RDTs were all PCR-
184 positive. Conversely, 37 (19.2%) of 193 specimens with negative results in all Ag-RDTs were
185 PCR positive. Figure 3 shows the distribution of Ag-RDT results in samples with discordant
186 results. The Ag-RDT that most often yielded a positive result in samples with negative results in

187 all other Ag-RDTs was the Lepu assay (n=23; 35.9%), followed by the Siemens assay (n=10,
188 15.6%). Table S2 summarizes the cycle threshold distribution across discordances.

189 To provide an estimate of misidentified cases—either false-positive or false-negative cases—that
190 can be used for making decisions in the public health setting, we modelled the positive and
191 negative predictive value for a prevalence range consistent with a mass screening of unexposed
192 asymptomatic individuals (Figure 4A). For the overall study sample, the estimated positive
193 predictive value (PPV) at a 1% prevalence ranged from 4.1% to 41.9%, with the Lepu assay and
194 the Abbott assay, respectively (Table S3). The estimated PPVs notably increased for the <30 Ct
195 subgroup of samples (Figure 4B), and when prevalence in the population was higher. The
196 estimated negative predictive value (NPV) at 1% prevalence ranged from 99.3% to 99.5%, with
197 the Surescreen assay and the Siemens assay, respectively.

198

199 Discussion

200 In this study, we compared head-to-head the sensitivity and specificity of five Ag-RDTs to screen
201 SARS-CoV-2 infected individuals with unknown exposure and no clinical suspicion of COVID-
202 19. Four of the tested Ag-RDTs (i.e., Abbott, Siemens, Roche, and Surescreen assays) showed a
203 specificity higher than 96%. Regarding sensitivity, despite it was low for the overall sample
204 (range 29% to 51%), the corresponding values for the subset of samples with a RT-qPCR value
205 Ct <30 were higher than 80% for the Siemens, Roche, and Lepu assays. This finding is of
206 particular interest for the proposed use of Ag-RDT as a reliable alternative to RT-qPCR for the
207 rapid detection of individuals with higher risk of infectivity in mass screening of asymptomatic
208 individuals. Pre-clinical studies have persistently reported a very low infectious capacity of
209 respiratory specimens with viral loads below 10^6 genome copies/mL, which usually correspond
210 to a Ct of approximately 29 – 31 [5,8,29]. These findings align with the significant increase of the
211 secondary attack rate for values of Ct <30 [30], indicating higher infectiousness among
212 individuals with viral loads below this Ct threshold.

213 Although sensitivity and specificity are important intrinsic characteristics of a test, the number of
214 expected errors when using the test for screening purposes strongly depends on the prevalence of
215 the infection in the screened sample. Hence, positive and negative predictive values are a mainstay
216 for making public health decisions regarding the use of a test. The reported prevalence of SARS-
217 CoV-2 infection in PCR-based untargeted screenings of the general population typically ranges
218 between 1% and 3%, depending on the virus transmission context [23,31]. In low prevalence
219 settings, Ag-RDTs will have a high NPV but a low PPV. According to our estimate, the NPV for
220 SARS-CoV-2 infections at 1% prevalence was higher than 99% for all test, suggesting that a
221 negative test may not require confirmation. In contrast, the PPV at 1% prevalence was lower than
222 50% in all tests, suggesting that a positive result will need immediate confirmation by RT-qPCR,
223 even for highly specific assays.

224 Our study has several strengths and limitations. We used the same fresh set of samples for
225 assessing five different Ag-RDTs and the sample size met the FIND recommendation for
226 retrospective assessments of the clinical performance of these tests. Furthermore, to our
227 knowledge, this is the first head-to-head comparison of Ag-RDT in asymptomatic screenings, an
228 intended use proposed by various authors [5,10,17,23]. On the other hand, our study was limited
229 by the small number of specimens with Ct <30, a threshold deemed of interest for the use of Ag-
230 RDT in screenings of the general population. In our sample, specimens below this threshold
231 accounted for 30%; however, other authors have reported proportions of nearly 60% in random
232 screenings of the general population [23]. Of note, we used specimens in transport medium. This
233 approach is convenient for mass screening strategies in which individuals with positive Ag-RDT
234 results may need further diagnostic confirmation by PCR. However, only one manufacturer (i.e.,
235 the Roche assay) provided instructions on how to process samples collected in virus transport
236 medium. The consistency of our results across assays, particularly regarding negative results,
237 suggests that the use of this media had a little or negligible impact on test performance. Finally,
238 it is worth mentioning that all nasopharyngeal swabs in our analysis were collected by trained

239 healthcare professionals. According to a recent report of lateral flow viral antigen detection
240 devices, the positivity rate might be lower in screenings performed by non-trained people [9].

241 Our results provide policymakers with evidence on the use of Ag-RDT for mass screening of
242 unexposed, asymptomatic individuals. Two commercial, widely available assays can be used for
243 SARS-CoV-2 antigen testing to achieve sensitivity in specimens with a $Ct < 30$ and specificity of
244 at least 80% and 96%, respectively. While these tests may overlook SARS-CoV-2 infection with
245 low viral loads, they accurately detect individuals with high viral loads and, therefore, at higher
246 risk of transmission. Our findings also support the idea that Ag-RDTs can be used for mass
247 screening in low prevalence settings and accurately rule out a highly infectious case in such
248 setting. In models according to population prevalence, all Ag-RDTs will have a NPV $> 99\%$ and
249 a PPV $< 50\%$ at 1% prevalence. Our results, together with the cumulative evidence on the limited
250 overlapping between PCR positivity and the presence of infectious viral particles in the
251 respiratory tract, encourage the design of public health interventions for containing viral COVID-
252 19 spread that shift from positivity testing to infectivity testing. The low cost and short turnaround
253 time of Ag-RDTs, which ease frequent testing, are additional advantages over assays better suited
254 for diagnostic use like NAATs. In low-income countries with limited laboratory resources, the
255 trade-off between targeted PCR analyses and massive screenings with Ag-RDTs should be
256 carefully considered.

257

258 *Funding*

259 This work was supported by Blueberry diagnostics, Fundació Institut d'Investigació en Ciències
260 de la Salut Germans Trias i Pujol, and #YoMeCorono.org crowdfunding campaign.

261 Bárbara Baro is a Beatriu de Pinós postdoctoral fellow granted by the Government of Catalonia's
262 Secretariat for Universities and Research, and by Marie Skłodowska-Curie Actions COFUND
263 Programme (BP3, 801370)

264 ISGlobal receives support from the Spanish Ministry of Science and Innovation through the
265 “Centro de Excelencia Severo Ochoa 2019-2023” Program (CEX2018-000806-S), and support
266 from the Generalitat de Catalunya through the CERCA Program.

267

268 *Acknowledgements*

269 The authors would like to thank Gerard Carot-Sans (PhD) for providing professional medical
270 writing support during the preparation of the manuscript. We thank Laia Comí, Josefa Gómez,
271 Maria Pilar Rodríguez and Aida Sanz for technical support with samples selection and storing.

272 *Contributors*

273 OM, BB, IB designed the study. PR, BB, AEB, SS, ESA, AA, MU, MCM, PMM performed the
274 laboratory procedures, and organized the data. BB, PR and DO verified the underlying data. DO
275 did statistical analysis. BB, OM wrote the first draft with revisions and input from IB, QB, CGB,
276 MVM, NP, JA, BC, MM. Funding acquisition by OM, IB, QB, JA, BC. All authors approved the
277 final version.

278 *Declaration of interests*

279 We declare no conflicts of interest.

280

281

282

283 **References**

- 284 1. Pavelka M, Van-Zandvoort K, Abbott S, et al. The effectiveness of population-wide, rapid
285 antigen test based screening in reducing SARS-CoV-2 infection prevalence in Slovakia.
286 medRxiv **2020**;
- 287 2. World Health Organization (WHO). Diagnostic testing for SARS-CoV-2. Available at:
288 <https://www.who.int/publications/i/item/diagnostic-testing-for-sars-cov-2>. Accessed 6
289 October 2020.
- 290 3. Gorzalski AJ, Tian H, Laverdure C, et al. High-Throughput Transcription-mediated
291 amplification on the Hologic Panther is a highly sensitive method of detection for SARS-
292 CoV-2. *J Clin Virol* **2020**; 129:104501.
- 293 4. Sethuraman N, Jeremiah SS, Ryo A. Interpreting Diagnostic Tests for SARS-CoV-2.
294 *JAMA - J Am Med Assoc* **2020**; 323:2249–2251.
- 295 5. Mina MJ, Parker R, Larremore DB. Rethinking COVID-19 test sensitivity—A strategy
296 for containment. *N Engl J Med* **2020**; 383:10.1056/NEJMp2025631.
- 297 6. Wölfel R, Corman VM, Guggemos W, et al. Virological assessment of hospitalized
298 patients with COVID-2019. *Nature* **2020**; 581:465–469.
- 299 7. La Scola B, Le Bideau M, Andreani J, et al. Viral RNA load as determined by cell culture
300 as a management tool for discharge of SARS-CoV-2 patients from infectious disease
301 wards. *Eur J Clin Microbiol Infect Dis* **2020**; 39:1059–1061.
- 302 8. Quicke K, Gallichote E, Sexton N, et al. Longitudinal Surveillance for SARS-CoV-2 RNA
303 Among Asymptomatic Staff in Five Colorado Skilled Nursing Facilities: Epidemiologic,
304 Virologic and Sequence Analysis. *medRxiv Prepr Serv Heal Sci* **2020**;
305 :2020.06.08.20125989.
- 306 9. PHE Porton Down & University of Oxford. Preliminary report from the Joint PHE Porton
307 Down & University of Oxford SARS-CoV-2 test development and validation cell: Rapid

- 308 evaluation of Lateral Flow Viral Antigen detection devices (LFDs) for mass community
309 testing. Available at: [https://www.ox.ac.uk/sites/files/oxford/media_wysiwyg/UK](https://www.ox.ac.uk/sites/files/oxford/media_wysiwyg/UK_evaluation_PHE_Porton_Down_University_of_Oxford_final.pdf)
310 [evaluation_PHE Porton Down University of Oxford_final.pdf](https://www.ox.ac.uk/sites/files/oxford/media_wysiwyg/UK_evaluation_PHE_Porton_Down_University_of_Oxford_final.pdf). Accessed 3 January 2020.
- 311 10. Larremore DB, Wilder B, Lester E, et al. Test sensitivity is secondary to frequency and
312 turnaround time for COVID-19 screening. *Sci Adv* **2020**; :eabd5393.
- 313 11. Porte L, Legarraga P, Vollrath V, et al. Evaluation of novel antigen-based rapid detection
314 test for the diagnosis of SARS-CoV-2 in respiratory samples. *Int J Infect Dis* **2020**;
315 99:328–333.
- 316 12. Lambert-Niclot S, Cuffel A, Le Pape S, et al. Evaluation of a rapid diagnostic assay for
317 detection of sars-cov-2 antigen in nasopharyngeal swabs. *J Clin Microbiol* **2020**; 58:0–3.
- 318 13. Lindner AK, Nikolai O, Kausch F, et al. Head-to-head comparison of SARS-CoV-2
319 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-
320 collected nasopharyngeal swab. *medRxiv*. 2020; :2003961.
- 321 14. Albert E, Torres I, Bueno F, et al. Field evaluation of a rapid antigen test (Panbio™
322 COVID-19 Ag Rapid Test Device) for COVID-19 diagnosis in primary healthcare centres.
323 *Clin Microbiol Infect* **2020**; 0.
- 324 15. Linares M, Pérez-Tanoira R, Carrero A, et al. Panbio antigen rapid test is reliable to
325 diagnose SARS-CoV-2 infection in the first 7 days after the onset of symptoms. *J Clin*
326 *Virol* **2020**; 133:104659.
- 327 16. Cerutti F, Burdino E, Milia MG, et al. Urgent need of rapid tests for SARS CoV-2 antigen
328 detection: Evaluation of the SD-Biosensor antigen test for SARS-CoV-2. *J Clin Virol*
329 **2020**; 132:104654.
- 330 17. Alemany A, Baró B, Ouchi D, et al. Analytical and clinical performance of the panbio
331 COVID-19 antigen-detecting rapid diagnostic test. *J Infect* **2021**; 18:16.
- 332 18. Torres I, Poujois S, Albert E, Colomina J, Navarro D. Evaluation of a rapid antigen test

- 333 (Panbio™ COVID-19 Ag rapid test device) for SARS-CoV-2 detection in asymptomatic
334 close contacts of COVID-19 patients. *Clin Microbiol Infect* **2021**; •:1.
- 335 19. Weitzel T, Legarraga P, Iruretagoyena M, et al. Comparative evaluation of four rapid
336 SARS-CoV-2 antigen detection tests using universal transport medium. *Travel Med Infect*
337 *Dis* **2021**; 39:101942.
- 338 20. Corman VM, Haage VC, Bleicker T, et al. Comparison of seven commercial SARS-CoV-
339 2 rapid Point-of-Care Antigen tests. *medRxiv*. 2020; :2020.11.12.20230292.
- 340 21. Foundation for Innovative New Diagnostics. Comparative evaluation of lateral flow assay
341 tests that directly detect antigens of SARS-CoV-2. Available at:
342 [https://www.finddx.org/wp-content/uploads/2020/04/20200421-COVID-Ag-RDT-](https://www.finddx.org/wp-content/uploads/2020/04/20200421-COVID-Ag-RDT-Evaluation-Synopsis.pdf)
343 [Evaluation-Synopsis.pdf](https://www.finddx.org/wp-content/uploads/2020/04/20200421-COVID-Ag-RDT-Evaluation-Synopsis.pdf). Accessed 20 December 2020.
- 344 22. Seegene Inc. Allplex™ 2019-nCoV Assay. Available at:
345 https://www.seegene.com/assays/allplex_2019_ncov_assay. Accessed 20 December
346 2020.
- 347 23. Pilarowski G, Lebel P, Sunshine S, et al. Performance characteristics of a rapid SARS-
348 CoV-2 antigen detection assay at a public plaza testing site in San Francisco. *medRxiv*
349 **2020**;
- 350 24. Altman DG, Bland JM. Statistics Notes: Diagnostic tests 1: Sensitivity and specificity.
351 *BMJ* **1994**; 308:1552.
- 352 25. Altman DG, Bland j. M. Statistics Notes: Diagnostic tests 2: Predictive values. *BMJ* **1994**;
353 309:102.
- 354 26. Collett D. *Modelling Binary Data*. Second. Chapman and Hall/CRC, 2002.
- 355 27. R Core Team. *R: A language and environment for statistical computing*. 2017. Available
356 at: <https://www.r-project.org>. Accessed 25 May 2020.
- 357 28. World Health Organization (WHO). Antigen-detection in the diagnosis of SARS-CoV-2

358 infection using rapid immunoassays. 2020. Available at:
359 [https://www.who.int/publications/i/item/antigen-detection-in-the-diagnosis-of-sars-cov-](https://www.who.int/publications/i/item/antigen-detection-in-the-diagnosis-of-sars-cov-2infection-using-rapid-immunoassays)
360 [2infection-using-rapid-immunoassays](https://www.who.int/publications/i/item/antigen-detection-in-the-diagnosis-of-sars-cov-2infection-using-rapid-immunoassays). Accessed 29 September 2020.

361 29. Singanayagam A, Patel M, Charlett A, et al. Duration of infectiousness and correlation
362 with RT-PCR cycle threshold values in cases of COVID-19, England, January to May
363 2020. *Eurosurveillance* **2020**; 25:2001483.

364 30. Marks M, Millat P, Ouchi D, et al. Transmission of Covid-19 in 282 clusters in Catalonia,
365 Spain: a cohort study. *Lancet Infect Dis* **2020**; :Manuscript accepted (in Press).

366 31. Oran DP, Topol EJ. Prevalence of Asymptomatic SARS-CoV-2 Infection : A Narrative
367 Review. *Ann Intern Med* **2020**; 173:362–367.

368

369

370 **Tables**

371 **Table 1. Sensitivity and specificity of the antigen-detecting rapid diagnostic tests for SARS-CoV-2.**

	<i>Abbott</i>	<i>Siemens</i>	<i>Roche</i>	<i>Lepu</i>	<i>Surescreen</i>
Overall Sensitivity	38·61% (29·09-48·82)	51·49% (41·33-61·55)	43·56% (33·72-53·8)	45·54% (35·6-55·76)	28·71% (20·15-38·57)
Detected	39	52	44	46	29
Not Detected	62	49	57	55	72
Total PCR+	101	101	101	101	101
Sensitivity in specimens with Ct<30	76·67% (57·72-90·07)	86·67% (69·28-96·24)	83·33% (65·28-94·36)	83·33% (65·28-94·36)	70% (50·6-85·27)
Detected	23	26	25	25	21
Not Detected	7	4	5	5	9
Total PCR+	30	30	30	30	30
Specificity	99·46% (97·03-99·99)	98·38% (95·33-99·66)	96·22% (92·36-98·47)	89·19% (83·8-93·27)	97·84% (94·56-99·41)
Detected	1	3	7	20	4
Not detected	184	182	178	165	181
Total PCR-	185	185	185	185	185

372

373 All samples were nasopharyngeal swabs collected from unexposed asymptomatic individuals during mass
374 screening campaigns. Sensitivity and specificity results are presented with the 95% confidence interval.

375

376

377

378

379 **Figure Legends**

380

381 **Figure 1. Flow-chart of sample inclusion.**

382 All samples were nasopharyngeal swabs collected from unexposed asymptomatic individuals
383 during screening campaigns.

384

385 **Figure 2. Sensitivity of the antigen-detecting rapid diagnostic tests according to the cycle**
386 **threshold value of the RT-qPCR analysis.**

387 Bars show the 95% confidence interval of the estimated sensitivity.

388

389 **Figure 3. Discordance analysis between Ag-RDTs.**


390 Bars show the number of samples for each discordance pattern. Black dots and grey dots indicate
391 the assays showing positive and negative results in each discordance pattern. Table S2
392 summarizes the cycle threshold distribution across discordances.


393


394 **Figure 4. Positive Predictive Value and Negative Predictive Value according to pre-test**
395 **probabilities.**

396 **A:** overall sample (n= 286). **B:** samples with cycle threshold <30 in the RT-qPCR assay. Table
397 S3 provides detailed values and confidence intervals for predicted false negative and false
398 positives in the investigated prevalence.

399











A



B





Click here to access/download

Supplementary file

Supplementaryfile__tables.docx

