

1 **Monoclonal antibody-based immunoassays for**
2 **cyprodinil residue analysis in QuEChERS-based fruit**
3 **extracts**

4 *Francesc A. Esteve-Turrillas^a, Antonio Abad-Somovilla^b, Guillermo Quiñones-*
5 *Reyes^b, Consuelo Agulló^b, Josep V. Mercader^{a,*}, Antonio Abad-Fuentes^{a,*}*

6
7 ^a *Department of Biotechnology, Institute of Agrochemistry and Food Technology, Consejo*
8 *Superior de Investigaciones Científicas (IATA–CSIC), Agustí Escardino 7, 46980 Paterna,*
9 *València, Spain*

10 ^b *Department of Organic Chemistry, Universitat de València, Doctor Moliner 50, 46100*
11 *Burjassot, València, Spain*

12

13

14 * Corresponding authors. Tel.: +34-963900022; fax: +34-963636301.

15 *E-mail addresses: jmercader@iata.csic.es (Josep V. Mercader), aabad@iata.csic.es (A.*
16 *Abad-Fuentes).*

17

18

19

20 **Abstract**

21 Cyprodinil is among the most common agrochemical residues found in highly perishable fruits,
22 such as strawberries. In the present study, high-affinity monoclonal antibodies to this
23 anilinopyrimidine fungicide were raised for the first time with the aim to produce valuable
24 immunochemical analytical assays. Cyprodinil bioconjugates and the generated novel
25 monoclonal antibodies were employed for sensitive competitive immunoassay development in
26 two different formats. The limits of detection of the optimized assays were 20 and 30 ng L⁻¹ for
27 the indirect and direct assay, respectively. Influence over assay parameters of different
28 physicochemical factors was studied. Strawberry samples were extracted following the
29 recommended QuEChERS procedure for pesticide residues in food, and analyzed by the
30 optimized immunoassays. Recoveries and coefficients of variation from fortified samples were
31 within standard values. In addition, the obtained immunochemical results with naturally
32 contaminated samples were statistically comparable, according to Deming regression analysis,
33 to those of a reference chromatographic method.

34

35 **Keywords**

36 ELISA, Anilinopyrimidine, Fungicide residues, Rapid methods, Food safety, Hapten.

37

38 **Chemical compounds studied in this article:**

39 Cyprodinil (PubChem CID: 86367); Pyrimethanil (PubChem CID: 91650); Mepanipyrim

40 (PubChem CID: 86296); Fluodioxonil (PubChem CID: 86398); Picoxystrobin (PubChem CID:

41 11285653).

42

43

44 1. Introduction

45 Cyprodinil is a systemic novel-generation fungicide of the anilinopyrimidine family that
46 was introduced by Ciba-Geigy in the mid-1990s. Its use was first approved by US authorities in
47 June 1998 and it was included in Annex I of Directive 91/414/EEC of the EU in May 2007.
48 Nowadays, cyprodinil is employed by foliar application in a variety of cereal, fruit, and
49 vegetable crops world-wide. It shows protective and curative effects against highly aggressive
50 pathogens such as *E. graminis*, *D. teres*, *V. inaequalis*, and *B. cinerea* (Knauf-Beiter, Dahmen,
51 Heye, & Staub, 1995). Unlike many fungicides that inhibit spore germination,
52 anilinopyrimidines show a different mode of action that consists of inhibiting germ-tube
53 elongation and initial mycelial growth (Rosslenbroich & Stuebler, 2000). Moreover, several
54 studies suggest that anilinopyrimidines also inhibit the secretion of extracellular proteins,
55 including fungal hydrolases associated with pathogenesis, and the biosynthesis of methionine
56 in *Botrytis* spp. (Masner, Muster, & Schmid, 1994; Milling & Richardson, 1995). Such distinct
57 biocide mechanisms make anilinopyrimidines very attractive plant protection products for
58 combined applications aiming to increase the pesticide efficiency and to prevent fungi
59 resistance. Thus, dual fungicide formulations containing cyprodinil have experienced great
60 success in the present decade, and due to extensive use in particular crops, residues are
61 commonly found in food samples. According to the 2010 EU Report on Pesticide Residues,
62 cyprodinil contents above the MRLs were found mainly in leek but also in apple, head cabbage,
63 lettuce, oat, peach, rye, strawberry, and tomato samples, and it was the most frequently found
64 pesticide in strawberries – almost 32% of the analyzed samples contained residues above the
65 detection limit (Scientific Report of EFSA, 2013). In 2012, cyprodinil residues were mainly
66 found in baby food, cherry tomatoes, plums, and winter squash in US food commodities (USDA
67 Pesticide Data Program, 2014). Despite of being moderately toxic to mammals, there is some
68 concern about a potential cyprodinil bioaccumulation effect – the FAO acceptable daily intake
69 (ADI) is 0.03 mg kg⁻¹ of body weight (DG SANCO Review Report, 2010).

70 Classical analytical methods for pesticide residue determination were based on gas
71 chromatography (GC) with electron capture or nitrogen–phosphorous detection, as well as
72 high-performance liquid chromatography (HPLC) with UV detectors. Nowadays, identification
73 and quantification is mainly performed with mass spectrometry (MS) or tandem MS detection.
74 Since 2007, the AOAC International and the European Committee for Standardization have
75 adopted the QuEChERS (Quick, Easy, Cheap, Effective, Rugged, and Safe) methodology
76 (Anastassiades, Lehotay, Stajnbaher, & Schenck, 2003) for pesticide residue extraction from
77 food matrices (Lehotay, 2007; European Committee for Standardisation, 2008). These
78 processing and measurement methods provide highly sensitive and robust analysis, and
79 possess the capacity to determine simultaneously a variety of residues. Recently, several
80 studies have been published describing chromatographic multiresidue analytical approaches
81 comprising cyprodinil analysis in wine (Rial-Otero, Yagüe-Ruiz, Cancho-Grande, & Simal-
82 Gándara, 2002; Vaquero-Fernández et al., 2008; Zhang et al., 2009; Fontana, Rodriguez, Ramil,
83 Altamirano, & Cela, 2011; Moeder, Bauer, Popp, van Pinxteren, & Reemtsma, 2012), vegetable
84 (Melo, Aguiar, Mansilha, Pinho, & Ferreira, 2012), and strawberry (Fernandes, Domingues,
85 Mateus, & Delerue-Matos, 2012) samples. Additionally, contemporary analysts requiring
86 sensitive, simple, rapid, specific, and high sample throughput approaches may recourse to
87 molecular receptor-based techniques, like the competitive enzyme-linked immunosorbent
88 assay (cELISA). Essentially, this method employs either rabbit polyclonal or mouse monoclonal
89 antibodies (mAbs) as receptors, together with a protein conjugate of a chemical derivative
90 mimicking the target compound. Both types of immunoglobulins are usually highly specific and
91 can reach very high affinities to pesticides – K_D values lower than 10^{-9} M. However, only
92 monoclonal antibodies guarantee a constant supply of well-defined and homogeneous reagent
93 for full immunoassay characterization, development, and validation, as well as wide
94 implementation in analytical laboratories through readily available commercial sources.

95 Structurally, cyprodinil is constituted by an anilinopyrimidine backbone with a methyl
96 and a cyclopropyl substituent at the *meta* positions of the pyrimidine ring. As prediction of the
97 optimum linker tethering site for immunizing haptens was not obvious, four cyprodinil
98 derivatives were prepared holding the spacer arm at alternative sites (Fig. 1). Previously, the
99 immunogenicity of the different haptens was assessed by raising rabbit polyclonal antibodies
100 (Esteve-Turrillas, Agulló, Abad-Fuentes, Abad-Somovilla, & Mercader, 2012; Esteve-Turrillas,
101 Mercader, Agulló, Abad-Somovilla, & Abad-Fuentes, 2015). The aim of the present study was
102 to generate high-affinity mAbs to cyprodinil and to develop and validate immunoassays for
103 rapid and sensitive analysis of residues of this fungicide in QuEChERS-based fruit extracts.
104 Strawberries are highly perishable fruits heavily consumed by children; moreover, they cannot
105 be peeled easily, and due to their rough surface, residues are hardly washed out, which makes
106 strawberries a complex and relevant commodity to evaluate the performance of the novel
107 methods. As far as we know, these are the first reported mAbs and immunoassays for the
108 determination of cyprodinil levels in this added-value foodstuff.

109 **2. Materials and methods**

110 *2.1. Reagents and instruments*

111 Cyprodinil (4-cyclopropyl-6-methyl-*N*-phenylpyrimidin-2-amine, CAS registry number
112 121552-61-2, Mw 225.29 g mol⁻¹) standard, Pestanal® grade, and other pesticide standards
113 were purchased from Fluka/Riedel-de-Haën (Seelze, Germany) or Dr. Ehrenstorfer (Augsburg,
114 Germany). Hybridoma fusion and cloning supplement (HFCS) was obtained from Roche Applied
115 Science (Mannheim, Germany). P3-X63-Ag 8.653 mouse plasmacytoma cell line was acquired
116 from the European Collection of Cell Cultures (Salisbury, UK). Cell culture media (high-glucose
117 Dulbecco's modified Eagle's medium, DMEM), gentamicin solution, and
118 hypoxanthine–thymidine (HT) and hypoxanthine–aminopterin–thymidine (HAT) supplements
119 were purchased from Gibco BRL (Paisley, UK). Poly(ethylene glycol) (PEG1500), fetal bovine

120 serum, 200 mM alanyl–glutamine solution, red blood cell lysing buffer Hybri-Max[®], MEM non-
121 essential amino acid solution, Freund’s adjuvants, triphenylphosphate (TPP), and *o*-
122 phenylenediamine were obtained from Sigma/Aldrich (Madrid, Spain). HiTrap™ protein G HP
123 columns for mouse IgG purification were procured from General Electric Healthcare (Uppsala,
124 Sweden). Polyclonal rabbit anti-mouse immunoglobulins peroxidase conjugate (RAM–HRP)
125 was acquired from Dako (Glostrup, Denmark). Primary/secondary amine and organic solvents
126 were from Scharlab (Barcelona, Spain). Culture plastic ware and Costar[®] flat-bottom high-
127 binding polystyrene ELISA plates were from Corning (Corning, NY, USA).

128 An ELx405 washer and a PowerWave HT reader, both from BioTek Instruments
129 (Winooski, VE, USA), were employed for microtiter analysis. Absorbance was read in dual
130 wavelength mode using 650 nm as reference. Hapten conjugation was verified with a 5800
131 matrix-assisted laser desorption ionization time-of-flight mass spectrometry (MALDI–TOF/TOF)
132 apparatus from AB Sciex (Framingham, MA, USA). Ultra-performance liquid chromatography
133 (UPLC) was performed on a Waters Acquity UPLC system from Waters (Milford, MA, USA),
134 equipped with a binary solvent delivery system, an autosampler, and a BEH C18 (1.7 μm, 2.1 ×
135 50 mm) column. Tandem mass acquisitions were carried out in an Acquity triple quadrupole
136 MS detector, also from Waters, equipped with a Z-spray electrospray ionization source, with
137 3.5 kV capillary voltage, 120 °C source temperature and 300 °C desolvation temperature.

138 2.2. *Protein conjugates*

139 Bovine serum albumin (BSA), ovalbumin (OVA), and horseradish peroxidase (HRP)
140 conjugates from previous studies were employed for immunization, plate coating, and direct
141 assays, respectively (Esteve-Turrillas et al., 2012; Esteve-Turrillas et al., 2015). The hapten-to-
142 protein molar ratio (MR) of each conjugate was determined by MALDI–TOF. For this purpose,
143 conjugate samples were extensively dialyzed in Milli-Q water and 1 μL at ca 0.3 μg μL⁻¹ of
144 every sample solution was spotted onto the MALDI plate. After the droplets were air-dried at

145 room temperature, 0.75 μL of matrix – 5 mg mL^{-1} sinapinic acid (Bruker) in 0.1%
146 TFA– $\text{CH}_3\text{CN}/\text{H}_2\text{O}$ (7:3, v/v) – was added and allowed to air-dry at room temperature. The
147 resulting mixtures were analyzed in positive linear mode – 1500 shots every position – with a
148 mass range of 10000–120000 m/z . Previously, the plate was calibrated with 1 μL of the
149 TOF/TOF calibration mixture (AB Sciex) in 13 positions. Every sample was calibrated by close
150 external calibration method with a BSA, OVA, or HRP spectrum acquired in a nearby position.

151 2.3. *Monoclonal antibodies*

152 Animal manipulation was carried out in compliance with the laws and guidelines of the
153 Spanish Ministry of Agriculture, Food, and Environment. Mice were immunized with BSA
154 conjugates of haptens CDb and CDm (Fig. 1), following previously published schedules (Suárez-
155 Pantaleón, Mercader, Agulló, Abad-Somovilla, & Abad-Fuentes, 2008). Briefly, each mouse
156 received an intraperitoneal injection with 200 μL of a 1:1 emulsion between complete Freund's
157 adjuvant and 100 mM phosphate buffer, pH 7.4, containing 100 μg of BSA–hapten conjugate.
158 Next, two additional boosts were given at 21-day intervals with an equivalent immunogen
159 prepared with incomplete Freund's adjuvant. After a resting period of at least 21 days, and
160 four days before cell fusion, a fourth injection was applied without adjuvant using the same
161 immunizing conjugate dose as before. Cell fusion and hybridoma cloning were performed as
162 described elsewhere (Mercader, Suárez-Pantaleón, Agulló, Abad-Somovilla, & Abad-Fuentes,
163 2008a). High-affinity antibody producing hybridomas were selected following a double-
164 screening technique using culture supernatants as published before [Mercader, Suárez-
165 Pantaleón, Agulló, Abad-Somovilla, & Abad-Fuentes, 2008b]. In brief, parallel and
166 simultaneous assays were run in a first screening procedure with 50 μL per well of culture
167 supernatant plus 50 μL per well of 10 mM phosphate buffered saline, pH 7.4, containing 140
168 mM NaCl (PBS) with 100 nM cyprodinil or without cyprodinil. Second, supernatants affording
169 saturated signals or noticeable inhibition in the first screening assay were reevaluated by

170 checkerboard indirect cELISA. Hybridomas were cloned twice by limiting dilution; then, they
171 were stabilized and expanded in culture plates. Finally, antibodies were purified from late
172 stationary-phase culture supernatants by salting out with saturated ammonium sulfate
173 solution and by affinity chromatography. Citrate buffer 100 mM, pH 2.5, was used for antibody
174 elution and fractions were immediately neutralized with Tris-HCl 1 M, pH 9.5.

175 2.4. *Competitive immunoassays*

176 Monoclonal antibody-coated direct cELISAs were carried out with plates coated by
177 overnight incubation at room temperature with 100 μL per well of a $1 \mu\text{g mL}^{-1}$ antibody
178 solution in 50 mM carbonate–bicarbonate buffer, pH 9.6. After each incubation step, plates
179 were washed four times with 150 mM NaCl containing 0.05% (v/v) Tween 20. Each washed
180 microwell received 50 μL of cyprodinil standard solution in PBS plus 50 μL of enzyme tracer
181 solution in PBST (PBS containing 0.05% (v/v) Tween 20). The competitive reaction was carried
182 out during 1 h at room temperature. Retained peroxidase activity was revealed with 100 μL
183 per well of a 2 mg mL^{-1} *o*-phenyldiamine freshly prepared solution in 25 mM citrate and 62
184 mM phosphate buffer, pH 5.4, containing 0.012% H_2O_2 . Signal generation was stopped after 10
185 min at room temperature with 100 μL per well of 1 M H_2SO_4 . Absorbance was immediately
186 read at 492 nm.

187 Conjugate-coated indirect cELISAs were performed with plates coated by overnight
188 incubation at room temperature with 100 μL per well of assay conjugate solution in 50 mM
189 carbonate–bicarbonate buffer, pH 9.6. Between incubation steps, plates were washed as
190 described above. The competitive reaction was done during 1 h at room temperature with 50
191 μL per well of cyprodinil standard solution in PBS plus 50 μL per well of monoclonal antibody
192 solution in PBST. The immunochemical reaction was detected by incubation at room
193 temperature during 1 h with 100 μL per well of RAM–HRP solution (1/2000 dilution) in PBST.
194 Finally, signals were generated and plates were read as indicated for direct assays.

195 Absorbance values were plotted versus the logarithm of standard concentration and
196 were mathematically fitted to a four-parameter logistic equation using the SigmaPlot software
197 package from Systat Software Inc. (Chicago, IL, USA). Sigmoidal curves were normalized and
198 average parameters were calculated from independent experiments. Assay sensitivity was
199 estimated as the analyte concentration reducing 50% (IC_{50}) the maximum absorbance (A_{max}),
200 and the limit of detection (LOD) was calculated as the IC_{10} value. Percentage cross-reactivity
201 (CR) was assessed from the quotient between the IC_{50} for cyprodinil and the IC_{50} for the
202 evaluated molecule, both in molar concentration units.

203 *2.5. Ionic strength, pH, and solvent tolerance*

204 To evaluate ionic strength and pH influence, the competitive reaction was carried out
205 with cyprodinil standards in Milli-Q water and tracer or antibody solutions, for direct and
206 indirect cELISAs, respectively, in the studied buffer. Twenty millimolar phosphate buffers, pH
207 7.4, were prepared containing 0.05% (v/v) Tween 20 and increasing concentrations of NaCl.
208 Regarding pH, a series of 20 mM phosphate buffers were arranged containing 0.05% (v/v)
209 Tween 20 and the adequate NaCl concentration so to keep the ionic strength constant ($I = 332$
210 mM). Solvent tolerance was assessed with cyprodinil standard curves prepared in Milli-Q
211 water containing increasing amounts of acetonitrile. The competitive reaction mixture was
212 accomplished with enzyme tracer or monoclonal antibody solutions prepared in 20 mM
213 phosphate, pH 7.4, containing 0.05% (v/v) Tween 20 and 280 mM NaCl.

214 *2.6. Sample extraction and analysis*

215 Strawberries were acquired from different local markets along a 1 month period, and
216 residues were extracted following the QuEChERS method (Lehotay, 2007). Briefly, 15 g of
217 sample homogenate was mixed with 1.5 g of sodium acetate and 6 g of anhydrous magnesium
218 sulfate in 15 mL of acetonitrile containing 1% (v/v) acetic acid. Fifty microliters of TPP was
219 added as internal control. The mixture was vigorously stirred and the organic phase was

220 separated by centrifugation at 2200×g during 5 min. Primary/secondary amine (50 mg) was
221 used to cleanup 1 mL of extract in the presence of 150 mg of anhydrous magnesium sulfate by
222 vortexing. After a second centrifugation step, clean extracts were filtrated through a Teflon
223 filter (0.2 μm) and stored at –20 °C in amber glass vials.

224 Extracts of cyprodinil-free strawberry samples, as checked by UPLC, were employed for
225 recovery studies. Moreover, samples were analyzed by the developed cELISAs and by a
226 reference chromatographic method. Deming regression analysis was applied for method
227 comparison using the SigmaPlot software (version 12.0). Immunochemical analysis was carried
228 out by direct and indirect cELISA. Strawberry extracts were diluted with Milli-Q water and
229 enzyme tracer or antibody solutions were prepared in 20 mM phosphate, pH 7.4, containing
230 0.05% (v/v) Tween 20 and 280 mM NaCl. A cyprodinil standard curve was run in each plate
231 under the same conditions. Standards were prepared in borosilicate glass tubes by 4-fold serial
232 dilution in Milli-Q water starting from a 10 μg L⁻¹ cyprodinil solution. UPLC with a triple
233 quadrupole MS detector was selected as reference methodology for cyprodinil sensitive
234 analysis. The injection volume was set at 5 μL and the mobile phase consisted of a mixture
235 between 0.5% (v/v) formic acid in water and acetonitrile. The gradient started at 50% of the
236 acidic aqueous solution, it was linearly increased to 95% in 4 min, and then to 98% in 2 min at
237 a flow rate of 400 μL min⁻¹. The retention times obtained under the aforementioned
238 conditions were 0.7 and 2.2 min for cyprodinil and TPP, respectively. The employed
239 parameters were: ESI+, parent ion 226.2 *m/z*, daughter ions 93.1 and 108.1 *m/z*, 33 eV collision
240 energy and 45 V cone energy for cyprodinil; and ESI+, parent ion 328.3 *m/z*, daughter ions 77.0
241 and 152.4 *m/z*, 30 eV collision energy and 30 V cone energy for TPP.

242 **3. Results and discussion**

243 *3.1. Immunoreagent evaluation*

244 Hapten-to-protein MR values of BSA, OVA, and HRP conjugates with cyprodinil
245 derivatives were determined by MALDI–TOF, and results were compared with those obtained
246 previously by UV determination (Table S1). For all conjugates, similar MRs were retrieved by
247 both methods, except for BSA and OVA conjugates with hapten CD*b*. In these cases, the
248 differential absorbance method seemed to underestimate the hapten density. On the
249 contrary, the UV method slightly overestimated the MRs of most HRP conjugates. In general, it
250 is considered that BSA contains 30 lysine ϵ -amine groups accessible for conjugation
251 (Hermanson, 1996). Thus, the prepared BSA conjugates had about half of the available lysine
252 residues coupled to a hapten molecule as expected, with coupling reaction yields between
253 54% for hapten CD*m* and 67% for haptens CD*b* and CD*n*. In the case of OVA conjugates,
254 coupling yields were slightly higher than those of BSA, reaching 69% for haptens CD*m* and CD*p*
255 and 77% for haptens CD*b* and CD*n*, based on 20 lysine residues in OVA. Regarding HRP
256 conjugates, all of the available ϵ -amine groups seemed to hold a cyprodinil hapten.

257 A battery of mAbs to cyprodinil has been raised in mice by immunization with BSA
258 conjugates of haptens CD*b* and CD*m*. These haptens were selected for immunization
259 considering the good immune responses obtained previously in rabbits [Esteve-Turrillas et al.,
260 2012; Esteve-Turrillas et al., 2015]. Moreover, these synthetic compounds hold the linker at
261 opposite sites of the cyprodinil molecule, thus displaying different immunodeterminant
262 chemical moieties at the distal end for antibody formation. Each novel monoclonal binder was
263 evaluated using conjugates of homologous and heterologous haptens (same or different
264 hapten compared to the hapten of the immunizing conjugate, respectively) by checkerboard
265 direct and indirect cELISA. Not surprisingly, none of the antibodies recognized the
266 corresponding conjugate of hapten CD*n* – holding the linker at a central position of the
267 cyprodinil molecule – in either of the two cELISA formats (results not shown). Regarding the
268 direct format (Table 1), most CD*b*-type antibodies – obtained from mice immunized with a
269 conjugate of hapten CD*b* – only bound the homologous tracer, and only one mAb (CD*b*#33)

270 could bind haptens *CDm* and *CDp* – with the linker at the aniline ring. On the contrary, a
271 significantly high number of *CDm*-type antibodies – generated from mice immunized with a
272 conjugate of hapten *CDm* – bound not only hapten *CDp* – with the linker at a contiguous
273 position in the same ring – but also hapten *Cdb* – holding the spacer at the opposite ring. This
274 result indicates that the methyl moiety at the pyrimidine ring in the immunizing hapten *CDm*
275 has probably increased the probability to find antibodies recognizing the heterologous hapten
276 *Cdb*, which holds the spacer in the same position of the methyl substituent of cyprodinil.
277 Concerning the indirect format (Table S2), most antibodies bound the two heterologous assay
278 haptens, independently of the aromatic ring at which the linker was located, thus showing a
279 more versatile behavior of this cELISA format to spacer tethering site heterology. Regarding
280 affinity, mAbs reaching IC_{50} values to cyprodinil as low as 0.2 nM were obtained. Furthermore,
281 the sensitivity of homologous assays could be improved with most mAbs in both cELISA
282 formats using linker-site heterologous haptens.

283 Antibody specificity was studied by homologous indirect cELISA. The anilinopyrimidine
284 family of fungicides is composed nowadays by three members: cyprodinil, pyrimethanil, and
285 mepanipyrim. Remarkably, *CDm*-type antibodies, derived from a hapten mainly exposing to
286 the immune system the unique moieties in cyprodinil, displayed a higher specificity than
287 antibodies from hapten *Cdb* (Table S3). Moreover, mepanipyrim was better recognized than
288 pyrimethanil by most binders. The most specific mAb was *CDm*#12, with CR values below 5%,
289 and the most generic mAb for this group of pesticides was *Cdb*#31, with CR values around
290 60%. Anyhow, no interferences should occur with pyrimethanil or mepanipyrim in food
291 analysis because cyprodinil is usually employed in combination with other fungicides with a
292 different mode of action, but not with another anilinopyrimidine, in order to avoid fungi
293 resistance and to increase biocide efficiency. Common active ingredients that are formulated
294 together with cyprodinil are picoxystrobin and fludioxonil. Studies with these compounds for

295 all of the generated mAbs did not cause signal inhibition when assayed at 10 μM . As a result,
296 no interferences can be expected from other pesticides when food samples will be analyzed.

297 3.2. Immunoassay characterization

298 Further studies were carried out with particular mAb and assay conjugate pairs
299 affording inhibition curves with A_{max} values higher than 1.0, low IC_{50} values, and moderate
300 slopes. Higher sensitivity was attained in both formats when a slight linker site heterology
301 (from *meta* to *para* position) was introduced. Thus, mAb CDm#31 by direct cELISA and mAb
302 CDm#21 by indirect cELISA together with the corresponding assay conjugate of hapten CDp
303 seemed to be the most suitable immunoreagent combinations (Tables 1 and S2). Assay
304 conditions and analytical features of the optimized immunoassays are listed in Table 2. In both
305 cELISA formats, the inhibition curves showed very low IC_{50} – between 0.1 and 0.2 $\mu\text{g L}^{-1}$ – and
306 LOD values – at or below 0.03 $\mu\text{g L}^{-1}$. Background signals were negligible for the two assays.
307 Moreover, inter and intra-day precision of A_{max} and IC_{50} values were satisfactory.

308 Influence of pH and ionic strength as well as solvent concentration over the two
309 selected immunoassays was examined (Fig. 2). Interestingly, for both assays, acidic pH values
310 decreased the maximum signal of the assay, probably due to protonation of the pyrimidine
311 moiety, whereas no substantial changes in A_{max} were observed at basic pHs. Besides, sensitivity
312 of the two immunoassays remained almost constant at any of the studied pH values.
313 Regarding ionic strength, high tolerance to salt concentration was found with the direct cELISA.
314 On the other hand, the indirect assay was more robust at low than at high ionic strength
315 values. Since the QuEChERS procedure employs acetonitrile for pesticide residue extraction,
316 the effect of that solvent over the developed cELISAs was assessed. With the indirect assay, a
317 decrease in A_{max} was found at low acetonitrile concentration and an increase of IC_{50} value was
318 seen particularly at solvent concentrations higher than 2% (Fig. 2). On the contrary, only a
319 slight decrease of maximum signal and sensitivity was observed with the direct assay.

320 3.3. Strawberry sample analysis

321 In order to evaluate the trueness and precision of the developed cELISAs, fortified
322 cyprodinil-free strawberry extracts were diluted at least 100 times in water. The higher matrix
323 effect over the indirect assay – required dilution was 300 fold – was probably due to a lower
324 tolerance to acetonitrile (Fig. 2). From this study, recoveries between 81% and 114% for the
325 direct assay and between 82% and 110% for the indirect format were obtained when a
326 cyprodinil concentration range from 10 to 5000 $\mu\text{g L}^{-1}$ was evaluated (Table 3). Relative
327 standard deviations (RSD) were mostly below 20%, except with the direct immunoassay at the
328 highest fungicide levels. Considering the assayed lowest fortified level that met the EU method
329 performance acceptability criteria – mean recovery value in the 70–120% range and a RSD \leq
330 20% (DG SANCO, 2007) –, the limit of quantification of both cELISAs could be set at 10 $\mu\text{g L}^{-1}$,
331 500 times lower than the EU and the US MRL for cyprodinil in strawberry samples (5 mg kg^{-1}).

332 3.4. Method assessment

333 Incidence of cyprodinil residues in 50 commercial strawberry samples was determined
334 by the optimized direct and indirect cELISA, and the positive samples were analyzed by UPLC–
335 MS/MS as a reference chromatographic method. No samples with cyprodinil levels above the
336 MRLs were encountered, yet 9 samples contained residue contents that were detectable by
337 the three applied analytical procedures. According to UPLC–MS/MS analysis, the highest
338 cyprodinil concentration was 907 ng mL^{-1} , the lowest was 41 ng mL^{-1} , and the average content
339 of the positive samples was 274 ng mL^{-1} (Table S4). Statistical data analysis was carried out by
340 Deming regression which considers the standard deviation values of positive samples in both
341 methods (Fig. 3). Comparison of the direct cELISA with the chromatographic method afforded
342 a regression with a slope of 0.90 (95% confidence interval was from 0.68 to 1.13) and an
343 intercept of 16 (95% confidence interval was from –25 to 57). The corresponding Deming
344 regression of the indirect cELISA showed a slope of 0.90 (95% confidence interval was from

345 0.78 to 1.03) and an intercept of -8 (95% confidence interval was from -26 to 11), indicating
346 that agreement of both immunoassays with the reference chromatographic method exists.

347 **4. Conclusions**

348 The first reported mAbs to cyprodinil were generated within the present study. High-
349 affinity antibodies were obtained from both of the immunizing haptens; however, the lowest
350 IC₅₀ values and the highest specificity were observed with CDM-type mAbs. Regarding linker
351 site heterology, sensitivity improvements were commonly achieved, though haptens with a
352 contiguous tethering site compared to the homologous hapten were better recognized than
353 those with opposite linker sites. Two immunoassays were selected with IC₅₀ values below 0.2
354 µg L⁻¹, and LODs of 0.02 and 0.03 µg L⁻¹ for the indirect and direct assay, respectively. The
355 direct assay was slightly more robust to physicochemical changes than the indirect format.
356 A_{max} of both cELISAs was sensitive to acidic pH values, whereas high salt concentrations could
357 increase the IC₅₀ value of the indirect assay. In order to extract strawberry samples, the
358 QuEChERS methodology using acetonitrile was applied. Recoveries of cyprodinil-fortified
359 samples were acceptable between 10 and 5000 µg L⁻¹ for both immunoassays at different
360 dilutions. When strawberry samples from local markets were studied, no violative samples
361 (over 5 mg L⁻¹) were found, although 18% contained detectable levels below 1 mg L⁻¹. The
362 results obtained by the developed cELISAs were in good agreement with those observed by
363 UPLC-MS/MS.

364

365 **Acknowledgements**

366 This work was supported by the Spanish Ministerio de Ciencia e Innovación (AGL2009-
367 12940-C02-01-02/ALI and AGL2012-39965-C02-01-02/ALI) and cofinanced by FEDER funds.
368 J.V.M. and F.A.E.-T. were hired by the Consejo Superior de Investigaciones Científicas (CSIC),
369 the former under a Ramón y Cajal contract and the latter under a JAE-doc contract, both of

370 them financed by Ministerio de Ciencia e Innovación and the European Social Fund. G.Q.-R. is a
371 recipient of a PROMEP predoctoral fellowship from Secretaría de Educación Pública of the
372 Mexican government. Proteomic analysis was carried out in the SCIE University of Valencia
373 Proteomics Unit, a member of the ISCIII ProteoRedProteomics Platform.

374

375 **Appendix A. Supplementary material**

376 Supplementary data associated with this article can be found, in the online version, at

377

378 **References**

- 379 Anastassiades, M., Lehotay, S. J., Stajnbaher, D., Schenck, F. J. (2003). Fast and easy
380 multiresidue method employing acetonitrile extraction/partitioning and “dispersive
381 solid-phase extraction” for the determination of pesticide residues in produce. *Journal*
382 *of AOAC International*, *86*, 412–431.
- 383 DG SANCO. (2007). Method validation and quality control procedures for pesticide residues
384 analysis in food and feed. Document No. SANCO/2007/3131.
- 385 DG SANCO Review Report. (2010). Review report for the active substance cyprodinil.
386 Document No. SANCO/10014/2006.
- 387 Esteve-Turrillas, F. A., Agulló, C., Abad-Fuentes, A., Abad-Somovilla, A., & Mercader, J. V.
388 (2012). Immunoreagent generation and competitive assay development for cyprodinil
389 analysis. *Journal of Agricultural and Food Chemistry*, *60*, 4803–4811.
- 390 Esteve-Turrillas, F. A., Mercader, J. V., Agulló, C., Abad-Somovilla, A., & Abad-Fuentes, A.
391 (2015). Moiety and linker site heterologies for highly sensitive immunoanalysis of
392 cyprodinil in fermented alcoholic drinks. *Food Control*, *50*, 393–400.
- 393 European Committee for Standardisation (CEN) EN 15662:2008. (2008). Foods of plant origin –
394 Determination of pesticide residues using GC–MS and/or LC–MS/MS following
395 acetonitrile extraction/partitioning and clean-up by dispersive SPE – QuEChERS-
396 method.
- 397 Fernandes, V. C., Domingues, V. F., Mateus, N., & Delerue-Matos, C. (2012). Pesticide residues
398 in Portuguese strawberries grown in 2009–2010 using integrated pest management
399 and organic farming. *Environmental Science and Pollution Research*, *19*, 4184–4192.
- 400 Fontana, A. R., Rodriguez, I., Ramil, M., Altamirano, J. C., & Cela, R. (2011). Solid-phase
401 extraction followed by liquid chromatography quadrupole time-of-flight tandem mass
402 spectrometry for the selective determination of fungicides in wine samples. *Journal of*
403 *Chromatography A*, *1218*, 2165–2175.

404 Hermanson, G. T. (1996). *Bioconjugate Techniques*. (1st ed.). San Diego: Academic Press, (Part
405 III).

406 Knauf-Beiter, G., Dahmen, H., Heye, U., & Staub, T. (1995). Activity of cyprodinil: Optimal
407 treatment timing and site of action. *Plant Disease*, *79*, 1098–1103.

408 Lehotay, S. J. (2007). Determination of pesticide residues in foods by acetonitrile extraction
409 and partitioning with magnesium sulfate: Collaborative study. *Journal of AOAC*
410 *International*, *90*, 485–520.

411 Masner, P., Muster, P., & Schmid, J. (1994). Possible methionine biosynthesis inhibition by
412 pyrimidinamine fungicides. *Pesticide Science*, *42*, 163–166.

413 Melo, A., Aguiar, A., Mansilha, C., Pinho, O., & Ferreira, I. M. P. L. V. O. (2012). Optimisation of
414 a solid-phase microextraction/HPLC/Diode Array method for multiple pesticide
415 screening in lettuce. *Food Chemistry*, *130*, 1090–1097.

416 Mercader, J. V., Suárez-Pantaleón, C., Agulló, C., Abad-Somovilla, A., & Abad-Fuentes, A.
417 (2008a). Hapten synthesis and monoclonal antibody-based immunoassay development
418 for detection of the fungicide trifloxystrobin. *Journal of Agricultural and Food*
419 *Chemistry*, *56*, 2581–2588.

420 Mercader, J. V., Suárez-Pantaleón, C., Agulló, C., Abad-Somovilla, A., & Abad-Fuentes, A.
421 (2008b). Production and characterization of monoclonal antibodies specific to the
422 strobilurin pesticide pyraclostrobin. *Journal of Agricultural and Food Chemistry*, *56*,
423 7682–7690.

424 Milling, R. J., & Richardson, C. J. (1995). Mode of action of the anilinopyrimidine fungicide
425 pyrimethanil. 2. Effects on enzyme secretion in *Botrytis cinerea*. *Pesticide Science*, *45*,
426 43–48.

427 Moeder, M., Bauer, C., Popp, P., van Pinxteren, M., & Reemtsma, T. (2012). Determination of
428 pesticide residues in wine by membrane-assisted solvent extraction and high-

429 performance liquid chromatography-tandem mass spectrometry. *Analytical and*
430 *Bioanalytical Chemistry*, 403, 1731–1741.

431 Rial-Otero, R., Yagüe-Ruiz, C., Cancho-Grande, B., & Simal-Gándara, J. (2002). Solid-phase
432 microextraction-gas chromatographic-mass spectrometric method for the
433 determination of the fungicides cyprodinil and fludioxonil in white wines. *Journal of*
434 *Chromatography A*, 942, 41–52.

435 Rosslenbroich, H. J. & Stuebler, D. (2000). *Botrytis cinerea* – history of chemical control and
436 novel fungicides for its management. *Crop Protection*, 19, 557–561.

437 Scientific Report of EFSA. (2013). The 2010 European Union report on pesticide residues in
438 food. European Food Safety Authority. *EFSA Journal*, 11, 3130.

439 Suárez-Pantaleón, C., Mercader, J. V., Agulló, C., Abad-Somovilla, A., & Abad-Fuentes, A.
440 (2008). Production and characterization of monoclonal and polyclonal antibodies to
441 forchlorfenuron. *Journal of Agricultural and Food Chemistry*, 56, 11122–11131.

442 USDA, Pesticide Data Program. (2014). Annual summary, calendar year 2012.

443 Vaquero-Fernández, L., Saenz-Hernández, A., Sanz-Asensio, J., Fernández-Zurbano, P., Sainz-
444 Ramírez, M., Pons-Jubera, B., López-Alonso, M., Epifanio-Fernández, S. I., & Martínez-
445 Soria, M. T. (2008). Determination of cyprodinil and fludioxonil in the fermentative
446 process of must by high-performance liquid chromatography-diode array detection.
447 *Journal of the Science of Food and Agriculture*, 88, 1943–1948.

448 Zhang, K., Wong, J. W., Hayward, D. G., Sheladia, P., Krynitsky, A. J., Schenk, F. J., Webster, M.
449 G., Ammann, J. A., & Ebeler, S. E. (2009). Multiresidue pesticide analysis of wines by
450 dispersive solid-phase extraction and ultrahigh-performance liquid chromatography-
451 tandem mass spectrometry. *Journal of Agricultural and Food Chemistry*, 57,
452 4019–4029.

453

454


455 **Figure legends**

456 **Fig. 1.** Chemical structure of cyprodinil haptens.

457 **Fig. 2.** Influence of pH, ionic strength, and acetonitrile over A_{\max} and IC_{50} values of a direct
458 cELISA with mAb CDM#31 and an indirect cELISA with mAb CDM#21.

459 **Fig. 3.** Deming regression between the developed cELISAs and UPLC-MS/MS with a 95%
460 confidence interval.

461




462

463

464 Figure 1

465




467

468

469 Figure 2

470

471



472

473 Figure 3

474

475

476

477

Table 1
Checkerboard assays by direct cELISA.

mAb	HRP-CD _b					HRP-CD _m					HRP-CD _p				
	[mAb] ^a	[HRP] ^b	A _{max}	Slope	IC ₅₀ ^c	[mAb]	[HRP]	A _{max}	Slope	IC ₅₀	[mAb]	[HRP]	A _{max}	Slope	IC ₅₀
CD _b #21	1000	100	- ^d	-	-	1000	100	-	-	-	1000	100	-	-	-
CD _b #31	1000	30	0.85	1.04	0.4	1000	100	-	-	-	1000	100	-	-	-
CD _b #32	1000	10	0.85	0.88	2.0	1000	100	-	-	-	1000	100	-	-	-
CD _b #33	1000	10	0.96	0.88	3.1	1000	10	1.31	0.87	3.4	1000	3	0.75	0.80	3.0
CD _m #11	1000	100	-	-	-	1000	100	1.02	0.98	4.5	1000	100	-	-	-
CD _m #12	1000	100	-	-	-	1000	100	1.17	1.30	5.7	1000	100	-	-	-
CD _m #13	1000	100	-	-	-	1000	10	1.01	1.14	1.9	1000	100	1.70	1.04	1.9
CD _m #14	1000	100	-	-	-	1000	100	0.85	0.96	2.1	1000	100	-	-	-
CD _m #21	1000	100	-	-	-	1000	100	-	-	-	1000	100	-	-	-
CD _m #22	1000	3	0.87	1.09	1.5	1000	3	1.24	1.06	6.4	1000	3	1.28	1.06	6.4
CD _m #23	1000	100	-	-	-	1000	30	1.05	0.94	1.2	1000	100	1.31	1.10	1.2
CD _m #31	1000	10	0.93	1.49	0.2	1000	3	1.07	1.61	0.7	1000	3	1.21	1.02	0.3
CD _m #32	1000	10	0.72	1.24	0.3	1000	3	1.40	1.10	0.8	1000	3	1.29	0.91	0.8
CD _m #33	1000	10	1.00	1.08	0.3	1000	3	1.16	1.24	0.7	1000	3	1.22	0.91	0.3
CD _m #34	1000	30	1.15	1.08	0.4	1000	3	1.20	0.92	1.3	1000	3	1.28	0.72	1.8
CD _m #35	1000	10	0.76	1.23	0.3	1000	3	1.31	1.13	0.9	1000	3	1.27	0.95	1.0
CD _m #36	1000	10	0.75	0.75	0.4	1000	3	1.03	0.98	1.1	1000	3	1.04	0.73	1.5

^a Antibody concentration in ng mL⁻¹.

^b Tracer concentration in ng mL⁻¹.

^c Values are in nM.

^d A_{max} values were lower than 0.7.


478

479

480

481

Table 2 482
Optimized assay conditions and analytical parameters of
cyprodinil standard curves. 183



	Immunoassay		
	Direct	Indirect	
<i>Immunoreagents</i>			489
mAb	CDm#31	CDm#21	
[mAb] (µg L ⁻¹)	500	100	490
Conjugate	HRP-CDp	OVA-CDp	
[Conjugate] (µg L ⁻¹)	2.5	100	491
<i>Curve parameters</i>			
A _{max}	1.23 ± 0.07	1.34 ± 0.16	
A _{min}	0.02 ± 0.01	0.03 ± 0.01	492
Slope	-1.25 ± 0.21	-1.12 ± 0.05	
IC ₅₀ (µg L ⁻¹)	0.17 ± 0.02	0.14 ± 0.01	493
R ²	0.998	0.999	
LOD (µg L ⁻¹)	0.03	0.02	
<i>Intra-day precision (% , n=5)</i>			494
A _{max}	6	3	
IC ₅₀	12	9	495
<i>Inter-day precision (% , n=5)</i>			
A _{max}	6	14	
IC ₅₀	10	20	496

497

498

499

500

Table 3
Recoveries from spiked strawberry samples
measured by the two developed monoclonal
cELISAs.^a

Spiked ($\mu\text{g L}^{-1}$)	Direct		Indirect	
	Recovery (%)	RSD (%)	Recovery (%)	RSD (%)
10	98	12	95	12
50	81	6	85	4
100	91	8	90	4
500	113	14	82	10
1000 ^b	114 ^b	22	92	5
5000 ^b	112 ^b	26	110	19





^a Results are the average of 4 independent measurements.
Extracts were diluted 100 fold for the direct cELISA and 300
fold for the indirect assay.

^b Extracts were diluted 1000 fold for both assays.

501


Figure(s)



[Click here to download high resolution image](#)



Figure(s)

[Click here to download high resolution image](#)





Supplementary Material

[Click here to download Supplementary Material: Supplementary Material_revised.docx](#)