



provided by Repository Open Access to Scientific Information.
24-29 OCTOBER, 2010

MAR DEL PLATA, ARGENTINA

Poster-S3g-06

ACID DIGESTION FOR CADMIUM DETERMINATION IN MINERAL SALT MIXING FOR ANIMAL NUTRITION BY THERMOSPRAY FLAME FURNACE ATOMIC ABSORPTION SPECTROMETRY

Marcelo Braga Bueno Guerra ¹, Ana Rita Araujo Nogueira ², <u>Patricia Helena Toniolo Silva</u> ³, <u>Edenir Rodrigues Pereira Filho</u>

*anarita@cppse.embrapa.br *patytoniolo@yahoo.com.br

¹UFSCar, Deparment of Chemistry, Rod. Washington Luiz, km 235. São Carlos-SP, Brasil
 ²EMBRAPA, CPPSE, Rodovia Washington Luiz, km 234, São Carlos, SP, Brasil
 ³Universidade Federal de São Carlos, Departamento de Química, Rod. Washington Luiz, km 235; P.O. Box 676, Brasil

Cadmium was determined in a sample of mineral salt mixing to milk cattle that has been preparing to be a reference material. A volume of 3 mL of HNO3 3.0 mol L-1 plus 2.0 of H2O2 30% (w w-1) was added to 250 mg of sample in a cavity microwave oven. After digestion performance the volume was adjusted to 10 mL with water, centrifuged and the supernatant determined by thermospray flame furnace atomic absorption spectrometry (TS-FF-AAS). The system consist of a peristaltic pump, a lab-made commutator injector, and the thermospray flame furnace unit, which includes the ceramic capillary (0.5 mm ID, 2.0 mm ED and 100 mm of length) and the Ni tube (10 cm) which is located on the AAS burner. For reaching a higher temperature inside the tube, six holes of 2 mm diameter were drilled in the bottom of the tube. Another orifice was drilled at 90° to the bottom holes for inserting the thermospray capillary. The capillary tip is about 1 mm inside the tube furnace. The manifold was assembled with 0.5 mm i.d. PTFE tubing, sample volume of 100 µL, 0.5 mL min-1 of flow rate, and 0.5 % (w v-1) EDTA plus 0.1 % (v v-1) Triton X 100 as carrier stream. Limits of detection (LOD) and quantification (LOQ) obtained were 0.31µg L-1 and 1.02µg L-1, respectively. The obtained results were 402 ± 35µg L-1. The accuracy was checked by performing GFAAS (465 ± 95µg L-1) and the results proved to be similar by using t test. The high rsd in the GFAAS was probably due the high chloride sample amount.