

Edinburgh Research Explorer

A Window in the Brain

Citation for published version:

Abdullateef, S, Jordan, B, McLellan, A, Nenadovic, V, Escudero, J & Lo, M 2022, A Window in the Brain: Applying data science to quantitatively detect seizures with minimal-density EEG montage. in 14th European Paediatric Neurology Society Congress., EPNS21-406, 14th European Paediatric Neurology Society Congress, Glasgow, United Kingdom, 28/04/22.

Link:

Link to publication record in Edinburgh Research Explorer

Document Version:

Publisher's PDF, also known as Version of record

Published In:

14th European Paediatric Neurology Society Congress

General rights

Copyright for the publications made accessible via the Edinburgh Research Explorer is retained by the author(s) and / or other copyright owners and it is a condition of accessing these publications that users recognise and abide by the legal requirements associated with these rights.

Take down policy
The University of Edinburgh has made every reasonable effort to ensure that Edinburgh Research Explorer content complies with UK legislation. If you believe that the public display of this file breaches copyright please contact openaccess@ed.ac.uk providing details, and we will remove access to the work immediately and investigate your claim.





A Window in the Brain: Applying data science to quantitatively detect seizures with minimal-density EEG montage

Shima Abdullateef ¹, Brian Jordan², Ailsa McLellan², Vera Nenadovic³, Javier Escudero¹, Tsz-Yan Lo^{2,4}





¹ School of Engineering, IDCOM, University of Edinburgh (Edinburgh), ² Royal Hospital for Children & Young Person (Edinburgh), ³ BrainsView (Ontario), ⁴ Usher Institute, University of Edinburgh, NINE BioQuarter (Edinburgh)


TAKE-HOME MESSAGE

Epileptic seizures can be detected with as few as **four** EEG channels with **above 70** % accuracy.



1. BACKGROUND

- Up to 74% of seizure in PCCU do not have clear clinical marker.
- Gold-standard seizure detection using multi-channels EEG needs:
 - Highly specialised neurologist
 - Trained clinical physiologists
 - Round-the clock service
- PCCU team can reliably apply a 4-channels montage shown here.




2. METHODS

 Forty fully anonymised routinely collected EEG provided by the Royal Hospital for Children and Young People


Age range: 1 -12 years old

Total duration: 38.5

No. Seizures: 236



■ Seizure Onset Duration
■ Background EEG




3. RESULTS


Both features capture the number of connections between the pairs of channels.

 Using two features improved the accuracy for both 8 and 4-channels seizure detection.

Seizure detection using Phase Synchrony only



Seizure detection using phase Synchrony and Coherence



5. FUTURE WORK

- Multi -centre recruitment of data
- Improving the accuracy and reducing the missrate
- Developing a user-friendly and accurate seizure detection tool for frontline clinicals