




# Developing Framework for Flexible and Robust Real-Time Release Testing

Natalie Maclean<sup>1,2</sup>, Blair Johnston<sup>1,2</sup>, James D. Litster<sup>3</sup> & Daniel Markl<sup>1,2</sup>

<sup>1</sup>Strathclyde Institute of Pharmacy & Biomedical Sciences, University of Strathclyde, Glasgow, UK  
<sup>2</sup>Centre for Continuous Manufacturing and Advanced Crystallisation (CMAC), University of Strathclyde, Glasgow, UK.  
<sup>3</sup>Department of Chemical and Biological Engineering, University of Sheffield, Sheffield, UK.

## Right First-Time Manufacturing of Pharmaceuticals (RiFTMaP)



Increased opportunity to implement RTTR in manufacturing lines due to advances in process analytical technology and digital manufacturing. [1]

- A three-year, £1.5m EPSRC collaboration between the University of Sheffield, University College London, University of Strathclyde and Purdue University.
- Developing right-first-time smart manufacturing systems incorporating Industry 4.0 concepts within a systematic framework for smart continuous manufacturing.

Development of a framework and computational tools for optimal design of pharmaceutical processes with real-time process management system and a flexible real-time release framework, all verified at pilot scale.




Fig 2. RiFTMaP WP4: Developing flexible and robust RTTR framework

## Performance-Controlling Mechanisms of Tablets





Fig 3. Workflow for the classification of the performance-controlling mechanisms of immediate-release tablets. [2]

- Mechanisms determined by raw material properties and the process conditions
- Allows the identification of critical formulation and manufacturing parameters which influence the performance of a specific formulation
- Use this knowledge to identify the appropriate sensors for RTTR

## Identifying the Critical Formulation & Manufacturing Properties



- Wettability controlled  
Disintegration time sensitive to changes in tensile strength & porosity
- Dissolution controlled  
Disintegration time sensitive to changes in contact angle (wettability)
- Swelling controlled  
Disintegration time is not strongly affected by tensile strength, porosity or contact angle.

## Predicting Performance

Raw Material Properties Performance of Tablets

- Surface energy
  - Particle size & shape
  - Moisture sorption
  - Solubility & intrinsic dissolution rate
  - True density & intra-particle porosity
  - Surface area
- Compare and model dissolution profiles for a range of different API, formulations & manufacturing settings

Development of data driven models to predict product performance

## Acknowledgements

This work is part of a collaborative EPSRC-funded project which focuses on the Right-First Time Manufacturing of Pharmaceuticals (RiFTMaP). The authors would like to thank Dr James Mann, Mr Alexander Abbott, Ms Helen Williams, Ms Heather Mead and Dr Ibrahim Khadra for their contributions. The work depicted in Figures 3 and 4 was funded by EPSRC and AstraZeneca.

## References

1. Markl, D. et al. (2020) Review of real-time release testing of pharmaceutical tablets: State-of-the-art, challenges, and future perspective. *Int J Pharm*, 582, 119353.
2. Maclean, N. et al. (2021) Exploring the performance-controlling tablet disintegration mechanisms for direct compression formulations. *Int J Pharm*, 599, 120221.
3. Maclean, N. et al. (2022) Investigating the role of excipients on the physical stability of directly compressed tablets. *Int J Pharm: X*, 4, 100106.

MCC - microcrystalline cellulose; DCPA - dibasic calcium phosphate anhydrous; CCS - croscarmellose sodium; XPVP - crospovidone; L-HPC - low-substituted hydroxypropyl cellulose; SSG - sodium starch glycolate;  $\theta_{c0}$  - initial contact angle.

Scan for more on this project

