



L-Aspartate oxidase from *Sulfolobus tokodaii*: immobilization studies

Andrea Fiorati,^a Chiara Allegretti,^a Lorenzo Cerioli,^a Luciano Piubelli,^{b,c} Loredano Pollegioni,^{b,c} Elena Rosini,^{b,c} Davide Tessaro,^{a,b} Mattia Valentino,^{b,d} Paola D'Arrigo^{a,b}

^a Politecnico di Milano, Department of Chemistry, Materials, and Chemical Engineering "Giulio Natta", Piazza Leonardo da Vinci 32, 20133 Milano, Italy

^b The Protein Factory, Centro Interuniversitario di Biotecnologie Proteiche, Politecnico di Milano, Università degli Studi dell'Insubria e ICNR-CNR, Via Mancinelli 7, 20131 Milano, Italy

^c Dipartimento di Biotecnologie e Scienze della Vita, Università degli Studi dell'Insubria, Via Dunant 3, 21100 Varese, Italy

^d ICRM, CNR Via Mancinelli 7, 20131 Milano, Italy


Introduction

L-Amino acid oxidases (LAOs, EC 1.4.3.2) catalyze the stereospecific oxidative deamination of L-amino acids to α -keto acids, ammonia and hydrogen peroxide.^{1,2} In our labs L-aspartate oxidase from the thermophilic archaea *Sulfolobus tokodaii* (StLASPO, EC 1.4.3.16) has been efficiently produced as recombinant protein in *E. coli* in the active form as holoenzyme³ and fully characterized.

The enzyme is active only on L-aspartate and on L-asparagine and it's absolutely stereoselective. StLASPO shows a remarkable stability: up to 80 °C in the 7 – 10 pH range. Plotting the activity (on L-aspartate) against temperature in the 25 – 95 °C range, an increase with no evidence for any plateau or decrease up to 80 °C was apparent. The activity of the enzyme as a function of pH shows a bell-shaped behavior: it is maximal at about pH 10.

Kinetic resolution of D,L-aspartate

In a previous study StLASPO was used for the resolution of a 50 mM solution of D,L-aspartate, and the full resolution was reached in 4 hours in the 37-70 °C and pH 10-11 range using 0.3 U StLASPO/mL.



The kinetic resolution of D,L-aspartate was now used for evaluating the performance of the immobilized enzyme and was carried out by adding 0.5 mL of a 50 mM D,L-aspartic solution (in water, adjusted to pH 10.0), 200 μ L of water and 2 μ L of catalase to the StLASPO immobilized on resin. This reaction was performed at 70 °C on a thermomixer set to 600 rpm.

For CLEA preparation the assay was performed by adding 0.5 mL of a 50 mM D,L-aspartic solution in "magic buffer" (160 mM Tris-HCl, 160 mM Na₂CO₃, 160 mM H₃PO₄, 0.65 mM potassium chloride and 1% glycerol) padjusted to the desired pH (8, 9 or 10) and 1 μ L of catalase to the CLEA-StLASPO preparation (1.12 U); the reaction was performed at 70 °C on a thermomixer set to 600 rpm.

Immobilization on Relizyme™ HA403/SR

Cycle number	Time (h)	e.e. (%)
1	1	100
2	1	100
3	1	100
4	1	87
	2	100
5	1	77
	2	100

The resin Relizyme™ HA 403/SR (250 mg), which has amino moieties, was treated with 10 mL of 0.125% glutaraldehyde solution in water for 2.5 h on a rotating device at 18 °C. The glutaraldehyde solution was then removed and the resin was washed with 0.1 M phosphate buffer at pH 7.5.

After that, 1 mL of StLASPO (2 U) was added to the resin, with 1 mL of 0.1 M phosphate buffer at pH 7.5 and the mixture was incubated on a rotatory shaker for 24 h at 18 °C.

StLASPO was immobilized with a full yield in terms of activity.

The resin was then washed with 50 mM phosphate buffer at pH 7.5, and stored at 4 °C in the same conditions.

The immobilized enzyme showed an activity of 8 U/g matrix

The influence of D,L-aspartate concentration on enzymatic (2 U/118 mg resin) resolution was also evaluated, by adding 0.5 mL of a D,L-aspartic solution at different concentration (from 50 to 500 mM) at pH 10.0.

Final D,L- Asp (mM)	L-aspartate conversion (%)									
	Time (h)									
	0.5	1	2	3	4	5	7	21	42	
35	100	-	-	-	-	-	-	-	-	
70	62	100	-	-	-	-	-	-	-	
110	12	36	100	-	-	-	-	-	-	
140	-	-	-	67	100	-	-	-	-	
180	-	-	-	-	53	69	84	-	-	
360	-	-	-	-	17	21	21	21	21	

Analytical

At predetermined intervals, 25 μ L aliquots were withdrawn from the reaction mixture and diluted in 100 μ L of distilled water. 10 μ L of this solution were derivatized with 25 μ L of OPA-NAC reagent, diluted with 25 μ L of 50 mM eluent and analyzed by HPLC chromatography (column Gemini 5 μ C18 Phenomenex; eluent: 50 mM sodium acetate buffer pH 5.2 / MeOH 90/10; flux of 0.8 mL/min; detection at 340 nm).

OPA-NAC reagent is a basic aqueous-methanolic solution of ortho-phthalaldehyde and N-acetylcysteine which are able to react with free amine giving a fluorescent derivate.

Other supports tested

Also other solid supports were evaluated:

Sepabeads EC-EP/A, Relizyme™ EP113, Relizyme™ HA403, Relizyme™ HA113, Eupergit® C, IDA-3 Sepharose® yielded poor results in term of immobilization except for IDA-3 resin: this matrix showed a high immobilization yield but the enzyme was leached out due to the non-covalent interaction with the matrix.

Immobilization on resin SEPABEADS EC-EP/S

Cycle number	Time (h)	e.e. (%)
1	3	68
	4	100
2	2	66
	4	100
3	2	65
	4	100
4	2	59
	4	81
	6	100
5	2	28
	4	45
	6	73
	8	80
	10	100

One unit of pure StLASPO was added to 60 mg of Sepabeads EC-EP/S, with oxirane moieties, and to 0.3 mL of 1.25 M KH₂PO₄ buffer at pH 8.0.

The mixture was incubated on a rotatory shaker for 18 h at 25 °C: approximately 50% of the enzymatic activity was present in solution suggesting that half of the starting enzyme has been immobilized.

This immobilized StLASPO preparation yielded full conversion of L-isomer of aspartate in 4 hours for 3 cycles. For the following cycles, the full resolution was obtained in a longer time, requiring up to 10 hours in the fifth cycle.

Immobilization as CLEA

The precipitation of StLASPO as a CLEA⁴ was conducted by adding to 0.5 mL of enzyme solution (2.24 U/mL), 4.5 mL of a 60% w/v (NH₄)₂SO₄ solution, in 0.1 M phosphate buffer at pH 7.5 and 20 μ L of glutaraldehyde.

The mixture was incubated on a rotatory shaker for 2.5 h at room temperature: no residual enzymatic activity was detected in the solution. CLEA-StLASPO was then washed twice with 0.1 M phosphate buffer at pH 7.5 and was stored at 4 °C in 0.05 M phosphate buffer at pH 7.5.

Before each cycle, the CLEA-StLASPO preparation was incubated with 1 μ L of FAD for 30 min in 0.1 M phosphate buffer at pH 7.5.

In the first cycle, the conversion was fully reached at all pH values tested, in a comparatively shorter time at highest pH. On the contrary, starting from the fourth cycle, the full oxidation of L-aspartate was obtained in a shorter time at pH 8 and 9.

Cycle number	Reaction time (h)	e.e. (%) ^a pH 8	e.e. (%) ^a pH 9	e.e. (%) ^a pH 10
1	2	49	32	53
	4	87	100	100
	6	100	-	-
2	2	49	63	66
	4	97	100	100
	6	100	-	-
3	2	69	69	64
	4	73	100	100
	6	100	-	-
4	2	88	54	42
	4	100	100	68
	6	-	-	100
5	2	91	91	40
	4	100	100	69
	6	-	-	95
	7	-	-	100
6	2	0	0	31
	4	0	0	56
	6	0	0	80
	8	0	0	100

Conclusions

Due to the very interesting properties of this enzyme, we propose the immobilized StLASPO as an attractive tool for biotechnological applications. We recognize that CLEA technique or Relizyme™ HA403/SR as solid support represent the best conditions for StLASPO immobilization: this step will allow to improve the stability as well as enzyme reusability.

[1] P. Macheroux, O. Seth, C. Bollschweiler, M. Schwarz, M. Kurfurst, LC Au, S. Ghisla, Eur J Biochem, (2001) 268, 1679-1686.

[2] L. Pollegioni, P. Motta, G. Molla, Appl. Microbiol. Biotechnol. (2013), Vol. 97, Issue 21, 9323-9341.

[3] D. Bifulco, L. Pollegioni, D. Tessaro, S. Servi, G. Molla, Appl Microbiol Biotechnol (2013), Vol. 97, Issue 16, 1-11.

[4] R. Schoevaart, M.W. Wolbers, M. Golubovic, M. Ottens, A.P.G. Kieboom, F. van Rantwijk, L.A.M. van der Wielen, R.A. Sheldon, Biotech. and bioeng. (2004) VOL. 87, NO. 6, 754 - 762