

A practical approach to render tuberculosis samples safe for application of tuberculosis molecular bacterial load assay in clinical settings without a biosafety level 3 laboratory[☆]

Bariki Mtafya^{a,b,*}, Paschal Qwaray^a, Joseph John^a, Emanuel Sichone^a, Alice Shoo^a, Stephen H. Gillespie^b, Nyanda Elias Ntinginya^{a,1}, Wilber Sabiiti^{b,1}

^a National Institute for Medical Research-Mbeya Medical Research Centre, Mbeya, Tanzania

^b University of St Andrews, School of Medicine, St Andrews, UK

ABSTRACT

Background: *Mycobacterium tuberculosis* is a category B infectious pathogen requiring level-3-containment laboratories for handling. We assessed the efficacy of heat and Guanidine thiocyanate (GTC) to inactivate *M. tuberculosis* prior to performance of tuberculosis Molecular Bacterial Load Assay (TB-MBLA).

Method: We performed *in vitro* experiments using *M.tb*, H37Rv reference strain and replicated in sputum specimens. A 0.5 MacFarland standard of *M. tuberculosis* was serially diluted to 1×10^1 CFU/mL and pooled sputum was homogenised prior to serial dilutions and Xpert MTB/RIF Ultra. Three replicates for each containing 1 mL for *M. tuberculosis* and sputum were inactivated at 80 °C for 20 min and with GTC for 15 min. Inactivated samples were processed for culture and TB-MBLA.

Results: No *M. tuberculosis* growth was observed in MGIT for GTC or heat treated H37Rv cultures. All untreated H37Rv dilutions were MGIT positive except the most diluted specimens. Heat and GTC treatment of H37Rv reduced TB-MBLA load by $2.1 \log_{10}$ ($P = 0.7$) and $1.8 \log_{10}$ ($P = 0.7$) respectively, compared to controls. In contrast, heat treated sputum had TB-MBLA bacterial load of $3.47 \pm 3.53 \log_{10}$ compared to $5.4 \pm 3.1 \log_{10}$ eCFU/mL for GTC ($p = 0.57$). All heat and GTC treated sputum were culture negative.

Conclusion: Heat or GTC renders *M. tuberculosis* non-viable and eliminates the need for BSL3 laboratory for performing TB-MBLA in routine healthcare settings.

1. Background

Tuberculosis (TB) is a global public health threat caused by a hazard group 3 infectious agent, *Mycobacterium tuberculosis* (*M.tb*) [1,2]. About 2 billion people are estimated to have latent *M.tb* infections of whom 10 million develop active TB and 1.5 million die each year [3]. Tanzania is among the 30 high burden countries that together contribute 86% of the total TB cases globally [4]. Despite the presence of effective diagnostic tools, still more than 40% of TB cases are treated clinically and effective tests to assess treatment response are not available [4,5]. Thus, the scale up of novel, effective tests that could monitor TB therapy in routine healthcare settings, obviate the need for complex and expensive biosafety level 3 laboratory infrastructure would increase and equalize access.

It has been recognized that TB laboratory staff have a higher risk of TB infection compared to staff in other healthcare occupations [6,7]. Laboratory acquired infections arise during the processing of specimens,

particularly sputum which is the main diagnostic specimen for active TB [8]. Procedures involving manipulation of TB specimens such as culture and extraction of nucleic acids that may generate aerosols are the main source of TB infections to laboratory staff [9]. The World Health Organization (WHO) recommends such procedures be carried in a biosafety level 3 (BSL3) laboratory infrastructure equipped with Class II biological safety cabinets and a negative pressure system. Such infrastructure are not available in most healthcare facilities in resource limited settings with high burden of TB and HIV co-infections, and need these diagnostics the most [10].

Microscopy and culture are the current reference standard for diagnosis and monitoring of TB treatment response [11,12]. Microscopy based tests are well established in most peripheral healthcare facilities while culture tests are available in few zonal and referral laboratories with BSL3 laboratory infrastructure. Unlike culture and standard smear microscopy, molecular tests such as Line Probe Assay (LPA) and GeneXpert MTB/RIF Assays are rapid and sensitive for diagnosis of TB and

[☆] On behalf of the St Andrews Africa Health Research (StAAHR) Network

* Corresponding author. NIMR-Mbeya Medical Research Centre, P.O. BOX 2410, Mbeya, Tanzania.

E-mail address: bmtafya@nimr-mmrc.org (B. Mtafya).

¹ Contributed equally.

accelerate rapid initiation of TB therapy [13–15]. The Xpert MTB/RIF Assay uses a closed cartridge system with chaotropic buffers and can be performed at a primary healthcare facility in absence of BSL3 laboratory infrastructure. Similar to culture, LPA requires manual extraction of DNA and therefore the assay is limited in referral or zonal laboratories equipped with infrastructure for manual extraction of nucleic acid from patient samples [16,17].

Lack of BSL3 laboratory infrastructure is the major bottleneck for the scale up of novel TB tests such as the TB molecular bacterial load assay (TB-MBLA), LPA and standard culture in routine healthcare settings. The WHO End TB strategy and Stop TB Partnership's Global Plan to End TB have set ambitious target to reduce TB incidences by 90% in 2035 [18, 19]. To achieve such milestone it will require, among other tools, rapid implementation of accurate TB diagnostics and monitoring tests in healthcare settings, ensuring their availability to the widest population of TB patients [20]. To this end, we investigated the effectiveness of Guanidine Thiocyanate (GTC) and heat to inactivate viable *M. tb* prior to RNA extraction and explore whether these innovative methods may be applied in routine healthcare laboratories for TB-MBLA test.

2. Methods

2.1. Ethical statements

This work was nested to the TB-MBLA translational study in routine healthcare practice conducted at Rungwe district hospital in Mbeya, Tanzania. The TB-MBLA study was approved by the Mbeya Medical Research and Ethics Committee (SCEC-2439/R. E/V.1/82) and Medical Ethics and Research Coordinating Committee of the National Institute for Medical Research (NIMR/HQ/R.8a/Vol.IX/3687). Additional external approval was granted by the University Teaching and Ethics Committee of the University of St Andrews, United Kingdom (MD15364). All participants provided a written consent or witnessed verbal consent for those who could not write or read in regard to the use

of samples for research purposes including validation of new techniques or assays.

2.2. Study samples

In this study, we used *M. tb*, H37Rv laboratory strain for *in vitro* experiments and replicated the experiments in real pulmonary TB sputum samples. Using these samples, we tested the effectiveness of heat and GTC for inactivation of TB bacilli as part of biosafety risk assessment for performing TB-MBLA in routine healthcare laboratory without BSL3 laboratory infrastructure.

2.3. *In vitro* experiments using *M. tb*, H37Rv laboratory strain

In vitro experiments were performed using the 14 days old culture of *M. tb* H37Rv strain. A 0.5 MacFarland standard equivalent to an approximate 1.0×10^8 CFU/mL of *M. tb* was prepared from pure culture propagated in Lowenstein Jensen media (LJ) and serially diluted in Middlebrook 7H9 media to 1.0×10^1 CFU/mL (Fig. 1). Three replicates containing 1 mL of each dilution were inactivated by boiling at 80°C on a water bath for 20 min and other batch of three aliquots were mixed with 4 mL of GTC (Sigma) containing 1% of β -mercaptoethanol (Sigma) at room temperature for 15 min. Thereafter, GTC treated cultures in 15 mL centrifuge tubes were centrifuged at 3000 g for 30 min whereas heat inactivated cultures in 1.5 mL tubes at 20,000 g for 20 min. The resulting cell pellet from each sample was re-suspended with 1 mL of phosphate buffer solution (PBS, pH 6.8, Sigma). To verify complete inactivation, 500 μL suspension per sample was inoculated into Mycobacterial Growth Indicator tube (MGIT), liquid culture supplemented with Oleic acid, Albumin, Dextrose and Catalase (OADC, Sigma) in absence of antibiotics and incubated in BD BACTEC MGIT 960 Systems following the manufacture instructions [21]. The remaining 500 μL of the suspension was mixed with 950 μL of lysis buffer (RNA Pro solution, MP biomedical) in a 2.0 mL lysing matrix (MP biomedical) and RNA extracted




Fig. 1. Flow of laboratory experiments for H37Rv, *M. tb* reference strain and real pulmonary TB sputum. A fourteen-day old culture of *M. tb* was serially diluted to obtain seven 10-fold dilutions of undiluted culture (10^0) to 1×10^{-8} CFU/ml, and real patients' sputum was diluted to obtain five 10-fold dilutions (fresh un-diluted neat sputum to 1×10^{-5}). Three technical replicates of each of the H37Rv dilutions or sputum was treated with heat by boiling at 80°C and GTC. The H37Rv, *M. tb* controls were not treated with either heat or GTC. For real pulmonary TB sputum, experimental sputum controls were treated with the standard NALC-NaOH procedure for culture and TB-MBLA.

Table 1

Xpert MTB/RIF Ultra results from pooled sputum compared to TB-MBLA and MGIT. Xpert MTB/RIF Ultra Assay detected low bacillary load specimens which were negative by both TB-MBLA and culture. Three technical replicates of each of the sputum aliquot (fresh un-diluted neat sputum to the 10⁻⁵ diluted sputum) were performed. Data presented as average CFU/mL with standard deviation.

| | Xpert MTB/RIF Ultra results | | TB-MBLA | | MGIT-TTP (days: hrs) |
|------------------------------|-----------------------------|-----------------|-----------------|-------------------------------|----------------------|
| | Semi-quantitative results | Ct values (±SD) | Ct values (±SD) | Bacterial load (Mean ± SD/mL) | (±SD) |
| Undiluted (10 ⁰) | Mtb-High | 16 ± 0.1 | 18.37 ± 0.21 | 2.26E+0± 1.41E+03 | 4.22 ± 0.001 |
| 10 ⁻¹ | Mtb-High | 16.3 ± 0.1 | 23.73 ± 0.5 | 6.36E+03± 1.95E+03 | 6.37 ± 0.014 |
| 10 ⁻² | Mtb-Medium | 16.5 ± 0.1 | 27.57 ± 0.85 | 4.23E+02± 8.51E+01 | 10.31 ± 0.018 |
| 10 ⁻³ | Mtb-Low | 18.2 ± 0.1 | 29.89 ± 0.67 | 9.34E+01± 2.81E+01 | 19.50 ± 0.026 |
| 10 ⁻⁴ | Mtb-Low | 20.6 ± 0.2 | Negative | Negative | Negative |
| 10 ⁻⁵ | Mtb-Very Low | 24.53 ± 0.25 | Negative | Negative | Negative |




Fig. 3. GTC and heat treatment of pulmonary TB sputum samples A) Comparison of TB-MBLA results between GTC and heat-treated TB sputum specimens; B) Comparison of the time to positivity in MGIT between heat and GTC treated TB sputum specimens and with un-treated controls.

4. Discussion

Presence of innovative methods with the capability to inactivate live *M. tb* in pulmonary patient’s specimens without compromising the integrity of nucleic acids will accelerate implementation and uptake of TB MBLA, opening the prospect of widespread availability of accurate treatment monitoring. We confirm that application of heat (80 °C for 20 min) or GTC inactivates *M. tb* and does not compromise integrity of the 16S rRNA of *M. tb* measured by TB-MBLA. This laboratory based study confirms results reported in early previous studies [24,25] and adds additional insight that TB-MBLA could be performed in routine health-care settings that lack BSL3 laboratory infrastructure thus, providing an opportunity for improved assessment of TB therapy.

Except for the Xpert MTB/RIF Assay based tests, most TB molecular tests require extraction of nucleic acids in complex clinical specimen such as sputum [17,26]. The use of manual or automated extraction can create aerosols and present biosafety risk to laboratory staff, justifying the need for BSL3 laboratory equipped with Class II biological safety cabinets [25]. We observed a higher yield of Xpert MTB/RIF Ultra for detection of *M. tb* than both culture and TB-MBLA in low bacillary load sputum samples. The higher yield in the Xpert MTB/RIF Ultra may be explained by the presence of an automated, closed extraction and PCR systems whilst TB-MBLA employ manual extraction that includes several purification steps that may compromise the RNA yield and quantity [22, 27]. Interestingly we were able to demonstrate similar detection threshold between TB-MBLA and MGIT liquid culture in low bacterial load samples which concurs with our previous reports with a Spearman’s rank correlation coefficient (r) of -0.81; 95% confidence interval (CI) of -0.86 to -0.74 [28,29]. This relationship further confirm that the two tests detect viable *M. tb* in clinical specimens [30].

Routinely, TB tests that require extraction of nucleic acids such as Line Probe Assay (LPA) or sputum processing for culture and phenotypic drug susceptibility test (DST) with potential to generate aerosols are not performed in routine healthcare laboratories. Samples for such tests are referred to the Zonal or TB reference laboratory equipped with BSL3 laboratory infrastructure for processing. However, it is well known that BSL3 laboratory infrastructures are very expensive to establish and maintain in resource limited countries overburdened with both TB and HIV diseases [16], therefore, such infrastructure are very limited. In this study, both heat and GTC treatment completely inactivated viable *M. tb* and results were confirmed by negative growth of *M.tb* in standard MGIT liquid culture [21]. Our result suggests that these innovative methods may obviate the need for performing TB-MBLA test in BSL3 laboratory infrastructure, procedures that are very much needed in high TB burden settings.

We note that both heat and GTC inactivate viable *M. tb* and offer an adequate RNA yield for TB-MBLA. However, we found a slightly higher yield of RNA for GTC treated H37Rv, *M. tb* culture and real pulmonary TB sputum samples compared to heat treatment by heating at 80 °C. It is well known that GTC is chaotropic salt which inactivates RNases enzymes responsible for degradation of RNA after cell death [25]. Our previous work showed a higher degradation of RNA after addition of RNases in heat inactivated cultures than those without RNases, which may explain reason for low RNA yield in heat treated H37Rv cultures or sputum compared to GTC treated specimens [24]. This means the use of GTC for inactivation of *M. tb* have a beneficial effect in preserving RNA required for TB-MBLA from degrading enzymes [25] and adds additional information that RNases may not be eliminated by heat inactivation [24].

In summary, the findings from the present study confirms that *M. tb* is nonviable after inactivation with either heat or GTC buffer and reduction in 16S rRNA quantity is not significant and thus may not compromise downstream quantification of bacillary load by TB-MBLA. This has implications on the methodologies for application in health-care settings without BSL3 laboratory infrastructure for scale up of novel TB-tests. These results will provide the guidance to scientists,

researchers, National TB programme managers and policy makers implementing TB molecular tests in routine healthcare laboratories.

Transparency and declaration

Authors have no conflicts of interest to declare.

CRediT authorship contribution statement

Bariki Mtafya: Conceived and designed the study, Developed the research tools, performed laboratory experiments, analysed the data, created figures, tables and drafted the manuscript. **Paschal Qwaray:** performed laboratory experiments. **Joseph John:** performed laboratory experiments. **Emanuel Sichone:** performed laboratory experiments. **Alice Shoo:** performed laboratory experiments. **Stephen H. Gillespie:** Conceived, designed the study and proofread the manuscript. **Nyanda Elias Ntinginya:** Conceived and designed the study, Developed the research tools and proofread the manuscript. **Wilber Sabiiti:** Conceived and designed the study, Developed the research tools, analysed the data and proofread the manuscript, All authors reviewed the manuscript before submission.

Acknowledgment

This work was supported by the United Kingdom Research and Innovation (UKRI) for the TB-MBLA translation study in routine healthcare settings through the Food Security and Health for East Africa Project (Grant number# EP/T01525X/1) as part of postdoctoral work for Dr Bariki Mtafya at the National Institute for Medical Research Centre-Mbeya Medical Research Centre, Mbeya, Tanzania with the University of St Andrews, UK. We thank the patients who provided the samples for research experiments and staff at the Rungwe district hospital for supervising recruitment of the study participants.

References

- [1] Lange C, Dheda K, Chesov D, et al. Management of drug-resistant tuberculosis. *Lancet* 2019;394:953–66.
- [2] Xu Y, Wang G, Xu M. Biohazard levels and biosafety protection for Mycobacterium tuberculosis strains with different virulence. 2020.
- [3] Centers for Disease Control and Prevention (CDC). CDC at the forefront of innovation in the global fight against TB, vols. 1–2; 2021.
- [4] Global Tuberculosis Report. Global tuberculosis report. Geneva, Switzerland. 2021.
- [5] Senkoro M, Mfinanga S, Egwaga S, et al. Prevalence of pulmonary tuberculosis in adult population of Tanzania: a national survey, 2012. *Int J Tubercul Lung Dis* 2016;20:1014–21.
- [6] Harrington JM, Shannon HS. Incidence of tuberculosis, hepatitis, brucellosis, and shigellosis in British medical laboratory workers. *Br Med J* 1976;1:759–62.
- [7] Reid DD. Incidence of tuberculosis among workers in medical laboratories. *Br Med J* 1957;2:10.
- [8] Sewell DL. Laboratory-associated infections and biosafety. *Clin Microbiol Rev* 1995;8:389–405.
- [9] World Health Organization. Tuberculosis laboratory biosafety manual. *World Heal Organ Publ*; 2012. p. 1–60. 978 92 41504638.
- [10] Knudsen RC, Mathews HM, Richmond JY, et al. Laboratory management of agents associated with hantavirus pulmonary syndrome; interim biosafety guidelines.
- [11] Nolan CM, Blumberg HM, Taylor Z, et al. American thoracic society/centers for disease control and prevention/infectious diseases society of America: controlling tuberculosis in the United States. *Am J Respir Crit Care Med* 2005;172:1169.
- [12] Zijenah LS. The world health organization recommended TB diagnostic tools. In: *Tuberculosis*. Geneva, Switzerland, pp. 72–90.
- [13] Boehme CC, Nicol MP, Nabeta P, et al. Feasibility, diagnostic accuracy, and effectiveness of decentralised use of the Xpert MTB/RIF test for diagnosis of tuberculosis and multidrug resistance: a multicentre implementation study. *Lancet* 2011;377:1495–505.
- [14] Yadav RN, Singh BK, Sharma SK, et al. Comparative evaluation of GenoType MTBDR plus line probe assay with solid culture method in early diagnosis of multidrug resistant tuberculosis (MDR-TB) at a tertiary care centre in India. *PLoS One* 2013;8:e72036.
- [15] Hai HT, Vinh DN, Thu DDA, et al. Comparison of the Mycobacterium tuberculosis molecular bacterial load assay, microscopy and GeneXpert versus liquid culture for viable bacterial load quantification before and after starting pulmonary tuberculosis treatment. *Tuberculosis* 2019;119:101864.
- [16] Parsons LM, Somoskovi A, Gutierrez C, et al. Laboratory diagnosis of tuberculosis in resource-poor Countries: challenges and opportunities. *Clin Microbiol Rev* 2011; 24:314–50.
- [17] Allen V, Nicol MP, Tow LA. Sputum processing prior to Mycobacterium tuberculosis detection by culture or nucleic acid amplification testing : a narrative review. *J Microbiol bitotechnology* 2016;5:96–108.
- [18] Stop TB. PARTnership: the paradigm shift. 2016-2020. Global plan to end TB. Geneva, Switz Stop TB Partnersh.
- [19] Floyd K, Glaziou P, Houben R, et al. Global tuberculosis targets and milestones set for 2016–2035: definition and rationale. *Int J Tubercul Lung Dis* 2018;22:723–30.
- [20] Bimba JS, Lawson L, Kontogianni K, et al. Primestore MTM and OMNIgene sputum for the preservation of sputum for xpert MTB/RIF testing in Nigeria. *J Clin Med* 2019;8:1–8.
- [21] Siddiqi SH, Rüsçh-Gerdes S. *MGIT procedure manual: for BACTEC™ MGIT 960™ TB system*. Geneva, Switzerland: Foundation for Innovative New Diagnostics; 2006. http://www.finndiagnostics.org/export/sites/default/resource-centre/find_documentation/pdfs/mgit_manual_nov_2007.pdf.
- [22] Gillespie SH, Sabiiti W, Oravcova K. Mycobacterial load assay. In: Bishop-Lilly Kimberly A, editor. *Diagnostic bacteriology: methods and protocols, methods in molecular biology*; 2017. p. 89–105.
- [23] Bowness R, Boeree MJ, Aarnoutse R, et al. The relationship between mycobacterium tuberculosis mgit time to positivity and cfu in sputum samples demonstrates changing bacterial phenotypes potentially reflecting the impact of chemotherapy on critical sub-populations. *J Antimicrob Chemother* 2015;70: 448–55.
- [24] Sabiiti W, Azam K, Esmeraldo E, et al. Heat inactivation renders sputum safe and preserves Mycobacterium tuberculosis RNA for downstream molecular tests. *J Clin Microbiol* 2019;57:1–8.
- [25] Clinghan R, Anderson TP, Everett V, et al. Viability of Mycobacterium tuberculosis after processing with commercial nucleic acid extraction kits. *J Clin Microbiol* 2013;51:2010.
- [26] Anahtar MN, Bowman BA, Kwon DS. Efficient nucleic acid extraction and 16S rRNA gene sequencing for bacterial community characterization. *JoVE* 2016;15: 1–11.
- [27] Honeyborne I, McHugh TD, Phillips PPJ, et al. Molecular bacterial load assay, a culture-free biomarker for rapid and accurate quantification of sputum Mycobacterium tuberculosis bacillary load during treatment. *J Clin Microbiol* 2011;49:3905–11.
- [28] Honeyborne I, Mtafya B, Phillips PPJ, et al. The molecular bacterial load assay replaces solid culture for measuring early bactericidal response to antituberculosis treatment. *J Clin Microbiol* 2014;52:3064–7.
- [29] Sabiiti W, Mtafya B, De Lima DA, et al. A tuberculosis molecular bacterial load assay (TB-MBLA) | protocol. *JoVE* 2020:e60460.
- [30] Mtafya B, Sabiiti W, Sabi I, et al. Molecular bacterial load assay concurs with culture on NaOH-induced loss of mycobacterium tuberculosis viability. *J Clin Microbiol* 2019;57. <https://doi.org/10.1128/JCM.01992-18>. Epub ahead of print.