




A service of the National Library of Medicine
and the National Institutes of Health

www.pubmed.gov

My NCBI [?] [Sign In] [Register]

All Databases PubMed Nucleotide Protein Genome Structure OMIM PMC Journals Books

Search PubMed for [] Go Clear

Limits Preview/Index History Clipboard Details

Display AbstractPlus Show 20 Sort by Send to

All: 1 Review: 0

1: [Ophthalmic Physiol Opt.](#) 2005 Sep;25(5):436-40.

 **Full text** [Links](#)
blackwell-synergy.com

Comparison of the ICare rebound tonometer with the Goldmann tonometer in a normal population.

Fernandes P, Diaz-Rey JA, Queiros A, Gonzalez-Meijome JM, Jorge J.

Department of Physics (Optometry), School of Sciences, University of Minho, Braga, Portugal.
pfernandes@fisica.uminho.pt

The aim of this study was to evaluate the accuracy of measurement of intraocular pressure (IOP) using a new induction/impact rebound tonometer (ICare) in comparison with the Goldmann applanation tonometer (AT). The left eyes of 46 university students were assessed with the two tonometers, with induction tonometry being performed first. The ICare was handled by an optometrist and the Goldmann tonometer by an ophthalmologist. In this study, statistically significant differences were found when comparing the ICare rebound tonometer with applanation tonometry (AT) ($p < 0.05$). The mean difference between the two tonometers was 1.34 ± 2.03 mmHg (mean \pm S.D.) and the 95% limits of agreement were ± 3.98 mmHg. A frequency distribution of the differences demonstrated that in more than 80% of cases the IOP readings differed by < 3 mmHg between the ICare and the AT. In the present population the ICare overestimates the IOP value by 1.34 mmHg on average when compared with Goldmann tonometer. Nevertheless, the ICare tonometer may be helpful as a screening tool when Goldmann applanation tonometry is not applicable or not recommended, as it is able to estimate IOP within a range of ± 3.00 mmHg in more than 80% of the population.

PMID: 16101950 [PubMed - indexed for MEDLINE]

Related Links

Accuracy of the new ICare rebound tonometer vs. other portable tonometers in healthy eyes. [Optom Vis Sci. 2006]

Clinical evaluation of rebound tonometer. [Acta Ophthalmol Scand. 2006]

Comparison of the rebound tonometer (ICare) to the applanation tonometer (Tonopen XL) in normotensive dogs. [Vet Ophthalmol. 2006]

Three portable tonometers, the TGDc-01, the ICARE and the Tonopen XL, compared with each other and with Goldmann applanation tonometry. [Ophthalmic Physiol Opt. 2005]

Comparison of ICare tonometer with Goldmann applanation tonometer in glaucoma patients. [Acta Ophthalmol Scand. 2006]

[See all Related Articles...](#)

Display Show Sort by Send to

[Write to the Help Desk](#)
[NCBI](#) | [NLM](#) | [NIH](#)
[Department of Health & Human Services](#)
[Privacy Statement](#) | [Freedom of Information Act](#) | [Disclaimer](#)

Oct 30 2006 07:27:22