

DETERMINATION OF β -LACTAMS, FLUOROQUINOLONES, SULFONAMIDES/TRIMETHOPRIM, TETRACYCLINES AND COLISTIN IN POULTRY EXCRETA BY UHPLC-MS/MS

B. Gorissen^{1,2}, T. Reyns¹, M. Devreese², P. De Backer², J. Van Loco¹, S. Croubels²


¹Scientific Institute of Public Health, Department of Food, Medicines and Consumer Safety, Service of Chemical Residues and Contaminants, Rue Juliette Wytsmanstraat 14, 1050 Brussels, Belgium

²Ghent University, Department of Pharmacology, Toxicology and Biochemistry, Faculty of Veterinary Medicine, Salisburylaan 133, 9820 Merelbeke, Belgium

Introduction and aims

- Due to the risk of antimicrobial resistance development, the Belgian Royal Decree concerning the eradication of *Salmonella* (C – 2007/22784) prohibits treatment of poultry with antimicrobials against zoonotic *Salmonella* spp. Nowadays, the unauthorized use of antimicrobials is determined through analysis of tissue samples of sacrificed animals. There is a need to develop a more animal-friendly method to detect illicit use of antimicrobials.
- In this study, an analytical method using UHPLC-MS/MS for the determination of antimicrobial residues in samples of poultry manure was developed and validated for β -lactams (amoxicillin (AMO) and phenoxymethylpenicillin), fluoroquinolones (enrofloxacin, difloxacin (DFX) and flumequine), sulfonamides in combination with trimethoprim (TMP) (sulfachloropyridazine (SCP), sulfadiazine and sulfaclozine), tetracyclines (chlortetracycline and doxycycline (DOX)) and polymyxins (colistin).
- Next, an animal experiment was conducted to gain insight into excretion of the selected compounds. The developed method was applied in the determination of target compounds in derived samples.

Animal experiment




Twelve laying hens were equally divided into one control group and five experimental groups.

For each antimicrobial class, one active compound was selected and administered by route of drinking water.

Between every sample point, floor bedding, feed (and, during administration, medicated drinking water) were refreshed.

Sample preparation



Results




Figure 1. MS/MS chromatograms for a blank poultry manure sample spiked at 50 ng/g (LOQ)

Antimicrobial concentrations in poultry excreta




Figure 2. Results of an excretion study performed in laying hens, after oral administration of either DFX, SCP/TMP combination, DOX or AMO

Conclusions

- A quick and simple UHPLC-MS/MS method for determination of antimicrobials in non-invasive samples of poultry manure was developed and validated.
- The applicability of the developed method was tested on biological samples in a preliminary animal experiment; all studied compounds were successfully extracted.
- The method is an animal-friendly alternative, suitable for high-throughput analysis in official control programmes to tackle the resistance development of *Salmonella*.