

Solid-state heat stability of peptides: kinetics and mechanisms


 Matthias D'Hondt¹, Chien-Yu Peng², Ralf Hoffmann³ and Bart De Spiegeleer^{1,*}
¹ DruQuaR (Drug Quality & Registration) group, Faculty of Pharmaceutical Sciences, Ghent University, Harelbekestraat 72, B-9000 Ghent, Belgium.

*bart.despiegeleer@ugent.be (Our Ref.: 2014-359b)

² Institute of Statistical Science, Academia Sinica, Academia Road 128, 11529 Taipei, Taiwan.


³ Center for Biotechnology and Biomedicine, Institute for Bioanalytical Chemistry, Leipzig University, Deutscher Platz 5, 04103 Leipzig, Germany.

INTRODUCTION

 Biologically active peptides
 Therapeutics, e.g. oncology (buserelin)


Processing hot melt extrusion (HME)

1. High temperature exposure ← OBJECTIVE
2. Mechanical shear stress
3. Polymer/matrix influence



EXPERIMENTAL

Stability indicating UPLC method:

 Acquity BEH300 C18 1.7µm (2.1 ×100 mm)
 MF A: 95/5 H₂O/ACN + formic acid
 MF B: 5/95 H₂O/ACN + formic acid

 1.5 min isocratic hold at 100% A
 linear gradient from 0 to 21% B in 9.5 min
 7 min isocratic hold

Dry heat conditions

T (°C)	150	157.5	165	172.5	180
Time (t) (min)	40	25	15	10	10
	80	50	30	20	20
	120	75	45	30	30
	160	100	60	40	40

Detection

- DAD-UV (kinetics via normalized areas)
- MS/MS (degradant identification)

Kinetic data evaluation

•Statistical evaluation of 17 solid-state kinetic models:

- ✓ Nucleation (7)
- ✓ Geometrical contraction (2)
- ✓ Diffusion (4)
- ✓ Reaction-order (4)

•Extrapolation to HME-related conditions


RESULTS and DISCUSSION

Kinetic data evaluation per temperature

Ginstling-Brounshtein (Diffusion model): minimal AIC values

$$1 - (2\alpha/3) - (1 - \alpha)^{2/3} = k \times t$$


 α = fraction degraded

 5 degradation constant $k \Leftrightarrow$ temperature T


Predicted degradation at HME-related conditions

	Polynomial regr.	Linear regr. solid state
5 min 100°C	<0.01%	0.10%
5 min 125°C	<0.01%	1.33%

Degradant identification



Degradation mechanism

 β -elimination + fragmentation

pGlu-His-Trp-NH₂ *
 pGlu-His-Trp-Ser-Tyr-NH₂ *
 pyruvoyl-Tyr-Ser(tBu)-Leu-Arg-Pro-NH-Et

Backbone hydrolysis

pGlu-His-Trp *
 Tyr-Ser(tBu)-Leu-Arg-Pro-NH-Et *
 Ala-Tyr-Ser(tBu)-Leu-Arg-Pro-NH-Et *

Isomerisation*

pGlu-His-Trp-Ser-Tyr-Ser(tBu)-Leu-Arg-Pro-NH-Et

CONCLUSIONS

(1) Kinetics: Ginstling-Brounshtein degradation model: HME ✓

 (2) Degradant profiling:

1. β -elimination
2. Backbone hydrolysis
3. Isomerisation

On-going:

- (1) Manufacturing HME implant
- (2) Characterisation, incl. dissolution
- (3) Stability
- (4) Pharmacokinetics *in-vivo*: (mice subdermal)

REFERENCES