

DEVELOPMENT OF A CYTOMETRIC BEAD ARRAY SCREENING TOOL

Simultaneous detection of pro-inflammatory cytokines in plasma of lipopolysaccharide-challenged pigs

H. WYNS, S. CROUBELS, K. DEMEYERE, A. WATTEYN, P. DE BACKER, E. MEYER

Department of Pharmacology, Toxicology and Biochemistry, Ghent University, Faculty of Veterinary Medicine, Salisburylaan 133, 9820 Merelbeke, Belgium, Heidi.Wyns@UGent.be

Introduction

Lipopolysaccharide (LPS) has been widely used as a model of immune challenge in pigs as this compound induces the immediate synthesis of pro-inflammatory cytokines such as tumor necrosis factor- α (TNF- α), interleukin-1 β (IL-1 β) and IL-6. In research, multiplex assays currently are a very popular tool for the simultaneous detection of biomarkers of infection and inflammation. Specific and sensitive Enzyme-Linked Immuno Sorbent Assays (ELISAs) are well-suited to perform single factor analysis, yet for multiparameter analyses, this approach is time-consuming and expensive. Cytometric bead array (CBA) is a flexible, bead-based flow cytometric application for the simultaneous detection of various soluble proteins of interest. The aim of the present study was to develop and validate a CBA 3-plex assay for the major pro-inflammatory cytokines TNF- α , IL-1 β and IL-6. The results were compared to commercial ELISA kits.

Materials and Methods

Experimental design

Four male pigs (Seghers Hybrid), with a mean (\pm SD) body weight (BW) of 24.9 ± 3.17 kg were intravenously challenged with $15 \mu\text{g}$ ultrapure LPS/kg BW (*Escherichia coli* serotype O111:B4). ELISAs were purchased from R&D Systems.

CBA 3-plex assay for TNF- α , IL-1 β and IL-6

An overview of the CBA assay is summarized and illustrated in Figure 1. More details of the protocol are described by Wyns *et al.* (2012).




Figure 1. CBA assay

Results

Table 1 reports the limits of detection (LODs), intra- and inter-assay variations (CVs) and dynamic ranges of the CBA 3-plex cytokine assay. Following an *in vivo* LPS challenge, similar plasma concentration-time profiles were observed for all cytokines with CBA and ELISA as shown in Figure 2.




Figure 2. Mean (\pm SD) plasma concentration-time curves of TNF- α , IL-1 β and IL-6 obtained from LPS-challenged pigs ($n = 4$). Adjusted p -values are reported as a ($p < 0.0005$), b ($p = 0.033$), c ($p = 0.004$), d ($p = 0.005$), e ($p = 0.039$), f ($p = 0.019$) and g ($p = 0.026$). Concentrations were measured using both ELISA (left panel) and CBA (right panel). Note different scales on the y-axes for TNF- α , IL-1 β and IL-6.

Table 1. LOD, mean intra- and inter-assay CVs and dynamic ranges of each cytokine in the CBA 3-plex assay and commercially available ELISAs.

Cytokine	LOD (ng/mL)		Intra-assay (CV%)		Inter-assay (CV%)		Dynamic range (ng/mL)	
	CBA	ELISA	CBA	ELISA	CBA	ELISA	CBA	ELISA
TNF- α	0.363	0.004	8.45	4.86	15.71	8.90	0.50 - 50	0.023 - 1.50
IL-1 β	0.109	0.007	2.26	5.77	15.48	6.27	0.50 - 50	0.039 - 2.50
IL-6	0.005	0.002	2.33	4.00	11.41	5.97	0.50 - 50	0.010 - 1.20

Discussion and Conclusions

CBA and ELISA show similar cytokine concentration-time profiles in plasma. Therefore, the optimised and validated CBA 3-plex cytokine protocol provides a fast, flexible and cost-effective screening tool for simultaneous measurement of the major porcine pro-inflammatory cytokines TNF- α , IL-1 β and IL-6. This technique will be applied in future research to study the immunomodulatory properties of drugs in a porcine LPS inflammation model.

Wyns, H., Croubels, S., Demeyere, K., Watteyn, A., De Backer, P., Meyer, E., Development of a cytometric bead array screening tool for the simultaneous detection of pro-inflammatory cytokines in porcine plasma. *Vet. Immunol. Immunopathol.* (2012), <http://dx.doi.org/10.1016/j.vetimm.2012.09.041> (in press)

