#### University of Massachusetts Medical School

### eScholarship@UMMS

**Open Access Articles** 

Open Access Publications by UMMS Authors

2016-04-06

## DNA immunization as a technology platform for monoclonal antibody induction

Shuying Liu SL Consulting

Et al.

## Let us know how access to this document benefits you.

Follow this and additional works at: https://escholarship.umassmed.edu/oapubs



Part of the Immunoprophylaxis and Therapy Commons

#### **Repository Citation**

Liu S, Wang S, Lu S. (2016). DNA immunization as a technology platform for monoclonal antibody induction. Open Access Articles. https://doi.org/10.1038/emi.2016.27. Retrieved from https://escholarship.umassmed.edu/oapubs/2864

Creative Commons License



This work is licensed under a Creative Commons Attribution 4.0 License.

This material is brought to you by eScholarship@UMMS. It has been accepted for inclusion in Open Access Articles by an authorized administrator of eScholarship@UMMS. For more information, please contact Lisa.Palmer@umassmed.edu.

## npg

#### **REVIEW**

# DNA immunization as a technology platform for monoclonal antibody induction

Shuying Liu<sup>1</sup>, Shixia Wang<sup>2</sup> and Shan Lu<sup>2</sup>

To combat the threat of many emerging infectious diseases, DNA immunization offers a unique and powerful approach to the production of high-quality monoclonal antibodies (mAbs) against various pathogens. Compared with traditional protein-based immunization approaches, DNA immunization is efficient for testing novel immunogen designs, does not require the production or purification of proteins from a pathogen or the use of recombinant protein technology and is effective at generating mAbs against conformation-sensitive targets. Although significant progress in the use of DNA immunization to generate mAbs has been made over the last two decades, the literature does not contain an updated summary of this experience. The current review provides a comprehensive analysis of the literature, including our own work, describing the use of DNA immunization to produce highly functional mAbs, in particular, those against emerging infectious diseases. Critical factors such as immunogen design, delivery approach, immunization schedule, use of immune modulators and the role of final boost immunization are discussed in detail. *Emerging Microbes and Infections* (2016) 5, e33; doi:10.1038/emi.2016.27; published online 6 April 2016

Keywords: antigen delivery; DNA immunization; hybridoma; immunogenicity; monoclonal antibody

#### INTRODUCTION

The recent successful use of protective monoclonal antibodies as a life-saving treatment in Ebola virus-infected humans<sup>1</sup> highlights the need for the development of new technologies that are fast and efficient in eliciting functional monoclonal antibodies (mAbs) to control emerging infectious diseases.

DNA immunization as it exists today was pioneered in the early 1990s. Its initial use as a vaccination platform generated great excitement due to the overall simplicity of using DNA plasmids to deliver immunogens.<sup>2–5</sup> One particularly attractive feature of DNA vaccines is that immunogens are produced *in vivo*, giving them the ability to induce T-cell immune responses through endogenous antigen processing and presentation pathways. However, the application of DNA immunization for human vaccine development has encountered challenges, specifically the low immunogenicity identified in early clinical studies, in which DNA vaccines were delivered via traditional needle injection without the use of adjuvants or other types of delivery instruments.

In recent years, significant progress has been made in the application of DNA vaccines in humans via two strategies: (i) the use of physical delivery methods, such as a gene gun or electroporation, which have greatly improved the immunogenicity of DNA vaccines in human volunteers, <sup>6,7</sup> and (ii) the development of a heterologous prime-boost scheme, <sup>8</sup> in which the hosts are first immunized with a DNA vaccine, followed by boost immunizations with either recombinant protein antigens or traditional killed or live attenuated vaccines, <sup>9–11</sup> which are more effective than homologous prime-boost immunizations using multiple doses of the same recombinant proteins or traditional vaccines, <sup>12–16</sup>

At the same time, it is generally agreed that DNA immunization is effective in small animal models, which can be useful for the production of high-quality mAbs. However, early animal model studies have focused mainly on T-cell immune responses, given the unique advantage of DNA immunization to induce such a response. <sup>17,18</sup> Little attention has been paid to the value of DNA immunization to elicit high-quality antibody responses, and much less has been given to the potential of DNA immunization to induce high-quality B-cell responses, which can be useful for the production of functional mAbs.

The current review provides a comprehensive summary of the literature that has accumulated over the past two decades, including our own work, in using DNA immunization to produce highly functional mAbs, in particular, against emerging infectious diseases. DNA immunization is more useful than traditional approaches to generating mAbs against more difficult targets, especially membrane proteins.

# DESIGN OF DNA VACCINES FOR MONOCLONAL ANTIBODY INDUCTION

DNA vaccines are constructed as mammalian expression vectors. Both the choice of expression vector and the design of immunogen inserts are important for the final antibody response, including the production of mAbs against the desired antigenic structures, but their roles are different in the process of DNA immunization.

#### Optimal design of immunogen inserts

The *in vivo* immunogen expression feature of DNA vaccines offers a number of benefits.



Traditional protein-based immunization approaches have difficulty producing full-length protein immunogens by the recombinant protein method if the proteins are naturally expressed in a membraneassociated format, such as the multi-transmembrane G-protein coupled receptors (GPCRs) and ion channels. The DNA immunization approach can circumvent these problems because full-length proteins can be expressed in vivo when they are delivered in the form of DNA vaccines. Furthermore, it is well known that the structural integrity of proteins is critical for the induction of functional mAbs, yet these sensitive structures tend to be lost during the in vitro protein production process, regardless of whether they are produced as recombinant proteins or are extracted directly from cells or other sources in which the proteins are naturally expressed. Production of functionally active mAbs is highly dependent on the conformation of the proteins. Expressing intact immunogens in vivo by DNA immunization appears to have the best chance of inducing mAbs with the desired biological activities.

DNA vaccines possess the unique advantage of immunogen design flexibility. Immunogen inserts expressing the full-length sequences of target proteins are commonly used for all types of proteins, especially transmembrane proteins, with good success. <sup>19–30</sup> For intracellular proteins, one may assume that it is necessary to re-direct intracellular proteins into secretory pathways by adding a signal peptide to elicit a better antibody response. However, in a number of mAb production studies, native proteins have been used as DNA vaccine immunogen

Table 1 Types of DNA vaccine immunogens used for mAb induction

•••	•	
Immunogen inserts	Original types of proteins	References
Full length	Transmembrane	19–30
	Extracellular matrix	42
	GPI anchor	43
	Intracellular	25,28,31,32
	Secretory	25,44-47
	Viral	25,48-53
	Bacteria	54
Extracellular domain	Single transmembrane	33–37
Fragment/subunit	Any type	25,38,39,55
Mini-gene insert	Any type	40
Novel format, vaccibody	Secretory	56
Immunogen-transmembrane domain fusion	Viral non-structure	41
Immunogen-Fc fusion	Intracellular	57,58

Abbreviation: glycophosphatidylinositol, GPI.

inserts without any sequence modifications.<sup>25,28,31,32</sup> On the other hand, the gene sequences of immunogen inserts for DNA vaccines can be easily edited to express 'designer proteins.' For a singletransmembrane protein, the extracellular domain of a secretory protein can be selectively cloned as the immunogen insert when the objective is to generate mAbs against epitopes on the extracellular domain. 33-37 Such an approach has also been applied to bacterial toxins. For example, a truncated fragment can be used for immunization in place of a full-length potentially lethal toxin protein, thereby avoiding the introduction of unwanted biological activity during the production of DNA vaccines, as well as during animal immunization. 38,39 Additional immunogen manipulations include the production of a 'mini-gene insert' to express a short peptide sequence to cover a receptor-binding domain. 40 In this case, antigenic determinants in the angiotensinconverting enzyme 2 binding domain of the severe acute respiratory syndrome spike protein, which does not closely match other coronaviruses, were predicted using software PROTEAN to induce anti-spike protein antibodies. Alternatively, a transmembrane anchor sequence can be added to non-membrane-associated antigens.<sup>41</sup> As a simple and flexible immunogen design approach, DNA immunization offers a wide range of options to produce novel immunogen inserts for the induction of mAbs against even the most challenging targets (Table 1).

One unique feature of DNA immunization is the convenience of using the same DNA vaccine constructs to express antigens for mAb screening. The choice of reagents and methods depends on the original type of protein (Table 2). Cell-associated antigen-based screening has been widely and successfully used for mAbs targeting transmembrane proteins, viral envelope proteins and intracellular proteins. In these cases, cells expressing the immunogens were used without the need for protein purification to screen the binding activity of mAbs by either fluorescence-activated cell sorting analysis, <sup>19–21,23,24,26,27,29,30,34,35</sup> whole cell enzyme-linked immunosorbent assay (ELISA), <sup>33,43,48,49</sup> or immunohistochemistry (IHC) methods. <sup>25,31,32</sup> In one study, a novel in-cell Western screening method was developed and optimized to generate monoclonal antibodies against nuclear, cytoplasmic and transmembrane proteins. <sup>28</sup>

#### Role of expression vectors

After two decades of effort by many research groups, the design of commonly used DNA vaccine vectors has been significantly optimized. It has been established that the promoter region of these DNA vaccine vectors is their most critical element and that the function of promoters can be further enhanced by other associated regulatory components. For example, the intron A sequence

Table 2 Screening methods

Screening reagent	Original type of protein	Screening method	References
Proteins	Secretory	ELISA	37-40,42,44-47,50,51,54-56
	Intracellular	ELISA	36,57,58
	Single transmembrane	ELISA	21,22
Cells expressing immunogens	Transmembrane	FACS	19-21,23,24,26,27,29,30,34,35
		Whole cell ELISA	33,43
	Viral envelope	Whole cell ELISA	48,49
		Immunostaining	53
	Intracellular	In-cell western	28
		IHC	25,31,32
Viral particles	Viral	Dot enzyme immunoassay	41
		ELISA	52

Abbreviation: enzyme-linked immunosorbent assay, ELISA; fluorescence-activated cell sorting, FACS; immunohistochemistry, IHC



associated with the cytomegalovirus (CMV) promoter can significantly enhance the function of CMV promoter.<sup>59</sup> The selection of poly-A tail is also important, as is optimized codon usage, which is particularly important for DNA vaccines expressing genes from infectious pathogens that often employ different codon usage than mammalian proteins.

One area that is less studied is whether sustained immunogen expression in vivo has a major impact on the induction of high-quality mAbs. In one study of antibody generation via intravenous delivery of plasmid DNA,<sup>25</sup> the relative efficacy of CMV promoter and the human ubiquitin C promoter was compared using luciferase as the model immunogen. Using the ubiquitin promoter, which can lead to sustained antigen expression in the liver, resulted in significant levels of antibody titers seven weeks after a single hydrodynamic tail vein (HTV) delivery. In contrast, a single HTV delivery of the CMV promoter, which results in only short-term antigen expression, produced very low antibody titers over the same time frame. Nevertheless, both promoters enabled high antibody titers when delivered repeatedly. Thus, the study authors concluded that the choice of expression vector may reduce the number of genetic immunizations while still promoting the induction of high-titer antigen-specific antibodies.

Similarly, immune responses were compared after employing DNA vaccine plasmids encoding multi-drug resistant protein 4 (MRP4), a 12-transmembrane transporter, but with different promoters: the CMV early enhancer/chicken beta actin promoter, which is a strong synthetic promoter frequently used to constitutively drive high levels of gene expression in mammalian expression vectors, and the classical CMV promoter for transient expression.<sup>30</sup> The immunogenicity results showed that the chicken beta actin promoter induced a higher antigen-specific immune response with HTV delivery, even when the plasmid was injected less frequently, compared with the use of a CMV promoter-containing plasmid.

#### **DELIVERY APPROACH AND SCHEDULE**

#### Physical versus chemical delivery approaches

Since the early discovery of DNA immunization, a wide range of delivery approaches has been studied. These approaches can be divided into two main categories. One is traditional needle injection of DNA plasmids in various solutions. Additional facilitating agents such as lipids and nanoparticles can be included in the solution to enhance delivery efficacy, and the composition of the chemical solution determines the uptake efficacy of the DNA vaccine plasmid. The second delivery type is based on physical forces. The most representative approach is the use of a gene gun, which uses a 'ballistic' force to deliver the DNA plasmids. First, the DNA plasmids are coupled with gold particles, which are then delivered by the ballistic force to penetrate the cells of the targeted tissues. Early generation gene guns created ballistic forces by adding high-voltage electricity to a drop of water.<sup>60</sup> More recent generation gene guns are based on the release of compressed gas. 61,62 Another physical method of DNA delivery is the use of electroporation. In this approach, DNA vaccines are first delivered by needle injection, followed by the application of an electrical current at the DNA injection site.

The relative immunogenicity of chemical and physical delivery approaches has been well analyzed. In one study, it was shown that electroporation delivery following intradermal needle injection was effective in delivering DNA vaccines to multiple intracellular compartments (that is, transmembrane, cytoplasm and nucleus), leading to the successful induction of mAbs.<sup>28</sup> Another study showed that electroporation followed by intramuscular needle injection generated higher antibody responses than intramuscular needle injection alone. This approach was also more effective than intramuscular needle injection alone when the DNA vaccine was formulated with a chemical polymer and protein immunogen.<sup>63</sup> In a more complete analysis comparing the delivery of DNA vaccines by intramuscular needle injection, electroporation following intramuscular needle injection, and gene gun alone, it was shown that both gene gun and electroporation delivery were more effective than the traditional intramuscular needle injection at eliciting higher antibody response levels.<sup>64</sup> Both the gene gun and electroporation approaches are effective, but the gene gun approach requires only a few micrograms of DNA plasmid to achieve the same level of immune response elicited by the electroporation approach, which requires at least 100 µg of DNA plasmid at the first step of intramuscular injection, even in mice.

One interesting finding from literature is that although the delivery approach may be critical for the induction of high-level immune responses for human vaccine development, different DNA delivery approaches have been similarly successful in producing mAbs against a wide range of target antigens. Table 3 lists the mAbs elicited by the gene gun approach;<sup>22,24,31-34,37,43,55,57,58</sup> injection, including intramuscular<sup>19,20,23,38,41,44–46,48–50</sup> intradermal<sup>21,35,42</sup> injection; and electroporation following intramuscular or intradermal injection. 27,28,39,52,53,56

One unique but less-studied approach is hydrodynamic intravenous delivery. One large study compared the relative immunogenicity of HTV and hydrodynamic limb vein delivery methods with 18 different antigens, including different types of target antigens (that is, intracellular, transmembrane, and secretory). Both methods were successful, but the hydrodynamic limb vein delivery method was especially potent for generating antibodies against a wide range of targets.<sup>25</sup> However, the hydrodynamic intravenous method may be more suitable for larger animals, such as rats and rabbits, given the larger vein size in these hosts compared with that in mice. The HTV method was also used successfully to generate a mAb against a 12-transmembrane transporter, which is a very challenging target for mAb induction.<sup>30</sup>

There have also been reports of producing mAbs with a single intrasplenic injection. 40,54 The author generated hybridomas by fusing spleen cells at 2, 3, 5, 10 and 25 days after a single intrasplenic injection of DNA vaccine plasmids. The highest number of specific hybridomas was generated at day 5 after a single initial injection.<sup>54</sup> However, these mAbs appeared to be useful only for immunoblotting, and no additional studies were conducted to characterize their affinities.

#### Immunization schedule

The optimal delivery schedule for the induction of high-quality mAbs by DNA immunization remains to be determined. The classical vaccine literature would indicate that an extended time period with long rest intervals may be more effective in eliciting high-quality antibody responses than immunizations in quick succession. There have been relatively few studies focusing solely on an optimal delivery schedule rather than the delivery approach itself. For example, it is not clear whether physical delivery approaches (that is, gene gun or electroporation) can be used more frequently than chemical delivery approaches (such as needle injection) because the former is more effective than the latter.

In most animal studies, the generation of mAbs requires multiple immunizations, usually every 2-3 weeks. However, a faster immunization procedure has been reported, which delivered 3-5 immunizations within 10–11 days at multiple sites by gene gun. 33,57,58 In these studies, lymph nodes were used for fusion 48 h after the last injection without the need for a final boost. The mAbs generated were used for



Table 3 DNA vaccine delivery approaches used for mAb induction

Delivery approaches	Original types of protein	Reference
Gene gun	Single transmembrane (Flt-3R)	33
	Intracellular (PED/PEA-15)	57
	Intracellular (annexin-V)	58
	Single transmembrane (CAR)	22
	Two- transmembrane (P2X7)	24
	GPI anchored enzyme	43
	Intracelluar (BCL-6)	31
	Intracelluar (MALT1)	32
	Single transmembrane (MHCI-related gene A)	34
	Parasite lipoprotein	55
	Viral envelop (HIV gp120)	37
IM	Bacteria toxin (Helicobacter pylori vacuolating cyto toxin)	38
	Seven transmembrane, GPCR (TSHR)	19
	Viral envelop (HGV E2)	48
	Seven transmembrane, GPCR (TSHR)	20
	Viral non-structure (Dengue NS1)	41
	Viral envelop (H5N1)	49
	Secretory protein, enzyme (prostate-specific antigen)	45
	Viral surface (HBV preS2/S)	50
	Seven transmembrane, GPCR (TSHR)	23
	Secretory protein, cytokine (CKLF1)	44
	Secretory protein, cytokine (Interferon beta)	46
ID	Single transmembrane (RET and CD30)	21
	Extracellular matrix and plasma glycoprotein (Fibulin-1)	42
	Single transmembrane (CD2)	35
IM followed by EP	Seven transmembrane, GPCR (CCX-CKR)	27
	Viral envelop (H5N1)	52
	Secretory (mCherry)	56
	Viral envelop (H1)	53
	Bacteria toxin (C. difficile toxin A)	39
ID followed by EP	Multiple targets (transmembrane, intracellular)	28

Abbreviations: B-cell lymphoma 6 protein, BCL-6; coxsackievirus and adenovirus receptor, CAR; cluster of differentiation 2, CD2; cluster of differentiation 30, CD30; chemokine-like factor 1, CKLF1; electroporation, EP; fetal liver tyrosine kinase 3 receptor, Flt-3R; hepatitis B virus, HBV; hepatitis G virus E2 protein, HGV E2; intradermal, ID; intramuscular, IM; mucosa-associated lymphoid tissue lymphoma translocation gene 1, MALT1; PZX purinoceptor 7, PZX7; phosphoprotein over expressed in diabetes/phosphoprotein enriched in astrocytes, PED/PEA-15; rearranged during transfection, RET; thyroid stimulating hormone receptor, TSHR.

fluorescence-activated cell sorting, western blot analysis and enzyme-linked immunosorbent assay.

More studies are needed to further optimize the DNA immunization schedule to elicit high-quality mAbs. It is important to determine whether different DNA delivery approaches are optimal with certain delivery schedules and whether the same delivery schedule can be applied to different animal species.

#### **IMMUNE MODULATION**

#### Molecular adjuvants

Due to the low immunogenicity of DNA vaccines in early human studies, great effort has been devoted to the inclusion of various adjuvants as part of DNA vaccine formulation or immunization. Such studies have been widely reported in the literature, and the following section will review the use of molecular adjuvants only in the context of mAb induction. Adjuvants identified for previous vaccine development are likely to be useful for the induction of mAbs, but additional

dedicated studies are needed to confirm the actual incremental benefits, particularly if DNA vaccine delivery is optimized, as discussed above.

Escherichia coli chaperone protein (GroEL) was demonstrated to act as a potent molecular adjuvant for DNA immunization and has been shown to enhance the Ab response against GPCRs, which are very difficult targets for mAb induction.<sup>27,65</sup> The authors reported that DNA immunization in mice with a plasmid encoding the full-length endothelin A receptor (ETAR) fused to GroEL at its C terminus induced strong, specific antibody responses to native ETAR. Co-injection of plasmids that expressed ETAR and GroEL (ETAR +GroEL) induced lower antibody responses than the ETAR-GroEL plasmid. In contrast, no specific antibody responses were produced in mice that were immunized with ETAR.65 The authors also mentioned that this strategy has been successfully applied to other GPCR targets and suggested that GroEL will be capable of producing antibodies against most GPCRs. One caution raised by the author is that in the case of unstable GPCRs, co-immunization with GroEL may be preferable to fusion with GroEL to induce an antibody response because fusion of some carrier proteins decreases their expression. Functional mAbs against a different GPCR target, the chemoCentryx chemokine receptor (CCX-CKR), were generated by the co-immunization of a plasmid encoding GroEL with a DNA vaccine plasmid-encoding target.<sup>27</sup>

In another study, a more complicated immunization strategy was tested, in which a plasmid encoding fetal liver tyrosine kinase 3 ligand was delivered as a priming dose, followed by the co-delivery of a plasmid encoding granulocyte-macrophage colony-stimulating factor and another plasmid encoding immunogen MRP4. This study demonstrated that the addition of fetal liver tyrosine kinase 3 ligand and granulocyte-macrophage colony-stimulating factor as immune modulators significantly improved not only the overall immune response in the mice but also the induction of antibodies capable of recognizing native extra-cellular epitopes.<sup>30</sup>

#### DNA prime-protein boost

One approach that presents a great advantage for the induction of high-titer and high-quality antibody responses is the heterologous prime-boost approach. In this approach, the DNA vaccine is delivered as the priming immunization, followed by a boost with protein antigens as recombinant proteins, peptides, or traditional inactivated or live attenuated vaccines. One unexpected finding regarding DNA priming immunization is their ability to induce higher-level antigenspecific B-cell responses. <sup>66</sup> Our research group has shown that DNA primer immunization was more effective than protein immunization in activating germinal center B cells. Higher levels of antigen-specific B cells set the stage for more robust antibody responses. Whether this higher-level activation of B cells leads to better mAb cloning remains to be determined.

The DNA prime-protein boost approach has been used to generate mAbs in both mouse and rabbit models. This approach was effective in generating a panel of mAbs that are protective against *Clostridium difficile* toxin A challenge, as well as mAbs for use as sensitive reagents to detect toxin A in various testing samples.<sup>39</sup> We also reported the use of DNA immunization to generate rabbit mAbs with a high affinity for and a diverse epitope-binding profile to the human immunodeficiency virus type 1 (HIV-1) envelope antigen.<sup>37</sup> One unique rabbit mAb targeted an area on the envelope protein of HIV-1, which blocks the binding of CD4 and the HIV-1 receptor.<sup>67</sup> Another research group used a similar DNA prime-protein boost approach to generate a higher antibody titer and higher quality mAbs than those observed with protein immunization alone.<sup>36</sup>



#### **ROLE OF FINAL BOOST IMMUNIZATION**

Production of mAbs by the traditional hybridoma method requires the availability of many activated antigen-specific B cells in lymphoid organs for fusion. With the traditional protein immunization approach, this is achieved by a final intravenous or intraperitoneal injection 3–5 days before fusion. This same procedure may also be needed for DNA immunization.

Hybridoma fusions using cells from DNA-immunized animal hosts without a final boost have been reported by several groups. <sup>38,41,44,46,55</sup> However, the overall fusion efficiency was low, and the resulting antibodies had low binding affinity, with IgM as the dominant isotype. One study compared the ability of a final protein boost with no final protein boost with respect to hybridoma generation and concluded that despite significant antibody responses in the immunized animals, the fusion of mouse spleen cells yielded a low number of and low-quality hybridomas unless the mice were given a boost 3–5 days before fusion. <sup>21</sup>

Other studies included an additional DNA plasmid immunization by intramuscular or intradermal injection 3–5 days before fusion as a final boost. 42,45,51 Although the numbers of mAbs generated were small (that is, only a few mAbs from each fusion), mAbs with good binding affinity and diversity were reported. 45,51 The final DNA plasmid boost could also be delivered by hydrodynamic injection five days before fusion, and specific mAbs were successfully generated, 25 including some against very difficult targets, such as multi-transmembrane proteins. 27,30

Proteins are commonly used as final boost reagents. Purified protein has been used successfully for the final protein boost for both secretory proteins<sup>47,56</sup> and intracellular proteins.<sup>54,57,58</sup> For single-transmembrane proteins or glycophosphatidylinositol-anchored proteins, a purified extracellular domain can be used for the final boost if it is confirmed that the proteins retain a native conformation.<sup>34,37,43</sup> Furthermore, cells expressing antigen proteins have been used directly as the final boost. This approach was successful for both membrane proteins<sup>19–21,23,24,35,48,49</sup> and intracellular proteins.<sup>25,28,31,32</sup> The use of cells as a boosting reagent can work for secretory proteins that are difficult to purify by adding a glycophosphatidylinositol anchor. In the case of viral antigens, inactive viral particles have been used successfully as the final boost.<sup>51–53</sup> Table 4 provides a summary of final boost immunization options.

DNA immunization can be used to generate mAbs by either a traditional hybridoma approach or the single B-cell cloning approach in different animal models. Our research group has recently produced mAbs from human volunteers who were immunized via an HIV vaccine DNA prime-protein boost regimen using a single B-cell cloning method (paper in preparation). These results show that DNA

Table 4 Options for the final boost

Type of boost	Original types of proteins	References
DNA by IM or ID	Any type	42,45,51
DNA by HTV	Any type	25,27,30
Protein	Secretory	39,47,56
	Intracellular	54,57,58
	Single-transmembrane/GPI anchored	34,37,43
Cells expressing immunogen	Transmembrane	19–21,23,24,35,48,49
•	Intracellular	25,28,31,32
Viral particle	Viral surface	51–53

Abbreviations: hydrodynamic tail vein, HTV; intradermal, ID; intramuscular, IM.

immunization can be used in a wide range of hosts to produce high-quality mAbs.

#### **SUMMARY**

DNA immunization is a powerful approach to producing high-quality mAbs and offers several unique advantages: (i) DNA immunization is an efficient method of testing different immunogen designs; (ii) DNA immunization does not require the production or purification of proteins from a pathogen, which avoids any issues related to biosafety; (iii) DNA immunization allows a rapid response to an emerging infectious agent once the pathogen gene sequence is known; (iv) DNA immunization is effective in generating mAbs against conformation-sensitive targets; (v) DNA immunization can be used for mAb induction in a wide range of hosts, including mouse, rabbit and human.

#### **ACKNOWLEDGEMENTS**

This work is partially supported by NIH grants HIVRAD U19 AI082274 and R21/R33 AI0879191.

- 1 Lyon GM, Mehta AK, Varkey JB et al. Clinical care of two patients with Ebola virus disease in the United States. N Engl J Med 2014; 371: 2402–2409.
- 2 Ulmer JB, Donnelly JJ, Parker SE et al. Heterologous protection against influenza by injection of DNA encoding a viral protein. Science 1993; 259: 1745–1749.
- 3 Fynan EF, Webster RG, Fuller DH et al. DNA vaccines: protective immunizations by parenteral, mucosal, and gene-gun inoculations. Proc Natl Acad Sci USA 1993; 90: 11478–11482.
- 4 Wang B, Ugen KE, Srikantan V et al. Gene inoculation generates immune responses against human immunodeficiency virus type 1. Proc Natl Acad Sci USA 1993; 90: 4156–4160.
- 5 Wolff JA, Malone RW, Williams P et al. Direct gene transfer into mouse muscle in vivo. Science 1990; 247: 1465–1468.
- 6 Jones S, Evans K, McElwaine-Johnn H et al. DNA vaccination protects against an influenza challenge in a double-blind randomised placebo-controlled phase 1b clinical trial. Vaccine 2009; 27: 2506–2512.
- 7 Kalams SA, Parker SD, Elizaga M et al. Safety and comparative immunogenicity of an HIV-1 DNA vaccine in combination with plasmid interleukin 12 and impact of intramuscular electroporation for delivery. J Infect Dis 2013; 208: 818–829.
- 8 Lu S. Heterologous prime-boost vaccination. Curr Opin Immunol 2009; 21: 346-351.
- 9 Wang S, Kennedy JS, West K et al. Cross-subtype antibody and cellular immune responses induced by a polyvalent DNA prime-protein boost HIV-1 vaccine in healthy human volunteers. Vaccine 2008; 26: 3947–3957.
- 10 Vaine M, Wang S, Liu Q et al. Profiles of human serum antibody responses elicited by three leading HIV vaccines focusing on the induction of Env-specific antibodies. PLoS One 2010: 5: e13916.
- 11 Ledgerwood JE, Wei CJ, Hu Z *et al.* DNA priming and influenza vaccine immunogenicity: two phase 1 open label randomised clinical trials. *Lancet Infect Dis* 2011; **11**: 916–924.
- 12 Vaine M, Wang S, Hackett A et al. Antibody responses elicited through homologous or heterologous prime-boost DNA and protein vaccinations differ in functional activity and avidity. Vaccine 2010; 28: 2999–3007.
- 13 Wang S, Parker C, Taaffe J et al. Heterologous HA DNA vaccine prime-inactivated influenza vaccine boost is more effective than using DNA or inactivated vaccine alone in eliciting antibody responses against H1 or H3 serotype influenza viruses. Vaccine 2008; 26: 3626–3633.
- 14 Suguitan Jr AL, Cheng X, Wang W *et al.* Influenza H5 hemagglutinin DNA primes the antibody response elicited by the live attenuated influenza A/Vietnam/1203/2004 vaccine in ferrets. *PLoS One* 2011; **6**: e21942.
- 15 Pal R, Wang S, Kalyanaraman VS et al. Polyvalent DNA prime and envelope protein boost HIV-1 vaccine elicits humoral and cellular responses and controls plasma viremia in rhesus macaques following rectal challenge with an R5 SHIV isolate. J Med Primatol 2005: 34: 226–236.
- 16 Wang S, Arthos J, Lawrence JM et al. Enhanced immunogenicity of gp120 protein when combined with recombinant DNA priming to generate antibodies that neutralize the JR-FL primary isolate of human immunodeficiency virus type 1. J Virol 2005; 79: 7933–7937.
- 17 Boyer JD, Chattergoon MA, Ugen KE *et al.* Enhancement of cellular immune response in HIV-1 seropositive individuals: A DNA-based trial. *Clin Immunol* 1999; **90**: 100–107.
- 18 Boyer JD, Chattergoon M, Shah A et al. HIV-1 DNA based vaccine induces a CD8 mediated cross-clade CTL response. Dev Biol Stand 1998; 95: 147–153.
- 19 Costagliola S, Rodien P, Many MC et al. Genetic immunization against the human thyrotropin receptor causes thyroiditis and allows production of monoclonal antibodies recognizing the native receptor. J Immunol 1998; 160: 1458–1465.



- 20 Costagliola S, Franssen JD, Bonomi M et al. Generation of a mouse monoclonal TSH receptor antibody with stimulating activity. Biochem Biophys Res Commun 2002; 299: 891–896
- 21 Nagata S, Salvatore G, Pastan I. DNA immunization followed by a single boost with cells: a protein-free immunization protocol for production of monoclonal antibodies against the native form of membrane proteins. J Immunol Methods 2003; 280: 59–72.
- 22 Carson SD, Switzer BL, Tracy SM et al. Monoclonal antibody against mouse CAR following genetic immunization. Hybrid Hybridomics 2004: 23: 19–22.
- 23 Costagliola S, Bonomi M, Morgenthaler NG et al. Delineation of the discontinuous-conformational epitope of a monoclonal antibody displaying full in vitro and in vivo thyrotropin activity. Mol Endocrinol 2004: 18: 3020–3034.
- 24 Adriouch S, Dubberke G, Diessenbacher P et al. Probing the expression and function of the P2X7 purinoceptor with antibodies raised by genetic immunization. Cell Immunol 2005: 236: 72–77.
- 25 Bates MK, Zhang G, Sebestyen MG et al. Genetic immunization for antibody generation in research animals by intravenous delivery of plasmid DNA. Biotechniques 2006; 40: 199–208
- 26 Fofana I, Krieger SE, Grunert F et al. Monoclonal anti-claudin 1 antibodies prevent hepatitis C virus infection of primary human hepatocytes. Gastroenterology 2010; 139: 953–964. 64 e1-4.
- 27 Takatsuka S, Sekiguchi A, Tokunaga M et al. Generation of a panel of monoclonal antibodies against atypical chemokine receptor CCX-CKR by DNA immunization. J Pharmacol Toxicol Methods 2011; 63: 250–257.
- 28 Daftarian P, Chowdhury R, Ames P *et al. In vivo* electroporation and non-protein based screening assays to identify antibodies against native protein conformations. *Hybridoma* (*Larchmt*) 2011; **30**: 409–418.
- 29 Fofana I, Xiao F, Thumann C et al. A novel monoclonal anti-CD81 antibody produced by genetic immunization efficiently inhibits Hepatitis C virus cell-cell transmission. PLoS One 2013; 8: e64221.
- 30 Hazen M, Bhakta S, Vij R et al. An improved and robust DNA immunization method to develop antibodies against extracellular loops of multi-transmembrane proteins. MAbs 2014: 6: 95–107.
- 31 Garcia JF, Maestre L, Lucas E et al. Genetic immunization: a new monoclonal antibody for the detection of BCL-6 protein in paraffin sections. J Histochem Cytochem 2006; 54: 31–38
- 32 Maestre L, Fontan L, Martinez-Climent JA *et al.* Generation of a new monoclonal antibody against MALT1 by genetic immunization. *Hybridoma (Larchmt)* 2007; **26**: 86–91.
- 33 Kilpatrick KE, Cutler T, Whitehorn E et al. Gene gun delivered DNA-based immunizations mediate rapid production of murine monoclonal antibodies to the Flt-3 receptor. Hybridoma 1998; 17: 569–576.
- 34 Wongsena W, Sconocchia G, Cho HS *et al.* Production and characterization of monoclonal antibodies against major histocompatibility complex class I chain-related gene A. *Tissue Antigens* 2008; **72**: 431–440.
- 35 Brady JL, Mannering SI, Kireta S et al. Monoclonal antibodies generated by DNA immunization recognize CD2 from a broad range of primates. *Immunol Cell Biol* 2009; 87: 413–418.
- 36 Xue T, Xue XP, Huang QS *et al.* Monoclonal antibodies against human aspartyl (asparaginyl) beta-hydroxylase developed by DNA immunization. *Hybridoma (Larchmt)* 2009; **28**: 251–257.
- 37 Chen Y, Vaine M, Wallace A et al. A novel rabbit monoclonal antibody platform to dissect the diverse repertoire of antibody epitopes for HIV-1 Env immunogen design. I Virol 2013: 87: 10232–10243
- 38 Ulivieri C, Burroni D, Telford JL *et al.* 1996. Generation of a monoclonal antibody to a defined portion of the Helicobacter pylori vacuolating cytotoxin by DNA immunization. *J Biotechnol* 1996; **51**: 191–194.
- 39 Zhang C, Jin K, Xiao Y et al. Potent monoclonal antibodies against Clostridium difficile toxin A elicited by DNA immunization. Hum Vaccin Immunother 2013; 9: 2157–2164.
- 40 Yu XF, Liang LH, She M et al. Production of a monoclonal antibody against SARS-CoV spike protein with single intrasplenic immunization of plasmid DNA. *Immunol Lett* 2005; 100: 177–181.
- 41 Puttikhunt C, Kasinrerk W, Srisa-ad S et al. Production of anti-dengue NS1 monoclonal antibodies by DNA immunization. J Virol Methods 2003; 109: 55–61.
- 42 Pupa SM, Forti S, Invernizzi AM *et al.* Monoclonal antibody to fibulin-1 generated by genetic immunization. *J Cell Biochem* 2003; **89**: 647–652.
- 43 Koch-Nolte F, Glowacki G, Bannas P et al. Use of genetic immunization to raise antibodies recognizing toxin-related cell surface ADP-ribosyltransferases in native conformation. Cell Immunol 2005; 236: 66–71.
- 44 Chen Y, Zhang T, Li T et al. Preparation and characterization of a monoclonal antibody against CKLF1 using DNA immunization with in vivo electroporation. Hybridoma (Larchmt) 2005; 24: 305–308.

- 45 Leinonen J, Niemela P, Lovgren J et al. Characterization of monoclonal antibodies against prostate specific antigen produced by genetic immunization. J Immunol Methods 2004; 289: 157–167.
- 46 Cantelli CP, da Gloria Martins Teixeira M, Santos EA *et al.* Generation of monoclonal antibodies against human recombinant interferon beta using genetic immunization with simultaneous expression of IgM and IgG isotypes. *Hybridoma (Larchmt)* 2009; **28**: 211–214.
- 47 Bian C, Zhang F, Wang F *et al.* Development of retinol-binding protein 4 immunocolloidal gold fast test strip using high-sensitivity monoclonal antibodies generated by DNA immunization. *Acta Biochim Biophys Sin (Shanghai)* 2010; **42**: 847–853.
- 48 Schmolke S, Tacke M, Schmitt U et al. Identification of hepatitis G virus particles in human serum by E2-specific monoclonal antibodies generated by DNA immunization. J Virol 1998; 72: 4541–4545.
- 49 Horimoto T, Fukuda N, Iwatsuki-Horimoto K et al. Antigenic differences between H5N1 human influenza viruses isolated in 1997 and 2003. J Vet Med Sci 2004; 66: 303–305
- 50 Ni Y, Ma K, Ni J et al. A rapid and simple approach to preparation of monoclonal antibody based on DNA immunization. Cell Mol Immunol 2004; 1: 295–299.
- 51 Martin-Folgar R, Lorenzo G, Boshra H et al. Development and characterization of monoclonal antibodies against Rift Valley fever virus nucleocapsid protein generated by DNA immunization. MAbs 2010; 2: 275–284.
- 52 Wei J, Yan B, Chen Z *et al.* Production and characterization of monoclonal antibodies against the hemagglutinin of H5N1 and antigenic investigation of avian influenza H5N1 viruses isolated from China. *Can J Microbiol* 2011; **57**: 42–48.
- 53 Tan GS, Krammer F, Eggink D *et al.* A pan-H1 anti-hemagglutinin monoclonal antibody with potent broad-spectrum efficacy in vivo. *J Virol* 2012; **86**: 6179–6188.
- 54 Velikovsky CA, Cassataro J, Sanchez M et al. Single-shot plasmid DNA intrasplenic immunization for the production of monoclonal antibodies. Persistent expression of DNA. J Immunol Methods 2000; 244: 1–7.
- 55 Al-Qaoud KM, Al-Omari MM, Al-Aghbar M et al. Production of monoclonal antibodies against the 8 kDa subunit of Echinococcus granulosus Antigen B (EgAgB8/2) using DNA immunization. Hybridoma (Larchmt) 2008; 27: 431–438.
- 56 Oynebraten I, Lovas TO, Thompson K et al. Generation of antibody-producing hybridomas following one single immunization with a targeted DNA vaccine. Scand J Immunol 2012; 75: 379–388.
- 57 Kilpatrick KE, Danger DP, Hull-Ryde EA *et al.* High-affinity monoclonal antibodies to PED/PEA-15 generated using 5 microg of DNA. *Hybridoma* 2000; **19**: 297–302.
- 58 Kilpatrick KE, Kerner S, Dixon EP *et al.* In vivo expression of a GST-fusion protein mediates the rapid generation of affinity matured monoclonal antibodies using DNA-based immunizations. *Hybrid Hybridomics* 2002; **21**: 237–243.
- 59 Wang S, Farfan-Arribas DJ, Shen S *et al.* Relative contributions of codon usage, promoter efficiency and leader sequence to the antigen expression and immunogenicity of HIV-1 Env DNA vaccine. *Vaccine* 2006; **24**: 4531–4540.
- 60 Groisman A, Simonnet C, Rinberg D Pneumatic capillary gun for ballistic delivery of microscopic particles into tissue. US Patent 7892836 B2, February 22, 2011.
- 61 Robinson HL, Torres CA. DNA vaccines. Semin Immunol 1997; 9: 271-283.
- 62 Wang S, Joshi S, Lu S. Delivery of DNA to skin by particle bombardment. *Methods Mol Biol* 2004; **245**: 185–196.
- 63 Alexandrenne C, Wijkhuisen A, Dkhissi F *et al.* Generating antibodies against the native form of the human prion protein (hPrP) in wild-type animals: a comparison between DNA and protein immunizations. *J Immunol Methods* 2009; **341**: 41–49.
- 64 Wang S, Zhang C, Zhang L et al. The relative immunogenicity of DNA vaccines delivered by the intramuscular needle injection, electroporation and gene gun methods. Vaccine 2008: 26: 2100–2010.
- 65 Fujimoto A, Kosaka N, Hasegawa H et al. Enhancement of antibody responses to native G protein-coupled receptors using E. coli GroEL as a molecular adjuvant in DNA immunization. J Immunol Methods 2012; 375: 243–251.
- 66 Li W, Wang S, Lu S. Pilot Study on the Use of DNA Priming Immunization to Enhance Y. pestis LcrV-Specific B Cell Responses Elicited by a Recombinant LcrV Protein Vaccine. *Vaccines (Basel)* 2013; **2:** 36–48.
- 67 Pan R, Chen Y, Vaine M et al. Structural analysis of a novel rabbit monoclonal antibody R53 targeting an epitope in HIV-1 gp120 C4 region critical for receptor and co-receptor binding. Emerg Microbes Infect 2015; 4: e44.



This work is licensed under a Creative Commons Attribution 4.0 International License. The images or other third party material in this

article are included in the article's Creative Commons license, unless indicated otherwise in the credit line; if the material is not included under the Creative Commons license, users will need to obtain permission from the license holder to reproduce the material. To view a copy of this license, visit http://creativecommons.org/licenses/by/4.0/