



Article LAMP Assay for the Detection of *Echinococcus multilocularis* Eggs Isolated from Canine Faeces by a Cost-Effective NaOH-Based DNA Extraction Method

Barbara J. Bucher¹, Gillian Muchaamba¹, Tim Kamber¹, Philipp A. Kronenberg^{1,2}, Kubanychbek K. Abdykerimov^{3,4}, Myktybek Isaev⁵, Peter Deplazes¹ and Cristian A. Alvarez Rojas^{1,*}

- ¹ Institute of Parasitology, Vetsuisse and Medical Faculty, University of Zurich, 8057 Zurich, Switzerland; barbara.bucher@uzh.ch (B.J.B.); gillian.muchaamba@uzh.ch (G.M.); timpeter.kamber@uzh.ch (T.K.); philipp.kronenberg2@uzh.ch (P.A.K.); deplazesp@access.uzh.ch (P.D.)
- ² Graduate School for Cellular and Biomedical Sciences, University of Bern, 3012 Bern, Switzerland
- ³ Section of Epidemiology, Vetsuisse Faculty, University of Zurich, 8057 Zurich, Switzerland; kubanychbek.abdykerimov@uzh.ch
- ⁴ Life Science Zurich Graduate School, University of Zurich, 8057 Zurich, Switzerland
- ⁵ Department of Parasitology, Kyrgyz Research Institute of Veterinary Medicine Arstanbek Duisheev, Ministry of Education and Science of the Kyrgyz Republic, Bishkek 720033, Kyrgyzstan; isaev-ww-1988@mail.ru
- Correspondence: cristian.alvarezrojas@uzh.ch

Abstract: The detection of Echinococcus multilocularis in infected canids and the environment is pivotal for a better understanding of the epidemiology of alveolar echinococcosis in endemic areas. Necropsy/sedimentation and counting technique remain the gold standard for the detection of canid infection. PCR-based detection methods have shown high sensitivity and specificity, but they have been hardly used in large scale prevalence studies. Loop-mediated isothermal amplification (LAMP) is a fast and simple method to detect DNA with a high sensitivity and specificity, having the potential for field-application. A specific LAMP assay for the detection of E. multilocularis was developed targeting the mitochondrial nad1 gene. A crucial step for amplification-based detection methods is DNA extraction, usually achieved utilising silica-gel membrane spin columns from commercial kits which are expensive. We propose two cost-effective and straightforward methods for DNA extraction, using NaOH (method 1A) and InstaGeneTM Matrix (method 1B), from isolated eggs circumventing the need for commercial kits. The sensitivity of both assays with fox samples was similar (72.7%) with multiplex-PCR using protocol 1A and LAMP using protocol 1B. Sensitivity increased up to 100% when testing faeces from 12 foxes infected with more than 100 intestinal stages of E. multilocularis. For dogs, sensitivity was similar (95.4%) for LAMP and multiplex-PCR using protocol 1B and for both methods when DNA was extracted using protocol 1A (90.9%). The DNA extraction methods used here are fast, cheap, and do not require a DNA purification step, making them suitable for field studies in low-income countries for the prevalence study of E. multilocularis.

Keywords: Echinococcus; NaOH; LAMP; PCR; DNA extraction; taeniid egg isolation

1. Introduction

Echinococcus multilocularis is the canine intestinal cestode responsible for alveolar echinococcosis (AE) in humans [1]. Primarily foxes, and to a lesser extent, other wild canids and domestic dogs play a role as definitive hosts for *E. multilocularis* [2]. Humans act as dead-end hosts, acquiring the infection by oral intake of viable eggs of *E. multilocularis* from infected food, soil, water, or hand-to-mouth contact after a direct interaction with infected dogs and contaminated matrices [3]. Establishing the exact link between AE and its origin (source attribution) is challenging due to the long incubation time of the disease of up to 15 years. Therefore, it is essential to have tools for the accurate determination of



Citation: Bucher, B.J.; Muchaamba, G.; Kamber, T.; Kronenberg, P.A.; Abdykerimov, K.K.; Isaev, M.; Deplazes, P.; Alvarez Rojas, C.A. LAMP Assay for the Detection of *Echinococcus multilocularis* Eggs Isolated from Canine Faeces by a Cost-Effective NaOH-Based DNA Extraction Method. *Pathogens* 2021, 10, 847. https://doi.org/10.3390/ pathogens10070847

Academic Editor: Vito Colella

Received: 20 May 2021 Accepted: 29 June 2021 Published: 5 July 2021

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2021 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). the environmental contamination with *E. multilocularis* eggs and the detection of canine intestinal infections.

Different methodologies for the diagnosis of *E. multilocularis* in canids have been standardised and published over the years. Post-mortem examination (i.e., sedimentation and counting technique, SCT) remains the gold standard [4]. However, a post-mortem examination has significant disadvantages, requiring either expensive infrastructure, strict biosafety measures and involves the killing of canids. Several methods based on the detection of E. multilocularis antigens [5–8] and DNA [9–11] in faeces have been developed (reviewed in [12,13]). Coproantigen detection (copro-ELISA) has the advantage of detecting prepatent infections [5,14]. Although some coproantigen tests are commercially available, they are not regularly used for massive screening in canid faeces from endemic areas. On the other hand, DNA of *E. multilocularis* can be potentially detected from the prepatent period as free DNA or from cells of the growing worms; and mostly from proglottids/eggs during patency. DNA can be isolated from total faeces and used, for example, in nested PCR (copro-DNA-PCR) [9] or from isolated taeniid eggs through floatation and sieving [10] and used in multiplex-PCR (egg-DNA-PCR) [15]. Although the egg-DNA-PCR focuses only on patent infections, it has the advantage of concentrating eggs and decreasing the presence of PCR inhibitors.

PCR-based methods for identifying parasites in faecal samples rely on isolation and DNA concentration using silica-gel membrane in spin columns from commercial kits removing PCR inhibitors (proteinases, bile salts, polyphenols, and acids). However, some publications report that the total removal of inhibitors from faeces is not always possible [16,17]. Furthermore, the usage of commercial kits requires trained staff performing several liquid handling steps, and the commercial kit itself, which can be a costly item when large numbers of samples need to be analysed. The cost of such kits ranges between 3 and 11 USD per sample, which is unreasonable for governments and researchers in low-income countries. Besides the PCR inhibitors, the successful disruption of the typical taeniid eggshell keratin-like layer by homogenisation or alkaline lysis is essential for DNA detection. Unlike polymerases used in PCR, the in silico designed Bst 2.0 DNA Polymerase used in LAMP (loop-mediated isothermal amplification) displays improved tolerance for inhibitors [18-20]. LAMP assays have been developed to detect E. multilocularis [21,22] and E. granulosus [23-26] in faecal samples and also several other parasites, including Clonorchis sinensis [27], Trichuris muris [28], Strongyloides spp. [29], Necator americanus [30], and Taenia spp. [31]. LAMP also offers a robust tool for DNA amplification on the presence of reagents like NaOH used in alternative DNA extraction methods suggesting that commercial kits for DNA isolation could be circumvented. Methods for DNA extraction without silica-gel membranes have been described in the literature since the 1990s [32–35]. A reappraisal of such methods for developing cost-effective DNA extraction (i.e., NaOH-based) has provided evidence to produce a template of sufficient quality for subsequent PCR from different tissues [36–38]. Furthermore, single *Echinococcus* eggs have been used for PCR- and LAMP-based methods directly after treatment with NaOH [26,39-42]. In the present study, we propose two cost-effective DNA extraction methods (without purification) from isolated taeniid eggs, based on NaOH and InstaGeneTM Matrix (Chelex 6%) to be used directly as a template for multiplex-PCR [15] and LAMP (developed in this study) for the detection of E. multilocularis. We also labelled two of the primers used in LAMP to combine them with lateral flow dipstick for the same diagnostic purpose. The proposed cost-effective DNA extractions can be used in large prevalence studies to detect E. multilocularis eggs and potentially *E. granulosus* in canine faecal samples in low-income endemic countries.

2. Results

2.1. Optimisation, Specificity, Analytical Sensitivity, Limit of Detection, and Stability of LAMP

After optimisation of the LAMP assay, the reaction was set up with 6 mM MgSO₄, 0.8 M betaine, 0.004% malachite green, and $1 \times$ primer mix (1.6 μ M of each FIP and BIP primer, 0.2 μ M of each F3 and B3 primer and 0.4 μ M of the loop primer); amplification occurred at 65 °C for 60 min. The selected LAMP primer set showed successful amplification

when using *E. multilocularis* DNA rendering a clear blue-green colour, a ladder-like structure on agarose gel after electrophoresis and a specific band on the lateral-flow-dipstick (when using the labelled BIP and FIP primers). Conversely, a negative result followed by testing DNA of other Echinococcus spp., Taenia spp., and other cestodes as a template (see material and methods). The limit of detection (analytical sensitivity) of LAMP with the selected primer set was 1 pg/ μ L of *E. multilocularis* DNA in all serial dilutions performed in triplicate. On the other hand, the limit of detection of the multiplex-PCR [15] was 10 fg/ μ L of E. multilocularis DNA. Serial dilutions of samples derived from protocol 1A were performed to exclude inhibition of LAMP reactions. The 1:50 dilution showed to perform best, as no inhibition could be detected. The lowest amount of eggs possible to detect with LAMP and multiplex-PCR was one egg treated with the protocol 1B in 4/4 replicates and two eggs in 3/4 replicates with protocol 1A. In order to investigate if free DNA of the parasite was present in the supernatant of the egg suspensions used for the limit of detection we used it as a template for LAMP and multiplex-PCR [15], all reactions were negative. In this way, we can confirm that DNA originated from the eggs in suspension. Aliquots of a prepared master mix for LAMP containing all the reagents except DNA were tested weekly over six weeks. A positive LAMP reaction was achieved in all the period when adding *E. multilocularis* DNA (2 ng).

2.2. Examination of Field Samples with Protocols 1A and 1B

A schematic representation of the sample processing is shown in Figure 1. In the case of negative samples, we collected the sediment from the 21 μ m filter, in the same way as in positive samples, to perform the DNA extraction.




Figure 1. Schematic representation of the process of field samples analysed in the present study. Step 7 shows the two methods used for DNA extraction in this study (protocols 1A and 1B). The visualisation of multiplex-PCR occurred after electrophoresis and in the case of LAMP results were assessed after change of colour of the reaction, electrophoresis, and lateral flow dipsticks as shown in step 9.

2.2.1. Foxes

Intestinal stages of *E. multilocularis* (between 2 to >100) were detected in 44/62 foxes after necropsy/SCT. Microscopic detection of taeniid eggs was successful for 30 samples from the 44 foxes documented to be infected with *E. multilocularis*. On the other hand, 18 faecal samples were collected from foxes without intestinal stages of *E. multilocularis* at necropsy/SCT; from them, it was possible to observe taeniid eggs in five samples (Table 1); 8/18 foxes mentioned before were infected with *Taenia* spp. When using protocol 1A, LAMP showed a positive reaction in 26/44 samples from animals harbouring *E. multilocularis* in their intestine and 1/18 samples from foxes without intestinal stages of *E. multilocularis*. Conversely, multiplex-PCR was positive for *E. multilocularis* in 32/44 samples from foxes infected with *E. multilocularis*. Using multiplex-PCR, other cestodes and *Taenia* spp. were identified, producing an amplicon (267 bp) in five samples from the 18 animals without intestinal stages of *E. multilocularis* (protocol 1A and 1B), suggesting that if PCR inhibitors were present, they did not impede the amplification of DNA (See detailed list of results in Table S1).

The second aliquot of faeces from the same foxes was used to isolate eggs for treatment with protocol 1B. In this case, taeniid eggs were observed in 31/44 faeces from animals with *E. multilocularis* in their intestine and in 5/18 faecal samples from negative animals. After treating eggs with protocol 1B, we found positive LAMP reactions in 32/44 foxes infected with *E. multilocularis* and in 1/18 foxes with no intestinal stages of *E. multilocularis*. Conversely, multiplex-PCR was positive for *E. multilocularis* in 30/44 positive and 1/18 negative foxes, respectively. A positive result for the amplification of 267 bp corresponding to other cestodes including *Taenia* spp. was found in 10 samples with a negative result for *E. multilocularis* using the DNA from protocol 1B (Table S1).

Furthermore, multiplex-PCR using protocols 1A and 1B showed a positive result for the amplification of 267 bp corresponding to other cestodes, including *Taenia* spp. in all the eight foxes identified to be infected with *Taenia* spp. at necropsy/SCT.

2.2.2. Dogs

The detection of *E. multilocularis* in dogs was performed by a multiplex-PCR in a previous prevalence study. In the case of the dog group (Table 1), taeniid eggs were observed in 20/22 samples from dogs previously identified to be infected with *E. multilocularis*. No eggs were observed in samples that were considered negative to *E. multilocularis*. When the isolated eggs were treated with protocol 1A and used for DNA amplification, we found a positive result for *E. multilocularis* in 20/22 samples using LAMP and multiplex-PCR. No positive reactions occurred with LAMP and multiplex-PCR in the ten samples previously identified as negative to *E. multilocularis*. When isolating taeniid eggs to be treated with protocol 1B, we observed eggs in 21/22 samples from dogs previously identified positive to *E. multilocularis*. In this group, LAMP and multiplex-PCR showed a positive result in 21/22 samples and a negative result in 10/10 samples identified as negative to *E. multilocularis*.

In the case of fox samples, sensitivity in LAMP and multiplex-PCR was 72.7% (C.I. 95%: 57.2–85) using protocol 1B and 1A, respectively (Table 2). Specificity values were not 100% for fox samples, the highest value was 94.4% (72.7–99.8) achieved using both protocols for LAMP and multiplex-PCR (Table 2). In the dog group, sensitivity for LAMP and multiplex-PCR was 95.4% (77.1–99.8) achieved with protocol 1B. The specificity in LAMP and multiplex-PCR with both protocols was 100% (Table 2).

Table 1. Results of LAMP and multiplex-PCR¹ amplification of *E. multilocularis* using the product of two DNA extraction methods (protocols 1A and 1B) as a template from isolated taeniid eggs from 62 fox faecal samples and 32 dog faecal samples. In the case of foxes, the status of *E. multilocularis* infection was established at necropsy with the Sedimentation and Counting Test (SCT); for dogs, *E. multilocularis* infections status was established based on egg isolation/multiplex-PCR [15].

	Number of Animals Confirmed	Protocol 1A: NaOH/Dilution 1:50 with Tris-HCl			Protocol 1B: NaOH + InstaGene TM Matrix		
Host	Positive (+) or Without (–) Detected <i>E. multilocularis</i> Infection (Method for Diagnosis)	Microscopic Detection of Taeniid Eggs ²	LAMP Positive	Multiplex-PCR ¹ Positive	Microscopic Detection of Taeniid Eggs ²	LAMP Positive	Multiplex-PCR ¹ Positive
	44 + (necropsy/SCT)	30 positive	25	29	31 positive	28	27
Easter		14 negative	1	3	13 negative	4	3
Foxes	18 – (necropsy/SCT)	5 positive	1	1	5 positive	1	1
		13 negative	0	0	13 negative	0	0
	22 + (multiplex-PCR ¹)	20 positive	20	20	21 positive	21	21
Dogs		2 negative	0	0	1 negative	0	0
0	$10 - (multiplex-PCR^{1})$	10 negative	0	0	10 negative	0	0

¹ as described in [15]; ² independent egg isolation from two aliquots of two grams of faeces performed as previously described [10].

Table 2. Sensitivity and specificity values including confidence intervals in brackets for LAMP and multiplex-PCR for the detection of *E. multilocularis* using as template material treated with protocols 1A and 1B (see Table 1 for details) from foxes and dogs. Foxes were considered truly infected if intestinal stages of *E. multilocularis* were found in the intestine at necropsy/SCT. Dogs were considered as truly infected based on a positive multiplex-PCR [15] performed previously as part of a prevalence study in Kyrgyzstan.

	Foxes Protocol 1A		Foxes Pr	exes Protocol 1B Dogs		otocol 1A	Dogs Protocol 1B	
	LAMP	Multiplex- PCR	LAMP	Multiplex- PCR	LAMP	Multiplex- PCR	LAMP	Multiplex- PCR
Sensitivity	59	72.7	72.7	68.1	90.9	90.9	95.4	95.4
(IC 95%)	(43.2–73.6)	(57.2–85)	(57.2–85)	(52.4–81.3)	(70.8 - 98.8)	(70.8–98.8)	(77.1–99.8)	(77.1–99.8)
Specificity	94.4	94.4	94.4	94.4	100	100	100	100
(IC 95%)	(72.7–99.8)	(72.7–99.8)	(72.7–99.8)	(72.7–99.8)	(69.1–100)	(69.1–100)	(69.1–100)	(69.1–100)

2.3. Examination of Field Samples with Protocols 1A and 1B According to Worm Burden

From 28 fox samples (positive at necropsy/SCT) collected in 2020, we assessed the worm burden for *E. multilocularis* after a thorough examination of the intestinal mucosa. Worm burden ranged between 2 and >100 *E. multilocularis* intestinal stages (Table 3). We divided the samples based on worm burden from foxes with two to 20 intestinal stages of *E. multilocularis* (9 foxes), 21–100 (7) and samples from foxes with >100 worms (12). We also included 17 samples from foxes that were negative at necropsy/SCT for *E. multilocularis*. LAMP and multiplex-PCR showed a positive result in 1/9 samples with worm burden between two and 20 parasites. When worm burden was between 21–100 worms, LAMP could detect 4/7 of the samples with protocol 1A and 6/7 were positive with multiplex-PCR. With protocol 1B, LAMP detected 6/7, and multiplex-PCR found 7/7 to be positive. Finally, when the worm burden was higher than 100 worms LAMP and multiplex-PCR detected 12/12 positive samples with protocol 1B. Multiplex-PCR detected 12/12 samples positives with protocol 1A (Table 3). Finally, from the 17 samples with no *E. multilocularis* found at necropsy/SCT, LAMP and multiplex-PCR rendered a positive result in one sample treated with both protocols.

Table 3. Number of positive results for *E. multilocularis* in LAMP and multiplex-PCR using protocols 1A and 1B in faecal samples from foxes necropsied in 2020 related to the total worm burden for *E. multilocularis*.

		Protocol 1A: NaOH/Dilution 1:50 Tris-HCl		Protocol 1B: NaOH + InstaGene TM Matrix		
# E. multilocularis at SCT	# Examined	LAMP Positive (Sensitivity)	Multiplex-PCR Positive (Sensitivity)	LAMP Positive (Sensitivity)	Multiplex-PCR Positive (Sensitivity)	
0	17	1	1	1	1	
2-20	9	1 (11.1%)	1 (11.1%)	2 (18.1%)	2 (18.1%)	
21-100	7	4 (57.1%)	6 (85.7%)	6 (85.7%)	7 (100%)	
>100	12	10 (83.3%)	12 (100%)	12 (100%)	12 (100%)	

2.4. DNA Extraction from Whole Faeces (Protocol 2)

We investigated if it was possible to use NaOH+InstaGeneTM Matrix (called protocol 2A) to treat a faecal sample to be used directly as a template for LAMP reaction (from this study) and a multiplex-PCR [15]. For comparison, we isolated DNA from the same samples using the Qiagen mini stool kit (named here protocol 2B). Using protocol 2A, we found a positive LAMP reaction for *E. multilocularis* in 10/30 faeces from foxes positive to *E. multilocularis* at necropsy/SCT; and a positive result in 17/30 in multiplex-PCR (Table 4). When an aliquot of the same samples was used for DNA isolation with the commercial kit, LAMP was positive in 16/30 and multiplex-PCR in 20/30 samples. Furthermore, a positive amplification corresponding to other cestodes, including *Taenia* spp. [15] was found in seven and five samples which were negative to *E. multilocularis* using a template from protocol 2A and the QIAGEN stool kit, respectively, with the multiplex-PCR. Detailed record of the multiplex-PCR results showing amplification of other cestodes, including *Taenia* spp. can be found in Table S1.

In the case of dogs, LAMP was positive in 2/18 samples and 7/18 samples in multiplex-PCR with protocol 2A. Using DNA isolated with a commercial kit, we found a positive LAMP reaction in 7/18 samples and 6/18 using the multiplex-PCR. Overall, the sensitivity achieved with the multiplex PCR using protocol 2B was 66% (CI: 47.1–82.7) followed by multiplex-PCR using protocol 2A (56.6%, 37.4–74.5) both in fox faeces. For dog faeces, the sensitivity achieved was 38.8% (17.3–64.2) with multiplex-PCR using protocol 2A and LAMP using protocol 2B.

3. Discussion

Loop-mediated isothermal amplification (LAMP) [43] has been portrayed as an affordable alternative to PCR for the detection of different pathogens, including parasites in faeces [18,19,44,45]. *Bst* polymerase used in LAMP offers high robustness being able to withstand harsh chemicals like NaOH and to overcome the presence of PCR inhibitors. In this study, we used the *Bst* 2.0 DNA polymerase which has improved performances regarding amplification speed, yield, tolerance to salt, and inhibitors present in different matrices including faeces [18–20]. Nevertheless, in the present study, a dilution of samples was necessary for protocol 1A, either to dilute NaOH or inhibitors present in samples. We developed a simple, cost-efficient methodology for DNA extraction involving lysis of isolated taeniid eggs with NaOH alone or adding Chelex resin (InstaGeneTM Matrix). Then, we coupled the DNA extraction with a LAMP test developed in this study. Additionally, we included the multiplex-PCR [15] to compare the results for the amplification of the *E. multilocularis* target.

Necropsy/SCT remains the gold standard for the detection of *Echinococcus* spp. in canids but it is not rational to perform. We had access to freshly hunted foxes from which a necropsy and SCT could be performed. Therefore, we were able to compare the performance of the methods for DNA extraction using LAMP and multiplex-PCR [15] to detect *E. multilocularis* in fox faecal samples against the gold standard. As seen in Table 2, the sensitivity of both assays with fox samples was 72.7 (57.2–85) with multiplex-PCR using protocol 1A and LAMP using protocol 1B. We performed a thorough worm count in the intestinal content of foxes necropsied in 2020, showing that the sensitivity of the tests increased up to 100% in multiplex-PCR using protocol 1A for DNA extraction in faeces from 12 foxes with more than 100 worms. Protocol 1B LAMP and multiplex-PCR also detected 12/12 positive samples. It is possible to suggest that sensitivity could be improved further using a higher amount of faeces as starting material to increase the chances of finding eggs in animals with low worm burden.

In samples from dogs, the infection status was assessed with multiplex-PCR from isolated eggs from faeces collected from the ground in a highly endemic area for AE in Kyrgyzstan. In this case, sensitivity was 95.4% for LAMP and multiplex-PCR using protocol 1B (Table 2) and 90.9% for both methods when DNA was extracted using protocol 1A. In the present study, we were able to detect other cestodes, including *Taenia* spp. and *Mesocestoides* spp. in some of the samples which were negative for *E. multilocularis* using the multiplex PCR developed by Trachsel et al. [15], suggesting that no inhibition was present in these samples. Furthermore, we detected Taenia in all the samples from naturally infected foxes harbouring adult Taenia specimens in their intestine. By using whole faeces (400–500 mg) as starting material, multiplex-PCR was more sensitive than LAMP detecting E. multilocularis infection in foxes and dogs that using DNA isolated with a commercial kit. We used a large amount of faecal sample (500 mg) as input for the QIA amp FAST DNA Stool Mini Kit. Using the same commercial kit, Skrzypek et al. [46] used 1 g of faeces to detect E. multilocularis reporting positive results in 45.7% and 48.6% of the faeces from infected foxes (diagnosed at necropsy/SCT) with nested and multiplex-PCR, respectively, which is similar to the results found by us.

II.e.t	Number of Animals Confirmed Positive (+)	Protocol 2A: Na	OH + InstaGene TM Matrix	Protocol 2B: QIAamp DNA Stool Mini Kit		
Host	or Without (-) Detected E. multilocularis infection (Method for Diagnosis)	LAMP Positive	Multiplex-PCR ¹ Positive	LAMP Positive	Multiplex-PCR ¹ Positive	
Foxes	30 + (necropsy/SCT)	10 ²	17	16 ²	20	
	5 - (necropsy/SCT)	0	0	0	0	
Dogs	18 + (multiplex-PCR)	2 ²	7	7	6	
	10 - (multiplex-PCR)	0	0	0	0	

Table 4. Results of LAMP for *E. multilocularis* (from this study) and multiplex-PCR [15] using as a template the supernatant of 500 mg of fox (from Switzerland) and dog faeces (from Kyrgyzstan) treated with protocol 2A and from the DNA isolated with the QIAGEN stool kit (protocol 2B).

¹ as described in [15]; ² all samples found to be positive in LAMP were also positive in multiplex-PCR.

9 of 16

Previously, a LAMP assay has been used to detect *E. multilocularis* in canine faeces in China, reporting higher sensitivity than PCR in experimentally and naturally infected dogs, detecting as low as 1 pg DNA [21] which is the same limit of detection in the LAMP from the present study. The LAMP assay developed by Ni et al. has been used for detecting *E. multilocularis* [21] in wastewater in China [47]. Interestingly, in the work by Ni et al. [48], DNA was isolated with the QIAamp DNA Stool Mini Kit which in our hands did not produce satisfactory results. However, the results using the commercial kit were better than using NaOH-based method directly in faeces. The same kit was used in another report of a LAMP assay for the detection of *E. multilocularis* and other taeniids in a Tibetan rural area [22]. The LAMP assay detected 1 pg of DNA from *E. multilocularis* and a minimum of two eggs from this parasite [22]. They found that multiplex-PCR was able to detect more samples infected with *E. granulosus* and *T. hydatigena* than LAMP.

Finally, LAMP has been portrayed as a simple and easy to implement method that could be used in field surveillance of pathogens, specifically for Echinococcus species. However, so far it has not been used in large epidemiological studies in endemic areas. One of the advantages of LAMP is that there is no need for a thermocycler; however, it still requires some equipment including a heating block or water bath, a centrifuge, gel electrophoresis system (if results need to be confirmed in case of ambiguous colour change), and a clean lab with trained staff. In LAMP reactions, large quantities of DNA are amplified and pose a high risk of contamination when tubes are opened. To minimise this risk, a typical lab where PCR is routinely used, needs different separate environments for DNA extraction (pre-LAMP), setup of LAMP reactions, and visualisation of LAMP results (post-LAMP). The irregular shedding of eggs by the definitive hosts hampers the correct diagnosis of patent *E. multilocularis* infections. Furthermore, there are differences in the egg shedding between hosts; the highest egg output occurs 37 to 42 days post experimental infection in foxes and between 43 and 45 days in dogs in the same conditions [49]. Additionally, the patent period lasts one month for 98% of the worm burden, with a residual worm burden that lasts several months [49]. Therefore, detecting eggs from naturally infected canids in a single sample from naturally infected animals is challenging, as we have shown here since the presence of eggs will depend on the stage of the infection. Furthermore, the taeniid egg isolation has some disadvantages which is important to keep in mind, it can be time consuming especially analysing large number of samples, it requires the purchase of nylon filters which need to be inserted in the lid of the tubes or PET bottles [50] for sieving and it involves the use of a floatation solution which in the case of zinc chloride can be toxic for the manipulator and the environment [51]. A sugar solution [52] could potentially replace the zinc chloride avoiding the toxicity. We also tested the feasibility of using lateral flow dipstick coupled with LAMP reaction, and the results are promising; however, the high cost of the dipsticks precluded us from suggesting them as a cost-effective tool for the diagnosis of E. multilocularis.

Methods based on DNA detection from *E. multilocularis* have been published and extensively reviewed [12,13,53]. In general, they can be divided into those using whole faeces and the ones based on isolated taeniid eggs as starting material for DNA extraction. We focus in this study on the latter method which allows the enrichment of helminth eggs from faecal samples; in doing so, we aim to decrease the presence of PCR inhibitors commonly found in faecal samples based on the thorough wash of the filters used in the sieving process [10,15]. The taeniid egg isolation based on floatation/sieving allows the use of a relatively large amount of faeces as starting material, increasing the chances of finding eggs in canines with low worm burden for example. Up to 20 mL of faeces suspended in ethanol were used in the original publication describing the method [10] reporting a positive PCR in 33 out of 35 samples from foxes infected with *E. multilocularis* (diagnosed at necropsy). On the other hand, protocols for DNA extraction from whole faeces using commercial kits generally accept from 100 up to 500 mg of faeces as starting material and have the advantage of acquiring high yield purified DNA, in theory without PCR inhibitors. But they are expensive, costing between 3.1 and 10.2 USD per sample, and

require between 40 to 90 min to be completed for 10 samples [16], making them unsuitable for use in large-scale prevalence studies. Therefore, there is a need to investigate and standardise cost-effective and straightforward methods for DNA extraction to be used in such studies. The use of NaOH-based DNA extraction methods from different tissues can produce a template of sufficient quality for PCR [36–38,54]. Furthermore, DNA has been extracted from Echinococcus spp. eggs, and also protoscoleces and cyst tissue, using NaOH for direct PCR and LAMP for genotyping studies without a purification step in several publications [26,39–41,55,56]. Similarly, the use of Chelex as a DNA extraction method has also been proposed since the 1990s [33] and used in recent publications as an alternative to producing input for PCR-based diagnostics [57,58] and mammal identification of scats [58] for example. NaOH and Chelex are components of the method used in seminal publications reporting on DNA isolation from taeniid eggs in which a commercial kit was also included [10,11,15,59]. The DNA extraction methods used in the present study offer a cost-effective alternative to the commercial kits for this purpose. If we consider the cost of a NaOH 0.2 M commercially available, the cost per sample is 0.014 USD for protocol 1A; in the case of protocol 1B the cost of InstaGeneTM Matrix per sample is 0.95 USD. Furthermore, the time of DNA extraction is reduced in protocol 1A to 20 and 40 min with protocol 1B for processing ten samples, offering advantages over the 90 min required for DNA extraction described [59]. However, the time needed to isolate eggs needs to be considered as 3 h are required for 10–20 samples (Figure 1).

To conclude, we show two cost-effective methods for DNA extraction from isolated taeniid eggs for direct use in multiplex-PCR and LAMP developed in this study. These procedures circumvent the use of commercial kits for DNA extraction. Compared with the gold standard, the sensitivity of the test remains lower mainly because our methodology is focusing on the detection of eggs, therefore being unable to detect prepatent or low worm burden infections. Nevertheless, considering the reduction in costs and time, we propose the methods for DNA extraction as a valuable tool that can be used in extensive prevalence studies investigating the presence of eggs in canine environmental faecal samples in endemic areas of *E. multilocularis*.

4. Materials and Methods

4.1. Parasites

Eggs of *E. multilocularis* were isolated from the faeces of a naturally infected fox with no *Taenia* or *Mesocestoides* spp., detected at necropsy. Faeces were stored at -80 °C for five days for biosafety reasons and then subjected to floatation with zinc chloride (specific gravity 1.45) and sieving protocol with different nylon filters [10]. Eggs were stored at 4 °C in PBS with penicillin-streptomycin until further use. Multiplex-PCR [15] confirmed that only eggs of *E. multilocularis* and no *Taenia* spp. or *Mesocestoides* spp. were present in the suspension.

DNA was isolated from *E. multilocularis* adult worms and metacestode tissue cultured in vitro [60] using the QIAamp DNA Mini Kit (Qiagen, Hilden, Germany). The same procedure was used to isolate DNA from metacestode or adult stages of different cestodes identified morphologically including *E. granulosus sensu stricto*, *E. equinus*, *E. ortleppi*, *E. intermedius* (G6 and G7), *E. vogeli*, *E. shiquicus*, *Taenia polyacantha*, *T. multiceps*, *T. ovis*, *T. saginata*, *T. solium*, *T. crassiceps*, *T. hydatigena*, *Hydatigera taeniaeformis*, *T. pisiformis*, *T. krabbei*, *Dipylidium caninum*, *Diphyllobothrium latum*, *Mesocestoides litteratus*, and *M. lineatus*. Parasite identification was confirmed with PCR/sequencing of a section of the *cox1* gene [61].

4.2. LAMP

4.2.1. Primer Design

Multiple alignments of the mitochondrial genome of *E. multilocularis* and other *Echinococcus* and *Taenia* species were performed to select a region for designing LAMP primers using Genious R10 V10.1.3. The mitochondrial *nad1* gene of *E. multilocularis* (accession number in GenBank AB668376) was selected as it showed sequence variation between species, and primers were created using Primer Explorer V5 (EikenChemicalCoLtd.) [62]

and Primer Designer 1.16 (Premier) [63]. Primers were tested in silico for specificity using BLAST [64] and Primer-BLAST [65]. In total, seven primer sets, four primers each, were synthesised by Microsynth (Balgach, Switzerland). Primer sets were tested, assessing their specificity and analytical sensitivity to detect DNA from *E. multilocularis*.

4.2.2. Specificity and Analytical Sensitivity

Initially, the primer sets were tested in LAMP reactions using DNA isolated from *E. multilocularis*. Serial dilutions (1:10) of *E. multilocularis* DNA were prepared in triplicate starting at 10 ng/ μ L to 0.1 fg/ μ L. The primer set showing the highest consistency, specifically amplifying only *E. multilocularis* DNA with the lowest concentration as a template, was chosen for standardisation of the LAMP reaction (Table 5). Additionally, we tested all primer sets with DNA from the different cestodes mentioned above. To allow the visualization of samples positive to *E. multilocularis* DNA with the HybriDetect-Universal Lateral Flow Assay Kit (Milenia Biotec, Giessen, Germany), the FIP and BIP primers were labelled at the 5' end with FAM (6-fluorescein amidite) and DIG (digoxigenin) as shown in Table 5.

Table 5. Primer set selected for the detection of *Echinococcus multilocularis* in a LAMP assay. For detection with the HybriDetect-Universal Lateral Flow Assay Kit (LFD), a FAM (6-fluorescein amidite) and DIG (digoxigenin) modifications were included in the Em-FIP and Em-BIP primers.

Primer	Sequence (5'-3')	5'-Modification (LFD)
Em-F3	GCTTGTTGTTGTTTCCATTGA	-
Em-B3	ACAAAACCACCAACC	-
Em-FIP	TCCCTTTCAGACTCCCCATAATCA-TTTTTGGTGTGTGTGTGCTATG	FAM
Em-BIP	AGCGGTATATACTTTACGTGTTTGT-TCATTACAACAATCAACCATGA	DIG
Em-LB	TTGCTTGTGAGTATATAGTTGTATATGTGT	-

4.2.3. LAMP Assay

For the optimisation of LAMP assays, different concentrations of MgSO₄ (4-8 mM); primer mix (2×, 1.5×, 1×, 0.75×, 0.5×, 0.25×); betaine (0.4, 0.8, 1 and 1.2 M); malachite green (0.004%, 0.008% and 0.016%) were tested. Bovine serum albumin (0.1% w/v) was added to the reaction mixture to improve the performance of the Bst Polymerase during amplification. Various amplification times (30–90 min) and temperatures (61–70 °C) were tested. The reagent concentrations showing the strongest amplification (see the visualisation of the LAMP reaction) without unspecific amplification were chosen for the field study. Finally, the LAMP reaction was set up in 25 μ L containing the isothermal amplification buffer from New England Biolabs (Ipswich, MA, USA) [20 mM Tris-HCl, 10 mM (NH₄)₂SO₄, 50 mM KCl, 2 mM MgSO₄, and 0.1% Tween 20], 6 mM MgSO₄, 1.4 mM of each dNTP, 0.8 M betaine, 0.004% malachite green, 8 U/mL Bst 2.0 DNA Polymerase Warmstart (New England Biolabs), 1.6 µM of each FIP(-FAM) and BIP(-DIG) primer, 0.2 µM of each F3 and B3 primer, and 0.4 µM of LB primer (loop primer). Finally, 2 µL DNA were added into each tube as a template. The amplification occurred at 65 °C for 60 min in a heating block. In every experiment, positive control (DNA of E. multilocularis) and negative control (water) were included.

4.2.4. Visualisation of the LAMP Reaction

Three different visualisation methods were used to assess the results of LAMP reactions. First, directly after the incubation time, visual judgment of the colour change was performed and documented. A positive reaction was indicated by different intensities of blue colour and in a negative reaction by a colour change from blue or green to colourless. Secondly, 3 μ L of each LAMP product were subjected to gel electrophoresis in 1.5% agarose gels stained with Gel Red (Biotium, Fremont, CA, USA) and visualised with a UV transilluminator (BioRad, Hercules, CA, USA). A DNA ladder-like pattern represented a positive reaction while no ladder was a negative result. Finally, we used a lateral flow dipstick (HybriDetect 2T, Milenia Biotec, Giessen, Germany) to assess the LAMP reactions with the FIP and BIP labelled primers (Table 5). Briefly, 100 μ L of citrate-phosphate-buffer (Milenia Biotec) were mixed with 10 μ L of the LAMP reaction in a 1.5 mL tube. Then, the dipstick was inserted into each tube and the results were assessed after five to ten minutes. LAMP amplicons labelled with FAM and DIG bind to the anti-digoxigenin antibodies in the strip. The test band will appear on the dipstick, along with the control band, whereas in a negative LAMP reaction, only the control band is visible.

4.2.5. Stability of the LAMP Master Mix

The stability of the LAMP master mix was assessed by preparing aliquots of 23 μ L of a LAMP master mix (as described above) without DNA template in 0.2 mL tubes stored at 4 °C, over a period of six weeks [66]. Every week, two tubes were taken from the fridge and 2 μ L of *E. multilocularis* DNA (1 ng/ μ L) were added to one tube and 2 μ L of water to the second tube as negative control. The amplification and visualisation of the reaction occurred as explained above.

4.2.6. Limit of Detection of LAMP

The limit of detection (No. of eggs) of the LAMP reaction with the primer set from Table 5 was established using eggs of *E. multilocularis* isolated from faeces of a naturally infected fox (as explained above). Eggs were aspirated individually with a micropipette under the microscope into 1.5 mL tubes. Four replicates of tubes containing ten, five, four, three, two, and one egg(s) were prepared in 30 μ L of water. Egg suspensions were subjected to protocols 1A and 1B as explained below for DNA extraction.

4.3. DNA Extraction (Protocols 1A and 1B)

The methods presented basically aim to break the eggs of *E. multilocularis*, and subsequently lyse cells of the embryo to release DNA to be used directly in LAMP and multiplex-PCR without a further purification step. We use here the term "DNA extraction" to refer to these methods for treating taeniid eggs. For protocol 1A, egg suspensions were treated with 0.2 M NaOH (ratio 1:1) and incubated at 95 °C for 10 min in a heating block [26,41]. Tubes were centrifuged (quick spin) and then dilutions with Tris-HCl pH 8.3 (100 mM) were prepared (ratios 1:2, 1:10, 1:50, and 1:100). Finally, 2 µL of each dilution were used as a template for LAMP (developed in this study using primers from Table 5) and multiplex-PCR [15]; the dilution offering the most consistent result was chosen for the DNA extraction in the field study. For protocol 1B, egg suspensions were treated as in protocol 1A without dilution with Tris-HCl. Subsequently, 100 µL of Instagene Matrix (Bio-Rad) were added to each tube and incubated for another 15 min at 56 °C. After vortexing and centrifugation (12,000 × g for 3 min), 2 µL of the supernatant were used directly as a template for LAMP (developed in this study) and multiplex-PCR [15].

4.4. Application of LAMP and Multiplex-PCR in Field Samples

4.4.1. Fox Faecal Samples

In total, 63 faecal samples from foxes shot during the official hunting season (January–February) in the surroundings of Zurich (Switzerland) from 2018 until 2020 were used in this study. At necropsy, fox intestines were examined for the presence of *E. multilocularis* and other parasites based on the SCT method [4]. Faecal samples were collected from the rectum, deposited in 50 mL tubes, and stored at -80 °C for five days to inactivate taeniid eggs. Worm burden was assessed in 45 foxes necropsied in 2020 including 28 positive foxes and 17 negative animals to the presence of *E. multilocularis*. Two aliquots of two grams of faeces per sample were taken for independent isolation of taeniid eggs as previously described [10]. The material retained in the 21 µm filter was collected and carefully screened with an inverted microscope for the presence of taeniid eggs in 10 mL tubes with a flat side. Then, all tubes were centrifuged at $200 \times g$ for 10 min, the supernatant discarded and the sediment was transferred to a 1.5 mL tube for treatments with protocols 1A or 1B and subsequently use as a template for LAMP (from this study) and multiplex-PCR [15].

4.4.2. Dog Faecal Samples

In total, we used 32 faecal samples from dogs collected as part of a large prevalence study for *Echinococcus* spp. in Kyrgyzstan between 2017 and 2018 [including a genetic characterisation of *E. multilocularis* [67] and identification by egg isolation + multiplex-PCR [15]]. In the study mentioned above, taeniid eggs were concentrated from dog faecal samples as previously described and DNA was isolated combining alkaline lysis and the QIA amp Kit [10,50]. Briefly, eggs were resuspended in 200 μ L of distilled water and $25 \ \mu L$ KOH (1M) and 7 μL of DTT (1M) were added and incubated at 65 °C for 15 min. Afterwards, 60 μ L 2M Tris-HCl pH 8.4 and 2 μ L HCl (12.4N/ \geq 37%) were pipetted into the tubes. Finally, 200 µL of Buffer AL (QIAamp Kit) and 20 µL of Proteinase K were added and incubated at 56 °C for 10 min. Then, 50 µL of Chelex solution (50%) were added and tubes were mixed in a rotator for 1 h at room temperature. After centrifugation, the supernatant (approximately 400 μ L) was transferred to a new tube with 200 μ L of ethanol (100%) and mixed; the content of each tube was then transferred into a Qiagen spin column. The protocol continues following the manufacturer instructions for the QIAamp DNA with two washes with buffer AW1 instead of one. DNA was eluted in 100 μ L of 10 mM Tris-HCl, pH 8.3, and stored at -20 °C until use as a template in multiplex-PCR [15].

4.5. Application of LAMP and Multiplex-PCR in Whole Faeces (Protocols 2A and 2B)

Two aliquots of 400–500 mg were taken from faecal samples of 30 foxes infected with *E. multilocularis* and five from negative foxes from the same group of animals used for protocols 1A and 1B (diagnosed at necropsy/SCT). Dog faecal samples included 18 *E. multilocularis* positive and ten negative (diagnosed using egg isolation + multiplex-PCR). One aliquot was treated with the protocol 2A: 400 μ L of NaOH solution (0.2 M) were added, mixed, and incubated at 95 °C for 10 min, then 1 mL of Instagene Matrix (Bio-Rad) was added and the tubes were kept at 56 °C for 15 min. The samples were vortexed, centrifuged at 12,000 × *g* for 3 min, and 500 μ L of the supernatant transferred to a new tube. From the supernatant, a 1:50 dilution with water was prepared and 2 μ L were used as a template for LAMP and multiplex-PCR. The second aliquot was used for DNA isolation following the QIAamp Fast DNA Stool Mini Kit[®] instructions for large stool volumes (protocol 2B).

4.6. Sensitivity and Specificity

The diagnostic sensitivities and specificities were calculated for protocol 1A and 1B by comparing the results of LAMP and multiplex-PCR with the SCT result in foxes and with the initial multiplex-PCR (from the prevalence study) result in the case of dogs.

Supplementary Materials: The following are available online at https://www.mdpi.com/article/10.339 0/pathogens10070847/s1, Table S1. List with all the results from foxes positive and negative to intestinal stages of *E. multilocularis* diagnosed at necropsy and used in the present study. Isolated eggs were treated with protocol 1A and 1B. Results for multiplex PCR also shows if an amplicon corresponding to other cestodes (oc) including *Taenia* spp. was also amplified. Em: *Echinococcus multilocularis*.

Author Contributions: Conceptualization, C.A.A.R., P.D.; methodology, B.J.B., T.K., C.A.A.R.; software, B.J.B., T.K.; validation, B.J.B., G.M., C.A.A.R.; formal analysis, B.J.B., C.A.A.R.; investigation, B.J.B., G.M., C.A.A.R.; resources, K.K.A., M.I., P.A.K., P.D., C.A.A.R.; writing—original draft preparation, B.J.B., P.D., C.A.A.R.; writing—review and editing B.J.B., G.M., T.K., P.A.K., P.D., C.A.A.R.; visualization B.J.B., C.A.A.R.; writing—review and editing B.J.B., G.M., T.K., P.A.K., P.D., C.A.A.R.; visualization B.J.B., C.A.A.R.; writing—review and editing B.J.B., G.M., T.K., P.A.K., P.D., C.A.A.R.; visualization B.J.B., C.A.A.R.; visualization of B.J.B. All authors have read and agreed to the published version of the manuscript.

Funding: Faecal samples from dogs were collected as part of a large prevalence study for *Echinococcus* spp. in Kyrgyzstan between 2017 and 2018 supported by the PT Swiss National Science Foundation (SNSF, grant agreement number 173131–Transmission modelling of emergent echinococcosis in Kyrgyzstan).

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Data Availability Statement: Not applicable.

Acknowledgments: Figure 1 was created using BioRender.com.

Conflicts of Interest: The authors declare no conflict of interest.

References

- 1. Eckert, J.; Deplazes, P. Biological, epidemiological, and clinical aspects of echinococcosis, a zoonosis of increasing concern. *Clin. Microbiol. Rev.* **2004**, *17*, 107–135. [CrossRef]
- 2. Romig, T.; Deplazes, P.; Jenkins, D.; Giraudoux, P.; Massolo, A.; Craig, P.S.; Wassermann, M.; Takahashi, K.; De La Rue, M. Ecology and life cycle patterns of *Echinococcus* species. *Adv. Parasitol.* **2017**, *95*, 213–314. [CrossRef]
- 3. Alvarez Rojas, C.A.; Mathis, A.; Deplazes, P. Assessing the contamination of food and the environment with *Taenia* and *Echinococcus* eggs and their zoonotic transmission. *Curr. Clin. Microbiol. Rep.* **2018**, *5*, 154–163. [CrossRef]
- 4. Eckert, J.; Deplazes, P.; Craig, P.S.; Gemmell, M.A.; Gottstein, B.; Heath, D.; Jenkins, D.J.; Kamiya, M.; Lightowlers, M. Chapter 3 Echinococcosis in animals: Clinical aspects, diagnosis and treatment. In WHO/OIE Manual on Echinococcosis in Humans and Animals: A Public Health Problem of Global Concern; Eckert, J., Gemmell, M.A., Meslin, F.X., Pawlowski, Z.S., Eds.; World Organisation for Animal Health: Paris, France, 2001.
- 5. Deplazes, P.; Gottstein, B.; Eckert, J.; Jenkins, D.J.; Ewald, D.; Jimenez-Palacios, S. Detection of *Echinococcus* coproantigens by enzyme-linked immunosorbent assay in dogs, dingoes and foxes. *Parasitol. Res.* **1992**, *78*, 303–308. [CrossRef]
- 6. Deplazes, P.; Alther, P.; Tanner, I.; Thompson, R.C.A.; Eckert, J. *Echinococcus multilocularis* coproantigen detection by enzyme-linked immunosorbent assay in fox, dog, and cat populations. *J. Parasitol.* **1999**, *85*, 115–121. [CrossRef]
- Sakai, H.; Nonaka, N.; Yagi, K.; Oku, Y.; Kamiya, M. Coproantigen detection in a survey of *Echinococcus multilocularis* infection among red foxes, *Vulpes vulpes schrencki*, in Hokkaido, Japan. J. Vet. Med. Sci. 1998, 60, 639–641. [CrossRef] [PubMed]
- 8. Machnicka, B.; Dziemian, E.; Rocki, B.; Kołodziej-Sobocińska, M. Detection of *Echinococcus multilocularis* antigens in faeces by ELISA. *Parasitol. Res.* 2003, 91, 491–496. [CrossRef] [PubMed]
- 9. Dinkel, A.; Von Nickisch-Rosenegk, M.; Bilger, B.; Merli, M.; Lucius, R.; Romig, T. Detection of *Echinococcus multilocularis* in the definitive host: Coprodiagnosis by PCR as an alternative to necropsy. *J. Clin. Microbiol.* **1998**, *36*, 1871–1876. [CrossRef]
- 10. Mathis, A.; Deplazes, P.; Eckert, J. An improved test system for PCR-based specific detection of *Echinococcus multilocularis* eggs. *J. Helminthol.* **1996**, *70*, 219–222. [CrossRef] [PubMed]
- 11. Bretagne, S.; Guillou, J.P.; Morand, M.; Houin, R. Detection of *Echinococcus multilocularis* DNA in fox faeces using DNA amplification. *Parasitology* **1993**, *106*, 193–199. [CrossRef] [PubMed]
- 12. Deplazes, P.; Dinkel, A.; Mathis, A. Molecular tools for studies on the transmission biology of *Echinococcus multilocularis*. *Parasitology* **2003**, 127, S53–S61. [CrossRef]
- 13. Siles-Lucas, M.; Casulli, A.; Conraths, F.J.; Müller, N. Chapter Three—Laboratory diagnosis of *Echinococcus* spp. in human patients and infected animals. *Adv. Parasitol.* **2017**, *96*, 159–257. [PubMed]
- 14. Yimam, A.E.; Nonaka, N.; Oku, Y.; Kamiya, M. Prevalence and intensity of *Echinococcus multilocularis* in red foxes (*Vulpes vulpes schrencki*) and raccoon dogs (*Nyctereutes procyonoides albus*) in Otaru City, Hokkaido, Japan. *Jpn. J. Vet. Res.* **2002**, 49, 287–296.
- 15. Trachsel, D.; Deplazes, P.; Mathis, A. Identification of taeniid eggs in the faeces from carnivores based on multiplex PCR using targets in mitochondrial DNA. *Parasitology* **2007**, *134*, 911–920. [CrossRef] [PubMed]
- Maksimov, P.; Schares, G.; Press, S.; Fröhlich, A.; Basso, W.; Herzig, M.; Conraths, F.J. Comparison of different commercial DNA extraction kits and PCR protocols for the detection of *Echinococcus multilocularis* eggs in faecal samples from foxes. *Vet. Parasitol.* 2017, 237, 83–93. [CrossRef] [PubMed]
- Videnska, P.; Smerkova, K.; Zwinsova, B.; Popovici, V.; Micenkova, L.; Sedlar, K.; Budinska, E. Stool sampling and DNA isolation kits affect DNA quality and bacterial composition following 16S rRNA gene sequencing using MiSeq Illumina platform. *Sci. Rep.* 2019, *9*, 3837. [CrossRef] [PubMed]
- 18. Kaneko, H.; Kawana, T.; Fukushima, E.; Suzutani, T. Tolerance of loop-mediated isothermal amplification to a culture medium and biological substances. *J. Biochem. Biophys. Methods* **2007**, *70*, 499–501. [CrossRef] [PubMed]
- 19. Grab, D.J.; Lonsdale-Eccles, J.; Inoue, N. Lamp for tadpoles. Nat. Methods 2005, 2, 635. [CrossRef] [PubMed]
- 20. Francois, P.; Tangomo, M.; Hibbs, J.; Bonetti, E.-J.; Boehme, C.C.; Notomi, T.; Perkins, M.D.; Schrenzel, J. Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunol. Med. Microbiol.* **2011**, *62*, 41–48. [CrossRef]
- Ni, X.; Mcmanus, D.P.; Yan, H.; Yang, J.; Lou, Z.; Li, H.; Li, L.; Lei, M.; Cai, J.; Fan, Y.; et al. Loop-mediated isothermal amplification (LAMP) assay for the identification of *Echinococcus multilocularis* infections in canine definitive hosts. *Parasites Vectors* 2014, 7, 254. [CrossRef]
- 22. Feng, K.; Li, W.; Guo, Z.; Duo, H.; Fu, Y.; Shen, X.; Tie, C.; Xiao, C.; Luo, Y.; Qi, G.; et al. Development of LAMP assays for the molecular detection of taeniid infection in canine in Tibetan rural area. *J. Vet. Med. Sci.* 2017, *79*, 1986–1993. [CrossRef]
- Avila, H.G.; Mozzoni, C.; Trangoni, M.D.; Cravero, S.L.P.; Pérez, V.M.; Valenzuela, F.; Gertiser, M.L.; Butti, M.J.; Kamenetzky, L.; Jensen, O.; et al. Development of a copro-LAMP assay for detection of several species of *Echinococcus granulosus sensu lato* complex. *Vet. Parasitol.* 2020, 277, 109017. [CrossRef] [PubMed]

- Ni, X.-W.; Mcmanus, D.P.; Lou, Z.-Z.; Yang, J.-F.; Yan, H.-B.; Li, L.; Li, H.-M.; Liu, Q.-Y.; Li, C.-H.; Shi, W.-G.; et al. A comparison of loop-mediated isothermal amplification (LAMP) with other surveillance tools for *Echinococcus granulosus* diagnosis in canine definitive hosts. *PLoS ONE* 2014, 9, e100877. [CrossRef] [PubMed]
- Salant, H.; Abbasi, I.; Hamburger, J. The development of a loop-mediated isothermal amplification method (LAMP) for *Echinococ-cus granulosus* coprodetection. *Am. J. Trop. Med. Hyg.* 2012, 87, 883–887. [CrossRef] [PubMed]
- Wassermann, M.; Mackenstedt, U.; Romig, T. A loop-mediated isothermal amplification (LAMP) method for the identification of species within the *Echinococcus granulosus* complex. *Vet. Parasitol.* 2014, 200, 97–103. [CrossRef] [PubMed]
- Rahman, S.M.M.; Song, H.B.; Jin, Y.; Oh, J.-K.; Lim, M.K.; Hong, S.-T.; Choi, M.-H. Application of a loop-mediated isothermal amplification (LAMP) assay targeting cox1 gene for the detection of *Clonorchis sinensis* in human fecal samples. *PLoS Negl. Trop. Dis.* 2017, *11*, e00059955. [CrossRef]
- 28. Fernández-Soto, P.; Fernández-Medina, C.; Cruz-Fernández, S.; Crego-Vicente, B.; Febrer-Sendra, B.; García-Bernalt Diego, J.; Gorgojo-Galindo, Ó.; López-Abán, J.; Vicente Santiago, B.; Muro Álvarez, A. Whip-LAMP: A novel LAMP assay for the detection of *Trichuris muris*-derived DNA in stool and urine samples in a murine experimental infection model. *Parasites Vectors* 2020, 13, 552. [CrossRef]
- Fernández-Soto, P.; Sánchez-Hernández, A.; Gandasegui, J.; Bajo Santos, C.; López-Abán, J.; Saugar, J.M.; Rodríguez, E.; Vicente, B.; Muro, A. Strong-LAMP: A LAMP Assay for *Strongyloides* spp. Detection in Stool and Urine Samples. Towards the Diagnosis of Human Strongyloidiasis Starting from a Rodent Model. *PLoS Negl. Trop. Dis.* 2016, *10*, e0004836. [CrossRef]
- Mugambi, R.M.; Agola, E.L.; Mwangi, I.N.; Kinyua, J.; Shiraho, E.A.; Mkoji, G.M. Development and evaluation of a Loop Mediated Isothermal Amplification (LAMP) technique for the detection of hookworm (*Necator americanus*) infection in fecal samples. *Parasites Vectors* 2015, *8*, 574. [CrossRef] [PubMed]
- 31. Nkouawa, A.; Sako, Y.; Nakao, M.; Nakaya, K.; Ito, A. Loop-mediated isothermal amplification method for differentiation and rapid detection of *Taenia* species. J. Clin. Microbiol. 2009, 47, 168–174. [CrossRef]
- 32. Steiner, J.J.; Polemba, C.J.; Fjellstrom, R.G.; Elliott, L.F. A rapid one-tube genomic DNA extraction process for PCR and RAPD analyses. *Nucleic Acids Res.* **1995**, *23*, 2569–2570. [CrossRef]
- 33. De Lamballerie, X.; Zandotti, C.; Vignoli, C.; Bollet, C.; De Micco, P. A one-step microbial DNA extraction method using "Chelex 100" suitable for gene amplification. *Res. Microbiol.* **1992**, *143*, 785–790. [CrossRef]
- Kawasaki, E.S. Sample preparation from blood, cells and other fluids. In PCR Protocols; Innis, M.A., Gelfand, D.H., Sninsky, J.J., White, T.J., Eds.; Academic Press: San Diego, CA, USA, 1990; pp. 146–152.
- 35. Wang, H.; Qi, M.; Cutler, A.J. A simple method of preparing plant samples for PCR. *Nucleic Acids Res.* **1993**, *21*, 4153–4154. [CrossRef] [PubMed]
- Osmundson, T.W.; Eyre, C.A.; Hayden, K.M.; Dhillon, J.; Garbelotto, M.M. Back to basics: An evaluation of NaOH and alternative rapid DNA extraction protocols for DNA barcoding, genotyping, and disease diagnostics from fungal and oomycete samples. *Mol. Ecol. Resour.* 2013, *13*, 66–74. [CrossRef]
- Satya, P.; Mitra, S.; Ray, D.P.; Mahapatra, B.S.; Karan, M.; Jana, S.; Sharma, A.K. Rapid and inexpensive NaOH based direct PCR for amplification of nuclear and organelle DNA from ramie (Boehmeria nivea), a bast fibre crop containing complex polysaccharides. *Ind. Crop. Prod.* 2013, 50, 532–536. [CrossRef]
- 38. Werner, O.; Ros, R.M.; Guerra, J. Direct amplification and NaOH extraction: Two rapid and simple methods for preparing bryophyte DNA for polymerase chain reaction (PCR). *J. Bryol.* **2002**, *24*, 127–131. [CrossRef]
- Hüttner, M.; Nakao, M.; Wassermann, T.; Siefert, L.; Boomker, J.D.F.; Dinkel, A.; Sako, Y.; Mackenstedt, U.; Romig, T.; Ito, A. Genetic characterization and phylogenetic position of *Echinococcus felidis* Ortlepp, 1937 (Cestoda: Taeniidae) from the African lion. *Int. J. Parasitol.* 2008, *38*, 861–868. [CrossRef]
- 40. Mulinge, E.; Magambo, J.; Odongo, D.; Njenga, S.; Zeyhle, E.; Mbae, C.; Kagendo, D.; Addy, F.; Ebi, D.; Wassermann, M.; et al. Molecular characterization of *Echinococcus* species in dogs from four regions of Kenya. *Vet. Parasitol.* **2018**, 255, 49–57. [CrossRef] [PubMed]
- 41. Nakao, M.; Sako, Y.; Ito, A. Isolation of polymorphic microsatellite loci from the tapeworm *Echinococcus multilocularis*. *Infect. Genet. Evol.* **2003**, *3*, 159–163. [CrossRef]
- 42. Wassermann, M.; Aschenborn, O.; Aschenborn, J.; Mackenstedt, U.; Romig, T. A sylvatic lifecycle of *Echinococcus equinus* in the Etosha National Park, Namibia. *Int. J. Parasitol. Parasites Wildl.* **2015**, *4*, 97–103. [CrossRef] [PubMed]
- 43. Notomi, T.; Okayama, H.; Masubuchi, H.; Yonekawa, T.; Watanabe, K.; Amino, N.; Hase, T. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* 2000, *28*, E63. [CrossRef] [PubMed]
- 44. Xu, J.; Rong, R.; Zhang, H.Q.; Shi, C.J.; Zhu, X.Q.; Xia, C.M. Sensitive and rapid detection of *Schistosoma japonicum* DNA by loop-mediated isothermal amplification (LAMP). *Int. J. Parasitol.* **2010**, *40*, 327–331. [CrossRef] [PubMed]
- Cai, X.Q.; Xu, M.J.; Wang, Y.H.; Qiu, D.Y.; Liu, G.X.; Lin, A.; Tang, J.D.; Zhang, R.L.; Zhu, X.Q. Sensitive and rapid detection of *Clonorchis sinensis* infection in fish by loop-mediated isothermal amplification (LAMP). *Parasitol. Res.* 2010, 106, 1379–1383. [CrossRef] [PubMed]
- Skrzypek, K.; Karamon, J.; Samorek-Pieróg, M.; Dąbrowska, J.; Kochanowski, M.; Sroka, J.; Bilska-Zając, E.; Cencek, T. Comparison of two DNA extraction methods and two PCRs for detection of *Echinococcus multilocularis* in the stool samples of naturally infected red foxes. *Animals* 2020, *10*, 2381. [CrossRef] [PubMed]
- 47. Lass, A.; Ma, L.; Kontogeorgos, I.; Xueyong, Z.; Li, X.; Karanis, P. Contamination of wastewater with *Echinococcus multilocularis* Possible implications for drinking water resources in the QTP China. *Water Res.* **2020**, *170*, 115334. [CrossRef]

- Chamai, M.; Omadang, L.; Erume, J.; Ocaido, M.; Oba, P.; Othieno, E.; Bonaventure, S.; Kitibwa, A. Identification of *Echinococcus granulosus* strains using polymerase chain reaction–restriction fragment length polymorphism amongst livestock in Moroto district, Uganda. *Onderstepoort J. Vet. Res.* 2016, *83*, 1–7. [CrossRef] [PubMed]
- 49. Kapel, C.M.O.; Torgerson, P.R.; Thompson, R.C.A.; Deplazes, P. Reproductive potential of *Echinococcus multilocularis* in experimentally infected foxes, dogs, raccoon dogs and cats. *Int. J. Parasitol.* **2006**, *36*, 79–86. [CrossRef]
- 50. Guggisberg, A.R.; Alvarez Rojas, C.A.; Kronenberg, P.A.; Miranda, N.; Deplazes, P. A sensitive, one-way sequential sieving method to isolate helminths' eggs and protozoal oocysts from lettuce for genetic identification. *Pathogens* 2020, *9*, 624. [CrossRef]
- 51. Bodar, C.W.M.; Pronk, M.E.J.; Sijm, D.T.H.M. The European Union risk assessment on zinc and zinc compounds: The process and the facts. *Integr. Environ. Assess. Manag.* 2005, *1*, 301–319. [CrossRef]
- 52. Matsuo, K.; Kamiya, H. Modified sugar centrifugal flotation technique for recovering *Echinococcus multilocularis* eggs from soil. *J. Parasitol.* **2005**, *91*, 208–209. [CrossRef]
- 53. Mathis, A.; Deplazes, P. Copro-DNA tests for diagnosis of animal taeniid cestodes. Parasitol. Int. 2006, 55, S87–S90. [CrossRef]
- 54. Meeker, N.D.; Hutchinson, S.A.; Ho, L.; Trede, N.S. Method for isolation of PCR-ready genomic DNA from zebrafish tissues. *Biotechniques* 2007, 43, 610–614. [CrossRef]
- 55. Boubaker, G.; Macchiaroli, N.; Prada, L.; Cucher, M.A.; Rosenzvit, M.C.; Ziadinov, I.; Deplazes, P.; Saarma, U.; Babba, H.; Gottstein, B.; et al. A Multiplex PCR for the simultaneous detection and genotyping of the *Echinococcus granulosus* complex. *PLoS Negl. Trop. Dis.* 2013, 7, e2017. [CrossRef] [PubMed]
- 56. Sean Walsh, P.; Metzger, D.A.; Higuchi, R. Chelex 100 as a medium for simple extraction of DNA for PCR-based typing from forensic material. *Biotechniques* **2013**, *54*, 134–139. [CrossRef]
- 57. Singh, U.A.; Kumari, M.; Iyengar, S. Method for improving the quality of genomic DNA obtained from minute quantities of tissue and blood samples using Chelex 100 resin. *Biol. Proced. Online* **2018**, *20*, 12. [CrossRef] [PubMed]
- 58. Paxinos, E.; Mcintosh, C.; Ralls, K.; Fleischer, R. A noninvasive method for distinguishing among canid species: Amplification and enzyme restriction of DNA from dung. *Mol. Ecol.* **1997**, *6*, 483–486. [CrossRef] [PubMed]
- Stefanic, S.; Shaikenov, B.S.; Deplazes, P.; Dinkel, A.; Torgerson, P.R.; Mathis, A. Polymerase chain reaction for detection of patent infections of *Echinococcus granulosus* ("sheep strain") in naturally infected dogs. *Parasitol. Res.* 2004, 92, 347–351. [CrossRef] [PubMed]
- 60. Laurimäe, T.; Kronenberg, P.A.; Alvarez Rojas, C.A.; Ramp, T.W.; Eckert, J.; Deplazes, P. Long-term (35 years) cryopreservation of *Echinococcus multilocularis* metacestodes. *Parasitology* **2020**, *147*, 1048–1054. [CrossRef]
- 61. Bowles, J.; Blair, D.; Mcmanus, D.P. Genetic variants within the genus *Echinococcus* identified by mitochondrial DNA sequencing. *Mol. Biochem. Parasitol.* **1992**, *54*, 165–173. [CrossRef]
- 62. ECC Ltd. PrimerExplorer V5. Available online: http://primerexplorer.jp/lampv5e/index.html (accessed on 24 August 2020).
- 63. PREMIER Biosoft. LAMP Designer. Available online: http://www.premierbiosoft.com/isothermal/index.html (accessed on 24 August 2020).
- 64. Altschul, S.F.; Gish, W.; Miller, W.; Myers, E.W.; Lipman, D.J. Basic local alignment search tool. *J. Mol. Biol.* **1990**, *215*, 403–410. [CrossRef]
- 65. Ye, J.; Coulouris, G.; Zaretskaya, I.; Cutcutache, I.; Rozen, S.; Madden, T.L. Primer-BLAST: A tool to design target-specific primers for polymerase chain reaction. *BMC Bioinform.* **2012**, *13*, 134. [CrossRef]
- 66. Kamber, T.; Koekemoer, L.L.; Mathis, A. Loop-mediated isothermal amplification (LAMP) assays for *Anopheles funestus* group and *Anopheles gambiae* complex species. *Med. Vet. Entomol.* 2020, 34, 295–301. [CrossRef] [PubMed]
- 67. Alvarez Rojas, C.A.; Kronenberg, P.A.; Aitbaev, S.; Omorov, R.A.; Abdykerimov, K.K.; Paternoster, G.; Mullhaupt, B.; Torgerson, P.; Deplazes, P. Genetic diversity of *Echinococcus multilocularis* and *Echinococcus granulosus sensu lato* in Kyrgyzstan: The A2 haplotype of *E. multilocularis* is the predominant variant infecting humans. *PLoS Negl. Trop. Dis.* **2020**, *14*, e0008242. [CrossRef] [PubMed]