#### **Ohio Wesleyan University**

## Digital Commons @ OWU

Student Symposium

2023

Apr 20th, 5:21 PM - 5:36 PM

## Ohio Wesleyan Bacillus Collection Student Research

Madeline Russell Ohio Wesleyan University

Nandini Arora Ohio Wesleyan University

**Emily Bode** Ohio Wesleyan University

Morgan Chmielenski Ohio Wesleyan University

Rebecca Coppenbarger Ohio Wesleyan University

See next page for additional authors

Follow this and additional works at: https://digitalcommons.owu.edu/studentsymposium



Part of the Biology Commons

Russell, Madeline; Arora, Nandini; Bode, Emily; Chmielenski, Morgan; Coppenbarger, Rebecca; Faeth, Ginny; Hutchison, Carianne; Klaus, Mindi; Mendenhall, Jonna; Molotievskiy, Natalia; Nichols, Jess; Paul, Jennifer; Schlesinger, Brynn; Smith, Olivia; Steed, Myles; Thierauf, Abby; and Vore, Ivan, "Ohio Wesleyan Bacillus Collection Student Research" (2023). Student Symposium. 3.

https://digitalcommons.owu.edu/studentsymposium/2023/presentations/3

This Presentation is brought to you for free and open access by the Student Scholarship at Digital Commons @ OWU. It has been accepted for inclusion in Student Symposium by an authorized administrator of Digital Commons @ OWU. For more information, please contact sbchaney@owu.edu.

resen	er Information							
adeline Russell, Nandini Arora, Emily Bode, Morgan Chmielenski, Rebecca Coppenbarger, Ginny Faet arianne Hutchison, Mindi Klaus, Jonna Mendenhall, Natalia Molotievskiy, Jess Nichols, Jennifer Paul, ynn Schlesinger, Olivia Smith, Myles Steed, Abby Thierauf, and Ivan Vore								

# BIOL 328: Bacterial Physiology

Dr. Laura Tuhela-Reuning

Maddy Russell, Teaching Laboratory Assistant

## **An Introduction**

- What is Bacterial Physiology?
  - The study of the functionality and life-process that keep bacteria alive
- What was the goal for students this semester?
  - Individually planned and executed research projects involving the OWU Bacillus Collection
- What is the Bacillus Collection?
  - An independent research project run by current and former faculty analyzing feather degrading bacteria isolated from songbird feathers across three states.

#### About me:

- Senior Microbiology Student
- Laboratory Teaching Assistant and Research Mentor





# Chief of Staphs

Mindi Klaus, Natalia Molotievskiy, Brynn Schlesinger, Abby Thierauf



What is the relationship between feather degradation and biofilm formation in bacteria?

### **Experiments:**

- 1. Screening of efficient feather-degrading strains
- 2. Screening of bacterial isolates for their ability to form biofilms
- 3. Quantitative analyze of feather degradation ability
- 4. Scanning Electron Microscopy analysis

## How did they address your question?

-We screened originally to find the best feather degraders and biofilm formers, and then subsequent experiments tested how feather degradation and biofilm formation influence one another

# Chief of Staphs











#### Conclusions:

- 5/10 strains showed substantial feather degradation
- 2/10 strains showed positive Biofilm formation ability
- 5154B was the overall best feather degrader!
- Biofilm data is still being collected, As well as the results from the dual inoculated flask

_									
-	1	2	3	4	5	6	7	8	9
Α	0.735	0.770	0.737	0.718	0.836	0.970	1.026	0.563	1.353
В	0.727	0.738	0.777	1.082	1.390	1.561	1.701	1.689	1.836
С	0.707	0.643	0.715	0.727	0.738	0.755	0.767	0.755	0.764
D	0.716	0.653	0.690	0.671	0.684	0.727	0.745	0.695	0.813
Е	0.732	1.060	0.921	0.998	1.152	1.149	1.176	0.968	1.304
F	0.696	0.718	0.687	0.708	0.732	0.743	0.760	0.782	0.837
G	0.363	0.359	0.369	0.362	0.360	0.371	0.358	0.363	0.388


# **Biochemical Warfare Agents**

Myles Steed, Jennifer Paul, Carianne Hutchison, & Rebecca Coppenbarger

Research Question: What feather bacterial isolates in the OWU Bacillus collection produce β-keratinase and what is their identity?



https://www.forbes.com/sites/neilyeoh/2018/05/17/startup-waste-feathers-sustainable-packaging/?sh=2f6038ce2dfb



https://en.wikipedia.org/wiki/Bacillus\_licheniformis

# **Biochemical Warfare Agents**

Myles Steed, Jennifer Paul, Carianne Hutchison, & Rebecca Coppenbarger

## Experimental Design:

- Isolated 7 bacteria from OWU Bacillus collection
- Isolated samples were screened for β-keratinase through:
  - Incolulating isolates in feather medium
  - Taking a sample of media every 12 hours
  - Centrifuging media samples to separate bacteria residue from feathers and further analyzing the residue
- Samples were run through a spectrophotometer to determine if there was any β-keratinase activity
  - Specifically were determining if there were any oligopeptides produced

# **Biochemical Warfare Agents**


Myles Steed, Jennifer Paul, Carianne Hutchison, & Rebecca Coppenbarger

#### Results:

- For our first set of data, 5 of the 7 isolates had significant oligopeptide production
- We are in the process of analyzing our second set of data of the

same isolates

Fig 1. Oligopeptide light absorbance in 5145V isolated strain




#### Conclusions

 We cannot make any final conclusions yet; however, it is likely that 5 of our 7 isolates produce β-keratinase

# Team AUGUR Research Project

Ivan Vore, Nandini Arora, Jonna Mendenhall, Ginny Faeth

- Research Question: do Bacillus spp. produce antibiotics?
  - Hypothesis: Bacillus spp. Produce antibiotics in the presence of competing microbes.
- Method: Cross Streaking a qualitative testing for antibiotic activity, measuring zones of inhibition (ZOI), indicating antibiotic production.




# **Results**

## **Bacillus spp. vs. Gram Negatives**









## Bacillus spp. vs. Bacillus spp.



## **Team Ur Mom**

Emily Bode, Olivia Smith, Jess Nichols, and Morgan Chmielenski

## Our Research Question:

"Are plasmids, like 5KfuGFR and tsPurple5K able to be inserted into gram positive bacteria, specifically the unknown *Bacillus* species isolated from birds."

## Our Experiments:

- Streaked 6 species of Bacillus onto slants
- Growth curve with spectrometer
- Plasmids were screened in unknown species
- Plasmid of interest were inserted into bacillus and transformed using a modified protoplast procedure

## Team Ur Mom


Emily Bode, Olivia Smith, Jess Nichols, and Morgan Chmielenski


#### Results

- Plasmid Extraction Bacillus
  - No liquid
- Plasmid Extraction E. coli
  - Successful
- Protoplasts
  - Successful
- Plating on Antibiotic plates
  - Purple = growth
  - pGlow = no growth

### Conclusion

- Conjugation = successful











Do you have any questions?

Special thanks to Melinda Endres, Kiley Lewin, and Kayce Tomcho for lab prep, experimental consultation, and emotional support!

**CREDITS**: This presentation template was created by **Slidesgo**, and includes icons by **Flaticon**, and infographics & images by **Freepik**