



University of Groningen

Ion-extraction from the CISe gas catcher

Cipagauta Mora, Jennifer; Brajkovic, Marko; Bouwman, Niels; Chen, X.; Even, Julia; Hartigan, Briain; Soylu, Arif; van der Werff, Lisa

IMPORTANT NOTE: You are advised to consult the publisher's version (publisher's PDF) if you wish to cite from it. Please check the document version below.

Publication date: 2023

Link to publication in University of Groningen/UMCG research database

Citation for published version (APA): Cipagauta Mora, J., Brajkovic, M., Bouwman, N., Chen, X., Even, J., Hartigan, B., Soylu, A., & van der Werff, L. (2023). *Ion-extraction from the CISe gas catcher*. Poster session presented at NWO Physics -2023, Veldhoven, Netherlands.

Copyright

Other than for strictly personal use, it is not permitted to download or to forward/distribute the text or part of it without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license (like Creative Commons).

The publication may also be distributed here under the terms of Article 25fa of the Dutch Copyright Act, indicated by the "Taverne" license. More information can be found on the University of Groningen website: https://www.rug.nl/library/open-access/self-archiving-pure/taverneamendment.

Take-down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Downloaded from the University of Groningen/UMCG research database (Pure): http://www.rug.nl/research/portal. For technical reasons the number of authors shown on this cover page is limited to 10 maximum.



2 CISe project

Our goal is to develop a new Chemical Isobaric Separation technique using a gas catcher. Gas catchers are used to slow down high-energy nuclear reaction products by interacting with a buffer gas. The ions lose energy in collisions with the buffer gas atoms through ionization and excitation. As a first step we are optimizing the ion extraction process within the gas catcher by testing different DC and RF settings.





The DC carpet's electric field gradient pushes ions towards the center, while the time-averaged repelling force of the RF voltage makes them hover. Applying a 90° phase difference generates a traveling potential wave, allowing the ions to surf over the carpet towards the exit hole.

Simulated DC and RF potential distribution



To study the extraction of ions we are using:

- A gas cell filled with 50 mbar of **Helium**. Helium has the advantage that it has a low density and high ionization potential, reducing the probability of ion neutralization.
- A ²²³Ra source that decays into ²¹⁹Rn. Radon is a noble gas with high ionization potential, making it one of the most challenging elements to extract.
- High electric field gradients to reduce the probability of recombination and to achieve a shorter extraction time.

Potential [V]

Anđelić, B. (2021). Direct mass measurements of No, Lr and Rf isotopes with SHIPTRAP and developments for chemical isobaric separation.



11	
4.77 m	

Parameter	Applied voltage
DC cage	300(1) V
DC target holder 1	240(1) V
DC Aluminium mesh	-2400(1) V
DC Carpet	40(1) V and 50(1) V
RF Carpet	63.4(1) Vpp to 88.35(1) Vpp

Extraction yields of 10-15% of 219 Rn have been achieved.

6 Summary

- In this study, alpha spectroscopy measurements showed that 219 Rn⁺, which has a half-life of 3.96 s, was extracted from the gas catcher.
- The results indicated a maximum extraction efficiency of 15%. Further studies will be conducted to optimize the extraction efficiency and unlock the full potential of this technique.

We thank the accelerator support staff, the mechanical workshop and the electronic workshop at UMCG-PARTREC.

This work is funded by the European Research Council through an ERC Starting Grant (No. 803740) and by the Dutch Research Council through a NWO Natuurkunde Projectruimte grant (No. 68091103) and by the Scholarship Program No. 906 of the Ministry of Science, Technology, and Innovation of Colombia.