

# A masking technique to perform high accuracy diffraction experiments on FIB-machined specimens

**Citation for published version (APA):**

Vermeij, T., Plancher, E., & Tasan, C. C. (2017). *A masking technique to perform high accuracy diffraction experiments on FIB-machined specimens*. Poster session presented at EBSD 2017, Oxford, United Kingdom. <https://doi.org/10.13140/RG.2.2.18529.20324>

**DOI:**

[10.13140/RG.2.2.18529.20324](https://doi.org/10.13140/RG.2.2.18529.20324)

**Document status and date:**

Published: 01/04/2017

**Document Version:**

Publisher's PDF, also known as Version of Record (includes final page, issue and volume numbers)

**Please check the document version of this publication:**

- A submitted manuscript is the version of the article upon submission and before peer-review. There can be important differences between the submitted version and the official published version of record. People interested in the research are advised to contact the author for the final version of the publication, or visit the DOI to the publisher's website.
- The final author version and the galley proof are versions of the publication after peer review.
- The final published version features the final layout of the paper including the volume, issue and page numbers.

[Link to publication](#)

**General rights**

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal.

If the publication is distributed under the terms of Article 25fa of the Dutch Copyright Act, indicated by the "Taverne" license above, please follow below link for the End User Agreement:

[www.tue.nl/taverne](http://www.tue.nl/taverne)

**Take down policy**

If you believe that this document breaches copyright please contact us at:

[openaccess@tue.nl](mailto:openaccess@tue.nl)

providing details and we will investigate your claim.


See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/322473826>

## A masking technique to perform high accuracy diffraction experiments on FIB-machined specimens


Poster · April 2017  
DOI: 10.13140/RG.2.2.18123.20174


CITATIONS  
0

3 authors, including:


 **Tijmen Vermeij**  
Technische Universiteit Eindhoven  
PUBLICATIONS 4 CITATIONS  
[SEE PROFILE](#)

Some of the authors of this publication are also working on these related projects:

 [Laue-DIC for high strain resolution](#) [View project](#)

 [Assessment and prevention of FIB damage](#) [View project](#)

READS  
27

 **Emeric Plancher**  
Massachusetts Institute of Technology  
10 PUBLICATIONS 86 CITATIONS  
[SEE PROFILE](#)




# A masking technique to perform high accuracy diffraction experiments on FIB-machined specimens

T. Vermeij<sup>1,2</sup>, E. Plancher<sup>1</sup>, C.C. Tasan<sup>1</sup>


<sup>1</sup> Massachusetts Institute of Technology, Department of Materials Science and Engineering, 77 Massachusetts Avenue, Cambridge, MA 02139, USA

<sup>2</sup> Eindhoven University of Technology, Department of Mechanical Engineering, Den Dolech 2, 5612 AZ Eindhoven, The Netherlands

T. Vermeij et al. / Ultramicroscopy 186 (2018) 35–41



## Goal



## Challenges

**FIB damage effects**

**Ion implantation:**


- Creates crystal defects by material-ion interactions
- Decreases crystallinity and EBSD quality
- Extra features on ECCI

**Redeposition:**


- Leaves extra amorphous layer in FIB areas
- Signal/noise ratio decrease in diffraction/channeling
- Diminution of contrast on EBSD and ECCI

• Can a protective mask avoid the damage?

## Method



## Without protection



## With protection

**Umbrella fabrication**

- Molded soft PDMS pillar, cut off with FIB
- Rough umbrella block on Silicon surface
- Inset: flat cutting of umbrella. Deposit layer of Pt
- Final umbrella attached to micro manipulator

**Umbrella assisted FIB & HR-EBSD (in situ)**

- HR-EBSD reference spots on Silicon substrate
- Umbrella placed on Silicon surface at spot locations
- FIB milling around umbrella
- Protected surface after umbrella removal, new HR-EBSD spots taken

**Correlate spots**

**HR-EBSD results**

• Using the umbrella mask, the protected surface retains a quality close to the original one!

## Discussion

- Without protection, the damaging effect of FIB on surface quality is devastating for EBSD and ECCI. This is mainly caused by amorphous redeposition layers
- Good protection not trivial, optimized umbrella is necessary
- Protection not perfect yet, redeposition can still come in under the block

## Conclusions

- Without protection, advanced diffraction techniques cannot be used near or on FIB machined specimens
- A soft rubber umbrella, with a flat Pt coated bottom is highly effective in protection
- The umbrella fabrication process allows freedom in size, materials and applications (e.g. protection of micro tensile bars)