

# Cold plasma for bacterial inactivation and its effects on keratinocytes

***Citation for published version (APA):***

Boekema, B. K. H. L., van Gils, K., Hofmann, S., Bruggeman, P. J., Middelkoop, E., & Kroesen, G. M. W. (2012). *Cold plasma for bacterial inactivation and its effects on keratinocytes*. Poster session presented at 4th International Conference for Plasma Medicine, ICPM 4, Orléans, France.

***Document status and date:***

Published: 01/01/2012

***Document Version:***

Publisher's PDF, also known as Version of Record (includes final page, issue and volume numbers)

***Please check the document version of this publication:***

- A submitted manuscript is the version of the article upon submission and before peer-review. There can be important differences between the submitted version and the official published version of record. People interested in the research are advised to contact the author for the final version of the publication, or visit the DOI to the publisher's website.
- The final author version and the galley proof are versions of the publication after peer review.
- The final published version features the final layout of the paper including the volume, issue and page numbers.

[Link to publication](#)

***General rights***

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal.

If the publication is distributed under the terms of Article 25fa of the Dutch Copyright Act, indicated by the "Taverne" license above, please follow below link for the End User Agreement:

[www.tue.nl/taverne](http://www.tue.nl/taverne)

***Take down policy***

If you believe that this document breaches copyright please contact us at:

[openaccess@tue.nl](mailto:openaccess@tue.nl)

providing details and we will investigate your claim.

# Cold plasma for bacterial inactivation and its effects on keratinocytes

Bouke Boekema<sup>1</sup>, Coen van Gils<sup>2</sup>, Sven Hofmann<sup>2</sup>, Peter Bruggeman<sup>2</sup>, Esther Middelkoop<sup>1,3</sup>, Gerrit Kroesen<sup>2</sup>

<sup>1</sup> Association Dutch Burn Centres, Beverwijk, The Netherlands

<sup>2</sup> Faculty of Applied Physics, Eindhoven University of Technology, The Netherlands

<sup>3</sup> Department of Plastic, Reconstructive and Hand Surgery, Research Institute MOVE, VU university medical center, Amsterdam, The Netherlands

E-mail: [bboekema@burns.nl](mailto:bboekema@burns.nl)

Plasmas have been used for many years for different applications. Atmospheric pressure plasmas deliver electrons, ions, UV radiation and an electric field, which together are effective in killing bacteria. Cold atmospheric pressure plasmas might provide additional means to reduce the bacterial load in a burn wound. It is however important to keep a balance between inactivating the bacteria and maintaining the wound healing potential. [1,2]

We studied the effects of cold argon plasma treatment on in vitro inactivation of *Pseudomonas aeruginosa*. For the treatment, a pulsed cold atmospheric plasma jet (13.56 MHz micro-jet) was used. Bacteria were diluted in culture broth (LB), phosphate buffered saline (PBS) or physiological salt (PS) and were treated with plasma for different times in microtiter plates. Growth or the absence thereof was recorded after overnight incubation. Alternatively, surviving bacteria were counted by plating dilutions.

The use of microtiter plates was a convenient way to test many different aspects in plasma treatment. However, it proved not to be useful when LB or PBS were used. Growth was absent after treatment but this appeared in part to be due to inhibition of growth rather than bacterial killing. In contrast, bacterial inactivation upon plasma treatment in PS was immediate and reached a log 6 reduction after 1 min. Our results indicated that a low pH during plasma treatment is critical. In addition, distance, duty cycle, plasma dissipated power and treatment time are important.

Because reactive radicals in cold plasma can interfere with the healing process, cell cultures of keratinocytes in PS were also treated with cold plasma. Two hours after treatment, the release of lactate dehydrogenase was determined as a measure of membrane leakage due to the treatment. In addition, activity of the cells was quantified with a tetrazolium based assay after 24 hours. Short treatments (10-60s) with argon did not result in extensive membrane damage or loss of activity when immediately after treatment culture medium was added. In conclusion, non-thermal argon plasma can be used to kill bacteria and yet preserve the viability of epidermal cells.

## Acknowledgment

This work is part of the research programme of the Foundation for Fundamental Research on Matter (FOM), which is part of the Netherlands Organisation for Scientific Research (NWO).

## References

- [1] Tupa RS and Kroesen GMW, IEEE Transactions on Plasma Science (2011), **39**, 2978-2979.
- [2] Kong MG, Kroesen GMW, Morfill G, Nosenko T, Shimizu T, van Dijk J and Zimmermann JL, New J. Phys. (2009), **11**, 115012.