

Cell Yield and Viability of Spray-On Skin Suspensions in Various Age Groups

Caitlin Stoudt

Jacob Orrico

Sigrid A. Blome-Eberwein MD

Follow this and additional works at: <https://scholarlyworks.lvhn.org/research-scholars>



Part of the [Medicine and Health Sciences Commons](#)

This Poster is brought to you for free and open access by LVHN Scholarly Works. It has been accepted for inclusion in LVHN Scholarly Works by an authorized administrator. For more information, please contact LibraryServices@lvhn.org.

Cell Yield and Viability of Spray-On Skin Suspensions in Various Age Groups

Caitlin Stoudt, Jacob Orrico, Sigrid Blome-Eberwein MD

Introduction

- RECELL® technology is used in burn surgery to create a spray-on skin solution which can be applied to burn wounds, using smaller donor sites than traditional split-thickness skin grafting (STSG).¹
- RECELL® is also believed to improve scar outcome compared to traditional meshed STSG.^{2,3}

Purpose

- Explore the relationship between patient age and percent cell viability of skin solution from RECELL® kit.

Methods

- Skin samples are taken from operating room and processed with RECELL® kit according to manufacturer instructions.
- Cells are stained with methylene blue viability assay (see Figure 3)
- Viable and nonviable cells are counted using hemocytometer (see Figure 4)
- Total cell viability counts (cells/mL) are calculated
- Percent viability of each spray-on skin suspension is determined




Figure 1
RECELL® Autologous Cell Harvesting Device:
A. Enzyme well
B. Buffer well
C. Cell suspension strainer
D. Scalpel for epidermal separation
E. Spray applicator




Figure 2
RECELL® spray applicator filled with skin suspension





Figure 3
Microscopic view of suspension stained with methylene blue on a hemocytometer. Epidermal skin cells which appear blue are nonviable, and those which appear white are viable.

Regional Burn Center
Lehigh Valley Health Network, Allentown, Pennsylvania


Methods

Calculating overall cell viability using hemocytometer data:



Results


Figure 5: Total Viable Cells in Skin Cell Suspension as a Function of Age



- So far, 11 cases have been studied.
- Study will continue until 100 skin samples from different donors have been processed.
- Our findings indicate that there is no plausible relationship between age and the overall number of viable cells.
- The total number of cells harvested from a piece of skin is largely dependent on preparation technique, thoroughness, and experience.

Results

Figure 6: Patient Age Versus Average Percent Viability of Prepared Skin Cell Suspension



- However, the results thus far indicate that there is a negative correlation between patient age and the percentage of cells which are viable in the spray-on skin suspension.
- There is a wide range of percent viability values observed so far, ranging from 32.5 % to 58.2 %.
- There is an important distinction between trends in the overall viable cell count vs. percent viability.

Discussion

- Results explain the clinical observation of slower healing of meshed skin grafts in elderly patients, which may be augmented with cell suspension treatment.³

- These results may be an important consideration in burn treatment.
- Patients of different age groups require different treatments to maximize their healing outcome.
- To account for the age trends in cell viability, providers should consider **applying a more dense coating of spray-on skin solution in older patients**, ensuring a sufficient quantity of viable epidermal cells cover the wound.




Figure 7
The expansion ratio of RECELL® is 1:80 (top) versus 1:4 in STSG (bottom)⁴

- References:**
- Holmes JH. A Brief History of RECELL® and Its Current Indications. *J Burn Care Res.* 2023;44(Suppl_1):S48-S49. doi:10.1093/jbcr/irac121
 - Holmes JH, Molnar JA, Carter JE, et al. A Comparative Study of the ReCell® Device and Autologous Split-Thickness Meshed Skin Graft in the Treatment of Acute Burn Injuries. *J Burn Care Res.* 2018;39(5):694-702. doi:10.1093/jbcr/iry029
 - Romanowski KS, Sen S. Wound healing in older adults with severe burns: Clinical treatment considerations and challenges. *Burns Open.* 2022;6(2):57-64. doi:10.1016/j.burnso.2022.01.002
 - Ren J, Liu J, Yu N, et al. The use of noncultured regenerative epithelial suspension for improving skin color and scars: A report of 8 cases and review of the literature. *J Cosmet Dermatol.* 2019;18(5):1487-1494. doi:10.1111/jocd.13071