

De La Salle University

## Animo Repository

---

Angelo King Institute for Economic and  
Business Studies

Units

---

2018

### Silver/Bismuth/Nafion Modified Pencil Graphite Electrode for Trace Heavy Metal Determination

Shirley Tiong Palisoc

Michelle Natividad

Marissa Noel

Follow this and additional works at: [https://animorepository.dlsu.edu.ph/res\\_aki](https://animorepository.dlsu.edu.ph/res_aki)



Part of the [Digestive System Diseases Commons](#)

---



# Silver/Bismuth/Nafion Modified Pencil Graphite Electrode for Trace Heavy Metal Determination

## Introduction

Frequent exposure to heavy metals can cause various health problems (Nagles Arancibia, Rojas, & Segura, 2012; Keawkim, Chuanuwatanakul, Chailapakul, & Motomizu, 2013; Xiao et al., 2014). Among the heavy metals, cadmium and lead have a great impact on biota. Cadmium can cause nephrotoxicity, and even lung cancer, when the amount of exposure to the said metal exceeds 26  $\mu\text{g}/\text{kg}$  body mass on a monthly basis (Li et al., 2011; Chen et al., 2014; Aragay & Meroçi, 2012). Lead can greatly affect the nervous system (encephalopathy) and the hormones, and it can trigger some cancers at more than 0.02-3  $\mu\text{g}/\text{kg}$  for adults and 0.03 to 9  $\mu\text{g}/\text{kg}$  for children on a daily basis (Flora, Gupta, & Tiwari, 2012).

Due to the adverse effects of heavy metals on human health, accurate and precise monitoring of these harmful metals in the environment is crucial. Thus, several techniques of detecting heavy metals have

Written by

Shirley Palisoc<sup>1,2</sup>  
Michelle Natividad<sup>1,2</sup>  
Marissa Noel<sup>3</sup>

<sup>1</sup>Condensed Matter Physics Laboratory,  
De La Salle University

<sup>2</sup>Condensed Matter Research Unit, CENSER,  
De La Salle University

<sup>3</sup>Chemistry Department,  
De La Salle University

been developed. Among these are spectroscopic and electrochemical techniques. Electrochemical techniques such as anodic stripping voltammetry (ASV) are more popular due to their low cost, high sensitivity, and portability. On the other hand, spectroscopic techniques such as atomic absorption spectroscopy, atomic emission spectroscopy, and mass spectroscopy are expensive, their availability is limited, they are not well suited for in situ measurements, and require complicated instrumentation (Willemse, Tlhomelang, Jahed, Baker, & Iwuoha, 2011).

This study developed a simple and relatively low-cost electrochemical sensor based on modified pencil graphite electrode (PGE) for fast and accurate detection of trace heavy metals in wastewater samples via anodic stripping voltammetry.

## Results and Discussion

Graphite rods (2HB) of 3 mm diameter from a commercial brand of pencils were modified by silver nanoparticles (AgNP), bismuth (Bi), and Nafion® via the drop coating technique. Characterization by field emission scanning electron microscopy confirmed the presence of the modifiers on the surface of the fabricated electrodes. The PGE modified by 3 mg AgNP and 2 mg Bi was deemed the best electrode as it yielded the highest anodic peaks as determined by ASV. The

limits of detection of the optimized electrode were found to be 0.19 parts per billion (ppb) for  $\text{Cd}^{2+}$  and 0.30 ppb for  $\text{Pb}^{2+}$ . To demonstrate the effectiveness of the fabricated electrode in sensing applications, wastewater samples were tested to determine their heavy metal content. The modified electrode was successful in determining  $\text{Cd}^{2+}$  and  $\text{Pb}^{2+}$  as well as  $\text{Cu}^{2+}$  and  $\text{Mn}^{2+}$  in the said real samples.

### **Policy Recommendation**

To protect the public from the possible adverse effects caused by heavy metal exposure, local government units and organizations involved in monitoring water pollution should carry out routine analysis using the low-cost sensor developed in this study to determine the presence of heavy metals in drinking water and aquatic systems. Moreover, government and private agencies must prioritize action plans to reduce the release of heavy metals into the environment and remedy polluted water resources.

## References

- Nagles, E., Arancibia, V., Rojas, C., & Segura, R.. (2012). Nafion-mercury coated film electrode for the adsorptive stripping voltammetric determination of lead and cadmium in the presence of pyrogallol red. *Talanta*, *99*, 119–124
- Keawkim, K., Chuanuwatanakul, S., Chailapakul, O., & Motomizu, S. (2013). Determination of lead and cadmium in rice samples by sequential injection/anodic stripping voltammetry using a bismuth film/crown ether/Nafion modified screen-printed carbon electrode. *Food Control*, *31*, 14-21.–
- Xiao, L., Xu, H., Zhou, S., Song, T., Wang, H., Li, S., . . . Yuan, Q. (2014). Simultaneous detection of Cd(II) and Pb(II) by differential pulse anodic stripping voltammetry at a nitrogen-doped microporous carbon/Nafion/bismuth-film electrode *Electrochimica Acta*, *143*, 143–151..
- Li, H., Li, J., Yang, Z., Xu, Q., Hou, C., Peng, J., Hu, X. (2011). Simultaneous determination of ultratrace lead and cadmium by square wave stripping voltammetry with in situ depositing bismuth at Nafion-medical stone doped disposable electrode. *Journal of Hazardous Materials*, *191*, 26–31..
- Chen, L., Li, Z., Meng, Y., Zhang, P., Su, Z., Liu, Y., . . . Yao, S. (2014). Sensitive square wave anodic stripping voltammetric determination of Cd<sup>2+</sup> and Pb<sup>2+</sup> ions at Bi/Nafion/overoxidized 2-mercaptoethanesulfonate-tethered polypyrrole/glassy carbon electrode. *Sensors and Actuators B: Chemical*, *191*, 94-101..
- Aragay, G., & Meroçi, A. (2012). Nanomaterials application in electrochemical detection of heavy metals. *Electrochimica Acta*, *84*, 49–61..
- Flora, G., Gupta, D., & Tiwari, A. (2012). Toxicity of lead: A review with recent updates. *Interdisciplinary Toxicology*, *5*, 47–58..
- Willemse, C., Tlhomelang, K., Jahed, N., Baker, P., & Iwuoha, E. (2011). Metallo-Graphene Nanocomposite Electrocatalytic Platform for the Determination of Toxic Metal Ions. *Sensors*, *11*, 3970–3987..

### CONTACT INFORMATION

**DLSU - Angelo King Institute for Economic and Business Studies (DLSU-AKI)**  
Room 223, St. La Salle Hall  
2401 Taft Avenue  
1004 Manila

Angelo King International Center  
Corner of Arellano Avenue and Estrada Street  
1004 Manila

+63-2-524-4611 loc. 287,  
+63-2-524-5333, +63-2-5245347 (Fax)  
<https://www.dlsu-aki.com>