

Multicenter performance evaluation of the Alinity m CMV assay for quantifying cytomegalovirus DNA in plasma samples

Lee, M.E.; Albert, E.; Wessels, E.; Kim, S.K.; Chung, H.S.; Gimenez, E.; ... ; Navarro, D.

Citation

Lee, M. E., Albert, E., Wessels, E., Kim, S. K., Chung, H. S., Gimenez, E., ... Navarro, D. (2023). Multicenter performance evaluation of the Alinity m CMV assay for quantifying cytomegalovirus DNA in plasma samples. *Journal Of Clinical Microbiology*, *61*(10). doi:10.1128/jcm.00415-23

Version:Publisher's VersionLicense:Creative Commons CC BY 4.0 licenseDownloaded from:https://hdl.handle.net/1887/3656233

Note: To cite this publication please use the final published version (if applicable).



8 | Virology | Full-Length Text

Check for updates

Multicenter performance evaluation of the Alinity m CMV assay for quantifying cytomegalovirus DNA in plasma samples

Miae Lee,¹ Eliseo Albert,² Els Wessels,³ Soo-Kyung Kim,¹ Hae-Sun Chung,¹ Estela Giménez,² Tom Vreeswijk,³ Eric C. J. Claas,³ Yan Chin Tai,⁴ Birgit Reinhardt,⁵ Mark M. Sasaki,⁶ David Navarro²

AUTHOR AFFILIATIONS See affiliation list on p. 8.

ABSTRACT Monitoring of cytomegalovirus (CMV) viral load is critical for informing treatment decisions in order to prevent the severe health consequences of CMV infection or reactivation of latent CMV in immunocompromised individuals. This first field evaluation examined the analytical and clinical performance of the Alinity m CMV assay. Analytical performance was assessed with a commercially available six-member panel, while the clinical performance evaluation compared the Alinity m CMV assay to the RealTime CMV assay and a laboratory-developed test (LDT) as the test of record at three large hospital-based clinical laboratories. Precision of the Alinity m CMV assay was demonstrated with total standard deviation (SD) between 0.08 and 0.28 Log IU/mL. A total of 457 plasma specimens were tested on the Alinity m CMV assay and compared to the test of record at each site (n = 304 with RealTime CMV and n = 153 with LDT CMV). The Alinity m CMV assay had excellent correlation (correlation coefficient $r \ge 0.942$) in comparison to the RealTime CMV or LDT CMV assays. The mean observed bias ranged from -0.03 to 0.34 Log IU/mL. Median onboard turnaround time of Alinity m CMV was less than 3 h. When the CMV assay is run on the Alinity m system, it has the capacity to shorten time to result and, therefore, to therapy.

KEYWORDS nucleic acid amplification test, DNA, immunocompromised, transplantation, high-throughput diagnostic assay

C ytomegalovirus (CMV) is a member of the herpesvirus family with a high seroprevalence worldwide (1). Primary infection can cause mild or subclinical disease in immunocompetent individuals and establishes life-long latency (2). Primary infection or reactivation of latent CMV has severe and life-threatening consequences in newborns and in patients who are immunocompromised, including the elderly, patients with AIDS, and stem cell and solid organ transplant recipients (2–4). Monitoring of CMV DNA load in plasma or blood of transplant recipients is critical to initiate preemptive therapy and, thus, prevent the occurrence of CMV end-organ disease (5) as CMV infection is associated with graft rejection and decreased patient survival (6).

Nucleic acid amplification tests (NAATs) have high sensitivity and specificity for target pathogens such as CMV and provide accurate quantitation to inform appropriate patient management (7). Many laboratories had established laboratory developed tests (LDTs) long time ago when commercial CMV assays were rare. However, facing a high testing burden and increased quality and documentation requirements, especially with regard to the new European *In vitro* Diagnostic Regulation, laboratories are looking for consolidating their CMV and other PCR assays on a single platform with high quality assays and low turnaround time (TAT). The Alinity m CMV assay (Abbott Molecular Inc., Des Plaines, IL, USA) (8) is a quantitative PCR assay that is run on the fully automated, continuous, random-access Alinity m analyzer with a processing capacity of 300 samples in approximately 8 h and a time to first result of <115 min. The purpose of this study

Editor Elitza S. Theel, Mayo Clinic, Rochester, Minnesota, USA

Address correspondence to Miae Lee, miae@ewha.ac.kr.

M.L. received grant funding from Abbott. E.A. received a travel grant from Abbott. B.R. is an employee of Abbott GmbH, Wiesbaden, Germany. M.M.S. is an employee of Abbott Molecular Inc., Des Plaines, IL, USA. Y.C.T. is an employee of Abbott Laboratories (Singapore) Pte. Ltd., Singapore. D.N. received grant funding and conference fees from Abbott. E.C.J.C., E.G., E.W., H.S.C., S.K.K., and T.V. have no conflicts of interest to disclose. This study was funded by Abbott Laboratories.

See the funding table on p. 8.

See the companion article at https://doi.org/ 10.1128/jcm.00472-23.

Received 4 April 2023 Accepted 30 June 2023 Published 20 September 2023

Copyright © 2023 Lee et al. This is an open-access article distributed under the terms of the Creative Commons Attribution 4.0 International license.

was to evaluate the analytical and clinical performance of the Alinity m CMV assay and compare it to the RealTime CMV assay (9) and an LDT CMV in plasma specimens (10).

MATERIALS AND METHODS

Assessment of the analytical performance

The analytical precision of the Alinity m CMV assay was assessed using a six-member CMV verification panel (Exact Diagnostics, Fort Worth, TX, USA) containing non-infectious intact whole virus at 2.3–6.6 Log IU/mL in ETDA plasma. For each site, three replicates per panel member were tested over 5 days (total n = 45 with 15 replicates per laboratory).

Reproducibility was assessed by evaluating the performance of the assay quality controls [QCs; low positive control (LPC) and high positive control (HPC)] from multiple reagent lots across testing sites. For each site, one replicate of HPC and LPC was tested over 5 days (total n = 15 with five replicates per laboratory).

Clinical performance evaluation and specimens

The performance of the Alinity m CMV assay (Abbott Molecular Inc., Des Plaines, IL, USA) was compared to two CMV assays in a multicenter, international study performed at three independent International Standard Organization accredited clinical laboratories: Ewha Womans University Mokdong Hospital, Seoul, South Korea; Hospital Clinico Universitario de Valencia, Valencia, Spain; and Leiden University Medical Center, Leiden, the Netherlands.

The study was performed in accordance with the principles of Good Clinical Practice and conducted in adherence with the Declaration of Helsinki. Only surplus remnant patient plasma specimens were used for this study. All specimens were anonymized before study initiation, and an identification number containing no patient identifiers was assigned to each remnant specimen. An approval by an ethics committee was obtained according to the institutional requirement.

Alinity m CMV assay performance was evaluated across 457 remnant clinical plasma specimens that were sent to the participating laboratories for routine CMV DNA load determination. They were collected from solid organ and hematopoietic stem cell transplant recipients and stored frozen at \leq -70°C for 0.5–6 years. The clinical specimens were retrospectively tested using the Alinity m CMV assay and the results were compared to the historical data of one of the two routine CMV assays based on real-time PCR, the RealTime CMV assay (Abbott Molecular Inc., Des Plaines, IL, USA) or an LDT CMV (Leiden, the Netherlands).

Alinity m CMV assay performance was also evaluated across the longitudinal course for 26 de-identified patients previously tested with LDT CMV. Patients were monitored between 4 and 296 days with a range of 2–10 times.

Molecular CMV assays and workflow evaluation

The Alinity m CMV assay is a dual target assay (UL34 and UL80.5) and standardized to the first WHO International Standard (11) with a quantitative range from 1.48 to 8.00 Log IU/mL.

The RealTime CMV assay (9) is also a dual target assay (UL34 and UL80.5) and standardized to the first WHO International Standard (11). It is run on the *m*2000 platform, including automated sample preparation on the *m*2000sp and PCR on *m*2000rt. The assay has a quantitative range from 1.49 to 8.19 Log IU/mL. Its analytical and clinical performance were evaluated in prior studies (12, 13).

The LDT CMV assay is based on the previously published PCR assay utilizing primers and probes to amplify a 126-bp fragment from the CMV immediate-early antigen region, with minor modifications. DNA isolation was performed using the MagNA Pure 96 platform with the DNA and Viral NA Small Volume Kit 2.0 using the Pathogen Universal 200 protocol (Roche Diagnostics, Almere, the Netherlands). The PCR reaction was carried out using the HotStar Taq master mix (Qiagen, Hilden, Germany) in the CFX96 real-time detection system (Bio-Rad, Veenendaal, the Netherlands) (10). The LDT CMV assay has been used in routine clinical practice for nearly 20 years, is standardized to the first WHO International Standard (11), and quantifies results between 1.7 and 7.0 Log IU/mL. In clinical routine, results below 250 IU/mL are reported as <2.4 Log IU/mL since quantitation of viral loads <250 IU/mL is not considered clinically relevant in our setting.

Onboard and processing TATs of the Alinity m system were evaluated based on the automatically documented timepoints of loading samples, sample aspiration and result reporting by Alinity m.

Statistical analysis

All analyses were performed using PC SAS (version 9.4) software (SAS, Cary, NC, USA). Relationships between quantitative variables were studied by means of Deming regression. Bland-Altman analysis was performed to evaluate the differences in quantification between the assays.

The following analysis was performed for each instrument and each panel member: The PROC MIXED procedure with the MIVQUE0 option in SAS was used to produce variance components for the model used in the analysis. The point estimates of the means and standard deviations (SDs) were reported. The SD was estimated for the within-day component, the between-day component, and the between-site component for each instrument and each panel member. All the effects were considered as random for the analyses. Any negative variance components were set to zero for these calculations. The total assay variability was defined as the sum of the within-day (residual error) component, the between-day component, and the between-site component estimates of variability. The following statistics were reported: *N*, mean, within-day SD, betweenday SD, between-site SD, and total SD. For the evaluation of the quality controls, a within-day component was not included in the analysis as only one replicate of each HPC and LPC was tested per day and site.

RESULTS

Analytical performance

Analytical precision of the Alinity m CMV assay was assessed across the three study sites by testing a commercially available plasma CMV verification panel. The observed quantitation was similar to the expected concentrations for all panels tested with a mean bias of 0.04 Log IU/mL. SD measure of each level tested was less than or equal to 0.28 Log IU/mL (Table 1). Reproducibility was characterized by a total SD of 0.12 Log IU/mL for either LPC and HPC (Table 2).

Difference Within-day Between-day Between-site Mean component component component Total Panel Target conc. conc. mean-target (Log IU/mL) (Log IU/mL) SD (Log IU/mL) SD (Log IU/mL) SD (Log IU/mL) SD (Log IU/mL) member Ν (Log IU/mL) 45 2.30 2.33 0.03 0.28 1 0.28 0.00 0.00 2 45 2.60 2.61 0.01 0.10 0.00 0.01 0.10 3 45 3.60 3.60 0.00 0.07 0.01 0.04 0.08 4 45 4.60 4.62 0.02 0.05 0.00 0.06 0.08 5 45 5.60 5.67 0.07 0.04 0.00 0.07 0.08 6 45 6.60 0.12 0.06 0.02 0.09 6.72 0.11

TABLE 1 Precision of Alinity m CMV assay using CMV verification panels tested across three laboratories (n = 45 with 15 replicates per laboratory)

^oOne replicate of panel member 1 was an outlier. If this outlier was removed from analysis, SD would be 0.17 Log IU/mL.

		Target conc. range	Mean conc.	Between-day component	Between-site component	Total
Control	Ν	(Log IU/mL)	(Log IU/mL)	SD (Log IU/mL)	SD (Log IU/mL)	SD (Log IU/mL)
HPC	15	5.99–6.00	6.01	0.09	0.08	0.12
LPC	15	3.14–3.30	3.24	0.04	0.11	0.12

TABLE 2 Reproducibility of testing Alinity m CMV quality controls across three laboratories using two lots of quality controls (n = 15 with five replicates per laboratory)^a

^aHPC, high-positive control; LPC, low-positive control.

Clinical performance

Overall, 457 plasma specimens were tested on the Alinity m CMV assay and the results were compared against one of the two comparator CMV assays: the RealTime CMV (304 specimens) or an LDT CMV assay (153 specimens).

A total of 304 remnant clinical specimens previously tested with the RealTime CMV assay were retested using the Alinity m CMV assay. The overall observed qualitative agreement, calculated as concordant negative and positive results (<LLOQ or quantitated results) between the two assays, was 89.1% (271/304) with a Cohen's kappa value of 0.66 representing substantial agreement (14) (Table 3). Of the 304 specimens tested on the Alinity m CMV assay and the RealTime CMV assay, 174 fell within the analytical measuring range (AMR) for both assays. The correlation coefficient was 0.942 (Deming regression equation, y = 1.06x - 0.20) and the mean observed bias -0.03 Log IU/mL (Bland-Altman analysis, Alinity m CMV – RealTime CMV, Fig. 1A and B).

The overall observed qualitative agreement between the Alinity m CMV and the LDT CMV assays was 96.7% (148/153) with a Cohen's kappa value of 0.87 representing almost perfect agreement (Table 4). Of the 153 specimens, 127 fell within the AMR of both assays. The correlation coefficient *r* was 0.959 (Deming regression equation, y = 1.19x –0.28) and the mean observed bias +0.34 Log IU/mL (Bland-Altman analysis, Alinity m CMV – LDT CMV, Fig. 1C and D).

Of the above 153 specimens, 91 tests were part of longitudinal analysis performed on 26 patients. After the baseline test, each patient underwent subsequent testing between 4 and 296 days with a range of 2–10 times. CMV viral load kinetics were very similar between Alinity m CMV and LDT CMV assays. Monitoring results for patients with at least three time points after baseline testing are shown in Fig. 2A through F.

Workflow analysis

The Alinity m systems allowed random and continuous loading of CMV study samples side by side with samples of routine assays that were processed simultaneously on the systems. In this real-world setting, observed median onboard TAT (from placement of the specimen on the analyzer to result reporting) for the Alinity m CMV assay was 2 h 54 min across the three study sites (ranging from 2 h 7 min to 4 h 41 min). Processing TAT (from sample aspiration to result reporting) was 112–116 min with 97% of all results being reported within 114 min, across all three study sites.

TABLE 3 Alinity m CMV clinical specimen agreement of results determined by the RealTime CMV assay

		RealTime CMV			
		Not detected	<lloq< th=""><th>Quantitated</th><th>Total</th></lloq<>	Quantitated	Total
	Not Detected	43	20	8 ^{<i>a</i>}	71
Alinity on CAAV Account	<lloq< td=""><td>5</td><td>25</td><td>26^b</td><td>56</td></lloq<>	5	25	26 ^b	56
Annity III Civiv Assay	Quantitated	0	3 ^c	174	177
	Total	48	48	208	304

^dEight specimens Not Detected by Alinity m CMV assay had a range of 1.49–1.93 Log IU/mL with the RealTime CMV assay.

^bTwenty-six specimens <LLOQ by Alinity m CMV had a range of 1.52–2.16 Log IU/mL with the RealTime CMV assay. Three specimens <LLOQ by RealTime CMV had a range of 1.57–1.81 Log IU/mL with the Alinity m CMV assay.

TABLE 4	Alinity m CMV clinical specime	n agreement of results determined b	y the LDT CMV assa
---------	--------------------------------	-------------------------------------	--------------------

		LI	LDT CMV		
		Not detected	Quantitated	Total	
	Not Detected	21	1 ^{<i>a</i>}	22	
Alinity m CMN/ Accov	<lloq< td=""><td>2</td><td>0</td><td>2</td></lloq<>	2	0	2	
Annity in Civiv Assay	Quantitated	2 ^{<i>b</i>}	127	129	
	Total	25	128	153	

^aOne specimen Not Detected by Alinity m CMV assay was quantitated at 1.81 Log IU/mL with the LDT CMV assay. ^bTwo specimens Not Detected by LDT CMV assay showed 1.59 and 1.61 Log IU/mL with the Alinity m CMV assay.

DISCUSSION

Evaluation of a CMV verification panel was performed across the three study sites using two lots of amplification reagents, three lots of lysis buffer, and three lots of sample preparation reagents. Despite testing with different instruments, operators, and reagent lots, the precision of the Alinity m CMV assay was high as shown by an SD of less than or equal to 0.28 Log IU/mL at each level tested across the AMR. One replicate of panel member 1 was an outlier which resulted in higher than expected SD. If this outlier was removed from analysis, SD would be 0.17 Log IU/mL. Unfortunately, no residual volume was left to perform any additional testing of this replicate. The observed means and total SDs for the Alinity m CMV assay HPC and LPC were 6.01 ± 0.12 and 3.24 ± 0.12 Log IU/mL, respectively.

This multicenter, international study demonstrated excellent correlation (correlation coefficient ≥ 0.942) between Alinity m CMV and currently used CMV assays in plasma. The Alinity m CMV assay demonstrated comparable quantitation to its predecessor, the RealTime CMV assay (bias -0.03 Log IU/mL). Overall bias between Alinity m CMV and the LDT CMV assays was 0.34 Log IU/mL, which is considered not significant in the context of biological variation and taking into account that residual stored specimens were used. This difference in quantitation could also be attributed to the intrinsic differences in the assay design such as sample extraction efficiency differences or amplicon size. Since the majority of the discordant results between the Alinity m and comparator CMV assays were at the low end of the AMR or below LLOQ of one or both assays, imprecision may also contribute to the difference in quantitation.

Of note, there was one sample with a difference of 1.75 Log IU/mL between LDT CMV and Alinity m CMV results (3.53 Log IU/mL with LDT and 5.28 Log IU/mL with Alinity m). Repeat testing with LDT CMV yielded a similar result as before (3.58 Log IU/mL), while resolution testing with the RealTime CMV assay (5.28 Log IU/mL) confirmed the Alinity m CMV assay result. Since the internal control that is present in every sample in the LDT CMV fulfills the criteria, inhibition of the sample is very unlikely. A possible explanation might be a mutation within the target region of LDT CMV, since this could impact the quantitation of this sample. The Alinity m CMV and RealTime CMV assays are dual target assays (UL34 and UL80.5), whereas the LDT CMV targets a single site in the CMV immediate-early antigen region. Another explanation could be a considerable amount of fragmented nucleic acid in this specimen that may be detected by the Alinity m CMV and RealTime CMV assays due to their slightly smaller amplicon size compared to LDT CMV (15).

Limitations of our study include the use of surplus samples for the assay comparison, which in some cases did not provide sufficient sample volume for additional analysis. Overall, longitudinal comparison of patient courses comparing LDT and Alinity m confirmed the high correlation between the two assays.

Current routine procedure of sample collection and batchwise testing with a frequency between twice a day and once a week and onboard TATs between 3.5 and 9 h (12, 16, 17) lead to an average total TAT of 1–8 days from sample receipt in the laboratory to result reporting across the three study sites. In contrast, random and continuous access capabilities of Alinity m without the need for batching combined with a median



FIG 1 Comparison of the Alinity m CMV assay and test of record CMV assay performance with clinical plasma specimens. Deming regression of CMV levels showing correlation between the Alinity m CMV assay and (A) the RealTime CMV assay on *m*2000 or (C) an LDT CMV assay. Bland-Altman analysis showing mean bias between the Alinity m CMV assay and (B) the RealTime CMV assay or (D) LDT CMV assay. The middle line indicates the mean bias, the lines above and below indicate ±1.96 × SD.



FIG 2 Comparison of longitudinal clinical plasma specimen results of Alinity m CMV and an LDT CMV assay. CMV viral load kinetics shown for six patients with \geq 3 tests. (A) Patient 1 tested at days 1, 47, 54, 74, 79, 92, 96, 110, 116, and 131. (B) Patient 2 tested at days 1, 29, 43, 57, 169, 204, and 260. (C) Patient 3 tested at days 1, 6, 28, 51, 122, and 128. (D) Patient 4 tested at days 1, 36, 50, 61, and 138. (E) Patient 5 tested at days 1, 110, 129, 144, and 235. (F) Patient 6 tested at days 1, 5, 8, and 33.

onboard TAT of less than 3 h enable same-day reporting of actionable results which is important for transplant patient management. Similar rapid and homogenous TAT results have already been reported for other Alinity m assays (16, 17).

In conclusion, this multicenter evaluation of the newly developed real-time PCRbased Alinity m CMV assay showed high precision, reproducibility, as well as accurate quantitation of CMV using verification panels and clinical plasma samples from solid organ or hematopoietic stem cell transplant recipients. The fully automated Alinity m platform enables same-day reporting of CMV test results allowing shortening the time between diagnosis and treatment and, thus, may improve patient care and outcome.

ACKNOWLEDGMENTS

The authors would like to thank Min-Kyung So (Ewha Womans University College of Medicine), Danijela Lucic (Abbott Molecular Inc.), and Jens Dhein (Abbott GmbH) for their assistance in the study.

This study was funded by Abbott Laboratories.

M.L.: Methodology, manuscript review, and editing, visualization, supervision, project administration, funding acquisition. E.A.: Validation, investigation, data curation, manuscript review and editing. E.W.: Methodology, validation, investigation, manuscript review and editing. S-K.K.: Data validation, review, and editing, visualization, IRB approval. H-S.C.: Investigation, resources, manuscript review, and editing. E.G.: Validation, investigation, data curation, manuscript review and editing. T.V.: Validation, investigation, manuscript review and editing. Y.C.T.: Data analysis, manuscript writing, review and editing. B.R.: Conceptualization, data analysis, manuscript writing, review and editing. M.M.S.: Data analysis, manuscript writing, review and editing. M.M.S.: pata analysis, manuscript review and editing. D.N.: Methodology, resources, supervision, manuscript review and editing.

AUTHOR AFFILIATIONS

¹Ewha Womans University College of Medicine, Seoul, Republic of Korea
²Hospital Clinico Universitario de Valencia, Valencia, Spain
³Leiden University Medical Center, Leiden, the Netherlands
⁴Abbott Laboratories (Singapore) Pte. Ltd., Singapore
⁵Abbott GmbH, Wiesbaden, Germany
⁶Abbott Molecular Inc., Des Plaines, Illinois, USA

AUTHOR ORCIDs

Miae Lee ⁽¹⁾ http://orcid.org/0000-0001-9140-3814 Yan Chin Tai ⁽¹⁾ http://orcid.org/0009-0007-2541-2617 Birgit Reinhardt ⁽¹⁾ http://orcid.org/0009-0007-7045-1105

FUNDING

Funder	Grant(s)	Author(s)
Abbott Laboratories (Abbott)		Miae Lee
Abbott Laboratories (Abbott)		David Navarro

AUTHOR CONTRIBUTIONS

Miae Lee, Funding acquisition, Methodology, Project administration, Supervision, Visualization, Writing – review and editing | Eliseo Albert, Data curation, Investigation, Validation, Writing – review and editing | Els Wessels, Investigation, Methodology, Validation, Writing – review and editing | Soo-Kyung Kim, Validation, Visualization, Writing – review and editing | Bae-Sun Chung, Investigation, Resources, Writing – review and editing | Estela Giménez, Data curation, Investigation, Validation, Writing – review and editing | Tom Vreeswijk, Investigation, Validation, Writing – review and editing | Eric C. J. Claas, Methodology, Resources, Supervision, Writing – review and editing | Birgit Reinhardt, Conceptualization, Formal analysis, Software, Writing – original draft, Writing – review and editing | Mark M. Sasaki, Formal analysis, Writing – original draft, Writing

- review and editing | David Navarro, Methodology, Resources, Supervision, Writing - review and editing

REFERENCES

- Zuhair M, Smit GSA, Wallis G, Jabbar F, Smith C, Devleesschauwer B, Griffiths P. 2019. Estimation of the worldwide seroprevalence of cytomegalovirus: a systematic review and meta-analysis. Rev Med Virol 29:e2034. https://doi.org/10.1002/rmv.2034
- Griffiths P, Baraniak I, Reeves M. 2015. The pathogenesis of human cytomegalovirus. J Pathol 235:288–297. https://doi.org/10.1002/path. 4437
- Boeckh M, Geballe AP. 2011. Cytomegalovirus: pathogen, paradigm, and puzzle. J Clin Invest 121:1673–1680. https://doi.org/10.1172/JCl45449
- Camargo JF, Komanduri KV. 2017. Emerging concepts in cytomegalovirus infection following hematopoietic stem cell transplantation. Hematol Oncol Stem Cell Ther 10:233–238. https://doi.org/10.1016/j. hemonc.2017.05.001
- Kotton CN, Kumar D, Caliendo AM, Huprikar S, Chou S, Danziger-Isakov L, Humar A, The Transplantation Society International CMV Consensus Group. 2018. The third international consensus guidelines on the management of cytomegalovirus in solid-organ transplantation. Transplantation 102:900–931. https://doi.org/10.1097/TP.-000000000002191
- Razonable RR, Humar A. 2019. Cytomegalovirus in solid organ transplant recipients - Guidelines of the American Society of Transplantation Infectious Diseases Community of Practice. Clin Transplant 33:e13512. https://doi.org/10.1111/ctr.13512
- Razonable RR, Inoue N, Pinninti SG, Boppana SB, Lazzarotto T, Gabrielli L, Simonazzi G, Pellett PE, Schmid DS. 2020. Clinical diagnostic testing for human cytomegalovirus infections. J Infect Dis 221:S74–S85. https://doi. org/10.1093/infdis/jiz601
- 8. Abbott Molecular Inc. 2021. Alinity m CMV AMP kit [package insert]
- 9. Abbott Molecular Inc. 2017. RealTime CMV assay [package insert]
- Kalpoe JS, Kroes ACM, de Jong MD, Schinkel J, de Brouwer CS, Beersma MFC, Claas ECJ. 2004. Validation of clinical application of cytomegalovirus plasma DNA load measurement and definition of treatment criteria by analysis of correlation to antigen detection. J Clin Microbiol 42:1498– 1504. https://doi.org/10.1128/JCM.42.4.1498-1504.2004

- Fryer JF, Heath AB, Wilkinson DE, Minor PD, Collaborative Study Group. 2016. A collaborative study to establish the 1st WHO International Standard for human cytomegalovirus for nucleic acid amplification technology. Biologicals 44:242–251. https://doi.org/10.1016/j. biologicals.2016.04.010
- Clari MÁ, Bravo D, Costa E, Muñoz-Cobo B, Solano C, Remigia MJ, Giménez E, Benmarzouk-Hidalgo OJ, Pérez-Romero P, Navarro D. 2013. Comparison of the new Abbott RealTime CMV assay and the Abbott CMV PCR Kit for the quantitation of plasma cytomegalovirus DNAemia. Diagn Microbiol Infect Dis 75:207–209. https://doi.org/10.1016/j. diagmicrobio.2012.10.010
- Tremblay M-A, Rodrigue M-A, Deschênes L, Boivin G, Longtin J. 2015. Cytomegalovirus quantification in plasma with Abbott RealTime CMV and Roche Cobas Amplicor CMV assays. J Virol Methods 225:1–3. https:// doi.org/10.1016/j.jviromet.2015.08.010
- Landis JR, Koch GG. 1977. The measurement of observer agreement for categorical data. Biometrics 33:159–174. https://doi.org/10.2307/ 2529310
- Preiksaitis JK, Hayden RT, Tong Y, Pang XL, Fryer JF, Heath AB, Cook L, Petrich AK, Yu B, Caliendo AM. 2016. Are we there yet? Impact of the first international standard for cytomegalovirus DNA on the harmonization of results reported on plasma samples. Clin Infect Dis 63:583–589. https:/ /doi.org/10.1093/cid/ciw370
- Galindo LT, Hristov AD, Gentil LG, Scarpelli L, Santiago J, Levi JE. 2021. Performance evaluation of the fully automated molecular system Alinity m in a high-throughput central laboratory. J Clin Virol 137:104786. https: //doi.org/10.1016/j.jcv.2021.104786
- Obermeier M, Pacenti M, Ehret R, Onelia F, Gunson R, Goldstein E, Chevaliez S, Vilas A, Glass A, Maree L, Krügel M, Knechten H, Braun P, Naeth G, Azzato F, Lucic D, Marlowe N, Palm MJ, Pfeifer K, Reinhardt B, Dhein J, Joseph AM, Martínez-García L, Galán J-C. 2020. Improved molecular laboratory productivity by consolidation of testing on the new random-access analyzer Alinity m. J Lab Med 44:319–328. https:// doi.org/10.1515/labmed-2020-0102