



Apr 25th, 10:30 AM - 4:30 PM

In Situ Hybridization Screening of a Lambda Library for Chromosome Specific DNA: The Isolation of a Rabbit Chromosome Number 9 Probe

Nicole Zook

Illinois Wesleyan University

Leonard Maroun, Faculty Advisor

Southern Illinois University Carbondale

Follow this and additional works at: <http://digitalcommons.iwu.edu/jwprc>

Nicole Zook and Leonard Maroun, Faculty Advisor, "In Situ Hybridization Screening of a Lambda Library for Chromosome Specific DNA: The Isolation of a Rabbit Chromosome Number 9 Probe" (April 25, 1992). *John Wesley Powell Student Research Conference*. Paper 1.

<http://digitalcommons.iwu.edu/jwprc/1992/posters/1>

This Event is brought to you for free and open access by The Ames Library, the Andrew W. Mellon Center for Curricular and Faculty Development, the Office of the Provost and the Office of the President. It has been accepted for inclusion in Digital Commons @ IWU by the faculty at Illinois Wesleyan University. For more information, please contact digitalcommons@iwu.edu.

©Copyright is owned by the author of this document.

1 IN SITU HYBRIDIZATION SCREENING OF A LAMBDA LIBRARY FOR CHROMOSOME SPECIFIC DNA: THE ISOLATION OF A RABBIT CHROMOSOME NUMBER 9 PROBE

2 Nicole Zook, Dept. of Biology, IWU/Dept. of Medical Microbiology & Immunology, Southern Illinois University School of Medicine, Leonard Maroun*

3 Down's syndrome, a disease which affects 1 of every 783 newborns, is the result of aneuploidy of human chromosome 21. The closely related rabbit chromosome 20 can be used as a model for the study of the cause of the misdistribution of this chromosome during meiosis. The procedure, called karyotyping, requires long hours behind the microscope; however, it can be replaced with the more accurate, time-efficient restriction fragment length polymorphism (RFLP) assay. The purpose of this study was to test the feasibility of finding a chromosome specific probe by searching a lambda library using IN SITU hybridization. IN SITU hybridization involves annealing a radioactive DNA probe directly to metaphase chromosomes fixed to a microscope slide. A photographic emulsion is spread over the slides. Photographic grains appear over the chromosome to which the DNA has annealed thus identifying the chromosome from which the probe DNA originated (FIG. 1). Probe DNA was phenol extracted from lambda phage virus clones that carry rabbit DNA inserts. Inserts were purified by gel electrophoresis (FIG. 2). Using this procedure, a chromosome 9 probe was discovered in the 86 phage DNA clones isolated (FIG. 3). This suggests that a larger search would be capable of finding the chromosome 20 probe needed to develop a rabbit model for Down's syndrome.