



Illinois Wesleyan University Digital Commons @ IWU

John Wesley Powell Student Research
Conference

2005, 16th Annual JWP Conference

Apr 16th, 9:00 AM - 10:00 AM

Quantification of Fluoroquinolone Antibiotics with High Performance Liquid Chromatography Coupled UV/VIS Detection

Nathan Smoter

Illinois Wesleyan University

Stephen Hoffmann, Faculty Advisor

Illinois Wesleyan University

Follow this and additional works at: <http://digitalcommons.iwu.edu/jwprc>

Nathan Smoter and Stephen Hoffmann, Faculty Advisor, "Quantification of Fluoroquinolone Antibiotics with High Performance Liquid Chromatography Coupled UV/VIS Detection" (April 16, 2005). *John Wesley Powell Student Research Conference*. Paper 30.

<http://digitalcommons.iwu.edu/jwprc/2005/posters/30>

This Event is brought to you for free and open access by The Ames Library, the Andrew W. Mellon Center for Curricular and Faculty Development, the Office of the Provost and the Office of the President. It has been accepted for inclusion in Digital Commons @ IWU by the faculty at Illinois Wesleyan University. For more information, please contact digitalcommons@iwu.edu.

©Copyright is owned by the author of this document.

Poster Presentation P51

**QUANTIFICATION OF FLUOROQUINOLONE ANTIBIOTICS WITH
HIGH PERFORMANCE LIQUID CHROMATOGRAPHY
COUPLED UV/VIS DETECTION**

Nathan Smoter and Stephen Hoffmann*
Chemistry Department, Illinois Wesleyan University

Fluoroquinolones such as Ciprofloxacin, Cipro[®], are a class of widely used antibiotics that inhibit DNA Gyrase of bacteria. These compounds are used in both human and veterinary medicines. Due to increasing number of antibiotic resistant bacterial strains, their presence in wastewaters and soil is important. The antibiotic, when administered, is not fully absorbed by the body; the excess is excreted. In the case of animals, the excess is excreted directly into the soil, while for humans the compounds enter the wastewater stream. If they are not removed in wastewater treatment, they may be released into the environment. The analytical method centered on reverse-phase high performance liquid chromatography on a C-18 column coupled with UV/Vis detection. The mobile phase was composed of acetonitrile and 25mM o-H₃PO₄. Fluoroquinolone compounds were visualized at 278nm. Solid phase extraction can be done using a mixed phase cation exchange cartridge and eluted with MeOH and o-H₃PO₄. These compounds are detected at such low levels, ng/L, that it becomes critical to avoid contamination; this study explores contamination possibilities.