

Illinois Wesleyan University Digital Commons @ IWU

John Wesley Powell Student Research Conference

2004, 15th Annual JWP Conference

Apr 17th, 9:00 AM - 10:00 AM

A Novel Method for Super-Resolution Optical Microscopy

Brett Parmenter
Illinois Wesleyan University

Gabriel C. Spalding, Facutly Advisor *Illinois Wesleyan University*

Follow this and additional works at: http://digitalcommons.iwu.edu/jwprc

Brett Parmenter and Gabriel C. Spalding, Facutly Advisor, "A Novel Method for Super-Resolution Optical Microscopy" (April 17, 2004). *John Wesley Powell Student Research Conference*. Paper 21. http://digitalcommons.iwu.edu/jwprc/2004/posters/21

This Event is brought to you for free and open access by The Ames Library, the Andrew W. Mellon Center for Curricular and Faculty Development, the Office of the Provost and the Office of the President. It has been accepted for inclusion in Digital Commons @ IWU by the faculty at Illinois Wesleyan University. For more information, please contact digitalcommons@iwu.edu. ©Copyright is owned by the author of this document.

THE JOHN WESLEY POWELL STUDENT RESEARCH CONFERENCE - APRIL 2004

Poster Presentation P37

A NOVEL METHOD FOR SUPER-RESOLUTION OPTICAL MICROSCOPY

Brett Parmenter and Gabriel C. Spalding*
Department of Physics, Illinois Wesleyan University
Department of Physics, Northwestern University

Even neglecting aberrations, conventional optical microscopes are limited (essentially by the uncertainty principle) to a resolution on the order of the wavelength of light used. The Near-Field Scanning Optical Microscope (NSOM) is able to overcome this fundamental limitation, by using evanescent fields. The most common NSOM design uses a fiber optic pulled down to a very fine tip with a sub-wavelength aperture at the end. The resolution of this type of microscope is limited only by the size of the aperture used. Unfortunately, as shown by Hans Bethe, the intensity of the light allowed through a sub-wavelength aperture falls off as the sixth power of the radius. Our alternative uses an optically trapped, resonant nanoparticle as a sub-wavelength, scanned probe. We expect to achieve super-resolution with improved signal-to-noise.