Digital health technologies and machine learning augment patient reported outcomes to remotely characterise rheumatoid arthritis

Andrew P. Creagh^{*1, 2}, Valentin Hamy³, Hang Yuan^{2,4}, Gert Mertes^{1,2,4}, Ryan Tomlinson⁵, Wen-Hung Chen⁵, Rachel Williams⁵, Christopher Llop⁶, Christopher Yee⁶, Mei Sheng Duh⁶, Aiden Doherty^{†,2,4}, Luis Garcia-Gancedo^{†, 3}, and David A. Clifton^{†, 1}

¹Institute of Biomedical Engineering, Department of Engineering Science, University of Oxford, UK;

²Big Data Institute, University of Oxford, UK;

³Value Evidence and Outcomes (VEO), GSK, UK;

⁴Nuffield Department of Population Health, University of Oxford, UK;

⁵Value Evidence and Outcomes (VEO), GSK, US;

⁶Analysis Group (AG), Boston, MA, USA.

ABSTRACT

Digital measures of health status captured during daily life could greatly augment current in-clinic assessments for rheumatoid arthritis (RA), to enable better assessment of disease progression and impact. This work presents results from weaRAble-PRO. a 14-day observational study, which aimed to investigate how digital health technologies (DHT), such as smartphones and wearables, could augment patient reported outcomes (PRO) to determine RA status and severity in a study of 30 moderate-tosevere RA patients, compared to 30 matched healthy controls (HC). Sensor-based measures of health status, mobility, dexterity, fatigue, and other RA specific symptoms were extracted from daily iPhone guided tests (GT), as well as actigraphy and heart rate sensor data, which was passively recorded from patients' Apple smartwatch continuously over the study duration. We subsequently developed a machine learning (ML) framework to distinguish RA status and to estimate RA severity. It was found that daily wearable sensor-outcomes robustly distinguished RA from HC participants (F1, 0.807). Furthermore, by day 7 of the study (half-way), a sufficient volume of data had been collected to reliably capture the characteristics of RA participants. In addition, we observed that the detection of RA severity levels could be improved by augmenting standard patient reported outcomes with sensor-based features (F1, 0.833) in comparison to using PRO assessments alone (F1, 0.759), and that the combination of modalities could reliability measure continuous RA severity, as determined by the clinician-assessed RAPID-3 score at baseline (r^2 , 0.692; RMSE, 1.33). The ability to measure the impact of disease during daily life—through objective and remote digital outcomes-paves the way forward to enable the development of more patient-centric and personalised measurements for use in RA clinical trials.

1 Introduction

Rheumatoid arthritis (RA) patients follow subtle and unpredictable disease courses, patient-to-patient, with a progressive decline in physical function and quality of life and over timeoften leading to disability and difficulty to perform many tasks of daily life¹. RA symptoms include joint pain or tenderness, joint swelling, morning stiffness, reduction in joint range of movement (ROM), muscle pain, and fatigue¹. Currently, the gold-standard methods to measure the impact of RA on daily life rely on infrequent clinical visits that may often occur every 3-4 months, with assessments depending on a combination of subjective clinician-determined scores² and patient-reported outcomes³. These have inherent limitations, however, in that they can be subjective and are prone to recall bias^{4,5}. As such, there is a need to objectively measure the impact of RA on daily life⁶, remotely over a continuous period, rather than restricting assessments to only intermittent physician visits. In recent years, consumer-grade mobile applications (app.)

However, there remains a lack of sufficient evidence for how DHT can provide objective insights into the impact of therapies for RA, despite progress made in other disease areas^{15–22}. Particularly, the benefit of sensor-outcomes generated from prescribed active assessments compared with passive monitoring has not yet been explored together. While digitised patient-reported outcomes (PROs) enable a patient the ability to regularly record their "subjective" experience of

and wearable devices have shown promise to objectively measure participants' symptoms during daily life⁷; these digital health technologies (DHT) tools⁸ have shown to increase study engagement, improve patient convenience, streamline collection of PROs⁹, and potentially generate more frequent and accurate data that can characterise disease¹⁰. DHT have been shown to measure RA symptoms and functions, such as range of motion (ROM) and gait-specific metrics during prescribed "active" assessments^{11,12}. Other studies have shown how "passive" wearable actigraphy sensor-outcome measurements capture differences in RA physical activity (PA) in daily life, compared to healthy controls (HC)¹³, as well as to detect flaring of RA symptoms¹⁴.

^{*}Corresponding author: andrew.creagh@eng.ox.ac.uk;

[†]These authors jointly supervised.



Figure 1. Illustration detailing the objectives of this study. The weaRAble-PRO 14-day trial aimed to investigate how digital health technologies (DHT)—a wrist-worn Apple smartwatch and an iPhone device, with bespoke mobile apps.—could augment patient reported outcomes (PRO) to characterise the impact of rheumatoid arthritis (RA) during the daily life of 30 moderate-to-severe RA patients, compared to 30 matched healthy controls (HC). We explore the ability of machine learning (ML) models to (1) estimate categorical RA outcomes, such as identifying RA participants from healthy controls and (2) estimate continuous RA outcomes, such as RA severity, using a combination of PRO and sensor-outcomes. Icon key: •th, triaxial accelerometer + gyroscope sensor; th, touch screen sensor; th, heart rate sensor; th, machine learning model.

disease activity in remote settings²³, it remains unclear how "objective" sensor-outcomes could provide additional insights that can augment PROs to better characterise the impact of RA on daily life. As part of this characterisation, the sensitivity of DHT to measure RA symptoms, such as the volume of remote data required and the number of sensor-outcome measurements needed, will also need to be determined. Finally, the application of DHT sensor-outcomes to monitor RA during daily life remains yet to be validated against standard in-clinic administered assessments of RA impact²⁴.

In this study, we therefore aimed to investigate how active and passive sensor-based measurements should be combined using machine learning (ML) to distinguish RA status from healthy controls, to augment traditional patient self-reported outcome (PRO) data, and to estimate standard in-clinic assessments of RA severity. Our work offers the first comprehensive evaluation of how sensor data captured during daily life can characterise RA status and severity, which represents an important first step towards the development of more sensitive and patient-centric measurements for use in RA clinical trials and real-world studies.

In order to investigate the objectives of this study, we performed the following set of analysis and experiments. We first illustrate the variety of sensor-based measurements that can be extracted from daily prescribed (active) smartphone-based assessments and (passive) smartwatch-based activity monitoring in an RA cohort. In this, we evaluate how smartwatch-based daily physical activity patterns can be remotely estimated using our bespoke deep convolutional neural, pre-trained using multi-task self-supervised learning (SSL) on a large-scale open-source cohort. We next assess the ability of our sensorbased measurements to identify RA status from healthy controls and to distinguish RA severity levels. As part of our analysis, we also explore the volume of days and number of sensor-outcomes required to remotely distinguish RA status. Finally, we investigated the power of active and passive sensoroutcomes to augment routinely collected patient self-reported outcome (PRO) data to estimate RA severity—as measured by standard in-clinic assessments of RA, such as the RAPID-3²⁵.

2 Results

The GSK weaRAble-PRO study (GSK212295) was a 14day observational study which investigated how DHT tools could objectively measure the impact of RA on participants' daily lives. Digital wearable devices—a wrist-worn Apple Watch for passive monitoring and an iPhone, integrated with a bespoke mobile app. which prescribed daily guided assessments—collected high-frequency, objective sensor data in 30 RA patients and 30 matched Healthy Controls (HCs). Figure 1 provides an illustrative overview of the objectives of this study. Sensor-based measures of physical function, mobility, dexterity, and other RA specific symptoms were extracted from daily prescribed (active) iPhone guided tests using a combination of bespoke algorithms and proprietary algorithms developed by Apple ResearchKit, for instance, a wrist-range of motion exercise, a walking assessment, a nine-



Figure 2. Ability of individual sensor-outcomes to distinguish between RA status and RA severity levels. Comparison of the average feature distributions per participants, between healthy controls (HC), RA (moderate) and RA (severe) groups for: (a–c) selection of passively collected smartwatch features; (d–f) selection of guided test collected smartphone features; and (g–i) selection of patient self-reported outcomes recorded on the smartphone application. For all examples shown, medians were significantly different between HC and RA groups: One-way ANOVA determined from the Kruskal–Wallis H-test, p<0.001. Abbreviations: deg, degrees; HAQ-DI, Health Assessment Questionnaire-Disability Index; mins, minutes; mg, mili-gravity acceleration units; MVPA, moderate-to-vigorous physical activity; RASIQ, GSK RA symptom and impact questionnaire; sed, sedentary; sec, seconds.

hole peg test, as well as two pose transition-based mobility exercises, lie-to-stand (LTS) and sit-to-stand (STS). In addition, continuous (passive) actigraphy was recorded from participants' Apple smartwatch over the study duration in order to characterise daily activity patterns and sleep. In order to illustrate the various characteristics of RA we are interested in assessing, we have grouped measurements in Fig. 1 into four domains: physical function, daytime activity, daily living, and sleep; denoting particular types of measurements which may attribute to each domain. Note: this manuscript details a substudy of weaRAble-PRO; trial design, feasibility, participant adherence, and other primary related study outcomes will be published as part of a complementary manuscript. Two RA participants withdrew immediately after enrolling in the study. Data from these participants were not collected, leaving 28 RA participants, 28 matched HCs, and 2 unmatched HCs for a total of 58 participant

2.1 Assessing smartwatch-based daily physical activity patterns

The daily physical activity of RA participants and healthy controls were estimated with a deep convolutional neural network (DCNN) that was first pre-trained on 100,000 participants in the publicly available UK Biobank, following a multi-task self-supervised learning (SSL) methodology²⁰, which was subsequently fine-tuned on the free-living Capture-24 dataset²⁷ of <150 participants to determine broad activity patterns of interest {sleep, sedentary, light physical activity, moderate-to-vigorous physical activity (MVPA)^{28,29} and fine-grained activity prediction labels {sleep, sitting/standing, mixed, vehicle, walking, bicycling $\frac{27}{1}$. In this study, we build upon our previous work by adding a temporal dependency to the DCNN (SSL) through a hidden markov model (HMM), which was appended to obtain a more accurate sequence of predicted activities over the continuous study period. It was found that the DCNN (SSL) + HMM improved broad activity estimation in Capture-24 (κ , 0.862 \pm 0.088; F1, 0.815 \pm 0.103) as compared to a baseline random forest (RF) + HMM

medRxiv preprint doi: https://doi.org/10.1101/2022.11.18.22282305; this version posted October 27, 2023. The copyright holder for this preprint (which was not certified by peer review) is the author/funder, who has granted medRxiv a license to display the preprint in perpetuity. It is made available under a CC-BY-ND 4.0 International license.



Figure 3. Ability of combined sensor-outcomes to distinguish between RA status and RA severity levels. Comparison of (a) RA identification (RA vs. HC) performance and (b) RA severity level estimation (RA (mod) vs RA (sev)), using patient reported outcomes (PRO) and combined PRO, active, and passive sensor-based outcomes in the weaRAble-PRO study. auroc: Area under the receiver operator curve; κ , Cohen's Kappa statistic; F₁, Macro-F1 score.

approach (κ , 0.813 ± 0.108; F1, 0.775 ± 0.117)²⁷. Next, the fine-tuned DCNN (SSL) + HMM model transformed the raw Apple smartwatch sensor data in weaRAble-PRO to determine participants' daily activity patterns over the 14-day study period, for example, the time spent walking, the frequency of exercise, the length and quality of sleep, and other RA-specific measures, such as morning stiffness. Activity predictions were qualitatively evaluated over the entire RA and HC study population and demonstrated excellent face validity (see section A and section B for additional details).

2.2 Analysis of sensor-outcomes to distinguish RA status and severity levels

The raw smartphone and smartwatch data recorded during the (active) guided test exercises, and passively during the participants' daily life, respectively, were summarised as sensoroutcome features. Univariate analysis demonstrated that a total of 153 (93%) sensor-based features (passive, n=131 (94%); active, n=22 (88%)) displayed significantly different medians (after post-hoc correction for multiple comparisons) between HC and RA severity groups (Kruskal-Wallis H test, p<0.05). A further 47 (34%) passive features, compared to 6 (24%) active features, were also significantly different (Mann-Whitney U test, p<0.05) between healthy and RA participants. Figure 2 compares the (fortnightly) average feature distributions between healthy controls (HC), RA (moderate) and RA (severe) participants for a selection of examples of passively collected smartwatch features (Fig. 2a-2c) and active guided test sensor features (Fig. 2d-2f) and a selection of patient selfreported outcomes recorded on the smartphone application (Fig. 2g-2i).

In order to explore the ability of many wearable sensoroutcomes to distinguish symptoms of RA from otherwise healthy individuals, and therefore measure the impact of RA during daily life, we devised a number of multivariate classification-based experiments. First, we investigated the performance of regularised logistic regression (LR) to differentiate RA participants from healthy controls using both passively collected activity monitoring features and guided test exercise features. Comparing model performance between sources (Fig. 3a), passive activity monitoring-based sensor features better distinguished RA participants using fortnightly averaged features (F1, 0.786) versus active (guided test) features (F1, 0.778). It was found that 12 subjects were misclassified using active-only models and 12 for passive-only, with just 4/12 (33%) of the same subjects incorrectly identified by both sources, 3 of which were the same HC participants. Combining active and passive wearable sensor features yielded in the highest performing models to distinguish RA participants overall, for example, using fortnightly averaged features from both sources (F1, 0.807) (for further expansion of results, see supplementary table A.I). It should be also be noted that linear logistic regression was found to perform comparatively to non-linear ensembles of decision trees, a Random Forest (RF) model and Extreme Gradient Boosted Trees (XGB)-as such this work subsequently opted to explore simple linear models for further analysis (see supplementary material, table A.II for expansion of results).

This study next investigated the ability of multiple sensorbased outcomes to augment PRO data in order to stratify RA severity levels. In weaRAble-PRO, participants were denoted as having moderate or severe RA based on baseline clinicianassessed RAPID-3 scores. Following similar procedure to RA identification, LR regularised models were investigated in order to distinguish RA (mod) and RA (sev) as binary classification tasks using fortnightly averaged study data. The benefit of incorporating additional sensor-based outcomes to patient (self-) reported outcomes is presented in Fig. 3b (expanded in supplementary table A.IV). It was observed that the linear combination of PRO assessments could accurately stratify RA symptom severity (F1, 0.759). The fusion of PRO data and sensor-based outcomes improved RA severity level estimation further with the addition of active (F1, 0.750) or passive (F1, 0.786) sources. Finally, the amalgamation of PRO



Figure 4. The number of days of sensor-data required to remotely characterise RA impact. Comparison of (a) the minimal amount of days of data needed distinguish RA status, as measured by the F1 score across 5-fold cross validation (CV), between active, passive, and combined feature sources; (b) the feature (test-retest) reliability, as measured by the intraclass correlation coefficient (ICC), between RA participants and HC across the study duration (14 days); F1 scores and ICCs suggest that model performance and feature reliability stabilises once more than 7 days of data are used per participant.

outcomes with both active and passive sensor-based outcomes resulted in the most accurate RA severity level estimation (F1, 0.833)—an improvement of 10% compared to PRO outcomes alone (Fig. 3b). For further information on the selected PRO + sensor-outcomes, we refer the reader to supplementary table A.V.

2.3 Estimating the volume of days and number of sensor-outcomes required to remotely distinguish RA status

In weaRAble-PRO, participants performed daily guided test exercises-resulting in daily sensor features-and continuously recorded Apple Watch sensor data were summarised as daily activity monitoring-based features, over the 14-day study period. In this work, we aimed to determine the minimal number of days of sensor data required build a stable and robust estimate of disease status in RA participants compared to HC over the 14-day study period. Figure 4a represents an experiment exploring the (observation-wise) out-of-sample RA classification performance as a function of varying the number of non-contiguous days of data that are averaged per participant. Evaluated over 500 randomly sampled permutations of non-contiguous days, results (median + IQR) indicated that RA prediction stabilised once more than 7 non-contiguous days of data were used per participant. Furthermore, we found that averaging daily feature values over weekly and fortnightly periods improved model performance. However, it was observed that model performance using weekly averaged features was often similar to fortnightly averaged (we refer the reader to supplementary table A.I for further details).

To investigate feature consistency and reproducibility, the intra-class correlation coefficient (ICC) for each feature was evaluated over the study duration (14 days). ICCs were calculated for each feature using n = [2, 3, ..., 14] days of data per participant, individually for HC and RA participants. Higher ICC's suggest a high degree of similarity on the performance of each task over the course of the study, and lower coefficients mean that participants tended to perform the task differently each day of the study. ICC's for HCs ranged from 0.582 to 0.854, while those for RA participants ranged from 0.424 to 0.897. Figure 4b depicts the median + inter-quartile range (IQR) of ICC values for the LR-elastic net retained active + passive features. Intra-rater reliability analyses suggest that feature reliability stabilises to good (ICC=0.75–0.9) and excellent (ICC>0.9) once more than 7 contiguous days of data were used per participant.

In order to evaluate the number of sensor-outcomes required to remotely distinguish RA status, we compared various feature regularisation techniques, lasso (ℓ_1), ridge (ℓ_2), elastic-net $(\ell_1 + \ell_2)$, and sparse-group lasso, using fortnightly (i.e., study duration) averaged features. It was found that introducing sparsity through regularisation improved classification performance. In addition, active and passively recorded sensor-based features could be grouped into domains, based on the guided test they were extracted from, or the perceived functional domain of daily activity they were assumed to assess. Introducing group-wise sparsity with the sparse-group lasso (SG-lasso), regularising on the number of groups (i.e., the feature domains) and the coefficients within each group, resulted in the highest RA participant identification performance (F1, 0.807), compared to lasso (ℓ_1 , F1, 0.772), ridge $(\ell_2, F1, 0.792)$, and elastic net $(\ell_1 + \ell_2, F1, 0.792)$ regularisation (for expansion of results, see supplementary material, table A.II). The features and groups selected by each regularisation technique are illustrated in Fig. 5, represented as



Figure 5. The number of sensor-outcomes required to remotely distinguish RA status. Comparison of features selected between regularised logistic regression (LR) models for: (a) elastic-net (F1, 0.79) and (b) SG-lasso (F1, 0.81). The SG-lasso promotes group-wise sparsity (i.e., regularising the number of feature domains) and within-group sparsity (i.e., regularising the number of features per domain), achieving a similar performance to LR elastic-net, while selecting a fewer number of domains and features. Feature importance, denoted as the mean LR coefficient value (w) over cross validation, are illustrated by colour intensity. Sensor-based feature domain abbreviations: AF: activity fragmentation; DEM: demographics; LTS: lie-to-stand assessment; MORN: morning stiffness; NTR: night-time restlessness; PEG, 9-hole peg test; STS: sit-to-stand assessment; TVDA: total volume of daytime activity; WLK: walking assessment; WRT: wrist assessment.

the mean LR coefficient value w over CV per each feature and feature domain (coefficient values have been normalised between 0-1 to benefit comparison between models). Examining the feature sparsity of elastic-net $(\ell_1 + \ell_2)$ (Fig. 5a), it was observed that features from multiple domains were selected. In contrast, the SG-lasso, as shown in Fig. 5b, selected mostly passive activity-based smartwatch features-TVDA with some morning stiffness measures-to distinguish RA status. Group sparsity penalised simultaneously selecting from multiple feature domains, where within group-sparsity regularised the feature coefficient values within the selected domains. Using fewer domains and features, the SG-lasso was able achieve similar performance to LR elastic-net, even marginally improving performance (F1, 0.807). For further details on the features extracted, and selected, we refer the reader to the supplementary material sections A and E respectively.

2.4 Estimating in-clinic RA severity scores from PRO and sensor-based outcomes

Rheumatoid arthritis severity levels were denoted by a clinician administered RAPID-3 assessment²⁵ at baseline in the weaRAble-PRO study. The RAPID-3—a "rapid" and easy to administer questionnaire—is also validated against more exhaustive assessments for RA, such as the disease activity score 28 (DAS28) and clinical disease activity index (CDAI) in clinical trials and clinical care²⁵. In this work, we aimed to establish how the combination of PRO and sensor-based outcomes could stratify continuous RAPID-3 RA severity. Note: HC subjects were assigned a RAPID-3 score of zero at baseline. Through multivariate modelling, using LR elastic-net, it was determined that PRO and sensor-based features could accurately estimate RAPID-3 scores to within 1 point (r^2 , 0.69; MAE, 0.94; RMSE, 1.33), an improvement compared to using PRO measures alone (r^2 , 0.63; MAE, 1.16; RMSE, 1.45). The association between actual and PRO + sensor-outcome estimated RAPID-3 scores was found to be good-to-excellent (r > 0.75), Pearson's r=0.60, p<0.001; Spearman's ρ =0.83, p<0.001.

Participants in weaRAble-PRO were also administered a twice-daily interactive Joint Pain Map (JMAP) questionnaire on their iPhone¹¹, in order to more precisely record and localise perceived pain. Participant model-estimated RAPID-3 scores were further interpreted through detailed inspection of the daily smartphone-based patient-reported joint pain map (JMAP) total scores-an external validation measure, which was not included as a predictor in the model-as expanded in Fig. 6. The JMAP score, defined as the sum of all individual joint pain scores per recording, was intended as a coarse measure to holistically capture participants' overall level of perceived pain, in addition to validated PRO assessments. Higher JMAP scores indicate higher levels of pain experienced. It was observed that RAPID-3 estimations were reliable and robust, in that they faithfully characterised RA participant's perceived level of symptoms, through the JMAP. For example, in Fig. 6, the RA (sev.) participant with consistently the largest reported degree of pain across the 14-day study exhibited the highest actual RAPID-3 score (6.7), which was closely estimated by the model as 7.1. JMAP scores further enabled additional explanation of model performance, especially with respect to RAPID-3 estimations that were not reflective of actual RAPID-3 scores. For instance, the RA (mod) participant with the lowest estimated RAPID-3 score (0.2) actually reported zero pain experienced over the

14-day study duration, despite a RAPID-3 assignment of 3.7 at baseline. Non-zero estimated RAPID-3 scores for some HC could also often be contextualised, due to these participants frequently self-reporting low-levels of pain in their JMAP (i.e., non-zero JMAP entries) over the study period, despite being healthy. As such, it was determined that PRO and sensor-based RAPID-3 estimates reliably reflected participant's RA symptoms over the study.

3 Discussion

Our findings in the weaRAble-PRO study demonstrate how digital health technology (DHT) captured sensor-outcomes, recorded from smartphone-based active tests, and continuously collected passive smartwatch-based monitoring, could characterise meaningful aspects of rheumatoid arthritis (RA) impairment and physical function impacting daily life. Remotely collected wearable sensor-outcomes could distinguish RA status from healthy controls-demonstrating further improved performance when combining the sensor-data from both devices-and how objective sensor-outcomes could augment patient (self-) reported outcomes to remotely estimate RA severity. Furthermore, by the half-way point of the weaRAble-PRO study (day 7), a sufficient volume of data had already been collected to reliably distinguish the characteristics of RA participants. This work provides the first comprehensive evaluation how remote and objective digital sensor-outcomes enrich our ability to understand the impact of RA on daily life between clinical visits.

In this work, we detailed how raw data collected from smartphone and smartwatch sensors can be transformed into sensor-based outcomes that are reflective of disease status. In concurrence with previous studies, many remotely collected smartphone sensor-outcomes distinguished RA participants and RA severity levels. For example, it was observed that joint ROM features differentiated HC and RA groups-a similar finding to our previous work¹²—and that RA participants were less mobile, taking longer to move between positions (as measured during the lie-to-stand exercise)-as previously shown by Andreu-Perez, et al.³⁰. Continuously collected smartwatch sensor data, known as passive monitoring, allowed the measurement of aspects of RA daily life, such as physical activity, sleep, and other RA specific symptoms, such as morning stiffness, or night-time restlessness. In this study we trained an activity recognition model on the free-living Capture-24 dataset to estimate daily activity patterns in the weaRAble-PRO population. Leveraging the latest advances in self-supervised learning (SSL) allowed our model to be pretrained on 100,000 participants with 700,000 days of diverse, unlabelled wearable sensor data in the UK Biobank²⁶, which combined with HMM temporal smoothing, significantly improved activity prediction compared to our previous established RF-HMM based methods^{27,29}. Our SSL DCNN+HMM model enabled a more robust and fine-grained estimation of daily activity patterns beyond traditional acceleration magnitude levels^{13,14}, which we proposed could allow a richer

characterisation of PA and sleep in RA. Activity monitoring revealed distinct differences distinguishing RA status, for example the daily percent of the day in moderate-to-vigorous physical activity, and similar features, were significantly lower in the RA population compared to healthy controls-a similar finding by Prioreschi, et al.¹³, and an observation people with RA regularly self-report³¹. Other specific RA symptom measurements, like morning stiffness or disrupted sleep, were evident in certain RA participants. For example, the mean acceleration value > 30 [mins] after wake-up were lower in RA—also a similar finding to Keogh, et al.³²—or that the number of movement episodes during night-time sleep distinguished some specific RA participants. We also observed that after collecting 7 days of sensor-data in the weaRAble-PRO study, a sufficient volume of data had already been recorded to reliably distinguish RA participants from a healthy population; participant feature reliability (as measured ICC values) stabilised at good-to-excellent levels, maximal identification performance of RA participants plateaued, and that there was no additional benefit to averaging over a fortnight's worth of data versus a week. Therefore it is recommended that considering at least one week's worth of sensor data is collected, it might be more beneficial to gather less data from a greater number of participants, rather than greater duration of sensor data from the same participants.

Our work is the first study to combine active smartphone and passive wearable measurements to distinguish RA status and measure variations in RA severity. While models trained on only passive features tended to marginally outperform models trained solely on active guided test features, combining both active + passive features led to the best performance in RA identification for all models investigated. Interestingly, it was found that different subjects were misclassified by active versus passive models. For example, 12 subjects were misclassified using active-only models and 12 for passiveonly, with just 4/12 (33%) of the same subjects incorrectly identified by both sources, 3 of which were the same HC participants. In addition, further experiments with the LR-SG-lasso determined that only activity monitoring domain features were mainly needed in order to distinguish RA participants from health controls. This indicates that we sometimes do not need to prescribe all guided test assessments, or to parse all activity feature domains, but that a small number of prescribed assessments can be sufficient to characterise RA status. For example, including only the lie-to-stand assessment rather than also prescribing the similar, and highly correlated, sit-to-stand assessment in future studies; or removing the prescribed walking assessment (shown to have little predictive value in the weaRAble-PRO study), and using passive daily life walking predictions generated from the activity recognition model instead, which could reduce patient burden. Finally, we also found that combining patient-reported outcomes (PRO) and objective sensor-outcomes could better capture RAPID-3-based RA severity at baseline than PROs alone; most estimated RAPID-3 scores correctly stratified



Figure 6. The ability of remote PRO + sensor-outcomes to estimate in-clinic determined RA severity scores. Scatter plot of baseline RAPID-3 scores *y* versus predicted \hat{y} scores per subject, using elastic net with PRO + sensor-outcomes, over cross-validation (CV). Participant model-estimated RAPID-3 scores can be further interpreted through detailed inspection of the daily smartphone-based patient-reported joint pain map (JMAP) total scores—which was not included as a predictor in the model. Higher JMAP scores indicate higher levels of pain experienced. Additional interpretability, through the JMAP, demonstrated that PRO + sensor-based outcome estimation of the RAPID-3 could reliably reflect patient's perceived daily RA symptoms. Note: Baseline JMAP total scores, recorded on the same day as the baseline RAPID-3, are denoted in grey; the JMAP y-axis scale is the same among all subplots. HC subjects were assigned a RAPID-3 score of zero at baseline. A black line represents perfect predictions (r^2 , 0.692; MAE, 0.938; RMSE, 1.333).

participants across severity levels from healthy to moderate to severe RA, suggesting that sufficient information to characterise RA disease severity could be reflected in the remote monitoring outcomes derived in the 14-day weaRAble-PRO study. To the best of the authors knowledge, this offers the first evaluation and insight how remote monitoring outcomes in daily life can estimate in-clinic administered assessments of RA impact.

There are a number of limitations that must be considered in the weaRAble-PRO study. Despite rich individual level measurements, the study recruited a relatively small sample size (HC, n=30; RA, n=30). As such, a degree of variability and uncertainty existed in constructing cross-validated models to distinguish RA participants, RA severity levels, or estimate the in-clinic RAPID-3 assessment. Extrapolation of results aimed at generalising RA is therefore not possible without the availability of larger cohorts and further external validation. In addition, this study only recruited RA patients with moderate-to-severe levels of disease activity; future studies should also aim to characterise patients with lower levels of disease activity or those in remission. There were also limitations associated with modelling a clinician-administered assessment, or clinical labels formulated from in-clinic assessments. For instance, the RAPID-3 was assessed at baseline, with participants recalling the prior week, yet the PRO and sensor-based features were calculated as averages over subsequent 14-day trial period from baseline. As such, the baseline RAPID-3 may not have precisely reflected the participant's

disease status recorded earlier, due to the underlying mutability and heterogeneity of RA symptoms over short periods of time. The subjectivity of PRO predictors should also considered, for instance, pain or perceived quality of sleep is relative, and some healthy participants recorded experiencing pain or affected sleep in PRO questionnaires. As a result, some PRO values influenced HC RAPID-3 predictions greater than zero, i.e., indicating the presence of RA symptoms—albeit nonzero estimated RAPID-3 predictions for HCs were generally low (<2).

The weaRAble-PRO study typifies how continuously collected patient self-reported and sensor-based outcomes may more closely reflect participant perceived and experienced symptoms that impact daily life. While in-clinic assessments are considered the gold-standard means of assessing disease severity in RA, it is clear that remotely collected, continuous, patient-centric measurements generated from PRO and sensorbased outcomes offer promising insights that can undoubtedly augment in-clinic assessments for RA. We believe that our work-the first comprehensive evaluation how remote sensor data can augment traditional PRO measures to estimate clinician-determined RA severity-helps informs future DHT study design to better characterise the impact of RA on daily life, ultimately to expand the use of DHT to develop more sensitive, and patient-centric, endpoints in RA clinical trials and real-world studies.

4 Methods

4.1 Dataset

Remotely collected smartphone and smartwatch sensor data was obtained from the GSK study title: Novel Digital Technologies for the Assessment of Objective Measures and Patient Reported Outcomes in Rheumatoid Arthritis Patients: A Pilot Study Using a Wrist-Worn Device and Bespoke Mobile App. (212295, weaRAble-PRO). This observational study followed 30 participants diagnosed with moderate-to-severe RA and 30 matched HCs over 14 days. The population demographics, in-clinic, and relevant patient self-reported outcomes, as assessed at baseline, are reported in table 1. RA participants were denoted as displaying moderate disability, RA (mod), or severe disability, RA (sev), as determined by their baseline RAPID-3 score. Note: Two RA participants withdrew immediately after enrolling in the study. Data from these participants were not collected, leaving 28 RA participants, 28 matched HCs, and 2 unmatched HCs for a total of 58 participants. Further study details, including participant requirement and data collection, are outlined in the accompanying supplementary material.

Sensor-based data collection The Apple Watch and iPhone were used to collect high frequency raw sensor data from predefined, (active) guided tests on a daily basis. Participants were prescribed daily to perform five iPhone-based assessments: WRT, a wrist range of motion (ROM) exercise¹²; WLK, a 30 second walking exercise¹²; PEG, a digital 9-hole peg test³³; STS, a sit-to-stand transition exercise^{30,34}; and LTS, a lie-to-stand transition exercise^{30,34}. For more details on the (active) guided test sensor-based features extracted, see supplementary material C. A brief overview of the guided tests prescribed in weaRAble-PRO are presented in supplementary material section C.1. In addition, the Apple Watch was used to continuously collect background sensor data (denoted passive data), as the participants went about their daily activities. Participants were asked to maintain a charge on both the Apple Watch and the iPhone, so that interruptions to monitoring and data transfer were kept to a minimum. Since night-time activity was also monitored, while participants were asleep, it was requested that charging should be done during the day, in a way that fit the participants' schedules (e.g., charging in the morning while getting ready for the day). For more details on the activity monitoring features, see supplementary material section B.5.

Patient-reported outcomes Patient-reported outcomes (PRO), most often self-report questionnaires, were administered to assess disease activity, symptoms, and health status and quality of life from the patients' perspective^{35,36}. The weaRAble-PRO study administered a selection of validated PRO measures for RA in compliment to bespoke digital PRO assessments—that are validated in clinical trials, where the questions, response options, and the general approach to assessment were standardised for all participants. PROs were recorded on days 1, 7, and 14 of data collection. The PRO

Table 1. Population demographics, in-clinic, and selected patient self-reported outcomes, as assessed at baseline, where the mean \pm standard deviation across the population are reported.

	HC ^a (n=28)	RA (mod) ^b (n=13)	RA (sev) ^c (n=15)	p^1
Demographics				
Age, years	58.4 ± 9.9	56.9 ± 11.4	60.4 ± 7.1	0.33
Female, n (%)	25 (89%)	11 (84%)	14 (93%)	0.92
BMI	25.8 ± 4.6	31.1 ± 5.9	31.7 ± 8.6	0.96
In-clinic Outcome(s	5)			
RAPID-3	0 ± 0	3.2 ± 0.7	5.3 ± 1.1	< 0.001
Patient Reported O	utcome(s)			
HAQ-DI	0 ± 0	0.63 ± 0.36	1.03 ± 0.42	< 0.01
RASIQ-pain	3.1 ± 6.7	32.1 ± 20.8	56.2 ± 11.6	< 0.01
RASIQ-stiffness	5.9 ± 9.5	33.9 ± 18.9	51.6 ± 10.2	< 0.05
RASIQ-impact	47.3 ± 5.0	53.9 ± 5.1	50.8 ± 7.6	0.33
FACIT	49.2 ± 2.9	38.9 ± 4.3	31.9 ± 7.6	< 0.05
PROMIS-sleep	49.6 ± 2.8	52.7 ± 4.2	52.4 ± 4.3	0.83
PROMIS-pain	42.2 ± 4.8	54.2 ± 7.29	58.8 ± 4.6	0.09
JMAP total pain ²	0.20 ± 0.5	13.5 ± 13.9	18.8 ± 13.7	0.23

¹ *p*, p-value calculated from Mann Whitney U-test comparing severe vs. moderate RA participants;

² Note: self-reported JMAP is not a validated PRO in RA;

^a Matched HC to RA participants only;

^b RA participants with baseline RAPID-3: 6.1–12.

^c RA participants with baseline RAPID-3: >12;

assessments administered to participants are outlined in supplementary material section D.3.

4.2 Smartwatch-based estimation of daily life patterns

In order to generate unobtrusive measures characterising physical activity and sleep in RA participants during daily life, the raw Apple Watch actigraphy (i.e., accelerometer) sensor data was transformed through a human activity recognition (HAR) sensor processing and deep convolutional neural network (DCNN) pipeline. Figure 7 illustrates how a deep convolutional neural network (DCNN) can transform raw Apple smartwatch sensor data to estimate a participant's daily activity patterns in the weaRAble-PRO study using self-supervised learning (SSL). The construction of this pipeline yielded unobtrusively measured summary features of physical activity and sleep for RA participants, computed daily during normal life.

A deep convolutional neural network (DCNN) with a ResNet-V2 architecture was first pre-trained following a multitask self-supervised learning (SSL) methodology on 100,000 participants—each participant contributing 7 days yielding roughly 700,000 person days of data—in the open-source UK biobank²⁶. The SSL pre-trained model was then fine-tuned to perform activity recognition as a downstream task in the Capture-24 dataset.

The Capture-24 study is a manually labelled, free-living dataset—that is reflective of real-world environments—and is available for training an activity recognition model to be



Figure 7. Self-supervised learning pipeline. Continuous (passive) actigraphy was recorded from patients' Apple smartwatch over the study duration. Deep convolutional neural networks (DCNN) were pre-trained on 700,000 person days in the publicly available UK Biobank using self-supervised learning—and fine-tuned with the Capture-24 dataset—to estimate participant's daily activity patterns in the weaRAble-PRO study. Physical activity (PA) metrics of daily-life, for example, the time spent walking, the frequency of exercise, or the length and quality of sleep were investigated as markers to characterise symptoms of disease in people with RA compared to HC.

applied to the weaRAble-PRO study. In Capture-24, actigraphy data was collected for 24-hours from 132 healthy volunteer participants with a Axivity AX3 wrist-worn device as they went their normal day. Activity labels provided by photographs automatically captured roughly every 30 seconds by a wearable camera for each participant. Capture-24 was labelled with 213 activity labels, standardised from the compendium of physical activities²⁸. Activity labels were then summarised into a small number of free-living behaviour labels, defining activity classes in Capture-24. There are two major labelling conventions used within Capture-24 that the model was trained to predict, defined as broad activity: {sleep, sedentary, light physical activity, moderate-to-vigorous physical activity (MVPA)^{28, 29}; and fine-grained activity: {sleep, sitting/standing, mixed, vehicle, walking, bicycling}²⁷.

HAR model predictions are essentially independent meaning that the sequence of activities over each 30 second epoch incorporates no temporal information epoch-to-epoch, for instance how the previous epoch prediction affects the current, or next, activity prediction. In order to add temporal dependency to the HAR SSL model, a Hidden Markov Model (HMM) was implemented in a post-processing step to obtain a more accurate sequence of predicted activities over the continuous 14-day data collection period as per Willetts, et al.²⁷.

The Capture-24 fine-tuned HAR SSL-HMM model was then implemented to estimate daily activities in weaRAble-PRO study data. For additional information of the HAR deep network, SSL, and other related information, we refer the reader to our previous work²⁶. Further results relating to the SSL-HMM model are outlined in the supplementary material. The sensor processing pipeline developed for the Apple Watch in the weaRAble-PRO study is outlined in supplementary Fig. B.I and within the accompanying supplementary material.

4.3 Extraction of sensor-based outcomes

Wearable sensor-based features were derived from the smartphone during the active guided tasks and passively from the smartwatch during daily life. "Active" features, extracted from smartphone sensor-based measurements during the prescribed guided tests, aimed to capture specific aspects of RA physical function, related to pain, dexterity, mobility and fatigue¹². In addition "passive" features were extracted from smartwatch sensor-based measurements, collected continuously in the background over the 14-day period. Daily activity predictions from the ML SSL model were summarised into general features measuring activity levels, period, duration and type of activity, as well as sleep detection and sleeping patterns. Furthermore, devised under the guidance of Rheumatologists, additional activity monitoring features specifically aimed at characterising well-known RA symptoms were also developed, such as morning stiffness and night-time restlessness.

The supplementary material sections B and C also detail algorithms used to extract active and passive features in the weaRAble-PRO study. For a full list of extracted sensorbased features in weaRAble-PRO, we refer the reader to the supplementary material table E.I.

4.4 Statistical Analysis

Univariate Testing Pair-wise differences groups between groups, for example HC vs. RA, or RA (mod) vs. RA (sev) were analysed for the equality in population median using the non-parametric Mann-Whitney U test (MWUT)^{37–39}. One-way analysis of variance (ANOVA) tests were also used to assess differences between medians of multiple groups, for example HC vs. RA (mod) vs. RA (sev) were assessed using the Kruskal-Wallis (KWt) test by ranks⁴⁰. The Brown-Forsythe (BF) test by (absolute deviation) of medians was used to investigate if various groups of data have been drawn with equal variances⁴¹.

Correlation Analysis Correlation analysis was utilised to determine the association or dependence between sets of random variables, such as the dependence between features, or for assessing a features' clinical utility by measuring the association to an established clinical metric. This study investigated the (linear) Pearson's *r* correlation and the (non-linear) Spearman's Rho ρ correlation between features, between features and PROs, and between clinical assessments to determine levels of association. The strengths of the correlations were classified as good-to-excellent (r>0.75), moderate-to-good (*r*=0.50–0.75), fair (*r*=0.25–0.49) or no correlation (*r*<0.25)⁴².

Feature Reliability Intra-rater (i.e., test-retest) reliability was determined using intra-class correlation coefficient (ICC) values⁴³, which were used to assesses the degree of similarity between repeated features over the course of the study for each patient. In this work, the *ICC*(3, *k*) was calculated⁴⁴–which considers the two-way random average measures with *k* repeated measurements—for the 14-day session across subjects, where the raters *k* are the study days. Reliability was categorised as either poor (ICC<0.5), moderate (ICC=0.5–0.75), good (ICC=0.75–0.9) or excellent (ICC>0.9)⁴⁵.

Correcting for multiple hypothesis testing Multiple hypothesis testing was performed due to the large volume of features by controlling the false discovery rate (FDR) at level α using the linear step-up procedure introduced by Benjamini and Hochberg (BH)^{46,47}.

4.5 Machine-learning estimation of RA status and severity

This work explored how state-of-the art machine learning (ML) models characterise the impact of RA during the daily life of participants in over the 14-day weaRAble-PRO study. Multivariate modelling aimed to explore the ability of active, passive, and PRO measures to (1) distinguish RA participants from healthy controls (HC), and (2) to estimate RA disease severity: between RA participants with moderate symptoms

(RA mod) and severe symptoms (RA sev) as binary classification tasks. Expansions of this analysis subsequently investigate how the in-clinic RAPID-3 assessment, a continuous measure of RA severity, could be estimated from the combination of PRO and sensor-based outcomes.

Overview of models: This analysis compared both linear and non-linear ML models to transform PRO and sensor-based outcomes to capture RA status and severity. Regularised linear regression (LR) models, with combinations of ℓ_1 and ℓ_2 priors, such as LR-lasso (ℓ_1), LR-ridge (ℓ_2), and LR-elasticnet ($\ell_1 + \ell_2$) were compared to yield predictive, yet sparse model solutions⁴⁸. Further regularisation extensions were also investigated using the sparse-group lasso (SG-lasso)–an extension of the lasso that promotes both group sparsity and within group parameter-wise (ℓ_2) sparsity, through a group lasso penalty and the lasso penalty—which aims to yield a sparse set of groups and also a sparse set of covariates in each selected group^{49,50}

Linear regression regularised models were also compared to decision tree (DT) based non-linear models, for instance the off-the-shelf Random Forest (RF)⁵¹ and Extreme Gradient Boosted Trees (XGB)⁵². Both LR- and DT-based models can intrinsically perform regression or classification depending on the task required. In the LR case, classification is denoted as logistic regression (though a logit-link function). NOTE: in this analysis LR can refer to both linear regression for continuous outputs or logistic regression for classification outputs. In the DT case, the mean prediction of the individual trees creates a continuous output for regression. For further details on the models employed in this study, we refer the reader to the supplementary material section F.2.

Model evaluation: To determine the generalisability of our models, a stratified leave-k-subjects out cross-validation (CV) was employed. This consisted of randomly partitioning the dataset into folds with k=5 subjects in each fold, which was stratified with equal class proportions where possible. Subject data remained independent between training, validation, and testing sets. One set was denoted the training set (in-sample), and the remaining 20% of the dataset was then denoted testing set (out-of-sample) on which predictions were made.

Feature-wise and prediction-wise aggregation: In this work, we experimented with feature-wise and prediction-wise aggregation. In feature-wise aggregation, features were computed either as: daily feature values over the 14-day study period; the average daily feature value over a 7-day period (weekly); the average daily feature value over a 14-day period (fortnightly). Predictions could then be evaluated for each day (denoted *observation-wise*) or aggregated over all days through majority voting each individual prediction per subject (denoted *subject-wise*). For example, daily and weekly averaged features result in daily, or weekly predictions (i.e., *observation-wise*), which were summarised into *subject-wise* outcomes by majority voting over the repeated predictions.

Evaluation metrics: Multi-class classification metrics were reported as the *observation-wise* median and interquartile (IQR) range over one CV, as well as the *subject-wise* outcome for that CV, using: auroc, area under the receiver operating characteristic curve; k, Cohen's kappa statistic^{53,54}; F₁, F1-score. The coefficient of determination, r^2 , the mean absolute error (MAE), and root-mean squared error (RMSE) were used to evaluate modelling the (continuous) in-clinic RAPID-3 scores.⁵⁵.

Acknowledgements

The authors are grateful to all the study participants and their families for their time and dedication to this study. The authors would also like to thank Priyanka Bobbili PhD, Julien Bendelac BSc, Jessica Landry MSc, Maral DerSarkissian PhD, Mihran Yenikomshian MBA, and Med Kouaici (MEng) from Analysis Group (MA, USA) for their support in app. design & development and data collection, and to Elinor Mody from Reliant Medical Group (MA, USA) for patient recruitment. The weaRAble-PRO study was funded and sponsored by GSK Plc. The research described in this paper was funded by GSK Plc. Authors conducting this research also acknowledge support by the National Institute for Health Research (NIHR) Oxford Biomedical Research Centre (BRC). Aiden Doherty is supported by the Wellcome Trust [223100/Z/21/Z].

Competing Interests

A.P.C, H.Y, G.M, A.D, D.A.C are employees of the University of Oxford. A.P.C is a GSK postdoctoral fellow and acknowledges the support of GSK. D.A.C received research funding from GSK to conduct this work. In addition, A.D., H.Y., and G.M. acknowledge the support of Novo Nordisk plc. A.D. AD is supported by the Wellcome Trust [223100/Z/21/Z]. V.H, W-H.C, R.T, R.W and L.G-G are employees of GSK and own stock and or shares. C.L, C.Y, M.S.D are employees of Analysis Group, which received research funding from GSK to conduct the study.

Data Availability Statement

Anonymised individual participant data that support the findings of this study are available from the corresponding author, upon reasonable request and subject to GSK's approval.

Code Availability

Apple Watch sensor processing was performed using a bespoke version of the biobankAccelerometerAnalysis toolkit, found at: https://github.com/OxWearables/biobankAccelerometerAnalysis. Deep networks were built using Python v3.7 through a Py-Torch v1.7 framework. Our self-supervised learning activity prediction code and trained models are publicly available at: https://github.com/OxWearables/ssl-wearables, including pre-trained models on 100K participants in the UK Biobank. Some guided test exercises and health metrics calculated are proprietary to

Apple ResearchKit (http://researchkit.org/) and Apple HealthKit (https://developer.apple.com/documentation/healthkit) which we refer the reader for more details. Statistical and machine learning analysis was developed using scikit-learn v1.1.1. Further analysis code can be made available from the corresponding author upon reasonable request.

Ethics Statement

All study information, informed consent, study questions and instructions for conducting the guided tests were first drafted in the form of a survey instrument. The survey instrument was then programmed into the mobile app. All documentation including the study protocol, any amendments, and informed consent procedures, were reviewed and approved by Reliant Medical Group's IRB. All participants provided written informed consent before any study procedures were undertaken. The study was conducted in accordance with the International Committee for Harmonisation principles of Good Clinical Practice and the Declaration of Helsinki.

Author Contributions

A.P.C conceptualised the data analysis, designed methodology, software and interpretation. V.H, H.Y, and G.M contributed software application for analysis. V.H, R.T, W-H.C, R.W, L.G-G contributed to the design of the study and towards the data analysis and interpretation. C.L, C.Y, M.S.D were involved in the design of the study, data collection, and software for data acquisition. A.D, L.G-G, D.A.C jointly supervised. A.P.C wrote the manuscript; all other authors: review & editing.

References

- 1. Grassi, W., De Angelis, R., Lamanna, G. & Cervini, C. The clinical features of rheumatoid arthritis. *Eur. journal radiology* 27, S18–S24 (1998).
- Banderas, B., Skup, M., Shields, A. L., Mazar, I. & Ganguli, A. Development of the rheumatoid arthritis symptom questionnaire (rasq): a patient reported outcome scale for measuring symptoms of rheumatoid arthritis. *Curr. Med. Res. Opin.* 33, 1643–1651 (2017).
- **3.** Lubeck, D. P. Patient-reported outcomes and their role in the assessment of rheumatoid arthritis. *Pharmacoeconomics* **22**, 27–38 (2004).
- **4.** Campbell, R., Ju, A., King, M. T. & Rutherford, C. Perceived benefits and limitations of using patient-reported outcome measures in clinical practice with individual patients: a systematic review of qualitative studies. *Qual. Life Res.* 1–24 (2021).
- **5.** Gossec, L., Dougados, M. & Dixon, W. Patient-reported outcomes as end points in clinical trials in rheumatoid arthritis. *RMD open* **1**, e000019 (2015).

- **6.** Flurey, C. A., Morris, M., Richards, P., Hughes, R. & Hewlett, S. It's like a juggling act: rheumatoid arthritis patient perspectives on daily life and flare while on current treatment regimes. *Rheumatology* **53**, 696–703 (2014).
- 7. Piga, M., Cangemi, I., Mathieu, A. & Cauli, A. Telemedicine for patients with rheumatic diseases: systematic review and proposal for research agenda. In *Seminars in Arthritis and Rheumatism*, vol. 47, 121–128 (Elsevier, 2017).
- **8.** Taylor, K. I., Staunton, H., Lipsmeier, F., Nobbs, D. & Lindemann, M. Outcome measures based on digital health technology sensor data: data-and patient-centric approaches. *NPJ digital medicine* **3**, 1–8 (2020).
- **9.** Yun, H. *et al.* Assessing rheumatoid arthritis disease activity with patient-reported outcomes measurement information system measures using digital technology. *Arthritis Care & Res.* **72**, 553–560 (2020).
- Munos, B. *et al.* Mobile health: the power of wearables, sensors, and apps to transform clinical trials. *Annals New York Acad. Sci.* 1375, 3–18 (2016).
- **11.** Crouthamel, M. *et al.* Using a researchkit smartphone app to collect rheumatoid arthritis symptoms from real-world participants: feasibility study. *JMIR mHealth uHealth* **6**, e9656 (2018).
- **12.** Hamy, V. *et al.* Developing smartphone-based objective assessments of physical function in rheumatoid arthritis patients: the PARADE study. *Digit. biomarkers* **4**, 26–44 (2020).
- **13.** Prioreschi, A., Hodkinson, B., Avidon, I., Tikly, M. & McVeigh, J. A. The clinical utility of accelerometry in patients with rheumatoid arthritis. *Rheumatology* **52**, 1721–1727 (2013).
- **14.** Gossec, L. *et al.* Detection of flares by decrease in physical activity, collected using wearable activity trackers in rheumatoid arthritis or axial spondyloarthritis: an application of machine learning analyses in rheumatology. *Arthritis care & research* **71**, 1336–1343 (2019).
- **15.** Pourahmadi, M. R. *et al.* Reliability and concurrent validity of a new iphone[®] goniometric application for measuring active wrist range of motion: a cross-sectional study in asymptomatic subjects. *J. anatomy* **230**, 484–495 (2017).
- **16.** Pratap, A. *et al.* Evaluating the utility of smartphonebased sensor assessments in persons with multiple sclerosis in the real-world using an app (elevateMS): observational, prospective pilot digital health study. *JMIR mHealth uHealth* **8**, e22108 (2020).
- **17.** Webster, D. E. *et al.* Clinical validation of digital biomarkers and machine learning models for remote measurement of psoriasis and psoriatic arthritis. *medRxiv* (2022).

- Omberg, L. *et al.* Remote smartphone monitoring of parkinson's disease and individual response to therapy. *Nat. Biotechnol.* 40, 480–487 (2022).
- Creagh, A. P. *et al.* Smartphone- and smartwatch-based remote characterisation of ambulation in multiple sclerosis during the two-minute walk test. *IEEE J. Biomed. Heal. Informatics* 25, 838–849, DOI: 10.1109/JBHI.2020. 2998187 (2021).
- **20.** Creagh, A. *et al.* Smartphone-based remote assessment of upper extremity function for multiple sclerosis using the draw a shape test. *Physiol. measurement* **41**, 054002 (2020).
- **21.** Lipsmeier, F. *et al.* Reliability and validity of the Roche PD mobile application for remote monitoring of early parkinson's disease. *Sci. Reports* **12**, 1–15 (2022).
- **22.** Lipsmeier, F. *et al.* A remote digital monitoring platform to assess cognitive and motor symptoms in huntington disease: Cross-sectional validation study. *J. Med. Internet Res.* **24**, e32997 (2022).
- **23.** El Miedany, Y. *et al.* Toward electronic health recording: evaluation of electronic patient-reported outcome measures system for remote monitoring of early rheumatoid arthritis. *The J. rheumatology* **43**, 2106–2112 (2016).
- 24. Coravos, A., Khozin, S. & Mandl, K. D. Developing and adopting safe and effective digital biomarkers to improve patient outcomes. *NPJ digital medicine* 2, 1–5 (2019).
- 25. Pincus, T., Yazici, Y. & Bergman, M. J. Rapid3, an index to assess and monitor patients with rheumatoid arthritis, without formal joint counts: similar results to das28 and cdai in clinical trials and clinical care. *Rheum. Dis. Clin.* 35, 773–778 (2009).
- **26.** Yuan, H. *et al.* Self-supervised learning for human activity recognition using 700,000 person-days of wearable data. *arXiv preprint arXiv:2206.02909* (2022).
- **27.** Willetts, M., Hollowell, S., Aslett, L., Holmes, C. & Doherty, A. Statistical machine learning of sleep and physical activity phenotypes from sensor data in 96,220 uk biobank participants. *Sci. reports* **8**, 1–10 (2018).
- Ainsworth, B. E. *et al.* 2011 compendium of physical activities: a second update of codes and met values. *Med Sci Sports Exerc.* 43, 1575–1581 (2011).
- **29.** Walmsley, R. *et al.* Reallocating time from devicemeasured sleep, sedentary behaviour or light physical activity to moderate-to-vigorous physical activity is associated with lower cardiovascular disease risk. *MedRxiv* (2020).
- **30.** Andreu-Perez, J. *et al.* Developing fine-grained actigraphies for rheumatoid arthritis patients from a single accelerometer using machine learning. *Sensors* **17**, 2113 (2017).

- **31.** Sokka, T. *et al.* Physical inactivity in patients with rheumatoid arthritis: data from twenty-one countries in a cross-sectional, international study. *Arthritis Care & Res. Off. J. Am. Coll. Rheumatol.* **59**, 42–50 (2008).
- **32.** Keogh, A. *et al.* A thorough examination of morning activity patterns in adults with arthritis and healthy controls using actigraphy data. *Digit. Biomarkers* **4**, 78–88 (2020).
- **33.** Mathiowetz, V., Weber, K., Kashman, N. & Volland, G. Adult norms for the nine hole peg test of finger dexterity. *The Occup. Ther. J. Res.* **5**, 24–38 (1985).
- Bohannon, R. W. Sit-to-stand test for measuring performance of lower extremity muscles. *Percept. motor skills* 80, 163–166 (1995).
- **35.** of Health, U. D. *et al.* Guidance for industry: patient-reported outcome measures: use in medical product development to support labeling claims: draft guidance. *Heal. Qual. Life Outcomes* **4**, 79 (2006).
- 36. Mercieca-Bebber, R., King, M. T., Calvert, M. J., Stockler, M. R. & Friedlander, M. The importance of patient-reported outcomes in clinical trials and strategies for future optimization. *Patient related outcome measures* 9, 353 (2018).
- **37.** Mann, H. B. & Whitney, D. R. On a test of whether one of two random variables is stochastically larger than the other. *The annals mathematical statistics* 50–60 (1947).
- Hollander, M., Wolfe, D. A. & Chicken, E. Nonparametric statistical methods, vol. 751 (John Wiley & Sons, 2013).
- **39.** Gibbons, J. D. & Chakraborti, S. *Nonparametric Statistical Inference: Revised and Expanded* (CRC press, 2014).
- Kruskal, W. H. & Wallis, W. A. Use of ranks in onecriterion variance analysis. J. Am. statistical Assoc. 47, 583–621 (1952).
- **41.** Brown, M. B. & Forsythe, A. B. Robust tests for the equality of variances. *J. Am. Stat. Assoc.* **69**, 364–367 (1974).
- **42.** Portney, L. G., Watkins, M. P. *et al. Foundations of clinical research: applications to practice*, vol. 892 (Pearson/Prentice Hall Upper Saddle River, NJ, 2009).
- **43.** Weir, J. P. Quantifying test-retest reliability using the intraclass correlation coefficient and the sem. *The J. Strength* & *Cond. Res.* **19**, 231–240 (2005).
- **44.** Shrout, P. E. & Fleiss, J. L. Intraclass correlations: uses in assessing rater reliability. *Psychol. bulletin* **86**, 420 (1979).
- **45.** Koo, T. K. & Li, M. Y. A guideline of selecting and reporting intraclass correlation coefficients for reliability research. *J. chiropractic medicine* **15**, 155–163 (2016).

- **46.** Shaffer, J. P. Multiple hypothesis testing. *Annu. review psychology* **46**, 561–584 (1995).
- Benjamini, Y. & Hochberg, Y. Controlling the false discovery rate: a practical and powerful approach to multiple testing. *J. Royal statistical society: series B (Methodolog-ical)* 57, 289–300 (1995).
- **48.** Hastie, T., Tibshirani, R. & Friedman, J. *The elements of statistical learning: data mining, inference, and prediction* (Springer Science & Business Media, 2009).
- **49.** Friedman, J., Hastie, T. & Tibshirani, R. A note on the group lasso and a sparse group lasso. *arXiv preprint arXiv:1001.0736* (2010).
- Simon, N., Friedman, J., Hastie, T. & Tibshirani, R. A sparse-group lasso. *J. computational graphical statistics* 22, 231–245 (2013).
- **51.** Breiman, L. Random forests. *Mach. learning* **45**, 5–32 (2001).
- **52.** Chen, T. & Guestrin, C. Xgboost: A scalable tree boosting system. In *Proceedings of the 22nd acm sigkdd international conference on knowledge discovery and data mining*, 785–794 (2016).
- **53.** He, H. & Garcia, E. A. Learning from imbalanced data. *IEEE Transactions on knowledge data engineering* **21**, 1263–1284 (2009).
- **54.** Cohen, J. A coefficient of agreement for nominal scales. *Educ. psychological measurement* **20**, 37–46 (1960).
- **55.** Rao, C. R. *Linear statistical inference and its applications*, vol. 2 (Wiley New York, 1973).
- **56.** Beckerman, H. *et al.* Smallest real difference, a link between reproducibility and responsiveness. *Qual. Life Res.* **10**, 571–578 (2001).
- **57.** He, K., Zhang, X., Ren, S. & Sun, J. Identity mappings in deep residual networks. In *European conference on computer vision*, 630–645 (Springer, 2016).
- Henchoz, Y. *et al.* Physical activity and energy expenditure in rheumatoid arthritis patients and matched controls. *Rheumatology* 51, 1500–1507 (2012).
- **59.** Harkness, J. *et al.* Circadian variation in disease activity in rheumatoid arthritis. *Br Med J (Clin Res Ed)* **284**, 551–554 (1982).
- **60.** Earhart, G. M. *et al.* The 9-hole peg test of upper extremity function: average values, test-retest reliability, and factors contributing to performance in people with parkinson disease. *J. Neurol. Phys. Ther.* **35**, 157–163 (2011).
- **61.** Zijlstra, W. & Hof, A. L. Displacement of the pelvis during human walking: experimental data and model predictions. *Gait & posture* **6**, 249–262 (1997).
- **62.** Godfrey, A., Del Din, S., Barry, G., Mathers, J. & Rochester, L. Instrumenting gait with an accelerometer:

a system and algorithm examination. *Med. engineering* & *physics* **37**, 400–407 (2015).

- **63.** Zhao, H. *et al.* Smartphone-based 3d indoor pedestrian positioning through multi-modal data fusion. *Sensors* **19**, 4554 (2019).
- 64. Schwickert, L. *et al.* Inertial sensor based analysis of lie-to-stand transfers in younger and older adults. *Sensors* 16, 1277 (2016).
- **65.** Becker, B. *et al.* Development, psychometric evaluation and cognitive debriefing of the rheumatoid arthritis symptom and impact questionnaire (rasiq). *J. Patient-Reported Outcomes* **5**, 1–15 (2021).
- **66.** Amtmann, D. *et al.* Development of a promis item bank to measure pain interference. *Pain* **150**, 173–182 (2010).
- **67.** Buysse, D. J. *et al.* Development and validation of patient-reported outcome measures for sleep disturbance and sleep-related impairments. *Sleep* **33**, 781–792 (2010).
- **68.** Maska, L., Anderson, J. & Michaud, K. Measures of functional status and quality of life in rheumatoid arthritis: health assessment questionnaire disability index (haq), modified health assessment questionnaire (mhaq), multidimensional health assessment questionnaire (mdhaq), health assessment questionnaire ii (haq-ii), improved health assessment questionnaire (improved haq), and rheumatoid arthritis quality of life (raqol). *Arthritis care* & *research* **63**, S4–S13 (2011).
- **69.** Hewlett, S., Dures, E. & Almeida, C. Measures of fatigue: Bristol rheumatoid arthritis fatigue multi-dimensional questionnaire (braf mdq), bristol rheumatoid arthritis fatigue numerical rating scales (braf nrs) for severity, effect, and coping, chalder fatigue questionnaire (cfq), checklist individual strength (cis20r and cis8r), fatigue severity scale (fss), functional assessment chronic illness therapy (fatigue)(facit-f), multi-dimensional assessment of fatigue (maf), multi-dimensional fatigue inventory (mfi), pediatric quality of life (pedsql) multi-dimensional fatigue scale, profile of fatigue (prof), short form 36 vitality subscale (sf-36 vt), and visual analog scales (vas). *Arthritis care & research* **63**, S263–S286 (2011).
- **70.** Cella, D. *et al.* Validation of the functional assessment of chronic illness therapy fatigue scale relative to other instrumentation in patients with rheumatoid arthritis. *The J. rheumatology* **32**, 811–819 (2005).
- **71.** Hewlett, S., Hehir, M. & Kirwan, J. R. Measuring fatigue in rheumatoid arthritis: a systematic review of scales in use. *Arthritis Care & Res.* **57**, 429–439 (2007).
- 72. Ware Jr, J. E. Sf-36 health survey update. *Spine* 25, 3130–3139 (2000).
- **73.** ten Klooster, P. M. *et al.* Performance of the dutch sf-36 version 2 as a measure of health-related quality of life in patients with rheumatoid arthritis. *Heal. quality life outcomes* **11**, 1–9 (2013).

- **74.** Matcham, F. *et al.* The impact of rheumatoid arthritis on quality-of-life assessed using the sf-36: a systematic review and meta-analysis. In *Seminars in arthritis and rheumatism*, vol. 44, 123–130 (Elsevier, 2014).
- **75.** Langley, G. & Sheppeard, H. The visual analogue scale: its use in pain measurement. *Rheumatol. international* **5**, 145–148 (1985).
- **76.** Di, J. *et al.* Patterns of sedentary and active time accumulation are associated with mortality in us adults: The nhanes study. *bioRxiv* 182337 (2017).
- 77. Little, R. J. & Rubin, D. B. Statistical analysis with missing data. john wiley & sons. *New York* (2002).
- **78.** Sakia, R. M. The box-cox transformation technique: a review. *J. Royal Stat. Soc. Ser. D (The Stat.* **41**, 169–178 (1992).
- **79.** McCullagh, P. *Generalized linear models* (Routledge, 2018).
- **80.** Hosmer Jr, D. W., Lemeshow, S. & Sturdivant, R. X. *Applied logistic regression*, vol. 398 (John Wiley & Sons, 2013).
- **81.** Hastie, T., Tibshirani, R. & Wainwright, M. *Statistical learning with sparsity: the lasso and generalizations* (CRC press, 2015).
- **82.** Tibshirani, R. Regression shrinkage and selection via the lasso. *J. Royal Stat. Soc. Ser. B (Methodological)* 267–288 (1996).
- **83.** Strobl, C., Boulesteix, A.-L., Kneib, T., Augustin, T. & Zeileis, A. Conditional variable importance for random forests. *BMC bioinformatics* **9**, 307 (2008).
- Genuer, R., Poggi, J.-M. & Tuleau-Malot, C. Variable selection using random forests. *Pattern recognition letters* 31, 2225–2236 (2010).

A Supplementary Results

A.1 Extended Results: Assessing smartwatch-based daily physical activity patterns

Figure A.Ia–A.Ib summarises the population-wide daily activity by time of the day for HC and RA groups. The probability of an activity being performed at a specific time can be computed as the number of instances detected for that activity across all data points (participants and days) at that time divided by the total number of data points. Representative examples of the predicted daily activity patterns for an individual healthy control (HC) and RA (moderate) participant are depicted in Fig. A.Ic–A.Id respectively. The times when the Apple Watch was left to charge can be clearly seen in each example, indicated by the white non-wear times, typically occurring after wake-up or before bedtime. Both participants demonstrated consistent wake-up and bed times, day-to-day—which the activity prediction model tended to correctly identify.



Figure A.I. Assessing smartwatch-based daily physical activity patterns. Variation in the average predicted daily-activity (probability) over time for all (a) HC participants and (b) RA participants in the 14-day weaRAble-PRO study. Predicted daily activity patterns for an individual (c) healthy control (HC) participant, (Female, 66 yrs.) and (d) a moderate Rheumatoid Arthritis (RA mod) participant (Female, 50 yrs.; RAPID-3, 3.7). Moving average acceleration values are overlaid in black. Participant self-reported sleep / wake times are indicated with long-dashed black lines. Non-wear times (expected daily for watch charging) are indicated by white areas. Note: the acceleration y-axis scaling between (c) and (d) is not the same due to difference in the magnitude of acceleration between participants. MVPA, moderate-to-vigorous physical activity.

A.2 Extended Results: Multivariate feature analysis

The relationship between the wearable sensor-based features extracted in this study, for both active (smartphone) and passive (smartwatch) data sources were investigated using pair-wise Spearman's ρ correlation. Correlation analysis indicated good-to-excellent relationships ($\rho > 0.75$) between many features within feature domains (intra-source); for example, most TVDA features were highly correlated with each other (positively and negatively). Analysis also revealed good-to-excellent ($\rho > 0.75$) correlation between domains of features sources (inter-source); for example, TVDA features were not only correlated with each other, but with other passive feature domains, such as AF or MORN features. Much of the inter-source correlation was between similar domains, such as within the activity monitoring-based feature domains, or within the guided test (active) feature domains—suggesting a high degree of multicollinearity and redundancy. However, mostly fair correlations ($\rho=0.25-0.50$) between active and passively extracted sensor features suggested that different information may be learned during activity monitoring versus guided test exercises. The resulting correlation matrix is depicted in Fig. A.II.



Figure A.II. Assessing correlation and collinearity between PRO and sensor-based features. Pairwise Spearman's ρ correlation matrix for PRO, active, and passive features, labelled by feature domain. Feature association is bounded between +1/-1 denoting positive and negative correlation. Feature domain abbreviations: FACIT: Functional Assessment of Chronic Illness Therapy Fatigue; HAQ-DI: Health Assessment Questionnaire-Disability Index; PROMIS: Patient-Reported Outcomes Measurement Information System; RASIQ: GSK RA symptom and impact questionnaire; SF-36: Short-Form 36 (SF-36); AF: activity fragmentation; MORN: morning stiffness; NTR: night-time restlessness; TVDA: total volume of daytime activity; LTS: lie-to-stand assessment; PEG: 9-hole peg test; STS: sit-to-stand assessment; WRT: wrist assessment; WLK: walking assessment; DEM: demographics.

A.3 Extended Results: Distinguishing RA participants

Table A.I. Comparison of RA vs. HC classification performance across different source and feature combinations with 5-fold cross-validation (CV). Results are presented as: (1) the posterior overall *subject-wise* outcome for one cross-validation (CV) run as well as (2) the *observation-wise* median and inter-quartile range (IQR) across that CV in brackets. The best performing model for each source combination are highlighted in **bold**. auroc: Area under the receiver operator curve; κ , Cohen's Kappa statistic; F₁, Macro-F1 score.

source	feature ¹	model	auroc	К	F ₁
	daily	LR-elastic-net LR-SG-lasso	0.673 (0.681, 0.616–0.748) 0.817 (0.725, 0.676–0.930)	0.183 (0.180, 0.090–0.421) 0.542 (0.433, 0.258–0.708)	0.531 (0.648, 0.464–0.676) 0.755 (0.680, 0.610–0.831)
active	weekly	LR-elastic-net LR-SG-lasso	0.699 (0.757, 0.629–0.792) 0.802 (0.771, 0.729–0.825)	0.362 (0.374, 0.118–0.459) 0.614 (0.471, 0.408–0.545)	0.640 (0.667, 0.444–0.700) 0.800 (0.667, 0.667–0.762)
	fortnightly	LR-elastic-net LR-SG-lasso	0.779 (0.743, 0.686–0.800) 0.795 (0.743, 0.714–0.833)	0.333 (0.267, 0.098–0.633) 0.578 (0.500, 0.471–0.633)	0.655 (0.667, 0.545–0.800) 0.778 (0.727, 0.667–0.800)
	daily	LR-elastic-net LR-SG-lasso	0.821 (0.691, 0.680–0.701) 0.852 (0.712, 0.677–0.758)	0.439 (0.347, 0.323–0.362) 0.609 (0.369, 0.335–0.380)	0.714 (0.689, 0.617–0.691) 0.776 (0.691, 0.667–0.698)
passive	weekly	LR-elastic-net LR-SG-lasso	0.837 (0.795, 0.722–0.859) 0.812 (0.833, 0.804–0.906)	0.579 (0.441, 0.438–0.507) 0.612 (0.588, 0.571–0.607)	0.786 (0.706, 0.667–0.733) 0.795 (0.750, 0.714–0.828)
	fortnightly	LR-elastic-net LR-SG-lasso	0.790 (0.800, 0.657–0.943) 0.831 (0.867, 0.657–0.943)	0.615 (0.500, 0.500–0.814) 0.579 (0.500, 0.314–0.814)	0.807 (0.727, 0.727–0.889) 0.786 (0.727, 0.714–0.889)
	daily	LR-elastic-net LR-SG-lasso	0.814 (0.729, 0.704–0.891) 0.837 (0.744, 0.704–0.908)	0.473 (0.312, 0.196–0.552) 0.504 (0.317, 0.275–0.625)	0.727 (0.686, 0.611–0.767) 0.720 (0.683, 0.605–0.808)
active + passive	weekly	LR-elastic-net LR-SG-lasso	0.848 (0.842, 0.786–0.850) 0.848 (0.833, 0.700–0.893)	0.614 (0.538, 0.486–0.577) 0.614 (0.589, 0.254–0.814)	0.800 (0.769, 0.762–0.800) 0.800 (0.783, 0.692–0.889)
	fortnightly	LR-elastic-net LR-SG-lasso	0.857 (0.914, 0.867–0.914) 0.842 (0.867, 0.714–0.943)	0.579 (0.471, 0.351–0.814) 0.615 (0.676, 0.500–0.814)	0.786 (0.714, 0.667–0.889) 0.807 (0.833, 0.727–0.889)

¹ daily: daily feature values over the 14-day study period; weekly: the average daily feature value over a 7-day period; fortnightly: the average daily feature value over a 14-day period;

Table A.II. Comparison of RA vs. HC classification performance for logistic regression (LR) based models and decision trees (DT) across with 5-fold cross-validation (CV) with fortnightly (i.e., study duration) averaged active + passive features. Results are presented as: (1) the posterior overall *subject-wise* outcome for one cross-validation (CV) run as well as (2) the *observation-wise* median and inter-quartile range (IQR) across that CV in brackets. The best performing model for each feature representation are highlighted in **bold**. auroc: Area under the receiver operator curve; κ , Cohen's Kappa statistic; F₁, Macro-F1 score.

	model	auroc	К	fl
LR	-	0.853 (0.867, 0.800–0.943)	0.542 (0.500, 0.459–0.814)	0.755 (0.727, 0.727–0.889)
	lasso	0.788 (0.857, 0.667–0.867)	0.545 (0.500, 0.241–0.814)	0.772 (0.727, 0.714–0.889)
	ridge	0.777 (0.833, 0.767–0.914)	0.612 (0.657, 0.441–0.814)	0.792 (0.800, 0.769–0.889)
	elastic-net	0.801 (0.867, 0.667–0.914)	0.612 (0.657, 0.441–0.814)	0.792 (0.800, 0.769–0.889)
	SG-lasso	0.842 (0.867, 0.714–0.943)	0.615 (0.676, 0.500–0.814)	0.807 (0.833, 0.727–0.889)
DT	RF	0.862 (0.933, 0.829–0.957)	0.646 (0.657, 0.471–0.814)	0.800 (0.800, 0.727–0.889)
	XGB	0.851 (0.833, 0.829–0.914)	0.614 (0.676, 0.471–0.814)	0.800 (0.833, 0.714–0.889)

Table A.III represents a selection of features that were retained by LR-elastic-net. The model tended to pick features from all domains, but consistently tended to select many different features between cross validation splits. Some features however, for example, the mean transition time [sec] from standing to lying, daily average total time in MVPA bouts [mins], the average hazard of non-MVPA to MVPA bouts were constantly chosen over all data splits. Other features were selected less often but, when chosen, weighted highly in the model, for instance: median ROM [deg] or the number of movement episodes during night-time sleep [count/hr].

In order to determine the test-retest reliability of the selected features, we calculated intra-class correlation coefficient (ICC) values⁴³, which were used to assesses the degree of similarity between repeated features over the course of the study for each patient. Here, we calculated $ICC(3,k)^{44}$ -which considers the two-way random average measures with *k* repeated measurements—for the 14-day session across subjects, where the raters *k* are the study days.

The minimal detectable change, with a 95% confidence interval (CI) (MDC₉₅), was also calculated to determine the minimal change in a feature which is greater than the within subject variability and measurement error, indicating how much a measured change is likely to reflect true change from repeated measurement. First, the standard error of measurement (SEM), which provides an absolute index of precision⁴² was calculated:

$$SEM = SD \times \sqrt{(1 - ICC)} \tag{A.1}$$

where *SD* and *ICC* are the variance and intra-class correlation coefficients of the feature, *x*. Next, the minimal detectable change, with a 95% confidence interval (CI) (MDC₉₅) was determined^{42,56}:

$$MDC_{95} = 1.96 \times SEM \times \sqrt{2} \tag{A.2}$$

where MDC_{95} was expressed as percentages that are independent of the units of measurement for each feature:

$$MDC_{95}\% = \frac{MDC_{95}}{\bar{x}}$$
 (A.3)

using the respective mean feature value, \bar{x} .

Table A.III. Selection of top performing active (smartphone) and passive (smartwatch) extracted features for RA identification, as determined by logistic regression (LR) elastic-net across 5-fold subject-wise cross validation (CV), with fortnightly (i.e., study duration) averaged features. Features were ranked per CV fold by increasing shrinkage regularisation parameter λ and recording the percentage (%) of time that feature is selected in the subset that minimises the CV error in the validation set. Feature domain abbreviations: AF: activity fragmentation; DEM: demographics; LTS: lie-to-stand assessment; MORN: morning stiffness; NTR: night-time restlessness; STS: sit-to-stand assessment; TVDA: total volume of daytime activity; WLK: walking assessment; WRT: wrist assessment.

	\mathbf{w}^1	p ^a	ICC	MDC ₉₅ %	selected	source	domain	metric
1	$0.235 {\pm} 0.140$	0.015	0.977 [0.98, 0.99]	13.6%	100%	phone	LTS	Mean transition time from standing to lying [sec].
2	$-0.134 {\pm} 0.063$	< 0.001	0.820 [0.74, 0.88]	30.5%	100%	watch	AF	Average hazard of non-MVPA to MVPA bouts.
3	-0.121 ± 0.095	< 0.001	0.818 [0.74, 0.88]	197.6%	100%	watch	TVDA	Daily average total time in MVPA bouts [mins].
4	$-0.117 {\pm} 0.058$	< 0.001	0.862 [0.80, 0.91]	147.7%	100%	watch	TVDA	Median acc. magnitude while in MVPA bouts [mg].
5	-0.096 ± 0.068	< 0.001	0.820 [0.74, 0.88]	194.6%	100%	watch	AF	Average consecutive duration in MVPA bouts [mins].
6	$0.237 {\pm} 0.180$	0.003	0.986 [0.98, 0.99]	13.6%	80%	phone	LTS	Mean transition time from lying to standing [sec].
7	-0.113 ± 0.104	< 0.001	0.902 [0.86, 0.94]	55.6%	80%	watch	TVDA	Daily percent of time spent walking.
8	-0.060 ± 0.055	< 0.001	0.872 [0.82, 0.92]	143.2%	80%	watch	TVDA	Average acc. while in MVPA [mg].
9	$0.187{\pm}0.196$	0.072	0.982 [0.97, 0.99]	3.0%	60%	phone	WRT	Median range-of-motion (ROM) [deg]
10	$0.092{\pm}0.116$	0.14	0.833 [0.76, 0.89]	69.3%	60%	watch	MORN	SD acc. value 30 mins after wake-up [mg].
11	$0.069 {\pm} 0.081$	0.109	0.618 [0.46, 0.75]	147.1%	60%	watch	NTR	Midpoint of night-time sleep window [hours]
12	$0.184{\pm}0.267$	0.154	0.735 [0.62, 0.83]	220.2%	40%	watch	NTR	# of movements during night-time sleep [count/hr].
13	-0.083 ± 0.177	0.741	0.795 [0.71, 0.87]	16.5%	40%	watch	NTR	Awake period during night-time sleep [mins].
14	$0.079 {\pm} 0.114$	0.273	0.987 [0.98, 0.99]	9.8%	40%	phone	STS	Mean transition time from sitting to standing [sec].
15	$-0.030 {\pm} 0.059$	0.015	0.985 [0.98, 0.99]	17.2%	40%	phone	WRT	Range-of-motion (ROM) median velocity [deg/sec].

Abbreviations (Abbrv.): acc., acceleration; sec, seconds; mins, minutes; h, hours; deg, degrees; deg/sec, degrees per second; ROM, range of motion; m, meters; m/s, meters per second; mg: mili-gravity units of acceleration; SD, standard deviation.

¹ Refers to the mean \pm standard deviation in LR feature coefficient values, w, over all CV folds;

^a Differences in feature distributions between RA / non-RA participants were investigated using a Mann-Whitney U Test.

p-values were post-hoc corrected using methods described by Benjamini and Hochberg⁴;

ICC, intraclass correlation coefficient [95% confidence interval]; MDC_{95%}, minimal detectable change (in percent, %).

A.4 Extended Results: Distinguishing RA severity levels

Table A.IV. Comparison of RA severity level prediction using patient reported outcomes (PRO), versus using PRO + sensor-outcomes, over 5-fold cross-validation (CV) with fortnightly (i.e., study duration) averaged active + passive features. Results are presented as: (1) the posterior overall *subject-wise* outcome for one cross-validation (CV) run as well as (2) the *observation-wise* median and inter-quartile range (IQR) across that CV in brackets. The best performing model for each feature representation are highlighted in **bold**. auroc: Area under the receiver operator curve; κ , Cohen's Kappa statistic; F₁, Macro-F1 score.

features	model	auroc	К	f1
PRO	LR-SG-lasso	0.736 (1.000, 0.750–1.000)	0.403 (0.286, 0.286–0.667)	0.733 (0.667, 0.667–0.857)
	LR-elastic-net	0.819 (1.000, 0.750–1.000)	0.479 (0.286, 0.286–0.615)	0.759 (0.750, 0.667–0.800)
PRO + active	LR-SG-lasso	0.747 (1.000, 0.833–1.000)	0.327 (0.286, 0.118–0.615)	0.710 (0.667, 0.667–0.800)
	LR-elastic-net	0.901 (1.000, 0.833–1.000)	0.560 (0.545, 0.333–1.000)	0.750 (0.667, 0.571–1.000)
PRO + passive	LR-SG-lasso	0.791 (1.000, 0.900–1.000)	0.479 (0.286, 0.286–1.000)	0.759 (0.667, 0.667–1.000)
	LR-elastic-net	0.879 (1.000, 0.925–1.000)	0.555 (0.545, 0.286–0.667)	0.786 (0.800, 0.750–0.857)
PRO + active + passive	LR-SG-lasso	0.885 (1.000, 0.906–1.000)	0.479 (0.286, 0.286–0.545)	0.759 (0.667, 0.667–0.750)
	LR-elastic-net	0.907 (1.000, 0.889–1.000)	0.707 (1.000, 0.667–1.000)	0.833 (1.000, 0.667–1.000)

Table A.V. Top 10 selected features from PRO + sensor-outcome based RA severity level estimation, as determined by LR-elastic-net across 5-fold subject-wise cross validation (CV), with fortnightly (i.e., study duration) averaged features.

	\mathbf{w}^1	p ^a	selected	source	domain	metric
1	$1.876 {\pm} 0.670$	0.002	100%	PRO	RASIQ	Joint pain
2	$1.696{\pm}1.066$	0.07	100%	phone	PEG	9HPT total time [sec]
3	$1.324{\pm}0.762$	0.45	100%	-	DEM	Age range [5 years)
4	$1.177 {\pm} 1.128$	0.13	100%	phone	WRT	Range-of-motion (ROM) median velocity [deg/s]
5	$0.971 {\pm} 0.774$	0.09	100%	watch	TVDA	# continuous periods of walking > 30 mins [count]
6	-0.410 ± 0.276	0.53	100%	-	DEM	Sex [M/F]
7	$-0.965 {\pm} 0.897$	0.10	80%	phone	WRT	Range-of-motion (ROM) [deg]
8	$0.543 {\pm} 0.825$	0.07	80%	PRO	RASIQ	Joint stiffness
9	$0.430{\pm}0.410$	0.12	80%	watch	AF	Average hazard of sedentary to non-sedentary bouts
10	$0.904{\pm}0.887$	0.02	60%	PRO	HAQ	HAQ-DI total score

Abbreviations (Abbrv.): PRO, patient-reported outcome; DEM, demographics information; acc., acceleration; s, seconds; mins, minutes; h, hours; deg, degrees; deg/s, degrees per second; ROM, range of motion; m, meters; m/s, meters per second; mg: mili-gravity units of acceleration. ¹ Refers to the mean \pm standard deviation in LR feature coefficient values, *w*, over all CV folds;

^a Differences in feature distributions between RA (mod) / RA (sev) participants were investigated using a Mann-Whitney U Test.

p-values were post-hoc corrected using methods described by Benjamini and Hochberg⁴⁷;

B Methodology: Smartwatch sensor feature extraction

B.1 Sensor Processing Pipeline

The sensor processing pipeline developed for the Apple Watch in the weaRAble-PRO study is outlined in Fig. B.I, yielding unobtrusively measured summary features of physical activity and sleep for RA participants, computed daily during normal life.



Figure B.I. Sensor processing pipeline developed for the Apple Watch in the weaRAble-PRO study. The raw 3-axis accelerometer sensor data, $(\mathbf{a}_x, \mathbf{a}_y, \mathbf{a}_z)$, denoted by $\mathbf{X} \in \mathbb{R}^{N \times 3}$, where *N* are the number of continuously collected accelerometer samples over the 14-day study period ($N = 50 \ [Hz] \times 3600 \ [sec] \times 24 \ [hr] \times 14 \ [days]$), can be transformed into measures characterising physical activity and sleep, $\mathbf{X}^* \in \mathbb{R}^{N \times 3}$, where *M* is the new sampling range, *daily* (M = 14), and *P* are the number of measures of daily life (i.e., features). In this case, a $N \gg M$ problem has been reduced into useful $M \times P$ features, unobtrusively measuring the physical activity and sleep of RA participants during daily life.

An overview of the pipeline is as follows:

1. Pre-processing:

- (a) Data chunking, memory optimisation. Convert raw 50 Hz accelerometer data to daily chunks;
- (b) Calibration to local gravity, local UTC timestamps;
- (c) Resampling, 30 Hz;
- (d) Butterworth, low-pass filtering at 17 Hz;
- (e) Non-wear detection and segmentation;

2. Daily-activity Estimation:

- (a) Epoch windowing, 30 [sec];
- (b) Activity classification per epoch;
- (c) Epoch aggregation, daily;
- (d) Posterior activity prediction with hidden Markov model (HMM) smoothing, see section B.3;

3. Characterising daily life:

(a) Physical activity and sleep feature extraction, see section B.5 for more details;

B.2 Deep Network-based Activity Recognition

In this work, a deep learning-based activity recognition model, known as a deep convolutional neural network (DCNN), was trained on Capture-24 and then used to directly estimate daily activity in the weaRAble-PRO study.

B.2.1 Multi-Task Self-Supervised Learning

Developing robust activity classification models is challenging in clinical studies due to the lack of labelled data for training. Deep networks, in particular, need a lot of training data in order to be robust and generaliseable. Open-source HAR-based datasets have small sample sizes, with generally n < 100 participants as annotating free-living wearable data for human activity recognition (HAR) requires a concurrent video stream, and the labelling process is resource-intensive²⁶.

There are however massive-scale unlabelled wearable datasets, such as the UK Biobank (UKB), which have collected data on roughly 100,000 participants with over 6 billion samples available. This study build upon our previous work demonstrating how advances in self-supervised learning (SSL) could help exploit the hidden information in these large-scale unlabelled datasets²⁶. SSL consists of training a model on a pretext task in an unlabelled dataset (often in a multi-task problem). The supervised task

is devised based on labels manually created in the unlabelled dataset, such as distinguishing transformed versus original data. The SSL model has to determine for each sample, if a transform has been applied or not, and what transform(s) have been applied as a multi-task problem. Essentially you create a robust deep feature extractor, built on a diverse and large amount of data. This pre-trained model can then be fine-tuned on a downstream task, such as activity recognition in the smaller datasets, such as Capture-24. In the main text, Fig. 7 illustrates a multi-task self-supervised approach for feature learning in HAR. We treated each of the tasks as a binary problem predicting whether a self-supervised transformation has been applied. Our multi-task SSL training relied on the *unlabelled UKB*, which contains roughly 700,000 person-days of free-living activity data (100,000 participants, 7 days of wear). For more information we refer the reader to our previous work²⁶.

B.2.2 Deep Network Architecture

We used a deep convolutional neural network (DCNN) with a ResNet-V2 architecture, consisting of 18 layers and 1D convolutions⁵⁷ as a feature extractor (10M parameters). The learned feature vector was of size 1024. All the tasks will share the same feature extractor. Then, we attached a softmax layer for each of the self-supervised tasks. In the downstream evaluation, we added a fully-connected (FC) layer of size 512 in between the feature extractor and softmax readout. The network structure was fixed for all the downstream evaluations. We computed the cross-entropy loss for each task and weighed all the tasks equally in loss calculation.

B.3 Hidden Markov Model (HMM) smoothing

Human activity recognition (HAR) model predictions are essentially independent—meaning that the sequence of activities over each 30 second epoch incorporates no temporal information epoch-to-epoch, for instance how the previous epoch prediction affects the current, or next, activity prediction. In order to add temporal dependency to the human activity recognition (HAR) model developed a Hidden Markov Model (HMM) was implemented in a post-processing step to obtain a more accurate sequence of predicted activities over the continuous 14-day data collection period.



Figure B.II. Diagram of a Hidden Markov Model (HMM). The sequence of discrete hidden states $\mathbf{h} = \{h_1, \dots, h_{(t-1)}, h_t, h_{(t+1)}, \dots, h_N\}$ form a Markov chain. At each time step an observation is obtained by a draw from a probability distribution that is conditional on the value of *h* at that time. This results in a sequence of observable values $\mathbf{x} = \{x_1, \dots, x_{(t-1)}, x_t, x_{(t+1)}, \dots, x_N\}$.

Hidden Markov Model (HMM) The Hidden Markov Model (HMM) defines a Markov chain on hidden (or "latent") variables $h_t = \{h_1, h_2, ..., h_H\}$, such that only the recent past influences the future:

$$p(h_t|h_{(1:t-1)},...,h_{(t-1)}) = p(h_t|h_{t-1})$$

The observed (or "visible") variables are dependent on the hidden variables through an emission $p(x_t|h_t)$. This defines a joint distribution:

$$p(\mathbf{h}|\mathbf{x}) = p(x_1|h_1)p(h_1)\prod_{t=2}^{N} p(x_t|h_t)p(h_t|h_{t-1})$$
(B.1)

where $p(x_t|h_t)$ defined the emission probability; $p(h_t|h_{t-1})$ defines the transition probability between hidden states; The transition distribution $p(h_{t+1}|h_t)$ is defined by a $H \times H$ transition matrix: $\mathbf{A}_{i',i} = p(h_{t+1} = i'|h_t = i)$. The emission distribution, $p(x_t|h_t)$, has discrete states $x_t \in \{1, ..., V\}$, we can define a $V \times H$ emission matrix: $\mathbf{B}_{i,j} = p(x_t = i|h_t = k)$. For continuous outputs, h_t selects one of H possible output distributions $p(x_t|h_t)$, $h_t \in \{1, ..., H\}$. The most likely hidden path, i.e., sequence of states, arg max{ $p(\mathbf{h}|\mathbf{x})$ }, is then found via the is found via the Viterbi algorithm. A diagram of a HMM is shown in Fig. B.II.

At each time step t, we define h_t as one of k activity classes, $\{c_1, c_2, ..., c_k\}$ as {sleep, sedentary, light physical activity, moderate-to-vigorous physical activity (MVPA)}. While h_t are not observed directly, at each t step there is an dependent

Table B.I. Comparison of activity recognition performance in the Capture-24 dataset between baseline random forest (RF) model and ResNet-based deep convoutional neural network (DCNN), pre-trained on 700,000 person days in the UK Biobank following a self-supervised learning (SSL) framework. κ , Cohen's kappa statistic; F₁, macro-F1 score.

model	к	f1
RF	0.705 ± 0.103	0.704 ± 0.102
RF + HMM	0.813 ± 0.108	0.775 ± 0.117
DCNN (SSL)	0.760 ± 0.087	0.735 ± 0.091
DCNN (SSL) + HMM	$\textbf{0.862} \pm \textbf{0.088}$	$\textbf{0.815} \pm \textbf{0.103}$

* HMM: hidden markov model posterior smoothing, see section B.3 for more details.

	LTS ·	0	0.01	0.56	0.27	0.16
ť	PEG ·	0.01	0.02	0.3	0.39	0.27
uided Tes	STS ·	0.01	0.01	0.54	0.31	0.14
פֿ	Walk ·	0.01	0.01	0.59	0.25	0.14
	Wrist ·	0.01	0.01	0.43	0.4	0.15
		imputed	MVPA Pred	light s	edentary	' sleep

Figure B.III. Validation of SSL-HMM activity predictions in the weaRAble-PRO study. Normalised confusion matrix evaluation of the SSL DCNN HAR model on the weaRAble-PRO study with guided test timings as pseudo-labels.

observed stochastic emission x_t . The hidden state sequence **h** is defined as the true activity labels and the emission distribution $p(x_t|h_t)$ is estimated directly using the predicted activity probabilities from the HAR model in the training set. As such we use the training predictions of activity from HAR model to infer the most likely sequence of true activity states that would have given rise to those predictions.

HMM smoothing helps to correct for erroneous predictions, such as when the transitions between those two classes of activity are rare, for instance sleeping to walking.

B.4 Evaluation of activity recognition model

The performance of the SSL model, compared to a feature-based Random Forest (RF) as a baseline, is reported in table B.I. It was observed that the SSL model improved activity recognition performance in Capture-24 beyond feature-based approaches and training a model end-to-end.

Furthermore, an insight into the performance of the SSL DCNN HAR model in the weaRAble-PRO dataset was obtained from a set of simple experiments. While the Capture-24 study uses a similar device (Axivity AX3) and placement (non-dominant wrist) to the weaRAble-PRO's Apple Watch, the evaluation of the HAR predictions are unknown due to the lack of activity labels in the weaRAble-PRO study. However, as participant's performed prescribed guided tests during the day, HAR predictions could be benchmarked against the timing of these assessments. As such, guided test timings can act as pseudo-labels in order to evaluate the SSL DCNN HAR model's robustness in applied to the weaRAble-PRO study, shown in Fig. B.III. As expected, during activity-based guided test assessments, such as walking or sit-to-stand and lie-to-stand, the HAR model more often predicts that participants are performing light activity. Guided assessments that require participants to be stationary while performing the task, such as the PEG test or wrist ROM test, are predicted more as sedentary activities. The percentage of sleep-based predictions during the guided test assessments (although incorrect) are roughly in line with the overall probability

of daytime sleep, as observed depicted in figures A.Ic–A.Id. Further work is needed to fully characterise and appraise the predictions of daytime sleep, which are assumed to be incorrect predictions of sedentary activity.

B.5 Passive Features

The section below details the passively extracted, activity monitoring features in the weaRAble-PRO study. Activity monitoring features were developed based on broad activity prediction labels {sleep, sedentary, light physical activity, moderate-to-vigorous physical activity (MVPA)}^{28,29} and fine-grained activity prediction labels {sleep, sitting/standing, mixed, vehicle, walking, bicycling}²⁷

Measures of physical activity and sleep were summarised based on smartwatch actigraphy sensor data magnitude data per epoch, or aggregated by intensity levels, activity classification and bouts of activity. These features could be broadly grouped into physical activity domains: total volume of daytime activity (TVDA), fragmentation of activity (AF), along with two RA symptom-specific domains focusing on morning stiffness (MORN) and night-time restlessness (NTR):

- 1. **Total volume of daytime activity (TVDA)**: captures information around the overall physical activity and the ability to perform physical activity at varying levels of intensity during daytime, which are known to be altered in patients with RA^{13,58};
- 2. Activity Fragmentation (AF): Metrics in this domain attempt to capture information related to the ability to perform sustained activity. Frequent interruptions of physical activity may reflect a worse health condition. For example, RA patients may need to interrupt some activity due to increased joint pain.
- 3. Morning stiffness (MORN) is a common symptom of RA. These measures include information related to timing of activity after getting up in the morning as estimated from the activity classification¹.
- 4. **Night-time Restlessness** (**NTR**): These measures include information related to timing, duration and quality of sleep. There is evidence that RA patients experience fluctuations in disease activity following a circadian rhythm with worsening of the illness during the night⁵⁹, thus measures of overnight movement serve as a proxy for estimating night-time restlessness, reflecting the impact of the disease on sleep quality.

C Methodology: Smartphone sensor feature extraction

C.1 Smartphone guided tests administered in weaRAble-PRO

Table C.I. Overview of the daily prescribed smartphone-recorded Active Assessments (denoted "Guided Tests") in the weaRAble-PRO study.

Exercise	Abbrv.	Brief Description
Wrist ROM	WRT	Participants were requested to sit down and place their forearm at the edge of a table, holding the iPhone horizontally facing up in their hand, and to flex and extend their wrist joint to its maximum angle, repeating the motion for 10 seconds. The test was carried out twice, once using each hand ¹² .
Gait	WLK	Participants were asked to affix the iPhone to their leg (on the right thigh) facing outward using a provided strap, then to walk in a straight line for 30 seconds (while being allowed to turn around at any point in the middle of the test) ¹² .
9-hole peg test	PEG	Digital touch-screen version of the standard clinical assessment where participants are asked to place, and subsequently remove, 9 pegs into and from a round hole, in the fastest time possible ³³ .
Sit-to-stand	STS	Participants were requested to perform a sit-to-stand transition from a chair with the iPhone attached to their upper right thigh with a strap, repeating the exercise 5 times at their own pace ^{$30, 34$} .
Lie-to-stand	LTS	Participants performed a lie-to-stand transition: from lying still with legs stretched on a bed, to standing up on the floor. Participants were requested to affix the iPhone to their right thigh with a strap during the exercise and repeat twice at their own pace ^{30, 34} .

Abbreviations (abbrv.): ROM, range of motion;

C.2 Peg Test Algorithm

The 9-hole peg test is a two-step test of hand dexterity to measure the MSFC score in Multiple Sclerosis, or signs of Parkinson's disease or stroke⁶⁰. The digital 9-hole peg test (9HPT), and subsequent metrics calculated, are proprietary to Apple ResearchKit, see http://researchkit.org/docs/docs/ActiveTasks/ActiveTasks.html for more details. This smartphone version of the 9HPT task requires participants to use two fingers to touch the on-screen peg and drag it into the an on-screen hole. Once all 9 pegs have been placed in each hole, the task is repeated by removing the pegs in the same manner. The 9HPT is repeated using both dominant and non-dominant hands. The total duration that the user spent taking the test for each hand repetition is recorded.

C.3 Wrist ROM Test Algorithm

The wrist range of motion (ROM) test algorithm is outlined previously as part of the GSK PARADE study^{11,12}. The iPhone accelerometer sensor data is converted to angular positions; ROM is then computed based on the differences between angular maxima; angular velocity is determined from the gyroscope sensors.

C.4 Walk Test Algorithm

The gait test algorithm is outlined previously as part of the GSK PARADE study¹². During walking, initial (IC) and final (FC) feet contact timings are calculated from the smartphone accelerometer data with an inverted pendulum model, as described in^{61,62}. Contact points are determined by integration of the vertical component of the accelerometer signal, \mathbf{a}_y , and subsequent differentiation of that signal with a continuous wavelet transform (CWT, convolution of the accelerometer data and an analysing function, i.e., mother wavelet). ICs and FCs are then denoted by minima and maxima timings respectively in this transformed signal. Detected peaks were used to estimate the number of steps, cadence, step length⁶³, and walk velocity for each assessment.

C.5 Sit-to-stand Test Algorithm

Both the iPhone accelerometer and gyroscope data were used to determined when participants were in a sitting or standing position during the sit-to-stand (STS) exercise³⁰. The vertical accelerometer component, \mathbf{a}_y , determined sit-to-stand transitions by estimating the phone axis orientation (and change thereof) relative to gravity, as the participants moved between sitting

and standing—given the phone's axis should be fixed as it is strapped to the participants thigh using a strap. Gyroscope axis sensor components, $(\mathbf{g}_x, \mathbf{g}_y, \mathbf{g}_z)$, helped determine whether a participant had completed valid transitions during the exercise. The start and end points of the STS transitions could therefore be determined by identifying peaks in \mathbf{a}_y that were within the given thresholds where participants could be feasibly standing or lying. The length of a sit-to-stand transition was then calculated as the difference between the point where the participant first reached a standing position and the final point where the participant was sitting prior to beginning a standing motion.

C.6 Lie-to-Stand Test Algorithm

Following a similar STS algorithm³⁰, during the lie-to-stand test, accelerometer and gyroscope (y-axis) measurements were used to determine a participant's standing and lying transition points following the algorithm introduced in⁶⁴. Gyroscopic y-axis gravity relates the phone's orientation relative to gravity at a given moment in time and the phone acceleration helped to improve accuracy in determining participants' lying positions—for example, the point of minimum acceleration generally corresponded to moments when participants held a lying position (causing a short plateau in y-gravity).

D Further Study Details

D.1 Inclusion and Exclusion Criteria

For participation in the full pilot study, 30 RA participants and 30 HCs matched on age, gender, and race were recruited. Due to the small sample size and potentially limited pool for recruitment, ages were matched within a window of \pm 3 years.

The overall ratio of moderate to severe participants was chosen not to exceed 2:1 in either direction. Inclusion/Exclusion criteria were the same for RA participants and HCs, unless otherwise noted. All participants must have been able and willing to perform the pre-defined guided tests at the start of the study, follow audio instructions from an iPhone, and have a sufficient level of English language to ensure ability to understand mobile app and questionnaires.

Rheumatoid Arthritis participants, at least 21 years of age at date of consent for study, were selected based on clinically verified diagnosis of moderate-to-severe RA, with severity assessed using Routine Assessment of Patient Index Data 3 (RAPID3). Healthy controls were selected based on no prior or current diagnosis of a rheumatological disorder, inflammatory disorder, malignancy, or other relevant diseases. Further exclusion criteria for all participants included history of other inflammatory rheumatologic or systemic autoimmune disorder (e.g., Hashimoto's thyroiditis or Sjogren's syndrome), history of movement disorders, other neurological disorders or conditions resulting in significant physical impairments that impact joint movements to be assessed, history of postural hypotension, unexplained syncope, or other conditions that make it difficult for participants to perform guided tests such as the lie-to-stand test and any history of severe skin allergy. Participants were also excluded if they required use of a wheelchair, walking aids, artificial limbs, or had any active implantable device, such as a pacemaker.

D.2 Data collection timelines



Figure D.I. Overview of frequency and timeline of data capture for each data source in the weaRAble-PRO Study.

D.3 Patient-reported outcomes administered in weaRAble-PRO

Table D.I. Overview of the PRO assessments administered to RA and HC participants on days 1, 7, and 14 in the weaRAble-PRO study

PRO	Abbrv.	Domain(s) assessed
GSK RA symptom and impact questionnaire	RASIQ	Generalised measure of the severity of RA symptoms and their impact on the $patient^{65}$.
Patient-Reported Outcomes Measurement In- formation Systemt ¹	PROMIS pain PROMIS sleep	Pain interference developed to assess the degree to which pain interferes with participants' physical, mental, and social activities ⁶⁶ . Sleep disturbance as assessed self reported perceptions of sleep quality, including perceived difficulty in falling asleep, difficulty experienced in staying asleep, sleep depth, and satisfaction with sleep quality ⁶⁷ .
Health Assessment Questionnaire-Disability Index	HAQ-DI	Self reported functional status measures. It is one of the most widely used measure of function in RA, with demonstrated reliability and validity in RA patients ⁶⁸ .
Functional Assessment of Chronic Illness Therapy Fatigue	FACIT	Assesses four domains of fatigue ⁶⁹ : physical fatigue, functional fatigue, emotional fatigue, social consequences. This questionnaire has been validated for use with RA patients ^{70,71} .
Short-Form 36	SF-36	36-item questionnaire to allow participants to self-assess functional health and well-being ⁷² . Scores are provided for eight domains: general health, mental health, physical functioning, social functioning, physical role, emotional role, bodily pain, vitality. This instrument has been validated and used in studies with RA patients ^{73,74} .
Interactive joint-pain map ²	JMAP	Records the number and severity of up to 55 pre-specified joints experienced by the participant at a given time. For the joints where patients are experiencing any pain at the given moment, they are asked to score the level of pain as 1 (mild pain), 2 (moderate pain), or 3 (severe pain) ¹² .
Visual analogue scale ²	VAS	Patient's assessment of arthritis pain is a single-item question that assesses the level of pain severity the participant is currently experiencing using a visual analogue scale ranging from 0 to 100^{75} .

Abbrv., abbreviations ¹ An reduced set of items from the PROMIS item bank for pain and sleep domains were used for this study;

² Administered daily over the 14-day study, within an hour of completing the predefined guided tests, once in the morning, and once in the afternoon.

E List of Extracted Features

	Feature	Source	Dom.	Description
0	AvgBoutLen MVPA			Avg length of MVPA bouts
1	AvgBoutLen_light			Avg length of active (light) bouts
2	AvgBoutLen_nght			Avg. length of sedentary houts
3	AvgHazard MVPAToany			Avg hazard ¹ of MVPA to non-MVPA
4	AvgHazard anyToMVPA			Avg. hazard of non MVPA to MVPA
5	AvgHazard anyTosodontory			Avg. hazard of non-sedentary to sedentary
6	AvgHazard_adoptery/Toopy			Avg. hazard of sodontary to non-sodontary
7	Avglazaru_sedentary toany			Avg. length of conceptive time in MVDA
,	AvgLenTimeActive_WVPA			Avg. length of consecutive time in wry PA
8	AvgLenTimeActive_light			Avg. length of consecutive time active
9	AvgLenTimeActive_sedentary			Avg. length of consecutive time sedentary activity
10	BoutsHazard_MVPAToany		AE	Avg. hazard of MVPA to non-MVPA bouts
11	BoutsHazard_any IoM VPA	watch	AF	Avg. hazard of non-MVPA to MVPA bouts
12	BoutsHazard_anyTosedentary			Avg. hazard of non-sedentary to sedentary bouts
13	BoutsHazard_sedentaryToany			Avg. hazard of sedentary to non-sedentary bouts
14	BoutsTransitionPr_MVPAToany			Transition probability of MVPA to non-MVPA bouts
15	BoutsTransitionPr_anyToMVPA			Transition probability of non-MVPA to MVPA bouts
16	BoutsTransitionPr_anyTosedentary			Transition probability of non-sedentary to sedentary bouts
17	BoutsTransitionPr_sedentaryToany			Transition probability of sedentary to non-sedentary bouts
18	RatioBoutsToActive			Ratio of time in active bouts to overall time active
19	TransitionPr_MVPAToany			Transition probability of active to sedentary (acc. intensity defined)
20	TransitionPr_anyToMVPA			Transition probability of active to sedentary (acc. intensity defined)
21	TransitionPr_anyTosedentary			Transition probability of active to sedentary (acc. intensity defined)
22	TransitionPr_sedentaryToany			Transition probability of active to sedentary (acc. intensity defined)
22				A
23	age_range		DEM	Age range [5 years)
24	56.7			Stx [W/1]
25	LTS_mean_lie2stand			Mean lie-to-stand transition time [s]
26	LTS_mean_lying	nhono	ITS	Mean lying time [s]
27	LTS_mean_stand2lie	phone	LIS	Mean stand-to-lie transition time [s]
28	LTS_mean_standing			Mean standing time [s]
20	Morning stiffness: aug. 0:15:00			
30	Morning stiffness: auc - 0:15:00			
31	Morning stiffness: auc - 0:45:00			Daily AUC of acc. vector magnitude during the first $n=115$ 30 45 60
32	Morning stiffness: auc 1:00:00			120, 2401 mine after getting up
32	Morning stiffness: auc - 2:00:00			120, 240J mins after getting up
34	Morning stiffness: auc. 4:00:00			
35	Morning stiffness: mean 0:15:00			
26	Morning stiffnass: mean 0:20:00			
27	Morning stiffness: mean - 0.50.00			Daily avg. acc. vector magnitude during the first n=[15, 30, 45, 60,
29	Morning stiffness: mean 1:00:00			120, 240] mins after getting up
20	Morning stiffness: mean - 2:00:00			
39	Morning summess: mean - 2:00:00			
40	Morning summess: mean - 4:00:00			
41	Morning summess: $q(0.5) = 0.15:00$			
42	Morning stiffness: $q(0.5) = 0.30:00$		MODN	
43	Morning stiffness: $q(0.5) = 0.45:00$	watch	MORN	Daily median of acc. vector magnitude during the first n=[15, 30, 45,
44	Morning stiffness: $q(0.5) - 1:00:00$			60, 120, 240] mins after getting up
45	Morning stiffness: $q(0.5) - 2:00:00$			
46	Morning stiffness: $q(0.5) - 4:00:00$			
4/	Morning stiffness: $q(0.95) - 0.15:00$			
48	Morning stiffness: $q(0.95) - 0.30:00$			
49	Morning stiffness: $q(0.95) - 0.45:00$			Daily 95^{m} centile of acc. vector magnitude during the first n=[15, 30,
50	Morning stiffness: $q(0.95) - 1:00:00$			45, 60, 120, 240] mins after getting up
51	Morning stiffness: q(0.95) - 2:00:00			
52	Morning stiffness: q(0.95) - 4:00:00			
53	Morning stiffness: stdev - 0:15:00			
54	Morning stiffness: stdev - 0:30:00			
55	Morning stiffness: stdev - 0:45:00			Daily avg. of Standard deviation of acc. vector magnitude during the
56	Morning stiffness: stdev - 1:00:00			first n=[15, 30, 45, 60, 120, 240] mins after getting up
57	Morning stiffness: stdev - 2:00:00			
58	Morning stiffness: stdev - 4:00:00			
59	MSleep			Daily avg. midpoint time of night-time sleep window [mins]
60	NSleepEpisodes			Avg. number of sleep episodes per night-time sleep window [count]
61	PercentSleep			Avg. percent time of sleep per night-time sleep window [%]
62	RestEfficiency			Avg. percent night-time rest efficiency per night-time sleep window
63	RestFragmentation			Avg. night-time rest fragmentation per night-time sleep window
64	SleepDur	watch	NTR	Sleep duration [h]
65	SleepMov			The number of movement episodes per night-time sleep window
-				[count]
		1	1	• •

Table E.I. Description of the sensor-based features extracted in the weaRAble-PRO study

Continued on Next Page...

¹ The average hazard reflects the probability for an individual to remain in a specified activity at minute *t*, or change to any other activity, given that the subject has been in a specified activity up to minute $t - 1^{76}$.

66 SleepMovPerHr The number of sleep movements per hour per n (counthr) 67 TotDaySleep PBG. totalFailures phone PBG 68 PEG.rotalFailures phone PBG The total # of per failures [count] 70 STS_mean_sit2Stand The total # of per failures [count] The total # of per failures [count] 71 STS_mean_stand2sit Pbone STS Mean sit-to-sit transition time [s] 73 STS_mean_stand2sit Mean sit-to-sit transition time [s] Mean sit-to-sit transition time [s] 74 95thConiteAccMag The study-arg.acc. in stochary bouts [ns.*1] The study-arg.acc. in stochary bouts [ns.*2] 74 PashConitbour.MVPA The study-arg.acc. in stochary bouts [ns.*2] The daily acc. in sechary bouts [ns.*2] 75 DailyAcclinbour.MVPA The daily acc. in sechary bouts [ms.*2] The daily acc. in sechary bouts [ms.*2] 76 DailyAcgAcclibour.MVPA DailyAcgAcclibour.MVPA DailyAcgAcc.240-200ng 76 DailyAcgAcclibour.MVPA DailyAcgAcclibour.MVPA DailyAcgAcclibour.MVPA 77 DailyAcgAcclibour.MVPA DailyAcgAcclibour.MVPA DailyAcgAcclibour.MVPA	
67TotDaySkepThe total # of peg fullures [count] The total # of peg fullures [count] The total # of peg fullures [count] The total # of peg fullures [count]68PEG_totalFailuresphonePEG70STS_mean_sitt2stand 	night-time sleep window
68 PEG_toalFinitures phone PEG The total 9 dPT time [s] 70 STS_mean_sit2stand The total 9HPT time [s] 71 STS_mean_stat2stand 73 STS_mean_stat2stand 74 95thCentileAccMag 75 ArgAcchBoul_MVPA 76 ArgAcchBoul_MVPA 77 ArgAcchBoul_right 78 ArgAcchBoul_right 79 DailyAcchBoul_right 79 DailyAcchBoul_right 71 Presconduct [ms] 78 ArgAcchBoul_weatherap 79 DailyAcchBoul_sedentary 79 DailyAcchBoul_sedentary 70 DailyAcchBoul_sedentary 71 DailyAcchBoul_sedentary 73 DailyArgAccuPang 74 DailyArgAccuPang 75 DailyArgAccuPang 76 DailyArgAccuPang 77 DailyArgAccuPang 78 DailyArgAccuPang 79 DailyArgAccuPang 74 DailyArgAccuPang 75 DailyArgAccuPang 76 DailyArgAccuPang 77 DailyArgAccuPang 78 DailyArgAccuPang 79 DailyArgAccuPang 70 DailyArgAccuPang	
70 STS_mean_sit2stand 71 STS_mean_sita_stand 72 STS_mean_standLing 73 STS_mean_standLing 74 95thCentileAccMag 75 ArgAcclnBout_MVPA 76 ArgAcclnBout_MVPA 77 ArgAcclnBout_MVPA 78 ArgAcclnBout_MVPA 79 DailyAcclnBout_edentary 79 DailyAcclnBout_edentary 70 DailyAcclnBout_edentary 71 ArgAcclnBout_MVPA 70 DailyAcclnBout_edentary 71 DailyAcclnBout_edentary 72 DailyAcclnBout_edentary 73 DailyArgAcc_20-040mg 74 DailyArgAcc_200-240mg 75 DailyArgAcc_200-240mg 76 DailyArgAcc_200-240mg 77 DailyArgAcc_200-240mg 78 DailyArgAcc_200-240mg 79 DailyArgAcc_200-240mg 70 DailyArgAcc_200-240mg 71 DailyArgArc_20-20-240mg 72 DailyArgArc_20-20-240mg 73 DailyArgTime_10-200mg 74 DailyArgTime_10-200mg 75 DailyArgTime_10-200mg 76 DailyArgTime_10-200mg 77 DailyArgTime_10-200mg 78<	
71 STS_mean_stining phone STS Mean stand-to-sit transition time [s] 73 STS_mean_standing mean stand-to-sit transition time [s] 74 95thCentileAccMag The 95 th centile of acc. vector magnitude [m.s. ²] 75 AvgAccInBout_MVPA The study-avg. acc. in NVPA bouts [m.s ²] 76 AvgAccInBout_dight The study-avg. acc. in scdentary bouts [m.s ²] 78 AvgAccInBout_scdentary The study-avg. acc. in scdentary bouts [m.s ²] 79 DailyAccInBout_scdentary The daily acc. in NVPA bouts [m.s ²] 71 DailyAccInBout_scdentary The daily acc. in scdentary bouts [m.s ²] 76 DailyAccCaloolut_scdentary The daily acc. in scdentary bouts [m.s ²] 77 DailyAccCaloolut_scdentary The daily acc. in scdentary bouts [m.s ²] 78 DailyAvgAcc_20-20-20mg DailyAvgAcc_20-20-20mg 79 DailyAvgAcc_20-20-20mg DailyAvgAcc_20-20-20mg 79 DailyAvgAcc_20-20-20mg DailyAvgTime_120-160mg 70 DailyAvgAcc_20-20-20mg DailyAvgTime_120-160mg 70 DailyAvgTime_120-160mg Daily avg. time spent in NPPA [ms ²] 70 DailyAvgTime_200-240mg Daily avg. time spent in ilght activity boats [mins] 71 DailyAvgTime_200-240mg Daily avg. time spent in ilght activity boats </td <td></td>	
73 STS_mean_standing Mean standing time [s] 74 95thCentileAccMag The standing time [s] 75 AvgAcchBout, MVPA The standing time [s] 76 AvgAcchBout, Jight The standing acc, in MVPA bouts [m,s ⁻¹] 77 AvgAcchBout, Jight The standing acc, in is edentary bouts [m,s ⁻¹] 78 DailyAcchBout_sedentary The standing acc, in is edentary bouts [m,s ⁻¹] 78 DailyAcchBout_sedentary The daily acc, in sedentary bouts [m,s ⁻¹] 78 DailyAcchBout_sedentary The daily acc, in sedentary bouts [m,s ⁻¹] 78 DailyAcchBout_sedentary The daily acc, in sedentary bouts [m,s ⁻¹] 79 DailyAvgAcc_2160-200mg The daily acc, in sedentary bouts [m,s ⁻¹] 79 DailyAvgAcc_20-200mg DailyAvgAcc_200-300mg DailyAvgAcc_200-300mg 79 DailyAvgAcc_200-400mg DailyAvgAcc_200-400mg DailyAvgAcc_200-400mg 70 DailyAvgAcc_200-200mg DailyAvgTime_10-40-40mg DailyAvgTime_10-40-40mg 70 DailyAvgTime_200-200mg DailyAvgTime_200-200mg DailyAvgTime_200-200mg 70 DailyAvgTime_200-200mg DailyAvgTime_10-40-400mg Daily avgt time spent in MVPA louts [%]	
7495hCentiLAccMagThe 95h ⁶ centile of acc. vector magnitude [ms. ²¹]75AvgAccInBout_MVPAThe study-avg. acc. in stelentary bouts [ms. ²¹]77AvgAccInBout_sedentaryThe study-avg. acc. in stelentary bouts [ms. ²¹]78DailyAccInBout_lightThe study-avg. acc. in stelentary bouts [ms. ²¹]79DailyAccInBout_lightThe daily acc. in MVPA bouts [ms. ²¹]78DailyAccInBout_gightThe daily acc. in sedentary bouts [ms. ²¹]79DailyAccInBout_gightThe daily acc. in sedentary bouts [ms. ²¹]71DailyAvgAcc.200-240mgThe daily acc. in sedentary bouts [ms. ²¹]72DailyAvgAcc.200-240mgThe daily acc. in sedentary bouts [ms. ²¹]70DailyAvgAcc.200-300mgDailyAvgAcc.200-300mg70DailyAvgAcc.200-300mgDailyAvgAcc.200-300mg71DailyAvgAcc.200-300mgDailyAvgAcc.200-300mg72DailyAvgAcc.200-200mgDailyAvgAcc.300-300mg73DailyAvgAcc.200-100mgDailyAvgTime.100-100mg74DailyAvgTime.100-200mgDailyAvgTime.200-100mg75DailyAvgTime.200-200mgDailyAvgTime.200-200mg76DailyAvgTime.200-200mgDailyAvgTime.200-200mg70DailyAvgTime.200-200mgDailyAvgTime.200-200mg71DailyAvgTime.200-200mgDailyAvgTime.200-200mg72DailyAvgTime.200-200mgDailyAvgTime.200-200mg73DailyAvgTime.200-200mgDailyAvgTime.200-200mg74DailyAvgTime.200-200mgDailyAvgTime.200-200mg75DailyAvgTime.200-200mgDa	
JAvgAcclinBout_MVPAThe study-avg_acc. in light bouts [ms^-1]76AvgAcclinBout_lightThe study-avg_acc. in light bouts [ms^-1]77AvgAcclinBout_lightThe study-avg_acc. in light bouts [ms^-1]78DailyAcclinBout_gedentaryThe study-avg_acc. in selentary bouts [ms^-1]79DailyAcclinBout_gedentaryThe daily acc. in MVPA bouts [ms^-1]78DailyAcclinBout_gedentaryThe daily acc. in MVPA bouts [ms^-1]79DailyAcclinBout_gedentaryThe daily acc. in MVPA bouts [ms^-1]71The daily acc. const [ms Points [ms^-1]72DailyAvgAcc_2-040mg73DailyAvgAcc_2-00-280mg74DailyAvgAcc_2-00-280mg75DailyAvgAcc_2-00-280mg76DailyAvgAcc_2-00-280mg77DailyAvgAcc_2-00-280mg78DailyAvgAcc_2-00-280mg79DailyAvgAcc_2-00-280mg70DailyAvgTime_1Bout_geht71DailyAvgTime_1Bout_geht72DailyAvgTime_200-240mg73DailyAvgTime_200-240mg74DailyAvgTime_200-240mg75DailyAvgTime_200-240mg76DailyAvgTime_30-400mg77DailyAvgTime_1ght78DailyAvgTime_1ght79DailyAvgTime_1do-200mg70DailyAvgTime_1do-200mg71DailyAvgTime_1do-200mg72DailyAvgTime_30-400mg73DailyAvgTime_1do-200mg74DailyAvgTime_1do-200mg75DailyPcTimeInBout_MVPA76DailyPcTimeInBout_geht76	.s ⁻²]
77AvgAcclnBout_sedentary AvgAccMagThe study-arg, acc. in sedentary bouts $[m,s^2]$ 78AvgAccMagThe study-arg, acc. in sedentary bouts $[m,s^2]$ 79DailyAcclnBout_lightThe daily acc. in light bouts $[m,s^2]$ 70DailyAcclnBout_sedentaryThe daily acc. in light bouts $[m,s^2]$ 78DailyAcclnBout_sedentaryThe daily acc. in sedentary bouts $[m,s^2]$ 78DailyAcclnBout_sedentaryThe daily acc. in sedentary bouts $[m,s^2]$ 78DailyAcclnBout_sedentaryThe daily acc. in sedentary bouts $[m,s^2]$ 79DailyAvgAcc_20-200mgDailyAvgAcc_20-200mg79DailyAvgAcc_20-200mgDailyAvgAcc_20-300mg70DailyAvgAcc_20-200mgDailyAvgAcc_20-40mg70DailyAvgAcc_20-200mgDailyAvgTimelnBout_sedentary70DailyAvgTimelnBout_sedentaryDaily avg. time spent in MVPA bouts [mins]70DailyAvgTime_20-20mgDailyAvgTime_20-20mg70DailyAvgTime_20-200mgDaily avgTime_sedentary bouts [mins]70DailyAvgTime_20-200mgDailyAvgTime_30-300mg70DailyAvgTime_30-300mgDailyAvgTime_sedentary71DailyAvgTime_40-80mgDaily avgTime spent in MVPA [mins]72DailyAvgTime_40-80mgDaily avgTime spent in MVPA [mins]73DailyAvgTime_40-80mgDaily avgTime spent in MVPA [mins]74DailyAvgTime_40-80mgDaily avgTime spent in fight activity fouts [not75DailyPet_10-00mgDaily percent of time spent in sedentary fouts [not74DailyPet_20-200mgDaily percent	
78AvgAccMagThe study-avg. acc. vector magnitude $[m.s^2]$ 79DailyAcclnBout_lightThe daily acc. in MVPA bouts $[m.s^2]$ 78DailyAcclnBout_lightThe daily acc. in light bouts $[m.s^2]$ 78DailyAygAcc_0-40mgThe daily acc. in sedentary bouts $[m.s^2]$ 78DailyAygAcc_0-40mgThe daily acc. in sedentary bouts $[m.s^2]$ 78DailyAygAcc_200-240mgDailyAygAcc_230-306mg79DailyAygAcc_230-306mgDailyAygAcc_230-306mg70DailyAygAcc_230-306mgDailyAygAcc_230-306mg71DailyAygAcc_230-400mgDailyAygAcc_230-306mg72DailyAygAcc_230-306mgDailyAygAcc_320-306mg79DailyAygAcc_1BhDaily AygAcc_320-40mg70DailyAygAcc_1BhDaily avg. acc. during MVPA $[m.s^2]$ 73DailyAygTime_1Bout_getnDaily avg. time spent in MVPA bouts $[m.s^2]$ 74DailyAygTime_120-160mgDaily avg. time spent in Ight activity bouts $[m.s^2]$ 76DailyAygTime_200-240mgDaily avg. time spent in ach 40 mg range $[m.s^2]$ 77DailyAygTime_200-240mgDaily avg. time spent in MVPA bouts $[mins]$ 78DailyAygTime_200-240mgDaily avg. time spent in MVPA [mins]79DailyAygTime_200-240mgDaily avg. time spent in MVPA [mins]70DailyAygTime_200-240mgDaily avg. time spent in MVPA [mins]71DailyAygTime_200-240mgDaily avg. time spent in MVPA [mins]72DailyAygTime_200-240mgDaily percent of time spent in MVPA [mins]73DailyAygTime_200-240mgDaily percent of ti	']
9DailyAcclinBout_Ight The daily acc. in light bouts $[m, s^2]$ The daily acc. in light bouts $[m, s^2]$ The daily acc. in light bouts $[m, s^2]$ The daily acc. in sedentary bouts $[m, s^2]$ DailyAvgAcc_200-240mg DailyAvgAcc_200-306mg DailyAvgAcc_200-306mg DailyAvgAcc_200-306mg DailyAvgAcc_200-306mg DailyAvgAcc_200-240mg DailyAvgTimeInBout_Sedentary DailyAvgTimeInBout_Sedentary DailyAvgTimeInBout_Sedentary DailyAvgTime_10-040mg DailyAvgTime_10-040mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_200-240mg DailyAvgTime_10-400mg DailyAvgTime_200-240mg 	
a b alijAcclnBout_sedentaryThe daily acc. in sedentary bouts $[m.s^2]$ B alijAvgAcc_D-040mgThe daily acc. in sedentary bouts $[m.s^2]$ B alijAvgAcc_120-160mgThe daily acc. in sedentary bouts $[m.s^2]$ B alijAvgAcc_20-200-240mgDailyAvgAcc_200-240mgB DailyAvgAcc_230-300mgDailyAvgAcc_230-300mgD D D DilyAvgAcc_230-300mgDailyAvgAcc_320-300mgD D<	
22 DailyAvgAccMag 33 DailyAvgAcc_140-00mg 44 DailyAvgAcc_210-160mg 55 DailyAvgAcc_200-240mg 56 DailyAvgAcc_200-240mg 57 DailyAvgAcc_200-240mg 58 DailyAvgAcc_230-360mg 59 DailyAvgAcc_30-300mg 50 DailyAvgAcc_30-360mg 50 DailyAvgAcc_40-80mg 50 DailyAvgAcc_230-400mg 51 DailyAvgTarce_160-200mg 50 DailyAvgTime_160ut_light 50 DailyAvgTime_100ut_light 51 DailyAvgTime_100-160mg 52 DailyAvgTime_160-200mg 53 DailyAvgTime_100-200mg 54 DailyAvgTime_100-200mg 55 DailyAvgTime_200-240mg 56 DailyAvgTime_200-240mg 56 DailyAvgTime_30-300mg 56 DailyAvgTime_30-300mg 57 DailyAvgTime_30-300mg 58 DailyAvgTime_30-300mg 59 DailyAvgTime_30-300mg 50 DailyAvgTime_30-300mg 50 DailyAvgTime_30-400mg 50 <	
 bailyAvgKac_10-0100mg bailyAvgKac_200-160mg bailyAvgKac_200-240mg bailyAvgKac_240-280mg bailyAvgKac_240-280mg bailyAvgKac_260-300mg bailyAvgKac_30-300mg bailyAvgKac_30-300mg bailyAvgKac_30-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_60-40-80mg bailyAvgKac_80-120mg bailyAvgKac_80-120mg bailyAvgTimeInBout_light bailyAvgTimeInBout_light bailyAvgTime_10-40mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_40-80mg bailyAvgTime_60-300mg bailyAvgTime_40-80mg bailyAvgTime_10-40-80mg bailyAvgTime_10-40-80mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyAvgTime_10-200mg bailyPet_10-200mg bailyPet_10-200mg bailyPet_10-200mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-200mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-200mg bailyPet_20-280mg bailyPet_20-200mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-200mg bailyPet_20-280mg bailyPet_20-280mg bailyPet_20-2	
 bailyAvgTic_100-200mg bailyAvgAc_200-240mg bailyAvgAc_200-280mg bailyAvgAc_200-400mg bailyAvgAc_200-400mg bailyAvgAc_20-4080mg bailyAvgAc_20-4080mg bailyAvgAc_20-4080mg bailyAvgAc_20-4080mg bailyAvgAc_20-4080mg bailyAvgAc_20-4080mg bailyAvgTimeInBout_MVPA bailyAvgTimeInBout_Sedentary bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_200-240mg bailyAvgTime_106-200mg bailyAvgTime_200-240mg bailyAvgTime_200-240mg bailyAvgTime_200-240mg bailyAvgTime_200-240mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_106-200mg bailyAvgTime_200-240mg bailyAvgTime_106-200mg bailyPetTimeInBout_light bailyPetTimeInBout_sedentary bailyPetTimeInBout_Sedentary bailyPetTimeInBout_Sedentary bailyPetTimeInBout_sedentary bailyPet_240-280mg bailyPet_	
86 DailyAvgAcc.200-240mg 87 DailyAvgAcc.200-280mg 88 DailyAvgAcc.200-320mg 90 DailyAvgAcc.230-300mg 90 DailyAvgAcc.320-300mg 91 DailyAvgAcc.300-400mg 92 DailyAvgAcc.30-400mg 93 DailyAvgAcc.30-400mg 94 DailyAvgAcc.30-120mg 95 DailyAvgTimeInBout_MVPA watch 96 DailyAvgTimeInBout_MVPA watch 97 DailyAvgTime.Bout_sedentary Daily avg. time spent in MVPA houts [mins] 98 DailyAvgTime.106-200mg Daily AvgTime.200-240mg 100 DailyAvgTime.200-240mg Daily AvgTime.200-240mg 101 DailyAvgTime.200-240mg DailyAvgTime.200-240mg 102 DailyAvgTime.200-240mg DailyAvgTime.306-400mg 103 DailyAvgTime.306-400mg DailyAvgTime.306-400mg 104 DailyAvgTime.306-400mg DailyAvgTime.106-200mg 103 DailyAvgTime.306-400mg DailyAvgTime.306-400mg 104 DailyAvgTime.106-200mg DailyAvgTime.106-200mg 111 DailyAvgTime.100ut_light Daily percent of time spent in MVPA [mins]	
88 DailyAvgAcc_240-280mg 88 DailyAvgAcc_280-320mg 90 DailyAvgAcc_320-360mg 91 DailyAvgAcc_360-80mg 92 DailyAvgAcc_360-80mg 93 DailyAvgAcc_360-120mg 94 DailyAvgAcc_S0-120mg 95 DailyAvgAcc_sedentary 96 DailyAvgTimeInBout_Light 97 DailyAvgTimeInBout_MVPA 98 DailyAvgTime.pather 99 DailyAvgTime.pather 90 DailyAvgTime.pather 91 DailyAvgTime.pather 92 DailyAvgTime.pather 93 DailyAvgTime.pather 94 DailyAvgTime.pather 95 DailyAvgTime.pather 96 DailyAvgTime.pather 97 DailyAvgTime.pather 98 DailyAvgTime.pather 99 DailyAvgTime.pather 90 DailyAvgTime.pather 910 DailyAvgTime.pather 910 DailyAvgTime.gather 9110 DailyAvgTime.gather 9110 DailyAvgTime.gather 9111 DailyAvgTime.gather <td></td>	
a) Daily AvgAc_230-300mg 90 Daily AvgAc_300-400mg 91 Daily AvgAc_30-80mg 92 Daily AvgAc_20-80mg 93 Daily AvgAc_200-80mg 94 Daily AvgAc_201ght 95 Daily AvgAc_201ght 96 Daily AvgTimelBout_light 97 Daily AvgTimelBout_light 98 Daily AvgTime_10-040mg 99 Daily AvgTime_100-100mg 90 Daily AvgTime_100-200mg 910 Daily AvgTime_100-200mg 910 Daily AvgTime_200-240mg 910 Daily AvgTime_200-240mg 910 Daily AvgTime_200-240mg 910 Daily AvgTime_30-400mg 910 Daily AvgTime_30-400mg 910 Daily AvgTime_30-300mg 910 Daily AvgTime_30-200mg 910 Daily AvgTime_30-20-300mg 911 Daily AvgTime_10+100mg 911 Daily AvgTime_10+100mg 911 Daily AvgTime_10+100mg 911 Daily AvgTime_10+100mg 911 Daily PortimelBout_light 911 Daily PortimelBout_light <	
90 DailyAvgAcc. 360-400mg 91 DailyAvgAcc. 40-80mg 92 DailyAvgAcc. 40-80mg 93 DailyAvgAcc. MVPA 94 DailyAvgAcc. [ight 95 DailyAvgTime.Bout_MVPA 96 DailyAvgTime.Bout_MVPA 97 DailyAvgTime.Bout_MVPA 98 DailyAvgTime.Bout_sclentary 99 DailyAvgTime.Light 90 DailyAvgTime.10-040mg 101 DailyAvgTime.200-240mg 102 DailyAvgTime.200-240mg 103 DailyAvgTime.200-240mg 104 DailyAvgTime.230-320mg 105 DailyAvgTime.200-240mg 106 DailyAvgTime.320-360mg 107 DailyAvgTime.30-400mg 108 DailyAvgTime.30-400mg 109 DailyAvgTime.30-400mg 108 DailyAvgTime.30-400mg 109 DailyAvgTime.30-400mg 101 DailyAvgTime.30-400mg 102 DailyAvgTime.40-80mg 103 DailyAvgTime.40-80mg 110 DailyAvgTime.dout.MVPA 120 DailyAvgTime.30-400mg <t< td=""><td></td></t<>	
91 DailyAvgAcc_40-80mg 92 DailyAvgAcc_40-80mg 93 DailyAvgAcc_MVPA 94 DailyAvgAcc_light 95 DailyAvgAcc_light 96 DailyAvgTimeInBout_MVPA watch 97 DailyAvgTimeInBout_Gedentary 98 DailyAvgTime.0-40mg 99 DailyAvgTime_0-40mg 90 DailyAvgTime_1000_sedentary 91 DailyAvgTime_200-240mg 101 DailyAvgTime_200-240mg 102 DailyAvgTime_200-240mg 103 DailyAvgTime_200-240mg 104 DailyAvgTime_200-240mg 105 DailyAvgTime_200-240mg 106 DailyAvgTime_200-240mg 107 DailyAvgTime_200-240mg 108 DailyAvgTime_200-240mg 109 DailyAvgTime_40-80mg 108 DailyAvgTime_40-80mg 109 DailyAvgTime_1ight 101 DailyAvgTime_sedentary 110 DailyAvgTime_sedentary 111 DailyAvgTime_1ight 112 DailyPetTimeInBout_ight 113 DailyPet_20-200-240mg	
Juily AugAcc_MVPA Juily AugAcc_MVPA Juily AugAcc_Sedentary Juily AugAcc_sedentary Juily AugTimeInBout_light Daily AvgTimeInBout_light Daily AvgTimeInBout_sedentary Daily AvgTime.160-200mg Daily AvgTime.20-160mg Daily AvgTime.20-240mg Daily AvgTime.30-400mg Daily AvgTime.30-400mg Daily AvgTime.light Daily AvgTime.light Daily PerTimeInBout_MVPA Daily PerTimeInBout_Sedentary Daily PerTimeInBout_sedentary Daily Pert_120-160mg Daily Pert_120-160mg Daily Pert_120-160mg Daily Pert_120-240mg Daily Pert_240-380mg Daily Pert_240-280mg Daily Pert_240-280	
94DailyAvgAcc_lightDailyAvgTime_InBout_MVPAwatchTVDA95DailyAvgTimeInBout_lightDailyAvgTimeInBout_lightDailyAvgTimesent in MVPA [m.s.²]97DailyAvgTimeInBout_lightDaily avg. time spent in light activity bouts [mins]98DailyAvgTime_10-40mgDailyAvgTime_20-240mg100DailyAvgTime_240-280mgDailyAvgTime_240-280mg103DailyAvgTime_240-280mgDailyAvgTime_280-320mg104DailyAvgTime_240-280mgDailyAvgTime_80-120mg105DailyAvgTime_20-240mgDailyAvgTime_30-360mg106DailyAvgTime_240-280mgDailyAvgTime_80-120mg107DailyAvgTime_80-120mgDailyavg. time spent in MVPA [mins]110DailyAvgTime_80-120mgDailyavg. time spent in MVPA [mins]111DailyAvgTime_lightDaily avg. time spent in light activity [mins]112DailyAvgTime_lout_lightDaily avg. time spent in MVPA [mins]113DailyPetTimeInBout_sedentaryDaily percent of time spent in light activity [mins]114DailyPet_120-160mgDaily percent of time spent sedentary [mins]118DailyPet_240-280mgDaily percent of time in each 40 mg range [%]121DailyPet_360-400mgDaily percent of time in each 40 mg range [%]122DailyPet_30-400mgDaily percent of time in each 40 mg range [%]124DailyPet_240-280mgDaily percent of time in each 40 mg range [%]124DailyPet_360-400mgDailyPet_30-400mg125DailyPet_360-400mgDailyPet_40-80mg124DailyPet_30-400m	
95 DailyAvgTime_IBout_MVPA watch TVDA 96 DailyAvgTimeIBout_MVPA watch TVDA 97 DailyAvgTimeIBout_sedentary Daily avg. time spent in MVPA bouts [mins] 98 DailyAvgTime.nBout_sedentary Daily avg. time spent in light activity bouts [mins] 99 DailyAvgTime_10-40mg Daily avg. time spent in light activity bouts [mins] 100 DailyAvgTime_10-40mg Daily avg. time spent in light activity bouts [mins] 101 DailyAvgTime_10-200mg DailyAvgTime_30-360mg 102 DailyAvgTime_30-360mg DailyAvgTime_40-80mg 103 DailyAvgTime_30-360mg DailyAvgTime_me_80-120mg 104 DailyAvgTime_d-0-80mg DailyAvgTime_me_80-120mg 105 DailyAvgTime_work Daily avg. time spent in light activity [mins] 106 DailyAvgTime_light Daily avg. time spent in light activity [mins] 103 DailyAvgTimelnBout_light Daily avg. time spent in sedentary [mins] 114 DailyPetTimelnBout_sedentary Daily percent of time spent in light activity [mins] 113 DailyPet_200-240mg Daily percent of time in each 40 mg range [%] 114 DailyPet_200-240mg Daily percent of time	
30Daily Avg Time Jinden Bout_lightWatchI VDA97Daily Avg Time Spent in light activity bouts [in98Daily Avg Time_0-40mg100Daily Avg Time_120-160mg101Daily Avg Time_240-280mg102Daily Avg Time_240-280mg103Daily Avg Time_30-360mg104Daily Avg Time_30-360mg105Daily Avg Time_30-360mg106Daily Avg Time_30-360mg107Daily Avg Time_80-120mg108Daily Avg Time_80-120mg109Daily Avg Time_80-120mg110Daily Avg Time_80-120mg111Daily Avg Time_80-120mg112Daily PortTimeInBout_light113Daily PctTimeInBout_light114DailyPct_160-00mg115DailyPct_120-160mg116DailyPct_200-240mg118DailyPct_230-320mg120DailyPct_230-320mg121DailyPct_360-400mg122DailyPct_30-300mg123DailyPct_40-80mg124DailyPct_40-80mg124DailyPct_MVPA125DailyPct_MVPA126DailyPct_MVPA127DailyPct_MVPA128DailyPct_MVPA129DailyPct_MVPA124DailyPct_MVPA125DailyPc	
98Daily AvgTime.InBout_sedentaryDaily avg. time spent sedentary bouts [mins]99Daily AvgTime_120-160mg101Daily AvgTime_160-200mg102Daily AvgTime_240-280mg103Daily AvgTime_280-320mg104Daily AvgTime_360-400mg105Daily AvgTime_80-120mg106Daily AvgTime_100-280mg107Daily AvgTime_360-400mg108Daily AvgTime_80-120mg109Daily AvgTime_s0-120mg109Daily AvgTime_s0-120mg101Daily AvgTime_s0-120mg102Daily AvgTime_s0-120mg103Daily AvgTime_s0-120mg104Daily AvgTime_s0-120mg105Daily AvgTime_s0-120mg106Daily AvgTime_s0-120mg110Daily AvgTime_s0-120mg111Daily AvgTime_s0-120mg112Daily PctTimeInBout_MVPA113DailyPctTimeInBout_g0114DailyPctTimeInBout_g0115DailyPct_120-160mg116DailyPct_20-280mg120DailyPct_20-280mg121DailyPct_280-320mg122DailyPct_280-320mg123DailyPct_40-80mg124DailyPct_40-80mg123DailyPct_40-80mg124DailyPct_40-80mg125DailyPct_MVPA126DailyPct_MVPA127DailyPct_40-80mg128DailyPct_MVPA129DailyPct_40-80mg124DailyPct_MVPA125DailyPct_MVPA126DailyPct_MVPA	mins]
99DailyAvgTime_0.40mg100DailyAvgTime_120-160mg101DailyAvgTime_160-200mg102DailyAvgTime_240-280mg103DailyAvgTime_280-320mg104DailyAvgTime_320-360mg105DailyAvgTime_30-400mg106DailyAvgTime_80-120mg107DailyAvgTime_80-120mg108DailyAvgTime_sedentary109DailyAvgTime_sedentary111DailyAvgTime_sedentary112DailyPerTimeInBout_MVPA113DailyPerTimeInBout_gedentary114DailyPerTimeInBout_gedentary115DailyPerTimeInBout_sedentary116DailyPet_100-200mg117DailyPet_100-200mg118DailyPet_100-200mg119DailyPet_240-280mg110DailyPet_240-280mg111DailyPet_240-280mg112DailyPet_240-280mg113DailyPet_40-80mg114DailyPet_40-80mg115DailyPet_40-80mg120DailyPet_40-80mg121DailyPet_40-80mg122DailyPet_40-80mg123DailyPet_40-80mg124DailyPet_40-80mg125DailyPet_MVPA126DailyPet_MVPA127DailyPet_MVPA128DailyPet_MVPA129DailyPet_MVPA121DailyPet_MVPA122DailyPet_MVPA123DailyPet_MVPA124DailyPet_MVPA125DailyPet_MVPA126DailyPet_MVPA127<	
100DailyAvgTime_160-200mg101DailyAvgTime_200-240mg102DailyAvgTime_240-280mg103DailyAvgTime_280-320mg104DailyAvgTime_320-360mg105DailyAvgTime_30-400mg106DailyAvgTime_80-120mg107DailyAvgTime_80-120mg108DailyAvgTime_sedentary109DailyAvgTime_sedentary111DailyAvgTime_sedentary112DailyPerTimeInBout_MVPA113DailyPerTimeInBout_gight114DailyPerTimeInBout_gedentary115DailyPerTimeInBout_sedentary116DailyPet_100-200mg117DailyPet_160-200mg118DailyPet_160-200mg119DailyPet_240-280mg120DailyPet_240-280mg121DailyPet_30-360mg122DailyPet_40-80mg123DailyPet_40-80mg124DailyPet_40-120mg125DailyPet_MVPA125DailyPet_MVPA	
102DailyAvgTime_200-240mg103DailyAvgTime_240-280mg104DailyAvgTime_280-320mg105DailyAvgTime_320-360mg106DailyAvgTime_30-400mg107DailyAvgTime_80-400mg108DailyAvgTime_80-120mg109DailyAvgTime_light101DailyAvgTime_sedentary111DailyPetTimelnBout_MVPA112DailyPetTimelnBout_MVPA113DailyPetTimelnBout_gedentary114DailyPetTimelnBout_gedentary115DailyPet_120-160mg116DailyPet_120-240mg117DailyPet_240-280mg118DailyPet_240-280mg120DailyPet_240-280mg121DailyPet_40-80mg122DailyPet_40-80mg123DailyPet_40-80mg124DailyPet_40-80mg125DailyPet_MVPA125DailyPet_MVPA	
103DailyAvgTime_240-280mg104DailyAvgTime_280-320mg105DailyAvgTime_320-360mg106DailyAvgTime_30-400mg107DailyAvgTime_80-120mg108DailyAvgTime_MVPA109DailyAvgTime_light101DailyAvgTime_sedentary111DailyPetTimeInBout_MVPA112DailyPetTimeInBout_gett113DailyPetTimeInBout_ight114DailyPetTimeInBout_ight115DailyPet_100-200mg116DailyPet_100-200mg117DailyPet_200-240mg118DailyPet_280-320mg120DailyPet_280-320mg121DailyPet_30-400mg122DailyPet_40-80mg123DailyPet_40-80mg124DailyPet_40-200mg125DailyPet_MVPA125DailyPet_MVPA	
105 DailyAvgTime_320-360mg 106 DailyAvgTime_30-400mg 107 DailyAvgTime_40-80mg 108 DailyAvgTime_80-120mg 109 DailyAvgTime_light 101 DailyAvgTime_light 102 DailyAvgTime_sedentary 103 DailyPetTimeInBout_MVPA 114 DailyPetTimeInBout_MVPA 115 DailyPetTimeInBout_sedentary 116 DailyPet_100-40mg 117 DailyPet_100-200mg 118 DailyPet_240-280mg 120 DailyPet_320-360mg 121 DailyPet_40-80mg 122 DailyPet_40-80mg 123 DailyPet_40-80mg 124 DailyPet_40-80mg 125 DailyPet_MVPA	nins]
106DailyAvgTime_360-400mg107DailyAvgTime_40-80mg108DailyAvgTime_80-120mg109DailyAvgTime_MVPA101DailyAvgTime_light112DailyAvgTime_sedentary113DailyPetTimeInBout_MVPA114DailyPetTimeInBout_MVPA115DailyPet_0-40mg116DailyPet_120-160mg117DailyPet_240-280mg118DailyPet_230-360mg120DailyPet_320-360mg121DailyPet_40-80mg122DailyPet_40-80mg123DailyPet_40-80mg124DailyPet_MVPA125DailyPet_MVPA126DailyPet_MVPA127DailyPet_MVPA128DailyPet_MVPA129DailyPet_MVPA120DailyPet_40-80mg121DailyPet_40-80mg122DailyPet_40-80mg123DailyPet_MVPA124DailyPet_MVPA125DailyPet_MVPA126DailyPet_MVPA127DailyPet_MVPA128DailyPet_MVPA129DailyPet_MVPA120DailyPet_MVPA121DailyPet_MVPA122DailyPet_MVPA123DailyPet_MVPA124DailyPet_MVPA125DailyPet_MVPA126DailyPet_MVPA127DailyPet_MVPA128DailyPet_MVPA129DailyPet_MVPA129DailyPet_MVPA129DailyPet_MVPA129DailyPet_MVPA<	
107 DailyAvgTime_40-80mg 108 DailyAvgTime_80-120mg 109 DailyAvgTime_MVPA 110 DailyAvgTime_light 111 DailyAvgTime_sedentary 112 DailyPetTimeInBout_MVPA 113 DailyPetTimeInBout_MVPA 114 DailyPetTimeInBout_sedentary 115 DailyPet_10-40mg 116 DailyPet_10-40mg 117 DailyPet_200-240mg 118 DailyPet_240-280mg 120 DailyPet_320-360mg 121 DailyPet_30-400mg 122 DailyPet_40-80mg 123 DailyPet_80-120mg 124 DailyPet_80-120mg 125 DailyPet_MVPA	
109DailyAvgTime_MVPA109DailyAvgTime_MVPA110DailyAvgTime_light111DailyAvgTime_sedentary112DailyPetTimeInBout_MVPA113DailyPetTimeInBout_MVPA114DailyPetTimeInBout_gedentary115DailyPet_10-40mg116DailyPet_120-160mg117DailyPet_240-280mg118DailyPet_320-360mg120DailyPet_320-360mg121DailyPet_40-80mg122DailyPet_40-80mg124DailyPet_MVPA125DailyPet_MVPA	
110 DailyAvgTime_light Daily avg. time spent in light activity [mins] 111 DailyAvgTime_sedentary Daily avg. time spent in sedentary [mins] 112 DailyPctTimeInBout_MVPA Daily percent of time spent in sedentary [mins] 113 DailyPctTimeInBout_light Daily percent of time spent in light botts [%] 114 DailyPct_O-40mg Daily percent of time spent in light botts [%] 115 DailyPct_120-160mg Daily percent of time spent sedentary bouts [%] 116 DailyPct_120-160mg DailyPct_200-240mg 117 DailyPct_200-240mg DailyPct_200-240mg 118 DailyPct_230-320mg Daily percent of time in each 40 mg range [%] 120 DailyPct_300-400mg DailyPct_40-80mg 121 DailyPct_40-80mg DailyPct_40-80mg 122 DailyPct_40-80mg DailyPct_MVPA 123 DailyPct_MVPA DailyPct_for time in MVPA [%]	
111 DailyAvg1ime_sedentary Daily avg. time spent in sedentary [mins] 112 DailyPetTimeInBout_MVPA Daily percent of time spent in MVPA bouts [%] 113 DailyPetTimeInBout_ight Daily percent of time spent in light bouts [%] 114 DailyPet_120-160mg Daily percent of time spent in light bouts [%] 115 DailyPet_120-160mg Daily percent of time spent in sedentary bouts [%] 116 DailyPet_200-240mg DailyPet_200-240mg 117 DailyPet_200-240mg DailyPet_200-240mg 119 DailyPet_230-320mg Daily percent of time in each 40 mg range [%] 121 DailyPet_300-400mg DailyPet_40-80mg 122 DailyPet_40-80mg DailyPet_40-80mg 123 DailyPet_80-120mg DailyPet_MVPA 124 DailyPet_MVPA Daily percent of time in MVPA [%]	
113 DailyPctTimeInBout_light Daily percent of time spent in light bouts [%] 114 DailyPctTimeInBout_sedentary Daily percent of time spent in light bouts [%] 115 DailyPct_0-40mg Daily percent of time spent in light bouts [%] 116 DailyPct_120-160mg Daily percent of time spent sedentary bouts [%] 117 DailyPct_200-240mg DailyPct_200-240mg 118 DailyPct_200-240mg DailyPct_200-240mg 119 DailyPct_200-240mg DailyPct_200-240mg 120 DailyPct_200-300mg Daily percent of time in each 40 mg range [%] 121 DailyPct_300-400mg DailyPct_40-80mg 122 DailyPct_300-400mg DailyPct_40-80mg 124 DailyPct_80-120mg Daily percent of time in MVPA [%]	%]
114 DailyPctTimeInBout_sedentary Daily percent of time spent sedentary bouts [%] 115 DailyPct_0-40mg Daily percent of time spent sedentary bouts [%] 116 DailyPct_120-160mg Daily Det_160-200mg 117 DailyPct_160-200mg DailyPct_200-240mg 118 DailyPct_200-240mg DailyPct_200-240mg 119 DailyPct_200-320mg Daily percent of time in each 40 mg range [%] 120 DailyPct_320-360mg Daily percent of time in each 40 mg range [%] 121 DailyPct_300-400mg DailyPct_40-80mg 122 DailyPct_40-80mg DailyPct_80-120mg 124 DailyPct_80-120mg Daily percent of time in MVPA [%]]
115 DailyPet_0-40mg 116 DailyPet_120-160mg 117 DailyPet_160-200mg 118 DailyPet_20-240mg 119 DailyPet_240-280mg 120 DailyPet_280-320mg 121 DailyPet_320-360mg 122 DailyPet_300-400mg 123 DailyPet_40-80mg 124 DailyPet_80-120mg 125 DailyPet_MVPA 125 DailyPet_MVPA	%]
117 DailyPct_160-200mg 118 DailyPct_200-240mg 119 DailyPct_240-280mg 120 DailyPct_280-320mg 121 DailyPct_320-360mg 122 DailyPct_360-400mg 123 DailyPct_40-80mg 124 DailyPct_80-120mg 125 DailyPct_MVPA 125 DailyPct_MVPA	
118 DailyPct_200-240mg 119 DailyPct_240-280mg 120 DailyPct_280-320mg 121 DailyPct_320-360mg 122 DailyPct_300-400mg 123 DailyPct_40-80mg 124 DailyPct_80-120mg 125 DailyPct_MVPA 125 DailyPct_MVPA	
120 DailyPct_280-320mg 121 DailyPct_280-320mg 122 DailyPct_300-360mg 123 DailyPct_360-400mg 124 DailyPct_080-120mg 125 DailyPct_MVPA 125 DailyPct_MVPA	
121 DailyPct_320-360mg 122 DailyPct_360-400mg 123 DailyPct_40-80mg 124 DailyPct_80-120mg 125 DailyPct_MVPA 125 DailyPct_MVPA	6]
122 DailyPct_360-400mg 123 DailyPct_40-80mg 124 DailyPct_80-120mg 125 DailyPct_MVPA Daily percent of time in MVPA [%]	
125 DailyPct_80-120mg 125 DailyPct_MVPA 125 DailyPct_MVPA 126 Daily percent of time in MVPA [%]	
125 DailyPct_MVPA Daily percent of time in MVPA [%]	
126 DailyPct_Walking 127 DailyPct_light DailyPct_light Daily percent of time in light activity [%]	
128 DailyPct_sedentary DailyPct_sedentary	
129 Median AccInBout_MVPA Median acc. while in MVPA	
130 Median AccInBout_light Median acc. while in light activity 131 Median AccInBout sedentary Median acc. while in cedentary	
132 Median AccMag Median acc. vector magnitude	
133 PctTime_0-40mg	
134 PctTime_120-160mg 135 PctTime_160-200mg watch TVDA	
136 PctTime_200-240mg Watch 17DA	
137 PctTime_240-280mg Study percent of time in each 40 mg range [%	%]
138 Pct1ime_280-320mg 139 PctTime 320-360mg	-

Table E.I. Description of the sensor-based features extracted in the weaRAble-PRO study

Continued on Next Page...

	Feature	Source	Dom.	Description
140	PctTime_360-400mg			
141	PctTime_40-80mg			
142	PctTime_80-120mg			
143	PctTime_MVPA			Study percent of time in MVPA [%]
144	PctTime_light			Study percent of time spent in light activity [%]
145	PctTime_sedentary			Study percent of time in sedentary activity [%]
146	PctTime_walking			Study percent of time spent walking [%]
147	StdAccMag			Standard deviation in acc. vector magnitude [m.s ⁻²]
148	nPeriods_Walking_120to600			Number of continuous walking periods with duration 2 to 10 minutes
149	nPeriods_Walking_1800toinf			(with up to 30-second rest period) [count] Number of continuous walking periods with duration >30 minutes (with up to 1-minute rest period) [count]
150	nPeriods_Walking_600to1800			Number of continuous walking periods with duration 10 to 30 minutes (with up to 1-minute rest period) [count]
151	WALK_numberOfSteps			The daily number of steps [count]
152	WALK_stepFrequency			Cadence [steps/min]
153	WALK_stepLength	phone	WLK	Step length [cm]
154	WALK_stepTime			Step time [sec]
155	WALK_stepVelocity			Walk velocity [m/s]
156	WRIST_ROM_global			Range of motion (ROM) [deg]
157	WRIST_ROM_max			Maximum ROM [deg]
158	WRIST_ROM_mean			Mean ROM [deg]
159	WRIST_ROM_median			Median ROM [deg]
160	WRIST_ROM_min	phone	WRT	Minimum ROM [deg]
161	WRIST_angvel_max			Maximum angular velocity [deg/s]
162	WRIST_angvel_mean			Mean angular velocity [deg/s]
163	WRIST_angvel_median			Median angular velocity [deg/s]
164	WRIST_angvel_min			Minimum angular velocity [deg/s]

Table E.I. Description of the sensor-based features extracted in the weaRAble-PRO study

F Methodology: Machine learning analysis for characterising RA

F.1 Feature Pre-processing

Missing data extracted from the GTs and passive monitoring was imputed using a carry-last value forward per participant⁷⁷. In cases where the last value was missing, mean imputation was used instead.

Features were assessed for non-normality by visual inspection. Those non-normal features were transformed using box-cox transformations⁷⁸.

Features were normalised using the z-score to have unit variance using their respective mean μ and standard deviation $\sqrt{\sigma}$.

F.2 Classification and Regression Models

F.2.1 Linear Regression

A regression model explicitly describes a relationship between predictor(s) $\mathbf{X} \in \mathbb{R}^{N \times P}$ and continuous response variables $\mathbf{y} \in \mathbb{R}^N$, the most basic of which is *linear regression* (LR)⁵⁵. For an *i*th observation row of \mathbf{X} , $\mathbf{x} \equiv \mathbf{x}_i \in \mathbb{R}^{1 \times P}$:

$$\hat{y} = w_0 + \sum_{j=1}^{P} w_j x_j + \varepsilon$$
(F.1)
$$= \mathbf{w}^\top \mathbf{x} + b$$
(F.2)

where w_j values denote the slope (weights, or regression coefficients) of the x_j features; w_0 is the intercept term; and ε denote the residual (model) errors term, which are assumed to be normally distributed with constant variance, $\varepsilon \sim \mathcal{N}(0, \sigma^2)^{48, 55}$. Often a linear model is described in vector notation (F.2), $\mathbf{w} = [w_0, w_1, w_2, ..., w_P]$, where w_0 is denoted as the "bias", *b* term, and the ε -term is often omitted.

F.2.2 Logistic Regression

Generalised linear models (GLMs) are extensions of linear regression models that can have non-linear outputs⁷⁹. GLMs utilise canonical link functions, ϕ , to transform the outputs of a linear regression: $\varphi = \mathbf{w}^{\top} \mathbf{x}$ to another distribution, such as with a logistic $\phi = \sigma(\varphi)$ link function (or inversely the logit, representing the log-odds) which will be used to form Logistic Regression for binary classification tasks^{48,80}, in this case ϕ is sigmoidal and is bounded between [0,1], therefore the output of σ can be interpreted as the probability of y = 1:

$$p(\mathbf{x}; \mathbf{w}) = \frac{1}{1 + e^{-(\mathbf{w}^{\top} \mathbf{x})}}$$
(F.3)

A threshold can be applied to the probabilistic output p to determine a classification prediction \hat{y} for a Logistic Regression model; threshold values are typically chosen as 0.5, but this can be altered based on the use case.

Regularisation as Feature Selection Many statistical and machine learning models can easily overfit to the training data, resulting in poorer estimations, models that are not generalisable or too complex. This is likely in the weaRAble-PRO dataset where we have $P \gg N$ problem: the number of predictors *P* is much larger than the number of observations *N* (participants), a case where standard models fail.

Regularisation can be introduced to mitigate against the $p \gg n$ problem. For example large coefficient values in a regression can be penalised by adding a regularisation term to a loss function, or through reducing the number of parameters or features used in a model. The most common regularisers use the ℓ_p -norm defined by^{48,81}:

$$||\mathbf{x}||_{p} = \left(\sum_{i=1}^{N} |x_{i}|^{p}\right)^{1/p}$$
(F.4)

for any $\mathbf{x} \in \mathbb{R}^{N \times 1}$, where the real number $p \ge 1$ defines the ℓ_p space. In this work, we experimented with a number of regularisation techniques in order to perform the classification tasks introduced in section 2.2.

F.2.3 Least Absolute Shrinkage and Selection Operator

The Least Absolute Shrinkage and Selection Operator (LASSO)^{81,82} is a technique that conversely solves the ℓ_1 -penalised sum of squares (F.4) in a linear regression such that:

$$\hat{\boldsymbol{w}} = \operatorname*{arg\,min}_{\mathbf{w}} \left\{ ||\boldsymbol{y} - \mathbf{w}^{\top} \mathbf{x}||_{2}^{2} + \lambda ||\mathbf{w}||_{1} \right\}$$
(F.5)

This is equivalent to minimising the sum of squares with a constraint of the form: $||w||_1 = \sum_{j=1}^{N} |w_j| \le t$. Because of the form of the ℓ_1 -penalty, LASSO both shrinks coefficients but also encourages sparsity in a model's parameters and thus inherently forms

feature selection, shrinking non-important features to zero. The LASSO can also be extended to perform feature selection for classification by substituting a canonical link function (such as the logistic $\phi = \sigma(\phi)$) and following the same procedure outlined in equation F.3, essentially performing regularised-logistic regression (denoted LR-lasso in this work)⁸¹.

F.2.4 Ridge Regression

Ridge Regression (Tikhonov Regularisaion)^{48,81} is a technique which utilises the ℓ_2 -norm to impose a penalty on the size of the coefficients in a linear regression (F.4) such that:

$$\hat{\boldsymbol{w}} = \operatorname*{arg\,min}_{\mathbf{w}} \left\{ ||\boldsymbol{y} - \mathbf{w}^{\top} \mathbf{x}||_{2}^{2} + \lambda ||\mathbf{w}||_{2} \right\}$$
(F.6)

This is equivalent to minimising the sum of squares with a constraint of the form: $||w||_2 = \sum_j^N w_j^2 \le t$ (where *t* controls the amounts of shrinkage; there is an exact relationship between *t* and corresponding λ (denoted LR-ridge in this work).

F.2.5 Elastic Net

Elastic net linearly combines the ℓ_1 and ℓ_2 penalties of the lasso and ridge methods such that:

$$\hat{\boldsymbol{w}} = \underset{\mathbf{w}}{\arg\min} \left\{ ||\boldsymbol{y} - \mathbf{w}^{\top} \mathbf{x}||_{2}^{2} + (1 - \alpha)\lambda ||\mathbf{w}||_{1} + \alpha\lambda ||\mathbf{w}||_{2} \right\}$$
(F.7)

given non-negative values λ , and α that is strictly between 0 and 1, which determines the trade off between ℓ_1 and ℓ_2 regularisation (denoted LR-elastic-net in this work).

F.2.6 Sparse-Group LASSO

The sparse-group lasso is an extension of the lasso that promotes both group sparsity and within group parameter-wise sparsity, through a group lasso penalty and the lasso penalty:

$$\hat{\boldsymbol{w}} = \underset{\mathbf{w}}{\operatorname{arg\,min}} \left\{ ||\boldsymbol{y} - \mathbf{w}^{\top} \mathbf{x}||_{2}^{2} + \lambda_{1} ||\mathbf{w}||_{1} + \lambda_{2} \sum_{l=1}^{M} \sqrt{p_{l}} ||\mathbf{w}^{(l)}||_{2} \right\}$$
(F.8)

where λ_1 is the parameter-wise regularisation penalty and λ_2 is the group-wise regularisation penalty. The data **x** contains sub-grouping, such that $\mathbf{x}^{(l)}$, denoting the features in group l, and corresponding learned weights **w** containing contains sub-grouping $\mathbf{w}^{(l)}$; p_l is the length of $\mathbf{w}^{(l)}$ (i.e., the number of features in each group) and M are the total number groups. Therefore, the sparse group lasso penalty will yield a sparse set of groups and also a sparse set of covariates in each selected group. As denoted in⁵⁰, we use the term "groupwise sparsity" to refer to the number of groups with at least one nonzero coefficient, and "within group sparsity" to refer to the number of nonzero coefficients within each nonzero group (denoted LR-SG-lasso in this work).

F.2.7 Random Forest

Classification and Regression Trees (CART), specifically Random Forests (RF), are a multi-functional, non-linear method capable of performing regression, classification and feature selection⁵¹. Unlike the linear filter-based methods of feature selection, for example, lasso, RFs incorporate non-linear feature selection as part of the model methodology.

Random Forests consist of a large ensemble of decision trees arranged in a hierarchical structure. To build an individual tree, we recursively descend through the hierarchy, performing binary splits (decisions) at each level in the structure (a node, *j*) using a single feature $x_j \in \mathcal{X}^p$ based on a threshold value (splitting criterion) s_j , sub-partitioning the feature space \mathcal{X}_j at each node. A tree is typically expanded until all leaves are pure (i.e each partition \mathcal{X}_j represents only one class) or until all leaves contain less than the minimum number of samples in a partition \mathcal{X}_j required to split a node. An individual tree selects m < N random subset of observations (with replacement), and each node considers a random subset of *p* features for each split. Once all trees have been grown, each of the 'weaker' decisions are aggregated (or *ensembled*) creating a robust final prediction. For example, consider an ensemble of two trees. The prediction scores of each individual tree are summed to obtain the final prediction:

$$\hat{y}_i = \sum_{k=1}^K f_k(x_i), f_k \in \mathcal{F}$$
(F.9)

where *K* is the number of trees, f_k is a function in the functional space \mathcal{F} , and \mathcal{F} is the set of all possible CARTs. The objective function optimised is then given by:

$$E(\boldsymbol{\theta}) = \sum_{i}^{n} \mathcal{L}(y_i, \hat{y}_i) + \sum_{k=1}^{K} w(f_k)$$
(F.10)

where $w(f_k)$ is the complexity of the tree f_k .

Classification predictions are deduced as the majority class label of the observations present in each final partition \mathcal{X}_j , whereas for continuous prediction (i.e., regression) the mean of the (continuous) responses would be calculated instead. To determine the optimal split criterion s_j for each x_j to create each \mathcal{X}_j we evaluate the *Gini* importance, which quantifies the average gain of purity (i.e., the presence of one class) caused by splits of a given variable. For regression the mean decrease in mean square error (MSE) is assessed instead.

F.2.8 Extreme Gradient Boosted Trees (XGB)

Extreme Gradient Boosted Trees (XGB) iterates on the CART through a regularising gradient boosting framework—the difference being in how trees are built and combined. Rather than bagging, like in CART, specifically RF, gradient boosting improves a single weak model by combining it with a number of other weak models in order to generate a collectively stronger model⁵². This boosting is formed as additive strategy during training where a gradient descent algorithm is used to minimise (or maximise) an objective function $E(\theta)$ for each new tree that is added at each time step, *t*:

$$E(\theta)^{t} = \sum_{i=1}^{n} \mathcal{L}(y_{i}, \hat{y}_{i}^{(t)}) + \sum_{i=1}^{t} w(f_{i})$$
(F.11)

where \hat{y}_i^t is the prediction value at step *t*; and f_i are the parameters of a tree, i.e., the tree structure and the leaf scores that are needed to be learned; and $w(f_i)$ is the complexity of the tree. Therefore \hat{y}_i^t is determined by the prediction of the previous tree at t-1, $\hat{y}_i^{(t-1)}$:

$$\hat{y}_i^{(t)} = \sum_{k=1}^{t} f_k(x_i) = \hat{y}_i^{(t-1)} + f_t(x_i)$$
(F.12)

To optimise the error function, XGBoost computes gradients (Jacobian) and Hessians of the error, denoted below as:

$$g_i = \partial_{\hat{y}_i^{(t-1)}} \mathcal{L}(y_i, \hat{y}_i^{(t-1)})$$
(F.13)

$$h_i = \partial_{\hat{y}_i^{(t-1)}}^2 \mathcal{L}(y_i, \hat{y}_i^{(t-1)})$$
(F.14)

which can be subbed into the objective function in equation (F.11), yielding:

$$E(\theta)^{t} = \sum_{i=1}^{n} \left[g_{i} f_{t}(x_{i}) + \frac{1}{2} h_{i} f_{i}^{2}(x_{i}) \right] + w(f_{t})$$
(F.15)

Next, to define the complexity of the tree w(f), we first refine the function for the definition of a tree f(x) as:

$$f_t(x) = w_{q(x)}, w \in \mathcal{R}^T, q : \mathcal{R}^d \to \{1, 2, \dots, T\}$$
(F.16)

where w is the vector of scores on leaves, q is a function assigning each data point to the corresponding leaf, and T is the number of leaves. Then, in XGBoost, we can define the complexity over the number of leaves, T:

$$w(f) = \gamma T + \frac{1}{2}\lambda \sum_{j=1}^{T} w_j^2$$
(F.17)

where γ and λ are tune-able regularisation parameters. Subbing equation F.17 into equation F.15, becomes:

$$E(\boldsymbol{\theta})^{t} = \sum_{j=1}^{T} \left[G_{j} w_{j} + \frac{1}{2} (H_{j} + \lambda) w_{j}^{2} \right] + \gamma T$$
(F.18)

where $I_j = \{i | q(x_i) = j\}$; $G_j = \sum_{i \in I_j} g_i$; $H_j = \sum_{i \in I_j} h_i$. Following the convention in⁵², equation F.18 can be reformulated as solving over:

$$w_j^* = -\frac{G_j}{H_j + \lambda} \tag{F.19}$$

$$E(\theta)^* = -\frac{1}{2} \sum_{j=1}^T \frac{G_j^2}{H_j + \lambda} + \gamma T$$
(F.20)

Summarising, XGB iteratively trains an ensemble of shallow decision trees, with each iteration using the error residuals of the previous model to fit the next model. The objective $E(\theta)^*$ measures how good a tree structure q(x) is, optimising one level of the tree at a time, while regularising on the model complexity. For example, to split a leaf into two leaves, and the score will be:

$$Gain = \frac{1}{2} \left[\frac{G_L^2}{H_L + \lambda} + \frac{G_R^2}{H_R + \lambda} - \frac{(G_L + G_R)^2}{H_L + H_R + \lambda} \right] - \lambda$$
(F.21)

if the gain acquired by adding a branch is smaller than γ , it would not be added.

F.3 Parameter Tuning

Optimal model parameters were determined via randomised grid-search over internal 5-fold (subject-wise) CV with 500 iterations. In the case of XGB, early stopping was determined using roughly 10% of the training data, proportionally, as validation, with for 10 boosting rounds.

For regularised logistic regression models, a parameter search was determined over ℓ_1 and ℓ_2 regularisation terms on weights, $\lambda \in \{10^{-5}, ..., 10^{-1}, ..., 0, ..., 10^0, 10^1, ..., 10^5\}$; and the elastic-net mixing parameter $\alpha \in \{0, 0.1, ..., 1\}$;

RFs have relatively little hyperparameter tuning: the number of trees to build $p \in \{500, 1000, 1500\}$; the number of input variables chosen at each node $p \in \{\sqrt{P}, 2\sqrt{P}, \sqrt{P}/2\}$, where *P* are the number of features, as suggested in^{83,84}.

For XGB, a parameter search was determined over: the boosting learning rate $p \in \{0.01, 0.05, 0.1, 1\}$; number of boosting rounds $p \in \{100, 500, 1000, 1500\}$; maximum tree depth for base learners $p \in \{3, 4, 5, 8, 10\}$; the subsample ratio of the training instance (selecting random training instances with higher probability when the gradient and hessian are larger) $p \in \{0, 0.5, 1\}$; the subsample ratio of features when constructing each tree $p \in \{0.2, 0.6, 0.8, 1.0\}$; minimum sum of instance weight (hessian) needed in a child $p \in \{1, 5, 10, 50, 100\}$; the ℓ_1 and ℓ_1 regularisation term on weights, $p \in \{10^{-5}, ..., 10^{-1}, ..., 0, ..., 10^0, 10^1, ..., 10^5\}$