Check for updates

OPEN ACCESS

EDITED BY Michael Hudecek, University Hospital Würzburg, Germany

REVIEWED BY Alessandra Merlini, University of Turin, Italy Leila Amini, Center for Advanced Therapies (BeCAT), Germany

*CORRESPONDENCE Anthony E. Zamora, 🛙 azamora@mcw.edu

[†]These authors have contributed equally to this work and share first authorship

RECEIVED 09 October 2023 ACCEPTED 08 January 2024 PUBLISHED 01 February 2024

CITATION

Colina AS, Shah V, Shah RK, Kozlik T, Dash RK, Terhune S and Zamora AE (2024), Current advances in experimental and computational approaches to enhance CAR T cell manufacturing protocols and improve clinical efficacy. *Front. Mol. Med.* 4:1310002. doi: 10.3389/fmmed.2024.1310002

COPYRIGHT

© 2024 Colina, Shah, Shah, Kozlik, Dash, Terhune and Zamora. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.

Current advances in experimental and computational approaches to enhance CAR T cell manufacturing protocols and improve clinical efficacy

Alfredo S. Colina^{1†}, Viren Shah^{2†}, Ravi K. Shah³, Tanya Kozlik¹, Ranjan K. Dash², Scott Terhune^{1,2} and Anthony E. Zamora^{1,3}*

¹Department of Microbiology & Immunology, Medical College of Wisconsin, Milwaukee, WI, United States, ²Department of Biomedical Engineering, Medical College of Wisconsin and Marquette University, Milwaukee, WI, United States, ³Department of Medicine, Medical College of Wisconsin, Milwaukee, WI, United States

Since the FDA's approval of chimeric antigen receptor (CAR) T cells in 2017, significant improvements have been made in the design of chimeric antigen receptor constructs and in the manufacturing of CAR T cell therapies resulting in increased in vivo CAR T cell persistence and improved clinical outcome in certain hematological malignancies. Despite the remarkable clinical response seen in some patients, challenges remain in achieving durable long-term tumor-free survival, reducing therapy associated malignancies and toxicities, and expanding on the types of cancers that can be treated with this therapeutic modality. Careful analysis of the biological factors demarcating efficacious from suboptimal CAR T cell responses will be of paramount importance to address these shortcomings. With the ever-expanding toolbox of experimental approaches, single-cell technologies, and computational resources, there is renowned interest in discovering new ways to streamline the development and validation of new CAR T cell products. Better and more accurate prognostic and predictive models can be developed to help guide and inform clinical decision making by incorporating these approaches into translational and clinical workflows. In this review, we provide a brief overview of recent advancements in CAR T cell manufacturing and describe the strategies used to selectively expand specific phenotypic subsets. Additionally, we review experimental approaches to assess CAR T cell functionality and summarize current in silico methods which have the potential to improve CAR T cell manufacturing and predict clinical outcomes.

KEYWORDS

chimeric antigen receptor T cells, CAR T cell, manufacturing, functional assays, computational modeling


Introduction

Over the last few decades, efforts exploiting the immune system to target and eliminate malignant cells have grown in popularity and clinical utility. Recently, immunotherapeutic approaches have focused on redirecting T cells to preferentially target specific antigens expressed on cancerous cells. One of the most widely adopted approaches consists of using T cells engineered to express surface-bound synthetic chimeric antigen receptors (CARs).

10.3389/fmmed.2024.1310002

Since the first commercial authorization by the United States Food and Drug Administration (FDA) and European Medicine Agency (EMA) in 2017 and 2018, respectively, advancements in the genetic modification of CAR T cells have continued to progress at an astonishing rate (European Medicines Agency, 2018b; European Medicines Agency, 2018c; FDA, 2021a; FDA, 2021b). Each new generation of CAR T cell has focused on incorporating increasingly sophisticated engineering strategies into the CAR construct with the hope of driving greater persistence of CAR T cells in patients and better clinical efficacy. Yet, despite the constant evolution in CAR T cell design, a sizeable fraction of cancer patients will ultimately relapse after treatment and succumb to their disease and major barriers exist before comparable successes will be realized in the solid tumor setting (Park et al., 2018; Byrne et al., 2019; Locke et al., 2019; Xu et al., 2019; Melenhorst et al., 2022). As a result, recent efforts in CAR T manufacturing have focused on addressing lingering obstacles related to enhancing CAR T cell persistence, decreasing CAR-mediated toxicity, and maintaining cells in a highly functional state post-infusion.

While opportunities to continue improving CAR T cell performance by means of modifying or introducing new domains to the CAR construct remain, a recent switch in focus towards revisiting how we manufacture CAR T cells to address ongoing challenges is beginning to take center stage. Growing efforts aim to target one or more of the key steps along the CAR T cell manufacturing continuum (Figure 1) in hopes of generating a superior product with favorable biological attributes. Fueled by retrospective studies from the clinical application of CAR T cells, it is becoming clear that differences in the manufacturing process can directly impact the composition of the resulting product including the percent of cells expressing CARs, CD4:CD8 T cell ratios, and the relative frequency of T cells at specific stages of differentiation (Zurko et al., 2021). Importantly, the differentiation status of CAR T cells has been shown to impact post-infusion expansion and persistence due to inherent biological differences between the subsets (Wherry et al., 2003; Louis et al., 2011; Schmueck-Henneresse et al., 2017; Fraietta et al., 2018; Philip and Schietinger, 2019). Several aspects of the manufacturing process can be coopted to direct the cellular products toward optimal composition for improving therapeutic response (Gargett et al., 2019; Stock et al., 2019). These include methods of introducing genetic constructs into T lymphocytes, the selection of culture media additives to support cell growth, the activation and expansion of CAR T cells (including new activation agents and considering the duration of expansion), and the selection or depletion of specific subgroups of cells (Stemberger et al., 2012; Riddell et al., 2014; Schmueck-Henneresse et al., 2017; Xu et al., 2018; K et al., 2022). While we acknowledge that the discovery of novel CAR constructs and antigenic targets play an important role in the continuous success of this therapy (Holzinger and Abken, 2019), detailed review of these topics are outside the scope of this review and have been covered elsewhere (Sadelain et al., 2009; van der Stegen et al., 2015; Holzinger and Abken, 2019; Abken, 2021). The goal of this review is to cover specific steps of the CAR T cell manufacturing process and describe how changes in these steps can help address some of the ongoing challenges with CAR T cell efficacy. Furthermore, given that each step is manipulatable, we anticipate that by altering, optimizing, testing, and modeling these steps using




or depletion of specific cell types, T cell activation, gene delivery, CAR T cell expansion, and either cryopreserva-tion prior to infusion or infusion of fresh product.

in vitro assays and computational models, we will glean important biological insight that will ultimately impact clinical outcome.

To gain insight into how CAR T cells function once administered to patients with cancer, we must continue to develop better methods to assess product functionality (cytolytic ability, proliferation capacity, and cytokine production) and determine whether the phenotypic and functional trajectories seen in vitro recapitulate what occurs in vivo. By linking how specific perturbations of cellular activation alter the phenotypic and functional signatures of CAR T cells, we may better establish how these correlate with therapeutic efficacy. In this review, we will overview specific assays that can be used to measure cytokine secretion and cytotoxicity and how they play a vital role in determining pre-infusion CAR T cell function (FDA, 2011; FDA, 2023b). We also introduce how computational methods are being used as tools to help predict CAR T cell function and clinical efficacy by incorporating variables such as CAR surface expression, binding affinity between CAR and antigen, and product differentiation status. Overall, through identifying the relevant manufacturing variables and then detailing current experimental and computational methods that can inform the optimization of these variables, our aim is to contribute to the discourse on how to best modulate CAR T cell manufacturing to improve therapeutic responses.

Manufacturing approaches for CAR T-cells

The major steps for all current CAR T cell manufacturing processes include: i) collection of apheresis product, ii) enrichment



and/or depletion of specific cell types, iii) T cell activation, iv) gene delivery, v) CAR T cell expansion, vi) CAR T cell formulation and filing, and vii) either cryopreservation prior to infusion or infusion of fresh product. Many of the guiding principles for this process have been shaped by the development of tumor-infiltrating lymphocytes (TIL) as therapeutic agents, which have historically implemented a similar isolation, culturing, and expansion strategy (Rosenberg et al., 1990). As early considerations for infusion primarily relied on ensuring sufficient numbers of CARexpressing cells were obtained for infusion, the duration and formulations used during the manufacturing process often varied and resulted in inconsistent clinical responses across clinical trials (Wang and Rivière, 2016; Gajra et al., 2022). Given the lack of activation and expansion of first-generation CAR T cells, initial emphasis was placed on developing newer generations of constructs that would circumvent these

shortcomings (Boyiadzis et al., 2018; Guedan et al., 2018; Sermer and Brentjens, 2019). While second and third generation CAR constructs have displayed increased cell expansion (Li et al., 2017; Weinkove et al., 2019), the clinical results have varied greatly between design iterations. Recent studies point to the variegated cellular differentiation status amongst products leading to heterogenous responses between individuals; however, more universal is the data supporting less differentiated cells having increased function and persistence (Fraietta et al., 2018; Hoffmann et al., 2018). As a result, optimizing the manufacturing process in CAR T cell product generation is now a significant goal for translational researchers and clinicians alike and product manufacturing has been identified as a key factor that drives functional differences between infused products (Gargett et al., 2019; Stock et al., 2019; Tyagarajan et al., 2020).

Media composition

The composition of the culturing media plays a critical role in determining the differentiation and functional status of transduced cells during the activation and expansion steps of CAR T cell manufacturing. Recent findings have emphasized the utility of selecting culture conditions that minimize the length of manufacturing and at the same time yield high absolute numbers of functionally relevant subpopulations responsible for mediating tumor control (Ghassemi et al., 2018; Ghassemi et al., 2022). Specifically, previous analyses looking at the composition of infused cells have determined that cells expressing surface proteins such as CD27⁺, CCR7+, and CD45RA + are less differentiated and belong to subgroups of T cells, such as central memory (T_{CM}) and naive (T_N) subsets, with greater proliferative capacity and persistent antitumor function (Gattinoni et al., 2012; Busch et al., 2016; Stock et al., 2019).

Recent comparative studies using commercially available T cell growth media formulations, such as RPMI-1640, Optimizer[™], X-VIVO 15[™], and TexMACS[™], have examined how the choice of base media can influence the expansion of CAR T cells. The findings suggest that media selection should be based on which T cell subsets are desired at the end of manufacturing as different media lead to varied phenotypes and functions of the T cells (Sato et al., 2009; Lu et al., 2016). In addition to choosing an appropriate base medium, media additives can also alter the end product. Specific sources of serum significantly impact the potency and proliferation of CAR T cells (Gardner et al., 2017; Fraietta et al., 2018; Ghassemi et al., 2020), and consideration of xeno- and serum-free media can lead to increased viability and reduced T cell exhaustion (Sanyanusin et al., 2023; Sartorius, 2023). Addition of specific supportive cytokines during cell activation and expansion can also drastically alter the relative frequency of T cells across the spectrum of differentiation states (Gargett and Brown, 2015; Gong et al., 2019; Zhang X. et al., 2022). While many commercial products rely on interleukin (IL)-2 as the sole supportive growth factor, emerging studies have highlighted the benefits of replacing IL-2 with IL-7 and IL-15 and have shown that the later cytokines result in a final cell product with a less differentiated phenotype and enhanced therapeutic potential (Cieri et al., 2013; Alizadeh et al., 2019; Zhou et al., 2019; Battram et al., 2021; Kim et al., 2022).

To reduce the frequency of terminally differentiated cells, others have sought to inhibit specific metabolic processes via the alteration of cell signaling processes (Zheng et al., 2018; Li W. et al., 2023). The addition of Protein kinase B (AKT) inhibitors during the manufacturing of CD19-specific CAR T cells has been shown to increase antitumor function by blocking phosphatidylinositol 3kinase (PI3K) signaling while maintaining mitogen-activated protein kinase (MAPK) signaling, leading to accumulation of T_{CM} associated transcription factors such as FOXO1 (Klebanoff et al., 2017; Urak et al., 2017; Coleman et al., 2021; Delpoux et al., 2021). Another strategy for improving the quantity and quality of memory T cells is to target the mammalian target of rapamycin (mTOR) and IL-2 receptor signaling pathways. Previous studies have shown that introduction of rapamycin during manufacturing of T cell products can enhance T cell functionality, increase T cell viability and resistance to apoptosis, and alter the metabolic state of activated T cells by partially inhibiting mTOR (Araki et al., 2010; Schmueck et al., 2012; Alizadeh et al., 2019; Amini et al., 2019). Taking a different approach, the addition of Wnt/ β -catenin signaling pathway modulators has been shown to impede the transition of naïve to effector T cell subsets, which may provide additional avenues to increase the frequency of less differentiated and highly functional CAR T cell products (Gattinoni et al., 2009; Muralidharan et al., 2011; Kondo et al., 2018). Overall, these studies highlight the importance of carefully considering and optimizing the composition of culturing media to expand functionally potent CAR T cells.

Isolation and depletion of select subsets

The isolation and/or depletion of select subsets play an important role in the manufacturing process of CAR T cells, particularly when the goal is to increase the representation of a predefined cell population within the final product. Enrichment of specific CD8⁺ and CD4⁺ subsets has been explored as a strategy to enhance the desired characteristics of CAR T cells (Fraietta et al., 2018; Lee et al., 2018; Shah et al., 2020a; Kim-Hoehamer et al., 2023). As each cellular subset within the CD8/CD4 lineage present unique functional characteristics, the identification of populations responsible for improved responses has been pursued as a mechanism of therapeutic improvement. For example, removal of CD4⁺ T cells with enhanced regulatory capabilities, such as regulatory T cells (T_{reg}), results in increased overall antitumor activity (Onda et al., 2019). By incorporating these isolation and depletion techniques into the manufacturing process, researchers aim to optimize the composition and functionality of CAR T cells. An approach that has been shown to dictate efficacy of CAR T cells is the infusion of products with pre-defined ratios of select T cells. Specifically, less differentiated subsets (T_N and T_{CM}) of CD4⁺ and CD8⁺ CAR T cells were more effective in target eradication than those from effector memory (T_{EM}) subpopulations, contributing to the correlation between antitumor activity and peak proliferation of CD4⁺ and CD8⁺ CAR T cells (Sommermeyer et al., 2016). As each transduced T cell subset offers unique and sometimes synergistic antitumor functions, the administration of CAR T cell products that consider the proportion of each subset display greater antitumor activity and improved progression-free survival (Turtle et al., 2016; Gardner et al., 2017).

Prior to enrichment of specific CD8+ and CD4+ T cell subsets, an initial characterization of these cells in the aphaeresis product is important to determine whether enrichment would have any added benefit. If the pre-manufactured product contains high frequencies of terminally exhausted T cells, other approaches, such as allogeneic PBMCs from a healthy donor (Graham et al., 2018; Young et al., 2022), may need to be considered. If an allogenic approach is necessary, additional considerations must be taken to minimize the adverse effects of highly reactive CAR T cells and potential allogenic responses. One approach that has shown promise in these scenarios is the depletion of CD45RA-expressing naïve-like cells from donor PBMCs. This depletion strategy removes T cells possessing a broader T cell receptor (TCR) repertoire that can cross-react with antigens expressed by the recipient (Bleakley et al., 2015), but it also preserves the long-lived functional capabilities of differentiated memory-like T cells that have more restricted TCR repertoires that minimize allogeneic responses (Wherry et al., 2003; Klebanoff et al., 2005). Recent advancements in developing good manufacturing practices (GMP) for CD45RA + cell depletion have been described, facilitating its potential application in CAR T cell manufacturing (Zheng et al., 2018). This has been further tested in the treatment of subjects with lymphocytic leukemia, where the graft-versus-leukemia (GVL) effect is optimized while reducing the risk of graft-versus-host disease (GVHD) (Dutt et al., 2007; Biernacki et al., 2020; Kim-Hoehamer et al., 2023; Naik and Triplett, 2023).

CAR T cell in-vitro activation

Effective expansion of T cells requires adequate, sustained, and sequential activation signals through a series of surface receptors. The primary signal for T cell activation results from the engagement of the TCR with its ligand (signal 1) and the magnitude of T cell activation and effector function is primarily driven by engagement of costimulatory receptors (signal 2), such as CD28, 4-1BB, ICOS (Bretscher, 1999; Chen and Flies, 2013). This activation triggers the internal proliferative program of T cells, leading to their clonal expansion (Noak et al., 2021).

As specific doses of CAR-expressing T cells within the manufactured product have been shown to have greater clinical efficacy, this activation step can help achieve sufficient numbers of CAR T cells. Various methods of expansion including the use of cell-, antibody-, bead-, and polymer-based activation have been tested with differing levels of success (Gee, 2018; Abou-el-Enein et al., 2021). Activation strategies using cell-based approaches tend to mirror normal immune synapse interactions and drive efficient expansion, but this approach requires greater technical skill, is very laborious, and can be problematic due to having to source sufficient numbers of cells (Sasawatari et al., 2006; Schmidts et al., 2020). An alternative strategy, and one of the most common in current clinical manufacturing workflows, involves the use of super-paramagnetic anti-CD3/CD28 antibody-coated beads, such as Dynabeads[™], which provide a reproducible source of T cell activation (Barrett et al., 2014; Gargett and Brown, 2015; Vormittag et al., 2018). One drawback of using Dynabeads[™] for activation is the need to remove the beads once sufficient numbers of CAR T cells are obtained, which extends the vein-to-vein infusion time. Newer alternatives, such as a polymeric with bio-degradable nanomatrix impregnated anti-CD3/ CD28 monoclonal antibodies, eliminate the added step of bead removal, expediting the manufacturing process of CAR T cells and resulting in a less differentiated composition of cells (Gargett et al., 2019). It should be noted that the use of such nanomatrices do however increase the overall cost of manufacturing (Gee, 2018). While T cell activation is often a necessary step for introducing constructs using retroviruses, newer protocols aim to eliminate cell activation altogether to reduce the manufacturing time and generate cells with less exhausted phenotypes (Yang et al., 2020; Ghassemi et al., 2022). A recent manufacturing protocol that performs same day T cell isolation using CD3/CD28 beads and next-day lentiviral transductions aims to reduce vein-to-vein time and address the ex-vivo culturing time required during CAR-T manufacturing. A benefit of this protocol is that it leads to greater T cell stemness in the CAR T cell product, which increases cellular persistence following infusion (Ghassemi et al., 2018). Clinically, CAR-T cells manufactured using the nextday manufacturing process displayed superior expansion and a greater proportion of patients achieved favorable responses compared to conventional CAR T cell manufacturing (Jiang et al., 2020; Zhang C. et al., 2022; Yang et al., 2022).

CAR construct selection considerations

Development of constructs with unique target specificities and distinct intracellular signaling domains also impacts the extent of T cell activation, expansion, and differentiation (Ghorashian et al., 2019; Mezősi-Csaplár et al., 2023). While first-generation CAR T cells consisted of an extracellular binding domain, transmembrane, and intracellular signaling domain, the lack of expansion in these early CAR T cells has lead researchers to invest significant effort towards incorporating and modifying specific domains in newer iterations with the hope of driving greater expansion and function (Sadelain et al., 2009). Second and third generation CARs incorporate either single or dual costimulatory signaling domains respectively and intend to recapitulate native co-stimulation in T cells, which leads to enhanced intracellular signaling and greater T cell expansion and potency (van der Stegen et al., 2015; Feins et al., 2019; Subklewe et al., 2019; Shah et al., 2020b). The latest (fourth) generation of CARs constructs takes the previous principles of having co-stimulatory domains and incorporates further genetic modifications such as additional co-stimulatory ligands or transgenes for cytokine secretion (Subklewe et al., 2019; Abken, 2021). These latest generation of CAR T cells seek to expand the utility of CARs beyond just target antigen recognition and provide additional avenues to engineer CAR T cells with additional capacities, such as altering the local tumor microenvironment. Clinically, each new iteration of CAR construct is met with challenges, with earlier generations presenting with increased relapsed rates and later generations resulting in greater toxicities (Park et al., 2016; Cappell and Kochenderfer, 2023).

Methods for introducing CAR constructs to T cells

Traditionally, viral-based techniques have been used for CAR T cell production, which involves delivering CAR genes into T cells using retroviral or lentiviral vectors (Rosenberg et al., 1990; Scholler et al., 2012; Gardner et al., 2017; Marcucci et al., 2018). These viruses randomly integrate their genetic material into the host cell genome, enabling long-term expression of the CAR. This approach has demonstrated remarkable success in generating CAR T cells with potent antitumor activity and has led to significant advancements in clinical settings. However, viral-based methods come with challenges. Safety concerns arise due to the potential for viral integration causing insertional mutagenesis or activating oncogenes (Scholler et al., 2012; Cavazza et al., 2013). Manufacturing viral vectors is complex, time-consuming, and expensive, hindering scalability and clinical accessibility. Moreover, regulatory requirements associated with viral vectors add further obstacles to the development and commercialization of CAR T cell therapies.

To address these limitations, non-viral synthetic biology methods have emerged as alternative approaches for CAR T cell production. RNA-based delivery systems have gained attention for achieving transient CAR expression without genomic integration (Rabinovich et al., 2009; Zhao et al., 2010; Soundara Rajan et al., 2020). Messenger RNA (mRNA) encoding the CAR construct is introduced into T cells, leading to efficient CAR expression. This method allows for fine-tuned control of CAR expression levels by adjusting RNA dosage and duration of expression, thereby reducing off-target effects. Unlike DNA integration, RNA offers reversible and transient control over gene expression, enabling CAR expression to subside over time and without constitutive expression in progeny cells after CAR T cell expansion. This adaptability may facilitate adjustments of CAR expression in response to patient conditions, optimizing therapy efficacy and minimizing risks for more tailored outcomes. Furthermore, it is important to highlight the potential impact of the controlled expression on mitigating cytokine release syndrome (CRS), immune effector cell-associated neurotoxicity syndrome (ICANS) (Gauthier and Turtle, 2018; Siegler and Kenderian, 2020), and tonic signaling leading to the exhaustion of T cells (Lamarche et al., 2023)-all of these detrimental effects have been associated with CAR overexpression resulting from viral-based technologies. This reduction in potential adverse effects could contribute significantly to the safety and efficacy of CAR T cell therapies. However, repeated dosing is necessary for sustained CAR expression, posing challenges for large-scale manufacturing and clinical implementation (Moretti et al., 2022).

DNA transposon systems offer a non-viral method for stable CAR integration into the T cell genome. DNA transposons are DNA segments flanked by inverted repeats that can "jump" into the host genome with the assistance of a transposase enzyme. Incorporating the CAR gene within a transposon allows for its stable integration, providing long-term CAR expression without relying on viral vectors. The transposon, containing the CAR gene, is flanked by inverted terminal repeats (ITRs) and delivered as a (nano-) plasmid or minicircle DNA. The transposase enzyme assists in integrating the transposon into the acceptor DNA at specific target sites, such as -TA- (Sleeping Beauty) (Ivics et al., 1997) or -TTAA- (PiggyBac) (Gogol-Döring et al., 2016) sequences. While the integration may not be precisely targeted, it can occur in preferred insertion sites known as "safe harbors" (Querques et al., 2019) and offers long-term CAR expression. Researchers are actively exploring different transposon systems such as the Sleeping Beauty and piggyBac systems (Kebriaei et al., 2016; Magn et al., 2020; Zhang X. et al., 2022) to enhance their performance and suitability for clinical applications. Although DNA transposon systems offer improved safety, optimizing their efficiency and minimizing off-target effects are ongoing challenges. The CARTELL trial, which employed PiggyBac technology for the treatment of relapsed/refractory B cell malignancies, encountered a concerning outcome as two out of ten patients developed lymphoma (Micklethwaite et al., 2021). This outcome has raised safety concerns surrounding the use of PiggyBac and transposon/transposase-based cell therapies for the treatment of B cell malignancies. The development of lymphoma may be related to the manufacturing process or an increase in global copy number changes observed in the products. In response to these safety concerns, there is a recognized need for the development of enhanced preclinical genotoxicity models and optimization strategies. Transcription activator-like effector nucleases (TALEN) and CRISPR/Cas9-mediated genome engineering provide avenues for site specific modification in T cells. TALENs, composed of hybrid molecules merging DNA recognition proteins (transcription factors) with an endonuclease, utilize TAL units of 33-35 amino acids to recognize a single base pair on genomic DNA. Linking several TALs with an endonuclease creates a site-specific TALEN, a crucial tool for precision engineering in the development of CAR-T cells (Jo et al., 2022). CRISPR/Cas9-mediated genome engineering, on the other hand, uses the CRISPR/Cas9 (Jinek et al., 2012) system to introduce precise modifications in the T cell genome (Dimitri et al., 2022). By designing a guide RNA (gRNA) specific to the target site, the Cas9 enzyme can cleave the DNA, triggering DNA repair mechanisms and allowing for gene modifications such as CAR integration or endogenous gene knockouts (Eyquem et al., 2017; Roth et al., 2018; Kamali et al., 2021; Shy et al., 2021). This method offers versatility and precision in modifying the T cell genome, enabling customized CAR designs and improved therapeutic efficacy. Nevertheless, the application of CRISPR technology is not without risks, and ongoing research is dedicated to addressing concerns such as reducing off-target effects, managing mosaicism, addressing potential chromosome translocations, and refining DNA repair processes. Current research efforts are focused on enhancing efficiency, mitigating risks, and ensuring the safety of CRISPR-edited cells, with the goal of providing effective gene engineering strategies to alter the T cell's genome for therapeutic purposes.

Recently, an RNA-guided endonuclease, Fanzor (Fz) (Saito et al., 2023) protein, has been identified from the eukaryotic system and can be used for genome editing. Although this newly identified protein has not yet been used for the generation of CAR T cells, its eukaryotic origin and relatively small size compared to Cas9/12 make it an attractive starting point for further development and use for cellular engineering. Engineering strategies such as systematic mutagenesis and guide RNA engineering combined with in-depth screening of more Fanzors could further improve their genome editing performance, highlighting the potential of Fzs for cellular engineering.

The traditional viral-based methods, while successful, face challenges related to safety, scalability, and manufacturing complexity. Non-viral methods like RNA-based delivery systems and DNA transposon systems offer advantages in terms of controlled expression and stable integration, respectively. TALENs and CRISPR/Cas9 provide greater precision in genome editing, but CRISPR/Cas9 can also present specific risks such as chromosomal translocations. Fanzor protein, although promising, requires further development as this is a relatively new cell engineering strategy.

Experimental approaches to test functionality for potency testing

As CAR T cell therapy manufacturing generates final products consisting of heterogeneous phenotypes, implementing standardized potency testing, and establishing thorough release criteria before infusion are imperative to ensure safety, efficacy, and consistency of care. Historically, quality control testing has primarily focused on the absence of contaminants like *mycoplasma* and endotoxins (Gee, 2018). However, the evolving field of biologics manufacturing has prompted the FDA and EMA to recognize the need for updated criteria encompassing assays that evaluate whether CAR T cells can mediate antitumor functions through cytotoxic- or cytokine-based mechanisms (Food and Drug Administration, 2019; European Medicines Agency, 2022; European Medicines Agency, 2018b; Gee, 2015; Si et al., 2022). By adopting standardized assays, the safety and efficacy of CAR T therapies can improve, ensuring administration of only high-quality products to patients (Aleksandrova et al., 2019; Roddie et al., 2019). Lastly, the emergence of new technologies enables more efficient preinfusion testing, facilitating the evaluation of the therapeutic potential of CAR T cells before infusion.

Secretory profiling

Proteomic characterization of secreted analytes has become a necessary criteria prior to infusion due to the potential risks following CAR T cell infusion is developing CRS and neurotoxicity (Maude et al., 2014; Brudno and Kochenderfer, 2016; Gardner et al., 2019; Hopfinger et al., 2019; Sterner and Sterner, 2021). To mitigate this risk, current release criteria rely on using multi-analyte profiling $(xMAP^{TM})$ which are bulk assays that measure soluble factors from cell cultures through the use of labeled microspheres that capture multiple analytes at a time (Maude et al., 2014; Graham et al., 2019; Melenhorst et al., 2022). Unfortunately, by characterizing cytokine secretion in bulk, these assays provide a global view of CAR T functionality and may not accurately depict the heterogeneous nature of the CAR T cell products and can pose challenges when trying to predict clinical outcome. Newer technologies incorporating single-cell approaches can capture the heterogenous nature of CAR T cells and provide a more granular view of how each cell behaves rather than taking an average measure across all cells. The development of instruments such as Bruker Cellular Analysis' IsoSpark allows interrogation of the secretory profile of CAR T cells at single cell resolution and has provided useful correlations to be made between the product's secretome and clinical response (Lu et al., 2013; Mackay et al., 2017; Rossi et al., 2018; Xu et al., 2020; Zurko et al., 2022). Using the Polyfunctional Strength Index (PSI) score, which describes the percentage of polyfunctional single cells in a sample secreting two or more analytes (IsoPlexis, 2023), standardized metrics could be incorporated as part of the release criteria for CAR T cells and may ultimately provide a better measure of clinical response and help minimize the toxicities observed in many patients.

Cytotoxic profiling

One of the most relevant functions that must be evaluated in the final product is the cells' ability to have sufficient cytotoxicity against antigen-expressing target cells. By using instruments that measure changes in target cell density with either electrical impedance or image-based approaches, one can evaluate the cytotoxic potential in bulk co-culture studies. Impedance-based assays are real-time cell analysis systems that monitor interactions without requiring labels or invasive techniques (Logun et al., 2023). Systems such as Agilent's xCelligence[™] and Axion's Maestro [™] can indirectly measure processes such as proliferation and cytotoxicity by using specialized microplates with integrated electrodes and measuring electrical impedance (Kiesgen et al., 2021; Lamarche, 2021; Axion Biosystems, 2024). Image-based approaches measure cytotoxicity in real-time, using automated high-resolution imaging and fluorescence analysis (Strietz and Chen, 2022). This technology implemented in instruments such as Sartorius' Incucyte[™], Enrich Biosystem's Trovo[™], and Nanolive's 3D Cell Explorer[™] can track and analyze additional cellular processes over time, including proliferation, cytokine secretion, and cytotoxicity (Zah et al., 2020; Nanolive, 2023; Sartorious, 2019). Both detection approaches provide useful readouts for functional testing and may provide clinically relevant data that could also be incorporated as part of the CAR T cell release criteria prior to infusion. The real-time monitoring capabilities and analytical power of these instruments highlight their value in understanding the dynamic cellular processes occurring between effector and target cells and can help generate data that more accurately reflects how CAR T cells behave in vivo.

Polyfunctional characterization

Most commercially available CAR T cell products implement a one-size-fits-all manufacturing process and often ignore interpatient differences that may dictate treatment outcome. By incorporating single-cell functional assays using technology capable of measuring multiple functional characteristics in parallel, researchers can comprehensively assess how the isolation and expansion of specific CAR T cell subsets will drive clinical efficacy (Mocciaro et al., 2018; Bronevetsky, 2020; Le et al., 2020; Romain et al., 2022). The recent explosion of various single cell platforms has revolutionized the field of functional biology and has allowed researchers to not only interrogate how cells function over time using single cells or groups of cells, but has also allowed for the isolation of specific subpopulations of cells with unique functional characteristics from a patient's blood or tissue on a per-patient basis (Bandey et al., 2021; Miwa et al., 2022; Urbani et al., 2022). Commercial technologies that have implemented this type of multiparameter testing, such as CellChorus' TIMING[™] and Celldom's $^{^{\rm TM}}$ CloneXplorer $^{^{\rm TM}}$ platforms, have utilized well-based technology to seed effector and target cells and record effector function using brightfield and fluorescent imaging (Liadi et al., 2015; Celldom, 2022). Other platforms such as Bruker Cellular Analysis' Beacon[™] and Lightning[™] systems utilize nondestructive light to sort cells into NanoPen® chambers distributed on $\mathsf{OptoSelect}^{^{\mathrm{TM}}}$ Chips and have the capability of longitudinally performing both cytokine secretion and cytotoxic assays at the single cell level allowing researchers to identify cells that display the best response to cancer (Bronevetsky, 2020; Berkeley Lights, 2024a). Functional cells of interest can then be isolated for additional expansion or to perform downstream assays, such as transcriptional profiling to identify immune signatures depicting unique functional states (Liadi et al., 2015; Berkeley Lights, 2024b).

One drawback is that due to the temporal nature of these assays, the current implementation of these characterizations as release criteria would lengthen the manufacturing time and delay infusion.

Challenges and next steps

Current advancements in technology and the implementation of in-depth characterization are opening avenues for improving the quality of infused CAR T cell products. The emergence of single-cell functional experiments provide a valuable opportunity to investigate and quantify phenotypic, functional, and transcriptional profiles of single cells. The knowledge gained will result in improved product design and manufacturing processes. Nevertheless, to fully unlock the potential of adoptive immunotherapies in treating refractory diseases, more robust and in-depth cellular characterizations uploaded to public databases will accelerate the improvement of future therapies. Despite the technological advances that help link functional variables and therapeutic efficacy, gaps remain in understanding the impact of previous treatment modalities on CAR T cell efficacy and how other factors, such as the temporal evolution in the tumor cells' ability to evade T-cell-mediated killing, result in suboptimal responses. To gain insight into these limitations, temporal transcriptional profiling of patient samples at single-cell resolution and longitudinal characterization of CAR T cells preand post-infusion has helped validate unique gene-expression profiles that gave rise to highly functional post infusion phenotypes (Wilson et al., 2022). Despite increased efforts to characterize and validate CAR T cell products, there are inherent risks that come with any cell-based engineering strategy. Recently, the FDA has begun evaluating whether greater regulatory action is needed for CAR T cell therapies as potential risks of developing T cell malignancies and secondary hematological cancers have been reported in patients that have received this line of therapy. While the current consensus appears to be that the benefits of CAR T cell therapy outweigh the potential risks for patients with cancer, these recent cases highlight the need for improved strategies aimed at preventing random integration of CAR constructs. More stringent characterization of CAR T cell products may also help identify cases where CAR-associated malignancies arise (FDA, 2023a; Harrison et al., 2023; Mullard, 2023).

By optimizing specific steps in the manufacturing of CAR T cells and incorporating some of the potency assays highlighted in this review, patients will begin receiving a more robust product with a greater proportion of cytotoxic CAR T cells that have longer persistence. The next steps would be to consider whether combination therapies using other therapeutic modalities or to overcome an immunosuppressive helping tumor microenvironment would increase success rates (Huang et al., 2022). An added benefit in building more comprehensive multimodal datasets is that this information can also help inform whether specific manufacturing steps should be altered and what the likely clinical outcomes would be. Finally, these data could also allow for a greater ability to identify individuals likely to respond to CAR T cell therapy using specific biomarkers or help physicians decide whether a different course of treatment is warranted.

In silico approaches to investigate CAR T-cell function and predict therapeutic responses

As outlined in previous sections, there are several possible variables introduced during CAR T cell manufacturing which will influence the CAR T cell product. It can be extremely difficult to interrogate how these decisions ultimately impact product functionality and patient responses using currently available experimental techniques alone. Recently, there has been an explosion of computational models aiming to tease apart relationships between the CAR T cell product manufacturing process, product functionality, and the resulting therapy responses to predict how manufacturing variables should be modulated to improve overall clinical response. We direct the reader to substantial prior works for more in-depth reviews on computational models and modeling methods in CAR T cell therapy (Nukala et al., 2021; Qi et al., 2022). In this section, we will review notable approaches and discuss their potential to inform CAR T cell manufacturing and improve the clinical efficacy of the therapy.

Computational models for characterizing CAR T cell function

Changes in CAR insertion methods can modulate CAR surface expression and avidity toward target antigens leading to changes in CAR T cell functionality. Mechanistic models based on kineticproofreading concepts have been developed to investigate intricate relationships between antigen densities, CAR surface expression and ligand-receptor binding rates. These models demonstrate that it is possible to control CAR T cell product functionality through tuning CAR receptor expression and antigen binding rates in relation to tumor cell antigen densities (Ha et al., 2018). Additionally, these models hypothesize that it may be possible to restrict CAR T cell activity toward tumor cells and mitigate off-target effects (cytolysis of healthy cells) by tuning avidities (Ha et al., 2018; Rajakaruna et al., 2023). These models offer an intriguing design tool for optimizing the performance of CAR T cell products in specific tumors, and potentially in specific patients, as a function of tumor antigen densities.

Additionally, as each domain in the CAR construct can be switched and modified (Mazinani and Rahbarizadeh, 2022), exactly how these changes modulate intracellular signaling and product functionality are important to consider. A mechanismbased (mechanistic) ordinary differential equation (ODE) model has been developed to study the intricate kinetics of CAR receptor and intracellular extracellular signal-regulated kinase (ERK) activation in the presence and absence of CD28 co-stimulation supported by protein phosphorylation kinetics derived using phosphoproteomics (Rohrs et al., 2018; Rohrs et al., 2020). Insights from the model can be used to understand the mechanisms of CAR signal propagation in CAR T cells and help develop hypotheses for optimizing CAR designs to improve product functionality through modulating the CD3 ζ and CD28 signaling domains. It is important to note, however, that models such as these are difficult to develop and are limited to study specific construct designs and signaling contexts which ultimately limit their utility. As the field is quickly advancing toward widely variable CAR constructs (Mazinani and Rahbarizadeh, 2022), more robust computational methods need to be developed to enable optimization techniques.

Computational models for CAR T cell therapeutic responses

Frequently described as a "living drug," therapy dynamics for CAR T products are exceptionally challenging to model using established pharmacokinetic (PK) and pharmacodynamic (PD) approaches used for small molecule drugs. One key complication is that CAR T cells continue to differentiate during manufacturing and after infusion leading to variable kinetics of cytotoxicity and proliferation. Further complicating model development are the inter-product and inter-patient variabilities in CAR design, differences in manufacturing conditions and patient-specific baseline T cell phenotype and function following pretreatment, tumor burden and antigen expression, and tumor microenvironments, which collectively impact therapeutic response. Nonetheless, considerable efforts have been directed at developing computational models for therapy response. Several proposed models use ODE-based approaches with simple predator-prey (Lotka-Volterra) mechanisms for modeling dynamics between CAR T cells and tumor cells (Hardiansyah and Ng, 2019; Sahoo et al., 2020). Additionally, mathematical relationships derived from enzyme kinetics and empirical relationships have also been used to describe these interactions (Stein et al., 2019; Singh et al., 2020; Liu et al., 2021). Once formulated, therapy response models can be used to explore how various facets of the therapy ultimately influence response. Additional variables such as a patient's age, disease burden, type and duration of bridging therapy, lymphodepletion regimen, and whether prior lines of therapy have been administered can be incorporated into these models. Additionally, models can give insight into the impact of dosing and lymphodepletion regimens on resulting therapy kinetics (Hardiansyah and Ng, 2019; Stein et al., 2019; Kimmel et al., 2021; Owens and Bozic, 2021). Parameters derived from these models can help elucidate kinetic rate constants for target cytolysis and CAR T cell proliferation (Kimmel et al., 2020; Liu et al., 2021; Amini et al., 2022; Lickefett et al., 2023). These parameters can potentially be used in conjunction with functional assays and other variables to help evaluate future products.

Importantly, models have also begun to describe how different phenotypic and functional states of the CAR T cell product can shape *in vivo* expansion and tumor clearance (Stein et al., 2019; Mueller-Schoell et al., 2021; Paixão et al., 2022; Kirouac et al., 2023). With further validation against clinical and experimental datasets, these models have the potential to be used to define additional biomarkers for predicting the survival of patients receiving CAR T cell therapy based on assessing product composition pre-infusion and tracking changes in CAR T cell compositions post-infusion (Mueller-Schoell et al., 2021; Paixão et al., 2022). It may also be possible to use this information to develop optimal CAR T cell product compositions to improve clinical outcome which can then be used to further optimize variables at different steps of manufacturing. Additional applications for therapy response models are to predict the occurrence of serious side-effects associated with therapy administration such as CRS and neurotoxicity. Physiologically based PK and PD (PBPK/PD) models are especially useful in understanding the possible CAR T cell distribution after they are infused into the body and help develop exposure-response relationships for neurotoxicity potentially by proxy of CAR T trafficking into non-tumor compartments (Singh et al., 2020). This information may provide additional insight into which subsets of CAR T cells to deplete to prevent most of the cells trafficking to non-tumor compartments. CRS can also be examined through determining the rate of production of key cytokines post therapy administration (either by CAR T cells or other sources in the body). Models can help predict the dynamics of cytokine release to direct timing of intervention measures and offer methods to demonstrate the interaction between prophylactic drugs and CAR T cells (Hardiansyah and Ng, 2019; Stein et al., 2019; Zhang Z. et al., 2022). These models in conjunction with secretory profiling methods discussed in previous sections may provide additional depletion strategies to avoid CRS.

One major limitation of the ODE-based models for CAR T cell therapy overviewed thus far is that they require preexisting knowledge of physiological mechanism and robust datasets to estimate parameters which inherently limits the number of variables the models can include. To sidestep some of these obstacles, alternative approaches utilizing logic-based and agent-based techniques are also being explored. These models are able to capture more of the complexities observed in CAR T cell therapy such as relationships between key cytokines, CAR T cell subsets (CD4, CD8), intracellular signaling, and environmental factors to develop more comprehensive models (Selvaggio et al., 2022; Prybutok et al., 2022; Shah et al., 2023). These approaches are attractive for exploratory studies for CAR T cell therapy to help quickly interrogate the influence of multiple variables and drive hypotheses for fine tuning manufacturing conditions. It remains to be seen whether these models can be validated and/or provide utility beyond ideation.

It should be noted that therapy response models developed to-date and reviewed here either aim to study the system on a theoretical level or are otherwise descriptive of specific patient or patient-cohort therapy responses considered for developing the model. Given the wide range of factors during CAR T cell manufacturing that control possible therapy responses, models currently offer limited predictive and prognostic value for future products and patients as critical parameters used for the models cannot be derived prior to therapy administration and will generally differ on a product and patient basis. Using in vitro assays to drive modeling can offer a potential avenue to both validate physiological mechanisms used in models and to derive parameters necessary for model simulations ahead of therapy administration (Sahoo et al., 2020; Kirouac et al., 2023; Li et al., 2023b). Experimental techniques reviewed in previous sections may offer modelers additional methods with which to support model development and help enable clinical applications for these models.

Challenges and next steps

One of the major challenges in developing clinically relevant computational models for CAR T cell therapy is the availability of

well characterized and robust datasets that can be used to formulate these models. Current experimental and clinical datasets fall short in routinely and comprehensively characterizing CAR T cell products, compositions, cytolysis and proliferation subset rates. biodistribution, and receptor expression. There are similar limitations for tumor cells due to the lack of data on proliferation rates and antigen expression. This poses challenges when trying to develop models to determine how key manufacturing variables will influence response. With larger and better curated datasets, the methods used to develop the models reviewed here for describing CAR T cell product functionality and therapy response have the potential to inform new manufacturing protocols and enable predictive and precision medicine in CAR T cell therapy. There is opportunity to utilize these preliminary works to provide rationale for the collection of appropriate datasets, especially where the impact on clinical outcome is not immediately clear. The goal is to advance and begin implementing computational techniques that may offer avenues to ultimately improve therapeutic effectiveness in the long-term. An overview of how computational models can be used to simulate CAR T cell therapy response and temporal dynamics is depicted in Figure 2. With the recent expansion and incorporation of artificial intelligence (AI) in systems biology and medicine, there is growing interest in determining how to best apply AI to improve cell-based manufacturing. Development of automated AI-driven CAR T cell manufacturing processes have been proposed and aim to decrease the cost, reduce the amount of manual labor, and improve the efficacy of CAR T cell products (Hort et al., 2022). These computational approaches provide additional opportunities to rapidly advance the next line of CAR T cell products.

Author contributions

AC: Conceptualization, Writing-original draft, Writing-review and editing. VS: Conceptualization, Writing-original draft, Writing-review and editing. RS: Conceptualization, Writing-original draft, Writing-review editing. TK: and Writing-review and editing. Conceptualization, RD: Writing-original draft, Writing-review and editing. ST: Conceptualization, Writing-original draft, Writing-review and

References

Abken, H. (2021). Building on synthetic immunology and T cell engineering: a brief journey through the history of chimeric antigen receptors. *Hum. Gene Ther.* 32 (19–20), 1011–1028. doi:10.1089/hum.2021.165

Abou-el-Enein, M., Elsallab, M., Feldman, S. A., Fesnak, A. D., Heslop, H. E., Marks, P., et al. (2021). Scalable manufacturing of CAR T cells for cancer immunotherapy. *Blood Cancer Discov.* 2 (5), 408–422. doi:10.1158/2643-3230. BCD-21-0084

Aleksandrova, K., Leise, J., Priesner, C., Melk, A., Kubaink, F., Abken, H., et al. (2019). Functionality and cell senescence of CD4/CD8-selected CD20 CAR T cells manufactured using the automated CliniMACS Prodigy[®] platform. *Transfus. Med. Hemotherapy* 46 (1), 47–54. doi:10.1159/000495772

Alizadeh, D., Wong, R. A., Yang, X., Wang, D., Pecoraro, J. R., Kuo, C. F., et al. (2019). IL15 enhances CAR-T cell antitumor activity by reducing mTORC1 activity and preserving their stem cell memory phenotype. *Cancer Immunol. Res.* 7 (5), 759–772. doi:10.1158/2326-6066.CIR-18-0466

Amini, L., Silbert, S. K., Maude, S. L., Nastoupil, L. J., Ramos, C. A., Brentjens, R. J., et al. (2022). Preparing for CAR T cell therapy: patient selection, bridging therapies and editing. AZ: Conceptualization, Funding acquisition, Visualization, Writing–original draft, Writing–review and editing.

Funding

The author(s) declare financial support was received for the research, authorship, and/or publication of this article. This work was supported by the Paula and Rodger Riney Foundation. AC and VS are supported in part by the Medical College of Wisconsin Cancer Center's Graduate Fellowship.

Acknowledgments

We thank our colleagues and collaborators at the Medical College of Wisconsin for their engaging discussions and feedback on the contents of this review. We would also like to thank the Paula and Rodger Riney Foundation for supporting our research in multiple myeloma and in developing novel CAR T cell therapies. We would also like to thank the Cell Therapy Lab (CTL) Shared Resource at the Medical College of Wisconsin Cancer Center for providing resources that allow us to perform experiments to test the functionality and potency of CAR T cells.

Conflict of interest

AZ reports research support from Bruker Cellular Analysis. The remaining authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Publisher's note

All claims expressed in this article are solely those of the authors and do not necessarily represent those of their affiliated organizations, or those of the publisher, the editors and the reviewers. Any product that may be evaluated in this article, or claim that may be made by its manufacturer, is not guaranteed or endorsed by the publisher.

lymphodepletion. Nat. Rev. Clin. Oncol. 19 (5), 342-355. doi:10.1038/s41571-022-00607-3

Amini, L., Vollmer, T., Wendering, D. J., Jurisch, A., Landwehr-Kenzel, S., Otto, N. M., et al. (2019). Comprehensive characterization of a next-generation antiviral T-cell product and feasibility for application in immunosuppressed transplant patients. *Front. Immunol.* 10, 1148. doi:10.3389/fimmu.2019.01148

Araki, K., Youngblood, B., and Ahmed, R. (2010). The role of mTOR in memory CD8 T-cell differentiation. *Immunol. Rev.* 235 (1), 234–243. doi:10.1111/j.0105-2896. 2010.00898.x

Axion Biosystems (2024). Tracking the dynamics of cytotoxicity in real-time with the Maestro Z impedance assay. Available at: https://www.axionbiosystems.com/resources/application-note/tracking-dynamics-cytotoxicity-real-time-maestro-z-impedance-assay (Accessed July 26, 2023).

Bandey, I. N., Adolacion, J. R. T., Romain, G., Paniagua, M. M., An, X., Saeedi, A., et al. (2021). Designed improvement to T-cell immunotherapy by multidimensional single cell profiling. *J. Immunother. Cancer* 9 (3), e001877. doi:10.1136/jitc-2020-001877

Barrett, D. M., Singh, N., Liu, X., Jiang, S., June, C. H., Grupp, S. A., et al. (2014). Relation of clinical culture method to T-cell memory status and efficacy in xenograft models of adoptive immunotherapy. *Cytotherapy* 16 (5), 619–630. doi:10.1016/j.jcyt. 2013.10.013

Battram, A. M., Bachiller, M., Lopez, V., Fernández de Larrea, C., Urbano-Ispizua, A., and Martín-Antonio, B. (2021). IL-15 enhances the persistence and function of BCMA-targeting CAR-T cells compared to IL-2 or IL-15/IL-7 by limiting CAR-T cell dysfunction and differentiation. *Cancers* 13 (14), 3534. doi:10.3390/cancers13143534

Berkeley Lights (2024a). BROCHURE: comprehensive single T cell profiling. Available at: https://www.berkeleylights.com/resources/brochure-comprehensivesingle-t-cell-profiling/ (Accessed August 16, 2023).

Berkeley Lights (2024b). Datasheet: OptoSeq[™] single cell 3' mRNA kit. Available at: https://www.berkeleylights.com/resources/optoseq-single-cell-3'-mrna-kit/ (Accessed August 16, 2023).

Biernacki, M. A., Sheth, V. S., and Bleakley, M. (2020). T cell optimization for graftversus-leukemia responses. *JCI Insight* 5 (9), e134939. doi:10.1172/jci.insight.134939

Bleakley, M., Heimfeld, S., Loeb, K. R., Jones, L. A., Chaney, C., Seropian, S., et al. (2015). Outcomes of acute leukemia patients transplanted with naive T cell-depleted stem cell grafts. *J. Clin. Invest.* 125 (7), 2677–2689. doi:10.1172/JCI81229

Boyiadzis, M. M., Dhodapkar, M. V., Brentjens, R. J., Kochenderfer, J. N., Neelapu, S. S., Maus, M. V., et al. (2018). Chimeric antigen receptor (CAR) T therapies for the treatment of hematologic malignancies: clinical perspective and significance. *J. Immunother. Cancer* 6 (1), 137. doi:10.1186/s40425-018-0460-5

Bretscher, P. A. (1999). A two-step, two-signal model for the primary activation of precursor helper T cells. *Proc. Natl. Acad. Sci.* 96 (1), 185–190. doi:10.1073/pnas.96.1.185

Bronevetsky, Y. (2020). Directly test individual T cell function with fewer cells on the berkeley lights lightning[™] platform. *Cytotherapy* 22 (5), S119. doi:10.1016/j.jcyt.2020. 03.221

Brudno, J. N., and Kochenderfer, J. N. (2016). Toxicities of chimeric antigen receptor T cells: recognition and management. *Blood* 127 (26), 3321–3330. doi:10.1182/blood-2016-04-703751

Busch, D. H., Fräßle, S. P., Sommermeyer, D., Buchholz, V. R., and Riddell, S. R. (2016). Role of memory T cell subsets for adoptive immunotherapy. *Semin. Immunol.* 28 (1), 28–34. doi:10.1016/j.smim.2016.02.001

Byrne, M., Oluwole, O. O., Savani, B., Majhail, N. S., Hill, B. T., and Locke, F. L. (2019). Understanding and managing large B cell lymphoma relapses after chimeric antigen receptor T cell therapy. *Biol. Blood Marrow Transpl.* 25 (11), e344–e351. doi:10. 1016/j.bbmt.2019.06.036

Cappell, K. M., and Kochenderfer, J. N. (2023). Long-term outcomes following CAR T cell therapy: what we know so far. *Nat. Rev. Clin. Oncol.* 20 (6), 359–371. doi:10.1038/ s41571-023-00754-1

Cavazza, A., Moiani, A., and Mavilio, F. (2013). Mechanisms of retroviral integration and mutagenesis. *Hum. Gene Ther.* 24 (2), 119–131. doi:10.1089/hum.2012.203

Celldom (2022). Integrated platform: next-gen live cell assays. Available at: https://static1. squarespace.com/static/645c24ad9669f64e4b2502a2/t/64d6a02242f4b93d2d8ca950/ 1691787299551/CloneXplorer+Brochure+Feb2023.pdf (Accessed August 16, 2023).

Chen, L., and Flies, D. B. (2013). Molecular mechanisms of T cell co-stimulation and co-inhibition. *Nat. Rev. Immunol.* 13 (4), 227–242. doi:10.1038/nri3405

Cieri, N., Camisa, B., Cocchiarella, F., Forcato, M., Oliveira, G., Provasi, E., et al. (2013). IL-7 and IL-15 instruct the generation of human memory stem T cells from naive precursors. *Blood* 121 (4), 573–584. doi:10.1182/blood-2012-05-431718

Coleman, N., Moyers, J. T., Harbery, A., Vivanco, I., and Yap, T. A. (2021). Clinical development of AKT inhibitors and associated predictive biomarkers to guide patient treatment in cancer medicine. *Pharmacogenomics Pers. Med.* 14, 1517–1535. doi:10. 2147/PGPM.S305068

Delpoux, A., Marcel, N., Hess Michelini, R., Katayama, C. D., Allison, K. A., Glass, C. K., et al. (2021). FOXO1 constrains activation and regulates senescence in CD8 T cells. *Cell Rep.* 34 (4), 108674. doi:10.1016/j.celrep.2020.108674

Dimitri, A., Herbst, F., and Fraietta, J. A. (2022). Engineering the next-generation of CAR T-cells with CRISPR-Cas9 gene editing. *Mol. Cancer* 21 (1), 78. doi:10.1186/s12943-022-01559-z

Dutt, S., Tseng, D., Ermann, J., George, T. I., Liu, Y. P., Davis, C. R., et al. (2007). Naive and memory T cells induce different types of graft-versus-host disease. *J. Immunol.* 179 (10), 6547–6554. doi:10.4049/jimmunol.179.10.6547

European Medicines Agency (2018a). Potency testing of cell-based immunotherapy medicinal products for the treatment cancer - scientific guideline. Available at: https://www.ema.europa.eu/en/potency-testing-cell-based-immunotherapy-medicinal-products-treatment-cancer-scientific-guideline (Accessed July 25, 2023).

European Medicines Agency (2018b). Kymriah. Available at: https://www.ema. europa.eu/en/medicines/human/EPAR/kymriah (Accessed June 22, 2023).

European Medicines Agency (2018c). Yescarta. Available at: https://www.ema. europa.eu/en/medicines/human/EPAR/yescarta (Accessed June 22, 2023).

European Medicines Agency (2022). Development and data requirements of potency tests for cell-based therapy products the relation to clinical efficacy - scientific guideline.

Available at: https://www.ema.europa.eu/en/development-data-requirements-potencytests-cell-based-therapy-products-relation-clinical-efficacy (Accessed July 25, 2023).

Eyquem, J., Mansilla-Soto, J., Giavridis, T., Van Der Stegen, S. J. C., Hamieh, M., Cunanan, K. M., et al. (2017). Targeting a CAR to the TRAC locus with CRISPR/ Cas9 enhances tumour rejection. *Nature* 543 (7643), 113–117. doi:10.1038/nature21405

FDA (2011). Guidance for industry potency tests for cellular and gene therapy products. Available at: https://www.fda.gov/media/79856/download.

FDA (2021a). KYMRIAH (tisagenlecleucel). Available at: https://www.fda.gov/ vaccines-blood-biologics/cellular-gene-therapy-products/kymriah-tisagenlecleucel (Accessed October 14, 2021).

FDA (2021b). YESCARTA (axicabtagene ciloleucel). Available at: https://www.fda. gov/vaccines-blood-biologics/cellular-gene-therapy-products/yescarta-axicabtageneciloleucel (Accessed October 14, 2021).

FDA (2023a). BCMA-directed or CD19-directed autologous chimeric antigen receptor (CAR) T cell immunotherapies: FDA safety communication - FDA investigating serious risk of T-cell malignancy. Available at: https://www.fda.gov/safety/medical-product-safety-information/bcma-directed-or-cd19-directed-autologous-chimeric-antigen-receptor-car-t-cell-immunotherapies-fda (Accessed December 23, 2023).

FDA (2023b). Manufacturing changes and comparability for human cellular and gene therapy products; draft guidance for industry. Available at: https://www.fda.gov/media/170198/download.

Feins, S., Kong, W., Williams, E. F., Milone, M. C., and Fraietta, J. A. (2019). An introduction to chimeric antigen receptor (CAR) T-cell immunotherapy for human cancer. *Am. J. Hematol.* 94 (S1), S3–S9. doi:10.1002/ajh.25418

Fraietta, J. A., Lacey, S. F., Orlando, E. J., Pruteanu-Malinici, I., Gohil, M., Lundh, S., et al. (2018). Determinants of response and resistance to CD19 chimeric antigen receptor (CAR) T cell therapy of chronic lymphocytic leukemia. *Nat. Med.* 24 (5), 563–571. doi:10.1038/s41591-018-0010-1

Food and Drug Administration (2019). Potency tests for cellular and gene therapy products. Available at: https://www.fda.gov/regulatory-information/search-fda-guidance-documents/potency-tests-cellular-and-gene-therapy-products (Accessed June 23, 2023).

Gajra, A., Zalenski, A., Sannareddy, A., Jeune-Smith, Y., Kapinos, K., and Kansagra, A. (2022). Barriers to chimeric antigen receptor T-cell (CAR-T) therapies in clinical practice. *Pharm. Med.* 36 (3), 163–171. doi:10.1007/s40290-022-00428-w

Gardner, R. A., Ceppi, F., Rivers, J., Annesley, C., Summers, C., Taraseviciute, A., et al. (2019). Preemptive mitigation of CD19 CAR T-cell cytokine release syndrome without attenuation of antileukemic efficacy. *Blood* 134 (24), 2149–2158. doi:10.1182/blood. 2019001463

Gardner, R. A., Finney, O., Annesley, C., Brakke, H., Summers, C., Leger, K., et al. (2017). Intent-to-treat leukemia remission by CD19 CAR T cells of defined formulation and dose in children and young adults. *Blood* 129 (25), 3322–3331. doi:10.1182/blood-2017-02-769208

Gargett, T., and Brown, M. P. (2015). Different cytokine and stimulation conditions influence the expansion and immune phenotype of third-generation chimeric antigen receptor T cells specific for tumor antigen GD2. *Cytotherapy* 17 (4), 487–495. doi:10. 1016/j.jcyt.2014.12.002

Gargett, T., Truong, N., Ebert, L. M., Yu, W., and Brown, M. P. (2019). Optimization of manufacturing conditions for chimeric antigen receptor T cells to favor cells with a central memory phenotype. *Cytotherapy* 21 (6), 593–602. doi:10. 1016/j.jcyt.2019.03.003

Gattinoni, L., Klebanoff, C. A., and Restifo, N. P. (2012). Paths to stemness: building the ultimate antitumour T cell. *Nat. Rev. Cancer* 12 (10), 671–684. doi:10.1038/nrc3322

Gattinoni, L., Zhong, X. S., Palmer, D. C., Ji, Y., Hinrichs, C. S., Yu, Z., et al. (2009). Wnt signaling arrests effector T cell differentiation and generates CD8+ memory stem cells. *Nat. Med.* 15 (7), 808–813. doi:10.1038/nm.1982

Gauthier, J., and Turtle, C. J. (2018). Insights into cytokine release syndrome and neurotoxicity after CD19-specific CAR-T cell therapy. *Curr. Res. Transl. Med.* 66 (2), 50–52. doi:10.1016/j.retram.2018.03.003

Gee, A. P. (2015). Manufacturing genetically modified T cells for clinical trials. *Cancer Gene Ther.* 22 (2), 67–71. doi:10.1038/cgt.2014.71

Gee, A. P. (2018). GMP CAR-T cell production. Best. Pract. Res. Clin. Haematol. 31 (2), 126–134. doi:10.1016/j.beha.2018.01.002

Ghassemi, S., Durgin, J. S., Nunez-Cruz, S., Patel, J., Leferovich, J., Pinzone, M., et al. (2022). Rapid manufacturing of non-activated potent CAR T cells. *Nat. Biomed. Eng.* 6 (2), 118–128. doi:10.1038/s41551-021-00842-6

Ghassemi, S., Martinez-Becerra, F. J., Master, A. M., Richman, S. A., Heo, D., Leferovich, J., et al. (2020). Enhancing chimeric antigen receptor T cell anti-tumor function through advanced media design. *Mol. Ther. - Methods Clin. Dev.* 18, 595–606. doi:10.1016/j.omtm.2020.07.008

Ghassemi, S., Nunez-Cruz, S., O'Connor, R. S., Fraietta, J. A., Patel, P. R., Scholler, J., et al. (2018). Reducing *ex vivo* culture improves the antileukemic activity of chimeric antigen receptor (CAR) T cells. *Cancer Immunol. Res.* 6, 1100–1109. doi:10.1158/2326-6066.CIR-17-0405 Ghorashian, S., Kramer, A. M., Onuoha, S., Wright, G., Bartram, J., Richardson, R., et al. (2019). Enhanced CAR T cell expansion and prolonged persistence in pediatric patients with ALL treated with a low-affinity CD19 CAR. *Nat. Med.* 25 (9), 1408–1414. doi:10.1038/s41591-019-0549-5

Gogol-Döring, A., Ammar, I., Gupta, S., Bunse, M., Miskey, C., Chen, W., et al. (2016). Genome-wide profiling reveals remarkable parallels between insertion site selection properties of the MLV retrovirus and the piggyBac transposon in primary human CD4+ T cells. *Mol. Ther.* 24 (3), 592–606. doi:10.1038/mt.2016.11

Gong, W., Hoffmann, J. M., Stock, S., Wang, L., Liu, Y., Schubert, M. L., et al. (2019). Comparison of IL-2 vs IL-7/IL-15 for the generation of NY-ESO-1-specific T cells. *Cancer Immunol. Immunother.* 68 (7), 1195–1209. doi:10.1007/s00262-019-02354-4

Graham, C., Jozwik, A., Pepper, A., and Benjamin, R. (2018). Allogeneic CAR-T cells: more than ease of access? *Cells* 7 (10), 155. doi:10.3390/cells7100155

Graham, H., Chandler, D. J., and Dunbar, S. A. (2019). The genesis and evolution of bead-based multiplexing. *Methods* 158, 2-11. doi:10.1016/j.ymeth.2019.01.007

Guedan, S., Posey, A. D., Shaw, C., Wing, A., Da, T., Patel, P. R., et al. (2018). Enhancing CAR T cell persistence through ICOS and 4-1BB costimulation. *JCI Insight* 3 (1), e96976. doi:10.1172/jci.insight.96976

Hardiansyah, D., and Ng, C. M. (2019). Quantitative systems pharmacology model of chimeric antigen receptor T-cell therapy. *Clin. Transl. Sci.* 12 (4), 343–349. doi:10.1111/ cts.12636

Harris, D. T., Hager, M. V., Smith, S. N., Cai, Q., Stone, J. D., Kruger, P., et al. (2018). Comparison of T Cell activities mediated by human TCRs and CARs that use the same recognition domains. *J. Immunol.* 200 (3), 1088–1100. doi:10.4049/jimmunol.1700236

Harrison, S. J., Nguyen, T., Rahman, M., Er, J., Li, J., Li, K., et al. (2023). CAR+ T-cell lymphoma post ciltacabtagene autoleucel therapy for relapsed refractory multiple myeloma. *Blood* 142, 6939. doi:10.1182/blood-2023-178806

Hoffmann, J. M., Schubert, M. L., Wang, L., Hückelhoven, A., Sellner, L., Stock, S., et al. (2018). Differences in expansion potential of naive chimeric antigen receptor T cells from healthy donors and untreated chronic lymphocytic leukemia patients. *Front. Immunol.* 8, 1956. doi:10.3389/fimmu.2017.01956

Holzinger, A., and Abken, H. (2019). CAR T cells: a snapshot on the growing options to design a CAR. *HemaSphere* 3 (1), e172. doi:10.1097/HS9.00000000000172

Hopfinger, G., Jäger, U., and Worel, N. (2019). CAR-T cell therapy in diffuse large B cell lymphoma: hype and hope. *HemaSphere* 3 (2), e185. doi:10.1097/HS9. 00000000000185

Hort, S., Herbst, L., Bäckel, N., Erkens, F., Niessing, B., Frye, M., et al. (2022). Toward rapid, widely available autologous CAR-T cell therapy – artificial intelligence and automation enabling the smart manufacturing hospital. *Front. Med.* 9, 913287. doi:10. 3389/fmed.2022.913287

Huang, J., Huang, X., and Huang, J. (2022). CAR-T cell therapy for hematological malignancies: limitations and optimization strategies. *Front. Immunol.* 13, 1019115. doi:10.3389/fimmu.2022.1019115

IsoPlexis (2023). Polyfunctional Strength offers a new approach to defining quality in cell therapy. Available at: https://isoplexis.com/literature/polyfunctional-strength-offers-a-new-approach-to-defining-quality-in-cell-therapy/ (Accessed August 16, 2023).

Ivics, Z., Hackett, P. B., Plasterk, R. H., and Izsvák, Z. (1997). Molecular reconstruction of sleeping beauty, a Tc1-like transposin from fish, and its transposition in human cells. *Cell* 91 (4), 501–510. doi:10.1016/s0092-8674(00) 80436-5

Jiang, H., Dong, B., Gao, L., Liu, L., Ge, J., He, A., et al. (2020). Clinical results of a multicenter study of the first-in-human dual BCMA and CD19 targeted novel platform fast CAR-T cell therapy for patients with relapsed/refractory multiple myeloma. *Blood* 136, 25–26. doi:10.1182/blood-2020-138614

Jinek, M., Chylinski, K., Fonfara, I., Hauer, M., Doudna, J. A., and Charpentier, E. (2012). A programmable dual-RNA-guided DNA endonuclease in adaptive bacterial immunity. *Science* 337 (6096), 816–821. doi:10.1126/science.1225829

Jo, S., Das, S., Williams, A., Chretien, A. S., Pagliardini, T., Le Roy, A., et al. (2022). Endowing universal CAR T-cell with immune-evasive properties using TALEN-gene editing. *Nat. Commun.* 13 (1), 3453. doi:10.1038/s41467-022-30896-2

Kamali, E., Rahbarizadeh, F., Hojati, Z., and Frödin, M. (2021). CRISPR/Cas9mediated knockout of clinically relevant alloantigenes in human primary T cells. *BMC Biotechnol.* 21 (1), 9. doi:10.1186/s12896-020-00665-4

Kath, J., Du, W., Pruene, A., Braun, T., Thommandru, B., Turk, R., et al. (2022). Pharmacological interventions enhance virus-free generation of TRAC-replaced CAR T cells. *Mol. Ther. Methods Clin. Dev.* 25, 311–330. doi:10.1016/j.omtm.2022.03.018

Kebriaei, P., Singh, H., Huls, M. H., Figliola, M. J., Bassett, R., Olivares, S., et al. (2016). Phase I trials using Sleeping Beauty to generate CD19-specific CAR T cells. *J. Clin. Invest.* 126 (9), 3363–3376. doi:10.1172/JCI86721

Kiesgen, S., Messinger, J. C., Chintala, N. K., Tano, Z., and Adusumilli, P. S. (2021). Comparative analysis of assays to measure CAR T-cell-mediated cytotoxicity. *Nat. Protoc.* 16 (3), 1331–1342. doi:10.1038/s41596-020-00467-0

Kim, M. Y., Jayasinghe, R., Devenport, J. M., Ritchey, J. K., Rettig, M. P., O'Neal, J., et al. (2022). A long-acting interleukin-7, rhIL-7-hyFc, enhances CAR T cell expansion, persistence, and anti-tumor activity. Nat. Commun. 13 (1), 3296. doi:10.1038/s41467-022-30860-0

Kim-Hoehamer, Y. I., Riberdy, J. M., Zheng, F., Park, J. J., Shang, N., Métais, J. Y., et al. (2023). Development of a cGMP-compliant process to manufacture donor-derived, CD45RA-depleted memory CD19-CAR T cells. *Gene Ther.* 30 (3), 222–231. doi:10. 1038/s41434-021-00307-0

Kimmel, G. J., Locke, F. L., and Altrock, P. M. (2020). Response to CAR T cell therapy can be explained by ecological cell dynamics and stochastic extinction events. *Cancer Biol.* doi:10.1101/717074

Kimmel, G. J., Locke, F. L., and Altrock, P. M. (2021). The roles of T cell competition and stochastic extinction events in chimeric antigen receptor T cell therapy. *Proc. R. Soc. B Biol. Sci.* 288 (1947), 20210229. doi:10.1098/rspb.2021.0229

Kirouac, D. C., Zmurchok, C., Deyati, A., Sicherman, J., Bond, C., and Zandstra, P. W. (2023). Deconvolution of clinical variance in CAR-T cell pharmacology and response. *Nat. Biotechnol.* 1–12. doi:10.1038/s41587-023-01687-x

Klebanoff, C. A., Crompton, J. G., Leonardi, A. J., Yamamoto, T. N., Chandran, S. S., Eil, R. L., et al. (2017). Inhibition of AKT signaling uncouples T cell differentiation from expansion for receptor-engineered adoptive immunotherapy. *JCI Insight* 2 (23), e95103. doi:10.1172/jci.insight.95103

Klebanoff, C. A., Gattinoni, L., Torabi-Parizi, P., Kerstann, K., Cardones, A. R., Finkelstein, S. E., et al. (2005). Central memory self/tumor-reactive CD8+ T cells confer superior antitumor immunity compared with effector memory T cells. *Proc. Natl. Acad. Sci.* 102 (27), 9571–9576. doi:10.1073/pnas.0503726102

Kondo, T., Imura, Y., Chikuma, S., Hibino, S., Omata-Mise, S., Ando, M., et al. (2018). Generation and application of human induced-stem cell memory T cells for adoptive immunotherapy. *Cancer Sci.* 109 (7), 2130–2140. doi:10.1111/cas.13648

Lamarche, C., Ward-Hartstonge, K., Mi, T., Lin, D. T. S., Huang, Q., Brown, A., et al. (2023). Tonic-signaling chimeric antigen receptors drive human regulatory T cell exhaustion. *Proc. Natl. Acad. Sci.* 120 (14), e2219086120. doi:10.1073/pnas.2219086120

Le, K., Tan, C., Le, H., Tat, J., Zasadzinska, E., Diep, J., et al. (2020). Assuring clonality on the Beacon digital cell line development platform. *Biotechnol. J.* 15 (1), 1900247. doi:10.1002/biot.201900247

Lee, D. H., Cervantes-Contreras, F., Lee, S. Y., Green, D. J., and Till, B. G. (2018). Improved expansion and function of CAR T cell products from cultures initiated at defined CD4:CD8 ratios. *Blood* 132, 3334. doi:10.1182/blood-2018-99-111576

Li, R., Sahoo, P., Wang, D., Wang, Q., Brown, C. E., Rockne, R. C., et al. (2023b). Modeling interaction of Glioma cells and CAR T-cells considering multiple CAR T-cells bindings. *Immunoinformatics Amst Neth* 9, 100022. doi:10.1016/j.immuno. 2023.100022

Li, W., Guo, L., Rathi, P., Marinova, E., Gao, X., Wu, M. F., et al. (2017). Redirecting T cells to glypican-3 with 4-1BB zeta chimeric antigen receptors results in Th1 polarization and potent antitumor activity. *Hum. Gene Ther.* 28 (5), 437–448. doi:10.1089/hum.2016.025

Li, W., Pan, X., Chen, L., Cui, H., Mo, S., Pan, Y., et al. (2023a). Cell metabolism-based optimization strategy of CAR-T cell function in cancer therapy. *Front. Immunol.* 14, 1186383. doi:10.3389/fimmu.2023.1186383

Liadi, I., Singh, H., Romain, G., Rey-Villamizar, N., Merouane, A., Adolacion, J. R. T., et al. (2015). Individual motile CD4+ T cells can participate in efficient multikilling through conjugation to multiple tumor cells. *Cancer Immunol. Res.* 3 (5), 473–482. doi:10.1158/2326-6066.CIR-14-0195

Lickefett, B., Chu, L., Ortiz-Maldonado, V., Warmuth, L., Barba, P., Doglio, M., et al. (2023). Lymphodepletion – an essential but undervalued part of the chimeric antigen receptor T-cell therapy cycle. *Front. Immunol.* 14, 1303935. doi:10.3389/fimmu.2023. 1303935

Liu, C., Ayyar, V. S., Zheng, X., Chen, W., Zheng, S., Mody, H., et al. (2021). Modelbased cellular kinetic analysis of chimeric antigen receptor-T cells in humans. *Clin. Pharmacol. Ther.* 109 (3), 716–727. doi:10.1002/cpt.2040

Locke, F. L., Ghobadi, A., Jacobson, C. A., Miklos, D. B., Lekakis, L. J., Oluwole, O. O., et al. (2019). Long-term safety and activity of axicabtagene ciloleucel in refractory large B-cell lymphoma (ZUMA-1): a single-arm, multicentre, phase 1–2 trial. *Lancet Oncol.* 20 (1), 31–42. doi:10.1016/S1470-2045(18)30864-7

Logun, M., Colonna, M. B., Mueller, K. P., Ventarapragada, D., Rodier, R., Tondepu, C., et al. (2023). Label-free *in vitro* assays predict the potency of anti-disialoganglioside chimeric antigen receptor T-cell products. *Cytotherapy* 25 (6), 670–682. doi:10.1016/j. jcyt.2023.01.008

Louis, C. U., Savoldo, B., Dotti, G., Pule, M., Yvon, E., Myers, G. D., et al. (2011). Antitumor activity and long-term fate of chimeric antigen receptor-positive T cells in patients with neuroblastoma. *Blood* 118 (23), 6050–6056. doi:10.1182/blood-2011-05-354449

Lu, T. L., Pugach, O., Somerville, R., Rosenberg, S. A., Kochenderfer, J. N., Better, M., et al. (2016). A rapid cell expansion process for production of engineered autologous CAR-T cell therapies. *Hum. Gene Ther. Methods* 27 (6), 209–218. doi:10.1089/hgtb. 2016.120

Lu, Y., Chen, J. J., Mu, L., Xue, Q., Wu, Y., Wu, P. H., et al. (2013). High-throughput secretomic analysis of single cells to assess functional cellular heterogeneity. *Anal. Chem.* 85 (4), 2548–2556. doi:10.1021/ac400082e

Mackay, S., Flynn, B., Morse, K., Ng, C., Gallatin, K., and Paczkowski, P. (2017). Single-cell multiplex proteomics reveals synergistic impact of antigen and rimiduciddependent stimulatory signals on promoting polyfunctional GoCAR-T cells targeting prostate stem cell antigen (PSCA). *Blood* 130, 2281. doi:10.1182/blood.V130.Suppl_1. 2281.2281

Magnani, C. F., Gaipa, G., Lussana, F., Belotti, D., Gritti, G., Napolitano, S., et al. (2020). Sleeping Beauty-engineered CAR T cells achieve antileukemic activity without severe toxicities. *J. Clin. Invest.* 130 (11), 6021–6033. doi:10.1172/JCI138473

Marcucci, K. T., Jadlowsky, J. K., Hwang, W. T., Suhoski-Davis, M., Gonzalez, V. E., Kulikovskaya, I., et al. (2018). Retroviral and lentiviral safety analysis of gene-modified T cell products and infused HIV and Oncology patients. *Mol. Ther.* 26 (1), 269–279. doi:10.1016/j.ymthe.2017.10.012

Maude, S. L., Frey, N., Shaw, P. A., Aplenc, R., Barrett, D. M., Bunin, N. J., et al. (2014). Chimeric antigen receptor T cells for sustained remissions in leukemia. *N. Engl. J. Med.* 371 (16), 1507–1517. doi:10.1056/NEJMoa1407222

Mazinani, M., and Rahbarizadeh, F. (2022). CAR-T cell potency: from structural elements to vector backbone components. *Biomark. Res.* 10 (1), 70. doi:10.1186/s40364-022-00417-w

Melenhorst, J. J., Chen, G. M., Wang, M., Porter, D. L., Chen, C., Collins, M. A., et al. (2022). Decade-long leukaemia remissions with persistence of CD4+ CAR T cells. *Nature* 602 (7897), 503–509. doi:10.1038/s41586-021-04390-6

Mezősi-Csaplár, M., Szöőr, Á., and Vereb, G. (2023). CD28 and 41BB costimulatory domains alone or in combination differentially influence cell surface dynamics and organization of chimeric antigen receptors and early activation of CAR T cells. *Cancers* 15 (12), 3081. doi:10.3390/cancers15123081

Micklethwaite, K. P., Gowrishankar, K., Gloss, B. S., Li, Z., Street, J. A., Moezzi, L., et al. (2021). Investigation of product-derived lymphoma following infusion of piggyBac-modified CD19 chimeric antigen receptor T cells. *Blood* 138 (16), 1391–1405. doi:10.1182/blood.2021010858

Miwa, H., Dimatteo, R., de Rutte, J., Ghosh, R., and Di Carlo, D. (2022). Single-cell sorting based on secreted products for functionally defined cell therapies. *Microsyst. Nanoeng.* 8 (1), 84–18. doi:10.1038/s41378-022-00422-x

Mocciaro, A., Roth, T. L., Bennett, H. M., Soumillon, M., Shah, A., Hiatt, J., et al. (2018). Light-activated cell identification and sorting (LACIS) for selection of edited clones on a nanofluidic device. *Commun. Biol.* 1 (1), 41–48. doi:10.1038/s42003-018-0034-6

Moretti, A., Ponzo, M., Nicolette, C. A., Tcherepanova, I. Y., Biondi, A., and Magnani, C. F. (2022). The past, present, and future of non-viral CAR T cells. *Front. Immunol.* 13, 867013. doi:10.3389/fimmu.2022.867013

Mueller-Schoell, A., Puebla-Osorio, N., Michelet, R., Green, M. R., Künkele, A., Huisinga, W., et al. (2021). Early survival prediction framework in CD19-specific CAR-T cell immunotherapy using a quantitative systems pharmacology model. *Cancers* 13 (11), 2782. doi:10.3390/cancers13112782

Mullard, A. (2023). Tumour-infiltrating lymphocyte cancer therapy nears FDA finish line. *Nat. Rev. Drug Discov.* 23, 3-7. doi:10.1038/d41573-023-00206-6

Muralidharan, S., Hanley, P. J., Liu, E., Chakraborty, R., Bollard, C., Shpall, E., et al. (2011). Activation of Wnt signaling arrests effector differentiation in human peripheral and cord blood-derived T lymphocytes. *J. Immunol.* 187 (10), 5221–5232. doi:10.4049/jimmunol.1101585

Naik, S., and Triplett, B. M. (2023). Selective depletion of naïve T cells by targeting CD45RA. *Front. Oncol.* 12, 1009143. doi:10.3389/fonc.2022.1009143

Nanolive (2023). Nanolive's LIVE T Cell Assay: label-free immuno-therapy analysis. Nanolive - a complete solution for your label-free live cell imaging. Available at: https:// www.nanolive.ch/products/live-cell-analytics/live-t-cell-assay/ (Accessed December 23, 2023).

Noaks, E., Peticone, C., Kotsopoulou, E., and Bracewell, D. G. (2021). Enriching leukapheresis improves T cell activation and transduction efficiency during CAR T processing. *Mol. Ther. - Methods Clin. Dev.* 20, 675–687. doi:10.1016/j.omtm.2021.02.002

Nukala, U., Rodriguez Messan, M., Yogurtcu, O. N., Wang, X., and Yang, H. (2021). A systematic review of the efforts and hindrances of modeling and simulation of CAR T-cell therapy. *AAPS J.* 23 (3), 52. doi:10.1208/s12248-021-00579-9

Onda, M., Kobayashi, K., and Pastan, I. (2019). Depletion of regulatory T cells in tumors with an anti-CD25 immunotoxin induces CD8 T cell-mediated systemic antitumor immunity. *Proc. Natl. Acad. Sci.* 116 (10), 4575–4582. doi:10.1073/pnas. 1820388116

Owens, K., and Bozic, I. (2021). Modeling CAR T-cell therapy with patient preconditioning. Bull. Math. Biol. 83 (5), 42. doi:10.1007/s11538-021-00869-5

Paixão, E. A., Barros, L. R. C., Fassoni, A. C., and Almeida, R. C. (2022). Modeling patient-specific CAR-T cell dynamics: multiphasic kinetics via phenotypic differentiation. *Cancers* 14 (22), 5576. doi:10.3390/cancers14225576

Park, J. H., Geyer, M. B., and Brentjens, R. J. (2016). CD19-targeted CAR T-cell therapeutics for hematologic malignancies: interpreting clinical outcomes to date. *Blood* 127 (26), 3312–3320. doi:10.1182/blood-2016-02-629063

Park, J. H., Rivière, I., Gonen, M., Wang, X., Sénéchal, B., Curran, K. J., et al. (2018). Long-term follow-up of CD19 CAR therapy in acute lymphoblastic leukemia. *N. Engl. J. Med.* 378 (5), 449–459. doi:10.1056/NEJMoa1709919 Philip, M., and Schietinger, A. (2019). Heterogeneity and fate choice: T cell exhaustion in cancer and chronic infections. *Curr. Opin. Immunol.* 58, 98–103. doi:10.1016/j.coi.2019.04.014

Prybutok, A. N., Yu, J. S., Leonard, J. N., and Bagheri, N. (2022). Mapping CAR T-cell design space using agent-based models. *Front. Mol. Biosci.* 9, 849363. doi:10.3389/fmolb.2022.849363

Qi, T., McGrath, K., Ranganathan, R., Dotti, G., and Cao, Y. (2022). Cellular kinetics: a clinical and computational review of CAR-T cell pharmacology. *Adv. Drug Deliv. Rev.* 188, 114421. doi:10.1016/j.addr.2022.114421

Querques, I., Mades, A., Zuliani, C., Miskey, C., Alb, M., Grueso, E., et al. (2019). A highly soluble Sleeping Beauty transposase improves control of gene insertion. *Nat. Biotechnol.* 37 (12), 1502–1512. doi:10.1038/s41587-019-0291-z

Rabinovich, P. M., Komarovskaya, M. E., Wrzesinski, S. H., Alderman, J. L., Budak-Alpdogan, T., Karpikov, A., et al. (2009). Chimeric receptor mRNA transfection as a tool to generate antineoplastic lymphocytes. *Hum. Gene Ther.* 20 (1), 51–61. doi:10.1089/ hum.2008.068

Rajakaruna, H., Desai, M., and Das, J. (2023). PASCAR: a multiscale framework to explore the design space of constitutive and inducible CAR T cells. *Life Sci. Alliance* 6, e202302171. doi:10.26508/lsa.202302171

Riddell, S. R., Sommermeyer, D., Berger, C., Liu, L., Steven), Balakrishnan, A., et al. (2014). Adoptive therapy with chimeric antigen receptor-modified T cells of defined subset composition. *Cancer J.* 20 (2), 141–144. doi:10.1097/PPO.00000000000036

Roddie, C., O'Reilly, M., Dias Alves Pinto, J., Vispute, K., and Lowdell, M. (2019). Manufacturing chimeric antigen receptor T cells: issues and challenges. *Cytotherapy* 21 (3), 327–340. doi:10.1016/j.jcyt.2018.11.009

Rohrs, J. A., Siegler, E. L., Wang, P., and Finley, S. D. (2020). ERK activation in CAR T cells is amplified by CD28-mediated increase in CD3ζ phosphorylation. *iScience* 23 (4), 101023. doi:10.1016/j.isci.2020.101023

Rohrs, J. A., Zheng, D., Graham, N. A., Wang, P., and Finley, S. D. (2018). Computational model of chimeric antigen receptors explains site-specific phosphorylation kinetics. *Biophys. J.* 115 (6), 1116–1129. doi:10.1016/j.bpj.2018.08.018

Romain, G., Strati, P., Rezvan, A., Fathi, M., Bandey, I. N., Adolacion, J. R. T., et al. (2022). Multidimensional single-cell analysis identifies a role for CD2-CD58 interactions in clinical antitumor T cell responses. *J. Clin. Invest.* 132 (17), e159402. doi:10.1172/JCI159402

Rosenberg, S. A., Aebersold, P., Cornetta, K., Kasid, A., Morgan, R. A., Moen, R., et al. (1990). Gene transfer into humans — immunotherapy of patients with advanced melanoma, using tumor-infiltrating lymphocytes modified by retroviral gene transduction. *N. Engl. J. Med.* 323 (9), 570–578. doi:10.1056/NEJM199008303230904

Rossi, J., Paczkowski, P., Shen, Y. W., Morse, K., Flynn, B., Kaiser, A., et al. (2018). Preinfusion polyfunctional anti-CD19 chimeric antigen receptor T cells are associated with clinical outcomes in NHL. *Blood* 132 (8), 804–814. doi:10.1182/blood-2018-01-828343

Roth, T. L., Puig-Saus, C., Yu, R., Shifrut, E., Carnevale, J., Li, P. J., et al. (2018). Reprogramming human T cell function and specificity with non-viral genome targeting. *Nature* 559 (7714), 405–409. doi:10.1038/s41586-018-0326-5

Sadelain, M., Brentjens, R., and Rivière, I. (2009). The promise and potential pitfalls of chimeric antigen receptors. *Curr. Opin. Immunol.* 21 (2), 215–223. doi:10.1016/j.coi. 2009.02.009

Sahoo, P., Yang, X., Abler, D., Maestrini, D., Adhikarla, V., Frankhouser, D., et al. (2020). Mathematical deconvolution of CAR T-cell proliferation and exhaustion from real-time killing assay data. J. R. Soc. Interface 17 (162), 20190734. doi:10.1098/rsif.2019.0734

Saito, M., Xu, P., Faure, G., Maguire, S., Kannan, S., Altae-Tran, H., et al. (2023). Fanzor is a eukaryotic programmable RNA-guided endonuclease. *Nature* 620, 660–668. doi:10.1038/s41586-023-06356-2

Sanyanusin, M., Tudsamran, S., Thaiwong, R., Tawinwung, S., Nishio, N., Takahashi, Y., et al. (2023). Novel xeno-free and serum-free culturing condition to improve piggyBac transposon-based CD19 chimeric antigen receptor T-cell production and characteristics. *Cytotherapy* 25 (4), 397–406. doi:10.1016/j.jcyt.2022.11.009

Sartorius (2023). Application note: sampling of human pathogenic viruses. Available at: https://www.sartorius.com/en/applications/life-science-research/infectious-diseases/environmental-testing/air-monitoring/sampling-of-human-pathogenic-viruses (Accessed August 15, 2023).

Sasawatari, S., Tadaki, T., Isogai, M., Takahara, M., Nieda, M., and Kakimi, K. (2006). Efficient priming and expansion of antigen-specific CD8+ T cells by a novel cell-based artificial APC. *Immunol. Cell Biol.* 84 (6), 512–521. doi:10.1111/j.1440-1711.2006. 01462.x

Sato, K., Kondo, M., Sakuta, K., Hosoi, A., Noji, S., Sugiura, M., et al. (2009). Impact of culture medium on the expansion of T cells for immunotherapy. *Cytotherapy* 11 (7), 936–946. doi:10.3109/14653240903219114

Schmidts, A., Marsh, L. C., Srivastava, A. A., Bouffard, A. A., Boroughs, A. C., Scarfò, I., et al. (2020). Cell-based artificial APC resistant to lentiviral transduction for efficient generation of CAR-T cells from various cell sources. *J. Immunother. Cancer* 8 (2), e000990. doi:10.1136/jitc-2020-000990

Schmueck, M., Fischer, A. M., Hammoud, B., Brestrich, G., Fuehrer, H., Luu, S. H., et al. (2012). Preferential expansion of human virus-specific multifunctional central memory T cells by partial targeting of the IL-2 receptor signaling pathway: the key role of CD4+ T cells. J. Immunol. 188 (10), 5189–5198. doi:10.4049/jimmunol.1103763

Schmueck-Henneresse, M., Omer, B., Shum, T., Tashiro, H., Mamonkin, M., Lapteva, N., et al. (2017). Comprehensive approach for identifying the T cell subset origin of CD3 and CD28 antibody-activated chimeric antigen receptor-modified T cells. *J. Immunol.* 199 (1), 348–362. doi:10.4049/jimmunol.1601494

Scholler, J., Brady, T. L., Binder-Scholl, G., Hwang, W. T., Plesa, G., Hege, K. M., et al. (2012). Decade-long safety and function of retroviral-modified chimeric antigen receptor T cells. *Sci. Transl. Med.* 4 (132), 132ra53. doi:10.1126/scitranslmed.3003761

Selvaggio, G., Parolo, S., Bora, P., Leonardelli, L., Harrold, J., Mehta, K., et al. (2022). Computational analysis of cytokine release following bispecific T-cell engager therapy: applications of a logic-based model. *Front. Oncol.* 12, 818641. doi:10.3389/fonc.2022.818641

Sermer, D., and Brentjens, R. (2019). CAR T-cell therapy: full speed ahead. Hematol. Oncol. 37 (S1), 95-100. doi:10.1002/hon.2591

Shah, N. N., Highfill, S. L., Shalabi, H., Yates, B., Jin, J., Wolters, P. L., et al. (2020a). CD4/CD8 T-Cell selection affects chimeric antigen receptor (CAR) T-cell potency and toxicity: updated results from a phase I anti-CD22 CAR T-cell trial. *J. Clin. Oncol.* 38, 1938–1950. doi:10.1200/JCO.19.03279

Shah, N. N., Johnson, B. D., Schneider, D., Zhu, F., Szabo, A., Keever-Taylor, C. A., et al. (2020b). Bispecific anti-CD20, anti-CD19 CAR T cells for relapsed B cell malignancies: a phase 1 dose escalation and expansion trial. *Nat. Med.* 26 (10), 1569–1575. doi:10.1038/s41591-020-1081-3

Shah, V., Womack, J., Zamora, A. E., Terhune, S. S., and Dash, R. K. (2023). Simulating the evolution of signaling signatures during CART-cell -- tumor cell interactions. ArXiv. 2023 Feb 8;arXiv:2302.04338v1.

Shy, B. R., Vykunta, V., Ha, A., Roth, T. L., Talbot, A., Nguyen, D. N., et al. (2021). Hybrid ssDNA repair templates enable high yield genome engineering in primary cells for disease modeling and cell therapy manufacturing. *Synth. Biol.* doi:10.1101/2021.09.02.458799

Si, X., Xiao, L., Brown, C. E., and Wang, D. (2022). Preclinical evaluation of CAR T cell function: *in vitro* and *in vivo* models. *Int. J. Mol. Sci.* 23 (6), 3154. doi:10.3390/ ijms23063154

Siegler, E. L., and Kenderian, S. S. (2020). Neurotoxicity and cytokine release syndrome after chimeric antigen receptor T cell therapy: insights into mechanisms and novel therapies. *Front. Immunol.* 11, 1973. doi:10.3389/fimmu.2020.01973

Singh, A. P., Zheng, X., Lin-Schmidt, X., Chen, W., Carpenter, T. J., Zong, A., et al. (2020). Development of a quantitative relationship between CAR-affinity, antigen abundance, tumor cell depletion and CAR-T cell expansion using a multiscale systems PK-PD model. *mAbs* 12 (1), 1688616. doi:10.1080/19420862.2019.1688616

Sommermeyer, D., Hudecek, M., Kosasih, P. L., Gogishvili, T., Maloney, D. G., Turtle, C. J., et al. (2016). Chimeric antigen receptor-modified T cells derived from defined CD8+ and CD4+ subsets confer superior antitumor reactivity *in vivo. Leukemia* 30 (2), 492–500. doi:10.1038/leu.2015.247

Soundara Rajan, T., Gugliandolo, A., Bramanti, P., and Mazzon, E. (2020). *In vitro*-Transcribed mRNA chimeric antigen receptor T cell (ivt mRNA CAR T) therapy in hematologic and solid tumor management: a preclinical update. *Int. J. Mol. Sci.* 21 (18), 6514. doi:10.3390/ijms21186514

Stein, A. M., Grupp, S. A., Levine, J. E., Laetsch, T. W., Pulsipher, M. A., Boyer, M. W., et al. (2019). Tisagenlecleucel model-based cellular kinetic analysis of chimeric antigen receptor-T cells. *CPT Pharmacomet. Syst. Pharmacol.* 8 (5), 285–295. doi:10.1002/psp4.12388

Stemberger, C., Dreher, S., Tschulik, C., Piossek, C., Bet, J., Yamamoto, T. N., et al. (2012). Novel serial positive enrichment technology enables clinical multiparameter cell sorting. *PLOS ONE* 7 (4), e35798. doi:10.1371/journal.pone.0035798

Sterner, R. C., and Sterner, R. M. (2021). CAR-T cell therapy: current limitations and potential strategies. *Blood Cancer J.* 11 (4), 69–11. doi:10.1038/s41408-021-00459-7

Stock, S., Schmitt, M., and Sellner, L. (2019). Optimizing manufacturing protocols of chimeric antigen receptor T cells for improved anticancer immunotherapy. *Int. J. Mol. Sci.* 20 (24), 6223. doi:10.3390/ijms20246223

Strietz, J., and Chen, Y. (2022). *CELLCYTE XTM* | *Quantitative live cell imaging assays* for immunology. CYTENA. Available at: https://www.cytena.com/wp-content/uploads/2022/09/2022_09_ImmunologyApp-Note.pdf (Accessed July 26, 2023).

Subklewe, M., von Bergwelt-Baildon, M., and Humpe, A. (2019). Chimeric antigen receptor T cells: a race to revolutionize cancer therapy. *Transfus. Med. Hemotherapy* 46 (1), 15–24. doi:10.1159/000496870

Sartorious (2019). The evolution of immune cell killing assays using live-cell analysis. Available at: https://www.sartorius.com/en/products/live-cell-imaging-analysis/live-cell-analysis-resources/evolution-of-immune-cell-killing-assays-white-paper (Accessed July 26, 2023).

Turtle, C. J., Hanafi, L. A., Berger, C., Hudecek, M., Pender, B., Robinson, E., et al. (2016). Immunotherapy of non-Hodgkin's lymphoma with a defined ratio of CD8+ and CD4+ CD19-specific chimeric antigen receptor-modified T cells. *Sci. Transl. Med.* 8 (355), 355ra116. doi:10.1126/scitranslmed.aaf8621

Tyagarajan, S., Spencer, T., and Smith, J. (2020). Optimizing CAR-T cell manufacturing processes during pivotal clinical trials. *Mol. Ther. - Methods Clin. Dev.* 16, 136–144. doi:10.1016/j.omtm.2019.11.018

Urak, R., Walter, M., Lim, L., Wong, C. W., Budde, L. E., Thomas, S., et al. (2017). *Exvivo* Akt inhibition promotes the generation of potent CD19CAR T cells for adoptive immunotherapy. *J. Immunother. Cancer* 5 (1), 26. doi:10.1186/s40425-017-0227-4

Urbani, J., Rotatori, S., Han, Y. C., Mohammed, C., Wang, D., Liu, X., et al. (2022). Suppressing single cell noise in T cell selection. Available at: https://enrichbiosystems. com/wp-content/uploads/2022/12/Tcell_selection_WhitePaper.pdf.

van der Stegen, S. J. C., Hamieh, M., and Sadelain, M. (2015). The pharmacology of second-generation chimeric antigen receptors. *Nat. Rev. Drug Discov.* 14 (7), 499–509. doi:10.1038/nrd4597

Vormittag, P., Gunn, R., Ghorashian, S., and Veraitch, F. S. (2018). A guide to manufacturing CAR T cell therapies. *Curr. Opin. Biotechnol.* 53, 164–181. doi:10.1016/j. copbio.2018.01.025

Wang, X., and Rivière, I. (2016). Clinical manufacturing of CAR T cells: foundation of a promising therapy. *Mol. Ther. - Oncolytics* 3, 16015. doi:10.1038/mto.2016.15

Weinkove, R., George, P., Dasyam, N., and McLellan, A. D. (2019). Selecting costimulatory domains for chimeric antigen receptors: functional and clinical considerations. *Clin. Transl. Immunol.* 8 (5), e1049. doi:10.1002/cti2.1049

Wherry, E. J., Teichgräber, V., Becker, T. C., Masopust, D., Kaech, S. M., Antia, R., et al. (2003). Lineage relationship and protective immunity of memory CD8 T cell subsets. *Nat. Immunol.* 4 (3), 225–234. doi:10.1038/ni889

Wilson, T. L., Kim, H., Chou, C. H., Langfitt, D., Mettelman, R. C., Minervina, A. A., et al. (2022). Common trajectories of highly effective CD19-specific CAR T cells identified by endogenous T-cell receptor lineages. *Cancer Discov.* 12 (9), 2098–2119. doi:10.1158/2159-8290.CD-21-1508

Xu, H., Chaney, K., Johnson, B., Fenske, T. S., Hamadani, M., Dropulic, B., et al. (2020). Single-cell cytokine analysis of LV20.19 bispecific CAR T-cell products from a phase I clinical trial. *Blood* 136, 22. doi:10.1182/blood-2020-136805

Xu, H., Wang, N., Cao, W., Huang, L., Zhou, J., and Sheng, L. (2018). Influence of various medium environment to *in vitro* human T cell culture. *Vitro Cell Dev. Biol.* - *Anim.* 54 (8), 559–566. doi:10.1007/s11626-018-0273-3

Xu, X., Sun, Q., Liang, X., Chen, Z., Zhang, X., Zhou, X., et al. (2019). Mechanisms of relapse after CD19 CAR T-cell therapy for acute lymphoblastic leukemia and its prevention and treatment strategies. *Front. Immunol.* 10, 2664. doi:10.3389/fimmu.2019.02664

Yang, J., He, J., Zhang, X., Li, J., Wang, Z., Zhang, Y., et al. (2022). Next-day manufacture of a novel anti-CD19 CAR-T therapy for B-cell acute lymphoblastic leukemia: first-in-human clinical study. *Blood Cancer J.* 12 (7), 104–109. doi:10.1038/ s41408-022-00694-6

Yang, J., Jiang, P., Zhang, X., Li, J., Wu, Y., Xu, L., et al. (2020). Successful 24-hours manufacture of anti-CD19/CD22 dual chimeric antigen receptor (CAR) T cell therapy for B-cell acute lymphoblastic leukemia (B-all). *Blood* 136, 2–3. doi:10.1182/blood-2020-136866

Young, R. M., Engel, N. W., Uslu, U., Wellhausen, N., and June, C. H. (2022). Nextgeneration CAR T-cell therapies. *Cancer Discov.* 12 (7), 1625–1633. doi:10.1158/2159-8290.CD-21-1683

Zah, E., Nam, E., Bhuvan, V., Tran, U., Ji, B. Y., Gosliner, S. B., et al. (2020). Systematically optimized BCMA/CS1 bispecific CAR-T cells robustly control heterogeneous multiple myeloma. *Nat. Commun.* 11 (1), 2283. doi:10.1038/s41467-020-16160-5

Zhang, C., He, J., Liu, L., Wang, J., Wang, S., Liu, L., et al. (2022b). Novel CD19 chimeric antigen receptor T cells manufactured next-day for acute lymphoblastic leukemia. *Blood Cancer J.* 12 (6), 96–99. doi:10.1038/s41408-022-00688-4

Zhang, J., Yang, G., Ye, P., Li, N., Abassi, Y., and Lamarche, B. (2021). *Real-time potency assay for CART cell killing of adherent cancer cells*. Agilent Technologies. Available at: https://www.agilent.com/cs/library/technicaloverviews/public/tn-car-t-potency-assay-5994-3555en-agilent.pdf (Accessed July 26, 2023).

Zhang, X., Zhu, L., Zhang, H., Chen, S., and Xiao, Y. (2022a). CAR-T cell therapy in hematological malignancies: current opportunities and challenges. *Front. Immunol.* 13, 927153. doi:10.3389/fimmu.2022.927153

Zhang, Z., Liu, L., Ma, C., and Chen, W. (2022c). A computational model of cytokine release syndrome during CAR T-cell therapy. *Adv. Ther.* 5 (10), 2200130. doi:10.1002/adtp.202200130

Zhao, Y., Moon, E., Carpenito, C., Paulos, C. M., Liu, X., Brennan, A. L., et al. (2010). Multiple injections of electroporated autologous T cells expressing a chimeric antigen receptor mediate regression of human disseminated tumor. *Cancer Res.* 70 (22), 9053–9061. doi:10.1158/0008-5472.CAN-10-2880

Zheng, W., O'Hear, C. E., Alli, R., Basham, J. H., Abdelsamed, H. A., Palmer, L. E., et al. (2018). PI3K orchestration of the *in vivo* persistence of chimeric antigen receptormodified T cells. *Leukemia* 32 (5), 1157–1167. doi:10.1038/s41375-017-0008-6

Zhou, J., Jin, L., Wang, F., Zhang, Y., Liu, B., and Zhao, T. (2019). Chimeric antigen receptor T (CAR-T) cells expanded with IL-7/IL-15 mediate superior antitumor effects. *Protein Cell* 10 (10), 764–769. doi:10.1007/s13238-019-0643-y

Zurko, J. C., Chaney, K., Xu, H., Schneider, D., Hari, P., Johnson, B., et al. (2021). Manufacturing bispecific LV20.19 CAR T-cells with IL-7 and IL-15 for a shorter duration improves CAR T-cell immunophenotype while maintaining target cell dose. *Blood* 138, 3883. doi:10.1182/blood-2021-147405

Zurko, J. C., Xu, H., Chaney, K., Schneider, D., Szabo, A., Hari, P., et al. (2022). Bispecific targeting of CD20 and CD19 increases polyfunctionality of chimeric antigen receptor T-cell products in B-cell malignancies. *Cytotherapy* 24 (8), 767–773. doi:10. 1016/j.jcyt.2022.03.011