

GOPEN ACCESS

Citation: Mól CG, Vieira AGdS, Garcia BMSP, dos Santos Pereira E, Eid RAC, Schultz MJ, et al. (2024) Closed-loop oxygen control for critically ill patients—A systematic review and meta-analysis. PLoS ONE 19(6): e0304745. https://doi.org/ 10.1371/journal.pone.0304745

Editor: Vijay Hadda, AIIMS: All India Institute of Medical Sciences, INDIA

Received: February 5, 2024

Accepted: May 16, 2024

Published: June 12, 2024

Peer Review History: PLOS recognizes the benefits of transparency in the peer review process; therefore, we enable the publication of all of the content of peer review and author responses alongside final, published articles. The editorial history of this article is available here: https://doi.org/10.1371/journal.pone.0304745

Copyright: © 2024 Mól et al. This is an open access article distributed under the terms of the <u>Creative</u> Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Data Availability Statement: All relevant data are within the manuscript and its <u>Supporting</u> Information files.

RESEARCH ARTICLE

Closed-loop oxygen control for critically ill patients—A systematic review and metaanalysis

Caroline Gomes Mól¹[•], Aléxia Gabriela da Silva Vieira¹, Bianca Maria Schneider Pereira Garcia¹, Emanuel dos Santos Pereira¹, Raquel Afonso Caserta Eid¹, Marcus J. Schultz^{2,3,4,5‡}, Ana Carolina Pereira Nunes Pinto^{6‡}, Ricardo Kenji Nawa¹[•]*

1 Department of Critical Care Medicine, Hospital Israelita Albert Einstein, São Paulo, SP, Brazil,

2 Department of Intensive Care, Amsterdam UMC, Location AMC, Amsterdam, The Netherlands,

3 Mahidol–Oxford Tropical Medicine Research Unit (MORU), Mahidol University, Bangkok, Thailand,

4 Nuffield Department of Medicine, University of Oxford, Oxford, United Kingdom, 5 Department of

Anesthesia, General Intensive Care and Pain Management, Medical University Wien, Vienna, Austria,

6 Iberoamerican Cochrane Center, Biomedical Research Institute Sant Pau, Barcelona, Spain

So These authors contributed equally to this work.

‡ These authors also contributed equally to this work.

caroline.rodrigues@eisntein.br (CGM); ricardo.nawa@einstein.br (RKN)

Abstract

Background

The impact of closed-loop control systems to titrate oxygen flow in critically ill patients, including their effectiveness, efficacy, workload and safety, remains unclear. This systematic review investigated the utilization of closed-loop oxygen systems for critically ill patients in comparison to manual oxygen titration systems focusing on these topics.

Methods and findings

A search was conducted across several databases including MEDLINE, CENTRAL, EMBASE, LILACS, CINAHL, LOVE, ClinicalTrials.gov, and the World Health Organization on March 3, 2022, with subsequent updates made on June 27, 2023. Evidence databases were searched for randomized clinical parallel or crossover studies investigating closedloop oxygen control systems for critically ill patients. This systematic review and meta-analysis was performed following the Preferred Reporting Items for Systematic Review and Meta-analysis guidelines. The analysis was conducted using Review Manager software, adopting the mean difference or standardized mean difference with a 95% confidence interval (95% CI) for continuous variables or risk ratio with 95% CI for dichotomous outcomes. The main outcome of interest was the percentage of time spent in the peripheral arterial oxygen saturation target. Secondary outcomes included time for supplemental oxygen weaning, length of stay, mortality, costs, adverse events, and workload of healthcare professional. A total of 37 records from 21 studies were included in this review with a total of 1,577 participants. Compared with manual oxygen titration, closed-loop oxygen control systems increased the percentage of time in the prescribed SpO₂ target, mean difference (MD) **Funding:** The author(s) received no specific funding for this work.

Competing interests: MJS was the team leader of Research and New Technologies at Hamilton Medical AG, Bonaduz, Switzerland, from January 2022 until January 2023. The other authors declare no conflicts of interest. This does not alter our adherence to PLOS ONE policies on sharing data and materials. This does not alter our adherence to PLOS ONE policies on sharing data and materials. 25.47; 95% CI 19.7, 30.0], with moderate certainty of evidence. Current evidence also shows that closed-loop oxygen control systems have the potential to reduce the percentage of time with hypoxemia (MD -0.98; 95% CI -1.68, -0.27) and healthcare workload (MD -4.94; 95% CI -7.28, -2.61) with low certainty of evidence.

Conclusion

Closed-loop oxygen control systems increase the percentage of time in the preferred SpO₂ targets and may reduce healthcare workload.

Trial registration

PROSPERO: CRD42022306033.

Introduction

Critically ill patients frequently require supplemental oxygen due to the clinical manifestation of inadequate gas exchange [1,2]. Oxygen administration can be considered a lifesaving treatment and may reduce the mortality and morbidity of hypoxemic patients [3,4]. Oxygen is commonly utilized for critically ill patients and its prescription is frequently observed for patients undergoing respiratory support, including mechanical ventilation to correct or prevent hypoxemia [3,5,6].

Hypoxemia can result from many factors, such as lung and cardiovascular diseases, or reduced oxygen levels at high altitudes and can potentially progress to hypoxia, leading to potential organ damage [7-10]. In addition, hyperoxemia, the presence of excess oxygen, also carries risks, including oxidative stress and potential harm to vital organs [11-14]. Thus, finding the right balance is crucial [15-20]. Therefore, healthcare professionals must carefully titrate oxygen to ensure that they are adequate for the patient's needs, minimizing the risks associated with both hypoxemia and hyperoxemia while optimizing patient outcomes.

The literature suggests safe and acceptable targets of peripheral capillary oxygen saturation (SpO_2) ranging between 92–98% for patients without lung diseases and 88–92% for patients with previous lung diseases [21]. However, manual adjustment of the fraction of inspired oxygen (FiO₂) to precisely deliver oxygen within the target is both challenging and time-consuming for healthcare professionals [22–24]. Therefore, the adoption of automated technologies may be considered to partially reduce workload of healthcare professionals.

Closed-loop oxygen control systems have been developed to provide continuing monitoring and adjustments of oxygen titration based on patients' SpO₂, avoiding and treating hypoxemia and hyperoxemia episodes [23,24]. These systems are based on the feedback principle to maintain a target saturation prescribed by continuing titrating oxygen levels [22–24]. It remains uncertain whether closed-loop oxygen control systems are effective, efficient, and safe. It is also uncertain whether these systems affect the workload of healthcare providers. We conducted a systematic review and meta-analysis focusing on these topics.

Methods

This review is reported in accordance with the Preferred Reporting Items for Systematic Reviews and Meta-Analyses (PRISMA) statement guidelines [25,26]. This systematic review

was registered on PROSPERO [CRD42022306033], and additional details are available in the published protocol [27].

Search strategy

The search was performed on March 3, 2022, and updated on June 27, 2023, in the MEDLINE, CENTRAL, EMBASE, LILACS, CINAHL, and LOVE evidence databases (**S1–S6 Tables in <u>S1</u> File**). Furthermore, a search was performed at the ClinicalTrials.gov website and World Health Organization (WHO) International Clinical Trials Registry Platform to find *'ongoing'* and *'unpublished'* studies. There were no restrictions to language, date, or type of publication.

Eligibility criteria

The eligibility criteria [28] were (P) population: adult ICU patients requiring supplemental oxygen; (I) intervention: any system or device that allows an automatic oxygen titration, including for use with invasive or noninvasive ventilation or with low- or high-flow oxygen therapy; (C) comparator: manual adjustments of oxygen; and (O) outcome: percentage of time in the SpO_2 target, time for weaning from oxygen, length of stay, mortality, costs, adverse events, and workload of healthcare professionals. Different records identified by the same registration number were grouped and considered as a single study. This way, we guarantee that studies with multiple publications were included just once.

Study selection

Two investigators (C.G.M. and A.G.V.) independently screened the studies retrieved by the searches. A third investigator (R.K.N.) was consulted to resolve potential disagreements, if necessary.

Outcome measures

The primary outcome was the time in the SpO_2 target, defined as the percentage of time in which the SpO_2 values remained in the predefined range (e.g., 92 to 96% or 88 to 92%). Secondary outcomes included time for weaning from supplemental oxygen, defined as the total time spent in oxygen support during hospitalization, length of stay, mortality, costs, adverse events, and workload of healthcare professionals, defined as the need for the professional to provide oxygen adjustments (e.g., number of manual adjustments; time spent providing adjustments).

Assessment of characteristics of studies

Characteristics and outcome data from the included studies were independently extracted by two investigators (C.G.M. and A.G.S.) and revised by a third investigator (A.C.P.) using a predefined data collection form.

Risk of bias

The risk of bias for the outcomes was determined by using the 'Cochrane Risk of Bias 2' (RoB2) tool for randomized and crossover trials [29,30]. The risk of bias arising from the randomization process, deviations from intended interventions, missing outcome data, measurement of the outcome and selection of the reported result were assessed. For crossover studies the risk of bias arising from period and carry-over effects was also assessed.

Data synthesis and analysis

The mean difference or standardized mean difference was adopted to analyze continuous variables, with a 95% confidence interval [95% CI]. Dichotomous outcomes are presented as risk ratios (RRs) with 95% CIs. When possible, skewed data were adjusted for mean and standard deviation using Wan's methods and the RevMan Calculator [31]. For figures with good resolution, we extracted data using the '*Webplotdigitizer*' website. For crossover studies, the first period was used to perform the analysis to avoid carry-over effects. When substantial heterogeneity was identified ($I^2 \ge 50\%$), we conducted a predefined subgroup analysis for the type of devices. All analyses were performed using Review Manager software (RevMan, Version 5.4.1. Copenhagen: The Nordic Cochrane Centre, The Cochrane Collaboration, 2020).

Assessment of certainty of evidence

The Grading of Recommendations Assessment, Development and Evaluation (GRADE) system was used to measure and summarize the overall certainty of the current evidence of each outcome [32] through the '*GRADEpro Guideline Development Tool*' software [33].

Post hoc analysis

Two additional outcomes were added after peer review of the published protocol [27] to provide additional information to analyze the data along with the percentage of time in the SpO_2 target, as follows: (a) the percentage of time with hyperoxemia and (b) the percentage of time with hypoxemia. The criteria for hypoxemia and hyperoxemia were considered based on the definitions provided in the studies. We were concerned with carefully verifying the direction of effect and the impact of these data in meta-analyses, describing the contribution of these trials when necessary. In addition, we conducted two posthoc analyses, one comparing the percentage of time in the SpO2 target for subgroups according to the duration of the intervention, and one in subgroups according to the reason for admission, i.e., medical or surgical.

Results

Studies included

A total of 14,256 records were identified, and 37 records from 21 studies [34–54] enrolling 1,577 participants were included (Fig 1 and Table 1 and S7–S9 Tables in S1 File). Of those, 13 studies [34–38,40,42,45,51,52] investigated closed-loop oxygen control systems for use with invasive mechanical ventilation, and 8 studies [41,43,44,46–50] investigated closed-loop oxygen control systems for use with noninvasive respiratory support. Nine studies [36,37,40,44,46–49,54] could be used in the meta-analysis. Further details can be found in the S10 and S11 Tables in S1 File.

Risk of bias for randomized clinical trials

Most of the outcomes presented some concerns [40,42,46,47,50,51,54] mainly due to insufficient information on the randomization process such as allocation concealment and baseline characteristics (Fig 2). Only for costs [48], all domains presented a low risk of bias due to the completed and desirable information provided. For adverse events, length of stay in the intensive care unit (ICU), mortality and time to oxygen weaning, in addition to concerns due to the lack of information from the randomization process, the lack of details on the pre-specified analysis defined a priori motivated the judgment of some concerns among the evaluators [47,50,51,54]. Of note, for adverse events from the seven studies that evaluated this outcome,



Fig 1. Flow of trials through the review. *Abbreviations*: CENTRAL = Cochrane Central Register of Controlled Trials; CINAHL = Cumulative Index to Nursing and Allied Health Literature; EMBASE = Excerpta Medica dataBASE; LILACS = Latin American the Caribbean Literature in Health Sciences; MEDLINE = Medical Literature Analysis and the Retrieval System Online; PICO = Patient Intervention Comparator Outcome.

https://doi.org/10.1371/journal.pone.0304745.g001

Table 1. Summary of included studies.

Author	Study type	Participants	Intervention					
(year) Country			Exp	Con				
Arnal et al.[34] (2010) France	RCT crossover	 n = 43 Age, mean (SD), y = 64 (15) Gender = NR Eligibility = NR Exclusion criteria = NR 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = 2-hours Follow-up = NR 	 IMV with manual control of oxygen (ASV mode) Duration = 2-hours Follow-up = NR 				
Arnal et al.[35] (2012) France	RCT crossover	 n = 50 Age, mean (SD), y = 65 (16) Gender, n (%) = M: 32 (64%) F: 18 (36%) Eligibility = invasively and passively ventilated patients with moderate severity of ARF. Exclusion criteria = increased intracranial pressure, severe ARDS requiring a permissive hypercapnia strategy, chronic respiratory failure receiving long-term oxygen therapy and/or home non-invasive ventilation, severe cardiac arrhythmia, therapeutic hypothermia, brain-dead patients, pregnancy and bronchopleural fistula. 	 IMV with closed-loop control of oxygen (INTELLIVENT[®]-ASV) Duration = 2-hours Follow-up = NR 	 IMV with manual control of oxygen (ASV mode) Duration = 2-hours Follow-up = NR 				
Arnal et al.[36] (2018) France	RCT	• n = 60 • Age, mean (SD), $y = \text{Exp } 67$ (12); Con 62 (16) • Gender, n (%) = Exp M: 20 (66.7%), F: 10 (33.3%); Con M: 22 (73.3%), F: 8 (26.7%) • Eligibility = MV patients (expected duration \geq 48 h) • Exclusion criteria = brain injured patients with GCS < 6, broncho-pleural fistula, chronic or acute dyshemoglobinemia, ventilation drive disorder such as Cheynes Stokes breathing, moribund patient, patient with "Do Not Resuscitate" order before inclusion, chronic respiratory disease requiring long term home ventilation, patient ventilated for planned surgery, pregnancy, lack of informed consent, and the participation in another clinical study.	 IMV with closed-loop control of oxygen (INTELLIVENT[®]-ASV) Duration = 4-days Follow-up = 4-days 	 IMV with manual control of oxygen Duration = 4-days Follow-up = 4-days 				
Bialais et al.[37] (2016) Belgium	RCT	 n = 70 Age, median [IQR], y = Exp 59 [47–69], Con 58 [44–71] Gender, n = Exp: M: 25; F: 17; Con M: 22; F: 16 Eligibility = MV patients (duration ≥ 48 h) and BMI < 40 kg/m². Exclusion criteria = ventricular assistance with intra-aortic balloon counterpulsation, bronchopleural fistula, CO₂ gradient higher than 15 mmHg, pregnancy, adults under guardianship, people deprived of freedom and inclusion in another study protocol. 	 IMV with closed-loop control of oxygen (INTELLIVENT[®]-ASV) Duration = 48-hours Follow-up = NR 	 IMV with manual control of oxygen Duration = 48-hours Follow-up = NR 				
Buiteman- Kruizinga et al. [<u>53]</u> (2022) Netherlands	RCT crossover	 n = 18 Age, y = NR Gender = NR Eligibility = patients within 24-hour of start of ventilation and expected to need MV > 24 hours. Exclusion criteria = MV through a tracheostomy cannula, BMI > 40, any contraindication for use of the automated ventilation mode. 	 IMV with closed-loop control of oxygen (INTELLiVENT-ASV) Duration = 3-hours Follow-up = NR 	 IMV with manual control of oxygen Duration = 3-hours Follow-up = NR 				

Author	Study type	Participants	Intervention				
(year) Country			Exp	Con			
Chelly et al. ³⁸ (2020) France	RCT crossover	 n = 265 Age, mean (SD), y = 64 (14) Gender, n (%) = M: 172 (65); F: 93 (35) Eligibility = MV for at least 48 h, with a FiO₂ ≤ 60%. Exclusion criteria = prone positioning, use of neuromuscular blocking agents, pregnant women, contraindication to automated ventilation mode (delirium, broncho-pleural fistula, respiratory drive disorder such as Cheyne-Stokes breathing), and low-quality measurement for SpO₂. 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = 6-hours Follow-up = NR 	 IMV with manual control of oxygen Duration = 6-hours Follow-up = NR 			
Clavieras et al. [<u>39]</u> (2013) France	RCT crossover	 n = 7 Age, median [IQR], y = 58 [50–64] Gender, n = M: 10; F: 4 Eligibility = MV patients (spontaneous mode), with an expected duration of MV > 48 h. Exclusion criteria = clinical instability, decision to withhold life-sustaining treatment, pregnant women. 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = 24-hours Follow-up = NR 	 IMV with manual control of oxygen Duration = 24-hours Follow-up = NR 			
De Bie et al.[<u>40</u>] (2020) Netherlands	RCT	 n = 220 Age, median [IQR], y = Exp 70 [62–76], Con 70 [63–76] Gender, n (%) = Exp M: 73 (67), F: 36 (33); Con M: 80 (72.1), F: 31 (27.9) Eligibility = elective cardiac surgery requiring IMV in the ICU. Exclusion criteria = before surgery: BMI >35 kg m², history of pneumonectomy or lobectomy, COPD (Global Initiative for Chronic Obstructive Lung Disease Class III or IV), already enrolled in another interventional trial. After cardiac surgery: ECMO, hemodynamic instability, fast-track cardiac surgery. 	 IMV with closed-loop control of oxygen (INTELLIVENT[®]-ASV) Duration = 3-hours Follow-up = 30-day 	 IMV with manual control of oxygen Duration = 3-hours Follow-up = 30-day 			
Denault et al. [41] (2020) Canada	RCT crossover	 n = 30 Age, mean (SD), y = 62 (9) Gender, n (%) = M: 23 (76.7), F: 7 (23.3) Eligibility = patients with mild and moderate obesity who underwent CABG as sole procedure. Exclusion criteria = COPD, restrictive syndrome other than one associated with obesity, or obstructive sleep apnea requiring CPAP to focus on obesity as a sole potential risk factor for hyperoxia-induced hypercapnia. 	 NIV closed-loop control of oxygen (FreeO2 system) with a SpO₂ target of 90 ±2% (conservative) Duration = 30-min Follow-up = NR 	 Noninvasive manual control of oxygen with a SpO₂ target of 95% (liberal) Duration = 30-min Follow-up = NR 			
Eremenko et al. [42] (2021) Russia	RCT	 n = 80 Age, mean (SD), y = Exp 59 (8), Con 59.2 (10.8) Gender, n (%) = M: 57 (71.2), F: 23 (28.8) Eligibility = aged 30–76 years; uncomplicated cardiac surgery; BMI: 18–35 kg/m² Exclusion criteria = before surgery: severe renal, hepatic or cardiac failure. After surgery: acute myocardial infarction, hemodynamic instability, PaO₂/FiO₂ < 150, allergic reaction, seizures, delirium, acute stroke. 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = total duration of MV support Follow-up = hospital length of stay 	 IMV with manual control of oxygen Duration = total duration of MV support Follow-up = hospital length of stay 			

Author	Study type	Participants	Intervention					
(year) Country			Exp	Con				
Hansen et al.[<u>43]</u> (2018) Denmark	RCT crossover	 n = 19 Age, mean (SD), y = 72.4 (8.4) Gender, n (%) = M: 12 (60), F: 8 (40) Eligibility = COPD exacerbation patients with an estimated LOS > 48 h and SpO₂ ≤ 88% on room air. Exclusion criteria = other significant respiratory or cardiac conditions causing hypoxemia, severe ongoing malignancy, high risk for need of mechanical ventilation. 	 Nasal cannula with closed-loop control of oxygen (O2matic system) Duration = 4-hours Follow-up = NR 	 Nasal cannula with manual control of oxygen Duration = 4-hours Follow-up = NR 				
Harper et al.[44] (2021) New Zealand	RCT	 n = 20 Age, mean (SD), y = Exp 68.7 (17), Con 69.2 (19) Gender, n (%) = Exp M: 6 (60) F: 4 (40), Con M: 6 (60) F: 4 (40) Eligibility = patients with HFNC with a prescribed SpO₂ target ranges from (88–92%) or (92–96%). Exclusion criteria = incapability of brief interruption in oxygen or high-flow therapy, domiciliary use of CPAP or NIV, obstructive sleep apnoea, respiratory infection or colonization with multidrug resistant bacteria, infection with SARS-CoV-2, haemodynamic instability, end of life care, risk of barotrauma, as assessed by the investigator, nasal or facial conditions precluding use of HFNC, intracranial trauma or trans-nasal neurosurgery (within 6 weeks), any condition which limits the feasibility of continuous SpO2 monitoring, pregnancy or breastfeeding, cognitive impairment or impaired consent, implanted electronic medical device, any other condition which, at the investigator's discretion, is believed may present a safety risk or impact the feasibility of the study or the study results. 	 HFNC with closed-loop control of oxygen (Airvo3 system) Duration = 24-hours Follow-up = NR 	 HFNC with manual control of oxygen (Airvo2 system) Duration = 24-hours Follow-up = NR 				
Huynh Ky et al. [<u>50]</u> (2017) Canada	RCT	 n = 60 Age, mean (SD), y = 63 (12) Gender, n = M: 44, F: 16 Eligibility = patients with acute coronary syndrome. Exclusion criteria = severe COPD diagnosis. 	 Exp 1: NIV closed-loop control of oxygen (FreeO₂ system) with 92% of SpO₂ target Exp 2: NIV closed-loop control of oxygen (FreeO2 system) with 96% of SpO₂ target Duration = 24-hours Follow-up = NR 	 Noninvasive manual control of oxygen (FreeO₂ system) Duration = 24-hours Follow-up = NR 				
Johannigman et al.[45] (2009) USA	RCT crossover	 n = 15 Age, mean (SD), y = 35.6 (13.2) Gender, n (%) = M: 13 (86.7), F: 2 (13.3) Eligibility = trauma patients requiring IMV; current FiO₂ >35%; indwelling arterial line; and ability to monitor SpO₂ accurately Exclusion criteria = NR 	 IMV with closed-loop control of oxygen (Impact 754, Eagle) Duration = 4-hours Follow-up = 8-hours 	 IMV with manual control of oxygen (Impact 754, Eagle) Duration = 4-hours Follow-up = 8-hours 				
Kobayashi et al. [51] (2017) Japan	RCT	 n = 48 Age, y = NR Gender = NR Eligibility = patients undergoing cardiac surgery. Exclusion criteria = NR 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = 6-hours Follow-up = NR 	 IMV with manual control of oxygen (ASV mode) Duration = 6-hours Follow-up = NR 				
Komnov et al. [54] (2023) Russian	RCT	 n = 32 Age, y = NR Gender = NR Eligibility: Adult patients (age >30 years), elective cardiac surgery, BMI > 35 kg/m² and postoperative treatment with MV Exclusion criteria = NR 	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = total duration of MV support Follow-up = hospital LOS 	 IMV with manual control of oxygen Duration = total duration of MV support Follow-up = hospital LOS 				

Author	Study type	Participants	Intervention				
(year) Country			Exp	Con			
L'Her et al.[46] (2017) France	RCT	• n = 187 • Age, mean (SD), $y = \text{Exp } 74.6$ (13.2), Con 77.7 (12.2) • Gender, n (%) = Exp M: 56 (60.9), F: 37 (39.1); Con M: 48 (51.1), F: 46 (48.9) • Eligibility = emergency department admission, within < 2 h for acute respiratory disorder requiring \ge 3 L/min of O ₂ (SpO ₂ target \ge 92%). • Exclusion criteria = life-threatening hypoxaemia, clinical signs of ventilatory assistance requirement, emergent surgery or coronary angiography, pregnant or breastfeeding women, patients under administrative protective measures.	 Mask to administer either low or high O₂ flow with closed-loop control of oxygen (FreeO2 system) Duration = 3-hours Follow-up = hospital length of stay 	 Mask to administer either low or high O₂ flow with manual control of oxygen Duration = 3-hours Follow-up = hospital length of stay 			
L'Her et al.[<u>47</u>] (2021) France	RCT	 n = 198 Age, mean (SD), y = Exp 62.2 (12.6), Con 63.1 (13.5) Gender, (%) = Exp M: 55.2, F: 44.8; Con M: 55.3, F: 44.7 Eligibility = patients undergoing abdominal or thoracic surgery; ARISCAT ≥ 26; indication of ventilatory support after surgery. Exclusion criteria = BMI ≥ 35; obstructive sleep apnea and pregnant patients. 	 Nasal cannula/prongs (low flow) or open face mask (low or medium flow) with closed-loop control of oxygen (FreeO2 system) Duration = 3-hours/3-days Follow-up = 28-days 	 Nasal cannula/prongs (low flow) or open face mask (low or medium flow) with manual control of oxygen Duration = 3-hours/3-days Follow-up = 28-days 			
Lellouche et al. [52] (2013) Canada	RCT	• n = 60 • Age, mean (SD), $y = \text{Exp } 67.7$ (8.8), Con 65.3 (9.7) • Gender, n (%) = Exp M: 23 (76.7), F: 7 (23.3); Con M: 22 (73.3), F: 8 (26.7) • Eligibility = pre inclusion criteria were: (1) elective cardiac surgery; (2) age 18–90 years; (3) BMI < 40 kg/m ² ; (4) baseline PaCO ₂ < 50 mmHg; and (5) serum creatinine < 200 µmol/L. After surgery, hemodynamically stable, and urine output > 50 mL/ h. • Exclusion criteria = bronchopleural fistula, study ventilator not available.	 IMV with closed-loop control of oxygen (INTELLiVENT[®]-ASV) Duration = 4-hours Follow-up = ICU length of stay 	 IMV with manual control of oxygen Duration = 4-hours Follow-up = ICU length of stay 			
Lellouche et al. [<u>48]</u> (2016) Canada	RCT	 n = 50 Age, mean (SD), y = Exp 71 (8), Con 73 (8) Gender, n (%) = Exp M: 15 (60), F: 10 (40); Con M: 12 (48), F: 13 (52) Eligibility = hospitalized COPD exacerbation patients with oxygen therapy prescription; age ≥ 40 years; and past or current smoking history of at least 10 pack-years. Exclusion criteria = admitted for > 24 hours, multidrug-resistant bacteria infection, intermittent noninvasive ventilation, cognitive impairment. 	 Nasal cannula or nonocclusive mask with closed-loop control of oxygen (FreeO2 system) Duration = 7-days Follow up = 180-days 	 Nasal cannula or nonocclusive mask with manual control of oxygen Duration = 7-days Follow up = 180-days 			

Author	Study type	Participants	Intervention				
(year) Country			Exp	Con			
Roca et al.[49] (2022) pain	RCT crossover	 n = 45 Age, mean (SD), y = 35.6 (13.2) Gender, n (%) = M: 24 (53.3), F: 21 (46.7) Eligibility = ICU admission receiving HFNC with flow rate ≥ 40 L/min with FiO₂ ≥ 30%, and expected to receive HFNC for at least 8 hours after randomization. Exclusion criteria = pregnant women, indication or high risk for immediate intubation, indication for non-invasive ventilation, haemodynamic instability, severe acidosis, poor SpO₂ signal, chronic or acute dyshaemoglobinaemia, tracheostomy, previously enrolled in this trial or enrolled in another interventional study could not participate. 	 HFNC with closed-loop control of oxygen Duration = 4-hours Follow-up = NR 	 HFNC with manual control of oxygen Duration = 4-hours Follow-up = NR 			

Abbreviations: ALI = acute lung injury; ARDS = acute respiratory distress syndrome; ARF = acute respiratory failure; ARISCAT = assess respiratory risk in surgical patients in Catalonia; ASV = adaptive support ventilation; BMI = body mass index; bpm = beats per minute; CABG = coronary artery bypass grafting; CO₂ = carbon dioxide; Con = control; COPD = chronic obstructive pulmonary disease; CPAP = continuous positive airway pressure; DNP = daily nursing procedures; EtCO₂ = end-tidal carbon dioxide; Exp = experimental; F = female; FiO₂ = inspired fraction of oxygen; GCS = Glasgow coma scale; h = hour; HFNC = high flow nasal cannula; ICU = intensive care unit; IMV = invasive mechanical ventilation; kg = kilogram M = male; Max = maximum; mL = milliliter; Min = minimum; MV = mechanical ventilation; NIV = noninvasive ventilation; NHF = nasal high flow; NR = not reported; O₂ = oxygen; PBW = predicted body weight; P_{INSP} = inspiratory pressure; P/F ratio = partial pressure of oxygen/inspired fraction of oxygen ratio; P_{PEAK} = peak pressure; P_{PLAT} = plateau pressure; postop = postoperative; RCT = randomized controlled trial; RR = respiratory rate; ROX = respiratory rate oxygenation; SD = standard deviation; SpO₂ = peripheral oxygen saturation; Vt = tidal volume.

https://doi.org/10.1371/journal.pone.0304745.t001

in addition to the concerns cited about the randomization process, one study [51] lacks information from blinding of outcome assessors.

Risk of bias for randomized cross-over trials

Adverse event is the only outcome that included five studies and only one of them [43] had an overall 'low risk' of bias (Fig 3). We were concerned about the actual impact of the insufficient information available to assess both proper randomization [34,38,39,45,49,53] and appropriate carry-over time in all outcomes [38,39,45,49,35]. As information from the first phase was not available for most of the crossovers, it was not possible to precisely assess whether there was a deviation from the intended interventions [35,38,43]. We did not find any evidence of missing outcome data or differences between groups for 'measure of outcome' for all outcomes planned in this review and thus considered it as low risk of bias. Furthermore, it was not possible to assess whether a selection of the reported result of the outcomes was due to the absence or insufficient reporting information from the registered protocol [35,45].

Time spent in predefined SpO₂ targets, percentage of time with hypoxemia or hyperoxemia, and weaning

The included studies had different SpO₂ targets (S1 Fig in S1 File). In 7 studies [37,40,44,46-49], closed-loop oxygen titration devices increased the percentage of time in the predefined SpO₂ target with substantial heterogeneity, not completely explained by the planned and post hoc subgroup analysis (Fig 4, S2 and S3 Figs and S11 Table in S1 File). The sensitivity analysis for the percentage of time in the SpO₂ target, included only trials with a low risk of bias (S4 Fig in S1 File) and showed similar estimates of the intervention and inconsistency. There's no

Study ID	Experimental	Comparator	Outcome	<u>D1</u>	<u>D2</u>	<u>D3</u>	<u>D4</u>	<u>D5</u>	<u>Overall</u>
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	Adverse Events	+	+	+	+	+	+
Eremenko (2021)	Intellivent-ASV®	Conventional ventilatory modes	Adverse Events	!	+	+	+	+	!
Harper (2021)	Airvo 3 System	Manual oxygen administration	Adverse Events	+	+	•	+	+	+
Kobayashi (2017)	Intellivent-ASV®	ASV	Adverse Events	!	•	+	•	•	
Lellouche (2016)	FreeO ₂	Manual oxygen titration	Adverse Events	+	+	+	+	+	+
L'Her (2017)	FreeO ₂	Manual oxygen titration	Adverse Events	!	+	•	+	+	!
L'Her (2021)	FreeO ₂	Manual oxygen titration	Adverse Events	!	•	+	+	•	!
Poder (2018)	FreeO ₂	Manual oxygen titration	Cost	+	+	+	+	+	+
Arnal (2018)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	+	+	+	+	+	+
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	+	+	+	+	+	+
Eremenko (2021)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	!	+	+	+	+	!
Harper (2021)	Airvo 3 System	Manual oxygen titration	Healthcare workload	+	+	+	+	+	+
Komnov (2023)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	!	+	+	+	+	!
Lellouche (2013)	Intellivent-ASV®	Protocolized ventilation group	Healthcare workload	+	+	+	+	+	+
Harper (2021)	Airvo 3 System	Manual oxygen administration	Hospital length of stay	+	+	+	+	+	+
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	Hospital length of stay	+	+	+	+	+	+
Arnal (2018)	Intellivent-ASV®	Conventional ventilatory modes	Hospital length of stay	+	•	+	+	+	+
Lellouche (2016)	FreeO2	Manual oxygen titration	Hospital length of stay	+	+	+	+	+	+
L'Her (2017)	FreeO ₂	Manual oxygen titration	Hospital length of stay	!	+	+	+	+	!
L'Her (2021)	FreeO ₂	Manual oxygen titration	Hospital length of stay	!	+	+	+	+	!
Eremenko (2021)	Intellivent-ASV®	Conventional ventilatory modes	ICU and Hospital length of stay	!	•	•	•	+	!
Arnal (2018)	Intellivent-ASV®	Conventional ventilatory modes	ICU length of stay	+	+	+	+	+	+
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	ICU length of stay	+	+	+	+	+	+
De Bie (2020)	Intellivent-ASV®	Conventional ventilatory modes	ICU length of stay	!	+	+	+	+	!
Komnov (2023)	Intellivent-ASV®	Conventional ventilatory modes	ICU length of stay	!	+	•	+	•	!
L'Her (2017)	FreeO2	Manual oxygen titration	ICU length of stay	!	•	+	•	•	!
Arnal (2018)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	+	+	+	+	+	+
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	+	•	+	+	+	+
De Bie (2020)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	!	+	+	+	+	!
Eremenko (2021)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	!	•	+	•	•	!
Huynh-Ky (2017)	FreeO ₂	Manual oxygen titration	Mortality	!	+	+	+	•	!
Komnov (2023)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	!	+	+	+	•	!
L'Her (2017)	FreeO ₂	Manual oxygen titration	Mortality	!	+	•	+	+	!
L'Her (2021)	FreeO ₂	Manual oxygen titration	Mortality (28 days)	!	+	•	•	+	!
Lellouche (2013)	Intellivent-ASV®	Protocolized ventilation group	Not acceptable zone	+	+	+	+	+	+
Lellouche (2016)	FreeO2	Manual oxygen titration	Time for O2 weaning	+	+	+	+	+	+
L'Her (2017)	FreeO ₂	Manual oxygen titration	Time for O2 weaning	!	+	+	+	+	!
L'Her (2021)	FreeO ₂	Manual oxygen titration	Time for O2 weaning	!	•	+	•	+	!
De Bie (2020)	Intellivent-ASV®	Conventional ventilatory modes	Time in the optimal zone until extubation	!	+	+	+	+	!
Bialais (2016)	Intellivent-ASV®	Conventional ventilatory modes	Time in the SpO2 target	+	+	+	+	+	+
Harper (2021)	Airvo 3 System	Manual oxygen titration	Time in the SpO2 target	+	+	+	+	+	+
Lellouche (2016)	FreeO ₂	Manual oxygen titration	Time in the SpO2 target	+	+	•	•	+	+
L'Her (2017)	FreeO ₂	Manual oxygen titration	Time in the SpO2 target	!	+	+	+	+	!
L'Her (2021)	FreeO2	Manual oxygen titration	Time in the SpO2 target	!	+	+	+	+	
Lellouche (2013)	Intellivent-ASV®	Protocolized ventilation group	Time in the SpO2 target	+	+	+	•	+	+

Fig 2. Risk of bias of RCT included studies. D1 = randomization process; D2 = deviations from the intended interventions; D3 = missing outcome data; D4 = measurement of the outcome; D5 = selection of the reported results. *Abbreviations*: ASV = adaptive support ventilation; ICU = intensive care unit; PSV = pressure support ventilation; O2 = Oxygen; RoB2 = cochrane risk of bias 2.0 tool for randomized clinical trials; SpO2 = peripheral oxygen saturation.

https://doi.org/10.1371/journal.pone.0304745.g002

+ Low risk ! Some concerns - High risk

Study ID	Experimental	Comparator	Outcome	<u>D1</u>	<u>DS</u>	<u>D2</u>	<u>D3</u>	<u>D4</u>	<u>D5</u>	Overall
Arnal (2012)	Intellivent-ASV®	Conventional ventilatory modes	Adverse Events	+		!	+	+	!	!
Chelly (2020)	Intellivent-ASV®	Conventional ventilatory modes	Adverse Events	!	+	!	+	+	+	!
Clavieras (2013)	Intellivent-ASV®	Conventional ventilatory modes	Adverse Events	!	+	+	+	+	•	!
Hansen (2018)	O2matic	Manual oxygen titration	Adverse Events	+	+	•	+	•	•	+
Johannigman (2009)	Autonomous control of FiO2	Manual oxygen titration	Adverse Events	!	!	+	+	+		!
Clavieras (2013)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	!	+	•	+	+	+	!
Roca (2022)	HFNC with automatic oxygen control	Manual oxygen titration	Healthcare workload	!	+	+	+	+	•	!
Arnal (2012)	Intellivent-ASV®	Conventional ventilatory modes	Healthcare workload	+	!	!	+	+	!	!
Clavieras (2013)	Intellivent-ASV®	Conventional ventilatory modes	Mortality	!	+	+	+	+	+	!
Johannigman (2009)	Autonomous control of FiO2	Manual oxygen titration	Mortality	!	!	•	+	+	!	!
Roca (2022)	HFNC with automatic oxygen control	Manual oxygen titration	Oxygen sonsumption	!	+	+	+	+	•	!
Johannigman (2009)	Autonomous control of FiO2	Manual oxygen titration	Oxygen usage	!	!	+	+	+	!	!
Roca (2022)	HFNC with automatic oxygen control	Manual oxygen titration	Time in the SpO2 target	!	•	+	+	+	+	!
Denault (2020)	FreeO ₂	Manual oxygen titration	Time in the SpO2 target	+	!	+	+	+	+	!
Johannigman (2009)	Autonomous control of FiO2	Manual oxygen titration	Time in the SpO2 target		•	•	•	•	!	!
Chelly (2020)	Intellivent-ASV®	Conventional ventilatory modes	Time in the SpO2 target		+	!	+	+	+	!
Hansen (2018)	O2matic	Manual oxygen titration	Time in the SpO2 target	+	+	!	+	+	+	!
Arnal (2010)	Intellivent-ASV®	ASV	Time in the SpO2 target		+	!	+	+	+	!
Buiteman-Kruizinga (2022)	Intellivent-ASV®	Conventional ventilatory modes	Time in the SpO2 target	!	!	!	+	+	•	!
						ow risk	L Some o	oncarns	Uiah a	ick

Fig 3. Risk of bias of crossover included studies. D1 = randomization process; DS = bias arising from period and carryover effects; D2 = deviations from the intended interventions; D3 = missing outcome data; D4 = measurement of the outcome; D5 = selection of the reported result. *Abbreviations*: ASV = adaptive support ventilation; ICU = intensive care unit; FiO2 = fraction of inspired oxygen; HFNC = high-flow nasal cannula; PSV = pressure support ventilation; SpO2 = peripheral oxygen saturation.

https://doi.org/10.1371/journal.pone.0304745.g003

clinically important difference in effect estimates between the meta-analysis encompassing all studies and the one focusing solely on studies with a low risk of bias, any inconsistency is likely to be insignificant [55]. Despite the observed asymmetry in the funnel plot (**S5 Fig in S1 File**), we cannot definitively attribute it to publication bias due to the limited number of studies included in this outcome. In 6 studies [36,37,40,46,48,49], closed-loop oxygen titration reduced the percentage of time with hypoxemia, with substantial unexplained heterogeneity, due to imprecise results of 2 studies [36,49] (**S6 Fig in S1 File**). In 4 studies [36,46,48,49] closed-loop oxygen titration did not reduce the percentage of time with hyperoxemia, with substantial heterogeneity and inconsistency in 2 studies [36,46] (**S7 Fig in S1 File**).

Four studies [38,43,45,53] were not included in the quantitative analyses due to the absence of reported data from the first crossover period. The mean percentage of time spent in the SpO2 target, ranging from 92 to 96%, was higher in the intervention group when compared to the control group, 83 (SD 21) versus 33 (SD 36), respectively [45]. These findings were similar to one study [43] published, in which the estimated effect presented, a mean difference (MD) of 38.5% [95% CI, 27.8 to 49.3]. One study [38] investigated 265 mechanically ventilated patients during daily nursing procedures. The closed-loop group spent a mean percentage of time in the SpO2 target of 48 (SD 37), compared with the control group of 43 (SD 37). One study [53] investigated the duration during which mechanically ventilated patients maintained their SpO2 within the optimal breath zone, considering the SpO2 target and FiO2. They found

	Closed-Loop			Manual				Mean Difference	Mean Difference	
Study or Subgroup	Mean	SD	Total	Mean	SD	Total	Weight	IV, Random, 95% CI	IV, Random, 95% CI	
1.19.1 Invasive mechanic	cal venti	lation								
Bialais 2016 (i-ASV)	84.15	19.0936	42	69.9	24.3389	38	13.4%	14.25 [4.59, 23.91]		
De Bie 2020 (i-ASV)	55.2	28	109	25.5	29.3	111	15.9%	29.70 [22.13, 37.27]		
Subtotal (95% CI)			151			149	29.3%	22.28 [7.15, 37.41]		
Heterogeneity: Tau ² = 99.7	75; Chi² =	= 6.09, df =	= 1 (P =	= 0.01);	² = 84%					
Test for overall effect: Z =	2.89 (P =	= 0.004)								
1.19.2 Low flow oxygen										
L'Her 2017 (FreeO2)	81.3	20.7	93	51.8	30	94	16.1%	29.50 [22.12, 36.88]		
L'Her 2021 (FreeO2)	94	11.3	83	62.1	23.3	86	18.5%	31.90 [26.41, 37.39]		
Lelouche 2016 (FreeO2)	81.2	15.9	25	51.3	19.7	25	13.1%	29.90 [19.98, 39.82]		
Subtotal (95% CI)			201			205	47.7%	30.86 [26.83, 34.88]	•	
Heterogeneity: Tau ² = 0.00); Chi² =	0.30, df =	2 (P =	0.86); l²	² = 0%					
Test for overall effect: Z =	15.02 (P	< 0.0000	1)							
4 40 0 11 1 1										
1.19.3 High flow oxygen										
Harper 2021 (HFNC)	96.4	2.236	10	72.9	24.854	10	8.2%	23.50 [8.03, 38.97]		
Roca 2022 (HFNC)	95.7	4.34	22	79.8	19.95	23	14.9%	15.90 [7.55, 24.25]		
Subtotal (95% CI)			32			33	23.1%	17.62 [10.27, 24.97]		
Heterogeneity: Tau ² = 0.00; Chi ² = 0.72, df = 1 (P = 0.40); l ² = 0%										
Test for overall effect: Z =	4.70 (P <	< 0.00001)							
Total (95% CI)			384			387	100.0%	25.47 [19.93. 31.00]	•	
Heterogeneity: $Tau^2 = 35$	33. Chi2 -	- 17 07 dt	F = 6 (P	- 0.006	(12 - 67)	4				
Test for overall effect: 7 =	9 02 (P <		- 0 (i)	- 0.000	,, - 017	U			-50 -25 0 25 50	
Test for subgroup difference	ces: Chi ²	= 10.07.0	, df = 2 (I	P = 0.00	$(6), l^2 = 80$	1%			Favours [Manual] Favours [Closed-Loop]	

Fig 4. Subgroup analyses of time spent in the SpO₂ target based on pooled data from five studies of noninvasive devices and two studies of mechanical ventilation. *Abbreviations*: CI = confidence interval; HFNC = high flow nasal cannula; i-ASV = INTELLiVENT-Adaptive Support Ventilation; SD = standard deviation.

https://doi.org/10.1371/journal.pone.0304745.g004

that patients in the closed-loop group spent 61% of their time within the optimal oxygenation range, a statistically significant improvement compared to the 52% observed with conventional ventilation. One study [52] evaluated the time spent in the optimal zone (i.e., SpO2 target) in 60 participants. The mean time spent in the closed-loop and control group remained at the optimal zone of 192 (SD 52), and 25 (SD 124) minutes, respectively. This study [52] was not included in the quantitative analysis due to the absence of information on all patients included to change mean and SD data from minutes to percentage of time.

Clinical outcomes

Of all included studies, 10 studies [36,37,40,42,44,46-49,54] evaluated length of stay, and 10 studies [36,37,39,40,42,45,47,49,50,54] evaluated mortality. In most studies, the intervention was not applied throughout the total duration of oxygen supplementation (S11 Table in S1 File). We considered it inappropriate to meta-analyze these details.

Costs

Of all studies, only one study [48] investigated hospitalization costs, and 2 studies [45,49] evaluated oxygen consumption. These studies suggest that closed-loop oxygen titration systems may result in a reduction in costs and oxygen consumption (S11 Table in S1 File).

Adverse events

Of all studies, 14 studies [34,35,37–39,42–48,51,52] reported adverse events. Quantitative analyses were not performed due to the wide variety and the absence of a standardized definition of adverse events (S11 and S12 Tables in S1 File).

	Closed-loop			Manual Me				Mean Difference	Mean Difference
Study or Subgroup	Mean	SD	Total	Mean	SD	Total	Weight	IV, Random, 95% CI	IV, Random, 95% CI
1.9.1 Invasive mechan	ical vent								
Arnal 2018 (i-ASV)	2	2	30	6	5	30	29.5%	-4.00 [-5.93, -2.07]	
Bialais 2016 (i-ASV)	3	4	42	13	11	38	19.2%	-10.00 [-13.70, -6.30]	
Komnov 2023 (i-ASV) Subtotal (95% Cl)	1	1	16 88	6	2	16 84	34.2% 82.8%	-5.00 [-6.10, -3.90] -5.72 [-8.06, -3.39]	*
Heterogeneity: Tau ² = 3	.02; Chi ²	= 8.02	2, df = :	2 (P = 0	.02);	$ ^2 = 75^{\circ}$	%		
Test for overall effect: Z	= 4.80 (F	P < 0.0	00001)						
1.9.2 High flow oxygen Harper 2021 (HFNC) Subtotal (95% CI) Heterogeneity: Not appl Test for overall effect: Z	1 4 icable = 0.38 (F	3.9 P = 0.7	10 10 70)	4.8	5.4	10 10	17.2% 17.2%	-0.80 [-4.93, 3.33] - 0.80 [-4.93, 3.33]	-
Total (95% CI)			98			94	100.0%	-4.94 [-7.28, -2.61]	\bullet
Heterogeneity: $Tau^2 = 3.85$; $Chi^2 = 11.95$, $df = 3$ (P = 0.008); $l^2 = 75\%$ -1 -20 -10 0 Test for overall effect: Z = 4.15 (P < 0.0001)									

Fig 5. Forest plot for workload of healthcare professionals. *Abbreviations*: CI = confidence interval; HFNC = high flow nasal cannula; i-ASV = INTELLIVENT-Adaptive Support Ventilation; SD = standard deviation.

https://doi.org/10.1371/journal.pone.0304745.g005

Healthcare workload

Of the 7 studies [35-37,42,44,49,54] that evaluated the workload of healthcare professionals, 4 studies [36,37,44,54] were included in a quantitative analysis, suggesting that closed-loop oxygen control may reduce healthcare workload (Fig 5). In the sensitivity analyses excluding trials with some concerns according to the risk of bias assessment (S8 Fig in <u>S1 File</u>), we observed a slight increase in the heterogeneity (I² = 83%) maintaining the direction of the effect of interventions (MD -4.94, 95% CI -9.43 to -0.46).

Discussion

The results of this systematic review and meta-analysis show that closed-loop oxygen control devices may increase the percentage of time in the SpO₂ target and reduce the percentage of time with hypoxemia. Closed-loop oxygen titration was found not to be associated with a lower percentage of time with hyperoxemia, time for weaning supplemental oxygen, adverse events, length of stay, and mortality. The evidence suggests that closed-loop oxygen titration may reduce costs, oxygen consumption and healthcare workload.

Although both are harmful, episodes of hypoxemia (i.e., $SpO_2 < 90$ to 92%) are frequently regarded as a *'red flag'* when compared with hyperoxemia episodes (i.e., $SpO_2 > 94-96$ to 100%). Consequently, supplemental oxygen is typically considered to be widely available and liberally administered, with little to no attention given to avoiding excessive use. Closed-loop systems for oxygen titration seem to increase the time spent in the SpO_2 target, contributing to the correct dose adjustment of oxygen and avoiding hypoxia episodes. However, it may not effectively reduce the occurrence of hyperoxemia, particularly in mechanical ventilation systems, or when high oxygen targets are being used. For instance, one trial [36] showed a high incidence of hyperoxemia with automatic oxygen control; however, this incidence was also high in the control group (35% vs. 33%). The included studies reported the use of different upper limits of SpO_2 target and timing of intervention, making it difficult to compare the results. In this scenario, the predefined width range could have influenced the results of the percentage of time in the SpO_2 target, not fully explained by

subgroup analysis. The duration of the intervention and the different patient characteristics, could explain the substantial heterogeneity observed among the studies included in the subgroup analysis of patients undergoing mechanical ventilation support.

Supplemental oxygen use is frequently observed in patients admitted to hospitals and is considered one of the most common treatments provided to critically ill patients. The underrecognition and delay in the correction of hypoxemia remains a major problem in clinical practice associated with an increase in mortality [56], especially in settings where resources are limited and in countries or regions [57]. In addition to socioeconomic factors, the health crisis triggered by the coronavirus disease (COVID-19) pandemic highlighted inequities related to an unprecedented need for oxygen supply and other resources, such as ICU beds, specialized professionals, invasive and noninvasive ventilators, and other medical supplies [58], and consequently, the indices of mortality in patients with and without COVID-19 significantly increased [59,60].

Investments in technology to develop new innovative devices for closed-loop oxygen titration that can efficiently manage resources have grown exponentially in recent years for both clinical and research uses. This was also observed in our findings, as most of the included trials were industry-sponsored investigations. Furthermore, unusual situations, such as the COVID-19 pandemic highlighted the importance and need to rationally use oxygen in clinical practice, as prolonged use of unnecessary supplemental oxygen increases health costs, ICU, and hospital stays, exposing patients to infections and psychosocial complications resulting from the hospitalization process. Despite the benefits of closed-loop devices for oxygen titration, the clinical use at bedside is still limited, due to the high initial investment necessary to acquire such technology. The direct cost related to closed-loop oxygen titration systems was evaluated in only one study [48], and oxygen consumption, an indirect measure of cost, was evaluated in two studies [45,49], showing inaccurate results due to the small number of included patients. Future studies should investigate 'cost outcome' to assist decision-makers in providing an assertive plan for closed-loop oxygen titration for clinical and research uses.

Patients admitted to the ICU are frequently exposed to different procedures, such as physical therapy and rehabilitation interventions (e.g., airway clearance, suctioning, and positioning), which can promote respiratory and hemodynamic changes according to their clinical condition. In addition, lung dynamics are directly affected by the patient's clinical presentation. Thus, it can require frequent parameter adjustments during invasive or noninvasive ventilation support to maintain the predefined target of SpO₂, increasing the workload of healthcare professionals. Closed-loop devices cannot replace healthcare professionals as they cannot handle complex clinical changes in oxygen adjustments. However, they can ease the workload of professionals by automating oxygen adjustments and freeing up time for patient care. These devices do not eliminate the need for clinical discussions and adjustments, especially for significant changes. They require alarms and data visualization for informed decision-making. While their safety is not clearly defined, they offer a potential way to reduce bedside tasks; however, more research on adverse events is needed due to oxygen-related risks. The identification and reporting processes of adverse events can be challenging. Researchers may consider the inclusion of adverse events as an outcome, reporting additional details of the definition, occurrence, and classification of events.

Strengths and limitations

This review has strengths that have helped reduce potential bias in the review process, such as a thorough literature search, rigorous and well-established methods to minimize bias, including multiple reviewers to independently screen abstracts, review studies, extract data, assess the risk of bias and the certainty of evidence outlined by the Cochrane Collaboration [30]. In addition, primary study authors were contacted to provide additional information, data, or clarifications when needed. For crossover randomized controlled trials (RCTs), we also contacted authors to request data from the first phase of interventions to describe baseline characteristics and results separately from patients randomized to each stage. We, therefore, believe that it is unlikely that we missed any relevant trials. Thus, this systematic review and meta-analysis provided important insights. Additional studies investigating the use of closed-loop oxygen systems through the entire duration of oxygen support should be developed to improve the findings for certainty of evidence.

However, this review has some limitations. First, due to the short intervention period in most studies it was difficult to establish a direct relationship between the device and the patient's length of stay, mortality and time for oxygen weaning. It is important to emphasize that most studies in this meta-analysis featured brief intervention periods, ranging from only a few hours to 24 hours. This limitation diminishes the generalizability of the results to patients necessitating prolonged oxygen therapy, typically those in the most critical condition. Second, despite our considerable efforts to mitigate publication bias, we were unable to reliably assess it through the funnel plot due to the limited number of available studies. Third, there was insufficient information about the randomization process of some studies, so we recommended the use of reliable methods (e.g., computer-generated random numbers) and detailed reporting to ensure homogeneity between groups. Fourth, the washout period between interventions of crossover RCTs was short and heterogeneous, and the possibility of the carry-over effect cannot be ruled out. Although the devices provide the same substance (oxygen), the type of device can influence the patient's stability. Fifth, the findings of the invasive mechanical ventilation subgroup should be interpreted with caution. Automatic mechanical ventilation modes adjust a group of ventilation settings such as tidal volume, respiratory rate, FiO_2 , and positive end-expiratory pressure. Thus, the findings of this review regarding the percentage of time spent in the SpO₂ target cannot be solely attributed to the oxygen titration. Additionally, the number of adjustments in ventilation parameters may also have influenced healthcare workload outcomes.

In summary, while closed-loop oxygen titration systems probably increase the percentage of time spent within desired SpO_2 ranges and reduce healthcare workload and costs, their safety remains uncertain. Adverse events are likely underreported in existing studies. Moreover, the extent to which these systems enhance efficiency, such as by reducing the duration of oxygen weaning, and even more importantly, by affecting length of stay and mortality rates, remains unclear. This uncertainty primarily stems from the limited duration of the intervention in the scarce, and small existing studies. Future research should focus on evaluating effectiveness, safety, and efficacy over longer periods, and prioritize patient-centered endpoints.

Supporting information

S1 File. Additional analysis, detailed search strategies, and supplementary information on the included studies.

(DOCX)

Acknowledgments

The authors acknowledge Dr Thiago Domingos Corrêa for providing administrative support. We would like to thank Helena Spalic for the valuable contribution proofreading this systematic review.

Author Contributions

- **Conceptualization:** Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Ricardo Kenji Nawa.
- **Data curation:** Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Bianca Maria Schneider Pereira Garcia, Emanuel dos Santos Pereira, Raquel Afonso Caserta Eid, Ricardo Kenji Nawa.
- Formal analysis: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.
- **Investigation:** Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.
- Methodology: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Bianca Maria Schneider Pereira Garcia, Emanuel dos Santos Pereira, Raquel Afonso Caserta Eid, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.
- Project administration: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Ricardo Kenji Nawa.
- Supervision: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Ricardo Kenji Nawa.
- Validation: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Bianca Maria Schneider Pereira Garcia, Emanuel dos Santos Pereira, Raquel Afonso Caserta Eid, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.
- Visualization: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Ricardo Kenji Nawa.
- Writing original draft: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Bianca Maria Schneider Pereira Garcia, Emanuel dos Santos Pereira, Raquel Afonso Caserta Eid, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.
- Writing review & editing: Caroline Gomes Mól, Aléxia Gabriela da Silva Vieira, Marcus J. Schultz, Ana Carolina Pereira Nunes Pinto, Ricardo Kenji Nawa.

References

- O'Driscoll BR, Howard LS, Earis J, Mak V, British Thoracic Society Emergency Oxygen Guideline Group, BTS Emergency Oxygen Guideline Development Group. BTS guideline for oxygen use in adults in healthcare and emergency settings. Thorax [Internet]. 2017; 72:ii1–90. Available from: http://dx.doi. org/10.1136/thoraxjnl-2016-209729.
- Petersson J, Glenny RW. Gas exchange and ventilation-perfusion relationships in the lung. Eur Respir J [Internet]. 2014 [cited 2022 Aug 24]; 44:1023–41. Available from: https://erj.ersjournals.com/content/ erj/44/4/1023. https://doi.org/10.1183/09031936.00037014 PMID: 25063240
- Gilbert-Kawai ET, Mitchell K, Martin D, Carlisle J, Grocott MPW. Permissive hypoxaemia versus normoxaemia for mechanically ventilated critically ill patients. Cochrane Database Syst Rev [Internet]. 2014;CD009931. Available from: https://doi.org/10.1002/14651858.CD009931.pub2 PMID: 24801519
- de Jonge E, Peelen L, Keijzers PJ, Joore H, de Lange D, van der Voort PHJ, et al. Association between administered oxygen, arterial partial oxygen pressure and mortality in mechanically ventilated intensive care unit patients. Crit Care [Internet]. 2008; 12:R156. Available from: https://doi.org/10.1186/cc7150 PMID: 19077208
- Barbateskovic M, Schjørring OL, Russo Krauss S, Jakobsen JC, Meyhoff CS, Dahl RM, et al. Higher versus lower fraction of inspired oxygen or targets of arterial oxygenation for adults admitted to the intensive care unit. Cochrane Database Syst Rev [Internet]. 2019;2019. Available from: https://doi.org/ 10.1002/14651858.CD012631.pub2 PMID: 31773728
- O'Driscoll BR, Smith R. Oxygen Use in Critical Illness. Respir Care [Internet]. 2019; 64:1293–307. Available from: https://doi.org/10.4187/respcare.07044 PMID: 31409632

- Barbateskovic M, Schjørring OL, Krauss SR, Meyhoff CS, Jakobsen JC, Rasmussen BS, et al. Higher vs Lower Oxygenation Strategies in Acutely III Adults. Chest [Internet]. 2021; 159:154–73. Available from: http://dx.doi.org/10.1016/j.chest.2020.07.015.
- Valencia Gallardo JM, Solé Violán J, Rodríguez de Castro F. [Translated article] Oxygen therapy. Considerations regarding its use in acute ill patients. Arch Bronconeumol [Internet]. 2022; 58:T102–3. Available from: http://dx.doi.org/10.1016/j.arbres.2021.03.024.
- Chu DK, Kim LH-Y, Young PJ, Zamiri N, Almenawer SA, Jaeschke R, et al. Mortality and morbidity in acutely ill adults treated with liberal versus conservative oxygen therapy (IOTA): a systematic review and meta-analysis. Lancet [Internet]. 2018; 391:1693–705. Available from: https://doi.org/10.1016/S0140-6736(18)30479-3 PMID: 29726345
- Li L, Zhang Y, Wang P, Chong W, Hai Y, Xu P, et al. Conservative versus liberal oxygen therapy for acutely ill medical patients: A systematic review and meta-analysis. Int J Nurs Stud [Internet]. 2021; 118:103924. Available from: http://dx.doi.org/10.1016/j.ijnurstu.2021.103924.
- Dahl RM, Grønlykke L, Haase N, Holst LB, Perner A, Wetterslev J, et al. Variability in targeted arterial oxygenation levels in patients with severe sepsis or septic shock. Acta Anaesthesiol Scand [Internet]. 2015; 59:859–69. Available from: https://doi.org/10.1111/aas.12528 PMID: 25914095
- Helmerhorst HJF, Arts DL, Schultz MJ, van der Voort PHJ, Abu-Hanna A, de Jonge E, et al. Metrics of Arterial Hyperoxia and Associated Outcomes in Critical Care. Crit Care Med [Internet]. 2017; 45:187– 95. Available from: https://doi.org/10.1097/CCM.0000000002084 PMID: 27763912
- Kilgannon JH, Jones AE, Shapiro NI, Angelos MG, Milcarek B, Hunter K, et al. Association between arterial hyperoxia following resuscitation from cardiac arrest and in-hospital mortality. JAMA [Internet]. 2010; 303:2165–71. Available from: https://doi.org/10.1001/jama.2010.707 PMID: 20516417
- Lilien TA, van Meenen DMP, Schultz MJ, Bos LDJ, Bem RA. Hyperoxia-induced lung injury in acute respiratory distress syndrome: what is its relative impact? Am J Physiol Lung Cell Mol Physiol [Internet]. 2023; 325:L9–16. Available from: https://doi.org/10.1152/ajplung.00443.2022 PMID: 37129255
- Benoît Z, Wicky S, Fischer J-F, Frascarolo P, Chapuis C, Spahn DR, et al. The Effect of Increased Fio2 Before Tracheal Extubation on Postoperative Atelectasis. Anesthesia & Analgesia [Internet]. 2002 [cited 2022 Aug 24]; 95:1777. Available from: https://journals.lww.com/anesthesia-analgesia/FullText/ 2002/12000/The_Effect_of_Increased_Fio2_Before_Tracheal.58.aspx. https://doi.org/10.1097/ 00000539-200212000-00058 PMID: 12456458
- Rothen HU, Sporre B, Engberg G, Wegenius G, Högman M, Hedenstierna G. Influence of gas composition on recurrence of atelectasis after a reexpansion maneuver during general anesthesia. Anesthesiology [Internet]. 1995; 82:832–42. Available from: https://doi.org/10.1097/00000542-199504000-00004 PMID: 7717553
- Rothen HU, Sporre B, Engberg G, Wegenius G, Reber A, Hedenstierna G. Prevention of atelectasis during general anaesthesia. Lancet [Internet]. 1995; 345:1387–91. Available from: <u>https://doi.org/10.1016/s0140-6736(95)92595-3</u> PMID: 7760608
- Munshi L, Ferguson ND. Evolving Issues in Oxygen Therapy in Acute Care Medicine. JAMA [Internet]. 2020; 323:607–8. Available from: https://doi.org/10.1001/jama.2019.22029 PMID: 31977030
- Martin DS, Grocott MPW. Oxygen therapy in critical illness: precise control of arterial oxygenation and permissive hypoxemia. Crit Care Med [Internet]. 2013; 41:423–32. Available from: http://dx.doi.org/10. 1097/CCM.0b013e31826a44f6.
- 20. Lumb AB, Walton LJ. Perioperative oxygen toxicity. Anesthesiol Clin [Internet]. 2012; 30:591–605. Available from: https://doi.org/10.1016/j.anclin.2012.07.009 PMID: 23089497
- Semler MW, Casey JD, Lloyd BD, Hastings PG, Hays MA, Stollings JL, et al. Oxygen-Saturation Targets for Critically III Adults Receiving Mechanical Ventilation. N Engl J Med [Internet]. 2022; 387:1759– 69. Available from: https://doi.org/10.1056/NEJMoa2208415 PMID: 36278971
- Rush B, Celi LA, Stone DJ. Applying machine learning to continuously monitored physiological data. J Clin Monit Comput [Internet]. 2019; 33:887–93. Available from: http://dx.doi.org/10.1007/s10877-018-0219-z.
- Denault M-H, Péloquin F, Lajoie A-C, Lacasse Y. Automatic versus Manual Oxygen Titration in Patients Requiring Supplemental Oxygen in the Hospital: A Systematic Review and Meta-Analysis. Respiration [Internet]. 2019; 98:178–88. Available from: https://doi.org/10.1159/000499119 PMID: 31129662
- Sanchez-Morillo D, Olaby O, Fernandez-Granero MA, Leon-Jimenez A. Physiological closed-loop control in intelligent oxygen therapy: A review. Comput Methods Programs Biomed [Internet]. 2017; 146:101–8. Available from: https://doi.org/10.1016/j.cmpb.2017.05.013 PMID: 28688479
- Liberati A, Altman DG, Tetzlaff J, Mulrow C, Gøtzsche PC, Ioannidis JPA, et al. The PRISMA statement for reporting systematic reviews and meta-analyses of studies that evaluate health care interventions: explanation and elaboration. J Clin Epidemiol [Internet]. 2009; 62:e1–34. Available from: https://doi.org/ 10.1016/j.jclinepi.2009.06.006 PMID: 19631507

- 26. Cumpston M, Li T, Page MJ, Chandler J, Welch VA, Higgins JPT, et al. Updated guidance for trusted systematic reviews: a new edition of the Cochrane Handbook for Systematic Reviews of Interventions. Cochrane Database Syst Rev [Internet]. 2019; 10:14651858. Available from: https://research-information.bris.ac.uk/files/223404871/Cumpston_et_al_2019_Cochrane_Database_of_Systematic_Reviews.pdf. https://doi.org/10.1002/14651858.ED000142 PMID: 31643080
- Mol CG, Vieira AG da S, Garcia BMSP, Pereira EDS, Eid RAC, Pinto P, et al. Closed-loop oxygen control for patients with hypoxaemia during hospitalisation: a living systematic review and meta-analysis protocol. BMJ Open [Internet]. 2022; 12:e062299. Available from: https://doi.org/10.1136/bmjopen-2022-062299 PMID: 36523244
- 28. Website [Internet]. Available from: Thomas J, Kneale D, Mckenzie JE, et al. Chapter 2: Determining the scope of the review and the questions it will address. In 2021;6.Available from <u>www.training.cochrane.org/handbook</u>.
- Higgins JPT, Altman DG, Gotzsche PC, Juni P, Moher D, Oxman AD, et al. The Cochrane Collaboration's tool for assessing risk of bias in randomised trials [Internet]. BMJ. 2011. p. d5928–d5928. Available from: http://dx.doi.org/10.1136/bmj.d5928.
- Higgins JPT, Thomas J. Cochrane Handbook for Systematic Reviews of Interventions [Internet]. John Wiley & Sons; 2019. Available from: https://books.google.com/books/about/Cochrane_Handbook_for_ Systematic_Reviews.html?hl=&id=G1SyDwAAQBAJ.
- Wan X, Wang W, Liu J, Tong T. Estimating the sample mean and standard deviation from the sample size, median, range and/or interquartile range [Internet]. BMC Medical Research Methodology. 2014. Available from: http://dx.doi.org/10.1186/1471-2288-14-135.
- GRADE handbook [Internet]. [cited 2022 Jan 31]. Available from: http://guidelinedevelopment.org/handbook.
- 33. GRADEpro [Internet]. [cited 2022 Jan 31]. Available from: http://gradepro.org.
- 34. Arnal J-M, Wysocki M, Demory D, Durish G, Laubscher T, Novotni D, et al. Prospective Randomized Cross-over Controlled Study Comparing Adaptive Support Ventilation (ASV) And A Fully Close Loop Control Solution (Intellivent®) In Adult ICU Patients With Acute Respiratory Failure [Internet]. B42. INVASIVE MECHANICAL VENTILATION: SOUP TO NUTS. 2010. Available from: https://doi.org/10. 1164/ajrccm-conference.2010.181.1%5Fmeetingabstracts.a3004
- **35.** Arnal J-M, Wysocki M, Novotni D, Demory D, Lopez R, Donati S, et al. Safety and efficacy of a fully closed-loop control ventilation (IntelliVent-ASV®) in sedated ICU patients with acute respiratory failure: a prospective randomized crossover study. Intensive Care Med [Internet]. 2012; 38:781–7. Available from: http://dx.doi.org/10.1007/s00134-012-2548-6.
- Arnal J-M, Garnero A, Novotni D, Corno G, Donati S-Y, Demory D, et al. Closed loop ventilation mode in Intensive Care Unit: a randomized controlled clinical trial comparing the numbers of manual ventilator setting changes. Minerva Anestesiol [Internet]. 2018; 84:58–67. Available from: <u>https://doi.org/10.</u> 23736/S0375-9393.17.11963-2 PMID: 28679200
- Bialais E, Wittebole X, Vignaux L, Roeseler J, Wysocki M, Meyer J, et al. Closed-loop ventilation mode (IntelliVent®-ASV) in intensive care unit: a randomized trial. Minerva Anestesiol [Internet]. 2016; 82:657–68. Available from: https://www.ncbi.nlm.nih.gov/pubmed/26957117.
- Chelly J, Mazerand S, Jochmans S, Weyer C-M, Pourcine F, Ellrodt O, et al. Automated vs. conventional ventilation in the ICU: a randomized controlled crossover trial comparing blood oxygen saturation during daily nursing procedures (I-NURSING). Crit Care [Internet]. 2020; 24:453. Available from: http://dx.doi.org/10.1186/s13054-020-03155-3.
- Clavieras N, Wysocki M, Coisel Y, Galia F, Conseil M, Chanques G, et al. Prospective Randomized Crossover Study of a New Closed-loop Control System versus Pressure Support during Weaning from Mechanical Ventilation. Anesthesiology [Internet]. 2013 [cited 2022 Aug 16]; 119:631–41. Available from: https://pubs.asahq.org/anesthesiology/article-pdf/119/3/631/262967/20130900_0-00023.pdf. https://doi.org/10.1097/ALN.0b013e3182952608 PMID: 23619172
- 40. De Bie AJR, Neto AS, van Meenen DM, Bouwman AR, Roos AN, Lameijer JR, et al. Fully automated postoperative ventilation in cardiac surgery patients: a randomised clinical trial. Br J Anaesth [Internet]. 2020; 125:739–49. Available from: https://doi.org/10.1016/j.bja.2020.06.037 PMID: 32739044
- Denault M-H, Ruel C, Simon M, Bouchard P-A, Simard S, Lellouche F. Evaluation of hyperoxia-induced hypercapnia in obese patients after cardiac surgery: a randomized crossover comparison of conservative and liberal oxygen administration. Canadian Journal of Anesthesia/Journal canadien d'anesthésie [Internet]. 2019 [cited 2022 Aug 16]; 67:194–202. Available from: https://link.springer.com/article/10. 1007/s12630-019-01500-x.
- 42. Eremenko A.A., Komnov R.D., Titov P.A., Gerasimenko S.A., Chakal D.A. Comparing the Intellivent-ASV® Mode with Conventional Ventilation Modes during Weaning after Uncomplicated Cardiac

Surgery. Messenger of Anesthesiology and Resuscitation. 2021; 18(3):36–45. (In Russ.) https://doi.org/10.21292/2078-5658-2021-18-3-36-45.

- Hansen EF, Hove JD, Bech CS, Jensen J-US, Kallemose T, Vestbo J. Automated oxygen control with O2matic during admission with exacerbation of COPD. Int J Chron Obstruct Pulmon Dis [Internet]. 2018; 13:3997–4003. Available from: http://dx.doi.org/10.2147/COPD.S183762.
- 44. Harper J, Kearns N, Bird G, Braithwaite I, Eathorne A, Shortt N, et al. Automatic versus manual oxygen titration using a novel nasal high-flow device in medical inpatients with an acute illness: a randomised controlled trial. BMJ Open Respir Res [Internet]. 2021;8. Available from: <u>https://doi.org/10.1136/bmjresp-2020-000843 PMID: 34362762</u>
- Johannigman JA, Branson R, Lecroy D, Beck G. Autonomous Control of Inspired Oxygen Concentration During Mechanical Ventilation of the Critically Injured Trauma Patient. J Trauma Acute Care Surg [Internet]. 2009 [cited 2022 Aug 16]; 66:386. Available from: https://journals.lww.com/jtrauma/Abstract/ 2009/02000/Autonomous_Control_of_Inspired_Oxygen.13.aspx. https://doi.org/10.1097/TA. 0b013e318197a4bb PMID: 19204511
- 46. L'Her E, Dias P, Gouillou M, Riou A, Souquiere L, Paleiron N, et al. Automatic versus manual oxygen administration in the emergency department. Eur Respir J [Internet]. 2017 [cited 2022 Aug 16]; 50. Available from: https://erj.ersjournals.com/content/50/1/1602552.abstract. https://doi.org/10.1183/ 13993003.02552-2016 PMID: 28729473
- 47. L'Her E, Jaber S, Verzilli D, Jacob C, Huiban B, Futier E, et al. Automated closed-loop versus standard manual oxygen administration after major abdominal or thoracic surgery: an international multicentre randomised controlled study. Eur Respir J [Internet]. 2021 [cited 2022 Aug 16]; 57. Available from: https://erj.ersjournals.com/content/57/1/2000182.abstract. https://doi.org/10.1183/13993003.00182-2020 PMID: 32855218
- Lellouche F, Bouchard P-A, Roberge M, Simard S, L'Her E, Maltais F, et al. Automated oxygen titration and weaning with FreeO₂ in patients with acute exacerbation of COPD: a pilot randomized trial. COPD [Internet]. 2016 [cited 2022 Aug 16]; 11:1983–90. Available from: https://www.dovepress.com/getfile. php?fileID=32049.
- 49. Roca O, Caritg O, Santafé M, Ramos FJ, Pacheco A, García-de-Acilu M, et al. Closed-loop oxygen control improves oxygen therapy in acute hypoxemic respiratory failure patients under high flow nasal oxygen: a randomized cross-over study (the HILOOP study). Crit Care [Internet]. 2022; 26:108. Available from: https://doi.org/10.1186/s13054-022-03970-w PMID: 35422002
- 50. Huynh Ky MT, Bouchard P-A, Morin J, L'Her E, Sarrazin J-F, Lellouche F. Closed-Loop Adjustment of Oxygen Flowrate with FreeO2 in Patients with Acute Coronary Syndrome: Comparison of Automated Titration with FreeO2 (Set at Two SpO2 Target) and of Manual Titration. A Randomized Controlled Study. Am J Respir Crit Care Med. 2017; 195:A3766.
- Kobayashi T, Onodera Y, Suzuki H, Nakane M, Kawamae K. Manual ASV vs INTELLiVENT-ASV for the patients after cardiac surgery—are automated ventilators better for the patients? Intensive Care Medicine Experimental. ESICM LIVES 2017. 2017; 5(Suppl 2):0749. <u>https://doi.org/10.1186/s40635-017-0151-4</u>.
- Lellouche F, Bouchard P-A, Simard S, L'Her E, Wysocki M. Evaluation of fully automated ventilation: a randomized controlled study in post-cardiac surgery patients. Intensive Care Med [Internet]. 2013; 39:463–71. Available from: https://doi.org/10.1007/s00134-012-2799-2 PMID: 23338569
- Buiteman-Kruizinga LA, Tsonas AM, Botta M, Guseva AA, Serpa NA, Schultz M, et al. The Effect of INTELLiVENT–ASV on Oxygenation, FiO2 and PEEP in Critically III Invasively Ventilated Patients. Intensive Care Medicine Experimental. ESICM LIVES 2022: part 1. 2022; 10(2):000489. https://doi.org/ 10.1186/s40635-022-00468-1.
- 54. 42nd International Symposium on Intensive Care & Emergency Medicine. Crit Care [Internet]. 2023 [cited 2023 Aug 21];27:1–156. Available from: https://ccforum.biomedcentral.com/articles/10.1186/ s13054-023-04377-x.
- 55. Guyatt G, Zhao Y, Mayer M, Briel M, Mustafa R, Izcovich A, et al. GRADE guidance 36: updates to GRADE's approach to addressing inconsistency. J Clin Epidemiol. 2023; 158: 70–83. <u>https://doi.org/10.1016/j.jclinepi.2023.03.003</u> PMID: 36898507
- 56. Sutherland T, Musafiri S, Twagirumugabe T, Talmor D, Riviello ED. Oxygen as an Essential Medicine: Under- and Over-Treatment of Hypoxemia in Low- and High-Income Nations. Critical care medicine. 2016. pp. e1015–6. https://doi.org/10.1097/CCM.00000000001912 PMID: 27635514
- Meara JG, Leather AJM, Hagander L, Alkire BC, Alonso N, Ameh EA, et al. Global Surgery 2030: evidence and solutions for achieving health, welfare, and economic development. Int J Obstet Anesth. 2016; 25: 75–78. https://doi.org/10.1016/j.ijoa.2015.09.006 PMID: 26597405

- Kaye AD, Okeagu CN, Pham AD, Silva RA, Hurley JJ, Arron BL, et al. Economic impact of COVID-19 pandemic on healthcare facilities and systems: International perspectives. Best Pract Res Clin Anaesthesiol. 2021; 35: 293–306. https://doi.org/10.1016/j.bpa.2020.11.009 PMID: 34511220
- 59. Usher AD. Medical oxygen crisis: a belated COVID-19 response. Lancet. 2021; 397: 868–869. https:// doi.org/10.1016/S0140-6736(21)00561-4 PMID: 33676621
- 60. Website. Available: Every Breath Counts Coalition. One year after the launch of the access to COVID19 tools accelerator (ACT-A) oxygen emergency taskforce, what has been achieved?. Available at: https://stoppneumonia.org/wp-content/uploads/2022/02/ACTAOxygenTaskforceAnniversaryStatement23February2022.pdf. Accessed March 2, 2022.