

Gas-Diffusion Sequential Injection Enzymatic Determination of Ethanol in Wines: Study of Matrix Effects Derived From Direct Sample Introduction



Teresa F. M. Pais, Susana S. M. P. Vidigal, Ildikó V. Tóth and António O. S. S. Rangel

Escola Superior de Biotecnologia, Universidade Católica Portuguesa

Rua Dr. António Bernardino de Almeida, 4200-072 Porto, Portugal

e-mail: aorangel@esb.ucp.pt

DETERMINATION OF ETHANOL

Determination of alcoholic degree

Precision of reference methods [1]:
on wine distillate : 0.04 -0.6% (v/v)
on non-distilled wine: 0.2% (v/v)

Fermentation monitoring

Precision of methods [2]:
on non-distilled wine: ≈ 2 % (v/v)

Enzymatic determination of ethanol

Alcohol dehydrogenase (ADH)



- Higher stability
- Better selectivity
- Need of cofactor

Detection method
•spectrophotometric
•fluorimetric
•amperometric

Alcohol oxidase (AOD)



- Lower stability
- Cofactor bounded to the enzyme

Detection method
•chemiluminometric
•amperometric
•spectrophotometric

DEVELOPED SYSTEM




Figure 1. Sequential injection manifold for the enzymatic determination of ethanol in beverages, P, peristaltic pump; C, carrier, H₂O; B, buffer, pH 10; S, sample; R, reactor; HC, holding coil; SV, selection valve; ADH, alcohol dehydrogenase enzyme solution; GDU, gas diffusion unit; λ, spectrophotometer, 340 nm; W, waste.

Flow protocol sequence

	SV position	Time (s)	Volume (μL)
Aspirate buffer solution	1	6	100
Aspirate enzyme solution	2	1.5	25
Aspirate NAD solution	3	1.5	25
Aspirate the buffer solution	4	6	100
Flow reversal and propel to GDU	7	12	192
Washing of HC	8	10	200
Aspirate the sample	5	2	30
Flow reversal and propelling sample through the donor channel of GDU	6	50	800
Propelling the reaction zone towards detection	7	80	2200

Analytical characteristics

Parameter	Value
Application zone	0-25% (v/v)
Repeatability RSD (n=10)	3.43% (0.5% (v/v)) 0.89% (12.5% (v/v))
LOD	0.2% (v/v)
LOQ	0.3% (v/v)
Determination rate	21 h ⁻¹
Reagent consumption	
ADH	0.45 U/assay
NAD ⁺	1.3 mg/assay
Sample	32 μL/assay
Waste produced	3.4 mL/assay

ANALYSIS OF WINE SAMPLES

Results obtained in the analysis of certified reference materials

Samples	Concentration, % (v/v)		RD, %	Recovery %
	Certified	Found		
LGC 5405	15.04 ± 0.05	15.33 ± 0.07	2.0	101.9
BCR 613	0.539 ± 0.0095	0.576 ± 0.020	6.9	106.9

Results obtained in the analysis of wine samples of different origin

Samples	Concentration of ethanol, % (v/v)		Dilution	R.D., %
	Ref. met.	SIA*		
White table wine	10.5 ± 0.1	10.3 ± 0.1	No	-1.9
White table wine	11.0 ± 0.1	11.1 ± 0.3	No	0.9
Red table wine	12.0 ± 0.1	10.9 ± 0.4	No	-9.2
Red table wine	12.5 ± 0.1	12.6 ± 0.3	No	0.8
White table wine	12.5 ± 0.1	12.5 ± 0.1	No	0
Red table wine	12.5 ± 0.1	12.5 ± 0.1	No	0
Port wine, White	19.0 ± 0.1	19.1 ± 0.5	3.5	0.5
Port wine, Tawny	20.2 ± 0.1	20.6 ± 0.4	3.5	2.0
Port wine, White	19.2 ± 0.1	19.1 ± 0.3	3.5	-0.5
Port wine, Ruby	19.0 ± 0.1	20.8 ± 0.5	3.5	9.5

* Mean and standard deviation (n=3); R.D., Relative deviation




Figure 2. Effect of the dilution on the apparent concentration of a port wine sample.

Conclusions

- 1 The complex matrix affects the gas diffusion process.
- 2 The difficulties can be overcome by appropriate sample dilution.
- 3 Achieved precision is more adequate for a screening method or fermentation monitoring.

References

- [1] OIV, Recueil des methodes internationales d'analyse des vins et des mouts, OIV (2005).
[2] M. D. Luque de Castro, J. González-Rodríguez, P. Pérez-Juan, *Food Reviews International*, 21:231-265, 2005.