

Determination of tocopherol and carotenoid contents in ST muscle of suckling lambs using fresh or lyophilised muscle

M. Blanco, F. Molino, M. García-Durillo, G. Ripoll, S. Lobón, A. Sanz, M. Joy
Centro de Investigación y Tecnología Agroalimentaria de Aragón (Zaragoza, Spain). e-mail: mblanco@aragon.es



A percentage of lamb consumers demand meat from grass-fed animals, however, they want guarantees of the animal's diet. Carotenoids can be used to trace forage-feeding but the analysis is expensive. Lyophilisation might reduce the time and quantities of reagents used in carotenoids and tocopherols analyses.

Objective: determine the content of carotenoids and tocopherols in meat of suckling lambs of grazing and indoors systems in lyophilised and fresh meat

Materials and Methods

39 Churra Tensina pairs of ewe-lamb
From lambs' birth to slaughter (10-12 kg LW)



Grazing



Hay

Ewes: + 300 g/d concentrate

Sampling:

▪ feedstuffs: weekly

▪ *Semitendinosus* (ST) muscles

Frozen at -80 °C

Lyophilised & frozen at -80 °C

Extraction:

▪ feedstuffs with acetone

▪ ST muscle



Lyophilised Fresh

	Lyophilised	Fresh
Muscle, g	0.1	0.4
Ethyl alcohol, ml	0.1	0.4
n-hexane, ml	1	1
n-hexane, ml	1	1

Analyses of carotenoids and tocopherols: HPLC

Results


Feedstuffs

Content, µg/g DM	Pasture	Hay	Concentrate
β-carotene	723	23.7	n.d.
Lutein	570	100	1.0
α-tocopherol	147	10.9	9.3
γ-tocopherol	6.6	1.1	6.2


Greater contents

n.d.: not detected

Lyophilisation increased α-tocopherol and lutein contents



Effect of the type of forage



Conclusions

Lutein and α-tocopherol contents in muscle could be used as markers of grazing in ovine. The determination in lyophilised muscle was more appropriate than in fresh muscle