

University of Arkansas, Fayetteville
ScholarWorks@UARK

Mechanical Engineering Undergraduate Honors
Theses

Mechanical Engineering

8-2007

Experimental analysis of plasma conditions in an inverted cylindrical magnetron sputtering system

Bradley Handloser

University of Arkansas, Fayetteville

Follow this and additional works at: <http://scholarworks.uark.edu/meeguht>

Recommended Citation

Handloser, Bradley, "Experimental analysis of plasma conditions in an inverted cylindrical magnetron sputtering system" (2007).
Mechanical Engineering Undergraduate Honors Theses. 24.
<http://scholarworks.uark.edu/meeguht/24>

This Thesis is brought to you for free and open access by the Mechanical Engineering at ScholarWorks@UARK. It has been accepted for inclusion in Mechanical Engineering Undergraduate Honors Theses by an authorized administrator of ScholarWorks@UARK. For more information, please contact scholar@uark.edu.

Experimental Analysis of Plasma Conditions in an Inverted Cylindrical Magnetron Sputtering System

Honors Thesis submitted in Partial Fulfillment of the
Requirements for the Honors Program for the
Degree of Bachelor of Science
In Mechanical Engineering

2007

Bradley Scott Handloser

Project Advisor – Dr. Matt Gordon

Mechanical Engineering
University of Arkansas, Fayetteville

Table of Contents

List of Tables.....	iii
List of Figures.....	iv
Chapter One	
1.0 Abstract.....	1
1.1 Introduction.....	1
Chapter Two	
2.0 Physical Vapor Deposition.....	3
2.1 Magnetron Sputtering System.....	3
2.2 The Inverted Cylindrical Magnetron System.....	4
Chapter Three	
3.0 Experiments.....	6
3.1 The Deposition Process.....	7
3.2 Mass Spectrometer.....	8
3.3 Langmuir Probe.....	12
3.4 Theory.....	13
3.5 The Top Plate Design.....	16
Chapter Four	
4.0 Results and Discussion.....	18
4.1 EQP Data.....	21
4.2 ESPION Data.....	29
Chapter Five	
5.0 Conclusions.....	38
5.1 Future Work.....	39
Works Cited.....	41

List of Tables

Table 1: Plasma recipe for alpha-alumina growth.....	7
Table 2: Low power oxygen plasma recipe.....	19
Table 3: Low power argon plasma recipe.....	17
Table 4: Oxygen plasma results.....	33
Table 5: Oxygen plasma EEDF results.....	36
Table 6: Argon plasma results.....	37
Table 7: Argon plasma EED function results.....	37

List of Figures

Figure 1: Plasma created in cylindrical magnetron.....	1
Figure 2: Diagram of magnetic fields in inverted cylindrical magnetron system using AC power.....	5
Figure 3: Experimental sputtering system, ICM-10 made by Isoflux.....	7
Figure 4: Schematic of EQP system by Hiden Analytical, Inc.....	9
Figure 1: Electron impact ionization by hot filament.....	10
Figure 6: Quadrupole mass filter, z-axis propagation.....	10
Figure 7: Qaudrupole mass filter Mathieu Stability Diagrams.....	11
Figure 2: Schematic of ESPION Langmuir probe made by Hiden Analytical, Inc.....	12
Figure 9: Typical I-V curve generated by Langmuir probe.....	13
Figure 10: Top plate design 1.....	17
Figure 11: Blank design.....	17
Figure 12: Top plane design 2.....	18
Figure 13: Experimental system, EQP on ICM-10.....	20
Figure 14: No plasma, RGA mass scan (0-60 amu).....	21
Figure 15: Oxygen plasma, SIMS positive ion mass scan (0-60 amu).....	22
Figure 16: Oxygen plasma, RGA mass scan (0-16 amu).....	23
Figure 17: Oxygen plasma, RGA mass scan (16-30 amu).....	24
Figure 18: Oxygen plasma, RGA mass scan (35-60 amu).....	24
Figure 19: Argon plasma, RGA mass scan (0-60 amu).....	25
Figure 20: Oxygen plasma, RGA energy scan.....	26
Figure 21: Oxygen plasma, SIMS positive ion mass scan (0-30 amu).....	27

Figure 22: Oxygen plasma, SIMS positive ion mass scan (35-60 amu).....	28
Figure 23: Oxygen plasma, SIMS positive ion energy scan.....	28
Figure 24: Oxygen plasma at 3 targets (I-V curve).....	30
Figure 25: Example of typical I-V, first, and second differential curves.....	31
Figure 26: Example of typical EED function.....	32
Figure 27: Oxygen plasma, 8 scan avg. (I-V curve).....	32
Figure 28: Oxygen plasma, smooth curve (I-V curve).....	33
Figure 29: Oxygen plasma, 1st differential.....	34
Figure 30: Oxygen plasma, 2nd differential.....	35
Figure 31: Oxygen plasma EED function.....	35
Figure 32: Argon plasma, smooth curve (I-V Curve).....	36
Figure 33: Argon plasma, EED function.....	37

Chapter One

1.0 Abstract

In this report, we will investigate Langmuir probe and mass spectrometer plasma diagnostic systems used to characterize conditions within a magnetron sputtering system. The specific systems used were the *ESPION* Langmuir probe and EQP mass spectrometer by Hiden Analytical and the ICM-10 inverted cylindrical magnetron sputtering system by Isoflux, Inc. This paper will show how these systems work, are operated, and produce data. Low power oxygen and argon plasmas were used through this research for diagnostics. With the EQP, two modes were used to collect data: RGA for neutrals and radicals and SIMS for positive ions. With both modes, mass and energy scans were produced. With the *ESPION*, current vs. voltage (I-V), first differential, second differential, and electron energy distribution functions (EEDF) were created. Also I-V curves at different locations within our sputtering system were made and compared. This preliminary research will help to lay the ground work for others to produce a true characterization of our system and so that in the future correlations can be made between the plasma conditions present and the quality of the films deposited.

1.1 Introduction

Magnetron sputtering has become a widely used technique for the deposition of high quality thin films and coatings. Some common applications include thin film solar cells, magnetic recording media, metallic interconnects, high wear resistive coatings, Aluminum coatings and functional glass coatings [1]. The process can deposit a wide variety of materials including metals, metal alloys, nitrides, carbides and oxides can be

deposited [13]. One reason for the success of magnetron sputtering is its unique ability to retain specific design requirements in a cost-effective manner. Another reason magnetron sputtering is among the most widely used physical vapor deposition (PVD) processes is because it is seen as environmentally friendly due to the fact that it avoids using harsh chemical reactants [5].

Typically, magnetron sputtering sources are operated between a frequency of 10 to 350 kHz and a duty cycle from 50% to 90%. A low pressure, around 0.5 Pa or lower, is also an essential part of the process [11]. In order for the



Figure 1: Plasma created in cylindrical magnetron [10]

sputtering to occur within the system, a plasma, or partially ionized gas, must be created. This is done by charging the low-pressure gas with an electrical field which causes the gaseous particles to lose an electron or become ionized. These plasmas consist of various neutral atoms, ions, radicals, and free electrons. The positively charged ions, along with energetic neutrals, are then accelerated via a strong magnetic field into the target consisting of a desired material. When these highly energized particles collide with the target, atoms are ejected or sputtered from its surface. These sputtered particles, consisting of a specific material, eventually condense and react with the substrate to become a thin film layer [4]. The external parameters of the sputtering system including frequency, power, pressure, system configuration, gas or gases used, and target material govern the characteristics of the plasma and ultimately the films created. In order to analyze the plasmas being created, diagnostics tools such as Langmuir probes and mass spectrometers are used.

In this thesis, ways by which the magnetron plasma would be characterized by our research group is considered and a set of preliminary data has been reported here to illustrate what set of parameters can be investigated. Though the data reported here is not ideal, they were not taken at deposition conditions due to temperature issues, it gives the group a fair understanding of what to expect at high powers.

Chapter Two

2.0 Physical Vapor Deposition

Physical vapor deposition (PVD) is a process that produces a certain atomic, molecular, or ionic species that are transported and condensed on a substrate. The three fundamental PVD processes include sputtering, evaporating, and ion plating. The evaporation process requires that a material be vaporized by thermal heat at an adequately low pressure. The vapor then can move freely and condense on the substrate. Sputtering occurs when a solid surface is bombarded by highly energized particle causing surface atoms to be knocked free. Ion plating uses aspects of both the evaporation and sputtering processes [12].

2.1 Magnetron Sputtering System

This is a fine technique of depositing thin-films onto substrates. Unlike evaporation, the materials to be sputtered do not have to be heated. The material is in the form of a thin plate called a target. The target is mounted in a chamber filled with Argon gas at low pressure. Positive argon ions are produced in the plasma (formed when the target is

connected to negative high voltage) by collisions of electrons with argon atoms that are accelerated toward the negatively-charged target. Upon impact with the target, free atoms of the target material are expelled and deposited on substrates mounted in front of the target, building a thin-film. Unlike ordinary sputtering processes, all the ions that attempt to escape the target, are captured behind/at the sides of the targets and are confined to the immediate vicinity of the target by strong permanent magnets. This increases the ion current resulting in a fast deposition. The target temperature is low thus resulting in a high quality film deposited.

2.2 The Inverted Cylindrical Magnetron System

The three major sputtering deposition techniques are DC (diode), RF (radio frequency), and AC (alternating current). The sputtering system used for our experiments, the ICM-10 made by Isoflux Inc., is an inverted cylindrical magnetron sputtering system that uses an AC power supply [12]. The system is made of two hollow cylindrical targets that completely surround the substrate. This configuration is advantageous over planar magnetrons because it has the ability to coat complex 3-D shapes. The system also implements a patented unbalanced magnetic field. Each target has a set of magnetic poles that is stronger than the other causing the magnetic field lines to focus towards the center. The unbalanced magnetic field produces extremely high plasma densities, much higher than conventional DC sputtering systems. The ICM-10's AC power supply enables its electrically isolated targets to each switch back and forth between being a cathode and anode. This allows the system to naturally produce high deposition rates without damaging arcing. In most other configurations, complex wiring is needed to

avoid this problem. The system is also capable of producing current densities much higher than other systems at similar target power and bias voltages (Glocker ICM10 paper).




Figure 2: Diagram of magnetic fields in inverted cylindrical magnetron system using AC power [6]

Our specific system uses three targets approximately 33 cm in diameter and 7.5 cm tall. The three targets are stacked on top of one another with small spacers separating each of them. The middle of the three targets is dead or inactive and also acts as a spacer. Like most systems, it operates at an extremely low pressure (2mTorr) which is essential for plasma production. The system is operated on a PLC platform and RSView 32 HMI interface. The system uses the PEII AC power supply made by Advanced Energy. This power supply operates at 41 kHz and has a maximum power of 10 kW. The DC bias power supply to the substrate can also reach a maximum of 5 kW. Argon, Oxygen, and Nitrogen are the main gases used in the ICM system. Ar is always used first to ignite the

plasma before the other gases are introduced. The two active targets are made of aluminum and the walls are water-cooled to control temperature.

One of the primary goals for our ICM-10 system is to deposit α -alumina thin films on steel and industrial tools. This is a highly sought after goal in the field of Surface Engineering because of the outstanding properties of aluminum oxide films in metal cutting. . The conventional alumina deposition process requires a substrate temperature of 750-1000°C which is harmful to temperature-sensitive steels (degrade at 500-700 °C). Low temperature ($\leq 500^\circ\text{C}$) α -alumina deposition has proven to be difficult; however, preliminary data from the inverted cylindrical magnetron sputtering system indicate possible success.

Chapter Three

3.0 Experiments

In this section, the alumina thin film growth process, new plate designs and the plasma diagnostics data acquisition process with the state-of-the-art tools: Langmuir probe and EQP are mentioned briefly.



Figure 3: Experimental sputtering system, ICM-10 made by Isoflux

3.1 The Deposition Process

The ICM system is computer controlled by using a spreadsheet-like program to run the system. The required parameters for any growth process are entered into the spreadsheet (table 1). All parameters listed in table one can be altered within the limits of the equipment. These parameters eventually determine the quality of the material deposited, thus careful adjustments will result in quality films. For instance below is of the recipe one of the recipes used in the attempted creation of α -alumina plasma.

Table 1: Plasma recipe for alpha-alumina growth

	Step 1	Step 2	Step 3	Step 4
Argon (0 - 200 sccm)	150	150	80	45
Oxygen (0 - 200 sccm)	0	0	0	45
MDX Mode (I/V/P)	P	V	V	V
MDX (I = 0-10, V = 0-1000, P = 0-5000)	100	200	50	35
PE II Mode (I/V/P)	P	P	P	P
PE II (I = 0-35, V = 0-1100, P = 0-10000)	500	1000	5000	5000
Shutter (Open/Closed)	Closed	Closed	Closed	Open
Step Time (Seconds)	60	300	300	7200

	Step 5	Step 6	Step 7	Step 8
Argon (0 - 200 sccm)	150	45	150	45
Oxygen (0 - 200 sccm)	0	45	0	45
MDX Mode (I/V/P)	V	V	V	V
MDX (I = 0-10, V = 0-1000, P = 0-5000)	0	35	0	35
PE II Mode (I/V/P)	P	P	P	P
PE II (I = 0-35, V = 0-1100, P = 0-10000)	5000	5000	5000	5000
Shutter (Open/Closed)	Closed	Open	Closed	Open
Step Time (Seconds)	300	7200	300	7200

The primary function of steps 1 and 2 are to ignite the plasma. Only Ar gas is present in this step and the shutter to the target is closed. The shutter is a circular wall that prohibits any sputtered particles from getting to the target for deposition. The MDX, which is the bias on the substrate, is also warming up. In step 3, there is still no sputtering but a significant increase in power is initiated to warm up the power supply. Step 4 marks the beginning of sputtering and deposition. The power is sufficiently high at 5000V, Ar and O₂ are now present, the shutter to the targets is open, and the time is increased to 7200s. This deposition step is then followed by a cleaning step. The shutter is now closed and there is no bias on the substrate to stop deposition. Again only Argon plasma is present which continues the sputtering process; this essentially cleans the Aluminum targets. Steps 4 and 5 are then repeated until the desired thickness of the film is achieved.

3.2 Mass Spectrometer

Mass spectrometers are used in plasma diagnostics to measure characteristics of individual molecules present in a system. Their analysis includes a combined Mass/Energy measurement for positive and negative ions, neutrals, and radicals [3]. Neutrals and radicals must be ionized before sampling while the positive and negative ions are directly sampled. All of the species must be charged so they can be moved and

filtered by external electric and magnetic fields for sampling [12]. With this data, an energy resolved mass spectra of the molecules present and created within the plasma can be produced.

The mass spectrometer used for plasma diagnostics in our experiments was the HAL Electrostatic Quadrupole Plasma (EQP) Analyser built by Hiden Analytical Limited. The EQP is a “high-transmission 45° sector field ion energy analyzer and quadrupole mass spectrometer” [3]. It measures the energy and mass-to-charge ratio distributions of positive and negative ions, neutrals, and radicals generated in the source plasma. The EQP system contains ion extractor optics, ioniser, sector field energy analyzer, triple filter mass analyser, pulse ion counting detector, and MASsoft PC data system for control and analysis (EQP/EQS Analyser Operator's Manual, 1-5).




Figure 4: Schematic of EQP system by Hiden Analytical, Inc. [3]

There are two modes of analysis available with the EPQ system: Secondary Ion Mass Spectrometry (SIMS) and Residual Gas Analysis (RGA). The SIMS mode analyzes the mass and energy of ions generated in the plasma while the RGA mode analyzes the mass and energy of neutrals and radicals. Positive and negative ions can be simply extracted by the EQP system because they already possess a charge. For RGA, the neutrals and radicals must be ionized by the ionization source located closely behind the extractor opening. Ionization is achieved by electron impact (EI) via thermionic emission from a hot filament. Molecules are normally converted to positively charged particles in this process [10].




Figure 5: Electron impact ionisation by hot filament [10]




Figure 6: Quadrupole mass filter, z-axis propagation [10]

All ions, both SIMS and RGA mode, are then focused into an electrostatic energy field analyser, a triple-section quadrupole mass filter, and finally a detector [3]. The electro static energy field analyser is a 45° sector field which measures the curves of ions to determine their kinetic energies [1].

The triple-section quadrupole is a mass filter that differentiates ions present and produced for detection. The quadrupole is comprised

of two pairs of parallel metal rods of equal length. They contain fixed direct current (DC) and alternating radio frequency (RF) components [10]. One pair consist of a DC component (V_{DC}) at an opposite potential compared to the other, while the other pair has an RF ($V_{RF} \cos \omega t$) component 180° out of phase of its adjacent rod [1]. All ions entering the quadrupole are deflected from their original trajectory and those without the exact mass : charge (m/e) ratio desired are deflected and neutralized upon impact on the rods. The change in RF enables only one certain m/e ratio at a time to oscillate down the z-axis.




Figure 7: Quadrupole mass filter Mathieu Stability Diagrams [10]

The desired ions that make it through produce an ion current on the detector. The detector used by the EQP is a secondary electron multiplier (SEM). When an ion impacts the surface of the SEM, 2 or 3 electrons are generated. These electrons are further multiplied in a cascade affect, via an applied voltage, when they collide with the surface. The multiplier effect normally produces a magnitude of 10^3 greater than the actual count

[10]. The output from the detector is given in arbitrary units of ‘counts per second.’ This output is normally paired in a graph versus either mass or energy because the m/e ratio of the desired ion is known [1].

3.3 Langmuir Probe

The Langmuir Probe diagnostic tool was first developed by Mott-Smith and Langmuir in 1926 [1].




Figure 8: Schematic of ESPION Langmuir probe made by Hidden Analytical, Inc. [7]

The probe is a small electrode consisting of a conductive material that is connected in series to a variable voltage source. This source is necessary to bias the probe positively or negatively against the plasma. When placed in the plasma, the electrons and ions present produce a current in the probe. The current drawn from the plasma is then measured as a function of the bias voltage. This is known as the current-voltage (I-V) characteristic [13]. From the I-V curve, several basic plasma properties can be obtained including ion density (n_i), electron density (n_e), electron temperature distribution (T_e), plasma potential (V_p), and electron energy distribution function (EEDF) [7].

The plasma potential (V_p) is the curve located between the two lines in the upper right region (knee) of the I-V curve. At this location, the probe potential (V_B) is equal to the plasma potential. Mostly electron current is collected here due to the fact that electrons are considerably lighter and more mobile than other species (i.e. ions, neutral, and radical). The region to the right of the plasma potential with the slightly positive slope is known as the electron saturation region.




Figure 9: Typical I-V curve generated by Langmuir probe [7]

The Langmuir probe used in our experiments was the *ESPION* made by Hidden Analytical, Inc. It can measure ion (n_i) and electron density (n_e) over the range of 10^{14} to 10^{19} m^{-3} and electron temperature (T_e) up to 10 eV. It has a voltage range of -200V to 100V and a current range of 20 μ A to 1A. Our system has the optional automatic z – translator capable of taking measurement

over 300mm in the z-direction. The probe is RF compensated and is designed for high frequency plasma. Our probe tip is made out of titanium rather than the standard tungsten to stop oxidation from O_2 plasma and the probe area is 4.72 mm. It can acquire, analyzes, and present data as graphs or tables with the use its software program ESPsoft made by Hidden. Analysis of the I-V curve in this program is based on the “Orbital Motion Limited” theory which effectively increases the diameter of the tip wire for ions [2].

3.4 Theory

$$I_{es} = -\frac{en\bar{c}_e}{4} = en_e \sqrt{\frac{kT_e}{2\pi m_e}}$$

where: I_{es} is electron current at saturation

e (charge of an electron) = $1.602 \text{ E-}19 \text{ C}$

n is the number of particles

\bar{c}_e is the average electron velocity

n_e is the number of electrons

k (Boltzman's constant) = 1.38 E-23 J/K

T_e is the electron temperature

m_e (mass of an electron) = 9.11 E-31 kg

The region to the right of the plasma potential with the slightly positive slope is known as the electron saturation region.

When the probe has a slightly negative potential compared to the plasma, some of the electrons are repelled. This causes the exponential curve below the plasma potential.

$$I_e = I_{es} \exp\left(\frac{e(V_p - V_B)}{kT_e}\right)$$

where: I_e is the electron current

V_p is the plasma potential in volts

V_B is the probe potential in volts

T_e is the electron temperature in Kelvin

From this, $\ln(I_p + I_i)$ can be plotted against V_p to produce a straight line whose slope represents the electron temperature.

The point on the curve that crosses the x-axis (V) is the floating potential ($V_B = V_f$). Here only a small portion of the electrons can overcome the bias potential of the probe. The net current is equal to zero because the electron current is compensated by the ion current.

Finally, when the probe possess a negative potential high enough that no electrons can reach the probe, the ion current (I_{is}) can be measured. Here, only the ion flux is measured by the probe. I_{is} is constant when the bias potential is less than or equal to zero with respect to the plasma [13].

$$I_{is} = eA_p n \exp\left(-\frac{1}{2} * \sqrt{\frac{kT_e}{M_e}}\right)$$

where: A_p is the area of the probe tip = 4.72mm

With this information, the plasma current (I_p) can now be found. This is the current corresponding to the plasma potential (V_p) on the I-V curve.

$$I_p = I_e - I_i = I_{es} \exp\left(\frac{e(V_p - V_B)}{kT_e}\right)$$

The ion flux is another characteristic of plasmas that can be determined with a Langmuir probe. Flux is defined as the number per area per time.

$$\Gamma_i = \frac{n\bar{c}_i}{4}$$

where: Γ_i is ion flux

$$\bar{c}_i \text{ (avg. ion velocity)} = \sqrt{\frac{8KT}{\pi m_i}}$$

The characteristic distance in which a charged particle in the plasma is shielded from external electric fields can also be found. This is known as the Debye length.

$$\lambda_D = \sqrt{\frac{\epsilon_0 KT}{n_e e^2}}$$

where: ϵ_0 (permittivity of free space) = 8.854 E -12 F/m

n_i is the ion density

[8]

The EEDF can also be calculated. This is done by substituting the second differential values of the I-V curve into the Druyvesteyn equation (ESPsoft User Manual, 8-5)

$$F(\varepsilon) = \frac{4}{e^2 A} \left(\frac{-mV}{2e} \right)^{\frac{1}{2}} \frac{d^2 I}{dV^2}$$

where: $F(\varepsilon)$ is the electron energy distribution function (EEDF)

V = potential (V)

A = current (A)

$\frac{d^2 I}{dV^2}$ are the second differential values

3.5 The Top Plate Design

The top plate of the inverted cylindrical magnetron sputtering system (ICM-10) used in our experiments was not equipped for the state-of-the-art diagnostic instruments (ESPION Langmuir probe and EQP mass spectrometer) used to analyze the plasma in the ICM-10. Some design work was required to enable these systems to be utilized. Both the ESPION and EQP require a probe to be inserted into the plasma for measurement. To make this possible, a new solid top plate had to be designed for the system.




Figure 10: Top plate design 1

The first top plate built and designed includes seven openings in which both the *ESPION* and *EQP* (via an adapter) can be inserted and attached. There is one hole in the center and six around the edge progressively one-fifth of an inch further from the outside wall (where the sputtering occurs).




Figure 11: Blank design

Blanks were also designed to occupy the holes not being used to keep the essential vacuum in the ICM.

Unfortunately, the outer holes were too close to the walls/shutter to be safely used with either of the expensive and fragile diagnostic systems. We are currently working on a new top plate with three holes 2", 3", and 4" from the center of the target to the center of the opening. This gives us a clearance of 2" from the shutter on the outermost hole.




Figure 12: Top plate design 2

Air cooling has also been added to the *ESPION* system. The first probe used in our system malfunctioned because the temperature inside exceeded the probe's limit. The inductors inside the probe melted at about 300°C. Air pressure at 3-10 psi now runs down the probe during operation to help control the temperature.

Chapter Four

4.0 Results and Discussion

The parameters of the ICM-10 magnetron sputtering system determine the plasmas created within it and ultimately the thin films deposited on the substrate. The specific parameters of the sputtering system are known as a recipe to create each plasma. Two recipes were used throughout this research for plasma characterization by both the EQP and *ESPION* diagnostic systems. The first was a recipe for low power oxygen plasma

and the second was for low power argon plasma. Low power plasmas had to be used to keep the temperatures sufficient low to prohibit damaging the diagnostic systems. Below are the recipes used to create these plasmas.

Table 2: Low power oxygen plasma recipe

	Step 1	Step 2
Argon (0 - 200 sccm)	80	0
Oxygen (0 - 200 sccm)	45	45
MDX Mode (I/V/P)	V	V
MDX (I = 0-10, V = 0-1000, P = 0-5000)	35	0
PE II Mode (I/V/P)	P	P
PE II (I = 0-35, V = 0-1100, P = 0-10000)	500	500
Shutter (Open/Closed)	Closed	Closed
Step Time (Seconds)	60	3600

Table 3: Low power argon plasma recipe

	Step 1	Step 2
Argon (0 - 200 sccm)	120	60
Oxygen (0 - 200 sccm)	0	0
MDX Mode (I/V/P)	P	P
MDX (I = 0-10, V = 0-1000, P = 0-5000)	0	0
PE II Mode (I/V/P)	P	P
PE II (I = 0-35, V = 0-1100, P = 0-10000)	500	500
Shutter (Open/Closed)	Closed	Closed
Step Time (Seconds)	60	3600

Step 1 for both the oxygen and argon plasma is an ignition step. In the oxygen plasma, both oxygen and argon gas are present because it is difficult to initiate oxygen plasma on its own. Throughout both processes, the power is low (0.5 kW) to prevent high temperatures and the shutter is closed to prevent sputtering. In this preliminary research, our concern is with the plasma characteristics, not the sputtered particles. In the second step, only the gas for the plasma desired (oxygen or argon) is present. The chamber and convectron pressures were roughly 7.3×10^{-3} torr, while the foreline pressure remained at

about 9.3×10^{-2} torr throughout both processes. All data for the oxygen and argon plasmas were taken during the second step.

4.1 EQP Data



Figure 13: Experimental system, EQP on ICM-10

One of the useful benefits of the EQP is that it can show exactly what particles are present within your system and if there are any leaks in the vacuum or gas lines. Even when our ICM is pumped down to 10×10^{-7} torr, some particles are still present. Below is an RGA mass scan with no plasma at vacuum.




Figure 14: No plasma, RGA mass scan (0-60 amu)

It is obvious that even at vacuum, there are particles present. Oxygen gas (amu = 32) is by far the dominant species here. There is also some O^{16} (amu = 16) and argon gas (amu = 40) in the chamber. The argon gas is probably there because of either a leak in the argon gas line or some is left from previous argon plasma experiments. Either way, there is only a minute quantity present. O_2 is there because it is readily present in our atmosphere.

As stated in the theory, there are two modes of scanning that can be taken with the EQP. The RGA scan is for sampling neutrals and radicals and the SIMS scan is for positive and negative ions. Within each of those modes, an energy and mass scan can be taken simultaneously. The first set of scans to be examined is with low-power oxygen plasma. Because oxygen gas is being used to create the plasma, O_2 is the most common element in the system. It was necessary on the mass scans to break the scan into segments in order

to miss the 30-35 amu range because O_2 has a mass of 32 amu. If this range were measured (as shown below), the rest of the data range is dwarfed by the high oxygen count and the results are not detailed enough to be useful.




Figure 15: Oxygen plasma, SIMS positive ion mass scan (0-60 amu)

The RGA mass scan was broken into three segments (0-16, 16-30, & 30-60 amu) for greater detail and to purposely miss the 32 amu section.




Figure 16: Oxygen plasma, RGA mass scan (0-16 amu)

In the 0-16 amu range, it appears that some hydrogen gas (H_2 , amu = 2), carbon (C^{12} , amu = 12), and nitrogen gas (N_2 , amu = 14) are present. Even if so, their presence is negligible because their counts per second (c/s) is so small (2×10^3) compared to some of the major elements and compounds present (7×10^6 c/s).




Figure 17: Oxygen plasma, RGA mass scan (16-30 amu)

In the 16-30 amu range, some water (H₂O, amu = 18) and nitrogen gas (N₂, amu = 28) are present. However, there is far more Ar²⁰ (amu = 20) than anything else in that range.




Figure 18: Oxygen plasma, RGA mass scan (35-60 amu)

The last range from 35-60 amu contains Argon (amu = 40) which has a c/s of a whole order of magnitude greater than any other. There also appears to be some carbon dioxide (CO_2 , amu = 44), but it is hard to see because it has such a smaller count compared to Ar. For comparison, here is a RGA scan of argon plasma created with the recipe described earlier.




Figure 19: Argon plasma, RGA mass scan (0-60 amu)

Argon gas is the dominant species as expected followed by oxygen gas and Ar^{20} . Argon had a similar count in both the oxygen plasma and argon plasma RGA scans at just under $7 \text{ E}+6 \text{ c/s}$. Obviously there is much less oxygen gas in the argon plasma as in the oxygen plasma (which was not shown). Both scans also showed particles present at a mass of 20, 18, and 16 amu.

The RGA oxygen plasma energy scan below shows the energy in Volts against the c/s so that one can determine the energy of an element or compound by matching its c/s in the mass scan.




Figure 20: Oxygen plasma, RGA energy scan

For example, the peak of this RGA is at about 6.3×10^5 c/s which corresponds to Ar^{20} .

We can now see that the charge on that compound is about 1.25 volts.

A SIMS positive ion scan was also done on this same oxygen plasma. The mass scan was broken into two parts, 0-30 and 35-60 amu, again to miss oxygen gas ions at 32 amu.




Figure 21: Oxygen plasma, SIMS positive ion mass scan (0-30 amu)

Obviously, oxygen (O^{16} , amu = 16) is the most abundant species in this range. This is followed by what appears to be aluminum (Al^{27} , amu = 27), which could have possibly sputtered of the aluminum target and made its way around the shutter. Still it has a very low c/s and therefore is relatively unimportant.




Figure 22: Oxygen plasma, SIMS positive ion mass scan (35-60 amu)

Again, as in the RGA mass scan, Argon is the most abundant species, other than oxygen of course. Because this is a positive ion scan, this shows that argon had the most negatively ionized atoms in this plasma.




Figure 23: Oxygen plasma, SIMS positive ion energy scan

It is easily recognized that the SIMS ion energy scan has a much higher energy than the RGA scan which makes sense because ions possess a charge while the neutrals don't. Also, from this graph one can see that the energy is coming from elements and compounds with low c/s in the 0-70,000 range. This data helps to determine what species and the amount present for a specific plasma.

4.2 ESPION Data

As discussed in the theory, a Langmuir probe measures current produced in a plasma as a function of the probes bias voltage. This produces a current-voltage (I-V) curve which is used to find the plasmas' potential (V_p), floating potential, electron temperature (T_e), electron density (N_e), ion density (N_i), Debye length (λ_D), ion flux, and electron energy distribution function (EEDF). Our ESPION Langmuir probe uses the ESPsoft program both made by Hiden to collect and process the data taken. This ESPION also has a translating z-drive that is capable of taking measurements over a range of 30cm, which covers the entire height of our sputtering system.

In this first set of data, oxygen plasma was sampled at three different locations within our system. The three locations represent the centers of the three targets (aluminum).




Figure 24: Oxygen plasma at 3 targets (I-V curve)

The center of the lower target is about 25 cm from the top, the middle target (not active) is roughly 16.5 cm from the top, and the upper target is about 5.1 cm from the top. The graph produced represents the three I-V curves of the oxygen plasma at each location.

While this data is odd because of the jagged spikes on the upper right portion of the graph, it is very interesting because of the obvious differences at each target center. It is believed that the plasma should be uniform through the cylindrical sputtering system.

The graph below shows what kind of results one would expect to see from an I-V curve and its first and second differential.




Figure 25: Example of typical I-V, first, and second differential curves [3]

The regular I-V curve has an “I” to its right and is a smooth curve that has a slightly positive slope just below the x-axis (voltage) and turns up into an s-shaped exponential curve as it crosses and normally flats back out again after the floating potential. The second differential has the “ dI/dV ” to its left and follows closely to the I-V curve until it spikes and sharply turns down. Lastly, the second differential curve with the “ d^2I/dV^2 ” to its left again follows closely to the I-V curve’s path but spikes lower than the first differential and has a minimum spike below the x-axis. The second differential is used to produce a EED function which is represented in example below.




Figure 26: Example of typical EED function [3]

To produce more accurate and consistent data, we took multiple scans at one time, averaged them, and produced a smoothed curve with the ESPION and ESPsoft system.




Figure 27: Oxygen plasma, 8 scan avg. (I-V curve)




Figure 28: Oxygen plasma, smooth curve (I-V curve)

The first graph is an average of 8 scans of oxygen plasma taken at one time and the second is of the same data with a smooth curve. The graph produced is fairly consistent with what one would expect to see on an I-V curve with the exception of the small minimum and maximum spikes near the origin. The cause of this has not yet been discovered and will be continued to be researched. The I-V curve produced from the oxygen plasma recipe gave us these experimental plasma conditions:

Table 4: Oxygen plasma results

Floating Potential	-14.88	V
Plasma Potential (Vp)	23.33	V
Electron Temperature (Te)	4.42E+04	K
Electron Temperature (Te)	3.81	eV
Electron Density (Ne)	1.41E+17	m ⁻³
Ion Density (Ni)	2.05E+17	m ⁻³
Debye Length (λ_D)	3.87E-05	m
Ion Flux	6.18E+20	m ² *s

From the smooth I-V curve, a graph of the first differential was produced.




Figure 29: Oxygen plasma, 1st differential

This graph follows a very unusual path. While the beginning and end of the graph display what one would expect to see, the middle is full of unexpected peaks and valleys that are most likely a result of improper data acquisition or damages sustained to our probe at high power.

A differential was taken again from this graph to produce the second differential curve which again displayed unexpected results.




Figure 30: Oxygen plasma, 2nd differential

There should be only one major maximum and minimum, not two or three. Hopefully, as research continues, better data will be collected to produce more clear results.

Nevertheless, an EEDF function and its results were produced from this data.




Figure 31: Oxygen plasma EED function

Table 5: Oxygen plasma EEDF results

Floating Potential	-14.88	V
EED (Ne)	7.43E+16	m ⁻²
EED (Avg. Energy)	19.7	eV
EED (Te)	13.2	eV

The same scans were done on argon plasma for comparison. Below are I-V and EED graphs of argon plasma with their respective results.




Figure 32: Argon plasma, smooth curve (I-V Curve)




Figure 33: Argon plasma, EED function

Table 6: Argon plasma results

Floating Potential	26.07	V
Plasma Potential	65.4	V
Electron Temperature (Te)	3.97E+04	K
Electron Temperature (Te)	3.43	eV
Electron Density (Ne)	2.74E+17	m ⁻³
Ion Density (Ni)	0	m ⁻³
Debye Length (λ_D)	2.63E-05	m
Ion Flux	0	m ² *s

Table 7: Argon plasma EED function results

Floating Potential	26.07	V
EED (Ne)	1.89E+17	m ⁻²
EED (Avg. Energy)	11.5	eV
EED (Te)	7.69	eV

The I-V curve produced by the argon plasma was much better and did not contain any unexpected spikes. This also resulted in the EED function being much smoother and more of what one would expect to produce.

Comparing the oxygen and argon plasma results shows us some similarities and differences between them. The plasma potential of the argon is more than double that of the oxygen plasma at 65.4 V compared to 23.33 V. Their electron temperatures were fairly similar with oxygen plasma having a high temperature of only 0.38 eV. The electron density of the argon plasma was again almost double that of the oxygen, but the argon plasma showed no ion density or ion flux which doesn't seem possible. The Debye length of the argon was also a third less than the oxygen's.

Chapter Five

5.0 Conclusions

The primary goals of my research were to create preliminary characterization of the plasma conditions in our inverted cylindrical magnetron sputtering system by learning how our plasma diagnostic systems work, how to operate them, and how to produce data. All of the data and results were experimental and have neither been proved nor disproved. Others that follow will continue with this research and follow this ground work. They will learn to perfect the data acquisition process, create useable data, interpret the data correctly, and produce a true characterization of our system.

Many possible factors could have contributed to results being skewed or false. One of the major possibilities is that both systems (ESPION and EQP) were not equipped to operate at the high temperatures produced in our system. All of the data for this report was produced at low power (.5 KW) and temperatures (<300°C) because both systems malfunctioned before high power (5 KW) could be achieved. Higher

temperature probes are being obtained and both systems will be incorporated with water cooling which should allow for higher power operation. Novelty and our unfamiliarity of these complex systems could have also played a role in imperfect data. This could have contributed to the noise that was seen within the data. This problem should decrease as more time is spent operating these plasma diagnostic systems.

Interpretation of the data is a final vital step in this research. Interpretation is very difficult and can be futile with preliminary data that is unproven. This will become more important as these diagnostic techniques are perfected through time and practice. As more data is collected consistently, interpretation will become more realistic as one can compare and contrast results with similar and different parameters. Finally, no one really knows what results should be seen within this system. That is the whole premise of this research. Rough assumptions can be made, but we will not know for sure until the diagnostic techniques are perfected and the data is interpreted properly. Once this is done, correlations can be made between the deposition parameters, the plasma characteristics present, and the quality of the material deposited. With this information, it will be possible to fully understand exactly what materials are being deposited and more importantly why. My research should help us to move towards this goal.

5.1 Future Work

Research will continue with all three systems (ICM-10, *ESPION* and EQP) discussed in this paper. The first thing to be done will be adding water-cooling to both the *ESPION* and EQP. Both systems are exceeding their maximum operating temperatures and can

not take measurements at high power. This is desired because the plasma needed to grow α -alumina is created at high power (5 kW). Up to this point, measurements have only been taken at lower power and temperatures. The *ESPION* probe was recently sent back to Hiden Analytical and the EQP is soon to follow.

Two new top plates will also be needed to accommodate the water-cooled mass spectrometer. Because the modified probe will possess a larger diameter, the new designs can include only one or two openings. The first new top plate will have only one large opening in center for the EQP. The second will include two openings near the outer wall, one for the EQP and one for the *ESPION*. The opening will be as close to the wall as possible (about one inch from the shutter).

Plasma diagnostic research will be continued on the ICM-10 with the *ESPION* and EQP. As more experiments are performed, a better understanding of how to perform the processes and collect good data will be gained. Once these diagnostic processes are perfected, the next step will be to learn to correctly interpret the data. With good data that can be properly interpreted, the final goal of correlating what specific parameters and plasma conditions result in what quality of deposition can be achieved.

Works Cited

- [1] Backer, Hermann. Time-Resolved Measurements of Pulsed Plasma for Thin Film Deposition. University of Manchester Institute of Science and Technology. 2003. 1-67.
- [2] ESPsoft Users Manual. Hiden Analytical. 2000. 8_5-8_8.
- [3] EQP/EQS Analyser Operator's Manual. Warrington: Hiden Analytical Limited. 1_5-5_7.
- [4] F, Debal, Bretagne J, and Jumet M. Analysis of DC Magnetron Discharges in Ar-N₂ Gas Mixtures. University of Mons-Hainaut, University of Paris-Sud. 1998.
- [5] F, Debal, Bretagne J, Wautelet M, Dauchot J P, and Hecq M. Analysis of DC Magnetron Discharges in Ar-N₂ Gas Mixtures. Comparison of a Collisional-Radiative Model with Optical Emission Spectroscopy. University of Mons-Hainaut, University of Paris-Sud. 1998. 219-220.
- [6] Glocker, David, D. Hofmann, A. Feuerstein, H. Schussler, M. M. Romach, and V. W. Lindberg. ICM10. 2001.
- [7] "Hiden ESPION Series Electrostatic Plasma Probes." Hiden Analytical. 5 June-July 2007 <www.hidenanalytical.com>.
- [8] Kruger, Charles H., and W. G. Vincenti. Physical Gas Dynamics. Malabar Florida: Krieger Co., 1967. 1-210.
- [9] "Mass Spectrometry of Processing Plasmas Using the Hiden EQP System." Hiden Analytical. 5 June-July 2007 <www.hidenanalytical.com>.
- [10] "Quadrupole Mass Spectrometry Concepts." Hiden Analytical. 5 June-July 2007 <www.hidenanalytical.com>.

- [11] Seo, Sang-Hun, and Jeon-Geon Han. Effect of Duty Cycle on Plasma Parameters in the Pulsed De Magnetron Argon Discharge. Korea Advanced Institute of Science and Technology. 2005.
- [12] Srivastava, Abhishek. Optical Emission Spectroscopy in an Inverted Cylindrical Magnetron Plasma. University of Arkansas. 2006. 3-12.
- [13] Vetushka, Alena. Study of the Plasma Parameters in a Bi-Polar Pulsed Magnetron Discharge. University of Manchester. 2006. 14-49.

Appendix A: Common atomic weights from Hiden Analytical, Inc.

**ATOMIC WEIGHTS OF THE ELEMENTS BASED
ON THE CARBON 12 STANDARD**

Symbol	Amu	Abundance	Symbol	Amu	Abundance
H ¹	1.0078	99.99	Ca ⁴⁰	39.9626	96.92
H ²	2.0141	0.02	Ca ⁴²	41.9586	0.64
He ³	3.0160	0.00013	Ca ⁴³	42.9588	0.13
He ⁴	4.0026	100.00	Ca ⁴⁴	43.9555	2.13
Li ⁶	6.0151	7.52	Ca ⁴⁶	45.9537	0.0032
Li	7.0160	92.48	Ca ⁴⁸	47.9525	0.18
Be ⁹	9.0123	100.00	Sc ⁴⁵	44.9559	100.00
B ¹⁰	10.0129	18.98	Ti ⁴⁶	45.9526	7.95
B ¹¹	11.0093	81.02	Ti ⁴⁷	46.9518	7.75
C ¹²	12.0000	98.89	Ti ⁴⁸	47.9479	73.45
C ¹³	13.0034	1.11	Ti ⁴⁹	48.9479	5.51
N ¹⁴	14.0031	99.64	Ti ⁵⁰	49.9448	5.34
N ¹⁵	15.0001	0.37	V ⁵⁰	49.9472	0.24
O ¹⁶	15.9949	95.77	V ⁵¹	50.9440	99.76
O ¹⁷	16.9991	0.04	Cr ⁵⁰	49.9460	4.31
O ¹⁸	17.9992	0.20	Cr ⁵²	51.9405	83.76
F ¹⁹	18.9984	100.00	Cr ⁵³	52.9406	9.55
Ne ²⁰	19.9924	90.92	Cr ⁵⁴	53.9389	2.38
Ne ²¹	20.9939	0.26	Mn ⁵⁵	54.9381	100.00
Ne ²²	21.9914	8.82	Fe ⁵⁴	53.9306	5.90
Na ²³	22.9898	100.00	Fe ⁵⁶	55.9349	91.52
Mg ²⁴	23.9350	78.60	Fe ⁵⁷	56.9354	2.25
Mg ²⁵	24.9858	10.11	Fe ⁵⁸	57.9333	0.33
Mg ²⁶	25.9826	11.29	Co ⁵⁹	58.9332	100.00
Al ²⁷	26.9815	100.00	Ni ⁵⁸	57.9354	67.76
Si ²⁸	27.9769	92.18	Ni ⁶⁰	59.9308	26.16
Si ²⁹	28.9765	4.71	Ni ⁶¹	60.9311	1.25
Si ³⁰	29.9738	3.12	Ni ⁶²	61.3283	3.66
P ³¹	30.9738	100.00	Ni ⁶⁴	63.9280	1.16
S ³²	31.9721	95.02	Cu ⁶³	62.9296	69.09
S ³³	32.9715	0.75	Cu ⁶⁵	64.9278	30.91
S ³⁴	33.9679	4.22	Zn ⁶⁴	63.9291	48.89
S ³⁶	35.9671	0.11	Zn ⁶⁶	65.9260	27.81
Cl ³⁵	34.9689	75.40	Zn ⁶⁷	66.9271	4.11
Cl ³⁷	36.9659	24.60	Zn ⁶⁸	67.9249	18.56
A ³⁶	35.9675	0.34	Zn ⁷⁰	69.9253	0.62
A ³⁸	37.9627	0.063	Ga ⁶⁹	68.9257	60.20
A ⁴⁰	39.9624	99.60	Ga ⁷¹	70.9248	39.80
K ³⁹	38.9637	93.08	Ge ⁷⁰	69.9240	20.52
K ⁴⁰	39.9640	0.012	Ge ⁷²	71.9216	27.43
K ⁴¹	40.9618	6.91	Ge ⁷³	72.9233	7.76

Symbol	Amu	Abundance	Symbol	Amu	Abundance
Ge ⁷⁴	73.9210	36.54	Pd ¹⁰⁵	104.9042	22.60
Ge ⁷⁶	75.9213	7.76	Pd ¹⁰⁶	105.9029	27.10
As ⁷⁵	74.9217	100.00	Pd ¹⁰⁸	107.9035	26.70
Se ⁷⁴	73.9225	0.96	Pd ¹¹⁰	109.9045	13.50
Se ⁷⁶	75.9193	9.12	Ag ¹⁰⁷	106.9048	51.35
Se ⁷⁷	76.9200	7.50	Ag ¹⁰⁹	108.9046	48.65
Se ⁷⁸	77.9174	23.61	Cd ¹⁰⁶	105.9059	1.22
Se ⁸⁰	79.9165	49.96	Cd ¹⁰⁸	107.9041	0.89
Se ⁸²	81.9166	8.94	Cd ¹¹⁰	109.9033	12.43
Br ⁷⁹	78.9184	50.57	Cd ¹¹¹	110.9043	12.86
Br ⁸¹	80.9164	49.43	Cd ¹¹²	111.9031	23.79
Kr ⁷⁸	77.9202	0.35	Cd ¹¹³	112.9045	12.34
Kr ⁸⁰	79.9166	2.27	Cd ¹¹⁴	113.9036	28.81
Kr ⁸²	81.9135	11.56	Cd ¹¹⁶	115.9050	7.66
Kr ⁸³	82.9141	11.65	In ¹¹³	112.9043	4.16
Kr ⁸⁴	83.9115	56.90	In ¹¹⁵	114.9036	95.84
Kr ⁸⁶	85.9108	17.37	Sn ¹¹²	111.9051	0.95
Rb ⁸⁵	84.9120	72.15	Sn ¹¹⁴	113.9039	0.65
Rb ⁸⁷	86.9093	27.85	Sn ¹¹⁵	114.9034	0.34
Sr ⁸⁴	83.9133	0.56	Sn ¹¹⁶	115.9022	14.24
Sr ⁸⁶	85.9094	9.86	Sn ¹¹⁷	116.9031	7.67
Sr ⁸⁷	86.9090	7.02	Sn ¹¹⁸	117.9021	24.01
Sr ⁸⁸	87.9060	82.56	Sn ¹¹⁹	118.9032	8.58
Y ⁸⁹	88.9057	100.00	Sn ¹²⁰	119.9022	32.97
Zr ⁹⁰	89.9043	51.46	Sn ¹²²	121.9035	4.71
Zr ⁹¹	90.9053	11.23	Sn ¹²⁴	123.9052	5.98
Zr ⁹²	91.9047	17.11	Sb ¹²¹	120.9037	57.25
Zr ⁹⁴	(98.9076)	17.40	Sb ¹²³	122.9042	42.75
Zr ⁹⁶	95.9093	2.80	Te ¹²⁰	119.9045	0.09
Nb ⁹³	92.9057	100.00	Te ¹²²	121.9029	2.46
Mo ⁹²	91.9086	15.05	Te ¹²³	122.9043	0.87
Mo ⁹⁴	93.9059	9.35	Te ¹²⁴	123.9001	4.61
Mo ⁹⁵	94.0944	14.78	Te ¹²⁵	124.9046	6.83
Mo ⁹⁶	95.9050	16.56	Te ¹²⁶	125.9039	18.71
Mo ⁹⁷	96.9062	9.60	Te ¹²⁸	127.9056	31.79
Mo ⁹⁸	97.9060	24.00	Te ¹³⁰	129.9070	34.49
Mo ¹⁰⁰	99.9065	9.68	I ¹²⁷	126.9047	100.00
Ru ⁹⁶	95.9087	5.68	Xe ¹²⁴	123.9061	0.10
Ru ⁹⁸	(97.9540)	2.22	Xe ¹²⁶	125.9045	0.09
Ru ⁹⁹	98.9061	12.81	Xe ¹²⁸	127.9035	1.92
Ru ¹⁰⁰	(99.9050)	12.70	Xe ¹²⁹	128.9047	26.44
Ru ¹⁰¹	(100.9050)	16.98	Xe ¹³⁰	129.9035	4.08
Ru ¹⁰²	101.9037	31.34	Xe ¹³¹	130.9051	21.18
Ru ¹⁰⁴	103.9042	18.27	Xe ¹³²	131.9042	26.89
Rh ¹⁰³	102.9046	100.00			
Pd ¹⁰²	101.9049	0.80			
Pd ¹⁰⁴	103.9033	9.30			

Symbol	Amu	Abundance	Symbol	Amu	Abundance
Xe ¹³⁴	133.9054	10.44	Dy ¹⁶¹	160.9257	18.88
Xe ¹³⁶	135.9072	8.87	Dy ¹⁶²	161.9256	25.53
Cs ¹³³	132.9051	100.00	Dy ¹⁶³	162.9278	24.87
Ba ¹³⁰	129.9062	0.10	Dy ¹⁶⁴	163.9285	28.18
Ba ¹³²	131.9051	0.10	Ho ¹⁶⁵	164.9291	100.00
Ba ¹³⁴	133.9043	2.42	Er ¹⁶²	161.9287	0.14
Ba ¹³⁵	134.9056	6.59	Er ¹⁶⁴	163.9298	1.56
Ba ¹³⁶	135.9044	7.81	Er ¹⁶⁶	165.9287	33.41
Ba ¹³⁷	136.0956	11.32	Er ¹⁶⁷	166.9305	22.94
Ba ¹³⁸	137.9049	71.66	Er ¹⁶⁸	167.9308	27.07
La ¹³⁸	137.9067	0.09	Er ¹⁷⁰	169.9356	14.88
La ¹³⁹	138.9060	99.91	Tm ¹⁶⁹	168.9343	100.00
Ce ¹³⁶	135.9071	0.19	Yb ¹⁶⁸	167.9339	0.14
Ce ¹³⁸	137.9060	0.25	Yb ¹⁷⁰	(169.9340)	3.03
Ce ¹⁴⁰	139.9053	88.48	Yb ¹⁷¹	(170.9360)	14.31
Ce ¹⁴²	141.9093	11.07	Yb ¹⁷²	(171.9252)	21.82
Pr ¹⁴¹	140.9075	100.00	Yb ¹⁷³	(172.9380)	16.13
Nd ¹⁴²	141.9075	27.09	Yb ¹⁷⁴	(173.9256)	31.84
Nd ¹⁴³	142.9096	12.14	Yb ¹⁷⁶	(175.9427)	12.73
Nd ¹⁴⁴	143.9098	23.83	Lu ¹⁷⁵	(174.9420)	97.40
Nd ¹⁴⁵	144.9121	8.29	Lu ¹⁷⁶	175.9419	2.60
Nd ¹⁴⁶	145.9127	17.26	Hf ¹⁷⁴	(173.9390)	0.20
Nd ¹⁴⁸	147.9165	5.74	Hf ¹⁷⁶	175.9408	5.23
Nd ¹⁵⁰	149.9208	5.63	Hf ¹⁷⁷	176.9425	18.55
Sm ¹⁴⁴	143.9117	3.16	Hf ¹⁷⁸	177.9430	27.23
Sm ¹⁴⁷	146.9145	15.07	Hf ¹⁷⁹	178.9452	13.79
Sm ¹⁴⁸	147.9144	11.27	Hf ¹⁸⁰	179.9461	35.07
Sm ¹⁴⁹	148.9168	13.84	Ta ¹⁸⁰	179.9450	0.01
Sm ¹⁵⁰	149.9169	7.47	Ta ¹⁸¹	180.9458	99.99
Sm ¹⁵²	151.9194	26.63	W ¹⁸⁰	179.9445	0.13
Sm ¹⁵⁴	153.9219	22.53	W ¹⁸²	181.9460	26.31
Eu ¹⁵¹	150.9195	47.77	W ¹⁸³	182.9481	14.28
Eu ¹⁵³	152.9206	52.23	W ¹⁸⁴	183.9489	30.64
Gd ¹⁵²	151.9194	0.20	W ¹⁸⁶	185.9507	28.64
Gd ¹⁵⁴	153.9209	2.15	Re ¹⁸⁵	(184.9480)	37.07
Gd ¹⁵⁵	154.9227	14.73	Re ¹⁸⁷	(186.9530)	62.93
Gd ¹⁵⁶	155.9223	20.47	Os ¹⁸⁴	(183.9460)	0.02
Gd ¹⁵⁷	156.9241	15.68	Os ¹⁸⁶	(185.9520)	1.59
Gd ¹⁵⁸	157.9244	24.87	Os ¹⁸⁷	(186.9560)	1.64
Gd ¹⁶⁰	159.9272	21.90	Os ¹⁸⁸	(187.9569)	13.20
Tb ¹⁵⁹	158.9249	100.00	Os ¹⁸⁹	(188.9585)	16.10
Dy ¹⁵⁶	(155.9140)	0.05	Os ¹⁹⁰	(189.9548)	26.40
Dy ¹⁵⁸	157.9240	0.09	Os ¹⁹²	(191.9622)	41.00
Dy ¹⁶⁰	159.9239	2.29	Ir ¹⁹¹	(190.9634)	38.50

Symbol	Amu	Abundance
Ir ¹⁹³	(192.9667)	61.50
Pt ¹⁹⁰	(189.9640)	0.012
Pt ¹⁹²	(191.9648)	0.78
Pt ¹⁹⁴	(193.9633)	32.80
Pt ¹⁹⁵	(194.9662)	33.70
Pt ¹⁹⁶	(195.9671)	25.40
Pt ¹⁹⁸	197.9646	7.23
Au ¹⁹⁷	196.9672	100.00
Hg ¹⁹⁶	195.9658	0.15
Hg ¹⁹⁸	197.9667	10.02
Hg ¹⁹⁹	198.9682	16.84
Hg ²⁰⁰	199.9683	23.13
Hg ²⁰¹	200.9703	13.22
Hg ²⁰²	201.9706	29.80
Hg ²⁰⁴	203.9735	6.85
Tl ²⁰³	202.9719	29.50
Tl ²⁰⁵	204.9738	70.50
Pb ²⁰⁴	203.9731	1.37
Pb ²⁰⁶	205.9744	25.15
Pb ²⁰⁷	206.9759	21.11
Pb ²⁰⁸	207.9766	52.38
Bi ²⁰⁹	208.9743	100.00
Th ²³²	232.0380	100.00
U ²³⁴	234.0403	0.01
U ²³⁵	235.0428	0.72
U ²³⁸	238.0486	99.28

H																	He
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	Lu	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	Lr	Unq	Unp	Unh	Uns	Uno	Une									
		La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb		
		Ac	Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No		



Hiden Analytical
 420 Europa Boulevard
 Warrington
 WA5 7UN
 England

Tel: +44 (0) 1925 445225
 Fax: +44 (0) 1925 416518
 Toll Free (USA): 1-888 96 HIDEN (44336)
 E-mail: info@hiden.demon.co.uk
 Web Site: <http://www.HidenAnalytical.com>

Appendix B: Atomic weights of common molecules

Name	Formula	Peak 1 m/e %	Peak 2 m/e %	Peak 3 m/e %	rel sens
Acetone	C ₃ H ₆ O	43 100	58 33	15 20	3.6
Air		28 100	32 27	14 6	1.0
Ammonia	NH ₃	17 100	16 80	15 8	1.3
Argon	Ar	40 100	20 16	36 0.3	1.2
Benzene	C ₆ H ₆	78 100	77 19	52 16	5.9
Boron Trichloride	BCl ₃	81 100	58 33	15 20	1.0
Carbon Dioxide	CO ₂	44 100	16 9	14 6	1.4
Carbon Monoxide	CO	28 100	12 5	16 2	1.05
Carbon Tetrafluoride	CCl ₄	69 100	50 12	19 7	1.0
Diborane	B ₂ H ₆	26 100	27 97	24 90	1.0
Ethane	C ₂ H ₆	28 100	27 33	30 26	2.6
Fomblin Oil		69 100	20 28	16 16	1.0
Freon 12	CCl ₂ F ₂	85 100	87 32	50 16	2.7
Helium	He	4 100			0.14
Hydrogen	H ₂	2 100	1 2		0.44

Name	Formula	Peak 1 m/e %	Peak 2 m/e %	Peak 3 m/e %	rel sens
Hydrogen Chloride	HCl	36 100	38 32	35 17	1.6
Hydrogen Sulfide	H ₂ S	34 100	32 44	33 42	2.2
Krypton	Kr	84 100	86 31	82 21	1.7
Methane	CH ₄	16 100	15 85	14 16	1.6
Methanol	CH ₃ OH	31 100	32 67	29 65	1.8
Neon	Ne	20 100	22 10	21 0.3	0.23
Nitrogen	N ₂	28 100	14 5	29 1	1.0
Oxygen	O ₂	32 100	16 9		0.86
Phosphine	PH ₃	34 100	33 33	31 32	2.6
Pump Oil		57 100	55 73	43 73	1.0
Silane	SiH ₄	30 100	31 78	29 29	1.0
Silicon Tetrafluoride	SiF ₄	85 100	86 5	28 4	1.0
Sulfur Dioxide	SO ₂	64 100	48 50	32 10	2.1
Water	H ₂ O	18 100	17 21	16 2	0.9
Xenon	Xe	132 100	129 98	131 79	3.0

Appendix C: Periodic Table

Representative (main group) elements														Representative (main group) elements																								
IA														IIA																VIII A								
1	H 1.0079																											2	He 4.003									
2	Li 6.941													Be 9.012																								
Periodic Table of the Elements																																						
Transition metals																																						
3	Na 22.990	Mg 24.305													Al 26.982	Si 28.086	P 30.974	S 32.066	Cl 35.453	Ar 39.948																		
4	K 39.098	Ca 40.078	Sc 44.956	Ti 47.88	V 50.942	Cr 51.996	Mn 54.938	Fe 55.845	Co 58.933	Ni 58.69	Cu 63.546	Zn 65.39	Ga 69.723	Ge 72.61	As 74.922	Se 78.96	Br 79.904	Kr 83.8																				
5	Rb 85.468	Sr 87.62	Y 88.906	Zr 91.224	Nb 92.906	Mo 95.94	Tc 98	Ru 101.07	Rh 102.906	Pd 106.42	Ag 107.868	Cd 112.411	In 114.82	Sn 118.71	Sb 121.76	Te 127.60	I 126.905	Xe 131.29																				
6	Cs 132.905	Ba 137.327	La 138.906	Hf 178.49	Ta 180.948	W 183.84	Re 186.207	Os 190.23	Ir 192.22	Pt 195.08	Au 196.967	Hg 200.59	Tl 204.383	Pb 207.2	Bi 208.980	Po 209	At 210	Rn 222																				
7	Fr 223	Ra 226.025	Ac 227.028	Rf 261	Db 262	Sg 263	Bh 262	Hs 265	Mt 266	Uun 269	Uuu 272	Uub 277																										
Rare earth elements																																						
Lanthanides														Actinides																								
58 Ce 140.115 59 Pr 140.908 60 Nd 144.24 61 Pm 145 62 Sm 150.36 63 Eu 151.964 64 Gd 157.25 65 Tb 158.925 66 Dy 162.5 67 Ho 164.93 68 Er 167.26 69 Tm 168.934 70 Yb 173.04 71 Lu 174.967														90 Th 232.038 91 Pa 231.036 92 U 238.029 93 Np 237.048 94 Pu 244 95 Am 243 96 Cm 247 97 Bk 247 98 Cf 251 99 Es 252 100 Fm 257 101 Md 258 102 No 259 103 Lr 262																								

Copyright © 2000 Benjamin/Cummings, an imprint of Addison Wesley Longman, Inc.