

Citation: Dhariwal J, Kitson J, Jones RE, Nicholson G, Tunstall T, Walton RP, et al. (2015) Nasal Lipopolysaccharide Challenge and Cytokine Measurement Reflects Innate Mucosal Immune Responsiveness. PLoS ONE 10(9): e0135363. doi:10.1371/journal.pone.0135363

Editor: T. Mark Doherty, Glaxo Smith Kline, DENMARK

Received: November 14, 2014

Accepted: July 18, 2015

Published: September 14, 2015

Copyright: © 2015 Dhariwal et al. This is an open access article distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Data Availability Statement: All relevant data are within the paper and its Supporting Information files.

Funding: This study was funded by Glaxo SmithKline and supported by National Institute for Health Research (NIHR) Biomedical Research Centre (BRC) based at Imperial College NHS Trust and Imperial College, London, United Kingdom. Glaxo SmithKline provided support in the form of salaries for authors JK, REJ, GF, JG and RM, and played a full and active role in the study design, data collection and analysis, decision to publish, or preparation of the manuscript. The specific roles of **RESEARCH ARTICLE**

Nasal Lipopolysaccharide Challenge and Cytokine Measurement Reflects Innate Mucosal Immune Responsiveness

Jaideep Dhariwal¹, Jeremy Kitson², Reema E. Jones², Grant Nicholson¹, Tanushree Tunstall¹, Ross P. Walton¹, Grace Francombe³, Jane Gilbert³, Andrew J. Tan¹, Robert Murdoch², Onn Min Kon¹, Peter J. Openshaw¹, Trevor T. Hansel¹*

 Centre for Respiratory Infection (CRI), National Heart and Lung Institute (NHLI) at Imperial College, St. Mary's Hospital: Mint Wing, Entrance C, Paddington, London, W2 1NY, United Kingdom,
GlaxoSmithKline, Stevenage, SG1 2NY, United Kingdom, 3 GlaxoSmithKline, Stockley Park West, Uxbridge, UB11 1BT, United Kingdom

* t.hansel@imperial.ac.uk

Abstract

Background

Practical methods of monitoring innate immune mucosal responsiveness are lacking. Lipopolysaccharide (LPS) is a component of the cell wall of Gram negative bacteria and a potent activator of Toll-like receptor (TLR)-4. To measure LPS responsiveness of the nasal mucosa, we administered LPS as a nasal spray and quantified chemokine and cytokine levels in mucosal lining fluid (MLF).

Methods

We performed a 5-way cross-over, single blind, placebo-controlled study in 15 healthy nonatopic subjects (n = 14 *per protocol*). Doses of ultrapure LPS (1, 10, 30 or 100µg/100µl) or placebo were administered by a single nasal spray to each nostril. Using the recently developed method of nasosorption with synthetic adsorptive matrices (SAM), a series of samples were taken. A panel of seven cytokines/chemokines were measured by multiplex immunoassay in MLF. mRNA for intercellular cell adhesion molecule-1 (ICAM-1) was quantified from nasal epithelial curettage samples taken before and after challenge.

Results

Topical nasal LPS was well tolerated, causing no symptoms and no visible changes to the nasal mucosa. LPS induced dose-related increases in MLF levels of IL-1 β , IL-6, CXCL8 (IL-8) and CCL3 (MIP-1 α) (AUC at 0.5 to 10h, compared to placebo, p<0.05 at 30 and 100 μ g LPS). At 100 μ g LPS, IL-10, IFN- α and TNF- α were also increased (p<0.05). Dose-related changes in mucosal ICAM-1 mRNA were also seen after challenge, and neutrophils appeared to peak in MLF at 8h. However, 2 subjects with high baseline cytokine levels showed prominent cytokine and chemokine responses to relatively low LPS doses (10 μ g and 30 μ g LPS).

these authors are articulated in the author contributions section.

PLOS ONE

Competing Interests: This study was partly funded by Glaxo SmithKline, the employer of JK, REJ, GF, JG and RM. In the past 5 years TTH has been Principal Investigator for respiratory clinical studies carried out at Imperial College, that have been sponsored by Pharmaceutical Companies: Novartis, Merck, GlaxoSmithKline (GSK) and Dainippon-Sumitomo. TTH has received lecture fees from Novartis, Merck, GSK and Dainippon-Sumitomo. TTH is involved in setting up a spin-off company with Imperial Innovations called Mucosal Diagnostics, and has received support to attend Respiratory Conferences from Boehringer Ingelheim. There are no patents, products in development or marketed products to declare. This does not alter the authors' adherence to all the PLOS ONE policies on sharing data and materials, as detailed online in the guide for authors.

Conclusions

Topical nasal LPS causes dose-dependent increases in cytokines, chemokines, mRNA and cells. However, responsiveness can show unpredictable variations, possibly because baseline innate tone is affected by environmental factors. We believe that this new technique will have wide application in the study of the innate immune responses of the respiratory mucosa.

Key Messages

Ultrapure LPS was used as innate immune stimulus in a human nasal challenge model, with serial sampling of nasal mucosal lining fluid (MLF) by nasosorption using a synthetic absorptive matrix (SAM), and nasal curettage of mucosal cells. A dose response could be demonstrated in terms of levels of IL-1 β , IL-6, CXCL8 and CCL3 in MLF, as well as ICAM-1 mRNA in nasal curettage specimens, and levels of neutrophils in nasal lavage. Depending on higher baseline levels of inflammation, there were occasional magnified innate inflammatory responses to LPS.

Trial Registration

Clinical Trials.gov NCT02284074

Introduction

A range of inhaled and nasal challenges have been fundamental to our understanding of the pathophysiology of allergy and asthma [1] [2]. Challenge tests give a controlled stimulus after which serial measurements can be made, and the time course of responses can be assessed. Experimental challenge studies with components of microbes (including pathogen associated molecular patterns, PAMPs) and live viruses [3] give researchers the opportunity to study mucosal host-microbial interactions in terms of the specific and innate immune systems [4,5].

The respiratory mucosa of the nose has an intricate sentinel role interacting with bacteria, viruses and aeroallergens [6]. This is a complex and dynamic interaction that is mediated at a molecular level by PAMPs interacting with host pattern recognition receptors (PRRs) [7,8]. Lipopolysaccharide (LPS) is contained in the wall of Gram negative bacteria and in a range of microbes, and is a powerful PAMP with potent immunostimulatory properties through interactions with Toll-like receptor (TLR)4. LPS represents the prototypic PAMP involved in certain bacterial infections [7], atopy [9], and in acute lung injury [10]. Interestingly, TLR4 is involved in some viral infections [11]. House dust mite allergen has been shown to mimic an LPS-binding component of the TLR4 signalling complex and causes activation of airway structural cells [12,13]. Environmental exposure to LPS in childhood may facilitate development of tolerance to allergens [14], and this is one of a range of environmental microbial influences for the development of allergy and asthma [15].

LPS has been employed extensively as a challenge agent by David Peden and coworkers: with systemic administration by intravenous infusion [16], by inhalation from a nebulizer [17], by segmental bronchial challenge [18] and by nasal challenge [19–22]. In particular, the group of Lars-Olaf Cardell have carried out a series of studies with nasal LPS challenge: measuring CCL3 (MIP1 α) in nasal lavage [23], evaluating a CXCR2 antagonist [24], demonstrating that

allergen primes LPS responses [25], and showing that LPS challenge has augmented nasal cytokine release in allergic rhinitis [26].

Recently, synthetic absorptive matrices (SAM) have been used to sample mucosal lining fluid (MLF) in the nose, followed by measurement of mediators. Nasosorption has been employed in neonates [27], children with rhinitis [28] and in relation to nasal allergen challenge (NAC) [29]. The method of nasosorption has advantages in being non-invasive and being able to sample neat MLF, which is eluted in a known volume of buffer [30]. Because there is less dilution, nasosorption allows detection of higher levels of chemokines and cyto-kines than after nasal lavage. An additional non-invasive sampling method employs nasal epithelial curettage that has been employed to assess gene expression after nasal allergen challenge [31] and experimental rhinovirus infection [32].

In this study we carry out topical nasal ultra-pure LPS challenge in healthy non-atopic volunteers, and investigate the inflammatory response following a series of different doses of LPS. We employed serial nasosorption sampling and measured levels of soluble chemokines and cytokines using a multiplex immunoassay system. In addition, nasal epithelial curettage samples were taken for analysis of mRNA (transcriptomics) for intercellular cell adhesion molecule-1 (ICAM-1) that is elevated in the nasal mucosa in allergic rhinitis [33].

Methods

Ethics and Consent

This study was approved by the South West London Research Ethics Committee. Since this study was a topical challenge study for research purposes, and not for therapeutic benefit, the study did not require authorization from the UK Medicines and Healthcare Products Regulatory Agency (MHRA). The study was carried out on the Imperial Clinical Respiratory Research Unit (ICRRU) of St Mary's Hospital in accordance with the Declaration of Helsinki and Good Clinical Practice (GCP) guidelines.

Study Design

All healthy volunteers were recruited through newspaper advertising (Fig 1). The dose response to LPS was assessed by a 5-way crossover, randomized, single-blind, placebo-con-trolled study design (Fig 2). The dose randomization for the 5 arms (placebo, 1, 10, 30, 100µg LPS/nostril) was performed by an independent statistician, employing a balanced Latin square design to generate the order of challenges for each subject. The subjects were "blinded" to the dosage being employed. Subjects were assigned to a sequence in accordance with a randomisation schedule generated prior to the start of the study. The randomisation was balanced over the 5 nasal challenge sequences and was generated using RandAll, the GSK web-server based clinical trials randomisation system.

Inclusion Criteria at Screening

Subjects were aged 18 to 60 years. They were current non-smokers for at least 6 months with a smoking history of <5 pack years. These individuals had no history of atopy, rhinitis, asthma, or other respiratory disease. Individuals with a history of malignancy or systemic disease were excluded, and use of any inhaled or oral anti-inflammatory medication was not permitted.

LPS Administration

LPS solutions at concentrations of 10, 100, 300 and 1000µg/ml were prepared using ultra-pure LPS from *Escherichia coli* 0111:B4 (Invivogen, San Diego, CA). This is a high quality source of



Fig 1. Consort Flow Chart.

doi:10.1371/journal.pone.0135363.g001

LPS, that is specific for TLR4 and does not contain bacterial lipoprotein contaminants that activate TLR2 (manufacturer's specification). LPS potency was confirmed within our laboratory using the E-Toxate (*Limulus* Amoebocyte Lysate) kit (Sigma-Aldrich, Gillingham, UK). A standard curve and dilutions of the test LPS were prepared and mixed with the lysate. The ultra-pure LPS used in the clinical study contained 166–333 EU/ μ g, with the unit activity at the higher end of the manufacturer's specification on the *Limulus* assay.



Fig 2. Nasal lipopolysaccharide challenge study design and procedures. 14 subjects completed the entire protocol (*per protocol*), while 1 subject only received nasal challenges at placebo and LPS. TNSS-total nasal symptom score. AE-adverse events

doi:10.1371/journal.pone.0135363.g002

LPS and placebo were administered from a Pfeiffer Bidose nasal delivery device (Thermofisher, Epsom, UK) by a physician. A volume of 100µl was administered as a spray to each nostril, with a total dose per nostril of 1, 10, 30 or 100µg of LPS. Weights were recorded before and after each actuation from the Pfeiffer Bidose, to check that a 100µg weight of fluid had been delivered per nostril. A washout time of 12–28 days was employed between challenges.

Symptom Scores

Modified total nasal symptom scores (TNSS) were recorded at the initial screening visit and before and at frequent intervals after nasal LPS challenge (-15 minutes, 30 minutes, then hourly from 1 to 10 hours and then at 24h post-LPS). Nasal congestion, rhinorrhea, sneezing, and nasal itch were scored from 0 to 3 (0: no, 1: mild, 2: moderate, 3: severe symptoms). The scores were then summed to give a final TNSS out of a maximum of 12.

We also measured safety parameters (temperature, heart rate, blood pressure and oxygen saturations), performed regular nasal mucosal inspections, and assessed peak nasal inspiratory flow (PNIF)

Nasal Lavage

A modified nasal pool technique was adapted from the method of Greiff et al [34]. The subjects were seated in a forward-flexed neck position (60 degrees from the upright position) to prevent fluid from reaching the nasopharynx. To ensure adequate washing, the lavage fluid (5ml of 0.9% normal saline) was passed via a 10ml syringe slowly into the nasal cavity via an olive consisting of an oval, hollow, stainless steel device that was used to obstruct the nostril. The fluid was withdrawn into the syringe and gently flushed back into the nasal cavity 20 times over 1 minute.

Nasal lavage samples were initially centrifuged (4°C, 10 min at 400g), and the cell pellet resuspended in 0.01% dithiothreitol in phosphate buffered saline (PBS). The samples were gently agitated on a rolling mixer for 10 minutes and centrifuged; the supernatant was

discarded; the pellet was resuspended in PBS at 1 to 2 million cells/mL, and cytospin slides were prepared. Differential leukocyte counts were determined by assessment of 400 leukocytes on the stained cytospin.

Nasosorption

Nasosorption was performed with strips of synthetic absorptive matrix (SAM; Fibrous Polyester; Hunt Developments, Midhurst, West Sussex, UK). The SAM was inserted under direct vision into the nasal cavity, the length being applied laterally against the anterior inferior turbinate. Nose clips were used to ensure good contact with the mucosal surface. Pre-chilled assay buffer (PBS, bovine serum albumin (BSA) 1%; Tween-20 0.05%, sodium azide 0.08% (Milliplex Assay Buffer, cat no. L-AB, Merck Millipore, Billerica, Mass., USA), was dispensed in 500µL aliquots into filter cups within Eppendorf tubes (Costar spin-X, cellulose acetate). After removal of the SAM strips from the nose, they were placed in the assay buffer, and spin filtration was performed (4°C for 5 min at 16,000 g). The eluate was collected and stored as aliquots at -80°C, prior to assay of chemokines and cytokines by a multiplex immunoassay (Mesoscale Diagnostics, MSD, Gaithersburg, Md, USA).

Cytokine levels in Nasosorption Samples

Two custom 7-spot MSD plates (MSD) were used for measuring the concentrations of CCL3 (MIP-1α), IL-1β, IL-6, CXCL8 (IL-8), TNF-α, IFN-α, IL-10 (Plate 1) and IFN-γ, GM-CSF, CCL2 (MCP-1), CXCL10 (IP-10), IL-12p70, IL-17 (Plate 2). The plates provided are pre-coated with antibody for the specific analytes. Diluent 2 (MSD) (25µl) was added to each well and the plate was incubated for 30 min while being vigorously shaken. All incubations in the MSD immunoassay were done on an orbital shaker at 1000 rpm at room temperature. 25μ L of the sample or calibrator was added to each well. A standard curve of calibrators (10000, 2500, 625, 156.3, 39.1, 9.8 and 2.4pg/mL) was run in duplicate. Controls of known concentrations (625pg/mL and 9.8pg/mL) were done in quadruplicate. Three replicates of each sample were added. After an incubation of 2 hours on a shaker with the sample/calibrator the plates were washed three times with PBS + 0.05% Tween and then 25μ L of 1X Detection Antibody Solution was added to each well. After incubation for 2 hours with vigorous shaking the plate was washed 3 times with PBS + 0.05% Tween. 150μ L of 2X Read Buffer was added to each well and the plate was read immediately using a MSD Sector Imager 6000. The raw data was analysed using MSD Discovery Workbench 3.0 software. The concentration of each cytokine was calculated from the standard curve of the Calibrators.

Nasal Epithelial Curettage

A single nasal curettage was performed at baseline and again at 3, 6, and 24hrs to obtain epithelial cells for mRNA gene expression studies using a plastic nasal curette (Rhinoprobe, Arlington Scientific). All nasal scrapes were performed after the collection of nasosorption and nasal lavage samples. Nasal scrapes were taken from the posterior aspect of the inferior turbinate, without local anaesthesia, collecting a single 5mm mucosal sample. Nasal scrapes were placed in labelled sterile tubes containing TRIzol reagent (500µl) and stored at -80°C.

Analysis of ICAM-1 mRNA expression from nasal Rhinoprobe samples by RT-qPCR

Nasal scrapes samples were placed in TRIzol Reagent (Invitrogen, Carlsbad, CA, USA) and stored at -80°C until the total RNA was purified. Total RNA was extracted using RNeasy Mini-

Statistical analysis

Biomarker data was available for 15 subjects although subject number 3102 withdrew after 2 periods. For all biomarkers, the lower limit of quantitation (LLOQ) was taken as 9.8pg/ml and a value of LLOQ/2 (= 4.9pg/ml) was used for all values less than the LLOQ. The median was taken of the 3 adjusted readings and used in the analysis. The area under curve (AUC) from 0-10h of the observed data were calculated for each dose of LPS, and responses after LPS challenge described as changes relative to placebo AUC 0-10h. The distribution of the AUC data was tested using a Shapiro-Wilks test and significant deviations from normal data distribution were seen. Therefore non-parametric statistics (Wilcoxon signed rank test) were applied throughout.

Statistical analysis of ICAM-1 mRNA data was performed in SAS v9.2 using a mixed model analysis of variance. All subjects received each of the treatments at a separate visit and within each visit a response across time was obtained. Due to the complex design of the study 'time' was included in the model as a repeated measures and 'subject/ visit' was included as a random effect to allow the appropriate error structure to be used for assessing LPS effects at each time point. The response used was logged gene abundance and log GAPDH (glyceraldehyde-3-phosphate dehydrogenase) was included as a covariate so that the resulting estimates of the mean fold changes from placebo were normalised.

Results

Baseline Demographics

Subjects were recruited through newspaper adverts and all visits took place between November and April (Fig 2). Following initial entry of 17 subjects, a total of 14 subjects completed all study visits and sampling procedures (*per protocol* group). A single individual received 2 doses (placebo and LPS at 100µg) and was included in the analysis, so that 15 subjects were studied. Subjects ranged in age from 21 to 57 years (mean 38 years), and there were 10 males and 5 females.

One individual withdrew from the study following development of an infective pyelonephritis 10 days after a low dose of LPS (1 μ g), this being the first nasal challenge, and the event was considered to be "not clinically related to LPS administration". A second individual was withdrawn from the study after 2 of the 5 planned nasal challenges, as they had withheld medical information that excluded their participation. These 2 subjects were excluded from the analysis.

Clinical Symptom Scores

None of the participants experienced any significant clinical symptoms during the course of their study participation. Following nasal LPS administration there were no visible signs of inflammation on direct examination of the nasal mucosa by a physician.

Levels of Cytokines and Chemokines in Nasal MLF

Following nasal LPS challenge there was a dose response for IL-1 β , IL-6, CXCL8, CCL3 and IL-10 levels in nasal SAM eluate. Median levels of these cytokines and chemokines over time are shown in Fig 3. Although there was an increase in the response to nasal LPS at doses up to 30µg, the response at 30 and 100µg was approximately equivalent.

The Area Under the Curve (AUC) for individual biomarkers measured at 0.5 to 10h was compared with placebo. Following LPS there were significant increases in AUC for IL-1 β , IL-6, CXCL8, and CCL3 (<u>Table 1</u>). IL-10, IFN- α and TNF- α were increased significantly at 100 μ g alone. There were particularly marked IL-1 β and IL-6 responses and AUC data is shown in response to nasal LPS dose as medians with quartiles, and individual AUC responses in Fig 4.

Idiosyncratic Responses (2)

Examining the levels of cytokines and chemokines in nasal SAM eluates from individual subjects, it was noted that 2 subjects had unexpectedly large responses to LPS on a single occasion during dose-escalation. This is illustrated in <u>Fig 4</u> (lower portion) that shows levels of cytokines and chemokines as AUC from 0.5h to 10.0h. These markedly elevated responses occurred in





Fig 3. Median levels of chemokines and cytokines (pg/ml) in nasal mucosal lining fluid (MLF) are shown after nasal spray challenge with 4 different doses of LPS and a placebo (see key below). Lower limits of detection (LLOD) for cytokines and chemokines were as follows: IL-1 β (5.0pg/ml), IL-6 (5.0pg/ml), CXCL8/IL-8 (10.0pg/ml), CCL3/MIP-1 α (5.0pg/ml), CCL2/MCP-1 (25.0pg/ml), IL-10 (5.0pg/ml). Time -0.5h refers to nasosorption performed 30 min prior to topical nasal challenge, and this varies considerably between individuals due to variations in individual's microbial flora and mucosal immune responses. Nasal lavage was performed after nasosorption at -0.5h, but prior to nasal challenge, in order to partially wash the nose free of baseline inflammatory microbes and mediators, but causes detectable levels of inflammatory mediators at 0.5h to decrease markedly from -0.5h. There is a tendency for nasosorption levels of mediators to gradually increase from 0.5h to 10h in the placebo arm, but levels of inflammatory cytokines and chemokines after LPS are generally higher.

doi:10.1371/journal.pone.0135363.g003



| Cytokine or chemokine | Nasal challenge LPS, 1µg | Nasal challenge LPS, 10µg | Nasal challenge LPS, 30µg | Nasal challenge LPS, 100µg |
|-----------------------|--------------------------|---------------------------|---------------------------|----------------------------|
| IL-1β | ns | <0.05 | <0.05 | <0.01 |
| IL-6 | ns | ns | <0.05 | <0.05 |
| CXCL8 (IL-8) | ns | ns | <0.05 | <0.01 |
| CCL3 (MIP-1α) | <0.05 | ns | <0.05 | <0.001 |
| CCL2 (MCP-1) | ns | <0.05 | ns | ns |
| IL-10 | ns | ns | ns | <0.05 |
| IFN-α | ns | ns | ns | <0.05 |
| TNF-α | ns | ns | ns | <0.05 |

Table 1. Statistical analysis of levels of cytokines and chemokines in nasal mucosal lining fluid¹.

¹Footnote: Area Under the Curve (AUC) for individual cytokines or chemokines measured at 0.5 to 10h was compared with placebo. The non-parametric Wilcoxon signed rank test was used for statistical analysis. Probability (p) values are presented, ns = non significant. The levels of the following cytokines and chemokines were not significantly altered by LPS: IFN-γ, GM-CSF, CXCL10 (IP-10), IL-12p70, IL-17.

doi:10.1371/journal.pone.0135363.t001

subject 2009 following 10µg LPS, and in subject 2101 following 30µg LPS. With regards to the order of challenges, subject 2009 received 10µg as their second dose (100µg, 10µg, 1µg, placebo, 30µg), while subject 3101 received 30µg LPS as their first dose (30µg, 100µg, 1µg, placebo, 10µg). This suggests that sensitisation to LPS is not the explanation for these large idiosyncratic responses.

Following the 2 idiosyncratic responses to nasal LPS, the raw data for nasal SAM eluate levels of cytokines and chemokines is presented in Fig 5. In subject 2009 the response was mainly in terms of IL-6 and IL-1 β , while in subject 3101 CCL2 (MCP-1) was also elevated. It was note-worthy that on the occasion of both these idiosyncratic responses, there were high baseline levels of nasal MLF IL-1 β and IL-6.

Levels of Nasal Lavage Neutrophils

There was an increase in the number of neutrophils in nasal lavage from baseline at 4h with doses of LPS at 1, 30 and $100\mu g$ (Fig 6). However, changes did not reach significance, and there was not dose-related. This is in keeping with our experience that nasal lavage is a somewhat unreliable sampling method, and that subsequent cell counts can be variable.

ICAM-I mRNA

ICAM mRNA expression was measured in nasal curettage samples by quantitative PCR (qPCR). Samples were taken 30 minutes before LPS challenge and at 3, 6 and 24 hours post challenge. A dose dependent increase in ICAM-1 mRNA expression was observed following LPS treatment (Fig 7). The maximal increase in ICAM-1 mRNA over placebo was observed at 6 hours post LPS challenge: 1µg LPS (p<0.05), 10µg LPS (p<0.01), 30µg LPS (p<0.01) and 100µg LPS (p<0.0001). A significant increase over placebo was also observed at 3 hours for the 30µg and 100µg groups (both at p<0.05). ICAM-1 mRNA expression levels returned to baseline at 24h post LPS treatment.

Discussion

This 5-way crossover study documents the nasal mucosal immune response to LPS challenge by performing serial non-invasive sampling (nasosorption and nasal epithelial curettage) and then measuring levels of cytokines and chemokines together with mRNA for ICAM-1. Original features of this report include the use of ultrapure LPS, serial nasal sampling by the novel 

Fig 4. Levels of interleukins (IL) in nasosorption eluates presented as area under the curve (AUC) for times 0.5-10h in relation to topical nasal challenge (n = 14 or 15). AUC is derived from from 0.5 since the baseline sample at -0.5h is pre-nasal lavage and is highly variable due to difference in individuals in terms of nasal microbes and inflammatory mediators. A) IL-1 β as medians with quartiles; B) IL-6 as medians with quartiles; C) IL-1 β AUC data for individual subjects; D) IL- AUC data for individual subjects. Note the idiosyncratic large responses of subject 2009 after nasal challenge with 10µg LPS, and subject 3101 after nasal challenge with 30µg LPS, these responses being too large to be presented on the y axis.

doi:10.1371/journal.pone.0135363.g004



Fig 5. Idiosyncratic large nasal responses to lipopolysaccharide (LPS) on 2 occasions. Raw data for levels of chemokines and cytokines (pg/ml) in nasal mucosal lining fluid (MLF). Data is presented for individual subjects: A. Subject 2009 after nasal challenge with 10µg LPS. B. Subject 3101 after nasal challenge with 30µg LPS. Note that baseline levels of cytokines and chemokines at -0.5h are markedly raised, and come down to only a small extent after nasal lavage at 0.5h. It is suggested that these individuals may have had a non-symptomatic alteration in their nasal microbial flora prior to LPS challenge on these opccasions: possibly the individuals had a subclinical viral infection or bacterial colonisation.

doi:10.1371/journal.pone.0135363.g005



Fig 6. Nasal lavage leukocyte counts (medians expressed as cells/ml) following nasal LPS challenge. Cells were counted in a modified Neubauer chamber, and a leukocyte differential performed on a stained cytospin observed under a light microscope.

doi:10.1371/journal.pone.0135363.g006



Fig 7. ICAM-1 mRNA abundance within nasal curettage samples after nasal LPS challenge, Data is shown as means with 95% confidence intervals. Relative abundance of expression with Reverse Transcriptase quantitative Polymerase Chain Reaction (RT-qPCR) using normalized data and (glyceraldehyde-3-phosphate dehydrogenase) (GAPDH) as a housekeeper gene (see <u>Methods</u>).

doi:10.1371/journal.pone.0135363.g007

PLOS ONE

method of nasosorption, employing LPS at a range of 4 doses, measuring a range of cytokines and chemokines using a multiplex immunoassay, and assessment of mRNA in epithelial curettage specimens. We were able to observe a dose response to LPS up to 30μ g in terms of significantly increased levels of IL-1 β , IL-6, CXCL8 and CCL3, which suggests that the dose of 50μ g per nostril employed by Cardell and co-workers is optimal in non-atopic subjects [23] [26,35]. Other nasal LPS challenge studies have demonstrated increases in IL-6 [20,36], CXCL8 [21,35] and CCL3 [23].

In the situation of a nasal topical challenge and wishing to take a series of repeated samples from the nasal mucosa, nasosorption of mucosal lining fluid has advantages over nasal lavage. In a study if repeated sampling in healthy volunteers, it was found that nasal lavage levels of a range of mediators were approximately 10-times lower than those from nasosorption eluates [<u>37</u>]. This was confirmed in a study of levels of IL-4, IL-5 and IL-13 after nasal allergen challenge (NAC), where levels of these interleukins were higher in nasosorption eluates, and frequently undetectable in nasal lavages [<u>38</u>]. In a recent study it was found that repeating nasal

lavage within 0.5 to 4h caused markedly decreased levels of mediators, and they concluded that nasal lavage can only be performed once daily to get comparable results [39].

Using a modified nasal lavage collection method, a recent nasal LPS challenge study failed to show significant cytokine responses in allergic rhinitis patients, but amplified responses occurred when nasal allergen challenge preceded nasal LPS challenge [26]. The reason our study was able to document a dose response was probably because we utilized nasosorption to obtain nasal MLF, while other studies have generally measured more dilute nasal secretions present in nasal lavage.

An interesting finding in the study was that in 2 subjects there were large idiosyncratic responses to submaximal doses of LPS. Whereas the median response following higher doses of LPS may be up to 5 times greater than that after placebo, the idiosyncratic reactions involved IL-6 and IL-1 β responses of greater than 1000 times those of placebo. Despite these large inflammatory reactions, subjects reported no symptoms, and there were no visible signs of nasal inflammation. When we assessed the sequence of LPS administrations, there was no evidence of a sensitization or tolerization process occurring with earlier doses. The phenomenon of tolerance to endotoxin has been demonstrated when low dose endotoxin is administered nasally up to 3 times on consecutive days [19,22],

The idiosyncratic exaggerated reactions occurred in one subject after 10µg LPS and in another subject at 30µg LPS. The upper respiratory tract is now recognised to contain a vast array of commensal and pathogenic bacteria included within the nasal microbiome, and these may interact with aeroallergens [40] [9]. In asymptomatic healthy individuals there is frequently carriage of the potential pathogens *Staphylococcus aureus* and *Streptococcus pneumoniae*. There is then extensive crosstalk among PAMPs, TLRs and other innate PRRs in the mucosa [7]. Interestingly the mucosal relationship with bacteria and viruses is altered in the asthmatic or allergic airway. In particular, TLR4 signalling in stromal cells has been shown in mice to be critical for initiation of allergic Th2 responses [41].

We are not currently able to explain these idiosyncratic responses, but each case there were high levels of IL-1 β and IL-6 in nasal mucosal lining fluid before challenge. This leads us to speculate that the idiosyncratic reactions relate to pre-activation by variations in nasal exposure to aeroallergens, viruses or bacteria.

In addition LPS caused a dose response in terms of induction of gene expression for ICAM-1, that is upregulated following inflammation, and acts as an epithelial receptor for human rhinoviruses (HRV). ICAM-1 is upregulated in the nasal mucosa of allergic rhinitis patients [33], and HRV causes activation of epithelial cells through this receptor. Measurement of soluble-ICAM-1 in nasal MLF would confirm whether the increased ICAM-1 gene expression relates to increased production of protein.

In conclusion, this study demonstrates a dose response to nasal LPS challenge in terms of levels of cytokines and chemokines measured in serial nasal MLF samples, and also with raised levels of ICAM-1 mRNA being found in serial nasal epithelial curettage samples. However, we also highlight the occasional phenomenon of an idiosyncratic large response to LPS, which is not dose-related, and that occurred in 2 subjects on single occasions.

In future studies this method, we intent do examine mucosal immune responses to a series of different challenge agents ranging from allergens, diverse microbial components and live viruses. We will also study the longitudinal responses to LPS and other agents over time and to document the influences of one nasal challenge agent on another (as has been done for nasal allergen challenge followed by LPS challenge [26]). Finally, in assessing the response to these nasal challenges in patients with allergic rhinitis and asthma, it may be possible to characterize mucosal immunology in individual patients as a basis for rational therapy.

Supporting Information

S1 CONSORT Checklist. CONSORT Checklist. (DOC)

S1 File. Data Files. (XLS)S1 Protocol. Study Protocol. (DOC)

Author Contributions

Conceived and designed the experiments: TTH RM. Performed the experiments: JD AT GN JK RJ GF TT. Analyzed the data: JD JK JG RPW OMK PJO RM TTH GF TT. Wrote the paper: JD TTH. Clinical procedures: JD AT. Laboratory work: JD GN RJ JK. Review of Manuscript: JD AT GN JK RJ RPW OMK PJO RM TTH.

References

- 1. O'Byrne PM, Gauvreau GM, Brannan JD (2009) Provoked models of asthma: what have we learnt? Clin Exp Allergy 39: 181–192. CEA3172 [pii];doi: 10.1111/j.1365-2222.2008.03172.x PMID: 19187330
- Holgate ST (2012) Innate and adaptive immune responses in asthma. Nat Med 18: 673–683. nm.2731 [pii];doi: 10.1038/nm.2731 PMID: 22561831
- Pollard AJ, Savulescu J, Oxford J, Hill AV, Levine MM, Lewis DJ, Read RC, Graham DY, Sun W, Openshaw P, Gordon SB (2012) Human microbial challenge: the ultimate animal model. Lancet Infect Dis 12: 903–905. S1473-3099(12)70292-X [pii];doi: <u>10.1016/S1473-3099(12)70292-X</u> PMID: <u>23174372</u>
- Yan N, Chen ZJ (2012) Intrinsic antiviral immunity. Nat Immunol 13: 214–222. ni.2229 [pii];doi: <u>10.</u> <u>1038/ni.2229</u> PMID: <u>22344284</u>
- Holt PG, Sly PD (2012) Viral infections and atopy in asthma pathogenesis: new rationales for asthma prevention and treatment. Nat Med 18: 726–735. nm.2768 [pii];doi: 10.1038/nm.2768 PMID: 22561836
- Hansel TT, Johnston SL, Openshaw PJ (2013) Microbes and mucosal immune responses in asthma. Lancet. S0140-6736(12)62202-8 [pii];doi: <u>10.1016/S0140-6736(12)62202-8</u>
- Kawai T, Akira S (2011) Toll-like receptors and their crosstalk with other innate receptors in infection and immunity. Immunity 34: 637–650. S1074-7613(11)00190-7 [pii];doi: <u>10.1016/j.immuni.2011.05</u>. <u>006</u> PMID: <u>21616434</u>
- Netea MG, Wijmenga C, O'Neill LA (2012) Genetic variation in Toll-like receptors and disease susceptibility. Nat Immunol 13: 535–542. ni.2284 [pii];doi: <u>10.1038/ni.2284</u> PMID: <u>22610250</u>
- Simpson A, Martinez FD (2010) The role of lipopolysaccharide in the development of atopy in humans. Clin Exp Allergy 40: 209–223. CEA3391 [pii];doi: <u>10.1111/j.1365-2222.2009.03391.x</u> PMID: <u>19968655</u>
- Imai Y, Kuba K, Neely GG, Yaghubian-Malhami R, Perkmann T, van LG, Ermolaeva M, Veldhuizen R, Leung YH, Wang H, Liu H, Sun Y, Pasparakis M, Kopf M, Mech C, Bavari S, Peiris JS, Slutsky AS, Akira S, Hultqvist M, Holmdahl R, Nicholls J, Jiang C, Binder CJ, Penninger JM (2008) Identification of oxidative stress and Toll-like receptor 4 signaling as a key pathway of acute lung injury. Cell 133: 235– 249. S0092-8674(08)00340-1 [pii];doi: 10.1016/j.cell.2008.02.043 PMID: 18423196
- Shirey KA, Lai W, Scott AJ, Lipsky M, Mistry P, Pletneva LM, Karp CL, McAlees J, Gioannini TL, Weiss J, Chen WH, Ernst RK, Rossignol DP, Gusovsky F, Blanco JC, Vogel SN (2013) The TLR4 antagonist Eritoran protects mice from lethal influenza infection. Nature. nature12118 [pii];doi: <u>10.1038/</u> <u>nature12118</u>
- Trompette A, Divanovic S, Visintin A, Blanchard C, Hegde RS, Madan R, Thorne PS, Wills-Karp M, Gioannini TL, Weiss JP, Karp CL (2009) Allergenicity resulting from functional mimicry of a Toll-like receptor complex protein. Nature 457: 585–588. nature07548 [pii];doi: <u>10.1038/nature07548</u> PMID: <u>19060881</u>
- Hammad H, Chieppa M, Perros F, Willart MA, Germain RN, Lambrecht BN (2009) House dust mite allergen induces asthma via Toll-like receptor 4 triggering of airway structural cells. Nat Med 15: 410– 416. nm.1946 [pii];doi: <u>10.1038/nm.1946</u> PMID: <u>19330007</u>
- 14. Braun-Fahrlander C, Riedler J, Herz U, Eder W, Waser M, Grize L, Maisch S, Carr D, Gerlach F, Bufe A, Lauener RP, Schierl R, Renz H, Nowak D, von ME (2002) Environmental exposure to endotoxin and

its relation to asthma in school-age children. N Engl J Med 347: 869–877. doi: <u>10.1056/</u> NEJMoa020057347/12/869 [pii]. PMID: <u>12239255</u>

- Heederik D, von ME (2012) Does diversity of environmental microbial exposure matter for the occurrence of allergy and asthma? J Allergy Clin Immunol 130: 44–50. S0091-6749(12)00255-2 [pii];doi: <u>10.</u> <u>1016/j.jaci.2012.01.067</u> PMID: <u>22502794</u>
- Pillay J, Ramakers BP, Kamp VM, Lo ATL, Lam SW, Hietbrink F, Leenen LP, Tool AT, Pickkers P, Koenderman L (2010) Functional heterogenecity and differential priming of circulating neutrophils in human experimental endotoxemia. J Leukoc Biol 88: doi: <u>10.1189/jlb.1209793</u>
- Hernandez ML, Herbst M, Lay JC, Alexis NE, Brickey WJ, Ting JP, Zhou H, Peden DB (2012) Atopic asthmatic patients have reduced airway inflammatory cell recruitment after inhaled endotoxin challenge compared with healthy volunteers. J Allergy Clin Immunol 130: 869–876. S0091-6749(12)00864-0 [pii];doi: 10.1016/j.jaci.2012.05.026 PMID: 22770265
- Schaumann F, Muller M, Braun A, Luettig B, Peden DB, Hohlfeld JM, Krug N (2008) Endotoxin augments myeloid dendritic cell influx into the airways in patients with allergic asthma. Am J Respir Crit Care Med 177: 1307–1313. 200706-870OC [pii];doi: 10.1164/rccm.200706-870OC PMID: 18388357
- Peden DB, Tucker K, Murphy P, Newlin-Clapp L, Boehlecke B, Hazucha M, Bromberg P, Reed W (1999) Eosinophil influx to the nasal airway after local, low-level LPS challenge in humans. J Allergy Clin Immunol 104: 388–394. PMID: <u>10452761</u>
- Danuser B, Rebsamen H, Weber C, Krueger H (2000) Lipopolysaccharide-induced nasal cytokine response: a dose-response evaluation. Eur Arch Otorhinolaryngol 257: 527–532. PMID: <u>11195030</u>
- Nita I, Hollander C, Westin U, Janciauskiene S-M (2005) Prolastin, a pharmaceutical preparation of purified human alpha1-antitrypsin, blocks endotoxin-mediated cytokine release. Respiratory Research 6: 1–11.
- Doreswamy V, Alexis NE, Zhou H, Peden DB (2011) Nasal PMN response to repeated challenge with endotoxin in healthy volunteers. Inhal Toxicol 23: 142–147. doi: <u>10.3109/08958378.2011.553247</u> PMID: <u>21391782</u>
- Ekman AK, Fransson M, Rydberg C, Adner M, Cardell LO (2009) Nasal challenge with LPS stimulates the release of macrophage inflammatory protein 1alpha. Int Arch Allergy Immunol 149: 154–160. 000189199 [pii];doi: 10.1159/000189199 PMID: 19127073
- Virtala R, Ekman AK, Jansson L, Westin U, Cardell LO (2012) Airway inflammation evaluated in a human nasal lipopolysaccharide challenge model by investigating the effect of a CXCR2 inhibitor. Clin Exp Allergy 42: 590–596. doi: <u>10.1111/j.1365-2222.2011.03921.x</u> PMID: <u>22192144</u>
- Bachar O, Gustafsson J, Jansson L, Adner M, Cardell LO (2007) Lipopolysaccharide administration to the allergic nose contributes to lower airway inflammation. Clin Exp Allergy 37: 1773–1780. CEA2842 [pii];doi: 10.1111/j.1365-2222.2007.02842.x PMID: 17941911
- Ekman AK, Virtala R, Fransson M, Adner M, Benson M, Jansson L, Cardell LO (2012) Systemic Up-Regulation of TLR4 Causes Lipopolysaccharide-Induced Augmentation of Nasal Cytokine Release in Allergic Rhinitis. Int Arch Allergy Immunol 159: 6–14. 000335196 [pii];doi: <u>10.1159/000335196</u> PMID: <u>22555057</u>
- Folsgaard NV, Chawes BL, Rasmussen MA, Bischoff AL, Carson CG, Stokholm J, Pedersen L, Hansel TT, Bonnelykke K, Brix S, Bisgaard H (2011) Maternal atopic skewing of the neonatal nasal cytokine signature. Am J Resp Crit Care Med.
- Chawes BL, Edwards MJ, Shamji B, Walker C, Nicholson GC, Tan AJ, Folsgaard NV, Bonnelykke K, Bisgaard H, Hansel TT (2010) A novel method for assessing unchallenged levels of mediators in nasal epithelial lining fluid. J Allergy Clin Immunol 125: 1387–1389. S0091-6749(10)00172-7 [pii];doi: <u>10.</u> <u>1016/j.jaci.2010.01.039</u> PMID: <u>20304470</u>
- Nicholson GC, Kariyawasam HH, Tan AJ, Hohlfeld JM, Quinn D, Walker C, Rodman D, Westwick J, Jurcevic S, Kon OM, Barnes PJ, Krug N, Hansel TT (2011) The effects of an anti-IL-13 mAb on cytokine levels and nasal symptoms following nasal allergen challenge. J Allergy Clin Immunol 128: 800–807. S0091-6749(11)00770-6 [pii];doi: <u>10.1016/j.jaci.2011.05.013</u> PMID: <u>21719078</u>
- Scadding GW, Calderon MA, Bellido V, Koed GK, Nielsen NC, Lund K, Togias A, Phippard D, Turka LA, Hansel TT, Durham SR, Wurtzen PA (2012) Optimisation of grass pollen nasal allergen challenge for assessment of clinical and immunological outcomes. J Immunol Methods 384: 25–32. S0022-1759 (12)00174-3 [pii];doi: 10.1016/j.jim.2012.06.013 PMID: 22759401
- Dreskin SC, Dale SN, Foster SM, Martin D, Buchmeier A, Nelson HS (2002) Measurement of changes in mRNA for IL-5 in noninvasive scrapings of nasal epithelium taken from patients undergoing nasal allergen challenge. J Immunol Methods 268: 189–195. S0022175902002065 [pii]. PMID: <u>12215387</u>
- 32. Proud D, Turner RB, Winther B, Wiehler S, Tiesman JP, Reichling TD, Juhlin KD, Fulmer AW, Ho BY, Walanski AA, Poore CL, Mizoguchi H, Jump L, Moore ML, Zukowski CK, Clymer JW (2008) Gene expression profiles during in vivo human rhinovirus infection: insights into the host response. Am J

Respir Crit Care Med 178: 962–968. 200805-670OC [pii];doi: <u>10.1164/rccm.200805-670OC</u> PMID: <u>18658112</u>

- 33. Wang SZ, Ma FM, Zhao JD (2012) Expressions of nuclear factor-kappa B p50 and p65 and their significance in the up-regulation of intercellular cell adhesion molecule-1 mRNA in the nasal mucosa of allergic rhinitis patients. Eur Arch Otorhinolaryngol. doi: 10.1007/s00405-012-2136-y
- Greiff L, Pipkorn U, Alkner U, Persson CG (1990) The 'nasal pool' device applies controlled concentrations of solutes on human nasal airway mucosa and samples its surface exudations/secretions. Clin Exp Allergy 20: 253–259. PMID: <u>2364306</u>
- **35.** Virtala R, Ekman AK, Jansson L, Westin U, Cardell LO (2012) Airway inflammation evaluated in a human nasal lipopolysaccharide challenge model by investigating the effect of a CXCR2 inhibitor. Clin Exp Allergy 42: 590–596. doi: <u>10.1111/j.1365-2222.2011.03921.x</u> PMID: <u>22192144</u>
- Fransson M, Adner M, Uddman R, Cardell LO (2007) Lipopolysaccharide-induced down-regulation of uteroglobin in the human nose. Acta Otolaryngol 127: 285–291. 771160953 [pii];doi: <u>10.1080/</u> 00016480600801340 PMID: <u>17364366</u>
- Riechelmann H, Deutschle T, Friemel E, Gross HJ, Bachem M (2003) Biological markers in nasal secretions. Eur Respir J 21: 600–605. PMID: <u>12762342</u>
- Erin EM, Zacharasiewicz AS, Nicholson GC, Higgins LA, Tan AJ, Williams TJ, Murdoch RD, Durham SR, Barnes PJ, Hansel TT (2005) Topical corticosteroids inhibit IL-4, IL-5, and IL-13 levels in nasal secretions following local allergen challenge. Clin Exp Allergy 35: 1608–1614. PMID: <u>16393327</u>
- Hentschel J, Muller U, Doht F, Fischer N, Boer K, Sonnemann J, Hipler C, Hunniger K, Kurzai O, Markert UR, Mainz JG (2014) Influences of nasal lavage collection-, processing- and storage methods on inflammatory markers—evaluation of a method for non-invasive sampling of epithelial lining fluid in cystic fibrosis and other respiratory diseases. J Immunol Methods 404: 41–51. S0022-1759(13)00367-0 [pii];doi: 10.1016/j.jim.2013.12.003 PMID: 24370751
- Habibzay M, Saldana JI, Goulding J, Lloyd CM, Hussell T (2012) Altered regulation of Toll-like receptor responses impairs antibacterial immunity in the allergic lung. Mucosal Immunol 5: 524–534. mi201228 [pii];doi: 10.1038/mi.2012.28 PMID: 22549744
- Tan AM, Chen HC, Pochard P, Eisenbarth SC, Herrick CA, Bottomly HK (2010) TLR4 signaling in stromal cells is critical for the initiation of allergic Th2 responses to inhaled antigen. J Immunol 184: 3535–3544. jimmunol.0900340 [pii];doi: <u>10.4049/jimmunol.0900340</u> PMID: <u>20194715</u>