

Bone-anchored prostheses from rehabilitation and beyond: is what you see is what you get?

Laurent Frossard (PhD)

Queensland University of Technology, Brisbane, QLD, Australia
University of the Sunshine Coast, Maroochydore, QLD, Australia
YourResearchProject Pty Ltd, Brisbane, QLD, Australia

Frossard L. Bone-anchored prostheses from rehabilitation and beyond: is what you see is what you get? In *Cutting Edge Developments in P & O*. 1st Annual Scientific Meeting of Rehabilitation Medicine Society of Australia and New Zealand (RMSANZ). 2016. Melbourne, Australia. p 2

Bio

Dr Laurent Frossard is currently an Adjunct Professor at the Queensland University of Technology (QUT) and University of Sunshine Coast (USC) as well as the Director/Chief Scientist Officer at YourResearchProject.

He is a Biomechanist focusing on the locomotion and rehabilitation of individuals with lower limb loss. He is one of the very few independent experts in the clinical benefits of bone-anchored prostheses.

His academic track record includes over 130 publications, several grants, supervisions of students and international collaborations (www.LaurentFrossard.com).

Presentation Synopsis

Individuals with limb amputation fitted with conventional socket-suspended prostheses often experience socket-related discomfort leading to a significant decrease in quality of life.^[1-3]

Bone-anchored prostheses are increasingly acknowledged as viable alternative method of attachment of artificial limb.^[4-8] In this case, the prosthesis is attached directly to the residual skeleton through a percutaneous fixation.^[5, 9-12]

To date, a few osseointegrated devices are

commercially available, relying either on screw-type fixations or press-fit implants.^[4, 5, 9, 10, 13-27]

Several devices are at different stages of development particularly in Europe and the US.^[12, 18, 21, 26-43]

Clearly, surgical procedures are currently blooming worldwide. Indeed, Australia is one of the fastest growing populations worldwide.

Previous studies about bone-anchored prostheses have focused on fragmented surgical and biomechanics aspects as well as specific clinical benefits and safety of the procedure.^[6, 7, 13-15, 17, 19, 40, 44-87]

However, very few publications have synthesized this information and provided an overview of the current international developments in bone-anchored prostheses.

The purpose of the presentation will be to provide an overview of the state-of-art developments in bone-anchored prostheses with as strong emphasis on various treatments (e.g., design of fixations, surgical procedures), rehabilitation programs and the outcomes (e.g., benefits, harms).

Connect with the RMSANZ on social media

- Twitter: @LaurentFrossard

- LinkedIn: <https://au.linkedin.com/in/LaurentFrossard>
- Facebook: www.facebook.com/YourResearchProject
- Google+: <https://www.google.com/+Yourresearchprojectgoogle>
- YouTube: <https://www.youtube.com/channel/UCpdOB8YxO6bgz9YKz2u1Zyg>
- ePrint: http://eprints.qut.edu.au/view/person/Frossard_Laurent.html
- Google Scholar: <http://scholar.google.ca/citations?user=KVmDLtIAAAAJ>
- ResearchGate: http://www.researchgate.net/profile/Laurent_Frossard
- ORCID: <http://orcid.org/0000-0002-0248-9589>

References

1. Hagberg, K. and R. Brånemark, Consequences of non-vascular trans-femoral amputation: A survey of quality of life, prosthetic use and problems. *Prosthetics and Orthotics International*, 2001. 25(3): p. 186-194.
2. Mak, A.F., M. Zhang, and D.A. Boone, State-of-the-art research in lower-limb prosthetic biomechanics-socket interface: a review. *J Rehabil Res Dev*, 2001. 38(2): p. 161-74.
3. Bui, K.M., G.J. Raugi, V.Q. Nguyen, and G.E. Reiber, Skin problems in individuals with lower-limb loss: literature review and proposed classification system. *J Rehabil Res Dev*, 2009. 46(9): p. 1085-90.
4. Pitkin, M., One lesson from arthroplasty to osseointegration in search for better fixation of in-bone implanted prosthesis. *Journal of Rehabilitation Research & Development*, 2008. 45(4): p. 6-14.
5. Brånemark, R., P.I. Brånemark, B. Rydevik, and R.R. Myers, Osseointegration in skeletal reconstruction and rehabilitation: a review. *J Rehabil Res Dev*, 2001. 38(2): p. 175-81.
6. van Eck, C.F. and R.L. McGough, Clinical outcome of osseointegrated prostheses for lower extremity amputations: a systematic review of the literature. *Current Orthopaedic Practice*, 2015. 26(4): p. 349-357.
7. Leijendekkers, R.A., G. van Hinte, J.P. Frolke, H. van de Meent, M.W. Nijhuis-van der Sanden, and J.B. Staal, Comparison of bone-anchored prostheses and socket prostheses for patients with a lower extremity amputation: a systematic review. *Disabil Rehabil*, 2016: p. 1-14.
8. Potter, B., From Bench to Bedside: A Perfect Fit? Osseointegration Can Improve Function for Patients with Amputations. *Clinical Orthopaedics and Related Research*, 2015: p. 1-3.
9. Pitkin, M., Design features of implants for direct skeletal attachment of limb prostheses. *J Biomed Mater Res A*, 2013. 101(11): p. 3339-48.
10. Pitkin, M., C. Cassidy, R. Muppavarapu, R. J. M. Shevtsov, O. Galibin, and S. Rousselle, New method of fixation of in-bone implanted prosthesis. *J Rehabil Res Dev*, 2013. 50(5): p. 709-722.
11. Monument, M.J., D.M. Lerman, and R.L. Randall, Novel Applications of Osseointegration in Orthopedic Limb Salvage Surgery. *Orthopedic Clinics of North America*, 2014(0).
12. Rubin, L., L. Kennon, J. Keggi, and H. Aschoff, Surgical management of trans-femoral amputation with a

- transcutaneous, press-fit distal femoral intra-medullary device: analysis with minimum 2 year follow-up. in *Journal of Bone & Joint Surgery, British Volume*. 2012.
13. Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. *Bone Joint J*, 2014. 96(1): p. 106-113.
14. Hagberg, K., E. Hansson, and R. Branemark, Outcome of Percutaneous Osseointegrated Prostheses for Patients With Unilateral Transfemoral Amputation at Two-Year Follow-Up. *Arch Phys Med Rehabil*, 2014. 95(11): p. 2120-2127.
15. Aschoff, H.H., R.E. Kennon, J.M. Keggi, and L.E. Rubin, Transcutaneous, distal femoral, intramedullary attachment for above-the-knee prostheses: an endo-exo device. *J Bone Joint Surg Am*, 2010. 92 Suppl 2(Supplement 2): p. 180-6.
16. Aschoff, H.-H. and R. McGough, The Endo-Exo Femoral Prosthesis: a new rehabilitation concept following above knee amputation. *Journal of Bone & Joint Surgery, British Volume*, 2012. 94-B(SUPP XXXIX): p. 77.
17. Van de Meent, H., M.T. Hopman, and J.P. Frolke, Walking ability and quality of life in subjects with transfemoral amputation: a comparison of osseointegration with socket prostheses. *Arch Phys Med Rehabil*, 2013. 94(11): p. 2174-2178.
18. Hillock, R., J. Keggi, R. Kennon, E. McPherson, T. Clyburn, D. Brazil, and T. McTighe, A Global Collaboration - Osteointegration Implant (OI) for Transfemoral Amputation - Case Report. *Reconstructive review - Joint Implant Surgery & Research Foundation*, 2013. 3(2): p. 50-54.
19. Juhnke, D.-L., J.P. Beck, S. Jeyapalina, and H.H. Aschoff, Fifteen years of experience with Integral-Leg-Prosthesis: Cohort study of artificial limb attachment system. *Journal of Rehabilitation Research and Development*, 2015. 52(4): p. 407-420.
20. Isaacson, B. and S. Jeyapalina, Osseointegration: A review of the fundamentals for assuring cementless skeletal fixation. *Orthopedic Research and Reviews*, 2014. 6: p. 55-65.
21. Pitkin, M., G. Raykhtsaum, O.V. Galibin, M.V. Protasov, J.V. Chihovskaya, and I.G. Belyaeva, Skin and bone integrated prosthetic pylon: a pilot animal study. *J Rehabil Res Dev*, 2006. 43(4): p. 573-80.
22. Pitkin, M., On the way to total integration of prosthetic pylon with residuum. *J Rehabil Res Dev*, 2009. 46(3): p. 345-60.
23. Pitkin, M., J. Pilling, and G. Raykhtsaum, Mechanical properties of totally permeable titanium composite pylon for direct skeletal attachment. *J Biomed Mater Res B Appl Biomater*, 2012. 100(4): p. 993-9.
24. Pitkin, M., C. Cassidy, R. Muppavarapu, and D. Edell, Recording of Electric Signal Passage Trough a Pylon in Direct Skeletal Attachment of Leg Prostheses. *IEEE Transaction and Biomedical Engineering*, 2012. 59(5): p. 1349-1353.
25. Shevtsov, M.A., O.V. Galibin, N.M. Yudintceva, M.I. Blinova, G.P.

- Pinaev, A.A. Ivanova, O.N. Savchenko, D.N. Suslov, I.L. Potokin, E. Pitkin, G. Raykhtsaum, and M.R. Pitkin, Two-stage implantation of the skin- and bone-integrated pylon seeded with autologous fibroblasts induced into osteoblast differentiation for direct skeletal attachment of limb prostheses. *J Biomed Mater Res A*, 2014. 102(9): p. 3033-48.
26. Farrell, B.J., B.I. Prilutsky, R.S. Kistenberg, J.F. Dalton, and M. Pitkin, An animal model to evaluate skin-implant-bone integration and gait with a prosthesis directly attached to the residual limb. *Clinical biomechanics*, 2014.
27. Piotr Prochor, S.P., Eugeniusz Sajewicz, Biomechanical evaluation of a novel Limb Prosthesis Osseointegrated Fixation System designed to combine the advantages of interference-fit and threaded solutions. 2016.
28. Pendegrass, C.J., A.E. Goodship, and G.W. Blunn, Development of a soft tissue seal around bone-anchored transcutaneous amputation prostheses. *Biomaterials*, 2006. 27(23): p. 4183-91.
29. Kang, N.V., C. Pendegrass, L. Marks, and G. Blunn, Osseocutaneous integration of an intraosseous transcutaneous amputation prosthesis implant used for reconstruction of a transhumeral amputee: Case report. *The Journal of Hand Surgery*, 2010. 35(7): p. 1130-1134.
30. Fitzpatrick, N., T.J. Smith, C.J. Pendegrass, R. Yeadon, M. Ring, A.E. Goodship, and G.W. Blunn, Intraosseous transcutaneous amputation prosthesis (ITAP) for limb salvage in 4 dogs. *Vet Surg*, 2011. 40(8): p. 909-25.
31. Guirao, L., C. Samitier, M. Costea, J. Camos, M. Majo, and E. Pleguezuelos, Improvement in walking abilities in transfemoral amputees with a distal weight bearing implant. *Prosthetics and Orthotics International*, 2016.
32. Saunders, M.M., J.S. Brecht, M.C. Verstraete, D.B. Kay, and G.O. Njus, Lower limb direct skeletal attachment. A Yucatan micropig pilot study. *J Invest Surg*, 2012. 25(6): p. 387-97.
33. Hugate, R., R. Clarke, T. Hoeman, and A. Friedman, Transcutaneous Implants in a Porcine Model: The Use of Highly Porous Tantalum. *International Journal of Advanced Materials Research*, 2015. 1(2): p. 32-40.
34. Tomaszewski, P.K., M. van Diest, S.K. Bulstra, N. Verdonschot, and G.J. Verkerke, Numerical analysis of an osseointegrated prosthesis fixation with reduced bone failure risk and periprosthetic bone loss. *J Biomech*, 2012. 45(11): p. 1875-80.
35. Tomaszewski, P.K., N. Verdonschot, S.K. Bulstra, J.S. Rietman, and G.J. Verkerke, Simulated bone remodeling around two types of osseointegrated implants for direct fixation of upper-leg prostheses. *J Mech Behav Biomed Mater*, 2012. 15(0): p. 167-75.
36. Tomaszewski, P.K., B. Lasnier, G. Hannink, G.J. Verkerke, and N. Verdonschot, Experimental assessment of a new direct fixation implant for artificial limbs. *J Mech Behav Biomed Mater*, 2013. 21: p. 77-85.
37. Tomaszewski, P.K., N. Verdonschot, S.K. Bulstra, and G.J. Verkerke, A New Osseointegrated Fixation Implant for Amputated Patients.

- Journal of Biomechanics, 2012. 45(0): p. S322.
38. Shelton, T.J., J.P. Beck, R.D. Bloebaum, and K.N. Bachus, Percutaneous osseointegrated prostheses for amputees: Limb compensation in a 12-month ovine model. *J Biomech*, 2011. 44(15): p. 2601-6.
39. Jeyapalina, S., J.P. Beck, K.N. Bachus, O. Chalayan, and R.D. Bloebaum, Radiographic Evaluation of Bone Adaptation Adjacent to Percutaneous Osseointegrated Prostheses in a Sheep Model. *Clin Orthop Relat Res*, 2014: p. 1-12.
40. Khemka, A., L. Frossard, S.J. Lord, B. Bosley, and M. Al Muderis, Osseointegrated total knee replacement connected to a lower limb prosthesis: 4 cases. *Acta Orthop*, 2015: p. 1-5.
41. Hillock, R., D. Tatum, and D. E, Osseointegration Implant Post Coupling With External Prosthetic Limb Continuation of Previous Case Reports "Stage III". *Reconstructive review - Joint Implant Surgery & Research Foundation*, 2014. 4(2): p. 13-16.
42. Zak, A. and R. Hillock, Anesthesia and Pain Control for Osteointegration Implantation Into the Femur - Stage II Osteointegration Implant (OI) Skin Coupling Procedure. *Joint Implant Surgery & Research Foundation's*, 2014. 3(4): p. 67-69.
43. Hillock, R., J. Keggi, R. Kennon, E. McPherson, T. Clyburn, D. Brazil, and T. McTighe, Stage II Osteointegration Implant (OI) Skin Coupling Procedure. *Reconstructive review - Joint Implant Surgery & Research Foundation*, 2013. 3(4): p. 62-66.
44. Helgason, B., H. Palsson, T.P. Runarsson, L. Frossard, and M. Viceconti, Risk of failure during gait for direct skeletal attachment of a femoral prosthesis: a finite element study. *Med Eng Phys*, 2009. 31(5): p. 595-600.
45. Lee, W.C., J.M. Doocey, R. Branemark, C.J. Adam, J.H. Evans, M.J. Percy, and L.A. Frossard, FE stress analysis of the interface between the bone and an osseointegrated implant for amputees--implications to refine the rehabilitation program. *Clin Biomech (Bristol, Avon)*, 2008. 23(10): p. 1243-50.
46. Hagberg, K., E. Haggstrom, M. Uden, and R. Branemark, Socket versus bone-anchored trans-femoral prostheses: hip range of motion and sitting comfort. *Prosthet Orthot Int*, 2005. 29(2): p. 153-163.
47. Frossard, L., N. Stevenson, J. Smeathers, E. Haggstrom, K. Hagberg, J. Sullivan, D. Ewins, D. Lee Gow, S. Gray, and R. Branemark, Monitoring of the load regime applied on the osseointegrated fixation of a trans-femoral amputee: a tool for evidence-based practice. *Prosthet Orthot Int*, 2008. 32(1): p. 68-78.
48. Lee, W., L. Frossard, K. Hagberg, E. Haggstrom, and R. Brånemark, Kinetics analysis of transfemoral amputees fitted with osseointegrated fixation performing common activities of daily living. *Clinical Biomechanics*, 2007. 22(6): p. 665-673.
49. Hagberg, K. and R. Branemark, One hundred patients treated with osseointegrated transfemoral amputation prostheses-rehabilitation perspective. *J Rehabil Res Dev*, 2009. 46(3): p. 331-44.


50. Hagberg, K., R. Branemark, B. Gunterberg, and B. Rydevik, Osseointegrated trans-femoral amputation prostheses: Prospective results of general and condition-specific quality of life in 18 patients at 2-year follow-up. *Prosthetics and Orthotics International*, 2008. 32(1): p. 29 - 41.
51. Frossard, L., K. Hagberg, E. Häggström, D. Lee Gow, R. Brånemark, and M. Pearcy, Functional Outcome of Transfemoral Amputees Fitted With an Osseointegrated Fixation: Temporal Gait Characteristics. *JPO Journal of Prosthetics and Orthotics*, 2010. 22(1): p. 11-20.
52. Jacobs, R., R. Brånemark, K. Olmarker, B. Rydevik, D. Van Steenberghe, and P. Branemark, Evaluation of the psychophysical detection threshold level for vibrotactile and pressure stimulation of prosthetic limbs using bone anchorage or soft tissue support. *Prosthetics and Orthotics International*, 2000. 24(2): p. 133-142.
53. Pantall, A., S. Durham, and D. Ewins, Variability of SEMG of five stump muscles during stance phase of gait in TF amputees with osseointegrated prostheses. *Gait & Posture*, 2006. 24(Supplement 2): p. S153-S154.
54. Tranberg, R., R. Zugner, and J. Karrholm, Improvements in hip- and pelvic motion for patients with osseointegrated trans-femoral prostheses. *Gait Posture*, 2011. 33(2): p. 165-8.
55. Hagberg, K., E. Haggstrom, and R. Branemark, Physiological cost index (PCI) and walking performance in individuals with transfemoral prostheses compared to healthy controls. *Disabil Rehabil*, 2007. 29(8): p. 643-9.
56. Pantall, A., S. Durham, and D. Ewins, Surface electromyographic activity of five residual limb muscles recorded during isometric contraction in transfemoral amputees with osseointegrated prostheses. *Clin Biomech (Bristol, Avon)*, 2011. 26(7): p. 760-5.
57. Lee, W., L. Frossard, K. Hagberg, E. Haggstrom, D. Lee Gow, S. Gray, and R. Branemark, Magnitude and variability of loading on the osseointegrated implant of transfemoral amputees during walking. *Med Eng Phys*, 2008. 30(7): p. 825-833.
58. Haggstrom, E.E., E. Hansson, and K. Hagberg, Comparison of prosthetic costs and service between osseointegrated and conventional suspended transfemoral prostheses. *Prosthet Orthot Int*, 2013. 37(2): p. 152-60.
59. Pantall, A. and D. Ewins, Muscle activity during stance phase of walking: Comparison of males with transfemoral amputation with osseointegrated fixations to nondisabled male volunteers. *Journal of Rehabilitation Research & Development*, 2013. 50(4): p. 499-514.
60. Haggstrom, E., K. Hagberg, B. Rydevik, and R. Branemark, Vibrotactile evaluation: osseointegrated versus socket-suspended transfemoral prostheses. *J Rehabil Res Dev*, 2013. 50(10): p. 1423-34.
61. Schalk, S., N. Jonkergouw, F. van der Meer, W.M. Swaan, H.-H. Aschoff, and P. van der Wurff, The Evaluation of Daily Life Activities after Application of an Osseointegrated Prosthesis Fixation

- in a Bilateral Transfemoral Amputee: A Case Study. *Medicine*, 2015. 94(36): p. e1416.
62. Al Muderis, M., H. Aschoff, B. Bosley, G. Raz, L. Gerdesmeyer, and B. Burkett, Direct skeletal attachment prosthesis for the amputee athlete: the unknown potential. *Sports Engineering*, 2016: p. 1-5.
63. Al Muderis, M., K. Tetsworth, A. Khemka, S. Wilmot, B. Bosley, S.J. Lord, and V. Glatt, The Osseointegration Group of Australia Accelerated Protocol (OGAAP-1) for two-stage osseointegrated reconstruction of amputated limbs. *Bone & Joint Journal*, 2016. 98-B(7): p. 952-960.
64. Frossard, L.A., R. Tranberg, E. Haggstrom, M. Percy, and R. Branemark, Load on osseointegrated fixation of a transfemoral amputee during a fall: loading, descent, impact and recovery analysis. *Prosthet Orthot Int*, 2010. 34(1): p. 85-97.
65. Tomaszewski, P.K., N. Verdonshot, S.K. Bulstra, and G.J. Verkerke, A comparative finite-element analysis of bone failure and load transfer of osseointegrated prostheses fixations. *Ann Biomed Eng*, 2010. 38(7): p. 2418-27.
66. Tillander, J., K. Hagberg, L. Hagberg, and R. Branemark, Osseointegrated titanium implants for limb prostheses attachments: infectious complications. *Clin Orthop Relat Res*, 2010. 468(10): p. 2781-8.
67. Thompson, M.L., D. Backman, R. Branemark, and C.K. Mechefske, Evaluating the bending response of two osseointegrated transfemoral implant systems using 3D digital image correlation. *Journal of Biomechanical Engineering*, 2011. 133(5): p. 051006.
68. Nebergall, A., C. Bragdon, A. Antonellis, J. Karrholm, R. Branemark, and H. Malchau, Stable fixation of an osseointegrated implant system for above-the-knee amputees: titel RSA and radiographic evaluation of migration and bone remodeling in 55 cases. *Acta Orthop*, 2012. 83(2): p. 121-8.
69. Berlin, Ö., P. Bergh, M. Dalen, S. Eriksson, K. Hagberg, S. Inerot, B. Gunterberg, and R. Brånemark, Osseointegration in transfemoral amputees: the gothenburg experience, in *Journal of Bone & Joint Surgery, British Volume* 2012. p. 55.
70. Jeyapalina, S., J.P. Beck, K.N. Bachus, and R.D. Bloebaum, Cortical bone response to the presence of load-bearing percutaneous osseointegrated prostheses. *Anat Rec (Hoboken)*, 2012. 295(9): p. 1437-45.
71. Frossard, L., R. Tranberg, E. Haggstrom, M. Percy, and R. Branemark, Fall of a transfemoral amputee fitted with osseointegrated fixation: loading impact on residuum. *Gait & Posture*, 2009. 30(Supplement 2): p. S151-S152.
72. Isaacson, B.M., R.E. Vance, T.G. Chou, R.D. Bloebaum, K.N. Bachus, and J.B. Webster, Effectiveness of resonance frequency in predicting orthopedic implant strength and stability in an in vitro osseointegration model. *J Rehabil Res Dev*, 2009. 46(9): p. 1109-20.
73. Schwarze, M., C. Hurschler, F. Seehaus, T. Correa, and B. Welke, Influence of transfemoral amputation length on resulting loads at the osseointegrated prosthesis fixation

- during walking and falling. *Clinical Biomechanics*, 2013. In press(0).
74. Palmquist, A., S.H. Windahl, B. Norlindh, R. Branemark, and P. Thomsen, Retrieved bone-anchored percutaneous amputation prosthesis showing maintained osseointegration after 11 years-a case report. *Acta Orthop*, 2014. 85(4): p. 442-5.
75. Tsikandylakis, G., O. Berlin, and R. Branemark, Implant Survival, Adverse Events, and Bone Remodeling of Osseointegrated Percutaneous Implants for Transhumeral Amputees. *Clin Orthop Relat Res*, 2014: p. 1-10.
76. Xu, W. and K. Robinson, X-ray image review of the bone remodeling around an osseointegrated trans-femoral implant and a finite element simulation case study. *Ann Biomed Eng*, 2008. 36(3): p. 435-43.
77. Mollaret, E. and M. Bertrand Marchand, Osseointegration of the lower limb: Reduction of mechanical strain on the abutment. *Annals of Physical and Rehabilitation Medicine*, 2015. 58, Supplement 1: p. e18.
78. Al Muderis, M., A. Khemka, S.J. Lord, H. Van de Meent, and J.P.M. Fro lke, Safety of Osseointegrated Implants for Transfemoral Amputees: A Two-Center Prospective Cohort Study. *The Journal of Bone & Joint Surgery*, 2016. 98(11): p. 900-909.
79. Haket, L.M., J.P. Frolke, N. Verdonschot, P.K. Tomaszewski, and H. van de Meent, Periprosthetic cortical bone remodeling in patients with an osseointegrated leg prosthesis. *J Orthop Res*, 2016.
80. Frossard, L., K. Hagberg, E. Haggstrom, and R. Branemark, Load-relief of walking aids on osseointegrated fixation: instrument for evidence-based practice. *IEEE Trans Neural Syst Rehabil Eng*, 2009. 17(1): p. 9-14.
81. Frossard, L., D.L. Gow, K. Hagberg, N. Cairns, B. Contoyannis, S. Gray, R. Branemark, and M. Pearcy, Apparatus for monitoring load bearing rehabilitation exercises of a transfemoral amputee fitted with an osseointegrated fixation: a proof-of-concept study. *Gait Posture*, 2010. 31(2): p. 223-8.
82. Frossard, L.A., Load on osseointegrated fixation of a transfemoral amputee during a fall: Determination of the time and duration of descent. *Prosthet Orthot Int*, 2010. 34(4): p. 472-87.
83. Frossard, L., N. Stevenson, J. Sullivan, M. Uden, and M. Pearcy, Categorization of Activities of Daily Living of Lower Limb Amputees During Short-Term Use of a Portable Kinetic Recording System: A Preliminary Study. *JPO Journal of Prosthetics and Orthotics*, 2011. 23(1): p. 2-11.
84. Frossard, L., E. Haggstrom, K. Hagberg, and P. Branemark, Load applied on a bone-anchored transfemoral prosthesis: characterisation of prosthetic components – A case study *Journal of Rehabilitation Research & Development*, 2013. 50(5): p. 619–634.
85. Vertriest, S., P. Coorevits, K. Hagberg, R. Branemark, E. Haggstrom, G. Vanderstraeten, and L. Frossard, Static load bearing exercises of individuals with transfemoral amputation fitted with an osseointegrated implant: reliability of kinetic data. *IEEE Trans Neural Syst Rehabil Eng*, 2015. 23(3): p. 423-30.

86. Dumas, R., R. Branemark, and L. Frossard, Gait analysis of transfemoral amputees: errors in inverse dynamics are substantial and depend on prosthetic design. *IEEE Trans Neural Syst Rehabil Eng*, 2016. PP(99): p. 1-1.
87. Vertriest, S., P. Coorevits, K. Hagberg, R. Branemark, E.E. Haggstrom, G. Vanderstraeten, and L. Frossard, Static load bearing exercises of individuals with transfemoral amputation fitted with an osseointegrated implant: Loading compliance. *Prosthet Orthot Int*, 2016.

Bone-anchored prostheses from rehabilitation and beyond: is what you see is what you get?



Adj/Prof Laurent Frossard

Queensland University of Technology,
Australia

University of the Sunshine Coast,
Australia

www.LaurentFrossard.com

Rehabilitation Medicine Society of
Australia and New Zealand

1st Annual Scientific Meeting

19 October 2016

Crown Promenade, Melbourne

Disclosure

Conflict

- No conflicts of interest
- No financial benefits associated with BAP
- No financial benefits associated with treatments
- Done some research work for:
 - Sahlgrenska University Hospital, Sweden
 - Ghent University Hospital, Belgium
- Done some consulting work for:
 - Munjed Al Muderis – Orthopaedic Surgery
 - University of the Sunshine Coast

Bone-anchored prostheses from rehabilitation and beyond: is what you see is what you get?

Background

Attachment

Prosthetic attachment options
Example: Young active trauma



<http://www.burgalaw.com/practice-areas/rancho-cucamonga-personal-injury-attorney/>

Background

Attachment

Socket




<http://www.laurentfrossard.com/biomechanics-of-individuals-with-limb-loss/population>

Background

Attachment

Socket

BAP



<http://www.laurentfrossard.com/biomechanics-of-individuals-with-limb-loss/population>

Background


Attachment

Socket

BAP

Screw-type

Osseointegrated Prosthesis for Rehabilitation of Amputees, (OPRA), Integrum Ab, Sweden



Dr Rickard
Branemark



<http://opraosseointegration.com/>

Background


Attachment

Socket

BAP

Screw-type

Osseointegrated Prosthesis for Rehabilitation of Amputees, (OPRA), Integrum Ab, Sweden



Dr Rickard
Branemark



Nebergall, A., C. Bragdon, A. Antonellis, J. Kärrholm, R. Brånemark, and H. Malchau, Stable fixation of an osseointegrated implant system for above-the-knee amputees. Acta Orthopaedica, 2012. 83(2): p. 121-128

Background

Attachment

Socket

BAP

Screw-type

Osseointegrated Prosthesis for Rehabilitation of Amputees, (OPRA), Integrum Ab, Sweden



Dr Rickard Branemark



<http://www.laurentfrossard.com/biomechanics-of-individuals-with-limb-loss/population>

Background

Attachment

Socket

BAP

Screw-type

Press-fit

Integral Leg Prosthesis (ILP), Orthodynamics Pty Ltd, UK



Dr Horst Aschoff



Lunow, C., K. Staubach, and H. Aschoff, [Endo-exo Femoral Prosthesis]. Zeitschrift Der Unfallchirurg, 2010.

Background

Attachment

Socket

BAP

Screw-type

Press-fit

Integral Leg Prosthesis
(ILP), Orthodynamics Pty Ltd, UK



Dr Horst Aschoff



Lunow, C., K. Staubach, and H. Aschoff, [Endo-exo Femoral Prosthesis]. Zeitschrift Der Unfallchirurg, 2010.

Background

Attachment


Socket

BAP


Screw-type

Press-fit

Integral Leg Prosthesis
(ILP), Orthodynamics Pty Ltd, UK



Dr Horst Aschoff



Source: Unknown

Background

Attachment

Socket

BAP

Screw-type

Press-fit

Osseointegrated Prosthesis Limb
(OPL) Permedica, Italy



Dr Munjed Al
Muderis



<http://www.osseointegrationaustralia.com.au/>

Background

Attachment

Socket

BAP

Screw-type

Press-fit

Osseointegrated Prosthesis Limb
(OPL) Permedica, Italy



Dr Munjed Al
Muderis



Al Muderis, M., K. Tetsworth, A. Khemka, S. Wilmot, B. Bosley, S.J. Lord, and V. Glatt, The Osseointegration Group of Australia Accelerated Protocol (OGAAP-1) for two-stage osseointegrated reconstruction of amputated limbs. Bone & Joint Journal, 2016. 98-B(7): p. 952-960.

Background

Attachment

Socket

	ILP	OPL	OPRA
Interface fixation - bone	Press-fit	Press-fit	Screw

BAP

Screw-type

Press-fit

Background

Attachment

Socket

	ILP	OPL	OPRA
Interface fixation - bone	Press-fit	Press-fit	Screw
Nb of surgeries	2	2	2

BAP


Screw-type

Press-fit


Stage 1

Stage 2

Background



Background



Background

Attachment

Socket

BAP

Screw-type

Press-fit

	ILP	OPL	OPRA
Interface fixation - bone	Press-fit	Press-fit	Screw
Nb of surgeries	1	2	2
Duration rehabilitation *	0	4 mth	9-12 mth
* Estimation			

Stage 1

Background

Attachment

Socket

BAP

Screw-type

Press-fit

	ILP	OPL	OPRA
Interface fixation - bone	Press-fit	Press-fit	Screw
Nb of surgeries	2	2	2
Duration rehabilitation *	4 mth	4 mth	9-12 mth
Nb of years since first S1 *	10	4	17
* Estimation			

Background

Attachment

Socket

BAP

Screw-type

Press-fit

	ILP	OPL	OPRA	
Interface fixation - bone	Press-fit	Press-fit	Screw	
Nb of surgeries	2	2	2	
Duration rehabilitation *	4 mth	4 mth	9-12 mth	
Nb of years since first S1 *	10	4	17	
Nb of patients *	130	250	300	
	200			
	650-750			

* Estimation based on 2013 data

Most published and acknowledged

Background

Attachment

Socket

“Hybrid”

BAP

Screw-type

Press-fit

Keep Walking

Benefits of the Keep Walking implant

- The ischium support is improved, alleviating or removing its load.
- Bone load improves the femur density.
- A distal load occurs.
- Greater control of the prosthesis.
- The patient needs less support.
- Less fatigue and more self reliance.

Implant

Stem

Spacer

Plug

Screw

Benefits of the Advanced extension

- Greater comfort: the use of a socket is no longer necessary.
- The Advanced extension unifies the prosthesis with the femur.
- The loads are directly transmitted to the femur.

Advanced extension

Intermediate device

Superior connecting piece

Inferior connecting piece

Prosthetic leg

http://implantekeepwalking.es

Background

Attachment

Socket


“Hybrid”

BAP

Screw-type

Press-fit

BANP



<http://www.chalmers.se/en/news/Pages/Thought-controlled-prosthesis-is-changing-the-lives-of-amputees.aspx>

Bone-anchored prostheses from
rehabilitation and beyond: is what you see
is what you get?


Rehabilitation

Pathways



Rehabilitation


Pathways



Hagberg, K. and R. Branemark, One hundred patients treated with osseointegrated transfemoral amputation prostheses--rehabilitation perspective. J Rehabil Res Dev, 2009. 46(3): p. 331-44.

Rehabilitation


Pathways



Hagberg, K. and R. Branemark, One hundred patients treated with osseointegrated transfemoral amputation prostheses--rehabilitation perspective. J Rehabil Res Dev, 2009. 46(3): p. 331-44.


Rehabilitation

Pathways




Rehabilitation

Pathways




Rehabilitation

Pathways




Rehabilitation

Pathways



Rehabilitation

Pathways




Bone remodelling
= **Right load**
+ **Right time**

Vertriest S, Coorevits P, Brånemark R, Hagberg K, Brånemark R, Vanderstraeten G, Frossard L. Static load bearing exercises of individuals with transfemoral amputation fitted with an osseointegrated implant: Reliability of kinetic data. IEEE Transactions on Neural Systems and Rehabilitation Engineering. 2015.

Rehabilitation


Pathways



Frossard, L., D.L. Gow, K. Hagberg, N. Cairns, B. Contoyannis, S. Gray, R. Brånemark, and M. Percy, Apparatus for monitoring load bearing rehabilitation exercises of a transfemoral amputee fitted with an osseointegrated fixation: A proof-of-concept study. *Gait and Posture*, 2010. 31(2): p. 223-228

Rehabilitation



Pathways



Hagberg, K. and R. Brånemark, One hundred patients treated with osseointegrated transfemoral amputation prostheses - the rehabilitation perspective. *Journal of Rehabilitation Research & Development*, 2009. 43(3): p. 331-344

Rehabilitation



Pathways



Frossard, L., K. Hagberg, E. Haggstrom, and R. Branemark, Load-relief of walking aids on osseointegrated fixation: instrument for evidence-based practice. NSRE, IEEE Transactions on, 2009. 17(1): p. 9-14

Rehabilitation


Pathways



Frossard, L., E. Haggstrom, K. Hagberg, and P. Branemark, Load applied on a bone-anchored transfemoral prosthesis: characterisation of prosthetic components – A case study Journal of Rehabilitation Research & Development, 2013. 50(5): p. 619–634.

Rehabilitation

Pathways



Frossard L, Stevenson N, Sullivan J, Uden M, Percy M. Categorisation of activities of daily living of lower limb amputees during short-term use of a portable kinetic recording system: a preliminary study. Journal of Prosthetics and Orthotics. 2011. 23 (1). p 2-11.

Bone-anchored prostheses from rehabilitation and beyond: **is what you see is what you get?**

Benefits

Benefits

HRQoL

SF-36

Health-related quality of life: SF-36

Table III. Questionnaire for Persons with a Transfemoral Amputation (Q-TFA) and Short-Form (SF)-36 scores at baseline and change from baseline to 12- and to 24-month follow-up, respectively. Three patients failed to complete the whole questionnaire at each visit

Variable	Mean score (range) [median; SD] (no. patients)		
	Baseline	Change from baseline to 12 mths	Change from baseline to 24 mths
SF-36			
Physical function	35 (0 to 85) [30; 22] (n = 51)	22 (-40 to 70) [20; 24] (n = 47) [§]	23 (-23 to 75) [25; 21] (n = 45) [§]
Role-physical	41 (0 to 100) [25; 42] (n = 50)	24 (-50 to 100) [25; 44] (n = 45) [§]	22 (-50 to 100) [13; 36] (n = 44) [§]
Bodily pain	55 (10 to 100) [51; 26] (n = 51)	7 (-52 to 74) [0; 26] (n = 47)	6 (-61 to 59) [9; 30] (n = 45)
General health	78 (37 to 100) [82; 18] (n = 51)	3 (-32 to 40) [0; 17] (n = 47)	-1 (-42 to 40) [0; 18] (n = 45)
Vitality	60 (15 to 90) [60; 20] (n = 51)	5 (-50 to 45) [5; 19] (n = 47)	3 (-70 to 45) [5; 23] (n = 45)
Social function	78 (13 to 100) [88; 25] (n = 51)	2 (-50 to 50) [0; 24] (n = 47)	1 (-100 to 63) [0; 30] (n = 45)
Role-emotional	75 (0 to 100) [100; 39] (n = 50)	5 (0 to 100) [0; 43] (n = 46)	0 (0 to 100) [0; 45] (n = 44)
Mental health	74 (4 to 100) [80; 21] (n = 51)	2 (-44 to 40) [0; 18] (n = 47)	2 (-76 to 40) [4; 24] (n = 45)
SF-36 Physical Component Summary [†]	32 (18 to 55) [31;8] (n = 50)	8 (-11 to 29) [8;10] (n = 45)	8 (-8 to 29) [9;9] (n = 44)
SF-36 Mental Component Summary [†]	53 (19 to 69) [57; 13] (n = 50)	-2 (-33 to 23) [-2; 11] (n = 45)	-3 (-44 to 22) [0; 15] (n = 44)

* a Prosthetic Use Score of 0 means the patient is not using a prosthesis and consequently the Prosthetic Mobility Score, Problem Score and Global Score could not be answered, hence results for lower numbers of patients in those scores¹
† the Problem Score is reversed, which means a lower figure indicates fewer problems related to amputation and prosthesis
‡ SF-36 Physical and Mental Component Summaries are normalised to the general population (mean 50 (SD 10))¹³
§ p < 0.001

Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113. Errata : Haddad, Bone Joint J: 2014, 96-B 106-113

Benefits

HRQoL

SF-36

Health-related quality of life: SF-36

Table III. Questionnaire for Persons with a Transfemoral Amputation (Q-TFA) and Short-Form (SF)-36 scores at baseline and change from baseline to 12- and to 24-month follow-up, respectively. Three patients failed to complete the whole questionnaire at each visit

Variable	Mean score (range) [median; SD] (no. patients)		
	Baseline	Change from baseline to 12 mths	Change from baseline to 24 mths
SF-36			
Physical function		↑	↑
Role-physical			
Bodily pain			
General health			
Vitality			
Social function			
Role-emotional			
Mental health			
SF-36 Physical Component Summary [‡]		↑	↑
SF-36 Mental Component Summary [‡]			

* a Prosthetic Use Score of 0 means the patient is not using a prosthesis and consequently the Prosthetic Mobility Score, Problem Score and Global Score could not be answered, hence results for lower numbers of patients in those scores¹
† the Problem Score is reversed, which means a lower figure indicates fewer problems related to amputation and prosthesis
‡ SF-36 Physical and Mental Component Summaries are normalised to the general population (mean 50 (SD 10))¹³
§ p < 0.001

Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113. Errata : Haddad, Bone Joint J: 2014, 96-B 106-113

Benefits

HRQoL

SF-36

Q-TFA

Health-related quality of life: Q-TFA

Table III. Questionnaire for Persons with a Transfemoral Amputation (Q-TFA) and Short-Form (SF)-36 scores at baseline and change from baseline to 12- and to 24-month follow-up, respectively. Three patients failed to complete the whole questionnaire at each visit

Variable	Mean score (range) [median; SD] (no. patients)		
	Baseline	Change from baseline to 12 mths	Change from baseline to 24 mths
Q-TFA*			
Prosthetic use score	47 (0 to 100) [52; 37] (n = 51)	34 (-23 to 100) [29; 29] (n = 44) [§]	32 (-100 to 100) [29; 41] (n = 45) [§]
Prosthetic mobility score	52 (0 to 82) [56; 20] (n = 42)	14 (-29 to 46) [15; 17] (n = 36) [§]	18 (-29 to 48) [17; 16] (n = 37) [§]
Problem score [†]	44 (5 to 77) [48; 19] (n = 42)	-28 (-57 to 2) [-33; 16] (n = 36) [§]	-27 (-59 to 7) [-30; 16] (n = 37) [§]
Global score	38 (0 to 92) [33; 19] (n = 42)	37 (-17 to 84) [34; 26] (n = 36) [§]	39 (0 to 92) [34; 24] (n = 37) [§]
Overall situation (n, %)			
Extremely poor	5 (10)	n = 42	n = 45
Poor	15 (29)	Declined: 2 (5)	Declined: 3 (7)
Average	17 (33)	No change: 11 (26)	No change: 11 (24)
Good	9 (18)	Improved: 29 (69) [§]	Improved: 31 (69) [§]
Extremely good	5 (10)		

Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Benefits

HRQoL

SF-36

Q-TFA

Health-related quality of life: Q-TFA

Table III. Questionnaire for Persons with a Transfemoral Amputation (Q-TFA) and Short-Form (SF)-36 scores at baseline and change from baseline to 12- and to 24-month follow-up, respectively. Three patients failed to complete the whole questionnaire at each visit

Variable	Mean score (range) [median; SD] (no. patients)		
	Baseline	Change from baseline to 12 mths	Change from baseline to 24 mths
Q-TFA*			
Prosthetic use score		↑	↑
Prosthetic mobility score		↑	↑
Problem score†			
Global score			
Overall situation (n, %)			
Extremely poor			
Poor			
Average			
Good			
Extremely good			

Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Sitting



<http://osseointeg.ning.com/profile/ErikAx>



<http://www.sahlgrenska.se/su/osseointegration>

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Image

Body representation

N=13

*“ The prosthesis (OI-prosthesis) is a part of me since it works so well, and you don’t have to think that it’s a problem and that it should be hard and so forth . . . **it’s more like a substitute, my “pretend leg”** ”*



<http://news.bme.com/tag/amputation/>

Lundberg, M., K. Hagberg, and J. Bullington, My prosthesis as a part of me: a qualitative analysis of living with an osseointegrated prosthetic limb. *Prosthetics and Orthotics International*, 2011. 35(2): p. 207-214

Benefits

HRQoL


SF-36

Q-TFA

Sitting

Image

Olperception



Hagberg, K., E. Häggström, S. Jönsson, B. Rydevik, and R. Brånemark, Osseoperception and Osseointegrated Prosthetic Limbs, P. Gallagher, D. Desmond, and M. MacLachlan, Editors. 2008, Springer London. p. 131-140

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Image

Olperception

Attachment

Doning and doffing



Hagberg, K., E. Häggström, S. Jönsson, B. Rydevik, and R. Brånemark, Osseoperception and Osseointegrated Prosthetic Limbs, P. Gallagher, D. Desmond, and M. MacLachlan, Editors. 2008, Springer London. p. 131-140

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Image

Olperception


Attachment

Hip ROM

Hip range of movement



<http://osseointeg.ning.com/profile/ErikAx>



Tranberg, R., R. Zügner, and J. Kärrholm, Improvements in hip- and pelvic motion for patients with osseointegrated trans-femoral prostheses. Gait & Posture, 2011. 33(2): p. 165-168

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Image

Olperception


Attachment

Hip ROM

Gait

Walking abilities

N=12



Frossard, L., K. Hagberg, E. Haggstrom, D. Lee Gow, R. Branemark, and M. Percy, Functional outcome of transfemoral amputees fitted with an osseointegrated fixation: Temporal gait characteristics. Journal of Prosthetics and Orthotics, 2010. 22(1): p. 11-20

Benefits

HRQoL

SF-36

Q-TFA

Sitting

Image

Olperception

Attachment

Hip ROM

Gait

Number of steps




Source: Facebook accessed 21/11/2015

Benefits

- HRQoL
- SF-36
- Q-TFA
- Sitting
- Image
- Olperception
- Attachment
- Hip ROM
- Gait

Activities of daily living




Frossard L, Stevenson N, Smeathers J, Häggström E, Hagberg K, Sullivan J, Ewins D, Lee Gow D, Gray S, Brånemark R. Monitoring of the load regime applied on the osseointegrated fixation of a trans-femoral amputee: A tool for evidence-based practice. Prosthetics and Orthotics International. 2008. 32 (1).

Safety

Safety

Infections

Overview – Superficial infections



[1] Al Muderis, M., K. Tetsworth, A. Khemka, S. Wilmot, B. Bosley, S.J. Lord, and V. Glatt, The Osseointegration Group of Australia Accelerated Protocol (OGAAP-1) for two-stage osseointegrated reconstruction of amputated limbs. Bone & Joint Journal, 2016. 98-B(7): p. 952-960.

Safety

Infections

Overview – Superficial infections

	Inclusion	Follow-up
		?
Reference		[1]
Number of participants in study		86
Low-grade soft-tissue infection = Cellulitis with signs of inflammation (redness, swelling, warmth, stinging pain, pain that increases on loading, tense)		28%

[1] Al Muderis, M., A. Khemka, S. Lord, H. Van de Meent, and J. Frolke, Safety of Osseointegrated Implants for Transfemoral Amputees: A Two-Center Prospective Cohort Study. The Journal of Bone & Joint Surgery, 2016. 98(11): p. 900-909.

Safety

Infections

Overview – Superficial infections

	Inclusion		Follow-up	
			(2-3 yrs)	(S2-2 yrs)
Reference	[1]	[2]	[1]	[2]
Number of participants in study	39	51	39	51
Local soft tissue infection in the skin penetration area / Superficial infection	17%	11%	29%	80%

[1] Tillander, J., K. Hagberg, L. Hagberg, and R. Branemark, Osseointegrated Titanium Implants for Limb Prostheses Attachments: Infectious Complications. Clinical Orthopaedic Related Research, 2010. 468(10): p. 2781-2788

[2] Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Safety

Infections

Overview – Superficial infections

	Inclusion		Follow-up	
			(2-3 yrs)	(S2-2 yrs)
Reference	[1]	[2]	[1]	[2]
Number of participants in study	39	51	39	51
Local soft tissue infection in the skin penetration area / Superficial infection	17%	11%	29%	80%

Cleaning

100%

[1] Tillander, J., K. Hagberg, L. Hagberg, and R. Branemark, Osseointegrated Titanium Implants for Limb Prostheses Attachments: Infectious Complications. Clinical Orthopaedic Related Research, 2010. 468(10): p. 2781-2788

[2] Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Safety

Infections

Overview - Deep infections



[1] Al Muderis, M., K. Tetsworth, A. Khemka, S. Wilmot, B. Bosley, S.J. Lord, and V. Glatt, The Osseointegration Group of Australia Accelerated Protocol (OGAAP-1) for two-stage osseointegrated reconstruction of amputated limbs. Bone & Joint Journal, 2016. 98-B(7): p. 952-960.

Safety

Infections

Overview – Superficial infections

	Inclusion	Follow-up
		?
Reference		[1]
Number of participants in study		86
High-grade soft-tissue infection = Pus collection, purulent discharge, raised level of C-reactive protein		5%

[1] Al Muderis, M., A. Khemka, S. Lord, H. Van de Meent, and J. Frolke, Safety of Osseointegrated Implants for Transfemoral Amputees: A Two-Center Prospective Cohort Study. The Journal of Bone & Joint Surgery, 2016. 98(11): p. 900-909.

Safety

Infections

Overview - Deep infections

	Inclusion		Follow-up	
			(2-3 yrs)	(S2-2 yrs)
Reference	[1]	[2]	[1]	[2]
Number of participants in study	39	51	39	51
Definite implant infection / Deep implant infection	5%	11%	15%	6%

[1] Tillander, J., K. Hagberg, L. Hagberg, and R. Branemark, Osseointegrated Titanium Implants for Limb Prostheses Attachments: Infectious Complications. Clinical Orthopaedic Related Research, 2010. 468(10): p. 2781-2788

[2] Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Safety

Infections

Overview - Deep infections

	Inclusion		Follow-up	
			(2-3 yrs)	(S2-2 yrs)
Reference	[1]	[2]	[1]	[2]
Number of participants in study	39	51	39	51
Definite implant infection / Deep implant infection	5%	11%	15%	6%

↓

Short course of antibiotics

[1] Tillander, J., K. Hagberg, L. Hagberg, and R. Branemark, Osseointegrated Titanium Implants for Limb Prostheses Attachments: Infectious Complications. Clinical Orthopaedic Related Research, 2010. 468(10): p. 2781-2788

[2] Branemark, R., O. Berlin, K. Hagberg, P. Bergh, B. Gunterberg, and B. Rydevik, A novel osseointegrated percutaneous prosthetic system for the treatment of patients with transfemoral amputation: A prospective study of 51 patients. Bone Joint J, 2014. 96(1): p. 106-113.

Safety


Infections

Breakage

High-impact activities / Falls



Source: Facebook accessed 12/04/2016




Source: Facebook accessed 21/11/2015

Safety

Infections

Breakage

High-impact activities / Falls



Thompson M. Mechanical analysis of osseointegrated transfemoral implant systems. 2009. Master Thesis. Queen's University Kingston, Ontario, Canada

Safety


Infections

Breakage

Loosening

Titel RSA and radiographic

N=55



Strong bonding

Nebergall, A., C. Bragdon, A. Antonellis, J. Kärrholm, R. Brånemark, and H. Malchau, Stable fixation of an osseointegrated implant system for above-the-knee amputees. Acta Orthopaedica, 2012. 83(2): p. 121-128

Safety

Infections

Breakage


Loosening

Fractures

XXX

XXX

Overview




Bone-anchored prostheses from
rehabilitation and **beyond**: is what you see
is what you get?

Future developments

Future developments

Fixation

ITAP, Stanmore Implant, UK



Kang, N.V., C. Pendegrass, L. Marks, and G. Blunn, Osseocutaneous integration of an intraosseous transcutaneous amputation prosthesis implant used for reconstruction of a transhumeral amputee: Case report. The Journal of Hand Surgery, 2010. 35(7): p. 1130-1134.

Future developments

Fixation

University of Utah - Orthopaedics

Orthotics/Prosthetics

Researcher announces plans for FDA study of osseointegrated implants

September 24, 2013

ORLANDO, Fla. — A Food and Drug Administration early feasibility study is expected to begin next year to test osseointegrated implants on humans, according to the presenter of the Keynote Address here at the O&P World Congress.

Roy Bloebaum, PhD, research professor at the University of Utah School of Medicine and co-director of the Department of Veterans Affairs Bone and Joint Research Lab in Salt Lake City, said a percutaneous, osseointegrated implant could avoid socket complications, such as pain, discomfort, skin breakdown, pressure sores, phantom pain and muscle weakness. The study aims to improve mobility, comfort and the ability to quickly don and doff a prostheses.

Bloebaum and his colleagues have successfully implanted osseointegrated devices in sheep models, and they expect to begin implanting the devices on human participants in April 2014.


"We want to inhibit periprosthetic infections, and we want to get successful load bearing," he said.

The study, collaboration between the University of Utah, Walter Reed National Military Medical Center and Brooke Army Medical Center, will include 10 transfemoral amputee participants. The evaluation process will include psychological and social prescreening tests and pre-implant and post-implant assessments, as well as a two-stage implementation procedure. Patients will then undergo follow-up testing for a minimum of 12 months, and two independent review boards will monitor the study.



Roy
Bloebaum


<http://www.healio.com/orthotics-prosthetics/prosthetics/news/online/%7Bbf5a0e16-eb8c-4e89-aa8b-0e2941bc31fb%7D/researcher-announces-plans-for-fda-study-of-osseointegrated-implants>



Future developments

Fixation


Tufts University



Pitkin, M., C. Cassidy, R. Muppavarapu, R. J. M. Shevtsov, O. Galibin, and S. Rousselle, New method of fixation of in-bone implanted prosthesis. J Rehabil Res Dev, 2013. 50(5): p. 709-722.

Future developments


Fixation



Juhnke, D.-L., J.P. Beck, S. Jeyapalina, and H.H. Aschoff, Fifteen years of experience with Integral-Leg-Prosthesis: Cohort study of artificial limb attachment system. Journal of Rehabilitation Research and Development, 2015. 52(4): p. 407-420.

Future developments


Fixation



Frossard L. Evaluation framework to assess benefits and harms of bone-anchored prosthesis. OSOS. 2015

Future developments


Fixation



Frossard L. Evaluation framework to assess benefits and harms of bone-anchored prosthesis. OSOS. 2015

Future developments


Fixation



Frossard L. Evaluation framework to assess benefits and harms of bone-anchored prosthesis. OSOS. 2015

Future developments

Fixation





Frossard L. Evaluation framework to assess benefits and harms of bone-anchored prosthesis. OSOS. 2015

Future developments

Fixation

Focus

Electrical field for osseointegration




Isaacson, B.M., J.G. Stinstra, R.D. Bloebaum, P.F. Pasquina, and R.S. MacLeod, Establishing multiscale models for simulating whole limb estimates of electric fields for osseointegrated implants. IEEE Trans Biomed Eng, 2011. 58(10): p. 2991-4.

Future developments

Fixation

Focus

Surgical methods




Khemka A, Frossard L, Lord S, Bosley B, Al Muderis M. Transcutaneous bone-anchoring prosthesis with total knee replacement. Acta Orthopædica. 2015. In press.

Future developments

Fixation

Focus

Neuromuscular control of prosthesis



<http://www.chalmers.se/en/news/Pages/Thought-controlled-prosthesis-is-changing-the-lives-of-amputees.aspx>

Future developments

Fixation

Focus

Challenges

Accessible to population with diabetes

ORIGINAL ARTICLE

Estimating the Prevalence of Limb Loss in the United States: 2005 to 2050

Kathryn Ziegler-Graham, PhD, Ellen J. MacKenzie, PhD, Patti L. Ephraim, MPH, Thomas G. Trivison, PhD, Ron Brookmeyer, PhD

Conclusions: One in 190 Americans is currently living with the loss of a limb. Unchecked, this number may double by the year 2050.

K. Ziegler-Graham, E. J. MacKenzie, P. L. Ephraim, T. G. Trivison, and R. Brookmeyer, "Estimating the prevalence of limb loss in the United States: 2005 to 2050," *Arch Phys Med Rehabil*, vol. 89, pp. 422-9, Mar 2008.

Future developments

Fixation

Focus

Challenges

Primary intervention



<http://www.smh.com.au/news/national/lplate-collision-kills-policewoman/2008/04/20/1208629736691.html>

Future developments

Fixation

Focus

Challenges

High impact activities



<http://www.tulsaworld.com/>

Future developments

Fixation

Focus

Challenges

Pediatric applications



26 Quantum Leap: Transitioning from Pediatric to Adult Care

By Miki Fairley

<http://online.publicationprinters.com/launch.aspx?eid=202f0de3-6dd9-4208-8677-4c273a37b5d3>

Future developments

Fixation

Focus

Challenges

Accessible after natural disasters



<http://www.msf.org.uk/teaching-resources-level-geography>

Future developments

Fixation

Focus

Challenges

Accessible to low income countries



http://projecthopeinthefield.blogspot.ca/2010_04_01_archive.html

Conclusions



Bone-anchored prostheses from rehabilitation and beyond: is what you see is what you get?



Adj/Prof Laurent Frossard

**Queensland University of Technology,
Australia**

**University of the Sunshine Coast,
Australia**

www.LaurentFrossard.com

**Rehabilitation Medicine Society of
Australia and New Zealand**

1st Annual Scientific Meeting

19 October 2016

Crown Promenade, Melbourne



Laurent Frossard

Researcher | Leader | Educator | Entrepreneur

Website

www.laurentfrossard.com

www.YourResearchProject.com



LinkedIn

www.ca.linkedin.com/pub/laurent-frossard/5/4b4/b59/

Google+

www.plus.google.com/#113083134851353167716/about



Facebook

www.facebook.com/YourResearchProject



ePrint