RF Instrumentation for Same-Breath Triple Nuclear Lung MR Imaging of ¹H and Hyperpolarized ³He and ¹²⁹Xe at 1.5T

Madhwesha Rao and Jim M. Wild*

Purpose: The hyperpolarized gases ³He and ¹²⁹Xe have distinct properties and provide unique and complementary functional information from the lungs. A triple-nuclear, same-breath imaging examination of the lungs with ¹H, ³He, and ¹²⁹Xe can therefore provide exclusive functional information from the gas images. In addition, the ¹H images provide complementary co-registered structural information in the same physiological time frame. The goal of this study was to design an RF system for triple nuclear lung MRI at 1.5T, consisting of a dual-tuned transceiver coil for ³He and ¹²⁹Xe, RF switches and a nested ¹H receiver array.

Methods: A dual-tuned transmit-receive dual-Helmholtz RF coil for ³He and ¹²⁹Xe was designed and constructed to work in unison with a nested ¹H receiver array.

Results: Triple-nuclear imaging (structural and ventilation) and apparent diffusion coefficient mapping of the human lungs was performed in the same breath-hold using the integrated RF system. B₁ maps and volumetric ventilation imaging using a three-dimensional, balanced steady-state free precession pulse sequence performed with both hyperpolarized ³He and ¹²⁹Xe indicate good stand-alone performance of the coil for the respective nucleus.

Conclusion: Triple-nuclear same-breath lung imaging with a dualtuned coil (³He and ¹²⁹Xe) and a nested ¹H array has been demonstrated with a custom RF system. **Magn Reson Med 75:1841– 1848, 2016.** © The Authors Magnetic Resonance in Medicine published by Wiley Periodicals, Inc. on behalf of International Society for Magnetic Resonance in Medicine.

Key words: dual-tuned RF coil; triple nuclear imaging; pulmonary; lungs; hyperpolarized gas

Academic Unit of Radiology, University of Sheffield, Sheffield, United Kingdom

Grant sponsor: Engineering and Physical Sciences Research Council, National Institute for Health Research (NIHR) Research Professorship, and University of Sheffield PhD Scholarship.

*Correspondence to: Jim M. Wild, Ph.D., Unit of Academic Radiology, C Floor, Royal Hallamshire Hospital, Glossop Road, Sheffield, S10 2JF, United Kingdom. E-mail: j.m.wild@sheffield.ac.uk

The copyright line for this article was changed on 12 May 2017 after original online publication.

The views expressed in this publication are those of the authors and not necessarily those of the National Health Service, the NIHR, or the UK Department of Health.

Received 23 December 2014; revised 10 February 2015; accepted 11 February 2015

DOI 10.1002/mrm.25680

INTRODUCTION

Imaging the lungs with inhaled hyperpolarized gases ³He and ¹²⁹Xe has been shown to provide functional information that cannot be accessed with proton (¹H) MRI or other imaging modalities (1-6). The two gases have distinct physical properties, which provide different but complementary functional information (7-9). The ability to image both nuclei in the same breath alongside the ¹H anatomical images adds further structural and functional sensitivity to the acquisition. ³He is highly diffusive when compared with ¹²⁹Xe (10–12), and the visualization and quantification of lung ventilation and diffusion with these two gases at the same time can help address important physiological questions such as the position of the diffusion-convection front in the lungs. The capability to measure the diffusivity of both ³He and ¹²⁹Xe gases in the same lung inflation level also provides added information for measuring and modeling lung microstructure based on their measured apparent diffusion coefficients (ADC) (13). ¹²⁹Xe is also denser and more viscous than ³He and as such has different fluid dynamic properties that define airflow in the airways, which can be measured with phase contrast MRI (14). ¹²⁹Xe has the added feature that it is soluble in blood and has a wide range of chemical shift, which enables quantification of perfusion and gas exchange in the lungs (8,9,15–17). In recent years, MRI of perfluorinated ¹⁹F gases has also gained interest (18-20) as another MR-sensitive gaseous tracer of regional lung function.

Therefore, same-breath, multinuclear lung imaging with ³He-¹²⁹Xe mixtures and ¹H MRI provides a unique combination of functional and structural information that is spatio-temporally coregistered in the same physiological time frame (21,22). Preliminary studies have used separate and spatially nested transmit-receive (T-R) coils for each nucleus (22). The reliance on the ¹H MR system's birdcage body coil for signal reception constrains the signal-tonoise ratio (SNR) in ¹H images of the lung, which is already limited by the low proton density of lung parenchyma. In a recent study (23), we showed that the ¹H lung SNR in same-breath imaging can be improved with a nested ¹H receive array, which is compatible with operation with either a ³He or a ¹²⁹Xe T-R coil.

The motivation of this study was the design and construction of an integrated radiofrequency (RF) coil and T-R switching system for triple nuclear lung imaging in the same breath. To achieve this, we developed a new dual-tuned flexible T-R RF coil to operate in quadrature for both ³He and ¹²⁹Xe at 1.5T. For ¹H imaging, we incorporated the ¹H array developed in our previous study

Published online 13 May 2015 in Wiley Online Library (wileyonlinelibrary.com). © 2016 The Authors Magnetic Resonance in Medicine published by Wiley Periodicals, Inc. on behalf of International Society for Magnetic Resonance in Medicine. This is an open access article under the terms of the Creative Commons Attribution License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited.



FIG. 1. **a**: Schematic of the dual-tuned flexible T-R coil for ¹²⁹Xe and ³He. The in-phase and quadrature-phase ports are marked with 0° and 90°, respectively. **b**: Picture of the dual-tuned flexible T-R coil for ¹²⁹Xe and ³He. The in-phase and quadrature-phase ports are marked with 0° and 90°, respectively. The plastic housing at the port (0°, 90°) consists of high pass matching circuits (90 pF, 90 pF, and 300 nH) marked in the schematic and the two loop capacitors (47 pF and 220 pF). **c**: Illustration of application of dual-tuned T-R coil on the subject for same-breath ADC measurement. **d**: Illustration of ¹H array and dual-tuned coil nested for triple nuclear imaging. **e**: Picture of the setup on the scanner. The picture indicates the mouthpiece and the two Tedlar bags affixed to it.

designed to nest within either ³He or ¹²⁹Xe T-R RF coils (23). With the developed RF instrumentation, we demonstrated triple nuclear same-breath lung imaging with hyperpolarized ³He and ¹²⁹Xe ventilation images and ¹H anatomical images. With the same system ADC measurement of mixtures of ³He and ¹²⁹Xe were performed in the same breath at a particular lung inflation state.

METHODS

³He and ¹²⁹Xe Dual-Tuned Coil Design

A dual-tuned (³He-¹²⁹Xe) flexible quadrature T-R coil was constructed in-house. The conducting elements were made from self-adhesive copper tape (FE-5100-5276-7; 3M, Bracknell, UK) of 66-µm thickness and 6-mm width, which was fixed on a substrate of 0.5-mm-thick polytetrafluoroethylene (Direct Plastics, Sheffield, UK) as shown in Figure 1b. The capacitors used on the resonant circuit were of 10C package (Dalian Dalicap Technology Co., Ltd, Dalian, China). The thickness of the array with the foam was 6 mm (3 mm each side). The dual-tuned flexible T-R coil was a dual Helmholtz-like pair of quadrature design in which the Helmholtz for the in-phase resonance of the quadrature spans the anterior right lung to posterior left lung, connected over left trapezius. Similarly, the Helmholtz for the quadrature-phase resonance spans the anterior left lung to posterior right lung, via the right trapezius. The cross-over of the copper strip for each of the Helmholtz pairs (which forms a "figure eight" topology) was positioned such that it was within the other resonant element (anterior) and was balanced on either side to minimize coupling as shown in Figure 1a and 1b. The schematic of the dual-tuned T-R coil circuit is shown in Figure 1a, and a photograph is shown in Figure 1b. The assembled topology of the flexible coil constitutes a bib design wrapped around the subject longitudinally, as shown in Figure 1c. Both the elements of the dual Helmholtz were fitted with two traps; one trap at the ¹H frequency to enable ¹H imaging with this coil in situ and the other trap to dual-tune the coil to the ¹²⁹Xe and ³He Larmor frequencies. The trap design was based on the formalism established in our earlier study for multituned resonators (23), the frequency of the trap for dual-tuning was 47.81 MHz. A high-pass matching circuit was used to match the coil at both resonant frequencies of ${}^{3}\mathrm{He}$ (48.62 MHz) and ${}^{129}\mathrm{Xe}$ (17.65 MHz) at 1.5T. The ¹H trap was tuned with a 47-pF capacitor and a seven-turn wire wound inductor with a diameter of 6 mm. The trap for dual-tuning the coil was tuned with a 56-pF capacitor and a nine-turn wire wound inductor with a diameter of 6 mm. Wire wound inductors were constructed from 21 AWG insulated copper wire. RF measurements were performed with an Agilent 5061B Network Analyzer (Keysight Technologies, Santa Rosa, California, USA). For the RF measurement, the dual-tuned coil was wrapped longitudinally around the thorax of the subject, as shown in Figure 1c. This coil was designed to work with full functionality when the four-channel ¹H chest receiver array from our earlier study (23) was nested inside for in situ high SNR ¹H lung imaging.

MR Imaging Methods for ³He, ¹²⁹Xe and ¹H

All in vivo imaging with ³He and ¹²⁹Xe was performed with approval from the National Research Ethics Committee. The imaging was performed on a healthy male volunteer (age, 31 years; height, 185 cm; weight, 89 kg). Lung MRI was performed on a GE whole body 1.5T Signa HDx system with ³He and ¹²⁹Xe gas polarized with spin exchange optical pumping (24). The gas dosage and the imaging and pulse sequence parameters used for all three nuclei are shown in Table 1. The hyperpolarized Table 1

Gas	Mixture Do	osade. I	maaina	Parameters.	and Pulse	Sequence	Parameters	Used in the	Study

	Lung structure and			
Measurement	ventilation	ADC	Whole lung ventilation	Flip angle map
Physiological details	Triple-nuclear same breath	Dual-nuclear same breath	Single-nuclear separate	Single-nuclear separate
Thysiological dotails	mple hubbell, earlie breath	Buar nuclear, came broath	breath	breath
PE coil	¹ H array dual tupod coil	Dual tupod coil	Dual tupod coil	Dual tupod coil
Nuclei				
Nuclei	H, He, Ke	-He,Xe	-He,Xe	-He,Xe
Dosage (mL)				
'H	—	—	—	—
³ He	350	300	200	50
¹²⁹ Xe	500	500	500	100
Flip angle				
1 ¹ H	50°	_	_	_
3 _{Но}	80	Ô٥	10°	_
129	0	9 10°	10	—
	9	10*	10	—
1E (MS)				
'H	0.9			
°He	1.1	4.8	0.6	1.1
¹²⁹ Xe	3.6	12.5	2.1	3.6
TR (ms)				
¹ H	2.9	_	_	_
³ He	3.6	10	1.9	3.6
¹²⁹ Xe	18.9	27	6.4	18 9
Matrix	10.0		0.4	10.0
	100			
Phase	192	—	—	—
Frequency	256	—	—	—
³ He				
Phase	104	48	82	52
Frequency	80	64	80	44
¹²⁹ Xe				
Phase	78	48	82	52
Frequency	64	64	80	44
Slice thickness (mm)	04	04	88	44
	15			
H 3	15	_	_	_
	15	15	4	200
¹²⁹ Xe	15	15	10	200
Number of slices				
¹ H	3			
³ He	3	2	46	1
¹²⁹ Xe	3	2	24	1
Field of view (cm)				
¹ н	40			
3 ₄ 0	40	4.4	40	40
129	40	44	40	40
Ae	40	44	40	40
Axis	2D, coronal	2D, coronal	Coronal	2D, coronal
Pulse sequence	bSSFP, FSGRE	FSGRE	3D bSSFP	FSGRE
Imaging time (s)				
¹ H	1	_	_	_
³ He	2	6	7	0.9
¹²⁹ Xe	4	8	13	2.4
Multiphase				
¹ н	_	_	_	_
311a	_	—	—	
129	—	—	—	6
	—	—	—	6
b Value (s · cm ⁻²)				
'H	—	—	—	—
³ He	_	1.6	_	_
¹²⁹ Xe	—	8	—	—
Corresponding figure				
1H	3a, 3d, 3e	_	_	_
³ He	3b. 3d	2e	4a	2c
¹²⁹ Xe	3c. 3e	2f	4b	2d
	,		-	

Abbreviations: ADC, apparent diffusion coefficient; bSSFP, balanced steady-state free precession; FSGRE, fast spoiled gradient echo.



FIG. 2. **a:** Matching of dual-tuned coil at 17.65 MHz and 48.62 MHz. **b:** Isolation between dual-tuned coil and ¹H array over frequency span of 10–80 MHz. **c:** Flip angle map of dual-tuned coil for ³He. **d:** Flip angle map of dual-tuned coil for ¹²⁹Xe. The color bars indicate the flip angle in degrees. **e, f:** ADC measurement was taken in the same breath for ³He (e) and ¹²⁹Xe (f). The color bars indicate the ADC in cm² · s⁻¹.

³He and ¹²⁹Xe gas was delivered in separate Tedlar bags and was mixed at the mouth piece at the time of inhalation, as illustrated in Figure 1e. ³He had polarization of 25% (\approx 100% of He is ³He). ¹²⁹Xe had a polarization of 40%–50% (87% of Xe is ¹²⁹Xe).

RF Signal Routing and Calibration

To route the transmit RF signal (³He-¹²⁹Xe) from the appropriate T-R switch on the scanner to the dual-tuned coil and to route the received RF signal (³He-¹²⁹Xe) from the dual-tuned coil back to the appropriate T-R switch on the scanner, a 2-kW rated coaxial antenna RF switch (CX-SW2PL; Watson, Essex, UK) was used. The RF power required for the desired flip angle for the ³He and ¹²⁹Xe sequences was calculated based on a standard calibration procedure, whereby the rate of depletion of polarization was calculated from the decay of signal resulting from a set of hard RF pulse-acquires of equal amplitude. The period to prescribe calibration values on the spectrometer between the end of imaging a particular nucleus and initiation of the sequence for imaging the next nucleus was less than 4 s. The time required to operate the RF switch manually between acquisitions was 3 s.

Same-Breath ADC (³He and ¹²⁹Xe) and Triple Nuclear (³He, ¹²⁹Xe, and ¹H) Structure and Ventilation Lung Imaging Methods

For same-breath ADC measurement, the dual-tuned ³He-¹²⁹Xe coil was wrapped longitudinally as shown in Figure 1c, without the ¹H array nested inside. To demonstrate same-breath ADC maps, two sets of ADC measurements were acquired back-to-back in a single breath, with ³He ADC measurement followed by ¹²⁹Xe measurement. The imaging parameters are shown in Table 1.

For triple-nuclear lung imaging, the dual-tuned ³He -¹²⁹Xe coil and the ¹H array from our earlier study (23) were nested as shown in Figure 1d. To demonstrate imaging of all three nuclei in the same breath, three sets of images were acquired back-to-back in a single breath in the order, with ³He imaging followed by ¹²⁹Xe imaging, in turn followed by ¹H imaging. The imaging parameters are shown in Table 1.

The T_1 of hyperpolarized gases when inhaled into the lungs is sensitive to the oxygen partial pressure in the lung during the breath-hold (25). ³He is more sensitive to this effect because the gyromagnetic ratio of ³He is approximately three times larger than that of ¹²⁹Xe, as



FIG. 3. **a**: ¹H images from lungs. **b**: Samebreath ³He images from lungs. **c**: Samebreath ¹²⁹Xe images from lungs. **d**: ³He images superimposed over ¹H images. **e**: ¹²⁹Xe images superimposed over ¹H images.

such the dipolar coupling to the electrons in the paramagnetic oxygen molecule is stronger. This rationale for the order of acquisition is ³He followed by 129 Xe, in turn followed by ¹H.

Flip Angle Mapping and High-Resolution Imaging Performance of the Coil as a Stand-Alone T-R Coil for 3 He and 129 Xe

Flip angle maps of the dual-tuned coil at the ³He and ¹²⁹Xe frequencies were calculated by measuring the depletion of polarization of the hyperpolarized gas ³He and ¹²⁹Xe at each voxel in the lungs by repeated imaging at breath-hold with a two-dimensional spoiled gradient echo sequence. The imaging parameters for this measurement are shown in Table 1, and T₁ relaxation was neglected when calculating the flip angle. In addition, to demonstrate the coil's performance as a stand-alone ³He or ¹²⁹Xe T-R coil (without the ¹H array in situ), high-resolution, three-dimensional (3D) imaging data sets were acquired with a 3D balanced steady state sequence (26) with imaging parameters as shown in Table 1.

RESULTS

Dual-Tuned Coil RF Performance

The two traps on the coil at 47.81 MHz and 63.86 MHz (¹H trap) generated three resonant modes at 17.65 MHz (¹²⁹Xe Larmor frequency), 48.62 MHz (³He Larmor fre-

quency), and 79.2 MHz. The isolation between the two ports of the Helmholtz was less than -15 dB. The guality (Q) factor of the dual-tuned coil at the ¹²⁹Xe Larmor frequency (17.65 MHz) was 61 in the unloaded condition and 17 in the loaded condition. The Q factor of the dualtuned coil at the ³He Larmor frequency (48.62 MHz) was 32 in the unloaded condition and 7 in the loaded condition. Thus, the ratio of Q factor unloaded to loaded condition was 3.5 at the ¹²⁹Xe Larmor frequency (17.65 MHz) and 4.5 at the ³He Larmor frequency (48.62 MHz). Under the loaded condition, the dual-tuned coil was matched to less than -20 dB at both ports at both the ¹²⁹Xe (17.65 MHz) and ³He (48.62 MHz) Larmor frequencies, as shown in Figure 2a. The isolation between the dual-tuned coil and the ¹H array was less than -15 dB, as shown in Figure 2b. Flip angle maps from the dualtuned coil for the transmit RF power prescribed for the nominal flip angles used for triple nuclear same-breath imaging and ADC measurement (Table 1) are shown in Figure 2c for ³He and Figure 2d for ¹²⁹Xe. The standard deviation of the flip angle map was calculated to be 0.7° (mean = 8°) for ³He and 0.3° (mean = 9°) for ¹²⁹Xe.

Multinuclear Lung Imaging

Same-breath ADC measurement of ³He and ¹²⁹Xe performed in the same lung-inflation state is shown in Figure 2e and 2f. The ³He ADC map shown in Figure 2e



FIG. 4. Hyperpolarized gas images of the lungs from 3D balanced steady-state free precession sequence with dual-tuned coil in a stand-alone configuration (without ¹H array nested). **a:** Hyperpolarized ³He gas. **b:** Hyperpolarized ¹²⁹Xe gas.

and the $^{129}\mathrm{Xe}$ ADC map shown in Figure 2f were acquired in the same breath.

Same-breath triple nuclear lung (structure and ventilation) images are shown in Figure 3. The ¹H images shown in Figure 3a, ³He images shown in Figure 3b, and ¹²⁹Xe images shown in Figure 3c, all of which were acquired in the same breath, are coregistered as shown in superimposed images in Figure 3d and 3e. Volumetric ventilation images from the 3D balanced steady-state free precession sequence for the coil in operation as a stand-alone transceiver for ³He and ¹²⁹Xe are shown in Figure 4a and Figure 4b, respectively.

DISCUSSION

The construction of the flexible dual-tuned coil is in the form of a bib, which enables a close fit to the subject's thorax irrespective of body type. The design was optimized to the typical subject size mentioned earlier. As the shape/form deviates from the optimal design with other body types, the distributed inductance and T-R efficiency of the dual-tuned coil changes accordingly. The B₁ field homogeneity of the dual Helmholtz design is inherently inferior to that of a birdcage design (27), and the flexibility of the dual-tuned coil adds some variability in this respect. Considering the typical anatomy of a torso, the distance between the RF coil and lung air spaces generally increases from superior (upper) to inferior (lower). This means that sensitivity in the lower lung is reduced for two reasons: first, due to proximity of the conducting elements to the lungs, and second, as the parallel condition for a Helmholtz pair is disrupted. Despite these factors, the observed B_1 transmit homogeneity (variation in flip angle, 3% for ^{129}Xe and 9% for ³He) is comparable to (28) or better than (29) studies reported previously using single-tuned flexible T-R coils for ³He and ¹²⁹Xe lung imaging.

Because the RF switches are currently manually operated, and the spectrometer has an inherent delay time for precalibration for each nucleus, the method is not currently compatible with repetition time resonant frequency interleaved imaging, as demonstrated in our earlier study with same-breath ³He-¹H lung imaging (22). It should be noted that this limitation is not due to the RF coil design or configuration; instead, it is due to the MR system, which supports only one spectrometer T-R switch (single-nucleus) to be actively connected at any given point in time (in addition to ¹H). Both the dualtuned RF coil and the nested ¹H array from the earlier study (23) are capable of operating simultaneously. If we consider the coil's operation as part of the system for triple nuclear imaging, 50%-60% of the time (ie, 18-20 s of the breath-hold) is consumed by switching the spectrometer between the nuclei. This can be reduced with the appropriate spectrometer software engineering and using electrically driven RF switches (eg, PIN diodes and Field Effect Transistor).

The free diffusion (in air) of ³He is 0.88 cm² · s⁻¹ (30,31); in this study, we report 0.85 cm² · s⁻¹ for ³He in the trachea (slightly lower than ³He free diffusion). The free diffusion of ¹²⁹Xe is 0.14 cm² · s⁻¹ (31); in this study, we report 0.22 cm² · s⁻¹ for ¹²⁹Xe in the trachea. The higher ADC value for ¹²⁹Xe in the trachea, we presume is due to its mixture with the highly diffusive ³He (as shown in Table 1). In the ventilation images, any observed asymmetry beyond what can be attributed to the measured variation/asymmetry in the flip angle was verified to be caused by the distribution of ³He and ¹²⁹Xe as a gas mixture in the lung (variation in the local concentration of the gases). These findings are currently being investigated in future work studying the physiology of gas mixing in the lung with the two gases.

In contrast to our previous triple nuclear same-breath lung imaging experiments demonstrated at 3T on a Philips system using the ¹H body T-R coil, a ³He birdcage T-R coil, and a nested ¹²⁹Xe T-R vest coil (22), the design used in this study at 1.5T has several potential benefits. First, from the coil perspective, the use of the dual-tuned ³He-¹²⁹Xe coil minimizes the number of individually tuned coils, and the nested ¹H array (23) improves the ¹H SNR by closer proximity to the lung. Applications of this triple nuclear RF system for lung MRI are manifold and allow the different physical and physiological properties of the two gases to be explored in the same time course with added provision of high-quality and coregistered ¹H structural images.

In conclusion, we have demonstrated a system for triple nuclear same-breath lung imaging of ¹H with hyperpolarized gases ³He and ¹²⁹Xe at 1.5T using a custom integrated RF system. This system incorporates a new design of dual-tuned RF coil for ³He and ¹²⁹Xe and RF switches, together with a nested receiver array for ¹H imaging. With this system, we have demonstrated high-quality, same-breath ¹H with ³He and ¹²⁹Xe ventilation imaging and the capability for ADC mapping of ³He and ¹²⁹Xe in the same lung-inflation state. In addition, the image quality on all three nuclei is comparable with those acquired with separate RF coils for the given nucleus.

ACKNOWLEDGMENTS

We thank Juan Parra-Robles, Neil Stewart, Helen Marshall, and Guilhem Collier for discussion and assistance with the imaging experiments.

REFERENCES

- Donnelly LF, MacFall JR, McAdams HP, Majure JM, Smith J, Frush DP, Bogonad P, Charles HC, Ravin CE. Cystic fibrosis: combined hyperpolarized 3He-enhanced and conventional proton MR imaging in the lung—preliminary observations. Radiology 1999;212:885–889.
- Woodhouse N, Wild JM, Paley MNJ, Fichele S, Said Z, Swift AJ, Van Beek EJR. Combined helium-3/proton magnetic resonance imaging measurement of ventilated lung volumes in smokers compared to never-smokers. J Magn Reson Imaging 2005;21:365–369.
- Zheng J, Leawoods JC, Nolte M, Yablonskiy DA, Woodard PK, Laub G, Gropler RJ, Conradi MS. Combined MR proton lung perfusion/ angiography and helium ventilation: potential for detecting pulmonary emboli and ventilation defects. Magn Reson Med 2002;47:433– 438.
- Ebert M, Grossmann T, Heil W, Otten WE, Surkau R, Leduc M, Bachert P, Knopp MV, Schad LR, Thelen M. Nuclear magnetic resonance imaging with hyperpolarised helium-3. Lancet 1996;347:1297– 1299.
- Van Beek EJR, Wild JM, Kauczor HU, Schreiber W, Mugler Iii JP, De Lange EE. Functional MRI of the lung using hyperpolarized 3-helium gas. J Magn Reson Imaging 2004;20:540–554.
- Swift AJ, Wild JM, Fichele S, Woodhouse N, Fleming S, Waterhouse J, Lawson RA, Paley MNJ, Van Beek EJR. Emphysematous changes and normal variation in smokers and COPD patients using diffusion He-3 MRI. Eur J Radiol 2005;54:352–358.
- Kauczor HU, Surkau R, Roberts T. MRI using hyperpolarized noble gases. Eur Radiol 1998;8:820–827.
- Möller HE, Chen XJ, Saam B, Hagspiel KD, Johnson GA, Altes TA, De Lange EE, Kauczor HU. MRI of the lungs using hyperpolarized noble gases. Magn Reson Med 2002;47:1029–1051.
- Fain SB, Korosec FR, Holmes JH, O'Halloran R, Sorkness RL, Grist TM. Functional lung imaging using hyperpolarized gas MRI. J Magn Reson Imaging 2007;25:910–923.
- Kirby M, Svenningsen S, Owrangi A, et al. Hyperpolarized He-3 and Xe-129 MR imaging in healthy volunteers and patients with chronic obstructive pulmonary disease. Radiology 2012;265:600–610.
- 11. Kirby M, Svenningsen S, Kanhere N, Owrangi A, Wheatley A, Coxson HO, Santyr GE, Paterson NAM, McCormack DG, Parraga G. Pulmonary ventilation visualized using hyperpolarized helium-3 and xenon-129 magnetic resonance imaging: differences in COPD and relationship to emphysema. J Appl Physiol 2013;114:707–715.

- Ouriadov A, Farag A, Kirby M, McCormack DG, Parraga G, Santyr GE. Lung morphometry using hyperpolarized Xe-129 apparent diffusion coefficient anisotropy in chronic obstructive pulmonary disease. Magn Reson Med 2013;70:1699–1706.
- 13. Parra-Robles J, Wild JM. The influence of lung airways branching structure and diffusion time on measurements and models of short-range 3He gas MR diffusion. J Magn Reson 2012;225:102–113.
- Collier GJ, Wild JM. In vivo measurement of gas flow in human airways with hyperpolarized gas MRI and compressed sensing. Magn Reson Med 2015;73:2255–2261.
- Mugler III JP, Driehuys B, Brookeman JR, et al. MR imaging and spectroscopy using hyperpolarized 129Xe gas: preliminary human results. Magn Reson Med 1997;37:809–815.
- Wagshul ME, Button TM, Li HFF, Liang ZR, Springer CS, Zhong K, Wishnia A. In vivo MR imaging and spectroscopy using hyperpolarized Xe-129. Magn Reson Med 1996;36:183–191.
- Ruppert K, Mata JF, Brookeman JR, Hagspiel KD, Mugler JP. Exploring lung function with hyperpolarized Xe-129 nuclear magnetic resonance. Magn Reson Med 2004;51:676–687.
- Couch MJ, Ball IK, Li T, Fox MS, Littlefield SL, Biman B, Albert MS. Pulmonary ultrashort echo time 19F MR imaging with inhaled fluorinated gas mixtures in healthy volunteers: feasibility. Radiology 2013; 269:903–909.
- Couch MJ, Ball IK, Li T, Fox MS, Ouriadov AV, Biman B, Albert MS. Inert fluorinated gas MRI: a new pulmonary imaging modality. NMR Biomed 2014;27:1525–1534.
- 20. Halaweish AF, Charles HC. Physiorack: an integrated MRI safe/conditional, gas delivery, respiratory gating, and subject monitoring solution for structural and functional assessments of pulmonary function. J Magn Reson Imaging 2014;39:735–741.
- 21. Wild JM, Ajraoui S, Deppe MH, Parnell SR, Marshall H, Parra-Robles J, Ireland RH. Synchronous acquisition of hyperpolarised 3He and 1H MR images of the lungs—maximising mutual anatomical and functional information. NMR Biomed 2011;24:130–134.

- 22. Wild JM, Marshall H, Xu X, Norquay G, Parnell SR, Clemence M, Griffiths PD, Parra-Robles J. Simultaneous imaging of lung structure and function with triple-nuclear hybrid MR imaging. Radiology 2013; 267:251–255.
- Rao M, Robb F, Wild JM. Dedicated receiver array coil for 1H lung imaging with same-breath acquisition of hyperpolarized 3He and 129Xe gas. Magn Reson Med 2015;74:291–299.
- 24. Norquay G, Parnell SR, Xu X, Parra-Robles J, Wild JM. Optimized production of hyperpolarized 129Xe at 2 bars for in vivo lung magnetic resonance imaging. J Appl Phys 2013;113.
- 25. Saam B, Happer W, Middleton H. Nuclear relaxation of He3 in the presence of O2. Phys Rev A 1995;52:862–865.
- Wild JM, Teh K, Woodhouse N, Paley MN, Fichele S, de Zanche N, Kasuboski L. Steady-state free precession with hyperpolarized 3He: experiments and theory. J Magn Reson 2006;183:13–24.
- Hayes CE, Edelstein WA, Schenck JF, Mueller OM, Eash M. An efficient, highly homogeneous radiofrequency coil for whole-body NMR imaging at 1.5 T. J Magn Reson (1969) 1985;63:622–628.
- De Zanche N, Chhina N, Teh K, Randell C, Pruessmann KP, Wild JM. Asymmetric quadrature split birdcage coil for hyperpolarized 3He lung MRI at 1.5T. Magn Reson Med 2008;60:431–438.
- 29. Miller GW, Altes TA, Brookeman JR, de Lange EE, Mugler III JP. Hyperpolarized 3He lung ventilation imaging with B 1inhomogeneity correction in a single breath-hold scan. Magn Reson Mater Phy 2004;16:218–226.
- 30. Yablonskiy DA, Sukstanskii AL, Leawoods JC, Gierada DS, Bretthorst GL, Lefrak SS, Cooper JD, Conradi MS. Quantitative in vivo assessment of lung microstructure at the alveolar level with hyperpolarized 3He diffusion MRI. Proc Natl Acad Sci U S A 2002; 99:3111–3116.
- 31. Chen XJ, Möller HE, Chawla MS, Cofer GP, Driehuys B, Hedlund LW, Johnson GA, Fitzsimons EG. Spatially resolved measurements of hyperpolarized gas properties in the lung in vivo. Part I: diffusion coefficient. Magn Reson Med 1999;42:721–728.