



**Tinning, Peter William and Franssen, Aimee and Hridi, Shehla Unaiza and Bushell, Trevor and McConnell, Gail (2017) A 340/380 nm light emitting diode illuminator for Fura-2 AM ratiometric Ca<sup>2+</sup> imaging of live cells with better than 5 nM precision. In: Photonex, 2017-06-14 - 2017-06-14. ,**

This version is available at <https://strathprints.strath.ac.uk/61504/>

**Strathprints** is designed to allow users to access the research output of the University of Strathclyde. Unless otherwise explicitly stated on the manuscript, Copyright © and Moral Rights for the papers on this site are retained by the individual authors and/or other copyright owners. Please check the manuscript for details of any other licences that may have been applied. You may not engage in further distribution of the material for any profitmaking activities or any commercial gain. You may freely distribute both the url (<https://strathprints.strath.ac.uk/>) and the content of this paper for research or private study, educational, or not-for-profit purposes without prior permission or charge.

Any correspondence concerning this service should be sent to the Strathprints administrator: [strathprints@strath.ac.uk](mailto:strathprints@strath.ac.uk)

Peter W. Tinning<sup>1</sup>, Aimee J.P.M. Franssen<sup>2</sup>, Shehla U. Hridi<sup>2</sup>, Trevor J. Bushell<sup>2</sup> and Gail McConnell<sup>1</sup>

<sup>1</sup> Centre for Biophotonics, Strathclyde Institute of Pharmacy and Biomedical Sciences, University of Strathclyde, Glasgow, United Kingdom.

<sup>2</sup> Strathclyde Institute of Pharmacy and Biomedical Sciences, University of Strathclyde, Glasgow, United Kingdom.

## Introduction

Cytosolic Ca<sup>2+</sup> plays an integral role in cells and the study of its dynamics can reveal much about biological processes [1]. Fura-2 can provide quantitative data on cytosolic Ca<sup>2+</sup> changes by exciting at 340 nm and 380 nm and taking the ratio of the emission at both wavelengths [2].

Traditionally for this type of imaging an arc lamp had to be used for illumination as LEDs of the appropriate wavelengths were not available [3]. LEDs hold advantages over arc lamps by exhibiting high amplitude stability and the ability to rapidly switch between wavelengths. We aimed to test a new 340/380 nm LED system for use in ratiometric Fura-2 AM Ca<sup>2+</sup> imaging and present results using tsA-201 cells and hippocampal neurons.

## Methods

- Specimens were washed three times with HEPES-buffered saline solution and loaded with Fura-2 AM for 60 minutes at 37 °C. They were then washed a further three times before imaging.
- HBS Control solution:** - 1 litre of distilled water containing (in mM): NaCl, 140; KCl, 5; MgCl<sub>2</sub>, 2; HEPES, 10; D-glucose, 10; CaCl<sub>2</sub>, 2 at a pH of 7.4.
- Power at the specimen plane:** - 340 nm: 1.35 mW  
380 nm: 1.40 – 3.08 mW
- LED exposure time:** - 0.5 Hz imaging: 100 ms  
24.39 Hz video-rate imaging: 20.5 ms
- Emission detection:** - Hamamatsu ORCA-Flash 4.0 CMOS camera with a binning n = 2




Figure 1: Schematic diagram of Olympus BX50 microscope and imaging apparatus

## 0.5 Hz ratiometric Fura-2 AM Ca<sup>2+</sup> imaging of drug-mediated responses in tsA-201 cells and hippocampal neurons




Figure 2: A) 340 and 380 images of 1 μM Fura-2 AM loaded tsA-201 cells  
Figure 2: B) 340 and 380 images of 1 μM Fura-2 AM loaded hippocampal neurons




Figure 3: A) Representative traces of pharmacologically induced Ca<sup>2+</sup> concentration changes in tsA-201 cells and C) hippocampal neurons. B) Average Ca<sup>2+</sup> increase for each stimuli in tsA-201 cells (n = 572) and D) hippocampal neuron (n = 388). Ca<sup>2+</sup> increases are in agreement with previous experiments in the same cells illuminated by arc lamps [4 - 8]

## Average fluctuations in basal Ca<sup>2+</sup> levels




Figure 4: Example plot of baseline Ca<sup>2+</sup> fluctuations in hippocampal neurons

The fluctuations are below the theoretical 5 – 10 nM precision of Fura-2 [9]

## Ca<sup>2+</sup> imaging of tsA-201s loaded with lower concentrations of Fura-2AM




Figure 5: Ca<sup>2+</sup> responses obtained from tsA-201 cells loaded with different concentrations of Fura-2 AM

No significant change in obtained trypsin (100 nM) Ca<sup>2+</sup> response

Advantages of using lower concentrations:

- Increasing the number of uses from the vial
- Reduced cost**
- Improving cell viability
- Cells live longer**

## Fura-2 ratiometric imaging of synaptically-driven Ca<sup>2+</sup> events in hippocampal neurons




Figure 6: Synaptically driven Ca<sup>2+</sup> events captured at a rate of A) 0.5 Hz in a single neuron and B) 24.39 Hz in two neurons

## Conclusions

- This is the first application of a truly 340/380 nm LED illuminator for Fura-2 ratiometric Ca<sup>2+</sup> imaging of live cell specimens with a precision that is only limited by the response of Fura-2.
- Using this illuminator it is now possible to use Fura-2 AM dye concentrations as low as 250 nM offering both an economical and cell viability advantage.
- Video rate imaging of synaptically-driven Ca<sup>2+</sup> events combines high temporal and spatial resolution to obtain higher throughput information than previously possible [10].
- This LED illuminator combines optimum excitation and high stability wavelength switching to free Fura-2 imaging from illumination problems experienced in the past.