



**Karolinska  
Institutet**

This is the peer reviewed version of the following article  
**J Chromatogr A. 2014 Aug 15;1355:100-6.** which has  
been published in final form at:

<http://dx.doi.org/10.1016/j.chroma.2014.05.082>

**Benzo[a]pyrene-specific online high-performance liquid  
chromatography fractionation of air particulate extracts  
: a tool for evaluating biological interactions.**

**Lim, H; Bergvall, Christoffer; Jarvis, Ian WH; Mattsson,  
Å; Dreij, Kristian; Stenius, Ulla; Westerholm, Roger**

Access to the published version may require subscription.  
Published with permission from: **Elsevier**

1 Benzo[*a*]pyrene-specific online HPLC fractionation of air particulate extracts – A tool for  
2 evaluating biological interactions

3 Hwanmi Lim<sup>a</sup>, Christoffer Bergvall<sup>a</sup>, Ian Jarvis<sup>b</sup>, Åse Mattsson<sup>b</sup>, Kristian Dreij<sup>b</sup>, Ulla Stenius<sup>b</sup>, Roger Westerholm<sup>a,\*</sup>

4 <sup>a</sup> *Department of Analytical Chemistry, Stockholm University, Svante Arrhenius väg 16 SE-106 91 Stockholm,*  
5 *Sweden*

6 <sup>b</sup> *Institute of Environmental Medicine, Karolinska Institutet, Box 210, SE-171 77 Stockholm, Sweden*

7 \* Corresponding author: Department of Analytical Chemistry, Stockholm University, Svante Arrhenius väg 16 SE-  
8 106 91 Stockholm, Sweden. E-mail: roger.westerholm@anchem.su.se. Tel: +46-8-16 24 40. Fax: +46-8-15 63 91.

9

10 ABSTRACT

11 Benzo[*a*]pyrene (B[*a*]P) is a known human carcinogen and is commonly used as a surrogate for assessing the  
12 carcinogenic risk posed by complex mixtures of polycyclic aromatic hydrocarbons (PAHs) present in air particulate  
13 matter (PM). However, studies have shown that using B[*a*]P as a surrogate may underestimate the carcinogenic  
14 potential of PAH mixtures, as the risk assessment approach does not consider interaction effects. Thus, toxicological  
15 studies using B[*a*]P to assess its carcinogenic potential in environmentally derived complex mixtures, as opposed to  
16 single compound experiments, could improve risk assessment. The intention of the present study was to develop an  
17 online HPLC fractionation system for the selective removal of B[*a*]P from air PM extracts. Two serial pyrenylethyl  
18 (PYE) columns enabled selective separation of B[*a*]P from its isomers and other PAHs as well as a short  
19 fractionation cycle of 30 minutes. One run consisted of three collection steps: the first fraction contained PAHs  
20 eluting earlier than B[*a*]P, the second contained B[*a*]P and the last contained later-eluting PAHs. The selectivity and  
21 recovery of the system was investigated using extracts of Stockholm air PM samples. The overall recovery for all  
22 PAHs was approximately 80%, and the system proved to be selective, as it removed 94% of B[*a*]P and less than 3%  
23 of benzo[*b*]fluoranthene from the complex PAH mixture. Exposing human cells to blanks generated by the  
24 fractionation system did not induce cytotoxicity or DNA damage signalling. In conclusion, the online HPLC system  
25 was selective for B[*a*]P fractionation whilst minimising run-to-run variation and allowing repeated fractionations for  
26 larger samples due to its relatively short run time

27 *Keywords:* automation, benzo[*a*]pyrene, fractionation, PAH, pyrenylethyl, synergism

28

29 1. Introduction

30 Benzo[a]pyrene (B[a]P) belongs to the group of chemical compounds referred to as polycyclic aromatic  
31 hydrocarbons (PAHs) [1]. Recent studies on B[a]P have revealed its potential role as a cancer inducer in human skin  
32 breast cells, and neurotoxic effects resulting from exposure to B[a]P have been investigated in Sprague-Dawley rat  
33 [2-4].

34 B[a]P is currently the only PAH classified as carcinogenic to humans, whilst other PAHs, such as  
35 dibenz[a,h]anthracene, several dibenzopyrenes and benzo[fluoranthenes], have been classified as “probably or  
36 possibly carcinogenic to humans” based on studies using various in vitro and in vivo experimental models [1]. B[a]P  
37 is the key compound used in the two main approaches for assessing the carcinogenic risk of human exposure to  
38 PAHs in ambient air [5,6]. In the epidemiological approach, B[a]P is used as a surrogate for the entire complex  
39 mixture of PAHs together with occupational epidemiological data to determine quantitative risk estimates. The  
40 World Health Organisation has derived a unit risk value based on the increased risk for lung cancer in coke-oven  
41 workers [7,8]. However, a major weakness of this methodology is that it assumes that the relative PAH  
42 concentration profiles are stable between different exposure situations. The other method for risk assessment utilises  
43 animal experiment data on the carcinogenic potency of individual PAHs relative to that of B[a]P. These relative  
44 potency values, referred to as toxic equivalency factors (TEFs) or potency equivalency factors (PEFs), are then  
45 multiplied with the concentration of the respective PAH, creating a product called a B[a]P equivalent [5]. Summing  
46 the B[a]P equivalents and multiplying by the potency of B[a]P gives a risk estimate [5]. One major assumption  
47 using the TEF/PEF scheme is the additivity of individual risks, implying that interactions among the different PAHs  
48 in complex exposure situations do not modulate the carcinogenic potential of the PAH mixture. Chemical  
49 interactions can be additive, synergistic, potentiative or antagonistic when two or more chemicals are combined.  
50 That is, the toxicity of a PAH mixture can increase or decrease depending on the corresponding interactions, which  
51 can lead to under- or overestimations of its toxicological effect [9]. Previous studies have reported synergistic [10-  
52 15] and antagonistic effects [16-18] of complex PAH mixtures in different biological systems, and the risk of  
53 underestimating the carcinogenic potential by the use of TEFs/PEFs has been demonstrated in several in vitro and in  
54 vivo studies [10-12,14,19,20]. Moreover, a newer approach to environmental cancer research, in which synergistic  
55 effects are included as one of the factors in the study of mixtures, has been proposed by the US President’s Cancer  
56 Panel [21]. As a result, testing the toxicity of chemical mixtures or environmental samples has been increasing in  
57 comparison to using a single compound [21].

58 Sample preparation methods for toxicity testing are crucial, and various fractionation methods have been established  
59 to obtain cleaner and more defined fractions. Most of the techniques involve initial crude and subsequent fine  
60 fractionations, followed by chemical analysis. These methods have several different names, including bioassay-  
61 directed chemical analysis [22], bioassay-directed fractionation [23-26], toxicity-based fractionation [27] and effect-  
62 directed fractionation [28]. Previous studies have been conducted using a similar procedure. First, sample extracts  
63 were treated with silica gel or cleaned up with solid phase extraction (SPE). Then the eluates were subjected to  
64 liquid chromatography (LC) to obtain refined fractions and sub-fractions, if further separation was necessary. Then,  
65 an aliquot was analysed for chemical content using high performance LC with ultra-violet and/or fluorescence

66 detection (HPLC-UV/FLD), or gas chromatography/mass spectrometry (GC/MS), while another aliquot was  
67 designated for toxicity tests [22-30]. Other studies have been conducted using solvent extraction methods with  
68 various combinations of organic solvents, or by SPE using different eluting conditions [31,32]. Additionally, an  
69 automated normal-phase LC (NPLC) fractionation system for PAHs in sediments using three LC columns was  
70 presented by Varel et al. [33].

71 Despite high-resolution fractionation steps, the final solution is often a mixture of several compounds with similar  
72 chemical properties. As a result, further investigations are required to narrow the candidates and identify the main  
73 compounds contributing to toxicological effects. Of the various chromatographic techniques used for selective  
74 fractionation, HPLC is considered to be a good tool, especially when appropriate HPLC columns are used to  
75 maximise the separation efficiency. Silica-based columns modified with amino, nitro, cyano, alkylhydroxy,  
76 alkylamino and alkylcyano moiety have been used for PAH analysis to remove aliphatic hydrocarbons and  
77 fractionate PAHs, based on the number of aromatic rings [34-37]. Polymeric C<sub>18</sub> columns were found to have unique  
78 selectivity for PAHs and a reversed-phase LC (RPLC) with FLD has replaced traditional NPLC for both preparative  
79 and quantitative PAH analysis [38-41]. Polymeric C<sub>18</sub> columns selectively retain isomeric PAHs, as explained in the  
80 “slot model” in which PAH retention is highly dependent on the molecular shape of the PAHs, that is, their ability to  
81 fit in the slit-like holes of the stationary phase [42]. However, fractions from an RPLC system normally contain  
82 water, which hinders fast evaporation during the preparation of samples for toxicity tests.

83 The ability to investigate the effects of interactions between B[a]P and other compounds in complex mixtures by  
84 selective inclusion/exclusion would be of significant benefit for future toxicity analyses. Thus, the aim of the present  
85 study was to develop an NPLC system for B[a]P-specific fractionation. The online system was equipped with an  
86 autosampler and switching valves controlled by a computer for precise and reproducible fractionation. The system  
87 was optimised using PAH standards and finally evaluated for separation efficiency and recoveries using a complex  
88 PAH mixture extracted from air particulate matter (PM) collected from the urban atmosphere of Stockholm.

89

## 90 2. Materials and methods

### 91 2.1 Chemicals and reagents

92 Solvents used for extraction, sample preparation and analysis were hexane, acetone and toluene (HPLC grade,  
93 Walkerburn Scotland) and dodecane (anhydrous, ≥99%, Sigma-Aldrich, St. Louis, MO, USA). All PAH standards  
94 and internal standards (ISs) used in the present study are specified with their abbreviation, CAS number, supplier,  
95 purity, molecular formula and molecular weight in the Supplementary data (S) **Table S1**.

### 96 2.2 Air sampling

97 Air PM was collected at the rooftop level (22 metres from the ground) of the Arrhenius building at Stockholm  
98 University [43]. Sampling was performed on a fluorocarbon coated glass fibre filter ( $\varnothing = 235$  mm, Fiberfilm Filters,  
99 Pallflex, Pall Corporation, Putnam, CT, USA) using an in-house constructed pump device equipped with a flow  
100 meter. The filter was desiccated for at least 24 hours before and after sampling, wrapped in aluminium foil and  
101 stored at  $-20$  °C until extraction. The total volume of sampled air was  $5141$  m<sup>3</sup>, and the average sampling rate was  
102  $72.5$  m<sup>3</sup>/hr. The sampling was performed in January 2013 (January 18-21) to avoid possible pollen introduction.

### 103 2.3 Extraction and SPE clean-up

104 The sampled filter was extracted with toluene for  $5 \times 30$  min using pressurised liquid extraction (ASE 200  
105 Accelerated Solvent Extraction System, Dionex Co., Sunnyvale, CA, USA). Then,  $0.6$  mL of the extract,  
106 corresponding to  $56$  m<sup>3</sup>, was applied to a silica SPE column ( $100$  mg, IST Isolute, Biotage, Cardiff, UK) and eluted  
107 with  $2$  mL of hexane to remove polar compounds and obtain a PAH enriched fraction. The final hexane eluate was  
108 gently evaporated to  $20$ - $50$   $\mu$ L under a nitrogen stream before online HPLC fractionation. The extraction and SPE  
109 clean-up procedures are described in detail elsewhere [44].

### 110 2.4 Online HPLC fractionation system

111 A schematic illustration of the automated system for B[a]P-specific fractionation is shown in **Fig. S1**.

112 The chromatographic system comprised an HPLC pump (Varian 9001, Varian Inc., Palo Alto, CA, USA), an  
113 autosampler (CMA/240, CMA Microdialysis AB, Stockholm, Sweden), two 2-(1-pyrenyl)ethyl (PYE) columns  
114 (Cosmosil PYE,  $150 \times 4.6$  mm,  $5$   $\mu$ m, Nacalai Tesque Inc., Kyoto, Japan) and a UV detector (Varian 9050, Varian  
115 Inc., Palo Alto, CA, USA). In addition, two position-switching valves with 3 (valve 2 and 3) and 4 (valve 1) ports  
116 (Valco Instruments Co. Inc., Houston, TX, USA) were used to change the column flow and collect the fractions into  
117 appointed test tubes. The autosampler was controlled with the CMA/200 v2.02 software (CMA Microdialysis AB,  
118 Stockholm, Sweden) and all hardware communications, including sample injection, valve switching and monitoring  
119 of the UV signal, were regulated by a personal computer equipped with the ELDS Win v1.1 software  
120 (Chromatography Data Systems AB, Svartsjö, Sweden). The temperature of the columns was set to  $30$  °C with a  
121 column oven (Croco-cil, Cluzeau Info Labo, Sainte-Foy-La-Grande, France) during operation to minimise run-to-  
122 run retention time variation.

### 123 2.5 LC-GC/MS analysis

124 The collected fractions were spiked with ISs and the volume was reduced to  $100$   $\mu$ L under a gentle stream of  
125 nitrogen. Then, the samples were transferred to micro vials for LC-GC/MS analysis. Analytical parameters are  
126 described in detail elsewhere [45].

127

## 128 3. Results and discussion

### 129 3.1 LC column selection

130 A test solution in hexane was prepared by mixing B[a]P with six structural isomers commonly encountered in  
131 environmental samples: benzo[a]fluoranthene (B[a]F), benzo[b]fluoranthene (B[b]F), benzo[j]fluoranthene (B[j]F),  
132 benzo[k]fluoranthene (B[k]F), benzo[e]pyrene (B[e]P) and perylene (Per). Initially, a nitrophenylpropyl (NO<sub>2</sub>) silica  
133 column (4.0 × 125 mm, 5 μm, Phenomenex, Torrance, CA, USA) was evaluated for the separation of the seven  
134 PAHs with molecular weights of 252 Da using hexane as the mobile phase. Hexane was selected as the mobile  
135 phase because it is a common solvent in NPLC systems and it has been recommended for increasing the resolution  
136 of PAHs [30]. However, insufficient separation was achieved using this set-up, and a PYE column was investigated  
137 to determine if it could improve the separation. Previously, this type of column has mainly been used for the  
138 separation of polychlorinated biphenyls (PCBs) and polychlorinated dibenzo-*p*-dioxins (PCDDs) [46-49]. Studies on  
139 PAH analysis using PYE columns are scarce, but this type of column was previously used in HPLC/UV analysis for  
140 the quantitation of PCBs and PAHs [50]. The retention of analytes on the PYE column is due to the donor-acceptor  
141 interaction in which the electron-rich pyrenyl group on the column function as electron donor to the analytes via π-π  
142 charge transfer, resulting in increased retention of planar analytes over non-planar ones owing to less steric  
143 hindrance [46]. The PYE column is also known to behave in the same manner as NO<sub>2</sub> or NH<sub>2</sub> substituted columns  
144 when used for PAH analysis in NPLC mode, that is, it separates PAHs based on their ring size but with less group  
145 separation efficiency than NO<sub>2</sub> or NH<sub>2</sub> substituted columns [51]. However, the PYE column is considered to be  
146 more suitable for isomer analysis [51]. Therefore, the PYE column was chosen in the present study because B[a]P  
147 separation from its isomers was the main issue rather than separation according to number of aromatic rings. The  
148 separation of the seven PAHs was improved and performed in a shorter time using the PYE column than with the  
149 NO<sub>2</sub> column. However, B[a]P partially co-eluted with B[a]F and B[b]F as shown in **Fig. 1, A**. Peaks were identified  
150 by comparing the retention times obtained from the injection of individual PAHs.

### 151 3.2 Serial column separation

152 Based on the enhanced separation observed on the PYE column compared to the NO<sub>2</sub> column, further exploration  
153 using the PYE column in the LC system was carried out to attempt to obtain baseline resolution of B[a]P from the  
154 other 252 Da molecular weight isomers. The peak capacity of an LC column is usually increased by the introduction  
155 of one or more columns into the existing set-up, either comprehensively or serially [52]. A serially coupled LC  
156 system is preferred compared to a comprehensive LC set-up because the peak capacity is greatly dependent on the  
157 orthogonality of the columns and requires relatively long separation times [52]. Previously reported studies have  
158 applied the serial column separation concept to complex samples of essential oils and digested proteins [53,54].  
159 Accordingly, two PYE columns with the same dimensions were connected in series to reduce the co-elution of  
160 B[a]P with B[a]F and B[b]F. The separation of the seven 252 Da PAHs was enhanced using the PYE-PYE serial  
161 column set-up when compared to the system using one PYE column (**Fig. 1, A and B**).

### 162 3.3 The B[a]P-specific online HPLC fractionation system

163 A UV chromatogram generated from injection of a standard solution containing 41 PAHs and 6 ISs in hexane on the  
164 automated B[a]P-specific online HPLC fractionation system is shown in **Fig. 2** and depicts the time points for  
165 collection of the different fractions. B[a]F and B[j]F were only used to prepare the test solution of the 252 Da PAHs  
166 and were not included in the standard PAH solution.

167 The large peak at the beginning of the chromatogram (5 min) is from toluene, which was used as the solvent for  
168 preparing the PAH stock solution. Each fractionation step is described in **Fig. S2, A-F**, in which the flow paths are  
169 differently coloured. The red line indicates the column flow before the backflush and the reversed column flow after  
170 the backflush is indicated with a blue line.

171 The system began in standby mode, with the mobile phase flowing through the columns and valves into the waste  
172 bottle connected to valve V3 (**Fig. S2, A**). After injection, the system was put in standby mode for approximately six  
173 minutes to direct mobile phase and residual toluene to waste (W1), then valve V3 was switched to collect PAHs  
174 eluting earlier than B[a]P in fraction 1 (F1) (**Fig. S2, B**). Just before the elution of B[a]P, the system began  
175 collecting fraction 2 (F2) by diverting the flow through valve V2 (**Fig. S2, C**).

176 Then, the column flow was reversed while valve V2 was kept at the same position to complete the collection of  
177 residual B[a]P in the tubing between the UV detector and the collecting tube (**Fig. S2, D**). Excess mobile phase was  
178 sent to the waste (W2) through valve V3 (**Fig. S2, E**), and then the backflush peak containing PAHs eluting later  
179 than B[a]P was collected in fraction 3 (F3) (**Fig. S2, F**). After finishing the collection, the system was returned to  
180 standby mode to prepare for the next run (**Fig. S2, A**). All valve switching times are presented in **Table S2** and the  
181 run time for one total cycle of fractionation was 30 minutes.

### 182 3.4 Fractionation of air PM extract

183 SPE cleaned-up extracts of Stockholm air PM were fractionated on the system in triplicate. A chromatogram  
184 generated from an injection of a Stockholm air PM extract is presented in **Fig. 3**.

185 Valve switching times for collecting F3 were adjusted to resolve peak broadening due to sample complexity.  
186 GC/MS chromatograms obtained in selected ion monitoring (SIM) mode, from which ions with  $m/z$  252 were  
187 extracted to assess the separation selectivity of the B[a]P-specific online HPLC fractionation system (displayed in  
188 **Fig. 4**), are shown in **Figs. S3, S4** and **S5**. The first fraction, F1 contained all four benzofluoranthenes but no  
189 detectable amount of B[a]P. A minor amount of B[b]F (2%) was present in fraction F2. B[e]P and Per peaks were  
190 observed in fraction F3, along with a small amount of B[a]P (6%).

191 Waste was collected before F1 (W1) and between F2 and F3 (W2) (**Fig. S6**) to confirm that no PAHs were lost in  
192 the waste. Chromatograms of fractions from injections of a blank sample and hexane are shown in **Figs. S7** and **S8**.  
193 Two system peaks were observed in GC/MS chromatograms from sample waste fractions and from the fractions

194 generated by blank and hexane injections. The system peaks did not interfere with other PAH peaks in the GC/MS  
195 chromatogram. GC/MS full scan spectra of the system peaks are shown in **Fig. S9**. A NIST library search found no  
196 relevant matches for these peaks, so another GC/MS system with different chromatographic conditions (**Table S3**)  
197 was used to analyse the eluate from the B[a]P-specific online HPLC fractionation system. Three major peaks were  
198 detected (**Fig. S10**) and matched with three different substituted pyrenes from the NIST library. The highest ranked  
199 hits for each peak were 1-hydroxypyrene, 1-pyrene-carboxaldehyde and 1-acetylpyrene, respectively. These  
200 compounds most likely originate from the pyrenyl moiety in the stationary phases of the HPLC columns (**Fig. S11**).  
201 In addition, no peaks were found in the GC/MS chromatograms when the eluate collected from the system without  
202 PYE columns was injected, which eliminates the injector, HPLC pump, valves and tubing as possible sources of the  
203 system peaks.

### 204 3.5 Toxicological testing of chromatographic blanks

205 The blank eluates generated from the fractionation system were applied to toxicity tests aimed at evaluating DNA  
206 damage signalling and cytotoxicity (methods are described in further detail in **S7**). Chk1 and H2AX are key signal  
207 transduction proteins activated by phosphorylation in response to DNA damage resulting from cellular exposure to  
208 PAHs [10]. The results of the present study showed no phosphorylation of Chk1 or H2AX in human HepG2 cells  
209 following exposure to the eluates whereas a 17- and 2-fold induction of Chk1 and H2AX, respectively, was  
210 observed in response to 3  $\mu\text{M}$  B[a]P (**Fig. S12, A-C**). No effects were observed on the levels of p53 protein in any  
211 of the exposures (**Fig. S12, A and D**). Furthermore, no effects were observed on cell viability (**Fig. S12, E**). Taken  
212 together, these data demonstrate that the blank eluates generated from the fractionation system are non-cytotoxic to  
213 cells and do not induce DNA damage related-signalling.

### 214 3.6 Recovery of PAHs and separation efficiency of B[a]P

215 Quantitative PAH results are summarised in **Fig. 5**. Non-fractionated samples were prepared by spiking ISs into the  
216 SPE cleaned-up samples, while the collected fractions were spiked with ISs after the fractionation process. The  
217 average recovery of PAHs in the fractionation system was approximately 80%, varying from 61 to 105%. Main  
218 losses were due to residue remaining in the micro vial and injection needle. Therefore, the PAH concentrations of  
219 fractions generated on the B[a]P-specific online HPLC fractionation system should be re-analysed prior to  
220 toxicological tests to establish more accurate dose-responses.

221 A relative comparison between different fractions was performed to investigate the separation efficiency of the  
222 B[a]P-specific online HPLC fractionation system. As shown in **Fig. 6**, 94% of B[a]P was found in F2 and 6% was  
223 found in F3, but no detectable amount was found in F1. The amount found in F3 could be explained by peak tailing  
224 of B[a]P on the LC columns, that is, a slight overlap with the following B[e]P peak. The standard deviations of the  
225 mean values of the relative amount were less than 1%. In the fractionation, the valve switching times were set such  
226 that greater than 90% of B[a]P could be collected while minimising the introduction of other PAHs in the same  
227 fraction. In addition, full scan GC/MS analysis of blank, non-fractionated and fractionated samples were conducted



228 in order to examine co-eluting compounds in the same region of chromatograms as the 252 Da PAHs (**Fig. S13**).  
229 Three unknown compounds (peak 5, 6 and 10) were detected in the chromatogram of the non-fractionated sample.  
230 Peak 6 and 10 could potentially be methylated isomers of 252 Da PAHs because their mass spectra features a  
231 molecular ion at  $m/z$  266. Peak 6 was also detected in F1, but no visible peaks were detected in F2 and F3 (**Fig. S13**,  
232 **F1-F3**). Methylated PAHs usually elute after the native PAHs on the PYE column in NPLC mode [55], but some  
233 methylated PAHs such as 1-methyl perylene elute earlier than their native forms [56]. Among possible methylated  
234 isomers, methylated B[a]Ps are of concern related to their mutagenicity [57-59] though they have been reported to  
235 be present at relatively low levels compared to that of B[a]P in coal tar [60]. On the assumption of peak 6 being one  
236 isomer of methylated B[a]P, its presence in F1 might contribute to the toxicity of the B[a]P-free fraction. However,  
237 the contribution to the toxicity is assumed to be low because the peak is relatively small.

238 The PAH concentrations determined in the initial toluene extracts (**Table S4**) were comparable to those from a  
239 previous study conducted at the same sampling site [10]. The results of the present study demonstrated that the  
240 B[a]P-specific online HPLC fractionation system can selectively remove B[a]P from a complex Stockholm air PM  
241 extract. In addition, this system set-up can easily be modified by changing columns, mobile phase and valve  
242 switching times to selectively remove other PAHs of interest.

243

#### 244 4. Conclusions

245 The present study established an online HPLC fractionation system for the selective removal of B[a]P from the  
246 complex PAH mixture of an air PM extract using two PYE columns in series. The system showed a high selectivity,  
247 removing 94% of B[a]P and less than 3% of B[b]F from the complex mixture. One fractionation cycle requires only  
248 30 minutes, which makes it possible to do repeated fractionations for larger samples. The system is simple and easy  
249 to re-construct in a normal chromatography lab, because it only requires analytical scale tubing, a pump, an injector,  
250 columns and a detector. More fractions can easily be obtained by adding additional switch valves or by introducing  
251 fraction collectors to the current system. In addition, automatic fractionation minimises the exposure of operators to  
252 toxicants and increases the precision and time efficiency. Sample overloading can be avoided by repeated  
253 fractionations and the system can be modified by changing mobile phase and columns, allowing for selective  
254 fractionation of other target compounds and matrices. The mixture as a whole or lacking B[a]P obtained from the  
255 system can assist in the evaluation of B[a]P toxicity by considering interactions between B[a]P and other co-existing  
256 compounds. Blanks generated on the fractionation system were not cytotoxic to cells and did not induce DNA  
257 damage-related signalling. In conclusion, this system is a useful tool for elucidating the toxicological effects of  
258 B[a]P in complex PAH mixtures. Furthermore, the mixture-based toxicity evaluation of not only B[a]P, but also  
259 other target compounds, could improve current methods for risk assessment.

260

261 Acknowledgements

262 The authors are grateful to Ioannis Sadiktis for assistance with GC/MS analyses. This work was supported by  
263 Stockholm University, Karolinska Institutet and the Swedish Research Council Formas (R&D project No. 2012-  
264 478).

265 Appendix A. Supplementary data

266 Supplementary data associated with this article can be found, in the online version, at  
267 <http://dx.doi.org/10.1016/j.chroma.xxxx.xx.xxx>.

268 References

269 [1] IARC, Monographs on the Evaluation of Carcinogenic Risks to Humans: Some Non-heterocyclic Polycyclic  
270 Aromatic Hydrocarbons and Some Related Exposures, International Agency for Research on Cancer (IARC), Lyon,  
271 2010.

272 [2] J. Brinkmann, K. Stolpmann, S. Trappe, T. Otter, D. Genkinger, U. Bock, M. Liebsch, F. Henkler, C. Hutzler, A.  
273 Luch, *Toxicol. Sci.* 131 (2013) 351.

274 [3] Y. Chen, C. Huang, C. Bai, H. Gao, R. Ma, X. Liu, Q. Dong, *Toxicology* 304 (2013) 167.

275 [4] S. Cheng, Y. Xia, J. He, X. Liu, X. Chen, Y. Ding, Y. Wang, B. Peng, B. Tu, *Environ. Toxicol. Phar.* 36 (2013)  
276 648.

277 [5] M. Pufulete, J. Battershill, A. Boobis, R. Fielder, *Regul. Toxicol. Pharm.* 40 (2004) 54.

278 [6] C.E. Boström, P. Gerde, A. Hanberg, B. Jernström, C. Johansson, T. Kyrklund, A. Rannug, M. Törnqvist, K.  
279 Victorin, R. Westerholm, *Environ. Health Perspect.* 110 (Suppl. 3) (2002) 451.

280 [7] WHO, Air Quality Guidelines for Europe, World Health Organization (WHO) Regional Publications,  
281 Copenhagen, 1987.

282 [8] WHO, Air Quality Guidelines for Europe, Second Edition, World Health Organization (WHO) Regional Office  
283 for Europe, Copenhagen, 2000.

284 [9] D. Delistraty, *Toxicol. Environ. Chem.* 64 (1997) 81.

285 [10] I.W.H. Jarvis, C. Bergvall, M. Bottai, R. Westerholm, U. Stenius, K. Dreij, *Toxicol. Appl. Pharmacol.* 266  
286 (2013) 408.

287 [11] Å Mattsson, S. Lundstedt, U. Stenius, *Environ. Mol. Mutagen.* 50 (2009) 337.

288 [12] M. Niziolek-Kierecka, K. Dreij, S. Lundstedt, U. Stenius, *Chem. Res. Toxicol.* 25 (2012) 862.

289 [13] Y.C. Staal, D.G. Hebels, M.H. van Herwijnen, R.W. Gottschalk, F.J. van Schooten, J.H. van Delft,  
290 *Carcinogenesis* 28 (2007) 2632.

291 [14] A. Tarantini, A. Maitre, E. Lefebvre, M. Marques, C. Marie, J.L. Ravanat, T. Douki, *Mutat. Res. – Fund. Mol.*  
292 *M.* 671 (2009) 67.

- 293 [15] A. Tarantini, A. Maître, E. Lefèbvre, M. Marques, A. Rajhi, T. Douki, *Toxicology* 279 (2011) 36.
- 294 [16] L.A. Courter, A. Luch, T. Musafia-Jeknic, V.M. Arlt, K. Fischer, R. Bildfell, C. Pereira, D.H. Phillips, M.C.  
295 Poirier, W.M. Baird, *Cancer Lett.* 265 (2008) 135.
- 296 [17] B. Mahadevan, C.P. Marston, A. Luch, W.M. Dashwood, E. Brooks, C. Pereira, J. Doehmer, W.M. Baird, *Int.*  
297 *J. Cancer* 120 (2007) 1161.
- 298 [18] C.P. Marston, C. Pereira, J. Ferguson, K. Fischer, O. Hedstrom, W.M. Dashwood, W.M. Baird, *Carcinogenesis*  
299 22 (2001) 1077.
- 300 [19] D.W. Gaylor, S.J. Culp, L.S. Goldstein, F.A. Beland, *Risk Anal.* 20 (2000) 81.
- 301 [20] L.K. Siddens, A. Larkin, S.K. Krueger, C.A. Bradfield, K.M. Waters, S.C. Tilton, C.B. Pereira, C.V. Löhr,  
302 V.M. Arlt, D.H. Phillips, D.E. Williams, W.M. Baird, *Toxicol. Appl. Pharmacol.* 264 (2012) 377.
- 303 [21] NCI, Reducing Environmental Cancer Risk: What We Can Do Now, U.S. Department of Health and Human  
304 Services (HHS), National Institutes of Health, National Cancer Institute (NCI) (2010)
- 305 [22] D. Schuetzle, J. Lewtas, *Anal. Chem.* 58 (1986) 1060A.
- 306 [23] M.G. Nishioka, C.C. Howard, D.A. Contos, L.M. Ball, J. Lewtas, *Environ. Sci. Technol.* 22 (1988) 908.
- 307 [24] J. Arey, W.P. Harger, D. Helmig, R. Atkinson, *Mutat. Res.* 281 (1992) 67.
- 308 [25] C.H. Marvin, B.E. McCarry, J.A. Lundrigan, K. Roberts, D.W. Bryant, *Sci. Total Environ.* 231 (1999) 135.
- 309 [26] S. Biselli, N. Reineke, N. Heinzl, U. Kammann, S. Franke, H. Hühnerfuss, N. Theobald, *J. Soils Sediments* 5  
310 (2005) 171.
- 311 [27] P.A. Kosian, E.A. Makynen, P.D. Monson, D.R. Mount, A. Spacie, O.G. Mekenyan, G.T. Ankley, *Environ.*  
312 *Toxicol. Chem.* 17 (1998) 1021.
- 313 [28] W. Brack, K. Schirmer, T. Kind, S. Schrader, G. Schüürmann, *Environ. Toxicol. Chem.* 21 (2002) 2654.
- 314 [29] T. Alsberg, M. Strandell, R. Westerholm, U. Stenberg, *Environ. Int.* 11 (1985) 249.
- 315 [30] W. Brack, T. Kind, H. Hollert, S. Schrader, M. Möder, *J. Chromatogr. A* 986 (2003) 55.
- 316 [31] H.F. Wesp, X. Tang, R. Edenharder, *Mutat. Res.* 472 (2000) 1.
- 317 [32] R.E. Cochran, N. Dongari, H. Jeong, J. Beránek, S. Haddadi, J. Shipp, A. Kubátová, *Anal. Chim. Acta* 740  
318 (2012) 93.
- 319 [33] U. Lübcke-von Varell, G. Streck, W. Brack, *J. Chromatogr. A* 1185 (2008) 31.
- 320 [34] S.A. Wise, S.N. Chesler, H.S. Hertz, L.R. Hilpert, W.E. May, *Anal. Chem.* 49 (1977) 2306.
- 321 [35] C. Östman, A. Colmsjö, *Chromatographia* 25 (1988) 25.
- 322 [36] E.P. Lankmayr, K. Müller, *J. Chromatogr. A* 170 (1979) 139.
- 323 [37] J. Chmielowiec, A.E. George, *Anal. Chem.* 52 (1980) 1154.
- 324 [38] W.E. May, S.A. Wise, *Anal. Chem.* 56 (1984) 225.

- 325 [39] L.C. Sander, S.A. Wise, *Anal. Chem.* 56 (1984) 504.
- 326 [40] L.C. Sander, S.A. Wise, *J. Chromatogr. A* 316 (1984) 163.
- 327 [41] L.C. Sander, S.A. Wise, *Anal. Chem.* 59 (1987) 2309.
- 328 [42] S.A. Wise, L.C. Sander, *J. High Res. Chromatog* 8 (1985) 248.
- 329 [43] C. Bergvall, R. Westerholm, *Environ. Sci. Technol.* 41 (2007) 731.
- 330 [44] C. Bergvall, R. Westerholm, *Anal. Bioanal. Chem.* 391 (2008) 2235.
- 331 [45] I. Sadiktsis, J.H. Koegler, T. Benham, C. Bergvall, R. Westerholm, *Fuel* 115 (2014) 573.
- 332 [46] P. Haglund, L. Asplund, U. Järnberg, B. Jansson, *J. Chromatogr. A* 507 (1990) 389.
- 333 [47] K. Kimata, K. Hosoya, T. Araki, N. Tanaka, E.R. Barnhart, L.R. Alexander, S. Sirimanne, P.C. McClure, J.  
334 Grainger, D.G. Patterson Jr., *Anal. Chem.* 65 (1993) 2502.
- 335 [48] M.M. Krahn, G.M. Ylitalo, J. Buzitis, C.A. Sloan, D.T. Boyd, S.-L. Chan, U. Varanasi, *Chemosphere* 29 (1994)  
336 117.
- 337 [49] K. Kimata, K. Hosoya, H. Kuroki, N. Tanaka, J.R. Barr, P.C. McClure, D.G. Patterson Jr., E. Jakobsson, A.  
338 Bergman, *J. Chromatogr. A* 786 (1997) 237.
- 339 [50] C. Nerín, C. Domeño, *Analyst* 124 (1999) 67.
- 340 [51] U.L. Nilsson, A.L. Colmsjö, *Chromatographia* 34 (1992) 115.
- 341 [52] E. Blahová, P. Jandera, F. Cacciola, L. Mondello, *J. Sep. Sci.* 29 (2006) 555.
- 342 [53] M. Herrero, F. Cacciola, P. Donato, D. Giuffrida, G. Dugo, P. Dugo, L. Mondello, *J. Chromatogr. A* 1188  
343 (2008) 208.
- 344 [54] P. Donato, P. Dugo, F. Cacciola, G. Dugo, L. Mondello, *J. Sep. Sci.* 32 (2009) 1129.
- 345 [55] L.C. Sander, R.M. Parris, S.A. Wise, *Anal. Chem.* 63 (1991) 2589.
- 346 [56] S.A. Wise, L.C. Sander, *J. Chromatogr. A* 514 (1990) 111.
- 347 [57] O.G. Fahmy, M.J. Fahmy, *Cancer Res.* 33 (1973) 302.
- 348 [58] R.P. Iyer, J.W. Lyga, J.A. Secrist, G.H. Daub, T.J. Slaga, *Cancer Res.* 40 (1980) 1073.
- 349 [59] R. Santella, T. Kinoshita, A.M. Jeffrey, *Mutat. Res.* 104 (1982) 209.
- 350 [60] P. Garrigues, J. Bellocq, S.A. Wise, *Fresen. J. Anal. Chem.* 336 (1990) 106.

351 Figure captions

352

353 **Fig. 1.** HPLC/UV chromatogram of a PAH standard solution injected on (A) a PYE column and (B) a serial PYE  
354 column set-up. (A) Conditions: hexane 100%, flow rate 0.7 mL/min, ambient column temperature, UV detection at  
355 254 nm, (B) Conditions: hexane 100%, flow rate 0.8 mL/min, column temperature at 30 °C, UV detection at 254  
356 nm. UV detector sensitivity of 0.02 absorbance units full scale (AUFS). AU: absorbance units. Ant\*: Anthracene  
357 was added to the mixture as a retention indicator.

358 **Fig. 2.** HPLC/UV chromatogram of the 41 PAH standard solution from the online HPLC fractionation system: (A)  
359 full-scale chromatogram, (B) enlarged chromatogram depicting fractionation points. Conditions: two PYE columns  
360 (150 × 4.6 mm, 5 μm), hexane 100%, flow rate 0.8 mL/min, column temperature at 30 °C, UV detection at 254 nm,  
361 UV detector sensitivity of 0.02 AUFS. AU: absorbance units, F: fraction, W: waste.


362 **Fig. 3.** HPLC/UV chromatogram of an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC  
363 fractionation system. Instrumental conditions are described in **Fig. 2.**

364 **Fig. 4.** GC/MS extracted ion chromatograms ( $m/z$  252) of the fractions generated from an injection of a Stockholm  
365 air PM extract on the B[a]P-specific online HPLC fractionation system.


366 **Fig. 5.** PAH concentrations determined in a Stockholm air PM extract and in the fractions collected from the B[a]P-  
367 specific online HPLC fractionation system. Error bars show one standard deviation from the mean value ( $n = 3$ ).  
368 Compound abbreviations are shown in **Table S1.**

369 **Fig. 6.** Relative abundance of 252 Da molecular weight PAHs in the fractions generated by the B[a]P-specific online  
370 HPLC fractionation system. Error bars show one standard deviation from the mean value ( $n = 3$ ).


## Figures




**Fig. 1.** HPLC/UV chromatogram of a PAH standard solution injected on (A) a PYE column and (B) a serial PYE column set-up. (A) Conditions: hexane 100%, flow rate 0.7 mL/min, ambient column temperature, UV detection at 254 nm, (B) Conditions: hexane 100%, flow rate 0.8 mL/min, column temperature at 30 °C, UV detection at 254 nm. UV detector sensitivity of 0.02 absorbance units full scale (AUFS). AU: absorbance units. Ant\*: Anthracene was added to the mixture as a retention indicator.



**Fig. 2.** HPLC/UV chromatogram of the 41 PAH standard solution from the online HPLC fractionation system: (A) full-scale chromatogram, (B) enlarged chromatogram depicting fractionation points. Conditions: two PYE columns (150 × 4.6 mm, 5 μm), hexane 100%, flow rate 0.8 mL/min, column temperature at 30 °C, UV detection at 254 nm, UV detector sensitivity of 0.02 AUFS. AU: absorbance units, F: fraction, W: waste.




**Fig. 3.** HPLC/UV chromatogram of an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC fractionation system. Instrumental conditions are described in **Fig. 2**.




**Fig. 4.** GC/MS extracted ion chromatograms ( $m/z$  252) of the fractions generated from an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC fractionation system.





**Fig. 5.** PAH concentrations determined in a Stockholm air PM extract and in the fractions collected from the B[a]P-specific online HPLC fractionation system. Error bars show one standard deviation from the mean value ( $n = 3$ ). Compound abbreviations are shown in **Table S1**.



**Fig. 6.** Relative abundance of 252 Da molecular weight PAHs in the fractions generated by the B[a]P-specific online HPLC fractionation system. Error bars show one standard deviation from the mean value ( $n = 3$ ).

## 1 **Supplementary data**

2 Benzo[*a*]pyrene-specific online HPLC fractionation of air particulate extracts – A tool for  
3 evaluating biological interactions

4 Hwanmi Lim<sup>a</sup>, Christoffer Bergvall<sup>a</sup>, Ian Jarvis<sup>b</sup>, Åse Mattsson<sup>b</sup>, Kristian Dreij<sup>b</sup>, Ulla Stenius<sup>b</sup>, Roger Westerholm<sup>a,\*</sup>

5 <sup>a</sup> *Department of Analytical Chemistry, Stockholm University, Svante Arrhenius väg 16 SE-106 91 Stockholm,*  
6 *Sweden*

7 <sup>b</sup> *Institute of Environmental Medicine, Karolinska Institutet, Box 210, SE-171 77 Stockholm, Sweden*

8

## 9 **Contents**

10 **S-1.** List of PAH standards and internal standards

11 **S-2.** Schematic illustration of the B[*a*]P-specific online HPLC fractionation system

12 **S-3.** Valve positions during a run on the B[*a*]P-specific online HPLC fractionation system

13 **S-4.** GC/MS chromatograms of Stockholm air PM extract from the B[*a*]P-specific online HPLC fractionation system

14 **S-5.** GC/MS chromatograms of blank and hexane injections on the B[*a*]P-specific online HPLC fractionation system

15 **S-6.** GC/MS full scan mass spectra of the system peaks from the B[*a*]P-specific online HPLC fractionation system

16 **S-7.** Biological tests of eluates from the B[*a*]P-specific online HPLC fractionation system

17 **S-8.** GC/MS full scan analysis of blank, non-fractionated and fractionated samples

18 **S-9.** PAH concentrations in Stockholm air PM sample

19 *S-I. List of PAH standards and internal standards*

20 **Table S1.** List of PAH standards and surrogate internal standards with abbreviation, CAS Registry number, supplier, purity, molecular formula (M.f.) and  
21 molecular weight (M.w. in g/mol)

Name	Abbreviation	CAS Registry No.	Supplier	Purity (%)	M.f.	M.w.
Phenanthrene	Phe	85-01-8	Merck, Germany	98.7	C <sub>14</sub> H <sub>10</sub>	178.23
Anthracene	Ant	120-12-7	Sigma-Aldrich, St. Louis, MO, USA	99.6	C <sub>14</sub> H <sub>10</sub>	178.23
3-Methylphenanthrene	3-MPhe	832-71-3	Larodan Fine Chemicals AB, Sweden	99.9	C <sub>15</sub> H <sub>12</sub>	192.26
2-Methylphenanthrene	2-MPhe	2531-84-2	Sigma-Aldrich, St. Louis, MO, USA	93.8	C <sub>15</sub> H <sub>12</sub>	192.26
2-Methylanthracene	2-MAnt	613-12-7	Koch-Light Laboratories, UK	100	C <sub>15</sub> H <sub>12</sub>	192.26
9-Methylphenanthrene	9-MPhe	883-20-5	Chiron AS, Trondheim, Norway	100	C <sub>15</sub> H <sub>12</sub>	192.26
1-Methylphenanthrene	1-MPhe	832-69-9	Larodan Fine Chemicals AB, Sweden	98.8	C <sub>15</sub> H <sub>12</sub>	192.26
4H-Cyclopenta[def]phenanthrene	4H-CPP	203-64-5	Sigma-Aldrich, St. Louis, MO, USA	99.1	C <sub>15</sub> H <sub>10</sub>	190.24
2-Phenylanthracene	2-PN	612-94-2	EGA-Chemie, Steinheim, Germany	94.2	C <sub>16</sub> H <sub>12</sub>	204.27
3,6-Dimethylphenanthrene	3,6-DMPhe	1576-67-6	Larodan Fine Chemicals AB, Sweden	96.1	C <sub>16</sub> H <sub>14</sub>	206.28
3,9-Dimethylphenanthrene	3,9-DMPhe	66291-32-5	Chiron AS, Trondheim, Norway	99.7	C <sub>16</sub> H <sub>14</sub>	206.28

Fluoranthene	Flu	206-44-0	Sigma-Aldrich, St. Louis, MO, USA	97.2	C <sub>16</sub> H <sub>10</sub>	202.25
Pyrene	Pyr	129-00-0	Janssen Chimica, Belgium	97.6	C <sub>16</sub> H <sub>10</sub>	202.25
1-Methylfluoranthene	1-MFlu	25889-60-5	Chiron AS, Trondheim, Norway	99.7	C <sub>17</sub> H <sub>12</sub>	216.28
Benzo[ <i>a</i> ]fluorene	B[ <i>a</i> ]f	238-84-6	Chiron AS, Trondheim, Norway	98.8	C <sub>17</sub> H <sub>12</sub>	216.28
Benzo[ <i>b</i> ]fluorene	B[ <i>b</i> ]f	243-17-4	Sigma-Aldrich, St. Louis, MO, USA	99.1	C <sub>17</sub> H <sub>12</sub>	216.28
2-Methylpyrene	2-MPyr	3442-78-2	Chiron AS, Trondheim, Norway	98.3	C <sub>17</sub> H <sub>12</sub>	216.28
4-Methylpyrene	4-MPyr	3353-12-6	Chiron AS, Trondheim, Norway	99.6	C <sub>17</sub> H <sub>12</sub>	216.28
1-Methylpyrene	1-MPyr	2381-21-7	Larodan Fine Chemicals AB, Sweden	99.1	C <sub>17</sub> H <sub>12</sub>	216.28
Benzo[ <i>ghi</i> ]fluoranthene	B[ <i>ghi</i> ]F	203-12-3	Larodan Fine Chemicals AB, Sweden	99.5	C <sub>18</sub> H <sub>10</sub>	226.27
Benzo[ <i>c</i> ]phenanthrene	B[ <i>c</i> ]Phe	195-19-7	Chiron AS, Trondheim, Norway	99.5	C <sub>18</sub> H <sub>12</sub>	228.29
Benz[ <i>a</i> ]anthracene	B[ <i>a</i> ]A	56-55-3	Fluka AG, Switzerland	98.4	C <sub>18</sub> H <sub>12</sub>	228.29
3-Methylchrysene	3-MChr	3351-31-3	Chiron AS, Trondheim, Norway	99.1	C <sub>19</sub> H <sub>14</sub>	242.31
2-Methylchrysene	2-MChr	3351-32-4	Chiron AS, Trondheim, Norway	99.7	C <sub>19</sub> H <sub>14</sub>	242.31
6-Methylchrysene	6-MChr	1705-85-7	Chiron AS, Trondheim, Norway	100	C <sub>19</sub> H <sub>14</sub>	242.31
1-Methylchrysene	1-MChr	3351-28-8	Chiron AS, Trondheim, Norway	99.3	C <sub>19</sub> H <sub>14</sub>	242.31
Benzo[ <i>b</i> ]fluoranthene	B[ <i>b</i> ]F	205-99-2	Chem Service, West Chester, PA, USA	100	C <sub>20</sub> H <sub>12</sub>	252.31

Benzo[ <i>k</i> ]fluoranthene	B[ <i>k</i> ]F	207-08-9	Chem Service, West Chester, PA, USA	98.3	C <sub>20</sub> H <sub>12</sub>	252.31
Benzo[ <i>e</i> ]pyrene	B[ <i>e</i> ]P	192-97-2	Sigma-Aldrich, St. Louis, MO, USA	99.7	C <sub>20</sub> H <sub>12</sub>	252.31
Benzo[ <i>a</i> ]pyrene	B[ <i>a</i> ]P	50-32-8	Sigma-Aldrich, St. Louis, MO, USA	97.6	C <sub>20</sub> H <sub>12</sub>	252.31
Perylene	Per	198-55-0	Sigma-Aldrich, St. Louis, MO, USA	99.5	C <sub>20</sub> H <sub>12</sub>	252.31
Indeno[1,2,3- <i>cd</i> ]fluoranthene	I[1,2,3- <i>cd</i> ]F	193-43-1	Radiant Dyes, Wermelskirchen, Germany	98.4	C <sub>22</sub> H <sub>12</sub>	276.33
Indeno[1,2,3- <i>cd</i> ]pyrene	I[1,2,3- <i>cd</i> ]P	193-39-5	AccuStandard Inc., New Haven, CT, USA	99.8	C <sub>22</sub> H <sub>12</sub>	276.33
Dibenz[ <i>a,h</i> ]anthracene	DB[ <i>a,h</i> ]A	53-70-3	Fluka AG, Switzerland	99.4	C <sub>22</sub> H <sub>14</sub>	278.35
Picene	Pic	213-46-7	Larodan Fine Chemicals AB, Sweden	96.1	C <sub>22</sub> H <sub>14</sub>	278.35
Benzo[ <i>ghi</i> ]perylene	B[ <i>ghi</i> ]p	191-24-2	Janssen Chimica, Belgium	98.8	C <sub>22</sub> H <sub>12</sub>	276.33
Dibenzo[ <i>a,l</i> ]pyrene	DB[ <i>a,l</i> ]P	191-30-0	AccuStandard Inc., New Haven, CT, USA	96	C <sub>24</sub> H <sub>14</sub>	302.37
Dibenzo[ <i>a,e</i> ]pyrene	DB[ <i>a,e</i> ]P	192-65-4	LGC Promochem, Sweden	98	C <sub>24</sub> H <sub>14</sub>	302.37
Coronene	Cor	191-07-1	Radiant Dyes, Wermelskirchen, Germany	100	C <sub>24</sub> H <sub>12</sub>	300.35
Dibenzo[ <i>a,i</i> ]pyrene	DB[ <i>a,i</i> ]P	189-55-9	Sigma-Aldrich, St. Louis, MO, USA	96.4	C <sub>24</sub> H <sub>14</sub>	302.37


Dibenzo[ <i>a,h</i> ]pyrene	DB[ <i>a,h</i> ]P	189-64-0	Koch-Light Laboratories, UK	100	C <sub>24</sub> H <sub>14</sub>	302.37
Phenanthrene-D <sub>10</sub>	Phe-D <sub>10</sub>	1517-22-2	Chiron AS, Trondheim, Norway <sup>a</sup> , Larodan Fine Chemicals AB, Sweden <sup>b</sup>	99.3 <sup>a</sup> , 97.2 <sup>b</sup>	C <sub>14</sub> D <sub>10</sub>	188.29
Pyrene-D <sub>10</sub>	Pyr-D <sub>10</sub>	1718-52-1	Chiron AS, Trondheim, Norway <sup>a</sup> , Larodan Fine Chemicals AB, Sweden <sup>b</sup>	99.8 <sup>a</sup> , 95.7 <sup>b</sup>	C <sub>16</sub> D <sub>10</sub>	212.31
Benz[ <i>a</i> ]anthracene-D <sub>12</sub>	B[ <i>a</i> ]A-D <sub>12</sub>	1718-53-2	Chiron AS, Trondheim, Norway	98.6	C <sub>18</sub> D <sub>12</sub>	240.36
Benzo[ <i>a</i> ]pyrene-D <sub>12</sub>	B[ <i>a</i> ]P-D <sub>12</sub>	63466-71-7	Chiron AS, Trondheim, Norway <sup>a</sup> , Larodan Fine Chemicals AB, Sweden <sup>b</sup>	98.7 <sup>a</sup> , 98.1 <sup>b</sup>	C <sub>20</sub> D <sub>12</sub>	264.38
Benzo[ <i>ghi</i> ]perylene-D <sub>12</sub>	B[ <i>ghi</i> ]p-D <sub>12</sub>	93951-66-7	Chiron AS, Trondheim, Norway	99.6	C <sub>22</sub> D <sub>12</sub>	288.4
Dibenzo[ <i>a,i</i> ]pyrene-D <sub>14</sub>	DB[ <i>a,i</i> ]P-D <sub>14</sub>	158776-07-9	LGC Promochem, Sweden	98	C <sub>24</sub> D <sub>14</sub>	316.45
Benzo[ <i>a</i> ]fluoranthene <sup>c</sup>	B[ <i>a</i> ]F	203-33-8	Gift from NIST	-	C <sub>20</sub> H <sub>12</sub>	252.31
Benzo[ <i>j</i> ]fluoranthene <sup>c</sup>	B[ <i>j</i> ]F	205-82-3	Larodan Fine Chemicals AB, Sweden	100	C <sub>20</sub> H <sub>12</sub>	252.31

22 <sup>a</sup> Surrogate internal standard working solution used for sample preparation

23 <sup>b</sup> Surrogate internal standard used for PAH calibration standard preparation

24 <sup>c</sup> Not used for quantitation


25 S-2. Schematic illustration of the B[a]P-specific online HPLC fractionation system



26

27 **Fig. S1.** Scheme of the online HPLC fractionation system.

28 S-3. Valve positions during a run on the B[a]P-specific online HPLC fractionation system



29

30 **Fig. S2.** Scheme illustrating the valve positions and flow paths during a run of the B[a]P-specific online HPLC fractionation system: (A)-(C) and (D)-(F) indicate  
31 column flow before and after backflush, respectively. AS: autosampler, V: two position-switching valve, UVD: UV detector, F: fraction, W: waste.



32 **Table S2.** Valve positions during the different steps of a run on the B[a]P-specific online HPLC fractionation system


Time (min)	V1	V2	V3	Shown in Figure
0.00	0	0	0	<b>Fig. S2, A</b>
6.05 <sup>a</sup>	0	0	1	<b>Fig. S2, B</b>
12.35 <sup>b</sup>	0	1	1	<b>Fig. S2, C</b>
12.85 <sup>c</sup>	1	1	1	<b>Fig. S2, D</b>
13.05 <sup>b</sup>	1	0	0	<b>Fig. S2, E</b>
22.20 <sup>a</sup>	1	0	1	<b>Fig. S2, F</b>
29.20 <sup>a</sup>	1	0	0	<b>Fig. S2, A</b>
30.00	0	0	0	<b>Fig. S2, A</b>

33 <sup>a</sup> Adjusted time corrected by adding 0.2 minutes from the observed retention time due to the dead volume after the UV detector

34 <sup>b</sup> Adjusted backflush time corrected by subtracting 0.05 minutes from the observed retention time due to the time lag of the switching valve


35 <sup>c</sup> Adjusted time corrected by adding 0.15 minutes from the observed retention time due to the dead volume after the UV detector

36 S-4. GC/MS chromatograms of Stockholm air PM extract from the B[a]P-specific online HPLC fractionation system



37


38 **Figure S3.** GC/MS chromatogram obtained in SIM mode of fraction F1 from an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC  
39 fractionation system. 1: Phe, 2: Ant, 3: 3-MPhe, 4: 2-MPhe, 5: 2-MAnt, 6: 9-MPhe, 7: 1-MPhe, 8: 4*H*-CPP, 9: 2-PN, 10: 3,6-DMPhe, 11: 3,9-DMPhe, 12: Flu,  
40 13: Pyr, 14: 1-MFlu, 15: B[*a*]f, 16: B[*b*]f, 17: 2-MPyr, 18: 4-MPyr, 19: 1-MPyr, 20: B[*ghi*]F, 21: B[*c*]Phe, 22: B[*a*]A, 23: 3-MChr, 24: 2-MChr, 25: 6-MChr, 26:  
41 1-MChr, 27: B[*b*]F, 28: B[*k*]F, IS-1: Phe-D<sub>10</sub>, IS-2: Pyr-D<sub>10</sub>, IS-3: B[*a*]A-D<sub>12</sub>, IS-4: B[*a*]P-D<sub>12</sub>, IS-5: B[*ghi*]p-D<sub>12</sub>, IS-6: DB[*a,i*]P-D<sub>14</sub>.



42


43 **Figure S4.** GC/MS chromatogram obtained in SIM mode of fraction F2 from an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC  
 44 fractionation system. 27: B[b]F, 29: B[a]P, IS-1: Phe-D<sub>10</sub>, IS-2: Pyr-D<sub>10</sub>, IS-3: B[a]A-D<sub>12</sub>, IS-4: B[a]P-D<sub>12</sub>, IS-5: B[ghi]p-D<sub>12</sub>, IS-6: DB[a,i]P-D<sub>14</sub>.

45



46


47 **Figure S5.** GC/MS chromatogram obtained in SIM mode of fraction F3 from an injection of a Stockholm air PM extract on the B[a]P-specific online HPLC  
 48 fractionation system. 30: B[e]P, 29: B[a]P, 31: Per, 32: I[1,2,3-cd]F, 33: I[1,2,3-cd]P, 34: DB[a,h]A, 35: Pic, 36: B[ghi]p, 37: DB[a,l]P, 38: DB[a,e]P, 39: Cor,  
 49 40: DB[a,i]P, 41: DB[a,h]P, IS-1: Phe-D<sub>10</sub>, IS-2: Pyr-D<sub>10</sub>, IS-3: B[a]A-D<sub>12</sub>, IS-4: B[a]P-D<sub>12</sub>, IS-5: B[ghi]p-D<sub>12</sub>, IS-6: DB[a,i]P-D<sub>14</sub>.




50

51 **Figure S6.** GC/MS chromatogram obtained in SIM mode of waste W1 and W2 from an injection of a Stockholm air PM extract on the B[a]P-specific online  
52 HPLC fractionation system.

53 S-5. GC/MS chromatograms of blank and hexane injections on the B[a]P-specific online HPLC fractionation system




54  
55 **Figure S7.** GC/MS chromatograms obtained in SIM mode of fractions F1, F2, F3, W1 and W2 from a blank  
56 injection on the B[a]P-specific online HPLC fractionation system.



57

58 **Figure S8.** GC/MS chromatograms obtained in SIM mode of fractions F1, F2, F3, W1 and W2 from a hexane  
59 injection on the B[a]P-specific online HPLC fractionation system.

60 **S-6.** GC/MS full scan mass spectra of the system peaks from the B[a]P-specific online HPLC fractionation system





61

62 **Figure S9.** Full scan mass spectra of system peaks at (A) 29.88 and (B) 30.15 min.

63 **Table S3.** Analytical method for the GC/MS analysis


Name	Condition
<b>GC</b>	Agilent 6890N (Agilent Technologies, Palo Alto, CA, USA)
Injection	1 $\mu$ L, splitless at 300 $^{\circ}$ C
Carrier gas	He
Column	J & W DB-5ms 30 m $\times$ 0.25 mm, 0.25 $\mu$ m phase
Oven	60 $^{\circ}$ C for 2 min, 20 $^{\circ}$ C/min to 200 $^{\circ}$ C, 10 $^{\circ}$ C /min to 300 $^{\circ}$ C for 10 min
<b>MS</b>	MSD 5975C (Agilent Technologies, Palo Alto, CA, USA)
	Electron ionization at 70 eV, full scan mode ( $m/z$ 50–450), solvent delay 6.00 min, quadrupole at 150 $^{\circ}$ C, ion source at 300 $^{\circ}$ C



64

65 **Figure S10.** (A) GC/MS chromatogram of the eluate from the B[a]P-specific online HPLC fractionation system  
 66 with (B) full scan mass spectra from each system peak.

67



68

69 **Figure S11.** Stationary phase structure of the PYE column [1].




70 *S-7. Biological tests of eluates from the B[a]P-specific online HPLC fractionation system*

71 **Cell culture and exposure** Human hepatocellular carcinoma cells (HepG2) were obtained from the American Type  
72 Culture Collection (Rockville, MD, USA) and cultured as described previously [2]. For cell viability experiments,  
73 cells were seeded at  $5 \times 10^4$  cells/mL in 24-well plates whereas for Western blot experiments, cells were seeded at  $2$   
74  $\times 10^5$  cells/mL in 6-well plates. Cells were cultured for 72 hours, then subsequently exposed to solvent control  
75 (0.1% DMSO), positive control (3  $\mu$ M B[a]P), or eluates from the HPLC system for 24 hours. Eluates collected  
76 from the B[a]P-specific online HPLC fractionation system were evaporated and re-dissolved in DMSO for the tests.  
77 The same amount of neat hexane was taken and prepared in the same way as method blanks.


78 **Cell viability** Cell viability was determined by MTT assay as described previously [2]. Briefly, after exposure  
79 human HepG2 cells were incubated with 0.5 mg/mL MTT for 4 hours, then washed and formazan crystals were  
80 solubilised in DMSO. Optical density was measured at a wavelength of 590 nm. Data are expressed as percentage of  
81 control.

82 **Western blotting** Protein levels in cells were determined by Western blotting as described previously [2]. Briefly,  
83 cell protein samples were subjected to SDS-PAGE and thereafter blotted onto PVDF membrane (Bio-Rad, Hercules,  
84 CA). Protein bands were subsequently probed using antibodies. Cell Signaling Technology (Beverly, MA, USA)  
85 provided the following antibodies: Chk1 phosphorylated at Ser317 and H2AX phosphorylated at Ser139. Antibodies  
86 against Cdk2 (M-2), p53 (DO-1) and secondary anti-rabbit and anti-mouse antibodies were obtained from Santa  
87 Cruz Biotechnology (Santa Cruz, CA, USA). Signals were detected using Luminol Western Blotting reagent (Santa  
88 Cruz Biotechnology, Santa Cruz, CA, USA). Densitometric analysis was performed using ImageJ software version  
89 1.48f (National Institute of Health, USA).

90




91



92

93 **Figure S12.** Eluates from the B[a]P-specific online HPLC fractionation system do not induce DNA damage  
 94 signalling and are not cytotoxic in HepG2 cells. Phosphorylation of Chk1 at Ser317 (A, B) and H2AX at Ser139 (A,  
 95 C) and levels of p53 protein (A, D) were assessed by Western blotting after 24 hours. Cell viability was assessed by  
 96 MTT assay after 24 hours (E). Dotted line represents control levels. Experiments were performed in triplicate and  
 97 data points represent means  $\pm$  standard deviations. Where not shown, error lies within the data points. \*  $p < 0.05$  as  
 98 compared with control levels by one-way repeated measures ANOVA.



100  
 101  
 102 **Figure S13.** GC/MS chromatograms obtained in full scan mode of blank, non-fractionated and fractionated (F1, F2  
 103 and F3) samples. 1: B[b]F, 2: B[k]F, 3: B[j]F, 4: B[a]F, 5: Unknown ( $m/z$  268), 6: Unknown ( $m/z$  266), 7: B[e]P, 8:  
 104 B[a]P-D<sub>12</sub>, 9: B[a]P, 10: Unknown ( $m/z$  266), 11: Per. Chromatograms were background subtracted and mass data  
 105 was read from the apex of each peak.

106 *S-9. PAH concentrations in Stockholm air PM sample*

107 **Table S4.** PAH concentrations (pg/m<sup>3</sup>) in Stockholm air PM sample. Data obtained in the present study is presented  
 108 as a mean of three replicates ± standard deviation.

Name	Mean ± SD (pg/m <sup>3</sup> )	Name	Mean ± SD (pg/m <sup>3</sup> )
Phe	226 ± 3	B[a]A	304 ± 5
Ant	27.5 ± 0.5	3-MChr	27.1 ± 0.7
3-MPhe	20.2 ± 0.1	2-MChr	48.2 ± 0.8
2-MPhe	29.0 ± 0.6	6-MChr	31.7 ± 0.7
2-MAnt	5.26 ± 0.05	1-MChr	51.2 ± 0.8
9-MPhe	15.3 ± 0.2	B[b]F	412 ± 7
1-MPhe	34.5 ± 0.5	B[k]F	204 ± 5
4 <i>H</i> -CPP	41.0 ± 0.1	B[e]P	320 ± 5
2-PN	26.1 ± 0.4	B[a]P	297 ± 8
3,6-DMPhe	1.16 ± 0.05	Per	47.6 ± 0.8
3,9-DMPhe	5.94 ± 0.20	I[1,2,3- <i>cd</i> ]F	31.7 ± 0.9
Flu	345 ± 7	I[1,2,3- <i>cd</i> ]P	230 ± 4
Pyr	390 ± 8	DB[ <i>a,h</i> ]A	36.7 ± 0.4
1-MFlu	49.9 ± 1.1	Pic	38.2 ± 0.8
B[ <i>a</i> ]f	48.8 ± 0.7	B[ <i>ghi</i> ]p	352 ± 8
B[ <i>b</i> ]f	29.2 ± 0.3	DB[ <i>a,l</i> ]P	3.07 ± 0.17
2-MPyr	24.5 ± 0.6	DB[ <i>a,e</i> ]P	37.5 ± 0.8
4-MPyr	34.5 ± 0.4	Cor	204 ± 1

1-MPyr	40.7 ± 0.3	DB[a,i]P	11.3 ± 0.6
B[ghi]F	180 ± 5	DB[a,h]P	5.54 ± 0.38
B[c]Phe	59.8 ± 2.0		

---

109

110

111 **Reference**

112 [1] K. Kimata, K. Hosoya, H. Kuroki, N. Tanaka, J.R. Barr, P.C. McClure, D.G. Patterson Jr., E. Jakobsson, A.  
113 Bergman, J. Chromatogr. A 786 (1997) 237.

114 [2] I.W.H. Jarvis, C. Bergvall, M. Bottai, R. Westerholm, U. Stenius, K. Dreij, Toxicol. Appl. Pharmacol. 266  
115 (2013) 408.