Provided by NASA Technical Reports Server

Brief 69-10319

NASA TECH BRIEF



NASA Tech Briefs are issued to summarize specific innovations derived from the U.S. space program, to encourage their commercial application. Copies are available to the public at 15 cents each from the Clearinghouse for Federal Scientific and Technical Information, Springfield, Virginia 22151.

Miniature Oxygen Resuscitator

The problem:

Development of a small, lightweight, portable resuscitator for use during evacuation of patients to medical facilities.

The solution:

A miniature resuscitation system weighing 4 lb with a 10-minute supply of oxygen or 8 lb with a 30minute supply. A carrying case contains a modified resuscitator head, a cylinder of oxygen, a two-stage oxygen regulator, a low-pressure tube, and a mask for mouth and nose.

How it's done:

The standard resuscitator has been modified by elimination of the mixing of air with the oxygen; 100% oxygen is now administered. The unit is cycled to give a positive patient pressure of 14 ± 2 mm-Hg and a negative patient pressure of 9 ± 2 mm-Hg when the constant source pressure is 22 ± 2 lb/inch² (gauge). Increase (or decrease) by 1 lb/inch2 in the supply pressure increases (or decreases) the patientpressure range by about 2 mm-Hg. A safety valve limits the patient's positive pressure to 21 ± 2 mm-Hg and his negative pressure to 17 ± 2 mm-Hg; a manual override permits a continuous flow at the valve's setting. The flow of oxygen is generally constant throughout the positive and negative phases of the resuscitator's cycle, falling between 11 and 13 liter/min as measured at the 22-lb/inch2 supply line.

The 50-inch3 oxygen cylinder accords with MIL-C-7905 and MS-26545.

The regulator consists of a two-stage pressure reducer and a shut-off valve directly mounted on the tank. With an input pressure ranging from 50 to 2,000 lb/inch² the reducer delivers a steady 22 ± 3 lb/inch2 (gauge). The quick-turn (maximum of 180°) valve, on the high pressure side, isolates the oxygen supply from the rest of the system. The high-pressure side has an L-2 type of pressure gauge, a standard aircraft-type filler port, and a high-pressure-blowout port set at 2,600 to 2,800 lb/inch2. Only special lubricants may be used in the system.

The resuscitator head is permanently attached to the pressure regulator by an 18-inch flexible tube made of poly(vinyl chloride) reinforced with nylon. The tube has an internal diameter of 3/16 inch and a burst pressure of $\geq 1,000 \text{ lb/inch}^2$.

The colorless and transparent mask, equipped with an inflatable cushion for airtight sealing around the patient's mouth and nose, is directly attached to the patient-outlet nipple of the resuscitator head.

The resuscitator should be useful at first-aid stations at beaches, swimming pools, factories, and hospitals, and for rescue crews and ambulances.

Note:

No further documentation is available. Inquiries may be directed to:

> Technology Utilization Officer Kennedy Space Center Kennedy Space Center, Florida 32899 Reference: B69-10319

> > (continued overleaf)

This document was prepared under the sponsorship of the National Aeronautics and Space Administration. Neither the United States Government nor any person acting on behalf of the United States Government assumes any liability resulting from the use of the information contained in this document, or warrants that such use will be free from privately owned rights.

Patent status:

No patent action is contemplated by NASA.

Source: (1) J. T. Teegen
Kennedy Space Center
(2) Henry Waddell of
North American Rockwell Corporation
(3) George Johnson of
The Boeing Company
both under contract to
Kennedy Space Center
(KSC-10398)