



Light stimulated healing of chronic wounds by use of optical waveguides and light management microstructures

David Kallweit, Nicolas Glaser, and Jan Mayer

CSEM SA, Switzerland





UNIVERSITÄT
HEIDELBERG
ZUKUNFT
SEIT 1386



Project Facts


- Duration: 02/15 – 01/18
- Project Costs: 3 Mio €
- Co-funded by the European Union as a H2020 Research & Innovation Action
- www.medilight-project.eu



The project MEDILIGHT receives funding from the European Union's Horizon 2020 research and innovation programme under Grant Agreement No. 644267.

Health Facts...

- Chronic wounds caused by diabetes
- 170 Mio people worldwide affected
- Costs are 40 billion € per year
- Difficult to treat
- Therapeutic effect of visible light has been proven
- Red light (620-750nm):
growth of keratinocytes and fibroblasts
in deeper layers of the skin
- Blue light (450-495nm):
antibacterial effects at the skin surface




<http://americanfoot.com>



Medilight Device

consists of two parts:

- A) electronic module
- B) wound dressing (**disposable**)




A: Electronic Module

- Light sources
- Beam shaping optics
- Controller & Data Acquisition
- Data transmission unit
- Rechargeable battery


B: Wound Dressing

- Flexible large area optical waveguide
- Light management structures:
 - diffraction gratings (incoupling blue and red)
 - diffusor (homogenization, outcoupling)
- Integrated sensors
 - temperature
 - blood oxygenation


Waveguide Design



Diffraction grating (light incoupling)



Diffusing microstructure (light outcoupling and homogenization)




Selection of Waveguide Material




The following properties are desired:

- high optical clearness
- low scattering
- good surface quality
- highly flexible / bendable
- not sensitive to humidity (swelling, etc.)
- easy to be embossed / injection molded



- Different materials were tested:
PET, PMMA, PC, COP (different types of each)
- **Cycloc olefin polymer (COP) has been selected**



Surface of a measured PMMA sample



Incoupling Grating (Blue) – Design by RCWA





- **Theoretical maximum diffraction efficiency in T-1 is 50% (at 440-465nm)**
- **25% diffraction efficiency were measured (fabricated profile deviated from simulation model)**

Measurement


Angle (deg)	T-1
11	24.6%
12	25.0%
13	25.3%
14	24.8%
15	24.5%

New Design with Increased ZnS Coating



- ✓ Theoretical maximum was increased up to 70% at 10-20°
- ✓ Angle dependence is reduced (higher acceptance angle)



Fabricated Grating (for Red)



Measurement


Angle (deg)	T1
-5	0.0%
0	0.6%
3	7.7%
7	20.1%
8	36.6%
9	60.9%
10	61.0%
11	53.3%
12	45.4%
13	38.2%
14	35.3%
15	30.4%

61% efficiency measured




Further Improvement: Angular Evaporation of ZnS

- Provides a maximum efficiency of T-1 diffraction order at **normal incidence** ($\alpha_{in} = 0^\circ$) with a **high angular tolerance**






Fabrication is ongoing and the compliance with the calculations still needs to be evaluated




Homogenizing Diffusor

Different diffusor designs were developed and tested:




Challenges & Outlook

1. **Waveguide efficiency:** Currently we reach 15-20%. To be improved.
2. **Homogenous illumination** of the wound area
 - Suitable patterning of the diffusor matching the used light sources and the waveguide geometry
3. **Waveguide cladding:** nanoporous layer under development
4. **Waveguide shape:** wholes, «fingers» allowing the wound to breath
5. **Integration of sensors** (temperature, blood oxygenation)



Thank you for your attention.



Basel, Switzerland