

Concentration techniques for genotoxicity testing of effluents

M.. Adahchour*, R. Vreuls* and B. van Hattum

Institute for Environmental Studies (IVM), Faculty of Earth and Life Sciences

* Department of Analytical Chemistry and Applied Spectroscopy (ACAS),
Faculty of Exact Sciences.

Report number E01/08

December 2001

Study commissioned by the National Institute for Inland Water Management and Waste Water (RIZA). Contract nr. 39801/EMP.

IVM

Institute for Environmental Studies
Fac. Earth and Life Sciences (FALW)
Vrije Universiteit
De Boelelaan 1115
1081 HV Amsterdam
The Netherlands

Tel. ++31-20-4449 555

Fax. ++31-20-4449 553

E-mail: secr@ivm.vu.nl

ACAS

Dept. Analytical Chemistry and Applied Spectroscopy
Fac. Exact Sciences (FEW)

Vrije Universiteit
De Boelelaan 1083

Tel. 020-4447544

Fax. 020-4447543

E-mail: aac@chem.vu.nl

Copyright © 2001, Institute for Environmental Studies, ACAS, RIZA

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise without the prior written permission of the copyright holders.

Contents

Preface	iii
Executive summary	v
1. Introduction	1
2. Concentration techniques	3
2.1 Introduction	3
2.2 General survey and methods classification	3
2.3 Motivation for the selection of SPE	5
2.4 Sorbents for SPE	5
2.5 Breakthrough volumes <i>versus</i> Log K_{ow}	8
2.6 Compound classification	10
3. Physico-chemical properties of responsive compounds	13
3.1 Introduction	13
3.2 Estimation of K_{ow}	13
3.3 Estimation of Henry's Law Constants	14
3.4 Distribution of hydrophobicity	15
3.5 Distribution of Henry's Law Constants	17
4. Selection of suitable SPE methods	21
4.1 Introduction	21
4.2 Application of selection criteria	22
4.3 Constraints from the UmuC-test protocol	24
4.4 Various SPE approaches	24
4.5 Proposal	27
5. Discussion and concluding statements	29
References	31
Annex section	41
Annex 1 Factsheets	43
Annex 2 List of responsive compounds	69
Annex 3 Predicted and measured values for K_{ow}	77
Annex 4 Predicted and experimentally determined Henry Law constants	87
Annex 5 Detailed description of the selected SPE approaches	97

Preface

The present study was part of a preparatory programme on the introduction of genotoxicity testing and other acute and chronic bioassays in the integrated framework for biological and chemical evaluation of wastewater effluents (WEER; whole effluent environmental risk procedure; Totaal-Effluentbeoordeling) in the Netherlands. Based on previously made choices for the application of the UmuC-test, the main objective of the present study was to evaluate suitable preconcentration techniques, to define proper selection criteria, and to develop a proposal for a cost-effective preconcentration technique.

The project was commissioned by the National Institute for Inland Water Management and Waste Water (RIZA) of Rijkswaterstaat (Contract nr. 39801/EMP) under supervision of Dr Gert-Jan de Maagd. The project executed in a 10 week period by Drs. Mohamed Adahchour and Dr René Vreuls of the Department of Analytical Chemistry and Applied Spectroscopy (ACAS), Faculty of Exact Sciences (FEW), and Dr Bert van Hattum of the Institute for Environmental Studies (IVM), Faculty of Earth and Life Sciences (FALW), who took care of the overall project management. The authors acknowledge the support of Serge Rotteveel Msc and Ineke Blanckesteijn Bsc (both from RIZA), Dr Freek Ariese (ACAS) and Dr Juliette Legler (IVM) for their part in the review of the project and usefull contributions.

Executive summary

In the present study, preconcentration methods were evaluated for the extraction and concentration of genotoxic compounds from effluents and surface water. The project was commissioned within the framework of the development of a whole effluent environmental risks procedure (WEER) by the National Institute for Inland Water Management and Waste Water Treatment (RIZA). Based on previously made choices for the application of the UmuC-test, the main objective of the present study was to evaluate suitable preconcentration techniques, to define proper selection criteria, and to develop a proposal for a cost-effective preconcentration technique.

To date, over 500 different organic and inorganic compounds have been tested for genotoxicity (mutagenicity, DNA damage, and carcinogenicity) with a large variety (appr. 200), of test systems, ranging from prokaryotic tests to mammalian species. Based on a selection of 156 compounds, responsive in the UmuC-test, physico-chemical properties, such as the *n*-octanol/water partitioning coefficient (K_{ow}) and the Henry Law Constant (H in $\text{Pa}\cdot\text{m}^3/\text{Mol}$), were evaluated and classified. These parameters are considered as key factors describing the fate during sample-handling procedures and were important criteria for the selection of potentially suitable methods.

Measured values were derived from the Environmental Fate Data Base (EFDB, Syracuse Comp., USA). Missing values were estimated from structural parameters with the KOWWIN (fragment method) and HENRYWIN (bond contribution) modules in the EPIWIN suite from US-EPA. Genotoxic compounds differ largely in hydrophobicity. Approximately 22 % was classified as hydrophilic ($\text{Log } K_{ow} < 0$), 16% as polar ($\text{Log } K_{ow} 0-1$), 36% as medium-polar ($\text{Log } K_{ow} 1-3$) and 26% as hydrophobic ($\text{Log } K_{ow} > 3$). The majority of the compounds (56%) was classified as relatively non-volatile from water ($H < 10^{-3} \text{ Pa}\cdot\text{m}^3/\text{Mol}$), approximately 29% as intermediate ($10^{-3} < H < 1$), and 14% as prone to volatile losses during sample-handling ($H > 1$), with 7% above $H = 10$.

An overview is presented of preconcentration techniques together with a more in-depth treatment of currently available SPE/SSPE techniques (solid phase extraction/ sequential solid phase extractions), some of the critical performance characteristics, such as breakthrough volume, and a provisional classification of compounds as used in this study. Basic properties of a selection of 22 SPE/SSPE procedures, derived from over 80 recent publications, are presented in the form of detailed factsheets in the Annex section. Especially the breakthrough volume appeared to be a critical factor. It was shown that most commonly used hydrophobic SPE materials exhibit breakthrough losses in the $\text{Log } K_{ow}$ interval of 0-3, which coincides with $\text{Log } K_{ow}$ of more than 50% of the UmuC-responsive compounds.

Based on the properties of responsive compounds and constraints from the UmuC-testing protocol a selection of 5 SPE/SSPE for surface water trace-enrichment methods was evaluated in more detail: C-18, Oasis HLB, Lichlolut EN, Isolut ENV+ and GCB based methods. For the four different hydrophobicity intervals the performance of the sorbents were evaluated with respect to: breakthrough volume, recovery, enrichment factor, and ease of drying.

A combined SSPE method, including C-18 followed by GCB was proposed as the most comprehensive method, expected to cover a wide range ($\text{Log}K_{ow}$:-2 to 6). Expected losses are likely to occur in the hydrophilic range ($\text{Log}K_{ow} < -2$; matching with 10% of the UmuC-responsive compounds) and approximately 5% (of the UmuC-responsive compounds) in the super-hydrophobic range ($\text{Log}K_{ow} > 6$). Further losses due to evaporation during sample-handling may be present for approximately 14% ($n=22$) of the compounds, if a conservative cut-off of $H=0$ is chosen. A conservative estimate shows that approximately 70% of the UmuC-responsive compounds is expected to be extracted quantitatively with the proposed method. If the volatile fraction is of importance, than the additional application of purge-and-trap concentration techniques prior to genotoxicity testing is recommended.

The proposed C18/GCB SSPE procedure can be easily adapted to yield separate fractions (with different hydrophobicity class intervals) for further toxicity identity evaluation (TIE) studies. Several of the discussed sorbent materials can also be incorporated in a HPLC setting, for on-line separation in combination with automated fraction-collecting systems, dosing to e.g. microwell-plates which allow high-throughput testing with the UmuC-protocol.

With the current state-of-the-art in SPE/SSPE techniques it is not possible to design a single generic method, suitable for routine (waste)water screening, which is fully comprehensive. If compounds with extremely high or low hydrophobicity are of importance, or the attention is dedicated to highly volatile compounds, special dedicated methods need to be applied. A comprehensive review of the literature revealed that three parameters are the most critical for the recovery in the final extract: volatility, breakthrough and desorption characteristics.

The category of inorganic compounds was not included in this study. Although in theory adding specific ion-exchange resins in sequence with the other SSPE may result in the trapping of metallic compounds of interest (e.g. selenium, chromium compounds), it is not possible to isolate the redox-state and chemical species, as applied in the original genotoxicity tests with dosages by far exceeding normal levels in effluents and surface water. We recommend not to include inorganics in the generic screening approach. In cases where genotoxicity of inorganics may be of importance (e.g. effluents of specific industrial sectors), dedicated studies are considered as the most cost-effective approach.

1. Introduction

The present study was commissioned within the framework of the development of a whole effluent environmental risk procedure (WEER) by the National Institute for Inland Water Management and Waste Water Treatment (RIZA). This integrated approach will combine ecotoxicological tests (acute- and chronic aquatic toxicity, bioaccumulation, genotoxicity and other *in-vitro* bioassays) in combination with traditional chemical approaches and dedicated toxicity identity evaluation (TIE) studies, to explain observed bioassay responses (Tonkes et al., 1998).

To date, over 500 different organic and inorganic compounds have been tested for genotoxicity (mutagenicity, DNA damage, carcinogenicity) with a large variety (appr. 200; de Maagd and Tonkes, 2000), of test systems, ranging from prokaryotic (e.g. bacterial) tests to mammalian species. The databases hosted by IARC (monographs with classification; www.iarc.lyon.fr) and the US-EPA/IARC (GAP database; www.epa.gov/gapdb) provide comprehensive and continuously updated information. Many of the genotoxic organic compounds are emitted from industrial and non-point sources and may end up in aquatic environments, depending on persistence, volatility and affinity to sediments and biota. Usually concentrations in surface waters and - to a lesser extent - effluents are below the no-effect ranges of genotoxicity tests, and some form of concentration may be required to elicit a response. A wealth of different preconcentration techniques has been applied for aquatic genotoxicity testing, and since long the need for validation and standardisation of sample-prep has been advocated (Stahl, 1991).

This study, addresses the selection of suitable concentration techniques for surface and wastewater, and covering a wide range of matrix and compound properties. From previous studies it is known that few methods exist that extract and concentrate all relevant compounds to an equal extent (Stahl, 1991; de Maagd, 2001) and that volatile or highly polar genotoxic compounds may be underrepresented in concentrated extracts (Vargas et al., 1995). Some authors have reported on variation in genotoxic response of wastewater extracted and concentrated with different methods (Gauthier, 1993). Filipic and Toman (1996) stated that the genotoxic response of concentrated samples cannot be extrapolated directly to the response of the original sample.

In previous studies, various genotoxicity tests were evaluated for their suitability for screening of wastewater and surface water (Klamer et al., 1997; De Maagd et al., 1999; De Maagd & Tonkes, 2000) and the UmuC-test (Reiferscheid et al., 1996; Hamer et al., 2000) was assigned as one of the candidate prokaryotic tests for primary DNA-damage in the WEER programme. Before implementation various issues need to be addressed, ranging from optimisation and standardisation of sample-handling procedures to protocol development, validation and intercalibration studies (de Maagd et al., 1999).

The main goal for the present satellite-project was to review existing preconcentration methods for (waste)water and to evaluate and select the most promising methods for application in the UmuC-test procedure. As the currently known set of genotoxic compounds exhibit order of magnitude variations in their physico-chemical properties, special attention was given to the comprehensiveness for the full range of compounds.

The following objectives were defined:

- to review the state-of-the-art of currently known concentration techniques suitable for combination with aquatic genotoxicity testing;
- to inventorise and evaluate some key physico-chemical properties of responsive compounds: the *n*-octanol/water partitioning coefficient (K_{ow}) and the Henry's Law Constant (H);
- to formulate and apply selection criteria to identify the most promising methods;
- to discuss gaps in knowledge and relevant developments in the field of sample-prep in the near future.

In Chapter-2 an overview of preconcentration techniques is presented, together with a more in-depth treatment of currently available SPE/SSPE techniques (solid phase extraction/ sequential solid phase extractions), some of the critical performance characteristics, such as breakthrough volume, and a provisional classification of compounds as used in this study. Basic properties of a selection of more than 30 SPE/SSPE procedures, derived from over 80 recent publications, is presented in the form of detailed factsheets in Annex-1. Estimated and experimentally determined properties (K_{ow} and H) of UmuC-responsive compounds are presented in Annex 2-4 and are evaluated and analysed in Chapter-3. Based on the properties of responsive compounds and constraints from the UmuC-testing protocol a selection of 5 SPE/SSPE methods was evaluated in more detail. The development and application of selection criteria, resulting in a proposed combination of SPE methods is presented in Chapter-4. After a brief discussion of gaps in knowledge and future development in the field (Chapter-5), conclusions and recommendations were summarized in Chapter-6.

2. Concentration techniques

2.1 Introduction

Many organic contaminants are present in surface waters in trace amounts at the $\mu\text{g/L}$ level and often below. As a result, very low detection limits are required for direct monitoring. To date, several hundred organic contaminants have been classified as genotoxic (Reifferscheid et al., 1996). Because of the low concentrations, detecting them without some sort of concentration step would be nearly impossible. For that purpose various concentration techniques have been used, including lyophilization, reverse osmosis, liquid-liquid extraction, solid-phase extraction, etc. In general, sample preparation is the most tedious and time-consuming step. It is also a potential contributor to the imprecision and inaccuracy of the overall analysis, especially when the sample (extract) is used for toxicity testing (Stahl et al., 1991; de Maagd et al., 1999). Therefore, the choice of the method(s) for the concentration and isolation of the various types and classes of contaminants (toxicants) from water samples is complicated and depends on the analyte characteristics (polarity, ionic character, stability) and other parameters related to the applied method.

In this section of the report, a general discussion and classification will be presented of the concentration techniques that were most widely used in the last decade, which will help to define the selection criteria presented in section 4.

Articles were selected from the following three main databases: chemical abstracts, web of science and science direct. The search strategy was based on the following keywords: concentration techniques, solid-phase extraction, liquid-liquid extraction, water samples, polar compounds, isolation techniques, toxicity test, toxic compounds, genotoxic compounds, polar organic contaminants, effluent, waste water, ground water, tap water, UmuC-test. Beside these keywords a search was conducted on specific groups of compounds, according to the selected list of genotoxic analytes presented in section 3.

2.2 General survey and methods classification

As mentioned above, there are still many problems encountered in the isolation and preconcentration of organic compounds in water samples. Concentrating the sample may produce artefacts, sometimes leading to both false positive/negative results. Different concentration techniques have been used in the past to overcome such problems:

Lyophilization. Lyophilization, the removal of water by vacuum sublimation of ice, is a method for concentration in-volatile organic compounds in aqueous solution and for preserving biological samples. The essential steps are that the sample is first frozen and then placed under high vacuum. The heat absorbed by the sample from its surroundings causes the ice to sublime, which is then re-condensed on a large surface condenser held at cryogenic temperatures. This technique, also known as freeze frying, has several drawbacks. First, salts, bacteria, and other materials in the sample are not removed and may interfere with the genotoxicity test. Second, only small samples can be processed

(generally < 100 ml) which reduces the concentration factor achievable. The advantage is that an aqueous solution is obtained, which is fully compatible with the UmuC-test.

Reverse osmosis and evaporation. Reverse osmosis is a separation process in which (large) molecules are separated from solution by filtration through membranes. The operating pressure must exceed the natural osmotic pressure for the system resulting in the movement of solvent; usually water, from the solution of high analyte concentration to that of low analyte concentration.

Next to the drawbacks mentioned for lyophilization, reverse osmosis and evaporation can be inefficient and time consuming. These approaches have the potential to markedly alter the composition of a sample through sorption, leaching, evaporation, or degradation of the contaminants.

Freeze concentration. Unlike lyophilization, freeze concentration has been used for the concentration of aqueous solutions of organic volatiles and substances that are heat labile. For successful results, the contact layers between the liquid and solid phases should be continuously disturbed by stirring or shaking and part of the solution should remain unfrozen at the end of the concentrating procedure.

The most important weakness of the freeze concentration method lies in its completeness. Increasing the total dissolved solids concentration of highly mineralised samples could potentially elicit toxic effects from ion imbalance. Conductivities greater than 2,000 $\mu\text{S}/\text{cm}$ may be high enough to adversely affect freshwater species.

Liquid-liquid extraction (LLE). LLE has remained the preferred technique for the preparation of liquid samples for several years, especially in the environmental field, and is still widely used in standard and official methods. Its process is based on distribution of analytes between two immiscible solvents in which analytes have different solubilities. The efficiency of an extracting solvent depends on the affinity of a solute for the extracting solvent, the phase ratio and the number of extractions.

LLE with traditional organic solvents is more time-consuming than the other methods discussed, but provides qualitative information about the classes of organic compounds causing the genotoxic activity, i.e., acid, neutral, or basic compounds. However, the introduction of survey lists (see Chapter 3) containing (very) polar analytes such as some degradation products of organic micro-pollutants has also pointed out the need for alternative methods to LLE. Many polar analytes are partly soluble in water and cannot be extracted with good recoveries, whatever the organic solvent selected. Other disadvantages of this technique, e.g. matrix interferences, emulsion formation, the use of large volumes of hazardous solvents, have troubled the analyst.

Solid-phase extraction. SPE is an extraction and enrichment technique and its process, to a first approximation, can be considered as a simple liquid chromatographic process. The sorbent is the stationary phase and the mobile phase is the water constituting the aqueous sample during the extraction step and the organic solvent during the elution step. Compounds that do not elute with the water are trapped on the sorbent during the percolation step. High enrichment factors are obtained when analytes are strongly retained by the sorbent in the presence of water, allowing the percolation of a large volume of sample and when there is a low retention when eluting by organic solvents.

Solid-phase extraction (SPE) has been developed as an alternative to LLE (and other extraction techniques) and can (partially) overcome its drawbacks. Motivation of this choice is presented in some detail in the next section.

2.3 Motivation for the selection of SPE

The increasing popularity of SPE is resulting not only from its obvious advantages over LLE (e.g., minimised consumption of organic solvents, no emulsion formation, reduced contact for the analyst with potentially toxic substances), but also from the gradual refinements of this procedure, which minimise its original drawbacks. A wide choice of newly developed material, belonging mostly to one of the three major groups of sorbents (i.e., bonded silicas, polymers and carbon materials), makes SPE a suitable tool to cope with an increasing variability of modern organic contaminants. The theory and method development of SPE for water analysis were recently reviewed (Poole et al., 2000; Liska 2000; Hennion 1999; Hennion 2000; Huck et al., 2000; Leon-Gonzalez et al., 2000 and Pichonm, 2000). The applicability of SPE for a specific group of analytes has also been reviewed, including phenols (Rodriguez et al., 2000 and Bruzzoniti et al., 2000), acidic herbicides (Wells et al., 2000) and pesticides, triazines and degradation products (Sabik et al., 2000).

Different laboratories have studied comparison of SPE to LLE extensively. Di Corcia *et al.* (1993) compared the recoveries of 20 polar pesticides obtained from 2 L of water using a 1-g graphitized carbon black (GCB) cartridge to those obtained with using a 1-g C₁₈ silica cartridge or by using LLE with three separate 120-ml portions of methylene chloride. SPE (GCB) was superior to LLE, with recoveries of 75-102% and 13-89% respectively. The recoveries obtained with 1-g C₁₈ were much lower (3-55%, except for isocarbamid, hexazinon and metoxuron with recoveries higher than 80%) because breakthrough (see paragraph 2.5 for more details) occurred for most of the pesticides with a sample volume far lower than 2 L. Tolgyessy et al. (1999) investigated the efficiency of pre-concentration of selected organic compounds from aqueous solutions on various SPE materials including bonded silica, polymers and carbon materials. Simultaneously, the potential of newly emerging SPE procedures was compared to the results of traditional LLE methods using various organic solvents. A group of 19 test analytes was selected so as to represent different classes of organic compounds, which may occur in waters. The results showed that most of the tested materials were suitable for sufficient preconcentration of a substantial part of the tested analytes. However, specific differences in the recovery of one or more analytes were found for almost each sorbent, even in the case when the materials had a similar composition. This behaviour clearly indicates the need of a thorough testing of the capabilities of any SPE material intended for the identification of (un)known organic micropollutants of a wide polarity range in waters.

2.4 Sorbents for SPE

The key problem when applying SPE remains the method development and the primary decision for the analysts is the choice of the sorbent that is able to solve their trace-analysis problem. Before selecting a sorbent for SPE, it is necessary to take into account some physico-chemical aspects (Chapter 3), such as the functional groups of the analytes, the nature of the bonded phase, the energetics of the interactions, the secondary in-

interactions between the sorbent and the components of the sample matrix, and the interactions between the analytes and the sample matrix. Sorbent-analyte interactions fall into three categories: non-polar, polar and ionic.

To date, typical SPE materials are modified silicas with C₈, C₁₈, CN and other groups, carbon blacks and styrene-divinylbenzene copolymers (PS-DVB). The most problematic compounds are the polar ones, for which breakthrough volumes are low with most of the sorbents mentioned. This has encouraged investigators to evaluate alternative sorbent materials, such as chemically modified polymeric resins, highly cross-linked polymers and others.

Modified silica. Since the retention mechanism is primarily governed by hydrophobic interactions between the analyte and the carbon moieties of the alkyl chain, a relation has been observed between the retention factor of the analyte and its *n*-octanol/water partition coefficient (K_{ow}). It was demonstrated that, with an amount of sorbent of 500 mg, a recovery in the range 90-100% for a percolated volume of 500 ml was found for compounds with a $\text{Log } K_w > 3$ (Hennion et al., 1998); K_w is the retention factor of the solute eluted by water (see Section 2.5). This guidance value is the basis of the choice of the sorbent. No special problems occur during the extraction of moderately or non-polar analytes. However, these sorbents are inadequate for solving the problems involved in isolating polar contaminants from large water volumes (Di Corcia et al., 1993). Apart from that, these sorbents suffer from chemical limitations, such as residual surface silanol groups, a narrow pH stability range and poor selectivity.

Graphitized carbon black (GCB). This material is characterised by a highly homogeneous and ordered structure and by a specific surface area of ca. 120 m²/g. It has been shown that organic compounds are adsorbed on a GCB surface more strongly than on a C₁₈ surface when analysing polar compounds in water (Di Corcia et al., 1993, Hennion et al., 1998). The retention mechanism differs from the one observed on C₁₈ silica or PS-DVB polymers due to its crystalline structure made of large sheets held together by weak Van der Waals forces (Hennion et al., 1995). Analytes are retained on GCB by both hydrophilic- and electronic-type interactions, so that non-polar but also very polar and water-soluble analytes have been shown to be strongly retained in aqueous mobile phases (Guenu et al., 1996a and b). A new GCB sorbent, Carbograph 4, was used to determine phenolic compounds in water and they yielded better recoveries for the most polar analytes than Carbograph 1 (Pocurull et al., 1996 and Di Corcia et al., 1996).

A drawback of the GCB sorbents are that they result in excessive retention (some analytes can even be irreversibly adsorbed) but this can be overcome by performing the elution in the back flush mode (Di Corcia et al., 1996). Its mechanical property is also poor; which results in low-pressure resistance to be used in liquid chromatography. Porous immobilised graphitic carbon (PGC) is a new carbon-based sorbent in which the graphite is immobilized on a silica structure, and this is why PGC is more stable than GCB. The use of PGC was investigated for the trace enrichment of polar phenolic compounds and compared with polymeric sorbents. Except for 4-chloro-2-aminophenol, the breakthrough volumes and detection limits were worse with PGC than those obtained using the polymeric sorbents (Puig et al., 1996).

-δ interac-

tions with unsaturated analytes. The efficiency of the polymeric sorbents depends on various physico-chemical parameters such as particle size, surface area, pore diameter, pore volume, degree of cross-linking and particle size distribution (Puig et al., 1996 and Hennion et al., 1998). Resins with high specific surface areas in the range 700 to 1200 m²/g, are now available in disposable cartridges and its effectiveness has been demonstrated by Hennion et al. (1998).

In recent years, chemically modified resins containing different polar functional groups have been developed and used in the SPE of polar compounds from environmental waters. These modified resins have excellent hydrophilicity and they also give higher recoveries than their unmodified analogues (Masque et al., 1997). This has been attributed to an increase in surface polarity, which enables the aqueous sample to make better contact with the resins' surface. They have also been compared with other SPE materials such as PLRP-S and GCB for the preconcentration of pesticides and phenolic compounds, and they yielded higher breakthrough volumes for the most polar compounds than PLRP-S and GCB (Masque et al., 1997 and Masque et al., 1997). In the last few years, several new highly cross-linked polymers, such as Envi-Chrom P, LiChrolut EN, Styrosorb, Isolut ENV and HYSphere-1, have become available. These highly cross-linked PS-DVB sorbents are the sorbents to be selected for the extraction of very polar analytes when large sample volumes are required. Many examples can be found in the recent literature (Castillo et al., 2001; Steen et al., 2001; Lopez de Alda et al., 2001 and la Farre et al., 2001a-d).

One limitation of both reversed-phase silica and polymer sorbents is that they must be conditioned with a wetting solvent and remain wetted until sample application. One new patented sorbent has been recently introduced in the market. It is the so-called hydrophilic-lipophilic balanced sorbent (Oasis HLB from Waters), which is a copolymer of divinylbenzene and *N*-vinylpyrrolidone (Bobeldijk et al., 2001a and Bobeldijk et al., 2001b). It is presented as the "universal" extraction sorbent since it is capable of extracting acidic, basic and neutral compounds, whether polar or non-polar.

This new hydrophilic-lipophilic polymer is (more) suitable for the handling of biological samples like serum. There is an "universal procedure" which consists of applying 1-ml of sample, washing with 1 ml of 5% methanol before eluting with 1 ml of methanol. Since the sorbent combines a high specific surface area with some polar hydrogen groups, it is easy to explain why there is no breakthrough of analytes with 1 ml of samples whatever the organic analyte, being polar, neutral or ionised, even when 5% methanol is added for the washing step. The new polymer is designed to extract an extensive spectrum of analytes, i.e., lipophilic, hydrophilic, acidic, basic and neutral with a single cartridge with a simplified procedure since no conditioning is required. But the universal procedure (apply 1 ml of sample, wash and elute) is limited to the handling of biological samples since the sample volume that is recommended is 1 ml. More work should be

done to look at the limitation in the extraction of polar analytes from relatively high (water) sample volumes as it was done for other highly cross-linked PS-DVB.

To guarantee the success of the SPE process for a particular application, it is of high importance to understand the effect of the most critical parameters on analyte recovery during the SPE procedure. In chapter 4, a selection criterion is proposed according the described goal in the Introduction to select (an) optimal method(s).

2.5 Breakthrough volumes versus Log K_{ow}

The breakthrough volume is the most important parameter for determining the suitability of a sampling device for isolating the analytes of interest. It represents the maximum sample volume, which can be percolated with a theoretical 100% recovery. In the initial sampling phase a sample of fixed concentration, and usually at a constant velocity, enters the sampling device. The analytes are quantitatively retained during the initial sampling phase by the sorbent, up to the point that the sample volume exceeds the retention capacity of the sorbent.

A typical example of the relation between breakthrough volume and Log K_{ow} of a group of pesticides (van Loon et al., 1996; Liska et al., 1992; Brouwer et al., 1991; Hennion et al., 1993) is shown in Figure 2.1. From this figure, it is clear that compounds with Log $K_{ow} < 1$ have a breakthrough volume smaller than 10 ml. For compounds with Log $K_{ow} > 3$, the breakthrough volume is higher than 100 ml. In the 1–3 Log K_{ow} range, breakthrough volume increases with Log K_{ow} . To recover the (very) polar compounds, a polar phase is needed i.e., GCB or PS-DVB with high specific surface areas.




Figure 2.1 Relation between Log K_{ow} and breakthrough volume for 65 pesticides and other compounds on polymeric phases (PRLP-S and PRP), L 10 mm x ID 2 mm. Adapted from: Van Loon et al. (1996), Liska et al. (1992), Brouwer et al. (1991) and Hennion et al. (1993).

The breakthrough volume can be calculated or estimated by the following well known empirical equation:

$$V_b = (V_{\text{geom}} \cdot \hat{\alpha}) (1 - 2/N^{0.5}) 10^{a \text{Log}K_{\text{ow}} - b}$$

where V_b is the breakthrough volume, V_{geom} is the geometric volume of the column or the cartridge, $\hat{\alpha}$ is the porosity of the sorbent, N is the number of theoretical plates, K_{ow} is the *n*-octanol/water coefficient of the solute, and a and b are the linear regression coefficient of the relationship between $\text{Log}K_{\text{ow}}$ and $\text{Log}K_w$ (retention factor of the solute eluted by water). With average values of density and porosity of the sorbent used in cartridges, the void volume V_m ($V_m = V_{\text{geom}} \cdot \hat{\alpha}$) of the cartridge can be calculated. For example, V_m for C_{18} silica is estimated as 0.12 ± 0.01 ml per 100 mg of sorbent. It is much more difficult to measure the efficiency of an SPE cartridge or that of an extraction disk, so that N has to be estimated.

The determination of breakthrough volumes, particularly by off-line methods, is time-consuming. It can be estimated from solute properties using models that require a minimal number of experimental measurements. The breakthrough curves can be easily be modeled according to the above equation, provided that values are known for porosity of the sorbent, number of plates and amount of sorbent (Hennion et al., 1998). These curves have been constructed with the hypothesis of an amount of 450 mg of sorbent, an average porosity of 0.70 and a number of plates of 20. The effect of the number of plates was shown to be negligible when modelling with the same conditions. The only strong effect is the amount of sorbent, but the relation between V_b and V_m is straightforward in the equation mentioned above.

This model has been validated by experimental measurements of recoveries for several analytes using extraction disk packed with C_{18} silica and PS-DVB sorbents (Hennion et al., 1998). A good agreement was obtained between experimental and modelled curves for oxamyl having $\text{Log} K_w$ values of 1.7 and 2.8 on C_{18} silica and PS-DVB sorbent and desethylatrazine with $\text{Log} K_w$ value of 3.5 on the PS-DVB sorbent. A great difference in sample volumes that can be percolated with a good recovery was observed when using a disk containing C_{18} silica or the PS-DVB copolymer. Taking a value of 90% for recovery, the corresponding sample volumes were, respectively, 30 and 300 ml for oxamyl on the two sorbents.

The *n*-octanol/water partition coefficient has become a recognised parameter to estimate compound hydrophobicity of general importance for modelling numerous environmental and biological properties associated with the use and disposal of organic compounds. Nakamura et al. (1996) proposed a general guide for the selection of sample processing conditions for the isolation and recovery of agricultural chemicals from water based on estimated $\text{Log} K_{\text{ow}}$ values. A group of empirical rules were presented in the form of a decision tree and said to be useful for selecting reversed-phase sorbents and organic solvents to aid experimental trial-and-error approaches to method development. However, for very polar analytes, $\text{Log} K_{\text{ow}}$ was shown to be of limited help for predicting the SPE recoveries, especially for very polar analytes with $\text{Log} K_{\text{ow}} < 1.5$ (Hennion et al., 1998). Since the retention mechanism is primarily governed by hydrophobic interactions between the analyte and the carbonaceous moieties of the alkyl chains grafted at the silica surface, a relation has been observed between the retention factors of the analytes and

their K_{ow} . Few data have been published with regards to polar analytes, especially for very polar analytes with $\text{Log } K_{ow}$ below 1.5. For these compounds, a large difference is observed between $\text{Log } K_{ow}$ and extrapolated $\text{Log } K_w$ values. For very polar analytes, a more rapid method is certainly to have in the laboratory a 10- or 5-cm long C_{18} column, and to extrapolate $\text{log } K_w$ from K measurement in a methanol–water mixture containing as high as possible water content. This is rapid and can be easily with autosampler and HPLC devices. Table 2.1 shows the extrapolated $\text{Log } K_{ow}$ values from the relationships $\text{Log } K$ –percentage of methanol of 14 different compounds ($\text{Log } K_{ow}$ –1.2–2.4 range). Five different cartridges and disks were compared including, C_{18} , PRP-1, PS-DVB (disk), high surface area/PS-DVB (HSA/PS-DVB) and PGC (carbon). It is clear that these HSA/PS-DVB (LiChrolut EN or Isolut NV+) are the sorbents to be selected for the extraction of (very) polar analytes when large sample volumes are required (Hennion et al., 1998; Pichon et al., 1998). The potential of PGC for extracting very polar compounds is shown in Table 2.1 for dealkylated and hydroxylated degradation products of atrazine down to cyanuric acid. Most of them have $\text{Log } K_{ow}$ values lower than 0. The limitation of other sorbents is clearly shown for the very polar ammeline, ammelide and cyanuric acid with $\text{Log } K_w$ values lower than 0.5 whereas they are higher than 2 with PGC. Using a 200-mg PGC cartridge, recoveries were above 90% with the handling of 250 ml of water sample for all the metabolites except the three more polar ones for which a 500-mg cartridge was required to obtain similar recoveries.

Table 2.1. Comparison of $\text{Log } K_w$ values obtained with C_{18} silicas, various PS-DVB copolymers and porous graphitic carbon (Hennion et al., 1998).

Compounds	$\text{Log } K_{ow}$	$\text{Log } K_w$				
		C_{18}	PRP-1	PS-DVB (disk)	HAS/PS-DVB	PGC
Cyanuric acid	-0.2	< 0.5	< 0.5	nd	< 0.5	2.6±0.1
Ammeline	-1.2	< 0.5	< 0.5	nd	< 0.5	2.4±0.2
Ammelide	-0.7	< 0.5	< 0.5	nd	< 0.5	2.5±0.2
Hydroxy-DIA	-0.1	1.0±		nd	1.8±0.1	3.0±0.2
Hydroxy-DEA	0.2	1.5±		nd	2.3±0.2	2.8±0.2
DEDIA	0	1.3±		nd	nd	2.8±0.1
Deisopropylatrazine (DIA)	1.2	2.3±		3.2±	4.4±0.3	> 3.5
Deethylatrazine (DEA)	1.4	2.7±		3.5±	4.8±0.3	3.2±0.2
Simazine	2.3	3.4±		4.1±	5.9±0.3	> 4
2-Chlorophenol	2.4	2.9±		3.6±		> 4
Oxamyl	-0.47	1.7±		2.8±	4.1±0.3	nd
Aldicarb	1.4	2.5±		4.0±	5.3±0.3	nd
Carbendazim	1.5			nd	5.7±0.3	> 4
Chloridazon	1.6	2.3±		3.8±		> 4

$\text{Log } K_w$ values extrapolated from the relationships $\text{Log } K$ –percentage of methanol.

Cyanuric acid: 2,4,6-trihydroxy-1,3,5-triazine; ammeline: 2,4-diamino-1,3,5-triazine; ammelide: 2-amino-4,6-dihydroxy-1,3,5-triazine.

2.6 Compound classification

The compounds of interest have been divided in groups based on their n -octanol/water partitioning and volatility from water. Compounds, which are preferentially in the n -

octanol phase (> 99.9 % in a 1/1 v/v *n*-octanol/water system) were designed as hydrophobic ($\text{Log } K_{ow} > 3$) and in the water phase (> 50%) as hydrophilic ($\text{log } K_{ow} < 0$). Compounds in the range $\text{Log } K_{ow}$ 1–3 were designated as medium polar and those in the range 0–1 as polar. For the volatility several parameters are available, i.e. vapour pressure (P_v) describing the volatility of the pure compound, the boiling point describing the potential to separate compounds in GC chromatography, and the Henry Law Constant (H in $\text{Pa}\cdot\text{m}^3/\text{Mol}$) describing the partitioning of a compound between air and water, and therefore taken as the best indicator for the evaluation of volatility losses from water.

For the classification we used the value proposed by Thomas (1990) for involatiles from water ($\text{Log } H < -3$), arbitrarily assigned categories of low ($\text{Log } H -3$ to -2) and intermediate volatility ($\text{Log } H -2$ to 1), and for the category of volatiles ($\text{Log } H > 1$) the properties of one of the least volatile compounds that can be analysed with purge and trap techniques, i.e. naphthalene ($\text{Log } H = 1.6$). This scheme is similar to the classification proposed by EPA (2000): non volatile ($\text{Log } H < -3$), slightly volatile ($\text{Log } H -3$ to -1), moderately volatile ($\text{Log } H -1$ to 1), volatile ($\text{Log } H 1$ to 3) and very volatile ($\text{Log } H > 3$). This classification has been used in the rest of the text and summarised in the following table.

Table 2.2 Operational classification of compounds for this study.

Compounds	Operational classification
Hydrophobic	$\text{Log } K_{ow}$ range: > 3
Medium polar	$\text{Log } K_{ow}$ range: 1-3
Polar	$\text{Log } K_{ow}$ range: 0-1
Hydrophilic	$\text{Log } K_{ow}$ range: < 0
Volatile	$\text{Log } H^*$ range: > 1
Intermediate	$\text{Log } H$ range: -2 to 1
Low-volatility	$\text{Log } H$ range: -3 to -2
Involatile from water	$\text{Log } H$ range: < -3

* H in $\text{Pa}\cdot\text{m}^3/\text{mol}$.

3. Physico-chemical properties of responsive compounds

3.1 Introduction

For a proper assessment of the suitability of concentration methods and sample-handling procedures it is necessary to evaluate a number of key properties of the target compounds. For this study, two main properties were selected: the *n*-octanol/water partitioning coefficient (K_{ow}), which is an important parameter for the prediction of the affinity to different sorbent materials, and the Henry Law coefficient (H), describing the partitioning between air and water and providing information on potential evaporative losses during sample handling.

The review of Reifferscheid et al. (1996), which contains a list over 450 compounds which have been tested with the UmuC-test, was taken as a starting point. The list of positive compounds (n=171) in the UmuC- test is included in Annex-2. The list comprises names and CAS-nrs of compounds which exhibited response without or after metabolic activation (S9; n= 32).

Other sources for databases of genotoxic properties of compounds, such as IARC (<http://www.iarc.fr>) and the EPA/IARC-GAP database (<http://www.epa.gov/gapdb>) containing comprehensive information on approximately 500 chemicals, were not used as it was difficult to isolate the UmuC-responsive compounds. For further exercises on compounds responsive in other test systems especially the IPA/IARC-GAP database seems very useful.

Based on information provided by RIZA, 37 compounds of the list in Annex-2, have been identified previously in XAD-extracts from River Rhine or Meuse water during monitoring studies conducted by RIZA (pers. comm. Dr. J. Staeb).

Measured values of K_{ow} and H were derived from the Environmental Fate database of Syracuse (SRC, 1998), which are also included in the EPIWIN software, available from US-EPA (Meylan and Howard, 1999a), and used for the prediction of K_{ow} and H for compounds for which no measured data were available. In the following sections a brief description is given of the estimation-software used. For more detailed descriptions of the methodological aspects and the comparison with other estimation algorithms we refer to Boethling and Mackay (eds., 2000), Meylan and Howard (1999a), and the documentation available from the EPIWIN website:

www.epa.gov/opptintr/exposure/docs/episuite.htm

3.2 Estimation of K_{ow}

The KOWWIN version 1.6 software (Meylan and Howard (1999c) was used to predict *n*-octanol/water partitioning coefficients (log K_{ow}) of organic compounds. The estimation procedure is based on the Atom/Fragment Contribution (AFC) method described by Meylan and Howard (1995) and differs from the fragment constant-based ClogP programme (Hansch and Leo, 1979; Leo 1993; Leo 1995). Both methods provide relatively good predictions of K_{ow} and are generally recommended (Leo, 2000).

The KOWWIN program requires input as a CASnr (or other compound database formats, such as MDL Mol files). For unknown compounds, not included in the database, structures are entered as SMILES notation (Simplified Molecular Input Line Entry System). The output of the calculated Log K_{ow} consists of the contributions of different structural fragments and correction factors. A comparison is made with a database of over 12000 experimentally determined Log K_{ow} values.

Currently, KOWWIN has been tested on a validation dataset of 10,338 compounds. Regression statistics for predicted versus experimentally Log K_{ow} values were: $r^2 = 0.94$; $sd = 0.47$; $me = 0.35$ (Meyan and Howard, 1999c).

Casnrns from the list of responsive compounds ($n=159$ excluding inorganic compounds) were entered in the KOWWIN programme. For approximately 89 % of the compounds matching CASnrns were found directly, and after correction this increased to 94%. For remaining compounds ($n=9$; 6%) Smiles notations needed to be entered. Smiles notations were either derived from Chemfinder (CambridgeSoft Corp., Massachusetts MA, USA; www.chemfinder.com) or created manually. For 100 (62%) of the compounds experimentally derived values were available and used in the further analysis.

3.3 Estimation of Henry's Law Constants

The volatility of compounds usually are described with: the boiling point, the vapour pressure (V_p in Pa) and the Henry Law Constant (H), which can be regarded as a measure for the air/water partitioning of a compound, and is defined as a function of the vapour pressure and the aqueous solubility S (Mol/m^3):

$$H = V_p / S \quad [\text{Pa}\cdot\text{m}^3/\text{Mol}]$$

The vapour pressure usually is regarded as a measure for the partitioning of the pure compound and the atmosphere; the boiling point provides an indication of the GC-amenity of a compound, and H can be regarded as a good indicator for the ability to extract a compound ($H > 40$) with purge-and-trap techniques (Van Loon, 1996)

The software program HENRYWIN version 3 estimates the Henry's Law Constant (H) based on the bond and the group contribution method (Hine and Mokerjee 1975; Meylan and Howard 1991 and 1999b; Mackay et al., 2000). The group method estimate is generally preferred (Howard, 1999b) when all fragment values are available (which was not the case for many of the compounds in this study). The bond contribution method is able to estimate many more types of structures and was designated as more accurate than the group method in an evaluation by Altschuh et al. (1999). In the review of Mackay et al. (2000) the bond contribution method was one of the two recommended methods.

The programme includes an experimental Henry's law constant database of 1650 compounds. The experimental database of 1650 compounds is comprised of values that were either (1) measured directly or (2) derived from reliable vapor pressure and water solubility data (584 were measured directly). Values measured directly are indicated in the output together with the reference source of the experimental data. Values derived from vapor pressure and water solubility are referred to a SRC and list both the vapor pressure and water solubility and the derived Henry's law constant. The reference source of the vapor pressure and water solubility are not reported; they are taken from SRC's

PHYSPROP Database. Reliable vapor pressure and water solubility data have been shown to derive very accurate Henry's law constants (Meylan and Howard, 1991).

Meylan and Howard (1999b) report that major estimation errors may occur when the bond estimation is less than 1×10^{-8} atm.m³/Mol (appr. 10^{-4} Pa.m³/Mol). Compounds with a $H < 3.10^{-3}$ Pa.m³/Mol are considered as non-volatile from water (Thomas, 1990).

Based on information provided in the manual (Meylan and Howard, 1999b) the statistics of the regression of predicted versus experimental values were: $r^2 = 0.85$; s.d. = 0.87; $me = 0.51$ (n=1293; compounds with $H < 10^{-4}$ Pa.m³/Mol not included).

The HENYWIN 3.0 programme requires a similar input format as the KOWWIN programme. The output consists of experimentally-derived values for H at 25°C (in atm.m³/Mol) together with results from the bond and (if available) the group contribution method. The options for different temperature corrections were not used.

Group contribution derived predictions were available for 73 compounds (46%). Based on the recommendations in the reviews of Mackay et al. (2000) and for reasons of consistency only bond-contribution estimations were used for compounds for which no measured values were available. For 39 (24%) of the compounds experimentally derived values were available and used in the further analysis.

3.4 Distribution of hydrophobicity

Estimated K_{ow} values and (if available) experimentally derived values (including source references) are listed in Annex-3. In Figure 3.1 and Table 3.1 frequency distributions are presented of the measured and predicted values for Log K_{ow} of the selected 159 UmuC-responsive compounds.

Table 3.1 Frequency distribution of Log K_{ow} of UmuC-responsive compounds.

Log K_{ow} class	All compounds			S-9 only		
	Frequency	%	Cumulative %	Frequency	%	Cumulative %
-4	1	.6%	.6%	1	3%	3%
-4 to -3	0	.0%	.6%	0	0%	3%
-2	5	3.1%	3.8%	0	0%	3%
-1	6	3.8%	7.5%	0	0%	3%
-1 to 0	23	14.5%	22.0%	2	6%	10%
0 to 1	25	15.7%	37.7%	3	10%	19%
1 to 2	35	22.0%	59.7%	10	32%	52%
2 to 3	23	14.5%	74.2%	4	13%	65%
4	21	13.2%	87.4%	5	16%	81%
5	11	6.9%	94.3%	1	3%	84%
6	4	2.5%	96.9%	3	10%	94%
7	3	1.9%	98.7%	1	3%	97%
7 to 8	2	1.3%	100.0%	1	3%	100%
>8	0		100.0%	0		100%
Total	159			39		

S9-only: selection of compounds responsive after metabolic activation (S9).




Figure 3.1 Frequency distribution of predicted ($n=59$) and experimentally derived ($n=100$) Log-K_{ow} values for UmuC-positive compounds.

Table 3.2 Frequency distribution over classes of Log K_{ow} .

LogK _{ow} class	All comp.		S-9 only		RIZA only	
	Frequency	%	Frequency	%	Frequency	%
<0	35	22%	3	10%	3	8%
0-1	25	16%	3	10%	6	17%
1 to 3	58	36%	14	45%	17	47%
3 to 6	36	23%	9	29%	8	22%
>6	5	3%	2	6%	2	6%
Total	159		31		36	

It is apparent that the UmuC-responsive compounds differ largely in hydrophobicity (Log K_{ow} $-4 <$ to >8). It should be noted that for both the hydrophilic compounds ($\text{Log K}_{ow} < 0$) and the super-hydrophobics ($\text{Log K}_{ow} > 8$) the predictions tend to have a high uncertainty.

Observations for the different hydrophobicity classes are given in Table 3.2. Several conclusions can be drawn from the frequency distributions:

- UmuC-responsive compounds show large differences in their *n*-octanol/water partitioning coefficient. Approximately 52% of the responsive compounds falls within the class-ranges of polar (16%) to semi-polar (36%) substances. Hydrophobic compounds constitute 26% (of which 3% with $\text{Log K}_{ow} > 6$) and hydrophilic compounds approximately 22%.

- The group of compounds that is responsive after S9-metabolic activation (n=31) seems to have a slightly higher distribution of K_{ow} , with 10% within the class-ranges of hydrophilics, 10% polars, 45% semipolar, and 35% hydrophobics.
- Responsive compounds observed in surface water studies by RIZA exhibit a similar distribution pattern, with slightly lower frequencies in the category of the hydrophilic compounds.
- Compared to the information from Chapter-2 it is important to note that approximately 50% of the UmuC-responsive compounds lies within the critical interval for break-through losses of many commonly used SPE phases.

3.5 Distribution of Henry's Law Constants

Estimated and experimentally derived values for H are given in Annex-4. Values are given in $\text{Pa}\cdot\text{m}^3/\text{Mol}$. Experimental values were available for 39 (24%) of the compounds, of which 27 (17%) were calculated from reliable vapour pressure and solubility data in the SRC PhysProp database. Less reliable predictions with H values below the cut-off value of 10^{-4} ($\text{Pa}\cdot\text{m}^3/\text{Mol}$) have been listed in Annex-4 and will not be treated further in the analysis.




Figure 3.2 Frequency distribution of predicted (n=120) and experimentally derived (n=39) Log H ($\text{Pa}\cdot\text{m}^3/\text{Mol}$) values for UmuC-positive compounds. Compounds below $\text{Log H} = -3$ are involatile from water. Cut-off value for volatiles: $\text{Log H} = 1.6$ (naphthalene)

Frequency distributions of the most important classes are given in Figure 3.2 and Table 3.3.

Table 3.3 Frequency distribution of Log H of UmuC responsive compounds.

LogH class	Frequency	%	Cumulative %
<-4	80	50%	50%
-4 to -3	10	6.3%	56.6%
-3 to -2	18	11%	67.9%
-2 to 0	29	18%	86%
0-1	12	8%	94%
1-2	6	4%	97%
>2	4	3%	100%
Total	159		

Approximately 50% of the responsive compounds have extremely low Log-H values, below the cut-off value of -4 . Approximately 56% of the values is below the level for virtually involatiles from water (Log H <-3 ; Thomas, 1990). The arbitrarily assigned category of low-volatility (Log H >-3 and <-2) comprises 11%. The category of intermediate volatility from water (Log H <-2 to 1) constitutes 24%. Compared to naphthalene (Log H=1.6), which is considered as a relatively volatile PAH, approximately 7% has a similar or higher volatility from water (Log H > 1).




Figure 3.3 Plots of Log H versus Log K_{ow} of UmuC responsive compounds. Only observations with Log H >-4 are included..

Compounds sensitive to volatilisation during sample-handling probably will have a Log H > 0 . From the plot in Figure 3.3 it can be derived that the hydrophobicity of these

compounds (n=22) ranges from $\text{LogK}_{\text{ow}} = -3$ to $\text{LogK}_{\text{ow}} > 3$. In table 3.4 these compounds have been listed individually.

Table 3.4 Selection of UmuC-positive compounds prone to evaporative losses during sample-prep.

Compound	CASnr	LogH	Reference exp. value
Benz(a)anthracene(+ S9)	56-55-3	0.07	Bamford, HA et al. (1999)
Azobenzene(+ S9)	103-33-3	0.12	SRC
Epoxy styrene	96-09-3	0.19	SRC
Styrene oxide	96-09-3	0.19	SRC
t-Butylhydroperoxide	75-91-2	0.20	
Nystatin	1400-61-9	0.21	
1,3-Dioxane	505-22-6	0.46	
Epichlorhydrin	106-89-8	0.47	SRC
4-Methyl-1,3-dioxane	1120-97-4	0.59	
2,4-Dimethyl-1,3-dioxane	766-20-1	0.71	
Propylene oxide	75-56-9	0.83	SRC
8-Proplotactone	57-57-8	0.88	
5-Butyrolactone	36536-46-6	1.00	
1,2-Epoxybutan	106-88-7	1.25	SRC
Benzylchloride	100-44-7	1.61	SRC
Ethylene dibromide	106-93-4	1.82	SRC
Methylene bromide	74-95-3	1.91	Moore, R.M. et al. (1995)
Trichloronitromethane	76-06-2	2.30	Kawamoto, K & Urano, K. (1989)
Bromobenzene	108-86-1	2.38	Shiu, W.Y. & Mackay, D. (1997)
Methylbromide	74-83-9	2.79	SRC
1-Bromo pentane	110-53-2	3.29	SRC

4. Selection of suitable SPE methods

4.1 Introduction

The key problem when applying SPE remains the method development and the primary decision for the analysts is the choice of the sorbent that is able to solve their trace-analysis problem. Method development being related to the properties of the analytes of interest, many aspects due to the various physico-chemical properties of compounds included in a multiresidue analysis have to be considered. A first approach for the method development is the process, which occurs during the extraction and this approach has remained rather empirical based on experimental trial-and-error procedures. Alternative approaches based on computer-aided strategies and simulation require an appropriate level of theory so that at decision steps in the method development process fast simulation or calculation procedures can be used in place of trial-and-error experiments.

For the extraction of compounds from a wide range of polarities, the analyst usually focuses on the low retention of polar analytes, which can be lost during the extraction step. The potential of the highly cross-linked polymers and the carbon-based sorbents for the extraction of (very) polar compounds was previously demonstrated. However, for this type of analysis, difficulties other than those related to the low retention of polar compounds can be expected (see below).

A typical SPE sequence consists of four general steps: (i) conditioning the sorbent, (ii) percolating the sample, (iii) drying the cartridge and (iv) desorption. These steps are described and presented in the following scheme.

In chapter 3, a study of physicochemical properties of the 159 target analytes covering broad ranges of polarity and volatility was presented. Polarity is presented in the form of $\text{Log } K_{ow}$ ($\text{Log } K_{ow}$ range from -4 to 6) and volatility in the form of $\text{Log } H$. In general, it is always a challenge to extract as much as possible in one run in order to decrease the cost and the time of the analysis in the environmental field (Pichon et al., 1998 and Nakamura et al., 1996). It depends on the polarity of the most polar analytes and one has to check if these ones are extracted with good recoveries with the required sample volume to be handled. For the hydrophobic ones, in theory, there is no problem of breakthrough.

However, if the extraction can be predicted, very often practical problems coming from the physico-chemistry are encountered. One problem may occur during the sample percolation, because recoveries of hydrophobic analytes with very low water solubility are low unless a certain percentage of organic solvent or surfactant is added in the sample (Pichon et al., 1998 and Nakamura et al., 1996). But if the addition of an organic solvent solves the problem of the hydrophobic ones, it decreases the breakthrough volumes of the more polar ones, which can then be poorly recovered using this procedure. Another problem is in the reconstitution of the extract. When very polar and non-polar analytes are together, complete solubilization of the extracts is often impossible: addition of water (or methanol) for the more polar ones, whereas very hydrophobic analytes can only be dissolved in a non-polar organic solvent.

*Typical steps in SPE sequence***Conditioning solvent** (typically 3-5 hold-up volumes)

- (a) Ensures reproducible retention and flow. Critical step for particle-loaded membranes
- (b) Helps to minimise contamination of extracts by solvent impurities
- (c) Replace by sample solvent before processing sample

Flow rates (typical range 0.2-1.5 mm/s)

- (a) More critical for cartridges than disks due to their variable and heterogeneous packing density (channelling)
- (b) More critical when the sample volume exceeds the breakthrough volume as typical sampling devices provide too few theoretical plates for flow independent retention

Sample properties

- (a) Remove excessive particle matter by filtration or centrifugation to maintain a constant sample-processing rate
- (b) Add small volume of organic solvent (1-10%, v/v) to large volume water samples to ensure sorbent remains solvated and to maintain a constant (fast) sample-processing rate. Critical for polar analytes (see text)
- (c) Adjust pH to reduce ionisation of weak acids and bases for reversed-phase sampling

Drying time (typically 1-5 min, but sometimes considerably longer)

- (a) Sufficient to remove all sample solvent trapped in the sorbent pores
- (b) Excessive drying may result in low recovery of analytes from evaporation or retention in poorly solvated regions of the sorbent

Eluting solvent (ideally 2-3 hold-up volumes but often larger)

- (a) Should be a strong solvent able to displace all analytes from the sorbent with a small volume
- (b) Normally should be volatile and miscible with sample solvent

4.2 Application of selection criteria

From chapter 2, five possible candidate sorbents have been selected, which have the potential of covering the polarity range mentioned above i.e., -4-6.

Four critical selection criteria were used: breakthrough volume, desorption mixture/volume, recovery and enrichment factor. Other parameters, such as costs, simplicity and ease of implementation were not discriminatory between different SPE methods.

For the sake of clarity, we divided the polarity range into four classes of compounds: A for hydrophilic ($\text{Log } K_{ow} < 0$), B for polar ($0 < \text{Log } K_{ow} < 1$), C for medium polar ($1 < \text{Log } K_{ow} < 3$) and D for hydrophobic ($\text{Log } K_{ow} > 3$).

The scores from ‘-’ to ‘++’ of each SPE method for each parameter along this compound classification are presented (see Table 4.1). Each parameter of the four critical ones shown as legend is divided into four categories from bad ‘-’ to excellent ‘++’. For instance, breakthrough volume higher than the minimum selected volume (500 ml) is considered as excellent and lower than 5 ml is bad. Some scores shown in the Table like ‘+/-’, means that the results vary between the two scores, and that is because either the $\text{Log } K_{ow}$ is relatively large and/or the behaviour may depend on other experimental conditions and physico-chemical properties of compounds.

Three parameters are important to achieve a desired enrichment factor, i.e., breakthrough volume, desorption characteristics (mixture/volume dimensions), and recovery. Other aspects such as simplicity, costs and implementation are not discriminative regarding the final selection of an SPE method (see Table 4.2). The basic steps of the SPE procedure

are simple and straightforward and can be trained in a one-day course. All equipment can be purchased from various sources and costs differ mainly in the degree of automation, which can vary from single-sample off-line procedures up to fully automated batch processing. Drying of cartridges is only discriminative when methylene chloride is used as desorption solvent. The problem is, especially for carbon cartridges, that water remains in the pores of the stationary phase and that methylene chloride is not able to reach the analytes, which leads to inefficient desorption.

Table 4.1. Criteria table for the five selected sorbents.

SPE method	Breakthrough volume ¹	Desorption ²	Recovery ³	Enrichment Factor ⁴	Drying (cartridge)
Lichrolut EN	A: - / ±	A: ++	A: - / ±	A: - / ±	+ / ±
	B: +	B: ++	B: +	B: +	
	C: ++	C: ++	C: ++	C: ++	
	D: ++	D: +/-	D: +/-	D: +/-	
Isolut ENV+	A: - / ±	A: ++	A: - / ±	A: - / ±	+ / ±
	B: +	B: ++	B: +	B: +	
	C: ++	C: ++	C: ++	C: ++	
	D: ++	D: +/-	D: +/-	D: +/-	
Oasis HLB	A: -	A: -	A: -	A: -	+
	B: + / ±	B: +	B: +	B: ±	
	C: +	C: +	C: +	C: ±	
	D: +	D: +	D: +	D: ±	
C18	A: -	A: ++	A: -	A: -	+
	B: -	B: +	B: -	B: -	
	C: ±	C: +	C: ±	C: ±	
	D: ++	D: ++	D: ++	D: ++	
GCB	A: +	A: ++	A: +	A: +	-
	B: ++	B: +	B: ++	B: ++	
	C: ++	C: ±	C: ++	C: ++	
	D: +	D: -	D: - / ±	D: +	

Explanation

A. $\text{Log } K_{ow} < 0$	++	= excellent
B. $0 < \text{Log } K_{ow} < 1$	+	= good
C. $1 < \text{Log } K_{ow} < 3$	±	= intermediate
D. $\text{Log } K_{ow} > 3$	-	= bad

Scores

¹ Breakthrough (ml)	² Desorption (ml)	³ Recovery (%)	⁴ Enrichment factor
++ = 500	++ = < 2	++ = > 85	++ = > 1000
+ = 100–500	+ = 2–5	+ = 50–85	+ = 100–1000
± = 5–100	± = 5–10	± = 5–50	± = 10–100
- = < 5	- = > 10	- = 0–5	- = < 10

The groups presented here are derived from Annex A1 (factsheets of selected SPE methods) directly or indirectly as follows: target compounds present in Annex A1 are directly included. If one or a group of target compounds are missing, then comparison and estimation are made according to their physico-chemical properties.

The following example should clarify the table. LiChrolut EN can be used in the $\text{Log } K_{ow}$ 0–3 range efficiently (B: +; C: ++). Outside this range, two problems appear. The first one is the desorption problem, especially for very hydrophobic analytes (D: +/-). The second one is recovering the more polar ones. To solve the first problem, a com-

promise is needed. Addition of a miscible organic solvent to the sample (see below) before sampling can solve this problem but there is risk to loss (very) polar analytes. That is why a '+ / -' scores is given for the desorption of hydrophobic compounds from LiChrolut EN sorbent. The same selection strategy was used for the rest of the Table.

In Table 4.2 the scores for the non-selective criteria are indicated. Ease of handling, handling time required for the procedure and potential for implementation in routine laboratory settings do not differ much between the different methods. Disk-options usually are more rapid compared to cartridges. As personal costs are dominant, the slight differences in price between disks and cartridges are of minor influence.

Table 4.2 Non-selective criteria.

SPE method	Simplicity	Costs	Implementation
Lichrolut EN	+	+	+
Isolut ENV+	+	+	+
Oasis HLB	+	+	+
C18	+	+	+
GCB	+	+	+

4.3 Constraints from the UmuC-test protocol

In order to be able carry out the UmuC-test successfully, the concentration technique has to meet some requirements. First, enrichment factor must be high enough (1000 or higher) to achieve the sensitivity needed. A sample volume of 500 ml or higher is needed to be handled and to get a final extract volume of 0.3–0.5 ml. Second, the final extract solvent must be suitable for the UmuC-test. Methanol, acetone and DMSO (Rotteveel, personal communication) are tolerable solvents at 3% volume for the test (water is the ideal solvent if possible). Finally, to improve the enrichment factor and to avoid loss of volatiles during evaporation, a high boiling point solvent (DMSO or water) can be added to the extract.

4.4 Various SPE approaches

As can be seen in Table 4.1, each sorbent shows its advantages and disadvantages along the four Log K_{ow} ranges. It is very difficult to find a single method to extract the polar and non-polar analytes together in one procedure. Highly cross-linked polymers are capable to extract polar and moderately polar analytes efficiently. For the very polar (Log $K_{ow} < 0$), GCB is recommended. However, other kind of problems appears; compounds with low water solubility cannot be recovered. This can be due to adsorption of these compounds on tubings and vessels or to an incomplete desorption of strongly retained compounds on the sorbent. The difficulty to simultaneously extract compounds from a broad range of polarities was well described by Norberg et al. (1995) in the case of the study of organophosphorus compounds. When plotting the recoveries vs. the sample volume, three types of behaviour were observed. For the most polar compounds, a quantitative recovery was found for low sample volumes (early breakthrough). For the moderately polar compounds, recoveries were quantitative over the whole 10-100 ml range. For the most hydrophobic analytes, recoveries slowly increased with the sample volumes

because a fraction of the analytes was adsorbed in the preconcentration system. To eliminate adsorption phenomena, a low percentage of organic solvent (1-10%, v/v) or a low amount of surfactant can be added directly to the sample (Pichon et al., 1998 and Nakamura et al., 1996). However, the presence of these modifiers decreases the retention of polar compounds on the sorbent during the percolation and causes loss of these compounds. Therefore, a compromise has to be found between adsorption of the low water-soluble compounds and the retention of the highly polar ones.

Concerning the problem of incomplete desorption, especially for GCB cartridge, it can be difficult to find the solvent that will allow a complete desorption in a small volume in order to obtain a high enrichment factor. In case of large fraction, it can be reduced by evaporation of the residue but this step needs to be well controlled because of the risk of loss of high vapour pressure compounds, i.e., phenols (Castillo et al., 1997) or organophosphorus compounds (Molina et al., 1996). In most cases, in order to limit the losses of volatile compounds, the evaporation step is carried out under mild conditions, i.e., under a stream of nitrogen instead of the use of a rotary evaporator until a final volume of 0.5 ml. Besides, a degradation of compounds in contact with some solvents can occur during this step as it was mentioned for aryloxyphenoxypropionic acids (Lagana et al., 1998).

Many studies reported the difficulty to select proper solvents (Castillo et al., 1997 and Di Corcia et al., 1993). As an example, methanol was preferred above acetonitril for the desorption of catechol from PS-DVB sorbent (LiChrolut EN) owing to a better solubilisation of this compound with methanol. The authors recommended a mixture of methylene chloride-methanol (4:1) according to their eluting strength on reversed-phase sorbent. It was also pointed out that residual water in the GCB cartridge had to be reduced to a minimum and, that, when this was not done, low and irreproducible recoveries were obtained because the water can hinder intimate contact between the desorption mixture and the GCB.

In the case, that the range of polarities is too broad, two separate procedures are recommended, one optimised for the (very) polar and moderately polar ones, and in a subsequent suited to the non-polar ones (Pichon et al., 1998 and Nakamura et al., 1996). When two subsequent procedures are required, it can also be interesting to use two different sorbents to facilitate the desorption. As an example, two off-line procedures were defined: (i) the polar and moderately polar pesticides were extracted on a PS-DVB disk and subsequently eluted with acetonitrile (methanol is also possible), (ii) after the addition of 10% of methanol to the sample, the non-polar pesticides were extracted on a C₁₈ disk and eluted with a mixture of methylene chloride-methanol (4:1) (Pichon et al., 1998). However, in this study, a Log K_{ow} 1-6 range was considered and very polar analytes were not included. Very polar analytes may still be difficult to extract with these sorbents, and graphitized carbon has been shown to solve this problem for some analytes (Pichon et al., 1995; Hennion et al., 1995; and Guenu et al., 1996). For hydrophobic compounds, a C₁₈ silica is usually preferred over PS-DVB sorbent since recoveries are good for a sample volume of 500 ml or higher even after the addition of methanol before the extraction step.

Castillo et al. (2001a; 2001b; 1999a-c); Thomas et al. (2001), La Farre et al., 2001a-d) and Remstma et al. (1999) used a sequential SPE (SSPE) method to extract polar and


non-polar analytes by two separate procedures (C_{18} silica and Isolut ENV+ or LiChrolut EN). This method was used in order to extract non-polar compounds on C_{18} silica and moderately polar and polar ones on PS-DVB polymer (Isolut ENV+) (for more details, see Figure 4.1).

Putschew et al. (2001) used the same SSPE approach but with LiChrolut EN (200 mg) and Envi-Carb (250 mg) cartridges in order to extract the very polar and partly ionic analytes. Recoveries of the eight selected very polar iodinated benzene derivatives were in the range of 80–100%. The enrichment factor was higher than 1000. It was not possible to achieve good recoveries for all compounds by just one extraction step. The only drawback of this method is that the Envi-Carb cartridge had to be eluted against the extraction direction (back flush) with a mixture (8 ml) of acetonitrile-ultra pure water (1:1, v/v) and trace of ammonium acetate.

It is clear from Table 4.1 and the above discussed parameters that the conditions for SPE must be carefully selected to achieve the goal described in the Introduction. Considering the physico-chemical of the components ($\text{Log } K_{ow}$) and the characteristics of the sorbent, the guidelines for the selection of SPE conditions shown in Figure 4.1 should be useful. Although, from this figure it might appear that the combined procedures are complicated, the contrary is true. These procedures consist of two simple SPE procedures, although carried out in a sequence.

When very polar ($\text{Log } K_{ow} < 0$, i.e., 35 from the 170-component list) are not included in the compound list, the SSPE procedure using C_{18} silica and Isolut ENV+ in series is recommended. That is because no desorption problems occur and the addition of methanol (up to 10%) is possible by C_{18} silica. Moreover, good recoveries (higher than 80%) and high enrichment factors (higher than 1000) can be obtained for all compounds over $\text{Log } K_{ow}$ 0–6 range. A compromise is needed when very polar ($\text{Log } K_{ow} < 0$) analytes are added to the list. Using C_{18} silica and GCB cartridges in series can recover the $\text{Log } K_{ow}$ – 2 to 6 range but some desorption problems by GCB can be expected and sometimes a backflush desorption is necessary. These desorption problems were already obviously encountered with the GCB cartridges: pure methanol, acetonitrile or dichloromethane were unable to desorb many organic pollutants. With GCB cartridges, the desorption of many analytes occurred with 6 ml of dichloromethane: methanol (80:20, v/v) with 300 mg of sorbent and the backflush desorption is even recommended for the 1-g cartridges. Cartridges allowing percolation and desorption in the opposite way are now available.

Oasis HLB sorbent is selected because it is capable of extracting acidic, basic and neutral compounds whether polar or non-polar. Retention of some analytes measured on a column packed with this polymer was compared with that obtained on C_{18} silicas, showing a large increase of K_w values. For instance the $\text{Log } K_w$ value for catechol is around 2.5 whereas it is 1.1 on C_{18} silica. It was measured to 1.6 using PRP-1. The increase in retention is explained by both the specific area of $800 \text{ m}^2/\text{g}$ and the occurrence of the pyrrolidone group in the polymer, which is a hydrogen acceptor (dipole interactions). It is a new sorbent so that more work should be done to look at the limitation in the extraction of (very) polar analytes as it was done for other highly cross-linked PS-DVB polymer. To date, Oasis HLB is more used to extract drugs from plasma and it is not recommended when high enrichment (higher than 200) factor is required.



Abbreviations: C18: silica cartridge; GCB: carbon cartridge; H2O: water; C6: *n*-hexane; DCM: dichloromethane; MeOH: methanol; N2: gas nitrogen.

Figure 4.1. Flow diagram of the selected SPE procedures. Parts presented in bold are the most critical steps in each SPE procedure (see text for more details).

4.5 Proposal

The following conclusions are made as a result of a number of carefully selected studies on this field by the authors and members of the board group. The major difficulty in the environmental sciences is the separation/detection of compounds that are present in the environment at low concentrations (typically $\mu\text{g/l}$ level). An efficient concentration technique is required. SPE is the method of choice over other discussed methods. Next to its simplicity, efficiency and flexibility, SPE meets the UmuC-test requirements.

For the extraction of compounds varying over a wide range of polarities ($\text{Log } K_{ow} -2 - 6$ range), a single SPE procedure should ideally be applied and all the compounds should be eluted in a single fraction prior the UmuC-test. However, difficulties other than those related to the low retention of polar compounds oppose this ideal. In general, extraction of polar and hydrophobic compounds simultaneously is rather difficult due to desorption problems. In case of quantitative extraction of polar compounds, desorption of the (most) hydrophobic is nearly impossible.

To overcome this dilemma, a sequential SPE procedure is recommended of which a part is optimised for the (very) polar and moderately polar compounds, and the other one suits to the non-polar ones. A combination of two different sorbents is preferred to facilitate the desorption. C₁₈ silica and PS-DVB (Isolut ENV+) can be used efficiently in series to recover Log K_{ow} 0 – 6 range. The hydrophobic compounds are extracted on the modified C₁₈ silica after addition of 10% of methanol to the sample and eluted with a mixture of methanol–methylene chloride (9:1). Next, the polar and moderately polar compounds present in the residue from the C₁₈ cartridge, are extracted on the PS-DVB sorbent after pH adjustment of the residual sample, and subsequently eluted with methanol.

When very polar analytes (Log K_{ow} < 0) are added in the list, GCB can solve the problem and can be combined with modified C₁₈ silica in a similar sequential SPE procedure. The only difference is that analytes are eluted with a mixture of methylene chloride–methanol (4:1) from the GCB material. The only drawback of this method is that sometimes desorption problems by GCB occur so that a backflush desorption mode is necessary.

In the final step, the desorption solvent is evaporated under a gentle stream of nitrogen to a volume of 100–500 µl depending on the required enrichment factor. In this step, all traces of methylene chloride are removed.

The detailed descriptions of the SPE approaches in Figure 4.1 are presented in Annex-5.

5. Discussion and concluding statements

From the currently available method for enrichment of contaminants SPE methods have been selected as the method of choice. They meet the main requirements for the final UmuC-testing, i.e. an enrichment by a factor of thousand and compatibility with the solvent requirements for this test. Although simple procedures such as lyophilization, reverse osmosis and freeze drying and with a fully compatible extract, viz water, their limitations exclude them for further investigations.

With the current state-of-the-art in SPE techniques it was not possible to select one single generic method, suitable for routine (waste)water screening by the UmuC-test, which is fully comprehensive. For the group of organics with extremely high or low hydrophobicity, special dedicated single-step methods have to be applied. In both cases, recovery at one end automatically leads to losses at the other end. If the polarity range is more limited, several newly developed SPE materials are available, especially for these compounds having $\log K_{ow}$ values in the 0 – 6 range. If the range of compounds is limited to the more polar end of the scale, GCB materials can be selected. Generic experimental conditions are described in a special fact sheet in Annex 5.

From the combination of figures 2.1 and 3.1 it becomes clear that most commonly used SPE materials exhibit breakthrough losses in $\log K_{ow}$ interval of 0-3, which coincides with $\log K_{ow}$ of more than 50% of the UmuC-responsive compounds.

The proposed SSPE method covers a wide range ($\log K_{ow}$:-2 to 6), but probably 10% of the compounds in the hydrophilic range (<-2) and approximately 5% in the superhydrophobic range (>6) will be missed. Further losses due to evaporation during sample-handling may be present for approximately 14% ($n=22$) of the compounds, if a conservative cut-off of $H=0$ is chosen. From the plots in Figure 3.3 it can be derived that these compounds do not coincide with the superhydrophobics or super hydrophilic compounds, and that losses up to 30% of the total range of compounds may occur. If the volatile fraction is of importance, than the additional application of purge-and-trap concentration techniques prior to genotoxicity testing, may reduce the number of 'missed' compound with more than 10%.

The category of inorganic compounds has been excluded in this study. Although in theory adding specific ion-exchange resins in sequence with the other SSPE may result in the trapping of metallic compounds of interest (e.g. selenium, chromium compounds), it is not possible to isolate the redox-state and chemical species, which were tested as genotoxic. Usually these reactive species are only a fraction of the total element concentration in (waste)water. Inorganic compounds which are considered as genotoxic, usually were tested as pure solid compounds, representative for exposure in occupational conditions (e.g. paint and semi-conductor factories), and to a lesser extent as dissolved compounds. Based on a rough comparison of reported LEDs (lowest effective dose) and HIDs (highest ineffective dose) for Se and Cr compounds in the IARC/EPA-GAP database, exposure levels (mg/L to g/L) by far exceed levels commonly encountered in effluents of wastewater treatment installations. Against this background, we recommend not to include inorganics in the generic screening approach. In cases where genotoxicity of

inorganics may be of importance (e.g. effluents of specific industrial sectors), dedicated studies are considered as more cost-effective.

In the scheme in Figure 4.1, two sequential-SPE procedures are presented, that can be used to extend the polarity range. Addition of some 10% of MeOH to the sample excludes losses of hydrophobic analytes. Compounds not retained in the first SPE cartridge (C18) are subjected to a second stationary (more tentative) phase (GCB or PS-DVB). After drying an elution, both extracts are recombined and concentrated further by a gently nitrogen purge. If the toxicity has to be evaluated by means of identification of analytes, scheme 4.1 also proposes procedures to yield separate fractions (with 4 intervals) for these TIE-studies

In previous work, commissioned by RIKZ, detailed procedures have been documented (Adahchour et al., 2000a and 200b). Several of the indicated sorbent materials can also be incorporated in a HPLC setting, for on-line separation in combination with automated fraction-collecting systems, dosing to e.g. microwell-plates suitable for further testing with the UmuC-protocol, similar as the Oasis-based system (Tox-Print) described by Bobbeldijk et al. (2001).

In conclusion, the present study revealed that genotoxic compounds exhibit orders of magnitude variations in volatility and hydrophobicity. From the currently available concentration techniques SPE/SSPE techniques were assigned as the method of choice for the screening (waste)water and combination with the UmuC-test. A wide range of SPE materials, documented in factsheets is commercially available. The most critical parameters for the proper selection of sorbent materials and dimensions are: volatility, breakthrough volume and desorption characteristics. A combine SSPE method, including C-18 followed by GCB was proposed as the most comprehensive method for contaminants with Log Kow ranging from -2 to 6.

References

Chapter 1

- Bobbeldijk, I., A. Brandt, B. Wullings and Th. Noij (2001). High-performance liquid chromatography-Toxprint: chromatographic analysis with a novel (geno)toxicity detection. *J. Chrom. A* 918: 277-291.
- Castillo, M., M.C. Alonso, J. Riu, M. Reinke, G. Kloter, H. Dizer, B. Fischer, P.D. Hansen, and D. Barcelo (2001). Identification of cytotoxic compounds in European wastewaters during a field experiment. *Analytica Chimica Acta* 426:265-277.
- De Maagd, P.G.J. (2000a). Bioaccumulation tests applied in whole effluent assessment: A review. *Environmental Toxicology and Chemistry* 19:25-35.
- De Maagd, P.G.J. and M. Tonkes (2000b). Selection of genotoxicity tests for risk assessment of effluents. *Environmental Toxicology* 15:81-90.
- De Maagd, P.G.J., M. Tonkes, J.L. Maas en C. van de Guchte (1999). *Genotoxiciteit als effect-parameter binnen totaal-effluentbeoordeling*. RIZA werkdocument nr. 99.110x. RIZA, Lelystad. (In Dutch).
- Filipic, M. and M.J. Toman (1996). *Water Sci. Technol.* 34: 1-7. Cited in De Maagd and Tonkes (2000).
- Gauthier, L., M.A. Van der Gaag, J. L'Haridon, V. Ferrier, and M. Fernandez (1993). In vivo detection of waste water and industrial effluent genotoxicity: use of the Newt Micronucleus Test (Jaylet Test). *The Science of the Total Environment* 138:249-269.
- Hamer, B., N. Bihari, G. Reifferscheid, R.K. Zahn, W.E. Muller, and R. Batel (2000). Evaluation of the SOS/umu-test post-treatment assay for the detection of genotoxic activities of pure compounds and complex environmental mixtures. *Mutation Research* 466:161-171.
- Klamer, H.J.C., L.A. Villerius, J. Roelsma, P.G.J. de Maagd, and A. Opperhuizen (1997). Genotoxicity testing using the Mutatox(TM) assay: Evaluation of benzo(a)pyrene as a positive control. *Environmental Toxicology and Chemistry* 16:857-861.
- Reifferscheid, G. and J. Heil (1996). Validation of the SOS/umu test using test results of 486 chemicals and comparison with the Ames test and carcinogenicity data. *Mutation Research* 369:129-145.
- Stahl, R.G.J. (1991). The genetic toxicology of organic compounds in natural waters and wastewaters. *Ecotoxicology and Environmental Safety* 22:94-125.
- Tonkes M., H.B. Pols, H. Warner en V. Bakker (1998). *Totaal-Effluentbeoordeling*. RIZA report 98.34. RIZA, Lelystad. (In Dutch).
- Vargas, V.M., R.R. Guidobono, C. Jordao, and J.A. Henriques (1995). Use of two short-term tests to evaluate the genotoxicity of river water treated with different concentration/extraction procedures. *Mutation Research* 343:31-52.

Chapter 2 REVIEWS

- Bruzzoniti, M.C., C. Sarzanini and E. Mentasti (2000). Preconcentration of contaminants in water analysis. *Journal of Chromatography A* 902:289-309.
- Hennion, M.C. (1999). Solid-phase extraction: method development, sorbents, and coupling with liquid chromatography. *Journal of chromatography A* 856:3-54.
- Hennion, M.C. (2000). Graphitized carbons for solid-phase extraction. *Journal of chromatography A* 885:73-95.

- Huck, C.W. and G.K. Bonn (2000). Recent developments in polymer-based sorbents for solid-phase extraction. *Journal of chromatography A* 885:1-72
- León-González, M.E. and L.V. Pérez-Arribas (2000). Chemically modified polymeric sorbents for sample preconcentration, *Journal of chromatography A* 902:3-16.
- Liska, I. (2000) Fifty years of solid-phase extraction in water analysis - historical development and overview. *Journal of chromatography A* 885:3-16.
- Martha J.M. Wells and Lan Zhou Yu (2000). Solid-phase extraction of acidic herbicides. *Journal of chromatography A* 885:237-250.
- Noble, A. (1993). Partition coefficients (n-octanol-water) for pesticides. *Journal of chromatography A* 642:3-14.
- Pichon, V. (2000). Solid-phase extraction for multiresidue analysis of organic contaminants in water. *Journal of chromatography A* 885:195-215.
- Rodríguez, I., M.P. Llompарт and R. Cela (2000). Solid-phase extraction of phenols. *Journal of chromatography A* 885:291-304.
- Sabik, H., R. Jeannot and B. Rondeau (2000). Multiresidue methods using solid-phase extraction techniques for monitoring priority pesticides, including triazines and degradation products, in ground and surface waters. *Journal of chromatography A* 85:217-236.

Chapter 2 ARTICLES

- Aguilar, C., F. Borrull and R.M. Marcé (1996). On-line and off-line solid-phase extraction with styrene-divinylbenzene-membrane extraction disks for determining pesticides in water by reversed-phase liquid chromatography-diodearray detection, *Journal of Chromatography. A.* 754:77-84.
- Baltussen, E., F. David, P. Sandra, H. G. Janssen, C. A. Cramer (1998). Retention model for sorptive extraction-thermal desorption of aqueous samples: application to the automated analysis of pesticides and polyaromatic hydrocarbons in water samples, *Journal of Chromatography A.* 805:237:247.
- Baun, A., N. Bussarawit and N. Nyholm (August 1998). Screening of pesticide toxicity in surface water from an agricultural area at Phuket Island (Thailand). *Environmental Pollution* 102:185-190.
- Bennett, D.A., A.C. Chung, S.M. Lee (October 1997). Multiresidue method for analysis of pesticides in liquid whole milk. *Journal of AOAC International* 80:1065-1077.
- Bobeldijk, I., J.P.C. Vissers, G. Kearney, H. Majorand, J.A. van Leerdam. (2001a). Screening and identification of unknown contaminants in water with liquid chromatography and quadrupole-orthogonal acceleration-time-of-flight tandem mass spectrometry. *Journal of Chromatography. A.* 929: 63-74
- Bobeldijk, I., A. Brandt, B. Wullings and T. Noij (May 25, 2001b). High-performance liquid chromatography--ToxPrint: chromatographic analysis with a novel (geno)toxicity detection. *Journal of Chromatography. A.* 918:277-291.
- Bouzige, M., G. Machtalère, P. Legeay, V. Pichon and M.C. Hennion (1999). New methodology for a selective on-line monitoring of some polar priority industrial chemicals in waste water. *Waste Management* 19:171-180.
- Brouwer, E.R., I. Liska, R.B. Geerdink, P.C.M. Frintrop, W.H. Mulder, H. Lingeman and U.A.Th. Brinkman (1991) Determination of polar pollutants in water using an on-line dual-precolum liquid chromatographic preconcentration system, *Chromatographia* 32:445-452.
- Carneiro, M.C., L. Puignou and M.T. Galceran. (2000). Comparison of silica and porous graphitic carbon as solid-phase extraction materials for the analysis of cationic herbicides in water

- by liquid chromatography and capillary electrophoresis, *Analytica Chimica Acta* 408: 263-269.
- Castillo, M. and D. Barceló (12 January 2001a). Characterisation of organic pollutants in textile wastewaters and landfill leachate by using toxicity-based fractionation methods followed by liquid and gas chromatography coupled to mass spectrometric detection. *Analytica Chimica Acta* 426:253-264.
- Castillo, M., D. Barceló (1 September 1999a). Identification of polar toxicants in industrial wastewaters using toxicity-based fractionation with liquid chromatography/mass spectrometry. *Analytical Chemistry* 71:3769-3776.
- Castillo, M., D. Barceló, A.S. Pereira and F.R. Aquino Neto (January 1999b). Characterization of organic pollutants in industrial effluents by high-temperature gas chromatography-mass spectrometry. *TrAC Trends in Analytical Chemistry* 18:26-36.
- Castillo, M., D. Puig and D. Barceló. (1997). Determination of priority phenolic compounds in water and industrial effluents by polymeric liquid-solid extraction cartridges using automated sample preparation with extraction columns and liquid chromatography; Use of liquid-solid extraction cartridges for stabilization of phenols, *Journal of Chromatography A*. 778:301-311.
- Castillo, M., E. Martínez, A. Ginebreda, L. Tirapu, D. Barceló (October 2000). Determination of non-ionic surfactants and polar degradation products in influent and effluent water samples and sludges of sewage treatment plants by a generic solid-phase extraction protocol. *The Analyst* 125:1733-1739.
- Castillo, M., F. Ventura and D. Barceló (April 1999c). Sequential solid phase extraction protocol followed by liquid chromatography-atmospheric pressure chemical ionization-mass spectrometry for the trace determination of non ionic polyethoxylated surfactants in tannery wastewaters. *Waste Management* 19:101-110.
- Castillo, M., M.C. Alonso, J. Riu, M. Reinke, G. Klöter, H. Dizer, B. Fischer, P. D. Hansen and D. Barceló (12 January 2001b). Identification of cytotoxic compounds in European wastewaters during a field experiment. *Analytica Chimica Acta* 426:265-277.
- Cheung, J. and R.J. Wells. (1997). Analysis of phenolic compounds in effluent by solid-phase extraction and gas chromatography-mass spectrometry with direct on-column benzylolation: sensitive negative ion chemical ionisation gas chromatography-mass spectrometry detection of phenyl benzyl ethers, *Journal of Chromatography A*. 771: 203-211.
- Crescenzi, C., A. Di Corcia, G. Passariello, R. Samperi, M.I.T. Carou. (1996). Evaluation of two new examples of graphitized carbon blacks for use in solid-phase extraction cartridges, *Journal of Chromatography A*. 733:41-55.
- De Bruyn, A.M.H. and J. B. Rasmussen (2001). Freeze concentration of ambient waters for toxicity testing. *Environmental toxicology and chemistry* 20:1733-1739.
- De Maagd. P.G.J., M. Tonkes (1999). Selection of genotoxicity tests for risk assessment of effluents. *Environmental Toxicology* 15:81-90.
- Di Corcia, A., A. Bellioni, M. Diab Madbouly and S. Marchese. (1996). Trace determination of phenols in natural waters Extraction by a new graphitized carbon black cartridge followed by liquid chromatography and re-analysis after phenol derivatization, *Journal of Chromatography A*. 733:383-393.
- Di Corcia, A., R. Samperi, A. Marcomini, S. Stelluto, (1993). *Analytical Chemistry* 65:907-
- Fernandez-Alba, A.R., A. Agüera, M. Contreras, G. Peñuela, I. Ferrer and D. Barceló (1998). Comparison of various sample handling and analytical procedures for the monitoring of pesticides and metabolites in ground waters, *Journal of Chromatography A*. 823:35-47.

- Ferrer, I. and D. Barceló. Validation of new solid-phase extraction materials for the selective enrichment of organic contaminants from environmental samples, Pages 180-192.
- Furlong, E.T., M.R. Burkhardt, P.M. Gates, S.L. Werner and A. W. Battaglin (5 April 2000). Routine determination of sulfonyleurea, imidazolinone, and sulfonamide herbicides at nanogram-per-liter concentrations by solid-phase extraction and liquid chromatography/mass spectrometry. *The Science of The Total Environment* 248:135-146.
- Galassi, S. and E. Benfenati (11 August 2000). Fractionation and toxicity evaluation of waste waters. *Journal of Chromatography A* 889:149-154.
- Guenu, S. and M.C. Hennion. (1996a). Evaluation of new polymeric sorbents with high specific surface areas using an on-line solid-phase extraction-liquid chromatographic system for the trace-level determination of polar pesticides, *Journal of Chromatography A*. 737: 15-24.
- Guenu, S. and M.C. Hennion. (1996b). Prediction from liquid chromatographic data of obligatory backflush desorption from solid-phase extraction cartridges packed with porous graphitic carbon, *Journal of Chromatography A*. 725:57:66.
- Heberer, T. and H.J. Stan. (1997). Detection of more than 50 substituted phenols as their t-butyl dimethylsilyl derivatives using gas chromatography-mass spectrometry, *Analytica Chimica Acta* 341: 21-34.
- Hennion, M.C. and V. Coquart (1993). Comparison of reversed-phase extraction sorbents for the on-line trace enrichment of polar organic compounds in environmental aqueous samples. *Journal of chromatography A* 642:211-224.
- Hennion, M.C., V. Coquart, S. Guenu and C. Sella. (1995). Retention behaviour of polar compounds using porous graphitic carbon with water-rich mobile phases, *Journal of Chromatography A*. 712:287-301.
- Hennion, M.C., C. Cau-Dit-Coumes and V. Pichon. (1998). Trace analysis of polar organic pollutants in aqueous samples; Tools for the rapid prediction and optimisation of the solid-phase extraction parameters, *Journal of Chromatography A*. 823:147-161.
- Hongjun J., Y. Xuan, Y. Daqiang and Y. Hongxia (1999). A Case Study on Identifying the Toxicant in Effluent Discharged from a Chemical Plant. *Marine Pollution Bulletin* 39:122-125.
- Ingelse, B.A., R.C.J. van Dam, R. J. Vreeken, H.G.J. Mol and O.M. Steijger. (2001). Determination of polar organophosphorus pesticides in aqueous samples by direct injection using liquid chromatography-tandem mass spectrometry, *Journal of Chromatography A*. 918: 67-78
- Junker-Buchheit, A. and M. Witznbacher. (1996). Pesticide monitoring of drinking water with the help of solid-phase extraction and high-performance liquid chromatography, *Journal of Chromatography A*. 737: 67-74.
- Klinkow, N., J. Oleksy-Frenzel and M. Jekel (September 1998). Toxicity-directed fractionation of organic compounds in tannery wastewater with regard to their molecular weight and polarity. *Water Research* 32:2583-2592.
- la Farré, M., I. Ferrer, A. Ginebreda, M. Figueras, L. Olivella, L. Tirapu, M. Vilanova and D. ble 17 August 2001a). Determination of drugs in surface-water and wastewater samples by liquid chromatography-mass spectrometry: methods and preliminary results including toxicity studies with *Vibrio fischeri*. *Journal of Chromatography A*, **In Press, Uncorrected Proof**.
- la Farré, M., M.J. García, L. Tirapu, A. Ginebreda and D. Barceló (26 January 2001b). Wastewater toxicity screening of non-ionic surfactants by Toxalert[®] and Microtox[®] bioluminescence inhibition assays. *Analytica Chimica Acta* 427:181-189.

- la Farré, M., M.J. García, M. Castillo, J. Riu and D. Barceló (February 2001c). Identification of surfactant degradation products as toxic organic compounds present in sewage sludge. *J. Environ Monit* 3:232-237.
- la Farré, M., O. Pasini, M. Carmen Alonso, M. Castillo and D. Barceló (12 January 2001d). Toxicity assessment of organic pollution in wastewaters using a bacterial biosensor. *Analytica Chimica Acta* 426:155-165.
- Lagana, A., G. Fago, A. Marino, (1998). Determination of aryloxyphenoxypropionic acid herbicides in water using different solid-phase extraction procedures and liquid chromatography-diode array detection, *Journal of Chromatography A*. 796:309-318.
- Ling, Y., Y. Hongxia, Y. Daqiang and J. Hongjun (June 1999). Application of the simplified toxicity identification evaluation procedures to a chemical works effluent. *Chemosphere* 38:3571-3577.
- Liska, I., A. Kuthan and J. Krupcik (1990). Comparison of sorbents for solid-phase extraction of polar compounds from water. *Journal of chromatography A* 509:123-134.
- Liska, I., E.R. Brouwer, A.G.L. Ostheimer, H. Lingeman, U.A.Th. Brinkman, R.B. Geerdink and W.H. Mulder (1992) Rapid screening of a large group of polar pesticides in river water by on-line trace enrichment and column liquid chromatography, *Intern. J. Environ. Anal. Chem.* 47:267-291.
- Loon van W., 1996. Een analytisch-chemisch classificatieschema voor opgeloste organische stoffen in water: selectie van testset voor een integral analytisch schema voor opgeloste organische stoffen in opperflaktewater en drinkwater.
- López de Alda, M.J. and D. Barceló (Available online 5 September 2001). Use of solid-phase extraction in various of its modalities for sample preparation in the determination of estrogens and progestogens in sediment and water. *Journal of Chromatography A*, **In Press, Uncorrected Proof**.
- Mallat, E. and D. Barceló. (1998). Analysis and degradation study of glyphosate and of aminomethylphosphonic acid in natural waters by means of polymeric and ion-exchange solid-phase extraction columns followed by ion chromatography- post-column derivatization with fluorescence detection, *Journal of Chromatography A*. 823:129-136.
- Martínez, R.C., E. Rodríguez Gonzalo, M^a E. Fernández Laespada and F. J. Sánchez San Román (2000). Evaluation of surface- and ground-water pollution due to herbicides in agricultural areas of Zamora and Salamanca (Spain). *Journal of Chromatography A* 869:471-480.
- Masqué, N., M. Galià, R.M. Marcé and F. Borrull (1998a). New chemically modified polymeric resin for solid-phase extraction of pesticides and phenolic compounds from water, *Journal of Chromatography A*. 803:147-155.
- Masqué, N., M. Galià, R.M. Marcé and F. Borrull (1997). Chemically modified polymeric resin used as sorbent in a solid-phase extraction process to determine phenolic compounds in water, *Journal of Chromatography A*. 771:55-61.
- Masqué, N., R.M. Marcé and F. Borrull (1998b). New polymeric and other types of sorbents for solid-phase extraction of polar organic micropollutants from environmental water. *Trends in analytical chemistry* 17(6):384-394.
- Masqué, N., R.M. Marcé and F. Borrull (1998c). Comparison of different sorbents for on-line solid-phase extraction of pesticides and phenolic compounds from natural water followed by liquid chromatography, *Journal of Chromatography A*. 793: 257-263.
- Meylan, W.M., P.H. Howard and R.S. Boethling (1996). Improved method for estimating water solubility from octanol/water partition coefficient, *Environmental Toxicology and Chemistry* 15:100-106.

- Molina, C., P. Grasso, E. Benfenati and D. Barceló. (1996). Automated sample preparation with extraction columns followed by liquid chromatography-ionspray mass spectrometry Interferences, determination and degradation of polar organophosphorus pesticides in water samples, *Journal of Chromatography A*. 737:47-58.
- Nakamura, M., Ma. Nakamura, and S. Yamada (1996). Conditions for solid-phase extraction of agricultural chemicals in waters by using n-octanol-water partition coefficients, *Analyst* 121:469-475.
- Norberg, J., J. Slobodnik, R.J.J. Vreuls and U.A.Th. Brinkman (1995). On-line solid-phase extraction-liquid chromatography for screening and quantification of organophosphorus pesticides in surface water, *Analytical Methods & Instrumentation* 2:266-276.
- Pichon, V., C.C. Dit Coumes, L. Chen, S. Guenu and M.C. Hennion. (1996). Simple removal of humic and fulvic acid interferences using polymeric sorbents for the simultaneous solid-phase extraction of polar acidic, neutral and basic pesticides, *Journal of Chromatography A*. 737:25-33.
- Pichon, V., L. Chen, S. Guenu, M.-C. Hennion (1995). Comparison of sorbents for the solid-phase extraction of the highly polar degradation products of atrazine (including ammeline, ammelide and cyanuric acid). (1995). *Journal of Chromatography A*. 711:257-267.
- Pichon, V., M. Charpak, M. -C. Hennion (1998). Multiresidue analysis of pesticides using new laminar extraction disks liquid chromatography and application to the French priority list, *Journal of Chromatography A*. 795:83-92.
- Pichon, V., M.-C. Hennion (1994). Determination of pesticides in environmental waters by automated on-line trace-enrichment and liquid chromatography, *Journal of Chromatography A*. 665:269-281.
- Pocurull, E., M. Calull, R.M. Marcé and F. Borrull. (1996). Determination of phenolic compounds at low $\mu\text{g l}^{-1}$ levels by various solid-phase extractions followed by liquid chromatography and diode-array detection, *Journal of Chromatography A*. 719:105-112.
- Puig, D. and D. Barceló. (1996). Comparison of different sorbent materials for on-line liquid-solid extraction followed by liquid chromatographic determination of priority phenolic compounds in environmental waters, *Journal of Chromatography A*. 733:371-381.
- Puig, D. and D. Barceló. (1997). Determination of polar priority phenols at parts per trillion levels in water using on-line liquid-solid extraction followed by liquid chromatography with coulometric detection, *Journal of Chromatography A* 778:313-319.
- Putschew, A., S. Schittko and M. Jekel (September 2001). Quantification of triiodinated benzene derivatives and X-ray contrast media in water samples by liquid chromatography-electrospray tandem mass spectrometry. *Journal of Chromatography A* 930:127-134.
- Reemtsma, T. (2001). Prospects of toxicity-directed wastewater analysis. *Analytica Chimica Acta* 426:279-287.
- Reemtsma, T., A. Putschew and M. Jekel (1999). Industrial wastewater analysis: a toxicity-directed approach. *Waste Management* 19:181-188.
- Reifferscheid, G., J. Heil. (1996). *Mutat. Res.* 369:129-145.
- Rodríguez, I., M.C. Mejuto, M.H. Bollaín and R. Cela. (1997). Evaluation of two solid-phase extraction procedures for the preconcentration of chlorophenols in drinking water, *Journal of Chromatography A*. 786:285-292.
- Sabik, H., R. Jeannot, (1998). Determination of organonitrogen pesticides in large volumes of surface water by liquid-liquid and solid-phase extraction using gas chromatography with nitrogen-phosphorus detection and liquid chromatography with atmospheric pressure chemical ionisation mass spectrometry, *Journal of Chromatography A*. 818:197-207.

- Sakai, M. (2001). Chronic toxicity tests with *Daphnia magna* for examination of river water quality. *Journal of Environmental Science and Health. Part. B, Pesticides, Food Contaminants, and Agricultural Wastes* 36:67-74.
- Sánchez, D.M., R. Martín, R. Morante, J. Marín and M.L. Munuera (2000). Preconcentration speciation method for mercury compounds in water samples using solid phase extraction followed by reversed phase high performance liquid chromatography. *Talanta* 52:671-679.
- Sicilia, D., S. Rubio, D. Pérez-Bendito, N. Maniasso and E.A.G. Zagatto. (1999). Anionic surfactants in acid media: a new cloud point extraction approach for the determination of polycyclic aromatic hydrocarbons in environmental samples, *Analytica Chimica Acta* 392:29-38.
- Stahl, R.G., (1991). Genotoxicity of organic compounds in water, *Ecotox. Environ. Saf.* 22:94-125.
- Steen, R. J. C. A., I. Bobeldijk and U.A. Th. Brinkman. (2001). Screening for transformation products of pesticides using tandem mass spectrometric scan modes, *Journal of Chromatography A*. 915:129-137
- Thomas, K.V., M.R. Hurst, P. Matthiessen, D. Sheahan and R. J. Williams (July 2001). Toxicity characterisation of organic contaminants in stormwaters from an agricultural headwater stream in South East England. *Water Research* 35:2411-2416.
- Thomas, K.V., R.E. Benstead, J.E. Thain and M.J. Waldock (October 1999). Toxicity Characterization of Organic Contaminants in Industrialized UK Estuaries and Coastal Waters. *Marine Pollution Bulletin* 38:925-932.
- Tölgyessy, P. and I. Liška (1999). Study of applicability of various solid-phase extraction materials for sample handling in screening analysis of organic micropollutants in water. *Journal of Chromatography A* 857:247-254.
- Tomkins, B.A., and W.H. Griest, (1996). Determination of N-nitrosodimethylamine at part-per-trillion concentrations in contaminated groundwaters and drinking waters featuring carbon-based membrane extraction disks, *Analytical Chemistry*, 68:2533-2540.
- Toribio, L., M&z. ausco., J. del Nozal, J.L. Bernal, J.J. Jiménez and M&z. ausco., L. Serna. (1998). Packed-column supercritical fluid chromatography coupled with solid-phase extraction for the determination of organic microcontaminants in water, *Journal of Chromatography A*. 823: 163-170.
- Vreeken, R.J., P. Speksnijder, I. Bobeldijk-Pastorova and Th. H. M. Noij. (1998). Selective analysis of the herbicides glyphosate and aminomethylphosphonic acid in water by on-line solid-phase extraction-high-performance liquid chromatography-electrospray ionization mass spectrometry, *Journal of Chromatography A*. 794:187-199.
- Wells, M.L. and W.K Bruland (1998). An improved method for rapid preconcentration and determination of bioactive trace metals in seawater using solid phase extraction and high resolution inductively coupled plasma mass spectrometry. *Marine Chemistry* 63:145-153.
- Wernersson, A., G. Dave (April 1997). Phototoxicity identification by solid phase extraction and photoinduced toxicity to *Daphnia magna*. *Archives of Environmental Contamination and Toxicology* 32:268-273.
- Zamperlini, G C., M. Santiago-Silva and W. Vilegas (2000). Solid-phase extraction of sugar cane soot extract for analysis by gas chromatography with flame ionisation and mass spectrometric detection. *Journal of Chromatography A* 889:281-286.

Chapter 3

- Altschuh, J., R. Bruggemann, H. Santl, G. Eichinger, and O.G. Piringner (1999). Henry's law constants for a diverse set of organic chemicals: experimental determination and comparison of estimation methods. *Chemosphere* 39: 1871-87.

- Anon. (1995). AQUAPOL, chemical database, Institute for Inland Water Management and Waste Water Treatment, Lelystad, The Netherlands and Institute for Coast and Sea, Den Haag, The Netherlands, 1995.
- Boethling, R.S and Mackay, D. (eds.) (2000). *Handbook of property estimation methods for chemicals – Environmental Health Sciences*. Lewis Publishers, New York (USA).
- Daylight (1995). CLOGP Program. Daylight Chemical Information Systems. Von Karman Ave., Irvine, CA 92715.
- EPA (2000). Pollution Prevention (P2) framework. Report nr. EPA-748-B-00-001. US Environmental Protection Agency, Office of Pollution Prevention and Toxics, Washington, DC (USA).
- Hansch, C and A.J. Leo (1979). *Substituent Constants for Correlation Analysis in Chemistry and Biology*; Wiley: New York, 1979.
- Hine, J. and P.K. Mookerjee (1975). The intrinsic hydrophilic character of organic compounds. Correlations in terms of structural contributions. *J. Org. Chem.* 40: 292-298.
- Howard, P. (1986). BIOLOG, BIODEG, AND FATE/EXPOS: New Files on Microbial Degradation and Toxicity as well as Environmental Fate/Exposure of Chemicals. *Environ. Toxicol. Chem.* 5:977-988.
- Leo, A. (2000). Octanol/water partition coefficient. In: Boethling, R.S and Mackay, D. (eds.). *Handbook of property estimation methods for chemicals – Environmental Health Sciences*. Lewis Publishers, New York (USA). pp. 89-114
- Leo, A. (2000). Octanol/water partition coefficient. In: Boethling, R.S and Mackay, D. (eds.). *Handbook of property estimation methods for chemicals – Environmental Health Sciences*. Lewis Publishers, New York (USA). pp. 89-114
- Leo, A.J. (1993). Calculating log Poct from structures. *Chemical Reviews*, 93:1281-1306.
- Leo, A.J. (1993). Calculating log Poct from structures. *Chemical Reviews*, 93:1281-1306.
- Leo, D. (1995). Computer Programme ClogP, Biobyte Corp., Claremont, CA.
- Leo, D. (1995). Computer Programme ClogP, Biobyte Corp., Claremont, CA.
- Mackay D. W.Y. Shiu, K.C. Ma (1997). *Illustrated handbook of physical-chemical properties of environmental fate for organic chemicals*. Springer-Verlag, Berlin/Heidelberg.
- Mackay, D., W.Y. Shiu and K.C. Ma (2000). Henry's Law constant. In: Boethling, R.S and Mackay, D. (eds.). *Handbook of property estimation methods for chemicals – Environmental Health Sciences*. Lewis Publishers, New York (USA). pp. 69-88.
- Meylan, W., P. Howard (1999a). *User's guide for EPIWIN^c*. EPI Suite - Estimation Programs Interface for Microsoft Windows. Syracuse Research Corporation, North Syracuse (NY).
- Meylan, W., P. Howard (1999b). *User's guide for HENRYWIN^c* Version 3, Henry's law constant program for Microsoft Windows. Syracuse Research Corporation, North Syracuse (NY).
- Meylan, W., P. Howard (1999c). *User's guide for KOWWIN^c* Version 1.6, Log Octanol-Water Partition Coefficient Estimation Program for Microsoft Windows. Syracuse Research Corporation, North Syracuse (NY).
- Meylan, W.M. and P.H. Howard.(1995). Atom/fragment contribution method for estimating octanol-water partition coefficients. *J. Pharm. Sci.* 84: 83-92.
- Meylan, W. and Howard, P.H.; *Environ. Toxicol. Chem.* 10: 1283_93, 1991
- Nirmalakhandan, N., R.A. Brennan and R.E. Speece (1997). Predicting Henry's law constant and the effect of temperature on Henry's law constant. *Wat. Res.* 31(6): 1471-1481.

- Reifferscheid, G. and J. Heil (1996) Validation of the SOS/umu test using test results of 486 chemicals and comparison with the Ames test and carcinogenicity data. *Mutation Research* 369:129-145.
- TGD (1996). *Technical guidance document in support of commission directive 93/67/EEC on risk assessment for new notified substances and commission regulation (EC) No 1488/94 on risk assessment for existing substances*. Office for Official Publications of the European Communities, Luxembourg.
- Thomas, R. (1990). Volatilization from water. In: *Handbook of Chemical Property Estimation Methods. Environmental Behavior of Organic Chemicals*. Lyman, W.J, Reehl, W.F., and Rosenblatt, D.H. (editors). Washington DC: Amer Chem Soc pp. 15-1 to 15-30.

Chapter 4

- Adahchour, M., P.E.G. Leonards, U.A.Th. Brinkman (2000a). *Procedures for fractionation and clean-up techniques of environmental samples; application for toxicity tests and analytical determinations*. Report no. ACAS-002; Department of Analytical Chemistry and Applied Spectroscopy, Vrije Universiteit, Amsterdam.
- Adahchour, M., P.E.G. Leonards, U.A.Th. Brinkman (2000b). *Procedures for fractionation and clean-up techniques of environmental samples; application for toxicity tests and analytical determination – Part 2*. Report no. ACAS-003; Department of Analytical Chemistry and Applied Spectroscopy, Vrije Universiteit, Amsterdam.
- Bobbeldijk, I., A. Brandt, B. Wullings and Th. Noij (2001). High-performance liquid chromatography-Toxprint: chromatographic analysis with a novel (geno)toxicity detection. *J. Chrom. A* 918: 277-291.
- Brinkman, U.A.T. and R.J.J. Vreuls (1996) Solid-phase extraction for on-line sample treatment in capillary gas chromatography. *Lc Gc-Magazine Of Separation Science* 14:581-585.
- EPA (2000). Pollution Prevention (P2) framework. Report nr. EPA-748-B-00-001. US Environmental Protection Agency, Office of Pollution Prevention and Toxics, Washington, DC (USA).
- Van Stee, L.L.P., P.E.G. LEONARDS, R.J.J. Vreuls, and U.A.T. Brinkman (1999) Identification of non-target compounds using gas chromatography with simultaneous atomic emission and mass spectrometric detection (GC-AED/MS): analysis of municipal wastewater. *Analyst* 124:1547-1552.
- Vreuls, J.J., A.J.H. Louter, and U.A.T. Brinkman (1999) On-line combination of aqueous-sample preparation and capillary gas chromatography. *Journal Of Chromatography A* 856:279-314.

Annex section

List of Annex Tables

- A1. Factsheets of SPE techniques
- A2. UmuC responsive compounds included in the survey, occurrence in surface water surveys by RIZA
- A3. Predicted and measured values for Kow
- A4. Predicted and measured values for H

Annex 1 Factsheets

Example and explanation of information provided in the factsheets.

0 <i>Material</i>	Extraction technique (SPE) material used (supplier)
<i>Amount material (mg) and size</i>	Amount of (SPE) sorbents and size of the used cartridges
<i>Deactivation</i>	Conditioning of the (sorbent) cartridge
<i>Kind of sample (ml)</i>	Amount and kind of the analysed water samples
<i>Breakthrough volumes (ml)</i>	Are there breakthrough volumes reported for analytes under the used method and conditions?
<i>(Amount) eluent (ml)</i>	Amount and kind of desorption solvent
<i>Fraction (nr.)</i>	Hoe many fractions are there after desorption?
<i>Drying method</i>	Hoe is the obtained eluate further evaporated (enriched)?
<i>Enrichment factor</i>	Which enrichment factors are obtained?
<i>Compounds</i>	Which (target) compounds are analysed?
<i>Recovery</i>	Under the above conditions what are the recoveries for these compounds?
<i>Repeatability</i>	Is the method repeatable?
<i>LOD's</i>	What are the Limit Of Detections (LOD) obtained by this method
<i>Reference</i>	Article
<i>Remarks</i>	Additional remarks, which can be relevant for the selection of (SPE) extraction method.

1	Sequential SPE (SSPE)
Material	C18 + Isolute ENV+/LiChrolut EN in serie
<i>Amount material (mg) and size</i>	C18 (500 mg, 6 ml) + LiChrolut EN/Isolut ENV+ (200 mg, 6 ml)
<i>Deactivation</i>	7 ml MeOH + 3 ml water (1 ml/min)
<i>Kind of sample (ml)</i>	all 200 ml (filtered if needed) of raw tanner wastewater three wastewater (effluent) wastewater industrial effluent drinking water
<i>Breakthrough volumes (ml)</i>	?
<i>(Amount) eluent (ml)</i>	1. 2x5 ml hexane 2. 2x5 ml CH ₂ Cl ₂ /hexane (4:1) 3. 2x5 ml methanol/CH ₂ Cl ₂ (9:1) 4. 1 ml water (5mM TEA + 5mM acetic acid, pH6.5) 5. 9 ml methanol
<i>Fraction (nr.)</i>	1, 2, 3, 4, 5
<i>Drying method</i>	Under gentle stream of N ₂
<i>Enrichment factor</i>	Up to 1000
<i>Compounds</i>	Linear alkyl benzene-sulfonate (LAS) + Non-ionic polyethoxylated surfactant (PEO) + benzene-and naphthalene-sulfonates (BS and NPS) C10LAS; C11LAS; C12LAS; C13LAS; 1,5-NPDS; 2,6-NPDS; 2,7-NPDS; 1-hydroxy-3,6-NPDS; 1-amino-5-NPS; 1-amino-6- NPS; 1-amino-7-NPS; 1-hydroxy-1-amino-3-NPS; 3-nitro-BS; 1- hydroxy-4-NBS; 4-methyl-BS; 4-chloro-BS; 1-NPS; 2-NPS; 1- amino-6-NPS; 1-amino-7-NPS; catechol; phenol; 4- methylphenol; 2,4-dinitrophenol; 2,2'-dihydroxybiphenyl; penta- chlorophenol; naphthol; poly(ethylene glycol); nonylphenol polyethoxylate; alcohol ethoxylate; dibutyl phthalate; dimethyl phthalate; bis-2-ethylexyl phthalate; tributylphosphate; ethylben- zoate; 1-methyl-2-pyrrolidinone; 2-methylbenzenesulfonamide; 2,2'-dimethyl-1,3-propanediol; 4-chlorophenol; 2-chlorophenol; 2,4-dichlorophenol; 2,4,6-trichlorophenol; pentachlorophenol; Isoproturon; diuron; simazine; atrazine; 4-chloro-m-cresol;
<i>Recovery</i>	72-103%
<i>Repeatability</i>	1-11%
<i>LOD's</i>	Low ppb level
<i>Reference</i>	M. Castillo et al./ Analytica chimica acta 426 (2001) 265 M. Castillo et al./ Waste Management 19 (1999) 101 M. Castillo et al./ Analytica chimica acta 426 (2001) 253 M. Castillo et al./ Trends in Analytical Chemistry 18 (1999) 26 M. Castillo et al./ Analytical Chemistry 71 (1999) 3769 K.V. Thomas et al./ Wat. Res. 35 (2001) 2411 M. la Farre et al./ Analytica Chimica Acta 426 (2001) 155 M. la Farre et al./ Analytica Chimica Acta 427 (2001) 181 Th. Reemtsma et al./ Waste Management 19 (1999) 181
<i>Remarks</i>	The developed methodology with toxicity measurements

allowed to detect different groups of pollutants responsible for the toxicity of the studied wastewaters. This method can cover broad log Kow range (0-6) efficiently and no problems are reported compared to other methods (GCB).

2	SPE
Material	Lichrolut EN (Merck)
<i>Amount material (mg) and size</i>	200 mg
<i>Deactivation</i>	5 ml methanol + 5 ml HPLC-grade water
<i>Kind of sample (ml)</i>	1000 for atrazine 100 for degradation experim.
<i>Breakthrough volumes (ml)</i>	?
<i>Amount eluent (ml)</i>	3x3 ml methanol
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	evaporated to 0.5 under N ₂
<i>Enrichment factor</i>	200–2000
<i>Compounds</i>	atrazine, hydroxyatrazine, desethylatrazine, desisopropylatrazine
<i>Recovery</i>	?
<i>Repeatability</i>	?
<i>LOD's</i>	?
<i>Reference</i>	R.J.C.A. Steen et.al./ J. chromatogr. A 915 (2001) 129
<i>Remarks</i>	The rather polar transformation products (TPs) are retained on the sorbent used (LiChrolut EN). Fortunately, this material has already been shown to be a good choice to retain relatively polar TPs and even partially ionized compounds (if necessary, after pH adjustment)
3	SPE
Material	Isolute ENV+ (International Sorbent Technology, Ltd)
<i>Amount material (mg) and size</i>	200 mg, 850 Å pore size, 1100 surface area, 1 ml cartridge
<i>Deactivation</i>	MeOH, C ₆ , acetone (9 ml for each solvent) and 30 ml of Milli-Q water
<i>Kind of sample (ml)</i>	4–5 l filtered surface water samples by pH 7 (first column) and then the same procedure by pH 2 (second column) (30–40 ml/min flow rate)
<i>Breakthrough volumes (ml)</i>	-
<i>Amount eluent (ml)</i>	9 ml of C ₆ /acetone (85/15,v/v) followed by 6 ml of MeOH
<i>Fraction (nr.)</i>	2
<i>Drying method</i>	The obtained extracts from the two columns (pH 2.0, 7.0) were combined and the volume was reduced using a Kuderna–Danish evaporator with atmospheric air and a temperature of 30°C.
<i>Enrichment factor</i>	300–1500
<i>Compounds</i>	Not reported (only toxicity is studied)
<i>Recovery</i>	-
<i>Repeatability</i>	-
<i>LOD's</i>	-
<i>Reference</i>	A. Baun et al./Environmental Pollution 102 (1998) 185
<i>Remarks</i>	

4	SPE
Material	Oasis (Waters)
<i>Amount material (mg) and size</i>	?
<i>Deactivation</i>	?
<i>Kind of sample (ml)</i>	500 surface water filtered over 0.2 um
<i>Breakthrough volumes (ml)</i>	?
<i>Amount eluent (ml)</i>	4x2 ml acetonitrile (after 1 ml of 5% acetonitrile in water and 30 min drying with N ₂ evaporated to 0.5 ml under N ₂)
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	-
<i>Enrichment factor</i>	1000
<i>Compounds</i>	metribuzine, pirimicarb, diuron, isoproturon, metamitron, atrazine
<i>Recovery</i>	-
<i>Repeatability</i>	-
<i>LOD's</i>	lower than 0.25 ul/l
<i>Reference</i>	I. Bobeldijk et. Al./ j. chromatgr. A 929 (2001) 63
<i>Remarks</i>	Further research concerning alternative solvents, isolation, overall sensitivity and other parameters is needed to overcome the matrix and other interference effects.

5	SPE + HPLC
Material	Oasis (Waters) + pellicular C18 guard column + inertsil ODS-2 analytical Column
<i>Amount material (mg) and size</i>	?
<i>Deactivation</i>	ultrapure water + acetonitril (gradient, 181 ml)
<i>Kind of sample (ml)</i>	all filtered over 0.2 um 100 ml surface water 100 ml municipal water 20 ml industrial water 100 ml hospital water
<i>Breakthrough volumes (ml)</i>	?
<i>Amount eluent (ml)</i>	58 ml linear gradient of acetonitrile and water (58 ml)
<i>Fraction (nr.)</i>	(1 min fraction) and fractions 11-46 are collected
<i>Drying method</i>	-
<i>Enrichment factor</i>	-
<i>Compounds</i>	4-nitroquinoline-N-oxide and 2-aminoanthracene (test compounds)
<i>Recovery</i>	75-125%
<i>Repeatability</i>	?
<i>LOD's</i>	low ug/l range
<i>Reference</i>	I. Bobeldijk et. Al./ j. chromatgr. A 918 (2001) 277
<i>Remarks</i>	Oasis sorbent with an adsorption capacity for both lipophilic and hydrophilic compounds was experimentally shown to be capable to extract compounds with a broad range of polarities.

6	SPE
<i>Material</i>	SDB (Baker) Oasis HLB (Waters) Envicarb (Supelco) Carbograph (Alltech)
<i>Amount material (mg) and size</i>	200 200 500 150
<i>Deactivation</i>	all cartridges 2x 5 ml desorption solvent, 2x 5ml methanol and 2x 5ml demi-water
<i>Kind of sample (ml)</i>	100 ml surface water
<i>Breakthrough volumes (ml)</i>	Yes!! (see recovery)
<i>Amount eluent (ml)</i>	2x5 ml ethyl acetate 2x5 ml ethyl acetate 2x5 ml ethyl acetate DCM/MeOH (8:2) DCM/MeOH (8:2)
<i>Fraction (nr.)</i>	1, 1, 1, 1
<i>Drying method</i>	?
<i>Enrichment factor</i>	1000
<i>Compounds</i>	Acephate, methamidophos, omethoate, oxydemeton-methyl vamidothion
<i>Recovery</i>	acceptable recoveries for last two, but the first three (more polar) low recoveries are obtained
<i>Repeatability</i>	better than 20%
<i>LOD's</i>	0.01-0.03 ug/l depending on the analyte
<i>Reference</i>	B.A. Ingelse et.al./j. chromatogr. A 918 (2001) 67-78
<i>Remarks</i>	Very polar (Log Kow < 0) OPs like acephate, metamidophos and omethoate can not be extracted from water using the currently commonly available SPE cartridges. But Lichrolut RP- 18, Lichrolut EN and Isolut cartridges were not tested.

7	SPE
<i>Material</i>	1. Lichrolut EN (Merck) and 2. Oasis (Waters)
<i>Amount material (mg) and size</i>	1. 200 2. 300
<i>Deactivation</i>	1. Hexane, acetone, HPLC water pH2 (6, 6, 6 ml res.) 2. Methanol, HPLC water (6, 6 ml resp)
<i>Kind of sample (ml)</i>	1. 1 l of filtered ground water (pH2) 2. 1 l of filtered ground water (washing with 1ml 5% MeOH in water)
<i>Breakthrough volumes (ml)</i>	?
<i>Amount eluent (ml)</i>	1. Acetone, methanol, acetone (1, 2, 2 ml resp.) 2. CH ₃ CN:MeOH(70:30), CH ₃ CN:MeOH(70:30) [3 and 2 ml resp.(solvents acidified to pH 3.7)]
<i>Fraction (nr.)</i>	1. 1 2. 1
<i>Drying method</i>	1. Gentle stream of N ₂ to 100 ul and then reconstituted with MeOH to 300 ul. 2. Gentle stream of N ₂ to 100 ul and then reconstituted with 500 ul formate buffer
<i>Enrichment factor</i>	
<i>Compounds</i>	salicylic acid, ketoprofen, naproxen, diclofenac-Na, ibuprofen and gemfibrozil
<i>Recovery</i>	1. 76, 69, 83, 85, 91, and 70 (pH2) 2. <10, 108, 125, 38, and 109 (pH7)
<i>Repeatability</i>	1. 13, 15, 10, 6, 11, and 6 (%SD) 2. Nd, 4.8, 11.6, 24, 12.1 and 2.54 (%SD)
<i>LOD's</i>	1. 15, 28, 29, 5, 43, and 56 (ng/l, SIM mode) 2. Nd, 53, 52, 37, 123 and 72 (ng/l, SIM mode)
<i>Reference</i>	M. la Farre et.al./ J. Chromatorgr. A 000 (2001)000
<i>Remarks</i>	The SPE methode (using Lichrolut-EN) was applied to 12 surface water sample and three effluent samples. While pharmaceuticals are present, and their concentrations increase with observed toxicity, their contribution to the whole toxicity cannot be determined because these substances are acting as tracers of toxicity.

8	SPE
<i>Material</i>	1. Isolute ENV 2. LiChrolut EN (Merck) 3. LiChrolut RP-18 (Merck) 4. RP-18 (Baker) 5. Oasis HLB (Waters) 6. HySphere-resin-GP (Spark Holland) 7. PLRP-S (Polymer Labs)
<i>Amount material (mg) and size</i>	1. 500 mg, 3 ml 2. 200 mg, 3 ml 3. 500 mg, 3 ml 4. 10 mm x 2 mm 5. 10 mm x 2 mm 6. 10 mm x 2 mm 7. 10 mm x 2 mm
<i>Deactivation</i>	1. CH ₃ CN, MeOH and LC-grade water (7, 5 and 5 resp.) 2. CH ₃ CN, MeOH and LC-grade water (7, 5 and 5 resp.) 3. CH ₃ CN and LC-grade water (4 and ml resp.) 4. CH ₃ CN and LC-grade water (4 and ml resp.) 5. CH ₃ CN and LC-grade water (4 and ml resp.) 6. CH ₃ CN and LC-grade water (4 and ml resp.) 7. CH ₃ CN and LC-grade water (4 and ml resp.)
<i>Kind of sample (ml)</i>	50 to 1000 ml filtered water samples
<i>Breakthrough volumes (ml)</i>	extraction of large volumes of sample (1000ml), HySphere-Resin-GP and Lichrolut cartridges are the best
<i>Amount eluent (ml)</i>	2x4 ml acetonitrile (alle cartridges)
<i>Fraction (nr.)</i>	1 fraction (each cartridge)
<i>Drying method</i>	all extracts were blown down to dryness under N ₂ and reconstituted in MeOH to 0.5 ml.
<i>Enrichment factor</i>	
<i>Compounds</i>	Estriol, estradiol, norethindrone, ethynyl estradiol, estrone, levonogestrel, diethylstilbestrol and progesterone.
<i>Recovery</i>	Lichrolut RP-18 and HySphere-Resin-GP give the best recoveries (87-101%)
<i>Repeatability</i>	< 25% (RSD)
<i>LOD's</i>	1-20 ng/l
<i>Reference</i>	M. J. Lopez de Alda et.al/ J. Chromatogr. A 000(2001) 000
<i>Remarks</i>	the first three cartridges (LiChrolut RP-18, 500 mg was better than LiChrolut EN, 200 mg and Isolut ENV, 500 mg) were used in off-line mode and the last four in on-line mode. In general, on-line mode gives the best results.

9	SPE
<i>Material</i>	1. Lichrolut EN (Merck) 2. Lichrolut RP-18 (Merck) 3. Mix. Of the two
<i>Amount material (mg) and size</i>	1. 200 mg, 3 ml 2. 500 and 1000 mg, 3 ml 3. 100 mg of EN and 100 mg of RP-18
<i>Deactivation</i>	1. 3 ml MeOH 2. 3 ml H ₂ O 3. 3 ml MeOH, 3 ml H ₂ O and 3 ml H ₂ O (pH 4)
<i>Kind of sample (ml)</i>	1 l tap water (pH 4-6) spiked with 33-multicomponent pesticide standard.
<i>Breakthrough volumes (ml)</i>	not observed by 1 l water
<i>Amount eluent (ml)</i>	1. And 2. 2x3 ml MeOH:EtOAc(1:1) 3. 2x3 ml MeOH
<i>Fraction (nr.)</i>	each cartridge gives 1 fraction
<i>Drying method</i>	extracts were blown down to dryness and reconstituted in 1 ml of a mixture of ACN and ammonium acetate (20:80)
<i>Enrichment factor</i>	> 1000
<i>Compounds</i>	Deisopropylatrazine, metamidron, chloridazon, deethylatrazine, crimidine carbetamide, bromacil, simasine, cynazine, deethylterbutylazine, karbutilate, methabenzthiazuron, chlortoluron, atrazine, monolinuron, isoproturon, diuron, metobromuron, metazachlor, methoprotrotryne, dimefuron, sebutylazine, propazine, terbutylazine, linuron, chloroxuron, prometryne, chlorpropham, terbutryne, metolachlor, pencycuron, bifenox, pendimethaline.
<i>Recovery</i>	The best recoveries were obtained by LiChrolut EN (200 mg, 3 ml) and ranging from 81 to 118%.
<i>Repeatability</i>	0.1-13.2% (RSD)
<i>LOD's</i>	0.007 mg/l (S/N ratio of 4:1)
<i>Reference</i>	A.Junker-Buchheit et. Al/ J. Chromatogr. A 737 (1996) 67
<i>Remarks</i>	LiChrolut (200 mg, 3ml) was the best. One gram of LiChrolut RP-18 is not sufficient to achieve quantitative recoveries of the most polar components (low Log K _{ow}). This problem was circumvented by increasing the mass of the sorbent (up to 2000 mg). The method was applied to the analysis of drinking water and surface water. Problems arising from high contents of humic acid compounds were circumvented by using a mixed RP-EN stationary phase.

10	SPE
<i>Material</i>	1. LiChrolut RP-18 (Merck) 2. LiChrolut EN (Merck) 3. Isolut Env+ (IST) 4. PLRP-S (Polymer Labs) 5. Envchrom P (Supelco)
<i>Amount material (mg) and size</i>	10x3 mm I.D 1. 40-63 um 2. 40-63 um 3. 70-100 um 4. 20 um 5. -
<i>Deactivation</i>	10 ml MeOH and 10 ml H ₂ O (deionised)
<i>Kind of sample (ml)</i>	2-40 ml of deionised water, tap water and river water.
<i>Breakthrough volumes (ml)</i>	Varying the deionised volume water from 2 to 40 ml, the results obtained were different depending on the compound, and in the cas of some PAHs (naphthalene, biphenyl, acenaphthylene and chrysene) there was a significant decrease in the recoveries when the volume of sample was higher than 10 ml.
<i>Amount eluent (ml)</i>	?
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	Not needed (necessary)--> on-line mode.
<i>Enrichment factor</i>	-
<i>Compounds</i>	naphthalene, biphenyl, acenaphthene, fluorene, acenaphthylene, phenanthrene, 2-chlorophenol, fluoranthrene, 2,4-dinitrophenol, propham, benz[a]anthracene, 2,4-dichlorophenol, chlorpropham, chrysene, 2,4-dimethylphenol, phenol, terbutylazine, dinitro-o-cresol, 4-chloro-3-methylphenol, napropamide, carbofurane, linuron, chlorbromuron, 4-nitrophenol, desethylatrazine, desisopropylatrazine, aldicarbsulphone, fenuron, desmidepham, phenmedipham, warfarin, chlortoluron, monuron, chloroxuron, metoxuron.
<i>Recovery</i>	The best recoveries were obtained by Isolut Env+ and ranging from 35 to 110%.
<i>Repeatability</i>	Repeatability (RSD %) from 2 to 15 % and reproducibility (RSD %) from 3-14%
<i>LOD's</i>	LOD,s (ug/l, and S/N=10) from 0.4 to 2.6
<i>Reference</i>	L. Toribio et al/ J Chromatogr. A 823 (1998) 163
<i>Remarks</i>	The recoveries of some PAHs were low, probably due to high retention on this kind of sorbent (Isolut Env+), which resulted in an increase of their detection limits.

11	SPE
<i>Material</i>	1. Carbo-pack B 120/400 (Supelco) 2. HYSphere-1 5um (Spark Holland) 3. Bond Elut PPL 125 um (Varian)
<i>Amount material (mg) and size</i>	10x3 mm I.D 1. 120/400 2. 5 um 3. 125 um
<i>Deactivation</i>	ACN and milli-Q water (pH 2.5)
<i>Kind of sample (ml)</i>	different volumes of Milli-Q water, tap water and river water (50-200 ml); pH 2.5.
<i>Breakthrough volumes (ml)</i>	Breakthrough volumes of phenol: 1. 2 ml 2. 22 ml 3. 14 ml For all compounds: HYSphere-1 gives better recoveries varying volumes from 50-200 ml.
<i>Amount eluent (ml)</i>	ACN
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	Not needed (necessary)--> on-line mode.
<i>Enrichment factor</i>	
<i>Compounds</i>	oxamyl, methomyl, phenol, 4-nitrophenol, 2,4-dinitrophenol, bentazol, simazine, chlor- phenoxy acid, atrazine.
<i>Recovery</i>	The best recoveries were obtained by HYSphere-1 and gave a higher breakthrough volume for phenol for 100 ml water sample. Recoveries were between 67% for phenol and 86% for eight polar analytes.
<i>Repeatability</i>	?
<i>LOD's</i>	LOD's were between 0.03 and 0.17 ug/l.
<i>Reference</i>	N. Masque et al/ J. Chromatogr. A 793 (1998) 257
<i>Remarks</i>	Problems arising from high contents of fulvic and humic acids were solved by adding 10% Na2SO3 solution to the sample

12	SPE
<i>Material</i>	1. SDB-1 PS-DVB (Baker) 2. SDB-1 (Baker) 3. PLRP-S (Polymer Labs)
<i>Amount material (mg) and size</i>	8x2 mm I.D 1. 200 mg, 43-123 μ m 2. - 3. -
<i>Deactivation</i>	5 ml methanol + 10 ml HPLC-grade water (pH 3)
<i>Kind of sample (ml)</i>	1-100 ml-1 L LC-grade water, 100 ml-1 L tap water and 100-200 ml river water (all filtered through cellulose ester filters (HA type, diameter 47 mm, pore size 0.45 μ m).
<i>Breakthrough volumes (ml)</i>	Breakthrough volumes are higher than 100 ml for each pesticides when the sample is at Ph 2 using 8x2 mm I.D. SDB-1 column, which, in general, gave the best recoveries.
<i>Amount eluent (ml)</i>	2 ml MeOH:ACN (50:50)
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	In the off-line mode, the extracts were evaporated to dryness at 30C with a gentle stream of N ₂ . The dry extracts were desolved in MeOH -0.005 M phosphate buffer (pH 7) (25:85, v/v). The test of volatility was performed by directly spiking the desorption solution with 50-100 ng of each analyte, then evaporating it to dryness and reconstituting the dry extract in the mixture as above. No loss of analytes was observed.
<i>Enrichment factor</i>	50
<i>Compounds</i>	Log Kow range: -0.5 to 1.7 oxamyl, methomyl, DIA (metabolite of atrazine), monocrotophos, fenuron, metamitron, DEA (metabolite of atrazine), chloridazon, carbendazim, aldicarb, aminocarb, metribuzin, methoxuron
<i>Recovery</i>	The best recoveries were obtained by SDB-1 sorbent and ranging from 75 and 105% for all the three types of water samples (pH 2)
<i>Repeatability</i>	< 10 % (RSD)
<i>LOD's</i>	LOD's were between 0.1 and 0.5 μ g/l depending on the analytes.
<i>Reference</i>	S. Guenu, M-C. Hennion/ J. Chromatogr. A 737 (1996) 15
<i>Remarks</i>	The effect of the acidification of the samples was examined because when using off-line cartridges packed with 200 mg of the SDB-1 sorbent, arecovery of 100% was measured when 500 ml of water at pH 7 was percolated. The advantage was that humic materials were not co-extracted at pH 7 and the amount of interfering substances in surface water was low.

13	SPE
<i>Material</i>	1. LiChrolut EN (Merck) 2. GCB (Alltech)
<i>Amount material (mg) and size</i>	1. 200 mg 2. 300 mg
<i>Deactivation</i>	10 ml acetone, 3 ml MeOH and 3 ml water (pH 2)
<i>Kind of sample (ml)</i>	50 ml Milli-Q water, 50 ml effluent and 50-400 ml waste water (pH 2)
<i>Breakthrough volumes (ml)</i>	Breakthrough volumes > 400 ml for all phenols using LiChrolut EN
<i>Amount eluent (ml)</i>	3x2.5 ml of ACN:MeOH (50:50)
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	The combined eluate was basified with 150 ul of a 0.5 mol/l methanolic solution of TMAOH to minimise evaporation losses of volatile analytes and was concentrated to 1 ml under gentle stream of N ₂ .
<i>Enrichment factor</i>	
<i>Compounds</i>	phenol, p-cresol, 2-chlorophenol, 2,4-dichlorophenol, 2,4,6-trichlorophenol, 4-nitrophenol, 2,4-dinitrophenol, pentachlorophenol.
<i>Recovery</i>	The best recoveries were obtained with LiChrolut EN and ranging between 80 and 110%
<i>Repeatability</i>	< 12 % (RSD)
<i>LOD's</i>	2-10 ppb (HPLC-UV), 0.8-60 ppb (GC-MS, EI mode) and 0.02-0.6 ppb (GC-MS, NICI mode)
<i>Reference</i>	J. Cheung, R.J. Wells/ J. Chromatogr. A 771 (1997) 203
<i>Remarks</i>	LiChrolut EN cartridges can be re-used to preconcentrate water samples without any problems. The recoveries of phenols from effluent samples remained relatively unchanged after three repeated extraction.

14	SPE
<i>Material</i>	1. SDB-1 (Baker) 2. C18 Empore 3M disks (Baker) 3. SDB Empore 3M disks (Baker) 4. PLRP-S (Polymer Labs) 5. PRP-1 (Hamilton)
<i>Amount material (mg) and size</i>	1. 3 ml, 200 mg 2. 450 mg, 11 um 3. 450 mg, 6.8 um 4. 10x2 mm ID, 15-25 um), 20-80 mg 5. 100x4.6 mm ID, 10 um), 20-80 mg
<i>Deactivation</i>	5 ml MeOH, 10 ml LC-grade water
<i>Kind of sample (ml)</i>	500-1000 ml LC-grade water, tap water and river water
<i>Breakthrough volumes (ml)</i>	> 500 ml
<i>Amount eluent (ml)</i>	4 ml MeOH
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	Evaporation under a stream of nitrogen: After the desorption with MeOH, 50 ul of a mixture containing MeOH and ammonia (4:1) were added before evaporation step. Under these basic conditions, acidic compounds are ionised and therefore cannot be volatilised during the evaporation step.
<i>Enrichment factor</i>	> 200
<i>Compounds</i>	Log Kow range: -0.5 to 4 oxamyl, DIA (metabolite of atrazine), DEA (metabolite of atrazine), chloridazon, carbendazim, aldicarb, methoxuron, simazine, 2-chlorophenol, dicamba, cyanazine, bentazone, atrazine, carbaryl, isoproturon, ioxynil, MCPP, difenoxuron, 2,4-DB, 2,4,5-T, metolachlor, dinoterb.
<i>Recovery</i>	The best recoveries were obtained with SDB-1 cartridges, ranging between 80-105% (pH 7)
<i>Repeatability</i>	?
<i>LOD's</i>	0.1 ug/l (river water)
<i>Reference</i>	V. Pichon et al/ J Chromatogr. A 737 (1996) 25
<i>Remarks</i>	At pH 7, humic and fulvic acids are not co-extracted using SDB-1. At pH 3, these two acids are co-extracted and interfere the more polar analytes.

15	SPE
Material	Lichrolut EN (Merck)
<i>Amount material (mg) and size</i>	250 mg, 6 ml
<i>Deactivation</i>	10 ml ACN, 10 ml MeOH and 10 ml distilled, dionised water (pH <2)
<i>Kind of sample (ml)</i>	1000 ml drinking water and other water samples
<i>Breakthrough volumes (ml)</i>	> 1000 ml
<i>Amount eluent (ml)</i>	2x3 ml ACN
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	The eluate was reduced to 1 ml under gentle stream of N ₂
<i>Enrichment factor</i>	> 500
<i>Compounds</i>	phenol, m-cresol, o-cresol, 2,4-dimethylphenol, 2-chlorophenol, 4-chloro-3-methylphenol, 2,4-dichlorophenol, 2,4,6-trichlorophenol, 2-nitrophenol, 2,3,4,6- tetrachlorophenol, 4-nitrophenol, pentachlorophenol, 2,4- dinitrophenol, 4,6-dinitro-ocresol
<i>Recovery</i>	between 70 and 127%
<i>Repeatability</i>	< 20%
<i>LOD's</i>	LOQ's (0.015-0.075 ug/l and LOD's (0.005-0.025 ug/l)
<i>Reference</i>	T. Heberer, H.-J. Stan/ <i>Analitica Chimica Acta</i> 341 (1997) 21
<i>Remarks</i>	Applying SPE with LiChrolut EN adsorbent, phenols can be detected at the ng/l level in environmental sam- ples independent of the origine of the sample and its matrix load.

16	SPE
Material	1.Sep-Pak silica (Waters) 2. PGC (Shandon)
<i>Amount material (mg) and size</i>	1. 650 mg 2. 190 mg
<i>Deactivation</i>	1. - 2. 2 ml MeOH, 2 ml H ₂ O, 2ml MeOH and 2 ml H ₂ O
<i>Kind of sample (ml)</i>	5-250 ml distilled water, tap water and river water
<i>Breakthrough volumes (ml)</i>	1. > 250 ml for the present compounds 2. 50-250 ml depending on the components
<i>Amount eluent (ml)</i>	1. 2 ml of 8% MeOH in 6.0 M HCl 2. 2 ml TFA:ACN (20:80)
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	The eluates were evaporated at 50C almost to dryness. The residus were dissolved with 1 ml H ₂ O
<i>Enrichment factor</i>	> 200
<i>Compounds</i>	paraquat (PQ), EQ, diquat (DQ) and difenzoquat (DF)
<i>Recovery</i>	between 65 and 98%
<i>Repeatability</i>	4-10% RSD's
<i>LOD's</i>	LOD's (0.2-4 ug/l)
<i>Reference</i>	M.C. Carneiro et al./ <i>Analitica Chimica Acta</i> 408 (2000) 263
<i>Remarks</i>	

17	SPE
<i>Material</i>	1. Carbograph 1 2. Carbograph 4
<i>Amount material (mg) and size</i>	6.5x1.3 cm ID cartridges 1. 500 mg, 120x400 mesh size 2. 500 mg
<i>Deactivation</i>	10 ml DCM:MeOH (80:20) containing TBACl, 2 ml MeOH, 14 ml HCl-acidified water (pH 2)
<i>Kind of sample (ml)</i>	4000 ml drinking water and 1000 ml river water: before spiking, 0.5 g/l Na ₂ SO ₃ .5H ₂ O was added to the sample to avoid oxidation of the analytes
<i>Breakthrough volumes (ml)</i>	> 4000 ml for tap water and > 1000 ml for river water
<i>Amount eluent (ml)</i>	6 ml DCM:MeOH (80:20) containing TBACl
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	The extract was divided in two equal portions which were dried in a water bath at 27C under a gentle stream of N ₂ . The first residue was reconstituted with 200 ul of H ₂ O/MeOH (80:20) acidified with 0.6% (v/v) TFA. The second extract with 150 ul of H ₂ O/ACN (80:20) basified with Na ₂ CO ₃ , 0.1 mol/l.
<i>Enrichment factor</i>	> 30 000
<i>Compounds</i>	phenol, 2-nitrophenol, 2-chlorophenol, 2,4-dimethylphenol, 4- chloro-3-methylphenol, 2,6-dinitro-2-methylphenol, 2,4-dichlorophenol, 2,4,6-trichlorophenol, 4-nitrophenol, 2,4- dinitrophenol, pentachlorophenol.
<i>Recovery</i>	between 91 and 101 % with Carbograph 4
<i>Repeatability</i>	1-7% RSD's
<i>LOD's</i>	< 1 ng/l
<i>Reference</i>	A. Di Corcia et al. /J. Chromatogr. A 733 (1996) 383
<i>Remarks</i>	If TBACl was not added to the solvent mixture, seven evaporative losses of most of the phenols were observed on drying the extract. Apparently, the presence of TBACl in the DCM:MeOH mixture makes the analytes less vola- tile.

18	SPE
Material	Lichrolut EN (Merck)
<i>Amount material (mg) and size</i>	200 mg, 3 ml
<i>Deactivation</i>	?
<i>Kind of sample (ml)</i>	50 ml waste water
<i>Breakthrough volumes (ml)</i>	Not observed by 50 ml water sample
<i>Amount eluent (ml)</i>	5 ml MeOH for HPLC and/or 5 ml Ethyl acetate for GC-MS
<i>Fraction (nr.)</i>	One fraction per cartridge
<i>Drying method</i>	Evaporation by nitrogen stripping to 1 ml
<i>Enrichment factor</i>	50
<i>Compounds</i>	Log Kow: -06 - 5.16 simasine; 1,2,3,4-tetrachlorobenzene; 2-hydroxydesethylatrazine; terbutylazine; 2-phenoxyethanol; triphenyl phosphate; 2- butoxyethanol acetate; 2-methylthiobenzothiazole; dibenzyl phthalate; dibutyl phthalate; dimethyl phthalate; bis(2-ethylhexyl)phthalate; benzylbutyl phthalate; 3,5- dichloroaniline; 4-nonylphenol
<i>Recovery</i>	> 85
<i>Repeatability</i>	Not reported
<i>LOD's</i>	?
<i>Reference</i>	S. Galassi et al./ J. Chromatogr. A 889 (2000) 149
<i>Remarks</i>	

19	SPE
Material	Strong anion exchange resin + RP-102 in series
<i>Amount material (mg) and size</i>	500 mg, 6 ml + 1 g, 6 ml
<i>Deactivation</i>	60 ml of 1% acetic acid in acetone+ 40 ml of 1% acetic acid in water
<i>Kind of sample (ml)</i>	1 l of sample waters (pH 2.8-3)
<i>Breakthrough volumes (ml)</i>	For more polar analytes (1 l) (see recoveries)
<i>Amount eluent (ml)</i>	3x4 ml of 1% acetic acid in acetone (3 ml/min)
<i>Fraction (nr.)</i>	One fraction per two cartridges in series
<i>Drying method</i>	2 times evaporation to dryness by Zymark solvent reduction workstation and reconstituted in 0.1 ml acetoni- tril
<i>Enrichment factor</i>	10 000
<i>Compounds</i>	imazapyr; imazythapyr; flumetsulam; nicosulfuron; imazaquin; thifensulfuron methyl; metsulfuron methyl; chlorsulfuron; sulfometuron methyl; triasulfuron methyl; bensulfuron methyl; halosulfuron methyl; prosulfuron; chlorimuron ethyl; triflusulfuron methyl; primisulfuron methyl
<i>Recovery</i>	39-92
<i>Repeatability</i>	14-26% (SD's)
<i>LOD's</i>	< 10 ng/l level
<i>Reference</i>	E.T Furlong et al./ The Science of the Total environment 248 (2000) 135
<i>Remarks</i>	

20	SSPE
<i>Material</i>	Layered column comprising of C2 (5 g) on Isolut ENV+ (2 g) layer
<i>Amount material (mg) and size</i>	C2 (5 g) + Isolut ENV+ (2 g)
<i>Deactivation</i>	?
<i>Kind of sample (ml)</i>	20 l surface water
<i>Breakthrough volumes (ml)</i>	Not observed
<i>Amount eluent (ml)</i>	sequential elution with 5 ml MeOH and 5 ml DCM
<i>Fraction (nr.)</i>	One fraction
<i>Drying method</i>	Evaporation under nitrogen (40C, for 2 h) to 1 ml MeOH
<i>Enrichment factor</i>	20 000
<i>Compounds</i>	Log Kow: 0.43 - 5.12 diethylnitrosamine; dipropylnitrosamine; desmetryn; diazinon; pentachlorophenol; pirimephosmethyl; phenanthrene; carbophenothion; nonylphenol
<i>Recovery</i>	?
<i>Repeatability</i>	?
<i>LOD's</i>	?
<i>Reference</i>	K.V. Thomas et al./ Marine Pollution Bulletin 38 (1999) 925
<i>Remarks</i>	

21	GPC + SPE
<i>Material</i>	TSK gel HW 40s + Envi Chrom P + LiChrolut EN
<i>Amount material (mg) and size</i>	120 x 5 cm + 250 mg + 200 mg
<i>Deactivation</i>	MeOH + phosphate-buffer
<i>Kind of sample (ml)</i>	200 ml water samples
<i>Breakthrough volumes (ml)</i>	?
<i>Amount eluent (ml)</i>	Sequential elution with 3 ml MeOH/H ₂ O (1:1) and 3 ml MeOH
<i>Fraction (nr.)</i>	GPC (5 fractions) and SPE one fraction per two cartridges
<i>Drying method</i>	Evaporation by Speed-Vac concentrator to 0.1 ml
<i>Enrichment factor</i>	2000
<i>Compounds</i>	-
<i>Recovery</i>	-
<i>Repeatability</i>	-
<i>LOD's</i>	-
<i>Reference</i>	N. Klinkow et al./ Wat. Res. 32 (1998) 2583
<i>Remarks</i>	This goal of this work was to develop a separation/fractionation scheme of organic compounds according to their molecular size (GPC) and their polarity (SPE). The method was combined with luminescence inhibition test to measure the toxicity of the fractions.

22	SPE
<i>Material</i>	Bakerbond Speedisk divinylbenzene (DVB) disk C18-silica disks (Bakerbond Speedisk C18)
<i>Amount material (mg) and size</i>	300 mg DVB polymers (50-mm diameter and 0.5-mm high bed) +
<i>Deactivation</i>	750 mg of C18-silica (50-mm diameter and 1-mm high bed)
<i>Kind of sample (ml)</i>	10 ml acetonitrile, 10 ml MeOH and 10 ml of LC-grade water 250–1000 ml of surface water and drinking water samples (200 ml/min flow rate). For C18 addition of 10% of MeOH to the sample before preconcentration
<i>Breakthrough volumes (ml)</i>	Not observed
<i>Amount eluent (ml)</i>	DVB disks: 9 ml acetonitrile C18 disks: 10 ml DCM:MeOH (4:1, v/v)
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	After addition of 120- μ l volume of MeOH:ammonia (4:1, v/v), evaporation was commenced at 40°C using a rotary evaporator to 1 ml and then further evaporated under a gentle stream of nitrogen to a final volume of about 100 μ l.
<i>Enrichment factor</i>	10 000
<i>Compounds</i>	Alachlor; Aldicarb; Aminotriazol; Atrazine; Chlorpyrifos; Dinoterb; Diuron; Endosulfan á; Endosulfan â; Fenpropimorph; Fluzilazole; Ioxynil; Isoproturon; Lindane; Linuron; Oxydemeton-methyl; Simasine; Terbutylazine; Trifluraline; Triallate.
<i>Recovery</i>	77–101%
<i>Repeatability</i>	3–12%
<i>LOD's</i>	0.01–0.05 μ l/l in dinking or ground water, whereas these detection limits are still in the low 0.1 μ l/l levels in 1000 ml of dirty surface water samples.
<i>Reference</i>	V. Pichon et al./J. Chromatogr. A 795 (1998) 83
<i>Remarks</i>	Using the two SPE procedures, it is possible to extract any compound, either polar or non-polar (Log Kow 1–6 range). Addition of 10% of MeOH to the sample was necessary to make extraction and desorption of apolar compound efficiently. Only very polar analytes may still be difficult to extract with these sorbents. In this case, a graphitized carbon is recommended.

23	SPE
<i>Material</i>	Sep-Pak Plus C18; tC18; C8; tC2; CN; Diol; NH2; and PS-1 (all from waters)
<i>Amount material (mg) and size</i>	?
<i>Deactivation</i>	5 ml MeOH and 20 ml of distilled water
<i>Kind of sample (ml)</i>	Distilled water, tap water, ground water and river water (50–1000 ml volumes)
<i>Breakthrough volumes (ml)</i>	Only for Asulam (Log Kow –0.5), breakthrough volume of 500 ml. For the rest there is no problem even for 1000 ml if the good condition are used.
<i>Amount eluent (ml)</i>	1–5 ml of 50% v/v MeOH–water or 1–5 ml of MeOH depending on the polarity of the analytes.
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	-
<i>Enrichment factor</i>	Up to 1000
<i>Compounds</i>	Log Kow –05–5.38 range Butan-2-one; Acetanilidine; Benzonitrile; Acetophenone; Nitrobenzene; Benzene; Methyl benzoate; Fluorobenzene; Ethyl benzoate; Toluene; Chlorobenzene; Brobobenzene; Ethylbenzene; Iodobenzene; p-Xylene; o-Xylene; 1,4-Dichlorobenzene; Phenyl benzoate; Naphthalene; Propylbenzene; Benzyl benzoate; 1,2,4-Trichlorobenzene; Phenanthrene; Butylbenzene; Fluoranthene; 1,2,3,5-Tetrachlorobenzene; Asulam; Metolcarb; Propoxur; Carbaryl; Xylylcarb; Macbal; Cyanophos; Dichlobinil; Phosmet; Fenobcarb; p-Dichlorobenzene; o-Dichlorobenzene; Fenitrothion; Mecoprop; Chlorothalonil; Fenthion; Phenthoate; Isixathion; Quintozene; Tolclophos-methyl; Dichlofenthion.
<i>Recovery</i>	Depend on the polarity of the analytes, sorbent, sample volume and elution solvent. For polar analytes the best recoveries were obtained with PS-1 (85-100%). For non-polar ones, C18 gives the best results (85-100%).
<i>Repeatability</i>	-
<i>LOD's</i>	-
<i>Reference</i>	M. Nakamura et al./ The Analyst, 121 (1996) 469
<i>Remarks</i>	Recoveries of the studied compounds were related to their Log Kow. Apolar compounds were easily sorbed, but were difficult to elute. The extent of elution could be improved by adding methanol to the sample solution.

24	SPE
<i>Material</i>	Bakerbond C18 or Polar Plus C18 (Baker or Supelco) Apolar PS-DVB; Isolute ENV+ (IST) Empore disks: C18 and PS-DVB (Baker)
<i>Amount material (mg) and size</i>	500 mg C18 and 200 mg PS-DVB or Isolute ENV+
<i>Deactivation</i>	Cartridges: 5 ml MeOH + 10 ml LC-grade water Disks: 10 ml MeOH:ACN (1:1) + 10 ml MeOH + 20 ml LC-grade water
<i>Kind of sample (ml)</i>	500–1000 ml water samples (10 ml/min percolation flow rate for cartridges and 25 ml/min for disks).
<i>Breakthrough volumes (ml)</i>	Depend on Log Kw values an sorbent
<i>Amount eluent (ml)</i>	Cartridges: C18 (4 ml MeOH) and PS-DVB (2 ml MeOH:ACN (1:1)). Disks: C18 (2x4 ml MeOH:ACN (1:1)) and PS-DVB (2x6 ml MeOH:ACN (1:1))
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	Under gentle stream of nitrogen
<i>Enrichment factor</i>	> 1000
<i>Compounds</i>	Log Kow -0.5– 3 Oxamyl; Methomyl; Aldicarb; Carbendazim; Carbofuran; Aminocarb; Carbaryl Captan; Methiocarb; Fenuron; Metoxuron; Monuron; Monolinuron; Chlortoloron; Isoproturon; Diuron; Desethylatrazine; Hydroxyatrazine; Simazine; Simetryne; Prometon; Propazine; Terbutylazine; Metamitron; Metribuzine; Hexazinone.
<i>Recovery</i>	The best recoveries were obtained with Isolute ENV+ and PGC (carbon). The recoveries ranging from 85% to 100%.
<i>Repeatability</i>	-
<i>LOD's</i>	-
<i>Reference</i>	M.-C. Hennion et al./J. Chromatogr. A 823 (1998) 147
<i>Remarks</i>	The selection of the sorbents for the extraction of polar analytes is based on Log Kw values.

25	SPE
<i>Material</i>	Carbopack B 120/400 (Supelco) ENVI Chrom P (Supelco)
<i>Amount material (mg) and size</i>	500 mg of Carbopack B and 500 mg ENVI Chrom P (80–160 µm particle size)
<i>Deactivation</i>	10 ml of MeOH (ENVI Chrom P) or DCM (Carbopack B) + 10 ml water and 2 ml of a 5 mM tetrabutylammonium bromide (TBA) solution.
<i>Kind of sample (ml)</i>	500 ml of Milli-Q-purified water, tap water samples (pH 9.0)
<i>Breakthrough volumes (ml)</i>	> 1000 ml except phenol.
<i>Amount eluent (ml)</i>	Carbopack B: 5 ml DCM + 1% of acetic acid ENVI Chrom P: 5 ml MeOH + 1% of acetic acid
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	Under vacuum with rotary evaporator to 1 ml.
<i>Enrichment factor</i>	500
<i>Compounds</i>	Phenol; 4-Nitrophenol; 2,4-Dinitrophenol; 2-Chlorophenol; 2-Nitrophenol; 2,6-Dimethylphenol; 2,4-Dimethylphenol; 2-Methyl-4,6-dinitrophenol; 4-Chloro-3-methylphenol; 2,4-Dichlorophenol, 2,4,6-trimethylphenol; 2,4,6-Trichlorophenol; pentachlorophenol.
<i>Recovery</i>	All higher than 90% except for phenol (54%, breakthrough volume of 1000 ml).
<i>Repeatability</i>	< 10%
<i>LOD's</i>	Direct injection: 30–100 µg/l Carbon cartridge: 2–6 µg/l Polymere cartridge: 60–100 ng/l
<i>Reference</i>	E. Pocerull et al./J. Chromatogr. A 719 (1996) 105
<i>Remarks</i>	To increase the retention of the most polar compounds, mainly phenol, TBA was used as an ion-pair reagent in both sorbents, an increase in the breakthrough volumes, especially that correspond to phenol, which was the most polar compound studied.

26	SPE
<i>Material</i>	Isolute ENV+ PLRP-S LiChrolut PGC (carbon)
<i>Amount material (mg) and size</i>	10 x 0.2 mm I.D. stainless-steel precolumns
<i>Deactivation</i>	5 ml MeOH + 1 ml water (pH 3) at 1 ml/min
<i>Kind of sample (ml)</i>	10–200 ml water samples (pH 2.5) at 4 ml/min
<i>Breakthrough volumes (ml)</i>	Isolut ENV+ and Lichrolut have the highest breakthrough volumes. By PGC there is desorption problem, so the majority of the analytes cannot be analyzed. PGC gave good results only for aminophenols.
<i>Amount eluent (ml)</i>	Backflush mode desorption
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	- (on-line mode)
<i>Enrichment factor</i>	> 500
<i>Compounds</i>	Catechol; Phenol; 4-Nitrophenol; 2,4-Dinitrophenol; 2-Chlorophenol; 3-Chlorophenol; 4-Chlorophenol; 2-Nitrophenol; 2,6-Dimethylphenol; 2,4-Dimethylphenol; 2-Methyl-4,6-dinitrophenol; 4-Chloro-3-methylphenol; 2,4-Dichlorophenol, 2,4,6-trimethylphenol; 2,3,5-Trichlorophenol; 2,3,4-Trichlorophenol; 3,4,5-Trichlorophenol; pentachlorophenol; 2-Amino-4-chlorophenol; 4-Chloro-2-aminophenol.
<i>Recovery</i>	Good recoveries were obtained for all sorbents for 100 ml sample, except for phenol, catechol and 4-chloro-2-aminophenol (V = 50 ml). (55–105%).
<i>Repeatability</i>	< 12%
<i>LOD's</i>	Ground water: 0.01–7 µg/l River water: 0.02–14 µg/l
<i>Reference</i>	D. Puig et al./ J. Chromatogr. A 733 (1996) 371
<i>Remarks</i>	A matrix effect study showed that acidification of the sample is necessary to avoid binding of some analytes to the humic substances and to prevent their partial deprotonation.

27	SPE
<i>Material</i>	Carbograph 1 (GCB) Carbograph 4 (GCB) Carbograph 5 (GCB) LiChrolut (PS-DVB) (Merck) Envi-Chrom P (PS-DVB) (Supelco)
<i>Amount material (mg) and size</i>	GCB's: 37–150 µm particle size. LiChrolut: 1200 m ² /g surface area and 40–120 µm particle size range. Envi-Chrom P: 900 m ² /g surface area and 80–160 µm particle size range.
<i>Deactivation</i>	-
<i>Kind of sample (ml)</i>	4000 ml Tap water and 500–1000 ml river water samples
<i>Breakthrough volumes (ml)</i>	No breakthrough observed.
<i>Amount eluent (ml)</i>	GCB's: 1 ml MeOH + 6 ml of DCM:MeOH (4:1) + 6 ml of DCM:MeOH (4:1) containing 10 mmol/l tetrabutylammonium chloride (TBACl) or trifluoroacetic acid (TFA). LiChrolut and Envichrom P: 2x5 ml MeOH:ACN (1:1)
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	Extracts were dried in a water bath at 30°C by allowing a gentle stream of nitrogen. Before concentrating extracts, addition of KOH was necessary to avoid evaporative losses of the most volatile phenols.
<i>Enrichment factor</i>	2000–8000
<i>Compounds</i>	Omethoate; Aldicarb sulfone; Butocarboxim sulfoxide; Aldicarb sulfoxide; Butoxycarboxim; Oxamyl; Methomyl; Monocrotophos; Atrazine, desethyl; Metamitron; Metribuzin; Chloridazon; 2,4-Dinitroph; 4,6-Dinitro-2-methylphenol; Pentachlorophenol; Linuron; Aldicarb; Dichlorprop; 2,4,5-Thrchlorophenoxy acetic acid; Ioxynil; 2,4-Dichlorophenoxy acetic acid; 2,4-Dichlorophenoxy butyric acid; Mecoprop; 4-Nitrophenol; 2-Chlorophenol; 2-Nitrophenol; Pichloram; Chloramben; Dicamba; Bentazone; Phenol; Sulfohenyl-4-propionic acid; Sulphonyl-4-butyric acid; Sulfohenyl-4-valeric acid; Nonylphenoxyacetic acid; Dimethoate; Monuron; Simasine; Propoxur; Isoproturon.
<i>Recovery</i>	Carbograph 4 gave the best results with less problems compared to other sorbent and especially Carbograph 5. Recoveries ranging from 85% to 100% including very polar analytes.
<i>Repeatability</i>	-
<i>LOD's</i>	Very low ng/l level
<i>Reference</i>	C. Crescenzi et al./ J. Chromatogr. A 733 (1996) 41
<i>Remarks</i>	When analyzing analyte mixtures both acidic and neutral in nature, a differential elution scheme was adopted. Base/neutral compounds were firstly eluted by 1 ml of methanol followed by 6 ml of CH ₂ Cl ₂ :CH ₃ OH (80:20,v/v). Acidic compounds were successively eluted and collected separately by passing 6 ml of the solvent mixture having the same composition as that mentioned above, but containing 10 mmol/l TBACl.

28	SPE
<i>Material</i>	Carbograph 1 (Alltech) LiChrolut-EN (Merck) PS-DVB empore disks
<i>Amount material (mg) and size</i>	Carbograph 1: 37–150 µm particle size LiChrolut-EN: 1200 m ² /g surface area and 40–120 µm particle size PS-DVB empore disks: 47 mm diameter and 0.5 mm thick and 500 mg sorbent in it.
<i>Deactivation</i>	Carbograph 1: 7 ml of DCM:MeOH (4:1) acidified with 50 mmol/l formic acid + 5 ml MeOH + 20 ml of LC-grade water acidified with 10 mmol/l HCl. LiChrolut-EN: 5 ml of MeOH:ACN (1:1) + 5 ml of LC-grade water
<i>Kind of sample (ml)</i>	1000–2000 ml water samples.
<i>Breakthrough volumes (ml)</i>	Occur when polymer sorbent was used by the above mentioned sample volumes.
<i>Amount eluent (ml)</i>	Carbograph 1: 8 ml DCM:MeOH (4:1) acidified with 50 mmol/l formic acid. Lichrolut EN: 2x3 ml of MeOH:ACN (1:1). PS-DVB Empore disks: 3x3 ml of MeOH:ACN (1:1).
<i>Fraction (nr.)</i>	1
<i>Drying method</i>	Under gentle nitrogen stream by 40°C after adding 0.4 mol/l of tetrabutylammonium fluoride (TBAF) in methanol.
<i>Enrichment factor</i>	2000–4000
<i>Compounds</i>	Fluazifop; Clodinafop; Quizalofop; Fenoxaprop; Haloxyfop; Diclofop.
<i>Recovery</i>	Carbograph 1 was superior: Recoveries 90–98% range
<i>Repeatability</i>	-
<i>LOD's</i>	7-20 ng/l for drinking water and 16–36 ng/l for spring water.
<i>Reference</i>	A. Lagana et al./ J Chromatogr. A 796 (1998) 309
<i>Remarks</i>	A good extraction efficiency for the six herbicides from acidified water was obtained with both LiChrolut-EN cartridge and PS-DVB Empore disk.

Annex 2 List of responsive compounds

Seq. nr	Compound name	S9	RIZA	CAS.Reg.No	Smiles Epiwin
1	1,1-Dichloroacetone		yes	513-88-2	<chem>ClC(Cl)C(=O)C</chem>
2	1,2:3,4-Dibenzoanthracene			215-58-7	<chem>c5ccc4cc3c1cccc1c2cccc2c3cc4c5</chem>
3	1,2:3,4-Diepoxybutane			298-18-0	<chem>C1OC1C2CO2</chem>
4	1,2-Diaminobenzene(+ S9)	yes		95-54-5	<chem>Nc(c(N)ccc1)c1</chem>
5	1,2-Dimethylhydrazine			540-73-8	<chem>CNNC</chem>
6	1,2-Epoxybutan		yes	106-88-7	<chem>O(C1CC)C1</chem>
7	1,3-Dinitropyrene			75321-20-9	<chem>O=N(=O)c4c2c3c(c(c4)N(=O)=O)ccc1c3c(ccc1)cc2</chem>
8	1,3-Dioxane			505-22-6	<chem>O1COCCC1</chem>
9	1,5-Dinitronaphthalene			605-71-0	<chem>c1cc2c(N(=O)=O)cccc2c(N(=O)=O)c1</chem>
10	1,6-Dinitropyrene			42397-64-8	<chem>O=N(=O)c1ccc2ccc3c(ccc4ccc1c2c34)N(=O)=O</chem>
11	1,8-Dinitropyrene			42397-65-9	<chem>O=N(=O)c1ccc2ccc3ccc(N(=O)=O)c4ccc1c2c34</chem>
12	1-Bromo pentane			110-53-2	<chem>BrCCCCC</chem>
13	1-Chloro-2,4-dinitrobenzene		yes	97-00-7	<chem>O=N(=O)c(ccc(c1N(=O)=O)Cl)c1</chem>
14	1-Naphthylamine		yes	134-32-7	<chem>c(c(c(N)cc1)ccc2)(c2)c1</chem>
15	1-Nitronaphthalene		yes	86-57-7	<chem>O=N(=O)c(c(c(ccc1)cc2)c1)c2</chem>
16	1-Nitropyrene			5522-43-0	<chem>O=N(=O)c(c(c(c(c1)ccc2)c2cc3)c3c4)c1)c4</chem>
17	2,4,7-Trinitro-9-fluorene			34263-36-0	<chem>O=[N+](O)C3=CC1=C(C=C3)C2=C(C=C([N+](O)O)C=C2[N+](O)O)C1</chem>
18	2,4-Diaminotoluene (+ S9)	yes		95-80-7	<chem>Nc(c(ccc1N)C)c1</chem>
19	2,4-Dimethyl-1,3-dioxane			766-20-1	<chem>O(CCC1C)C(O1)C</chem>
20	2,4-Dinitrotoluene		yes	121-14-2	<chem>O=N(=O)c(ccc(c1N(=O)=O)C)c1</chem>
21	2,6-Diaminotoluene (+ S9)	yes		823-40-5	<chem>Nc(c(c(N)cc1)C)c1</chem>
22	2,7-Dinitro-8-fluorenone			31551-45-8	<chem>O=C(c(c(c1ccc(N(=O)=O)c2)ccc3N(=O)=O)c3)c12</chem>
23	2-Acetylaminofluorene (+S9)	yes	yes	53-96-3	<chem>O=C(Nc(ccc(c1Cc2cccc3)c23)c1)C</chem>
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MelQ (+ S9)	yes		77094-11-2	<chem>CC2=CC1=NC=CC=C1C3=C2N(C)C(N)=N3</chem>

25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	yes		76180-96-6	<chem>Cc3cc1c(ccc2ncnc12)nc3N</chem>
26	2-Amino-3-methyl-9H-pyrido(2,3-b)indole; MeA alpha C (+ S9)	yes		68006-83-7	<chem>Cc1cc2c(nc1N)nc3ccccc23</chem>
27	2-Amino-6-methyl-2-pyridopyrido(1,2-a:3',2'd)imidazole (Glu-P-1; + S9)	yes		67730-11-4	<chem>CC1=CC=Cn2c1nc3ccc(N)nc23</chem>
28	2-Amino-9H-pyrido(2,3-b)indole (+S9)(A alpha C)	yes		26148-68-5	<chem>Nc1ccc2c(nc3c2ccccc3)n1</chem>
29	2-Aminoanthracene(+ S9)	yes		613-13-8	<chem>Nc3ccc2cc1ccccc1cc2c3</chem>
30	2-Amino-dipyrido(1,2-a:3',2'-d)imidazole-2-amine (Glu-P-2; + S9)	yes		67730-1 0-3	<chem>Nc3ccc2nc1C=CC=Cn1c2n3</chem>
31	2-Aminofluorene(+S9)	yes	yes	153-78-6	<chem>Nc3ccc1c(Cc2ccccc12)c3</chem>
32	2-Amino-1-methyl-6-phenylimidazo[4,5-b]pyridine (PhIP)			105650-23-5	<chem>CN1C(N)=NC2=C1C=C(C3=CC=CC=C3)C=N2</chem>
33	2-Nitrofluorene			607-57-8	<chem>O=N(=O)c(ccc(c1Cc2ccccc12)c23)c1</chem>
34	2-Nitronaphthalene			581-89-5	<chem>O=N(=O)c2ccc1ccccc1c2</chem>
35	2-Nitro-p-phenylenediamine			5307-14-2	<chem>O=N(=O)c(c(N)ccc1N)c1</chem>
36	3,7-Dinitrofluoranthene			105735-71-5	<chem>O=N(=O)c1ccc2c4ccccc(c4c3c2c1ccc3)N(=O)=O</chem>
37	3,9-Dinitrofluoranthene			22506-53-2	<chem>[O-][N+](C2=CC=C3C1=C(C4=C3C=C([N+])([O-])=O)C=C4)C=CC=C12)=O</chem>
38	3-Amino-1,4-dimethyl-5H-pyrido[3,4-b]indole; Trp-P-1 (+ S9)	yes		62450-06-0	<chem>Cc2nc(N)c(C)c3c1ccccc1nc23</chem>
39	3-Amino-1-methyl-5H-pyrido[4,3-b]indole; Trp-P 2 (+ S9)	yes		62450-07-1	<chem>CC1=NC(N)=CC2=C1C(C=CC=C3)=C3N2</chem>
40	3-Methoxy-4-aminoazobenzene(+S9)	yes		3544-23-8	<chem>COc1c(N)ccc(N=Nc2ccccc2)c1</chem>
41	3-Methylcholanthrene (+ S9)	yes	yes	56-49-5	<chem>c(c(ccc1C)cc(c2ccc3ccccc4c34)(c1CC5)c25</chem>
42	3-Nitrofluoranthene			892-21-7	<chem>O=N(=O)c2ccc3c1ccccc1c4ccccc2c34</chem>
43	4,4'-Dinitrobiphenyl			1528-74-1	<chem>O=N(=O)c(ccc(c(ccc(N(=O)=O)c1)c1)c2)c2</chem>
44	4-Aminobiphenyl		yes	92-67-1	<chem>Nc(ccc(c(ccc1)c1)c2)c2</chem>
45	4-Methyl-1,3-dioxane			1120-97-4	<chem>O1COC(C)CC1</chem>
46	4-Nitro-o-phenylenediamine			99-56-9	<chem>O=N(=O)c(ccc(N)c1N)c1</chem>
47	4-Nitroquinoline-N-oxide		yes	56-57-5	<chem>O=N(=O)c1ccn(=O)c2ccccc12</chem>

48	4-Nitrosodiphenylamine			156-10-5	<chem>O=Nc(ccc(Nc(cccc1)c1)c2)c2</chem>
49	4-Nitroso-N,N-dimethylaniline			138-89-6	<chem>O=Nc(ccc(N(C)C)c1)c1</chem>
50	5-Butyrolactone			36536-46-6	<chem>O1C(=O)CC1C</chem>
51	5-Fluorouracil			51-21-8	<chem>N1C(=O)NC(=O)C(F)=C1</chem>
52	5-MOP (+ UV)			484-20-8	<chem>O(C=C1)c(cc(OC(=O)C=2)c3C2)c1c3OC</chem>
53	5-Nitro-2-furaldehyde			698-63-5	<chem>N(=O)(=O)c1ccc(C(=O))o1</chem>
54	5-Nitro-2-furylacrylacid			710-25-8	<chem>O1C(N(=O)=O)=CC=C1C=CC(=O)N</chem>
55	5-Nitroacenaphthene			602-87-9	<chem>O=N(=O)c(c(c(c(cc1)CC2)c2c3)c1)c3</chem>
56	6-Aminochrysene(+ S9)	yes		2642-98-0	<chem>Nc3cc2c1cccc1ccc2c4cccc34</chem>
57	6-chloro-9-(3-(2-chloroethylamino)propyl)amino-2-methoxyacridine*2HCl(ICR-191)			17070-45-0	<chem>COC3=CC=C2N=C1C=C(Cl)C=CC1=C(NCCCCNCCCl)C2=C3</chem>
58	6-Nitrochrysene			2-89-6	<chem>CC1CC2OC(=O)C(=C)C2C(OC(=O)CCCC(=O)OC3C4C(CC(C)C5C=CC(=O)C53C)OC(=O)C4=C)C6(C)C1C=CC6=O</chem>
59	7,12-Dimethylbanzanthracene (+ S9)	yes	yes	57-97-6	<chem>c(c(c(c(c1)ccc2)c2)c(c(c3ccc4)c4)C)(c3C)c1</chem>
60	8-MOP (+ UV)			298-81-7	<chem>O=C1C=Cc2cc3ccoc3c(OC)c2O1</chem>
61	8-Proplotactone			57-57-8	<chem>O=C(OC1)C1</chem>
62	9 Aminoacridine			90-45-9	<chem>Nc2c1cccc1nc3cccc23</chem>
63	Acridine Orange			494-38-2	<chem>CN(c1cc2c(cc1)cc3c(n2)cc(cc3)N(C)C)C</chem>
64	Acrinol(+S9)	yes		1837-57-6	<chem>NC1=C(C=C(OCC)C=C3)C3=NC2=C1C=CC(N)=C2</chem>
65	Adriamycin(Adr)			25316-40-9	<chem>COC1=C(C(C(C(O)=C([C@@H](O[C@H]5CC(N)C(O)[C@H](C)O5)C[C@@](C(CO)=O)(O)C4)C4=C3O)=C3C2=O)=O)C2=CC=C1</chem>
66	Afiatoxin G1			1165-39-5	<chem>C12C=COC1Oc3cc(OC)c(C(C5)=C(C(=O)OC5)C(=O)O4)c4c23</chem>
67	Aflatoxin B1			1162-65-8	<chem>C12C=COC1Oc3cc(OC)c(C(C5)=C(C(=O)C5)C(=O)O4)c4c23</chem>
68	Amphotericin B			1397-89-3	<chem>O=CCC(CC(CCC(C(CC(CC1(OC(CC(OC2C(C(N)C(C(O2)C)O)O)O)O)O)O)O)O)O)O)O</chem>
69	Amsacrine			51264-14-3	<chem>COc1cc(NS(C)(=O)=O)ccc1Nc3c2cccc2nc4cccc34</chem>
70	Auramine			2465-27-2	<chem>CN(C)c1ccc(cc1)C(=N)c2ccc(cc2)N(C)C</chem>
71	Avarol(+ S9)	yes		55303-98-5	<chem>OC(C(C[C@]([C@H](CC[C@@]23C)C)(C2CCC=C3C)C)=C1)=CC=C1O</chem>
72	Azaserine			115-02-6	<chem>NC(COC(=O)C=N#N)C(O)=O</chem>

73	Azinphos-Methyl		yes	86-50-0	<chem>S=P(OC)(OC)SCN1N=Nc2ccccc2C1(=O)</chem>
74	Azobenzene(+ S9)	yes	yes	103-33-3	<chem>N(=Nc(cccc1)c1)c(cccc2)c2</chem>
75	Benz(a)anthracene(+ S9)	yes	yes	56-55-3	<chem>c(c(c(c1)ccc2)c2)cc(c3ccc4)c4)(c1)c3</chem>
76	Benzidine (+ S9)	yes	yes	92-87-5	<chem>Nc(ccc(c(ccc(N)c1)c1)c2)c2</chem>
77	Benzo(a)pyrene(+ S9)		yes	50-32-8	<chem>c(c(c(cc1)ccc2)c2cc3)(c3cc(c4ccc5)c5)c14</chem>
78	Benzylchloride		yes	100-44-7	<chem>c(cccc1)(c1)CCl</chem>
79	Bleomycin			11056-06-7	<chem>O=C(N)C[C@@H](C1=NC(N)=C(C)C(C(N[C@@]([H])([C@@](C5=CNC=N5)([H])OC4C(OC6C(O)C(OC(N)=O)C(O)C(CO)O6)C(O)C(O)C(CO)O4)C(N[C@@]([C@H]([C@H](C)C(N[C@@]([C@](O)(C)[H])([H])C(NCCC2=NC(C3=NC(C(NCCC(SH)(C)C)=O)=CS3)=CS3)=O)=O)O)([H])C)=O)=O)=N1)NC[C@H](N)C(N)=O</chem>
80	Bromo acetic acid			79-08-3	<chem>O=C(O)CBr</chem>
81	Bromobenzene		yes	108-86-1	<chem>c(cccc1)(c1)Br</chem>
82	Captan		yes	133-06-2	<chem>O=C(N(SC(Cl)(Cl)Cl)C(=O)C1CC=CC2)C12</chem>
83	Chloramine T			127-65-1	<chem>ClN([Na])S(=O)(=O)c1ccc(C)cc1</chem>
84	Chromium(VI)oxide			1333-82-0	<chem>[Cr](=O)(=O)=O</chem>
85	Chrysene (+ S9)	yes	yes	218-01-9	<chem>c1ccc2ccc3c4ccccc4ccc3c2c1</chem>
86	Ciprofloxacin			85721-33-1	<chem>C1CNCCN1c2cc3N(C4CC4)C=C(C(=O)O)C(=O)c3cc2F</chem>
87	cis-Platinum			15663-27-1	<chem>Cl[Pt](Cl)(N(H)(H)H)N(H)(H)H</chem>
88	Citrinin			518-75-2	<chem>CC2OC=C1C(=C(C(O)=O)C(=O)C(=C1C2C)C)O</chem>
89	Cumenehydroperoxide			80-15-9	<chem>O(O)C(c(cccc1)c1)(C)C</chem>
90	Cupferron			135-20-6	<chem>NON(N=O)c1ccccc1</chem>
91	Danthron			117-10-2	<chem>O=C(c(c(c(O)cc1)C(=O)c2c(O)ccc3)c1)c23</chem>
92	Daunomycin			23541-50-6	<chem>COC1=C(C(C(C(O)=C([C@@H](O[C@H]5CC(N)C(O)[C@H](C)O5)C[C@@](C(C)=O)(O)C4)C4=C3O)=C3C2=O)=O)C2=CC=C1</chem>
93	Daunorubicin			20830-81-3	<chem>COc4ccccc5C(=O)c3c(O)c2CC(O)(CC(OC1CC(N)C(O)C(C)O1)c2c(O)c3C(=O)c45)C(C)=O</chem>
94	Dibromoacetoneitrile		yes	3252-43-5	<chem>N#CC(Br)Br</chem>
95	Dichlofluanid			1085-98-9	<chem>CN(C)S(=O)(=O)N(SC(F)(Cl)Cl)c1ccccc1</chem>
96	Diethylnitrosamine (+ S9)	yes	yes	55-18-5	<chem>O=NN(CC)CC</chem>
97	Diethylsulfate (+ S9)	yes	yes	64-67-5	<chem>O=S(=O)(OCC)OCC</chem>
98	Dimethylnitrosamine (+ S9)	yes	yes	62-75-9	<chem>O=NN(C)C</chem>

99	Dimethylsulfate		yes	77-78-1	<chem>O=S(=O)(OC)OC</chem>
100	Dimethylsulfoxide			67-68-5	<chem>O=S(C)C</chem>
101	Dithianone			3347-22-6	<chem>c1ccc2C(=O)C(S3)=C(SC(C#N)=C3C#N)C(=O)c2c1</chem>
102	Doxorubicinhydrochloride			25316-40-9	
103	Ellipticine			519-23-3	<chem>c1ccc2nc3c(C)c4ccncc4c(C)c3c2c1</chem>
104	Enoxacin			74011-58-8	<chem>C1CNCCN1c2nc3N(CC)C=C(C(=O)O)C(=O)c3cc2F</chem>
105	Epichlorhydrin		yes	106-89-8	<chem>O(C1CCl)C1</chem>
106	Epoxyystyrene			96-09-3	<chem>O(C1c(cccc2)c2)C1</chem>
107	Ethidium bromide (+ S9)	yes		1239-45-8	<chem>CCn3(Br)c(c1ccccc1)c2cc(N)ccc2c4ccc(N)cc34</chem>
108	Ethyl methanesulfonate (EMS)		yes	62-50-0	<chem>O=S(=O)(OCC)C</chem>
109	Ethylene dibromide		yes	106-93-4	<chem>BrCCBr</chem>
110	Folpet			133-07-3	<chem>O=C(N(SC(Cl)(Cl)Cl)C(=O)c1ccccc2)c12</chem>
111	Formaldehyde			50-00-0	<chem>O=C</chem>
112	Furazolidone			67-45-8	<chem>O=C(OCC1)N1N=CC(OC(N(=O)=O)=C2)=C2</chem>
113	Furylfuramide (AF-2)			3688-53-7	<chem>NC(=O)C(=Cc1ccc(o1)N(=O)=O)c2ccco2</chem>
114	Glutaraldehyde			111 30-8	<chem>O=CCCCC=O</chem>
115	Glyoxal			4405-13-4	<chem>C12OC(O)C(O)OC1OC(O)C(O)O2</chem>
116	Harmane			486-84-0	<chem>Cc1nccc2c3ccccc3nc12</chem>
117	Hydrazine sulfate (+ S9)	yes		10034-93-2	<chem>NNOS(O)(=O)=O</chem>
118	Hydrogen peroxide			7722-84-1	<chem>OO</chem>
119	Hydroxy urea			127-07-1	<chem>NC(=O)NO</chem>
120	m-Dinitrobenzene			99-65-0	<chem>O=N(=O)c(cccc1N(=O)=O)c1</chem>
121	Methapyrilene (+ S9)	yes	yes	91-80-5	<chem>CN(C)CCN(Cc1cccs1)c2cccn2</chem>
122	Methyl methanesulfonate (MMS)			66-27-3	<chem>O=S(=O)(OC)C</chem>
123	Methylbromide			74-83-9	<chem>BrC</chem>
124	Methylene bromide		yes	74-95-3	<chem>BrCBr</chem>
125	Methylenebromide	yes(?)	yes	74-95-3	<chem>BrCBr</chem>
126	Metronidazol			443-48-1	<chem>Cc1ncc(N(=O)=O)n1CCO</chem>
127	Mitomycin C			50-07-7	<chem>O=C1C(N)=C(C)C(=O)C2=C1C(COC(=O)N)C3(OC)C(N4)C4CN23</chem>
128	Nalidixic acid (free acid)			389-08-2	<chem>Cc1ccc2C(=O)C(C(=O)O)=CN(CC)c2n1</chem>
129	Nalidixic acid (sodium-salt)			3374-05-8	<chem>Cc1ccc2C(=O)C(C(=O)O[Na])=CN(CC)c2n1</chem>
130	N-Ethyl-N*-nitro-N-nitrosoguanidine			4245-77-6	<chem>CCN(N=O)C(=N)NN(=O)=O</chem>

131	N-Ethyl-N-nitrosoarea		759-73-9	<chem>O=C(N(N=O)CC)N</chem>
132	Nifuroxazide		965-52-6	<chem>O=C(NN=Cc(oc(N(=O)(=O))c1)c1)c(ccc(O)c2)c2</chem>
133	Nitrofurantoin		67-20-9	<chem>O=C(N(N=CC(OC(N(=O)=O)=C1)=C1)CC2=O)N2</chem>
134	Nitrofurazone		59-87-0	<chem>O=C(N)NN=CC(OC(N(=O)=O)=C1)=C1</chem>
135	Nitrogen dioxide		10102-44-0	<chem>O=N=O</chem>
136	N-Methyl-N-nitro-N-nitrosoguanidine		70-25-7	<chem>O=N(=O)NC(=N)N(N=O)C</chem>
137	N-Methyl-N-nitrosoarea		684-93-5	<chem>O=C(N(N=O)C)N</chem>
138	N-Nitrosobutylurea		869-01-2	<chem>CCCCN(N=O)C(N)=O</chem>
139	N-Nitrosocimetidine		73785-40-7	
140	N-Nitrosodiethanolamine		1116-54-7	<chem>O=NN(CCO)CCO</chem>
141	N-Nitrosodiphenylamine	yes	86-30-6	<chem>O=NN(c(cccc1)c1)c(cccc2)c2</chem>
142	N-Nitroso-N-butyl-N-propylamine		25413-64-3	<chem>CCCCN(N=O)CCC</chem>
143	Norharmane		244-63-3	<chem>c1ccc2c(c1)nc3cnc23</chem>
144	Nystatin		1400-61-9	<chem>O=C(O)CCCCCCCC=CC=CC=CC=CCC</chem>
145	o-Aminoazotoluene		97-56-3	<chem>N(=Nc(c(ccc1)C)c1)c(ccc(N)c2C)c2</chem>
146	Ofloxacin		82419-36-1	<chem>C1CN(C)CCN1c2c(F)cc3C(=O)C(C(=O)O)=CN4c3c2OCC4C</chem>
147	o-Nitroanisole		91-23-6	<chem>O=N(=O)c(c(OC)ccc1)c1</chem>
148	o-Tolidine	yes	119-93-7	<chem>Nc(c(cc(c(ccc(N)c1C)c1)c2)C)c2</chem>
149	Paraquat		4685-14-7	<chem>Cn1(Cl)ccc(cc1)c2ccn(Cl)(C)cc2</chem>
150	Phenobarbital		50-06-6	<chem>CCC1(C(=O)NC(=O)NC1=O)c2ccccc2</chem>
151	Phenylhydrazine		100-63-0	<chem>N(N)c(cccc1)c1</chem>
152	Pikrinic acid (+ S9)	yes	88-89-1	<chem>O=N(=O)c(cc(N(=O)=O)c(O)c1N(=O)=O)c1</chem>
153	Potassium chromate		11073-34-0	
154	Potassium dichromate (Cr207 2-)		7778-00-9	
155	p-Phenyldiamine	yes	106-50-3	<chem>Nc(ccc(N)c1)c1</chem>
156	Propane sultone		1120-71-4	<chem>O=S(=O)(OCC1)C1</chem>
157	Propylene oxide		75-56-9	<chem>O(C1C)C1</chem>
158	Pyrogallol		87-66-1	<chem>Oc(c(O)ccc1)c1O</chem>
159	Quercetin(+ S9)	yes	117-39-5	<chem>Oc1cc(O)c2C(=O)C(O)=C(c3cc(O)c(O)cc3)Oc2c1</chem>
160	Saccharin		81-07-2	<chem>O=C(NS(=O)(=O)c1cccc2)c12</chem>
161	Selene dioxide		7446-08-4	
162	Sodium nitrite(NO2-)		7632-00-0	<chem>[Na]ON=O</chem>

163	Sterigmatocystin		10048-13-2	<chem>COc4cc2OC1OC=CC1c2c5oc3cccc(O)c3C(=O)c45</chem>
164	Streptonigrin		3930-19-6	<chem>c42C(=O)C(OC)=C(N)C(=O)c4nc(c1nccc(c3c(O)c(OC)c(OC)cc3)c1N)cc2</chem>
165	Streptozotocin		18883-66-4	<chem>CN(N=O)C(=O)NC1C(O)OC(CO)C(O)C1O</chem>
166	Styrene oxide		96-09-3	<chem>O(C1c(cccc2)c2)C1</chem>
167	t-Butylhydroperoxide	yes	75-91-2	<chem>O(O)C(C)(C)C</chem>
168	Trichloroacetone	yes	918-00-3	<chem>ClC(Cl)(Cl)C(=O)C</chem>
169	Trichloronitromethane		76-06-2	<chem>O=N(=O)C(Cl)(Cl)Cl</chem>
170	Tris(2,3-dibromopropyl)phosphate		126-72-7	<chem>O=P(OCC(Br)CBr)(OCC(Br)CBr)OCC(Br)CBr</chem>

Annex 3 Predicted and measured values for Kow

Nr.	Compound	LogKow predicted	LogKow Exp	LogKow selected	References Kow exp	Remarks Kow
1	1,1-Dichloroacetone	0.2		0.2		
2	1,2:3,4-Dibenzoanthracene	6.7	6.41	6.41	Helweg,C et al. (1997a)	
3	1,2:3,4-Diepoxybutane	-0.58	-0.28	-0.28	Deneer,JW et al. (1988)	
4	1,2-Diaminobenzene(+ S9)	0.16	0.15	0.15	Hansch,C et al. (1995)	
5	1,2-Dimethylhydrazine	-0.54		-0.54		
6	1,2-Epoxybutan	0.86		0.86		
7	1,3-Dinitropyrene	4.57		4.57		
8	1,3-Dioxane	0.18		0.18		
9	1,5-Dinitronaphthalene	2.8	2.58	2.58	Debnath,AK et al. (1992)	
10	1,6-Dinitropyrene	4.57		4.57		
11	1,8-Dinitropyrene	4.57		4.57		
12	1-Bromo pentane	3.14	3.37	3.37	Hansch,C et al. (1995)	
13	1-Chloro-2,4-dinitrobenzene	2.27	2.17	2.17	Debnath,AK et al. (1991)	
14	1-Naphthylamine	2.25	2.25	2.25	Hansch,C et al. (1995)	
15	1-Nitronaphthalene	2.99	3.19	3.19	Hansch,C et al. (1995)	
16	1-Nitropyrene	4.75	5.06	5.06	BioByte (1995)	
17	2,4,7-Trinitro-9-fluorene	3.47		3.47		
18	2,4-Diaminotoluene (+ S9)	0.16	0.14	0.14	Debnath,AK et al. (1992)	
19	2,4-Dimethyl-1,3-dioxane	1.02		1.02		
20	2,4-Dinitrotoluene	2.18	1.98	1.98	Hansch,C et al. (1995)	
21	2,6-Diaminotoluene (+ S9)	0.16		0.16		
22	2,7-Dinitro-8-fluorenone	3.19	2.84	2.84	Debnath,AK & Hansch,C (1992)	
23	2-Acetylaminofluorene (+S9)	3.12		3.12		
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	2.12	1.98	1.98	Hansch,C et al. (1995)	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	1.66		1.66		

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
26	2-Amino-3-methyl-9H-pyrido(2,3-b)indole; MeA alpha C (+ S9)	2.48	2.9	2.9	BioByte (1995)
27	2-Amino-6-methyldipyrdo(1,2-a:3',2'd)imidazole (Glu-P-1;+ S9)	2.18	1.75	1.75	Hansch,C et al. (1995)
28	2-Amino-9H-pyrido(2,3-b)indole (+S9)(A alpha C)	1.93	2.6	2.6	Hansch,C et al. (1995)
29	2-Aminoanthracene(+ S9)	3.43		3.43	
30	2-Amino-dipyrdo(1,2-a:3',2'-d)imidazole-2-amine (Glu-P-2; + S9)	1.64	1.38	1.38	Hansch,C et al. (1995)
31	2-Aminofluorene(+S9)	3.1	3.14	3.14	Debnath,AK et al. (1992)
32	2-Amino-1-methyl-6-phenylimidazo[4,5-b]pyridine (PhIP)	2.16	2.23	2.23	Hansch,C et al. (1995)
33	2-Nitrofluorene	3.83	3.37	3.37	Debnath,AK & Hansch,C (1992)
34	2-Nitronaphthalene	2.99	3.24	3.24	Debnath,AK et al. (1992)
35	2-Nitro-p-phenylenediamine	0.55	0.53	0.53	Bronaugh,RL & Congdon,ER (1984)
36	3,7-Dinitrofluoranthene	4.57		4.57	
37	3,9-0initrofluoranthene	4.57		4.57	
38	3-Amino-1,4-dimethyl.SH-pyrido[3,4-b]indole; Trp-P-1 (+ S9)	3.02	1.97	1.97	BioByte (1995)
39	3-Amino-1-methyl-5H-pyrido[4,3-b]indole; Trp-P 2 (+ S9)	2.48	1.97	1.97	Hansch,C et al. (1995)
40	3-Methoxy-4-aminoazobenzene(+S9)	3.27		3.27	
41	3-Methylcholanthrene (+ S9)	7.05	6.42	6.42	Hansch,C et al. (1995)
42	3-Nitrofluoranthene	4.75		4.75	
43	4,4'-Dinitrobiphenyl	3.39		3.39	
44	4-Aminobiphenyl	2.84	2.86	2.86	Martin-Villodre,A et al. (1986)
45	4-Methyl-1,3-dioxane	0.6		0.6	
46	4-Nitro-o-phenylenediamine	0.55	0.88	0.88	Hansch,C et al. (1995)
47	4-Nitroquinoline-N-oxide	0.82	1.09	1.09	Hansch,C et al. (1995)
48	4-Nitrosodiphenylamine	3.16		3.16	
49	4-Nitroso-N,N-dimethylaniline	2.04		2.04	
50	5-Butyrolactone	-0.38		-0.38	

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
51	5-Fluorouracil	-0.81	-0.89	-0.89	Hansch,C et al. (1995)
52	5-MOP (+ UV)	2.14	1.93	1.93	(in press)
53	5-Nitro-2-furaldehyde	0.65	1.01	1.01	Hansch,C et al. (1995)
54	5-Nitro-2-furylacrylacid	0.01	0.65	0.65	Balaz,S et al. (1985)
55	5-Nitroacenaphthene	3.97	3.85	3.85	Debnath,AK et al. (1992)
56	6-Aminochrysene(+ S9)	4.6	4.99	4.99	Hansch,C et al. (1995)
57	6-chloro-9-(3-(2-chloroethylamino)propyl)amino-2-rmethoxyacridine*2HCl(ICR-191)	4.4		4.4	
58	6-Nitrochrysene	2.54	1.37	1.37	BioByte (1995)
59	7,12-Dimethylbanzanthracene (+ S9)	6.62	5.8	5.8	Hansch,C et al. (1995)
60	8-MOP (+ UV)	2.14	2	2	(in press)
61	8-Proplotactone	-0.8		-0.8	
62	9 Aminoacridine	2.4	2.74	2.74	Hansch,C et al. (1995)
63	Acridine Orange	3.67		3.67	
64	Acrinol(+S9)	2.06		2.06	
65	Adriamycin(Adr)	1.85	1.27	1.27	
66	Afiatoxin G1	0.5		0.5	
67	Aflatoxin B1	1.23		1.23	
68	Amphotericin B	-2.8		-2.8	
69	Amsacrine	3.89		3.89	
70	Auramine	2.98		2.98	
71	Avarol(+ S9)	7.02		7.02	
72	Azaserine	-2.36	-2	-2	Ellington,JJ & Stancil,FE (1988)
73	Azinphos-Methyl	2.53	2.75	2.75	Hansch,C et al. (1995)
74	Azobenzene(+ S9)	4.11	3.82	3.82	Hansch,C et al. (1995)
75	Benz(a)anthracene(+ S9)	5.52	5.76	5.76	Wang,L et al. (1986)
76	Benzidine (+ S9)	1.92	1.34	1.34	Hansch,C et al. (1995)
77	Benzo(a)pyrene(+ S9)	6.11	6.13	6.13	De Maagd,PG et al. (1998)

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
78	Benzylchloride	2.79	2.3	2.3	Hansch,C et al. (1995)
79	Bleomycin	-10.1			
80	Bromo acetic acid	0.43	0.41	0.41	Hansch,C et al. (1995)
81	Bromobenzene	2.88	2.99	2.99	Hansch,C et al. (1995)
82	Captan	2.74	2.8	2.8	Tomlin,C (1997)
83	Chloramine T	-0.5		-0.5	
84	Chromium(VI)oxide				
85	Chrysene (+ S9)	5.52	5.81	5.81	De Maagd,PG et al. (1998)
86	Ciprofloxacin	0	0.28	0.28	Takacs-Novak,K et al. (1992)
87	cis-Platinum	-2.75	-2.19	-2.19	Hansch,C et al. (1995)
88	Citrinin	0.45		0.45	
89	Cumenehydroperoxide	2.16		2.16	
90	Cupferron	-1.73		-1.73	
91	Danthron	3.94		3.94	
92	Daunomycin	2.19	1.83	1.83	Sangster (1993)
93	Daunorubicin	2.19	1.83	1.83	Sangster (1993)
94	Dibromoacetone	0.47		0.47	
95	Dichlofluanid	2.72	3.7	3.7	Tomlin,C (1997)
96	Diethylnitrosamine (+ S9)	0.34	0.48	0.48	Hansch,C et al. (1995)
97	Diethylsulfate (+ S9)	1.14	1.14	1.14	Hansch,C et al. (1995)
98	Dimethylnitrosamine (+ S9)	-0.64	-0.57	-0.57	Hansch,C et al. (1995)
99	Dimethylsulfate	0.16		0.16	
100	Dimethylsulfoxide	-1.22	-1.35	-1.35	Hansch,C et al. (1995)
101	Dithianone	2.98	2.84	2.84	Hansch,C et al. (1995)
102	Doxorubicinhydrochloride				
103	Ellipticine	4.47	4.8	4.8	Hansch,C et al. (1995)
104	Enoxacin	-0.21	-0.2	-0.2	Sangster (1994) (ion-correct)
105	Epichlorhydrin	0.63	0.45	0.45	Deneer,JW et al. (1988)

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
106	Epoxy styrene	1.59	1.61	1.61	Hansch,C et al. (1995)
107	Ethidium bromide (+ S9)	-0.38		-0.38	
108	Ethyl methanesulfonate (EMS)	-0.17		-0.17	
109	Ethylene dibromide	2.01	1.96	1.96	Hansch,C et al. (1995)
110	Folpet	2.84	2.85	2.85	Hansch,C et al. (1995)
111	Formaldehyde	0.35	0.35	0.35	Hansch,C et al. (1995)
112	Furazolidone	1.02	-0.04	-0.04	Debnath,AK et al. (1991)
113	Furylfuramide (AF-2)	1.14	0.15	0.15	BioByte (1995)
114	Glutaraldehyde	-0.18		-0.18	
115	Glyoxal	-2.89		-2.89	
116	Harmane	2.75	3.1	3.1	BioByte (1995)
117	Hydrazine sulfate (+ S9)	-4.05		-4.05	
118	Hydrogen peroxide	-1.57		-1.57	
119	Hydroxy urea	-1.68	-1.8	-1.8	Hansch,C et al. (1995)
120	m-Dinitrobenzene	1.63	1.49	1.49	Hansch,C et al. (1995)
121	Methapyrilene (+ S9)	2.55	2.87	2.87	Sangster (1994)
122	Methyl methanesulfonate (MMS)	-0.66		-0.66	
123	Methylbromide	1.18	1.19	1.19	Hansch,C et al. (1995)
124	Methylene bromide	1.52	1.7	1.7	Martiska,A & Bekarek,V (1990)
125	Methylenebromide	1.52	1.7	1.7	Martiska,A & Bekarek,V (1990)
126	Metronidazol	0	-0.02	-0.02	Hansch,C et al. (1995)
127	Mitomycin C	-1.18	-0.4	-0.4	Hansch,C et al. (1995)
128	Nalidixic acid (free acid)	1.64	1.59	1.59	BioByte (1995)
129	Nalidixic acid (sodium-salt)	-2.16			
130	N-Ethyl-N*-nitro-N-nitrosoguanidine	-0.43		-0.43	
131	N-Ethyl-N-nitrosourea	-0.02	0.23	0.23	Hansch,C et al. (1995)
132	Nifuroxazide	1.49		1.49	
133	Nitrofurantoin	-0.17	-0.47	-0.47	Hansch,C et al. (1995)

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
134	Nitrofurazone	0.23	0.23	0.23	Hansch,C et al. (1995)
135	Nitrogen dioxide				
136	N-Methyl-N-nitro-N-nitrosoguanidine	-0.92		-0.92	
137	N-Methyl-N-nitrosoourea	-0.52	-0.03	-0.03	Hansch,C et al. (1995)
138	N-Nitrosobutylurea	0.96	1.04	1.04	Hansch,C et al. (1995)
139	N-Nitrosocimetidine				
140	N-Nitrosodiethanolamine	-1.28		-1.28	
141	N-Nitrosodiphenylamine	3.16	3.13	3.13	Hansch,C et al. (1995)
142	N-Nitroso-N-butyl-N-propylamine	1.82	2.09	2.09	Vera,A et al. (1992)
143	Norharmane	2.2	3.17	3.17	Biagi,GL et al. (1989)
144	Nystatin	7.08		7.08	
145	o-Aminoazotoluene	4.29		4.29	
146	Ofloxacin	-0.2	-0.39	-0.39	Hansch,C et al. (1995)
147	o-Nitroanisol	1.89	1.73	1.73	Hansch,C et al. (1995)
148	o-Tolidine	3.02	2.34	2.34	Hansch,C et al. (1995)
149	Paraquat	-2.71		-2.71	
150	Phenobarbital	1.33	1.47	1.47	Hansch,C et al. (1995)
151	Phenylhydrazine	0.79	1.25	1.25	Hansch,C et al. (1995)
152	Pikrinic acid (+ S9)	1.54	1.33	1.33	Sangster (1994)
153	Potassium chromate				
154	Potassium dichromate (Cr207 2-)				
155	p-Phenylendiamine	-0.39	-0.3	-0.3	Hansch,C et al. (1995)
156	Propane sultone	-0.28		-0.28	
157	Propylene oxide	0.37	0.03	0.03	Hansch,C et al. (1995)
158	Pyrogallol	0.97		0.97	
159	Quercetin(+ S9)	1.48		1.48	
160	Saccharin	0.45	0.91	0.91	Hansch,C et al. (1995)
161	Selene dioxide				

Nr.	Compound	LogKow pre- dicted	LogKow Exp	LogKow selected	References Kow exp
162	Sodium nitrite(NO ₂ -)				
163	Sterigmatocystin	3.81		3.81	
164	Streptonigrin	-0.56		-0.56	
165	Streptozotocin	-1.61	-1.45	-1.45	Hansch,C et al. (1995)
166	Styrene oxide	1.59	1.61	1.61	Hansch,C et al. (1995)
167	t-Butylhydroperoxide	0.94		0.94	
168	Trichloroacetone	1.12		1.12	
169	Trichloronitromethane	1.32	2.09	2.09	Hansch,C et al. (1995)
170	Tris(2,3-dibromopropyl)phosphate	4.19	4.29	4.29	Sangster (1994) (avg)

References experimentally determined K_{ow} values
as cited in Meylan and Howard (1999c)

- Akacs-Novak, K., M. Jozan, I. Hermeicz and G. Szasz (1992). Lipophilicity antibacterial fluoroquinolones, *Int. J. Pharm.* 79(2-3) 89-96.
- Balaz, S., A. Kuchar, E. Sturdik, M. Rosenberg, L. Stibranyi and D. Ilavsky (1985). Interphase distributions of 2-furylethylenes, *Coll. Czech. Chem. Commun.*, 50(8), 1642-1647.
- Balaz, S., D. Ilavsky, E. Sturdik and J. Kovac (1985). Antimicrobial activity of methyl esters and nitriles of 2-cyano-3-(5'-R-2'-furyl)propenoic acid, *Folia Microbiol.* (Prague), 30, 34-41.
- Biagi, G.L., M.C. Pietrogrande, A.M. Barbaro, M.C. Guerra, P.A. Borea and G.C. Forti (1989). Study of the lipophilic character of a series of beta-carbolines, *J. Chromatogr.* 469, 121-126.
- BioByte (1995). ClogP for Windows Program (v1.0). BioByte Corp., Claremont, CA.
- Bronaugh, R.L. and E.R. Congdon (1984). Percutaneous absorption of hair dyes: correlation with partition coefficients, *J. Invest. Dermatol.* 83(2), 124-127.
- De Maagd, P.G. (1998). Physicochemical properties of polycyclic aromatic hydrocarbons: aqueous solubilities, n-octanol/water partition coefficients, and Henry's law constants. *Environ. Toxicol. Chem.* 17: 251-257.
- Debnath, A.K. and C. Hansch (1992). Structure-activity relationship of genotoxic polycyclic aromatic nitro compounds. Further evidence for the importance of hydrophobicity and molecular orbital energies in genetic toxicity, *Environ. Mol. Mutagen.* 20(2), 140-144.
- Debnath, A.K., G. Debnath, A.J. Shusterman and C. Hansch (1992). A QSAR investigation of the role of hydrophobicity in regulating mutagenicity in the Ames test.1., *Environ. Mol. Mutagen.* 19(1), 37-52.
- Debnath, A.K., R.L. Lopez de Compadre, G. Debnath, A.J. Shusterman and C. Hansch (1991). Structure-activity relationship of mutagenic aromatic and heteroaromatic nitro compounds. Correlation with molecular orbital energies and hydrophobicity, *J. Med. Chem.* 34(2), 786-797.
- Deneer, J.W., T.L. Sinnige, W. Seinen and J.L.M. Hermens (1988). A quantitative structure-activity relationship for the acute toxicity of some epoxy compounds to the guppy, *Aquatic Toxicol.* 13(3), 195-204.
- Ellington, J.J. and F.E. Stancil Jr. (1988). *Octanol/water coefficients for evaluation of hazardous waste land disposal: selected chemicals*. EPA/600/M-88-010, Aug 1988; Athens, GA: USEPA (Environ. Res. Lab.), NTIS PB89-120760.
- Ellington, J.J. and T.L. Floyd (1996). *Octanol/water partition coefficients for eight phthalate esters*. EPA/600/S-96/006, Sept. 1996; Athens, GA: USEPA (National Exposure Research Lab).
- Hansch, C., A. Leo and D. Hoekman (1995). *Exploring QSAR. Hydrophobic, Electronic, and Steric Constants*. ACS Professional Reference Book. Washington, DC: American Chemical Society.
- Helweg, C., T. Nielsen and P.E. Hansen (1997). Determination of octanol-water partition coefficients of polar polycyclic aromatic compounds (N-PAC) by high performance liquid chromatography. *Chemosphere* 34(8): 1673-1684.
- Helweg, C., T. Nielsen and P.E. Hansen (1997a). QSPRs compared with HPLC experiments for estimation of K_{ow} for N-PACs. (Chapter 12) In: *Quantitative Structure-Activity Relationships in Environmental Sciences - VII*. editors F. Chen and G. Schuurman (Proceedings of QSAR 96 24-28 June Elsinore, Denmark. SETAC Special Publication Series, Pensacola, FL

- Martin-Villodre, A., J.M. Pla-Delfina, J. Moreno, M.D. Perez-Buendia, J. Miralles, E.F. Collado, E. Sanchez-Moyano and A. Del Pozo (1986). Studies on the reliability of a bihyperbolic functional absorption model. I. Ring-substituted anilines, *J. Pharmacokinet. Biopharm.*, 14(6), 615-633.
- Martiska, A. and V. Bekarek (1990). Application of the effective Born's relative permittivity functions for the evaluation of polarity effects of solutes on the octanol-water partition coefficients. *Acta Univ. Palacki Olomuc., Fac. Rerum Nat.* 97(Chem. 29): 63-7.
- Sangster, J. (1993). LOGKOW Databank. A databank of evaluated octanol-water partition coefficients (Log P) on microcomputer diskette. Montreal, Quebec, Canada: Sangster Research Laboratories.
- Sangster, J. (1994). LOGKOW Databank. A databank of evaluated octanol-water partition coefficients (Log P) on microcomputer diskette. Montreal, Quebec, Canada: Sangster Research Laboratories.
- Tomlin, C. (1997). *The Pesticide Manual*, Eleventh edition, Crop Protection Publications; British Crop Protection Council, 49 Downing St, Farnham, Surrey GU9 7PH, United Kingdom.
- Vera, A. et al. (1992). Quantitative structure-activity relationship study of the biophysicochemical behaviour of nitrosamine. *J. Pharm. Sci.* 81: 791-6.
- Wang, L., X. Wang, O. Xu and L. Tian (1986). Determination of the n-octanol/water partition coefficients of polycyclic aromatic hydrocarbons and estimation of their aqueous solubilities, *Huanjing Kexue Xuebao*, 6(4), 491-497 (1986).

Annex 4 Predicted and experimentally determined Henry Law constants

Nr.	Compound	CAS.Reg.No.	H bond Pa.m ³ /Mol	H group	H exp	H selected	LogH selected	Ref HLC
1	1,1-Dichloroacetone	513-88-2	6.03E-01			6.03E-01	-0.22	
2	1,2:3,4-Dibenzoanthracene	215-58-7	4.80E-02	1.21E-02		4.80E-02	-1.32	
3	1,2:3,4-Diepoxybutane	298-18-0	8.08E-02	3.47E-03		8.08E-02	-1.09	
4	1,2-Diaminobenzene(+ S9)	95-54-5	6.60E-05	8.71E-05	7.06E-04	7.06E-04	-3.15	SRC
5	1,2-Dimethylhydrazine	540-73-8	6.82E-03			6.82E-03	-2.17	
6	1,2-Epoxybutan	106-88-7	2.08E+01	1.70E+01	1.77E+01	1.77E+01	1.25	SRC
7	1,3-Dinitropyrene	75321-20-9	1.27E-05	1.99E-05		1.27E-05	-4.90	
8	1,3-Dioxane	505-22-6	2.90E+00	2.95E+00		2.90E+00	0.46	
9	1,5-Dinitronaphthalene	605-71-0	8.03E-04	2.51E-03		8.03E-04	-3.10	
10	1,6-Dinitropyrene	42397-64-8	1.27E-05	1.99E-05		1.27E-05	-4.90	
11	1,8-Dinitropyrene	42397-65-9	1.27E-05	1.99E-05		1.27E-05	-4.90	
12	1-Bromo pentane	110-53-2	2.59E+03	2.24E+03	1.93E+03	1.93E+03	3.29	SRC
13	1-Chloro-2,4-dinitrobenzene	97-00-7	6.10E-03	3.09E-02		6.10E-03	-2.21	
14	1-Naphthylamine	134-32-7	1.82E-02	1.48E-02	1.09E-02	1.09E-02	-1.96	ABRAHAM,MH ET AL. (1994)
15	1-Nitronaphthalene	86-57-7	2.04E-01	3.02E-01	1.73E-01	1.73E-01	-0.76	ALTSCHUH,J ET AL. (1999)
16	1-Nitropyrene	5522-43-0	3.21E-03	2.40E-03		3.21E-03	-2.49	
17	2,4,7-Trinitro-9-fluorene	34263-36-0	1.00E-06	1.95E-06		1.00E-06	-6.00	
18	2,4-Diaminotoluene (+ S9)	95-80-7	7.29E-05	9.34E-05		7.29E-05	-4.14	
19	2,4-Dimethyl-1,3-dioxane	766-20-1	5.12E+00	1.86E+00		5.12E+00	0.71	
20	2,4-Dinitrotoluene	121-14-2	9.08E-03	3.89E-02	5.30E-03	5.30E-03	-2.28	ALTSCHUH,J ET AL. (1999)
21	2,6-Diaminotoluene (+ S9)	823-40-5	7.29E-05	9.34E-05		7.29E-05	-4.14	
22	2,7-Dinitro-8-fluorenone	31551-45-8	1.03E-06			1.03E-06	-5.99	
23	2-Acetylaminofluorene (+S9)	53-96-3	1.88E-05			1.88E-05	-4.73	

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MelQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
26	2-Amino-3-methyl-9H-pyrido(2,3-b)indole; MeA alpha C (+ S9)	68006-83-7	4.32E-09			4.32E-09	-8.36	
27	2-Amino-6-methyldipyrido(1,2-a:3',2'd)imidazole (Glu-P-1;+ S9)	67730-11-4	1.11E-10			1.11E-10	-9.96	
28	2-Amino-9H-pyrido(2,3-b)indole (+S9)(A alpha C)	26148-68-5	3.92E-09			3.92E-09	-8.41	
29	2-Aminoanthracene(+ S9)	613-13-8	1.78E-03	1.02E-03		1.78E-03	-2.75	
30	2-Amino-dipyrido(1,2-a:3',2'-d)imidazole-2-amine (Glu-P-2; + S9)	67730-1 0-3	1.00E-10			1.00E-10	-10.00	
31	2-Aminofluorene(+S9)	153-78-6	4.45E-03	1.38E-03		4.45E-03	-2.35	
32	2-Amino-1-methyl-6-phenylimidazo[4,5-b]pyridine (PhIP)	105650-23-5	2.73E-08			2.73E-08	-7.56	
33	2-Nitrofluorene	607-57-8	6.47E-02	2.81E-02		6.47E-02	-1.19	
34	2-Nitronaphthalene	581-89-5	2.04E-01	3.02E-01		2.04E-01	-0.69	
35	2-Nitro-p-phenylenediamine	5307-14-2	5.70E-06	7.25E-07		5.70E-06	-5.24	
36	3,7-Dinitrofluoranthene	105735-71-5	1.27E-05	1.99E-05		1.27E-05	-4.90	
37	3,9-Dinitrofluoranthene	22506-53-2	1.27E-05	1.99E-05		1.27E-05	-4.90	
38	3-Amino-1,4-dimethyl-SH-pyrido[3,4-b]indole; Trp-P-1 (+ S9)	62450-06-0	4.78E-09			4.78E-09	-8.32	
39	3-Amino-1-methyl-5H-pyrido[4,3-b]indole; Trp-P 2 (+ S9)	62450-07-1	4.32E-09			4.32E-09	-8.36	

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
40	3-Methoxy-4-aminoazobenzene(+S9)	3544-23-8	3.01E-05	5.01E-04		3.01E-05	-4.52	
41	3-Methylcholanthrene (+ S9)	56-49-5	2.91E-01	2.95E-02		2.91E-01	-0.54	
42	3-Nitrofluoranthene	892-21-7	3.21E-03	2.40E-03		3.21E-03	-2.49	
43	4,4'-Dinitrobiphenyl	1528-74-1	6.32E-04	2.88E-03		6.32E-04	-3.20	
44	4-Aminobiphenyl	92-67-1	1.43E-02	1.70E-02		1.43E-02	-1.84	
45	4-Methyl-1,3-dioxane	1120-97-4	3.85E+00	6.91E+00		3.85E+00	0.59	
46	4-Nitro-o-phenylenediamine	99-56-9	2.60E-07	7.25E-07		2.60E-07	-6.59	
47	4-Nitroquinoline-N-oxide	56-57-5	2.67E-09			2.67E-09	-8.57	
48	4-Nitrosodiphenylamine	156-10-5	1.08E-03			1.08E-03	-2.97	
49	4-Nitroso-N,N-dimethylaniline	138-89-6	8.79E-02			8.79E-02	-1.06	
50	5-Butyrolactone	36536-46-6	1.01E+01	2.95E+00		1.01E+01	1.00	
51	5-Fluorouracil	51-21-8	1.63E-05			1.63E-05	-4.79	
52	5-MOP (+ UV)	484-20-8	3.92E-03			3.92E-03	-2.41	
53	5-Nitro-2-furaldehyde	698-63-5	5.19E-03			5.19E-03	-2.29	
54	5-Nitro-2-furylacrylacid	710-25-8	1.02E-07			1.02E-07	-6.99	
55	5-Nitroacenaphthene	602-87-9	1.09E-01	4.79E-02		1.09E-01	-0.96	
56	6-Aminochrysene(+ S9)	2642-98-0	1.74E-04	7.08E-05		1.74E-04	-3.76	
57	6-chloro-9-(3-(2-chloroethylamino)propyl)amino-2-methoxyacridine*2HCl(ICR-191)	17070-45-0	4.17E-11			4.17E-11	-10.38	
58	6-Nitrochrysene	2-89-6	1.90E-14			1.90E-14	-13.72	
59	7,12-Dimethylbanzanthracene (+ S9)	57-97-6	5.98E-01	1.99E-01		5.98E-01	-0.22	
60	8-MOP (+ UV)	298-81-7	3.92E-03			3.92E-03	-2.41	

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
61	8-Propylotactone	57-57-8	7.58E+00	1.26E+00		7.58E+00	0.88	
62	9 Aminoacridine	90-45-9	2.32E-06	1.59E-05		2.32E-06	-5.63	
63	Acridine Orange	494-38-2	1.67E-06			1.67E-06	-5.78	
64	Acrinol(+S9)	1837-57-6	6.46E-11	7.08E-09		6.46E-11	-10.19	
65	Adriamycin(Adr)	25316-40-9	2.19E-18			2.19E-18	-17.66	
66	Afiatoxin G1	1165-39-5	4.86E-08			4.86E-08	-7.31	
67	Aflatoxin B1	1162-65-8	1.37E-08			1.37E-08	-7.86	
68	Amphotericin B	1397-89-3	2.07E-30			2.07E-30	-29.68	
69	Amsacrine	51264-14-3	1.65E-11			1.65E-11	-10.78	
70	Auramine	2465-27-2	3.57E-04			3.57E-04	-3.45	
71	Avarol(+ S9)	55303-98-5	6.73E-05			6.73E-05	-4.17	
72	Azaserine	115-02-6	3.19E-06			3.19E-06	-5.50	
73	Azinphos-Methyl	86-50-0	2.80E-05		2.34E-03	2.34E-03	-2.63	SRC
74	Azobenzene(+ S9)	103-33-3	1.44E+00	1.51E+00	1.32E+00	1.32E+00	0.12	SRC
75	Benz(a)anthracene(+ S9)	56-55-3	4.91E-01	1.74E-01	1.18E+00	1.18E+00	0.07	BAMFORD,HA ET AL. (1999)
76	Benzidine (+ S9)	92-87-5	5.07E-06	6.91E-06		5.07E-06	-5.29	
77	Benzo(a)pyrene(+ S9)	50-32-8	7.94E-02	1.99E-02	4.48E-02	4.48E-02	-1.35	TEN HULSCHER,TEM ET AL. (1992)
78	Benzylchloride	100-44-7	2.05E+02	3.89E+01	4.04E+01	4.04E+01	1.61	SRC
79	Bleomycin	11056-06-7						
80	Bromo acetic acid	79-08-3	6.19E-03	8.71E-03		6.19E-03	-2.21	
81	Bromobenzene	108-86-1	2.11E+02	2.19E+02	2.42E+02	2.42E+02	2.38	SHIU,WY & MACKAY,D (1997)
82	Captan	133-06-2	4.50E-04		6.86E-04	6.86E-04	-3.16	SRC

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
83	Chloramine T	127-65-1	8.45E-01			8.45E-01	-0.07	
84	Chromium(VI)oxide	1333-82-0	0.00E+00	0.00E+00				
85	Chrysene (+ S9)	218-01-9	4.91E-01	1.74E-01	5.13E-01	5.13E-01	-0.29	BAMFORD,HA ET AL. (1999)
86	Ciprofloxacin	85721-33-1	4.99E-14			4.99E-14	-13.30	
87	cis-Platinum	15663-27-1	#VALUE!			0.00E+00		
88	Citrinin	518-75-2	2.06E-10			2.06E-10	-9.69	
89	Cumenehydroperoxide	80-15-9	9.51E-02		2.12E-02	2.12E-02	-1.67	SRC
90	Cupferron	135-20-6	3.55E-04			3.55E-04	-3.45	
91	Danthron	117-10-2	5.34E-06			5.34E-06	-5.27	
92	Daunomycin	23541-50-6	1.40E-20			1.40E-20	-19.85	
93	Daunorubicin	20830-81-3	1.40E-20			1.40E-20	-19.85	
94	Dibromoacetonitrile	3252-43-5	3.98E-02			3.98E-02	-1.40	
95	Dichlofluanid	1085-98-9	6.61E-02		3.71E-03	3.71E-03	-2.43	SRC
96	Diethylnitrosamine (+ S9)	55-18-5	3.56E-01	1.70E-01	3.56E-01	3.56E-01	-0.45	MIRVISH,SS ET AL. (1976)
97	Diethylsulfate (+ S9)	64-67-5	4.50E-01		8.24E-01	8.24E-01	-0.08	SRC
98	Dimethylnitrosamine (+ S9)	62-75-9	2.02E-01	1.18E-01	1.78E-01	1.78E-01	-0.75	MIRVISH,SS ET AL. (1976)
99	Dimethylsulfate	77-78-1	2.55E-01		3.92E-01	3.92E-01	-0.41	SRC
100	Dimethylsulfoxide	67-68-5	4.86E-03	1.86E+02	1.48E-04	1.48E-04	-3.83	TAFT,RW ET AL. (1985)
101	Dithianone	3347-22-6	3.81E-12		5.54E-06	5.54E-06	-5.26	SRC
102	Doxorubicinhydrochloride	25316-40-9	0.00E+00	0.00E+00				
103	Ellipticine	519-23-3	1.32E-06	4.46E-06		1.32E-06	-5.88	
104	Enoxacin	74011-58-8	1.12E-16			1.12E-16	-15.95	
105	Epichlorhydrin	106-89-8	5.51E+00	2.57E-01	2.98E+00	2.98E+00	0.47	SRC
106	Epoxystyrene	96-09-3	9.51E-01		1.55E+00	1.55E+00	0.19	SRC

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
107	Ethidium bromide (+ S9)	1239-45-8	1.78E-15			1.78E-15	-14.75	
108	Ethyl methanesulfonate (EMS)	62-50-0	5.25E-01		2.54E-02	2.54E-02	-1.60	SRC
109	Ethylene dibromide	106-93-4	1.27E+02	1.51E+01	6.54E+01	6.54E+01	1.82	SRC
110	Folpet	133-07-3	1.51E-04		7.51E-03	7.51E-03	-2.12	SRC
111	Formaldehyde	50-00-0	9.11E+00	6.02E+00	3.30E-02	3.30E-02	-1.48	BETTERTON,EA & HOFFMAN,MR (1988)
112	Furazolidone	67-45-8	3.20E-06			3.20E-06	-5.50	
113	Furylfuramide (AF-2)	3688-53-7	4.56E-09			4.56E-09	-8.34	
114	Glutaraldehyde	111 30-8	1.08E-02	2.34E-03		1.08E-02	-1.97	
115	Glyoxal	4405-13-4	3.14E-13	9.12E-19		3.14E-13	-12.50	
116	Harmane	486-84-0	1.23E-05	7.77E-05		1.23E-05	-4.91	
117	Hydrazine sulfate (+ S9)	10034-93-2	8.55E-09			8.55E-09	-8.07	
118	Hydrogen peroxide	7722-84-1	2.09E-03			2.09E-03	-2.68	
119	Hydroxy urea	127-07-1	5.32E-06			5.32E-06	-5.27	
120	m-Dinitrobenzene	99-65-0	8.23E-03	3.63E-02	4.81E-03	4.81E-03	-2.32	ALTSCHUH,J ET AL. (1999)
121	Methapyrilene (+ S9)	91-80-5	3.18E-07			3.18E-07	-6.50	
122	Methyl methanesulfonate (MMS)	66-27-3	3.95E-01			3.95E-01	-0.40	
123	Methylbromide	74-83-9	8.33E+02	6.17E+02	6.12E+02	6.12E+02	2.79	SRC
124	Methylene bromide	74-95-3	9.59E+01	8.51E+01	8.06E+01	8.06E+01	1.91	MOORE,RM ET AL. (1995)
125	Methylenebromide	74-95-3	9.59E+01	8.51E+01	8.06E+01	8.06E+01	1.91	MOORE,RM ET AL. (1995)
126	Metronidazol	443-48-1	1.66E-06			1.66E-06	-5.78	
127	Mitomycin C	50-07-7	1.12E-19			1.12E-19	-18.95	
128	Nalidixic acid (free acid)	389-08-2	5.02E-11			5.02E-11	-10.30	
129	Nalidixic acid (sodium-salt)	3374-05-8	5.04E-11			5.04E-11	-10.30	
130	N-Ethyl-N*-nitro-N-	4245-77-6	1.59E-07			1.59E-07	-6.80	

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MeIQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9) nitrosoguanidine	76180-96-6	1.79E-09			1.79E-09	-8.75	
131	N-Ethyl-N-nitrosoourea	759-73-9	1.29E-05			1.29E-05	-4.89	
132	Nifuroxazide	965-52-6	1.92E-11			1.92E-11	-10.72	
133	Nitrofurantoin	67-20-9	1.30E-07			1.30E-07	-6.88	
134	Nitrofurazone	59-87-0	3.04E-08			3.04E-08	-7.52	
135	Nitrogen dioxide	10102-44-0						
136	N-Methyl-N-nitro-N-nitrosoguanidine	70-25-7	1.20E-07			1.20E-07	-6.92	
137	N-Methyl-N-nitrosoourea	684-93-5	9.72E-06			9.72E-06	-5.01	
138	N-Nitrosobutylurea	869-01-2	2.28E-05			2.28E-05	-4.64	
139	N-Nitrosocimetidine	73785-40-7	0.00E+00	0.00E+00				
140	N-Nitrosodiethanolamine	1116-54-7	4.76E-07	2.24E-11		4.76E-07	-6.32	
141	N-Nitrosodiphenylamine	86-30-6	1.19E-01			1.19E-01	-0.93	
142	N-Nitroso-N-butyl-N-propylamine	25413-64-3	8.32E-01	4.79E-01		8.32E-01	-0.08	
143	Norharmane	244-63-3	1.11E-05	5.62E-05		1.11E-05	-4.96	
144	Nystatin	1400-61-9	1.62E+00	1.99E-02		1.62E+00	0.21	
145	o-Aminoazotoluene	97-56-3	6.21E-04	7.08E-04	3.12E-03	3.12E-03	-2.51	SRC
146	Ofloxacin	82419-36-1	4.88E-15		0.00E+00	4.88E-15	-14.31	
147	o-Nitroanisol	91-23-6	1.24E-01	3.55E+00	4.21E-02	4.21E-02	-1.38	SRC
148	o-Tolidine	119-93-7	6.17E-06	7.94E-06		6.17E-06	-5.21	
149	Paraquat	4685-14-7	3.16E-08			3.16E-08	-7.50	
150	Phenobarbital	50-06-6	1.62E-09			1.62E-09	-8.79	
151	Phenylhydrazine	100-63-0	7.80E-04		4.33E-04	4.33E-04	-3.36	SRC
152	Pikrinic acid (+ S9)	88-89-1	1.07E-05	3.72E-08	1.67E-06	1.67E-06	-5.78	SRC
153	Potassium chromate	11073-34-0						

Nr.	Compound	CAS.Reg.No.	H bond	H group	H exp	H selected	LogH selected	Ref HLC
24	2-Amino-3,4-dimethyl-3H-imidazol(4,5-f)quinoline; MelQ (+ S9)	77094-11-2	3.82E-08			3.82E-08	-7.42	
25	2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline; IQ (+ S9)	76180-96-6	1.79E-09			1.79E-09	-8.75	
154	Potassium dichromate (Cr207 2-)	7778-00-9						
155	p-Phenylendiamine	106-50-3	6.60E-05	8.71E-05		6.60E-05	-4.18	
156	Propane sultone	1120-71-4	2.31E-01			2.31E-01	-0.64	
157	Propylene oxide	75-56-9	1.57E+01	1.21E+01	6.83E+00	6.83E+00	0.83	SRC
158	Pyrogallol	87-66-1	5.95E-10	9.77E-10		5.95E-10	-9.23	
159	Quercetin(+ S9)	117-39-5	6.47E-16			6.47E-16	-15.19	
160	Saccharin	81-07-2	1.21E-04			1.21E-04	-3.92	
161	Selene dioxide	7446-08-4						
162	Sodium nitrite(NO2-)	7632-00-0						
163	Sterigmatocystin	10048-13-2	2.78E-07			2.78E-07	-6.56	
164	Streptonigrin	3930-19-6	7.08E-25			7.08E-25	-24.15	
165	Streptozotocin	18883-66-4	7.70E-17			7.70E-17	-16.11	
166	Styrene oxide	96-09-3	9.51E-01		1.55E+00	1.55E+00	0.19	SRC
167	t-Butylhydroperoxide	75-91-2	1.57E+00			1.57E+00	0.20	
168	Trichloroacetone	918-00-3	2.13E-01			2.13E-01	-0.67	
169	Trichloronitromethane	76-06-2	1.80E-01		2.01E+02	2.01E+02	2.30	KAWAMOTO,K & URANO,K (1989)
170	Tris(2,3-dibromopropyl)phosphate	126-72-7	3.14E-07			3.14E-07	-6.50	

References experimentally determined HLCs
as cited in Meylan and Howard (1999b)


Sources containing experimentally measured Henry's law constants included the following selected references. Source data cited as SRC relate to HLCs calculated from available experimentally determined vapour pressure and solubility data in the PhysProp or the EFDB databases.

- Abraham, Mh Et Al. (1994). Abraham, M.H.; Chadha, H.S.; Whiting, G.S.; Mitchell, R.C.; Hydrogen Bonding. 32. An Analysis Of Water-Octanol And Water-Alkane Partitioning And The Log P Parameter Of Seiler.; *J. Pharm. Sci.*; 83:1085-100.
- Altschuh, J Et Al. (1999). *Chemosphere* 39(11): 1871-1887.
- Bamford, Ha Et Al. (1999). *Environ. Toxicol. Chem.* 18(9): 1905-1912.
- Betterton, E.A., M.R. Hoffmann (1988). Henry's Law Constants Of Some Environmentally Important Aldehydes, *Environ Sci Technol* 22:1415-8.
- Kawamoto, K., K. Urano (1989). Parameters For Predicting Fate Of Organochlorine Pesticides In The Environment. (I) Octanol-Water And Air-Water Partition Coefficients, *Chemosphere* 18:1987-96.
- Mirvish, S.S., P. Issenberg, H.C. Sornson (1976). Air-Water And Ether-Water Distribution Of N-Nitroso Compounds: Implications For Laboratory Safety, Analytic Methodology And Carcinogenicity For The Rat Esophagus, Nose And Liver, *J Natl Cancer Inst* 56:1125-9.
- Moore, R.M., C.E. Geen, V.K. Tait (1995). Determination Of Henry's Law Constants For A Suite Of Naturally Occurring Halogenated Methanes In Seawater, *Chemosphere* 30:1183-91.
- Shiu, Wy, D. Mackay (1997). Henry's Law Constants Of Selected Aromatic Hydrocarbons, Alcohols, And Ketones; *J Chem Eng Data* 42: 27-30.
- Taft, R.W., M.H. Abraham, R.M. Dougherty, M.J. Kamlet, (1985). The Molecular Properties Governing Solubilities On Nonelectrolytes In Water.; *Nature* 313:384-6.
- Ten Hulscher (1992). *Tem; Environ Toxicol Chem* 11: 1595-603.

Annex 5 Detailed description of the selected SPE approaches

C18 + Isolut ENV+	Description
<i>Water samples</i>	500 ml or more: filtration if needed with 0.45 mm membrane filters and pH adjustment (pH 7.0). Add some 10% of methanol to the sample to improve the desorption of hydrophobic compounds.
<i>C18 silica cartridge</i>	500 mg sorbent in 6 ml cartridge.
<i>Conditioning</i>	7 ml methanol and then 3 ml of HPLC-water at 1 ml/min: the sorbent is not allowed to become dry before the percolation of the sample.
<i>Percolation</i>	Percolation of the filtered sample (pH 7.0) into the C18 cartridge at 5 ml/min.
<i>Drying</i>	After sample percolation, the sorbent is completely dried under vacuum for 20–30 min.
<i>Desorption</i>	Differential elution can be applied to C18 cartridge to obtain three fractions (F) containing contaminants with different polarities and functional groups: F1: 2x5 ml hexane; F2: 2x5 ml dichloromethane:hexane (4:1, v/v) and F3: methanol:dichloromethane (9:1, v/v). If this fractionation is not needed, the desorption can be done only with F3. Wait 1–5 min between the two aliquots (2x5 ml) to allow sufficient contact time between the solvent and the trapped analytes.
<i>Isolut ENV+ or LiChrolut</i>	Same conditioning as by C18 cartridge. The residual (the C18-preconcentrated) water is acidified to pH 3–3.5 and loaded onto the Isolut ENV+ cartridge (200 mg, 6 ml) at 15 ml/min
<i>Drying</i>	After sample percolation, the sorbent is completely dried under vacuum for 20–30 min.
<i>Desorption</i>	Desorption by 2x5 ml of methanol or 2x5 ml of 1 ml (5 mM triethylamine and 5 mM acetic acid, pH=6.5) and 9 ml of methanol at 1 ml/min waiting 5 min between the two aliquots to allow sufficient contact time between the solvent and the trapped analytes.
<i>Evaporation</i>	Combine all the extracts and evaporate under a gentle stream of nitrogen to the desired volume (0.5 ml). Before starting the evaporation, add some dimethyl sulfoxide (DMSO) to avoid loss of volatiles during evaporation.
<i>UmuC-test</i>	The extract is ready for the UmuC-test.

C18 + GCB	Description
Water samples	500 ml or more: filtration if needed with 0.45 mm membrane filters and pH adjustment (pH 7.0). Add some 10% of methanol to the sample to improve the desorption of hydrophobic compounds.
<i>C18 silica cartridge</i>	500 mg sorbent in 6 ml cartridge.
<i>Conditioning</i>	10 ml methanol and then 10 ml of HPLC-water at 1 ml/min: the sorbent is not allowed to become dry before the percolation of the sample.
<i>Percolation</i>	Percolation of the filtered sample (pH 7.0) into the C18 cartridge at 5 ml/min.
<i>Drying</i>	After sample percolation, the sorbent is completely dried under vacuum for 20–30 min.
<i>Desorption</i>	Differential elution can be applied to C18 cartridge to obtain three fractions (F) containing contaminants with different polarities and functional groups: F1: 2x5 ml hexane; F2: 2x5 ml dichloromethane:hexane (4:1, v/v) and F3: methanol:dichloromethane (9:1, v/v). If this fractionation is not needed, the desorption can be done only with F3. Wait 1–5 min between the two aliquots (2x5 ml) to allow sufficient contact time between the solvent and the trapped analytes.
<i>GCB (carbograph 4)</i>	Same conditioning as by C18 cartridge. The only difference is that HPLC-water was acidified to pH 2-3. The residual (the C18-preconcentrated) water is acidified to pH 3–3.5 and loaded onto the carbon cartridge (500 mg, 6 ml) at 10 ml/min
<i>Drying</i>	After sample percolation, the sorbent is completely dried under vacuum for 20–30 min.
<i>Desorption</i>	Desorption by 2x5 ml of dichloromethane:methanol (80:20, v/v) mixture at 1 ml/min waiting 5 min between the two aliquots to allow sufficient contact time between the solvent and the trapped analytes.
<i>Evaporation</i>	Combine all the extracts and evaporate under a gentle stream of nitrogen to the desired volume (0.5 ml). Before starting the evaporation, add some dimethyl sulfoxide (DMSO) to avoid loss of volatiles during evaporation.
<i>UmuC-test</i>	The extract is ready for the UmuC-test.



Examples of SPE and SSPE approaches.