National Aeronautics and Space Administration

Electrochemical Detection of the Molecules of Life

Introduction

All forms of life on Earth contain cellular machinery that can transform and regulate chemical energy through metabolic pathways. These processes are oxidationreduction reactions that are performed by four key classes of molecules: flavins, nicotinamaides, porphyrins, and quinones. By detecting the electrochemical interaction of these redox-active molecules with an electrode, a method of differentiating them by their class could be established and incorporated into future life-detecting missions. This body of work investigates the electrochemistry of ubiquitous molecules found in life and how they may be detected. Molecules can oxidise or reduce the surface of an electrode - giving or receiving electrons - and these interactions are represented by changes in current with respect to an applied voltage. This relationship varies with: electrolyte type and concentration, working electrode material, the redox-active molecule itself, and scan rate. Flavin adenine dinucleotide (FAD), riboflavin, nicotinamide adenine dinucleotide (NADH), and anthraquinone are all molecules found intracellularly in almost all living organisms. An organism-synthesised extracellular redox-active molecule, Plumbagin, was also selected as part of this study. The goal of this work is to detect these molecules in seawater and assess its application in searching for life on Ocean Worlds.



Seamus Thomson^{1,2}, Richard Quinn², Jessica Koehne²

¹ University of Sydney; ² NASA Ames seamusdthomson@gmail.com

Methods

Redox-active molecules were independently added Marche The chemically-dynamic environments of Ocean Worlds to Synthetic Sea Water solutions (SSW) with render them suitable candidates for future life-detection increasing concentrations from 1nM to 1mM. SSW missions. In context of our solar system, Europa and solution was degassed using nitrogen gas for 10 Enceladus are of significant interest for future spacecraft minutes before electrochemical missions because of their favourable conditions that could running experiments. An AgIAgCI reference electrode and sencourage life to exist, and - perhaps - sustain itself. The Pt wire counter electrode were used. The working introduction of electrochemical sensors into a science electrode was a Glassy Carbon Electrode that was payload for these missions would generate valuable data pretreated with mechanical polishing, as to whether life-created redox-active molecules are electrochemical etching at 1.2V in 1M NaOH for present in a sample and aid in identifying their class. This 45s, and sonication for 5min in deionised water. work demonstrates a useful limit of detection for glassy Differential Pulse Voltammetry (DPV) was used to carbon electrodes in screening for molecules of life, and similarly how to differentiate them (Fig 2, Table 1). The determine the electrochemical features of the molecules of interest with a scan rate of 10mV/s. successful integration of these sensors into a microfluidic Peak analysis and baseline measurements were channel has shown a relatively simple yet effective performed in OriginLab software. application of this technology (Fig 3).

Applications

Limit of Detection	
AD	10nM
oflavin	30nM
nbagin	100nM
quinone	30nM
ADH	10uM

Conclusion and Future Work

This project demonstrate that it is possible to detect naturally-synthesised redox-active molecules in synthetic sea water. The glassy carbon electrodes detected a relatively low concertation of these molecules, and could be readily incorporated into a 3D-printed microfluidic channel. A limitation for further investigation of these ubiquitous molecules is that most of them exist intracellularly, and so future studies could be directed at the lysis of cells so that the sensors can directly access these molecules. Similarly, the scaling-down of the glassy carbon macro electrode to a micron scale may prove valuable in decreasing the limit of detection. Investigating other carbon-based electrode materials such as diamond-like-carbon and boron-doped diamond are also of future interest for their wide electrochemical range and robust properties. This is important for detecting porphyrins, where a low applied voltage is needed.

